Gebrauchsanweisung/Instructions for use/Mode d'emploi/
Instrucciones de uso/Istruzioni per l'uso/Instrukcja obsługi

AXION STIM-PRO Eco

2-Kanal Nervenstimulator/2-channel nerve stimulator/
Neurostimulateur à 2 canaux/Estimulador nervioso de 2 canales/
Elettrostimolatore nervoso a 2 canali/2-kanalowy stymulator nerwów

INHALTSVERZEICHNIS / CONTENT

2 – 14	Gebrauchsanweisung
15 – 25	Instructions for use
26 – 38	Mode d'emploi
39 – 51	Instrucciones de uso
52 – 64	Istruzioni per l'uso
65 – 77	Instrukcja obsługi

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ALLGEMEINES

Was ist TENS?

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GENERAL INFORMATION

What is TENS?

TENS (transcutaneous electrical nerve stimulation) is a proven procedure in pain therapy with few side effects. It is used for symptomatic relief of acute, chronic and intractable pain.

Functional principle of TENS

Transcutaneous electrical nerve stimulation is a non-invasive, drug-free method of pain management. TENS uses tiny electrical pulses sent through the skin to the nerves to alter your perception of pain. TENS does not cure a physiological problem; it only helps control pain. TENS doesn't work for everyone, but for most patients it helps reduce or eliminate pain.

IMPORTANT NOTICES

- ▶ Read the instructions for use before use. Make sure that you observe all warnings and safety instructions. Ignoring the instructions may result in user injury or device damage.
- ▶ Report all serious incidents that have occurred in connection with the product to the manufacturer and the competent authority of the Member State (where the user is located).
- ▶ TENS unit supplies energy to the human body by means of electrode pads. Electrode pads are classified separately as Class I products. CE 2460 is only valid for the therapy device, and CE is valid for electrodes (class I products).
- ▶ Consult your physician before use.

INTENDED USE

The purpose of the STIM-PRO Eco is the symptomatic relief of acute, chronic and intractable pain. The device is used and operated by the patient himself. Use is permitted only from the age of 18 years.

INDICATIONS/CONTRAINDICATIONS

Indications

TENS is used to treat the following pain patterns: Muscle, joint or bone pain, pain associated with diseases such as arthrosis or fibromyalgia, conditions such as low back pain, neck pain or tendinitis.

EN

Contraindications

- ▶ Patients with an implanted electronic device, such as a pacemaker or defibrillator, or any other metallic or electronic implant should not undergo TENS treatment without first consulting a physician.
- ▶ Patients with heart disease, epilepsy, cancer, fever or other illnesses must not undergo TENS treatment without prior medical consultation.
- ▶ Use the electrodes only on healthy, uninjured skin.
- ▶ If heart disease is suspected, TENS treatment should be performed only under medical supervision.
- ▶ Do not stimulate via the carotid sinus nerve, especially in patients with a known sensitivity to the carotid sinus reflex.
- ▶ The stimulation must not be applied over swollen, infected or inflamed areas or skin eruptions, such as phlebitis, thrombophlebitis, varicose veins, etc.
- ▶ The stimulation must not be applied over or in the vicinity of cancer.

Side effects

The side effects of TENS units are generally mild, even with prolonged use. If disturbing effects occur, a physician must evaluate the further use of the device.

Skin irritation and minor burns are possible side effects. If these occur, discontinue use and consult your doctor.

WARNINGS & SAFETY INSTRUCTIONS

Warnings

- ▶ The stimulation emitted by the device may be sufficient to cause an electric shock. Electric current of this magnitude must not flow through the rib cage (thorax) or through the chest and the upper back, as it can cause cardiac arrhythmia. **Do not place electrodes over the heart!**
- ▶ Do not place electrodes on the front of the throat as spasm of the laryngeal and pharyngeal muscles may occur. Stimulation over the cervical cavity (throat area) may close the airway, make breathing difficult, and have negative effects on heart rhythm or blood pressure.
- ▶ Do not place electrodes on your head or in places where the electric current can flow transcerebrally (through the head).
- ▶ If TENS therapy is ineffective or becomes uncomfortable, stimulation should be discontinued until its use is re-evaluated by a physician.
- ▶ The device does not have AP/APG protection. Do not use it in the presence of explosive atmosphere and flammable mixture.
- ▶ The device may only be used according to medical advice.
- ▶ Electrodes must not be placed over the eyes, the mouth or inside the body.
- ▶ The device is to be used only when awake.
- ▶ Do not use the device during pregnancy unless prescribed by your physician.
- ▶ Simultaneous connection of a patient to a high frequency surgical equipment and the stimulator may result in burns at the site of the stimulator electrodes and possible damage to the stimulator.
- ▶ Do not operate the product within 1 m of shortwave or microwave therapy equipment as it may produce instability in the stimulator output.

Safety instructions

- ▶ Use only original accessories from the manufacturer with your device. (see scope of delivery)
- ▶ Do not use this device for undiagnosed pain symptoms before consulting a physician.
- ▶ Do not use this device while driving, operating machinery, near water, or during activities where involuntary muscle contractions may expose the user to excessive risk of injury.
- ▶ Always turn off the unit before applying or removing TENS accessories to the skin.
- ▶ Keep this device out of the reach of children.
- ▶ The device has no curative effect.
- ▶ Electronic monitoring devices (such as ECG monitors and ECG alarms) may not function properly when electrical stimulation is used.
- ▶ The device should be cleaned with soapy water on a clean cloth after each use.
- ▶ The device is an electronic product and must therefore not be disposed of with residual waste. Further information at: <https://axion.shop/en/environment-and-disposal/>

MAINTENANCE, TRANSPORT, STORAGE


- ▶ The device must not be exposed to large quantities of liquids/water.
- ▶ After use, place the device back into the foam insert of the storage case for protection.
- ▶ Remove the battery and store the device in the storage case in a dry and cool place if you are not going to use it for a longer period of time.

SCOPE OF DELIVERY

Each STIM-PRO Eco includes the following accessories, which can be reordered at **www.axion.shop**.

QTY	Product name	Order no.
1x	TENS Device	1014
1x	4x Self-adhesive electrodes 50x50mm	20007
1x	2x Connection cable	20032
1x	Storage case	101409
1x	9.0 V Block battery	00012
1x	Instructions for use	101401

OPERATING THE DEVICE



Slide down the cover of the device to reach the controls.

1 ON / OFF, Regulating the intensity

The intensity can be adjusted stepless from 0–80mA by turning the intensity regulators on top of the device. The application starts automatically when the intensity is increased.

The numbers 1 to 8 on the regulators indicate the current intensity.

2 Change mode switch

There are three TENS modes available: **(B)** burst, **(C)** normal and **(M)** pulse width modulation.

TENS Burst Mode (B)	Burst frequency fixed at 2 Hz Pulse width adjustable from 30–260 μs
TENS Normal Mode (C)	Frequency adjustable from 2–150 Hz Pulse width adjustable from 30–260 μs
TENS Pulse Width Mod- ulation Mode (M)	This mode contains pulse width modulation. The pulse width is decreased down to 60% of the set value within 3 seconds. In the next 3 seconds the pulse width is increased to the set value again. The parameters can be set from 2–150 Hz and from 30–260 μs .

3 Regulating the frequency

The frequency is adjustable from 2 Hz to 150 Hz and determines how many electrical impulses are transmitted through the skin every second. If no instructions regarding the frequency are given in therapy, set the control to 70–120 Hz (Burst Mode 2Hz).

4 Regulating the pulse width

The pulse width can be adjusted from 30 μs to 260 μs . It regulates the time for which the electrical signal is applied. If no instructions regarding the frequency are given in therapy, set the control to 150–200 μs .

5 Timer ⌚


The device comes with a timer with three possible settings: 30 minutes, 60 minutes and continuous mode. The application ends automatically when the timer runs out.

EN

BATTERY CHANGE



The battery needs to be replaced from time to time to ensure the correct functionality of the device. Remove the battery if the device is not likely to be used for some time.

- ▶ Make sure the power is off.
- ▶ Open the battery compartment by sliding down and removing the front cover.
- ▶ Remove the battery.
- ▶ Insert the new batteries.
Please check that the polarity is correct!
- ▶ Close the battery compartment.

LABEL

backside

EN

	Manufacturer	REF	Article Number
CE 2460	Conformity marking with number of the notified body (device complies with the requirements of (EU)2017/745)	MD	Medical device
	Date of manufacture YYYY-MM-DD	SN	Serial number
	Attention! Check the instructions for use for important safety-related information such as warnings and precautions.		Follow instructions for use
			Follow operating instructions
IP22	IP protection class		Temperature limit
	Type BF device		Temperature limit Storage/Transportation
	Protect from sunlight		Protect from moisture
	Do not dispose of with residual waste		2D Datamatrix Barcode: EAN, date of manufacture and serial number encoded

FUNCTION AND SAFETY TEST

For safety reasons, the following precautions should be carried out once a week.

- ▶ Check the device for any external damage:
Check for deformation of the housing or damage to the electrode output sockets.
- ▶ Check that no type plates or stickers are damaged.
- ▶ Check the indication of the LEDs when switching on.
- ▶ Check cables and electrodes for damages.
- ▶ A safety check must be performed by authorized personnel at least every two years by the end of the month in which it is due.
- ▶ The instructions for use must always be carried with the unit. If there are any problems or defects, please contact your retailer.

REPAIR MINOR FAULTS YOURSELF

If the device does not work properly, the instructions below may help identify the cause of the fault. If the problem cannot be solved with any of these tips, contact the customer service.

- ▶ Check cables and connections! Have all plugs been connected correctly?
- ▶ Does the display not switch on? Change batteries.
- ▶ If a cable is externally damaged, replace it immediately with a good one. Only use original accessories.
- ▶ No current is flowing? Check that the electrodes are correctly connected to the same channel and that the current is strong enough.

CONFORMITY TO SAFETY STANDARDS

The STIM-PRO Eco device is in compliance with the safety standards EN 60601-1 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance & EN 60601-1-2 Medical electrical

equipment - Part 1-2: General requirements for basic safety and essential performance – Collateral Standard: Electro-magnetic disturbances - Requirements and tests.

TECHNICAL DESCRIPTION

01	Channel	2 channels with adjustable intensity
02	Intensity	Adjustable 0–80 mA with a load of 500 Ohm on each channel
03	Pulse amplitude	Asymmetrical, bi-phasic square pulse
04	Output Voltage	0–40V (Vpp), 500Ω Load
05	Power Source	1 pc. 9,0V battery
06	Size	9.5 cm (L) x 6.5 cm (W) x 2.35 cm (H)
07	Weight	115 g incl. battery
08	Frequency	Adjustable from 2–150 Hz
09	Pulse width	Adjustable from 30–260 μs
10	Modes	3 TENS modes with adjustable parameters: B (Burst) C (Normal) M (Modulation)
11	Operating condition	Temperature: 10°C to 40°C Relative humidity: 40 %–90 % Air pressure: 700 hPa–1060 hPa
12	Transport and storage conditions	Temperature -10°C to 60°C Relative humidity: 30 %–95 % Air pressure: 700 hPa–1060 hPa
13	Notice	All technical values include a tolerance of ± 5 %

WARRANTY

All devices have a warranty of 24 months from the date of purchase. The warranty applies to the stimulator only and covers both parts and labor relating thereto. The warranty does not apply to damage resulting from improper handling, the failure to follow the instructions for use, loss or dropping.

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axion GmbH