Koncentrator Tlenu

Cronos TP-B1 (Medical oxygen concentrators)

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2 Features

1. Description of operation

This device contains a physical mechanism for generating and delivering rich in Oxygen 29+/-2% concentrate using a vacuum pump and micromolecules.

2. Usage

Ambient temperature: 5~40 C Humidity: less than 90%RH

3. Storage

Ambient temperature: 0~40 C Humidity: less than 80%RH

4. Control button

CAUTION: The device should be kept away from dust and any oily fumes. When the device is used in an area where the temperature and humidity are not within the recommended scale or you are at a high altitude above sea level, it may malfunction.

3 Security measures

1. Do not turn off the device by pulling the plug from the outlet.

2. Children and people with disabilities should use the device under supervision.

3. Ensure that the room is well ventilated when the device is operating.

4. Do not clog the air inlets and ensure that the device does not stand close to elements that can interfere with airflow, and clean the filter regularly, at least every 6 months.

5. A damaged power cord must be replaced with a new one by the manufacturer.

6. Protect the device from getting any liquids into it.


7. Do not shake the device while it is running, and avoid bumping its elements.

NOTE: The power cord is an integral part of the device and we do not recommend its Replacement with other replacements.

4

Koncentrator tlenu

Product description


1. Oxygen breathing apparatus
2. Concentrate generator
3. Oxygen outlet
4. Oxygen tube
5. Capsule
6. Control button

5

Koncentrator tlenu

Installation and use

1. After unpacking, disconnect all components from the temporary mounting frame.

2. Connect all parts to the oxygen line and then to the generator concentrate.
3. Plug the power supply into an outlet.
4. Press the "Power" button to initiate concentrate generation. A lit "20" LED indicates that the device will generate concentrate for 20 minutes. Press the "SET" button to change the device's operation time between 20/40/60 minutes.
5. To turn off the device, first press the power button and then unplug the power cord from the outlet.

NOTE: We strongly recommend that you choose the 40-minute mode and use the device daily.

6 Membrane separation technology

Koncentrator tlenu

Since the velocity and solubility at the penetration of particles in the membrane varies, this technology can easily separate oxygen and nitrogen molecules producing oxygen concentrations of 28%+/-2%.

1. Applying a state-of-the-art biological separation method to produce oxygen-rich concentrate based on the theory of different degrees of solubility and rates at particle penetration in membranes.
2. The device generates a pure oxygen-rich concentrate.
3. The macromolecular membrane can fully filter out harmful substances in the air such as dust, viruses and bacteria, and refresh the air to make it health and environmentally friendly.

7 5 groups of potential consumers

Koncentrator tlenu

The oxygen-rich concentrate carries broad health benefits.

Regular breathing of oxygen-rich concentrate increases physical abilities and vital activity.

Office women: skin care and improved skin elasticity and nutrition. Businessmen: physical and mental relaxation, improved physical fitness and productivity at work.

Students: improving memory, refreshing and reducing mental fatigue.

Manual workers: reduction of physical and mental fatigue

Lovers of high altitude sports: improvement of aerobic capacity which helps on the high altitudes.

8 Occurrence of problems

1. During the initial period of operation, the genset and oxygen tube may generate unpleasant odors.
2. If the device is not placed on a hard surface, it may generate annoying sounds.
3. Tilting and lifting during operation can generate excessive noise.
4. The concentrate generator can be stopped by the power supply if an unforeseen event occurs. This is to protect the device from short circuits.
5. When the ambient temperature drops below 20C, the device will automatically start to heat up which can increase the airflow for heat removal while increasing noise emissions.

9 Technical parameters

1. **Model: TP-B1**
2. **Oxygen Concentrate: 28% +/- 2%**
3. **Airflow: ~3L/m**
4. **Power supply: AC100V-240V / DC12V-16.8V**
5. **Power consumption: 32W**
6. **Dimensions: 230x110x160(mm)**
7. **Volume: <48dB**
8. **Net weight: 2kg**
9. **Package dimensions: 30x20x24(cm)**
10. **Material of construction: ABS**

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