MD

Leaflet approval date: 15.02.2022
Data zatwierdzenia ulotki: 15.02.2022

info@reatape.com
www.reatape.com
tel.: +48 22 230 2057

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absorb inflammation

muscle can easily contract

- to remove, pull with direction of hair and move finger along to hold skin down, avoid lifting the skin,
- Keep out of reach of children,


What is RT Method?

It is the new method of taping developed by RT R&D team through long period of research and experience, which increases the effect of general taping methods. The basic principal of RT taping is to maintain broad type of 'I' shaped RT to effectively cover problematic areas of body.

Key Point

- Apply with correct application method,
- Test patch is recommended for those with high allergic tendencies,
- Stop usage immediately in the result of skin irritation,
- Apply with little or no stretch,
- Elongate the muscles or joints prior to applying tapes,
- Remove any excess oils, sweat, and/or dirt prior to applying the tape,
- Avoid storing product in direct sunlight and/or high temperatures,
- To remove, pull with direction of hair and move finger along to hold skin down, avoid lifting the skin,
- Keep out of reach of children

Application G

Step 3. TAPE 1
Apply the base of the first tape next to the ear.
Apply the lower strip of the tape towards your nose.
Apply the upper strip of the tape towards your forehead.

Step 4. TAPE 2
Apply the base of the second tape next to the ear.
Apply the lower strip of tape towards your nose.
Apply the upper strip of tape towards your forehead (above the brow).

Step 3: TAPE 1
Apply the base of the first tape in the middle of the lower jaw. Apply the lower part of the tape towards the lower part of the nasal lobe. Apply the middle part of the tape towards the middle of the nasal lobe. Apply the upper part of the tape towards the corner of the eye.

Step 4: TAPE 2
Apply the second tape to the ear area. Apply the lower part of the tape towards the bottom of the eye. Apply the middle part of the tape towards the outer corner of the eye. Guide the upper part of the tape under the eyebrow (eyelid).

Step 3. TAPE 1
Apply the base of the first tape to the inner corner of the eye and guide it downwards.

Step 4. TAPE 2
Apply the base of second tape approx. in the outer corner of the eye and guide it down.

Step 5. TAPE 3
Apply the base of the third tape below the second tape and guide it down.

Step 3, TAPE 1
Apply the base of the first tape in the area of the nasal lobe and attach it.

Step 4, TAPE 2
Apply the second tape to the temple area and guide it towards the cheekbone.

Step 5, TAPE 3
Apply the base of the third tape approx. in the ear area; apply the first piece towards the lower jaw, the second towards the cheek, and the third towards the cheekbone.

Step 3. TAPE 1
Apply the base of the first tape to the center of the collarbone.
Apply one piece of tape to the outside of the breast and the other to the inside of the breast.

Step 4. TAPE 2
Repeat on the other breast.

Round the corners,

Step 3, TAPE 1
Apply the base of the first tape to the lower part of the sternum (xiphoid process).
Apply the ends of the tape to form a 'fan'.

Step 4, TAPE 2
Repeat for the second tap.

Step 3. TAPE 1
Apply the base of the first tape to the hip spine.
Apply the first part of the tape over the navel,
the second under the navel.

Step 4. TAPE 2
Repeat for the other tape.

 **Storage and handling condition:**
Keep away from sunlight

② Single use product / Do not re-use

Caution:

- Stop usage immediately in the result of skin irritation
- Remove any excess oils, sweat, and/or dirt prior to applying the tape
- Elongate the muscles or joints prior to applying tape
- Applying tape with little or no stretch
- To remove, pull with direction of hair and move finger along to hold skin down, avoid lifting the skin
- Do not use the device over long time
- Avoid storing product in direct sunlight and/or high temperatures,

Patient target group:
Who have a muscle pain, joint pain

Device performance:
Reduce pain

Occurred serious incident, report it to the relevant authorities.