Operating instructions for pH and Redox electrodes

EN

Preparing the electrode

Carefully remove watering cap.
Flush electrode with water.

pH: Check inside of pH glass membrane for unwanted air bubbles. Shaking the electrode gently will cause any bubbles to rise to the top.

Storing the Electrode

Electrodes should be stored with the watering cap on, containing 3M KCl solution or Storage Solution. Electrodes stored dry will exhibit unstable values for a short period. If the electrode inadvertently dries out, it can be regenerated by placing it overnight in Storage Solution.

Cleaning

In general, you can use acids, alkaline solutions and commonly used solvents for brief periods to clean electrodes with a glass shaft. Flush with water after cleaning. After cleaning, electrodes are likely to exhibit sluggish response times for a certain period, so place it in storage solution for 15 minutes after cleaning.

Regeneration

pH: Immerse electrode for 10 minutes in 0.1 - 1.0 M NaOH, then 10 minutes in 0.1 - 1.0 M HCl. After regeneration, place it in storage solution for 15 minutes. **Redox:** clean metal surfaces with a mildly abrasive medium, e.g. toothpaste or very fine scouring powder.

Disposal

 The design of this sensor minimizes environmental impact. According to the EU directive 2012/19/UE this sensor should be disposed as waste of electrical and electronic equipment, and not in municipal waste, or it can be sent back to your supplier for disposal.

General

The life span of electrodes depends on the customer requirements of response time, zero point and slope. Measurement conditions - especially high temperatures and aggressive samples - may shorten life span. In favorable conditions, life expectancy can be 1 - 3 years at ambient temperature, and 1-3 months at 90°C. Some ageing also occurs during storage.

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