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INSTRUCTIONS FOR USE

Issue number and date: 1 / 11.2024

Anti-bedsores bubble mattress with alternating pressure

Model M100AM

Introduction

Thank you for choosing our product. We are confident that it will meet your expectations. We encourage you to read this document before using the product for the first time. These instructions are intended to provide the necessary information on the safe use of the medical device. The manual provides detailed guidance on the device's functions and capabilities. Retaining the document allows this information to be easily found in the future, facilitating correct use.

Description of the medical device

The anti-bedsores variable pressure bubble mattress (190x82 cm; with cushion) consists of two parts: the mattress and the pump. The pump uses a small compressor, making it quiet and energy-efficient. The control panel is simple and very easy to use. The mattress reduces pressure by sequentially deflating and inflating alternating air cells at intervals of approximately 6 minutes. It is widely accepted that constant pressure on bony protuberances constitutes a major cause of skin damage. The continuous movement provided by the device reduces areas of constant pressure and improves circulation.

Intended use

A medical device designed for people who are immobile or have limited mobility and are at risk of developing first - and second-degree pressure sores. Its main purpose is to prevent pressure sores by continuously varying the pressure in the individual cells of the mattress, which prevents prolonged pressure on the skin and improves blood circulation, reducing the risk of tissue damage. This mattress is used in hospitals, hospices, and in homes with people requiring long-term care.

Contraindications

The product should not be used by patients with allergies to the materials used, sleep or orientation disorders and without a basic mattress in medical beds, which may increase the risk of skin damage.

The user or their carer should consult a healthcare professional if they are unsure whether the medical device is appropriate for the patient in question.

Maximum permissible user weight – 135kg.

Any use not in accordance with the intended use is considered inappropriate and DANGEROUS. The manufacturer and/or distributor shall not be liable for damage caused by improper, incorrect and/or unreasonable use.

Set description

The anti-bedsore variable pressure bubble mattress (190x82 cm; with cushion) consists of two parts: the mattress and the pump.


Figure 1. Variable pressure pump


Figure 2. Electrical circuit diagram

Pump and mattress operation

1. Fix the pump in a horizontal position or suspend it from special hooks.
2. Spread the mattress on the bed surface and secure it by tucking the trim under the bed frame. Ensure that the air intake is located at the foot end of the bed.
3. Connect the air supply tubes to the mattress outlet and the pump inlet. Ensure that the tubes are not creased or twisted to ensure free airflow. Allow the inner compartments of the mattress to lay freely.
4. Connect the pump to the power source and start it using the appropriate switch. The indicator light should illuminate on start-up.
5. Move the pressure adjustment knob to the centre position. Leave the mattress for approximately 20-30 minutes to fully inflate all cells.
6. Cover the mattress with a clean sheet, avoiding stretching the fabric too much.
7. In order to adjust the firmness of the mattress to the individual patient's needs, the air pressure can be adjusted using an adjustment knob.

Mattress operation phases

Phase 1: All cells of the mattress are fully inflated.

Phase 2: The air is removed from the first part of the cells.

Phase 3: All cells are again filled with air.

Phase 4: The air is removed from the second part of the cells.

Safety considerations

- 1 Do not smoke on or near the pump.
- 2 Keep the pump away from hot surfaces, other heat sources.
- 3) There is a risk of explosion if used in the presence of flammable anaesthetics.
- 4 Replace the fuse as indicated: T1A 250V.
Power source: ~ 220V 50Hz, 0.1A.
- 5 The mattress can be inflated immediately after purchase, but can only be used 8 hours after the first inflation (ideally 12 hours), as the internal pull strap and air bed connectors require a buffering process.

Maintenance

Disconnection of the pump from the power supply is required before maintenance operations can be carried out.

- Regular cleaning of the product should be carried out with water, a sponge and mild detergents intended for plastics, such as soap or a mild cleaning agent approved by the National Institute of Hygiene. It is important to dry the product thoroughly after each use.
- Avoid the use of: strong detergents, abrasives, bleach, petrol and all solvents, as these may cause damage to the product, visible defects and lead to a loss of warranty.

Storage and transport conditions

The product should be stored in a dry and well-ventilated place, protected from moisture and extreme temperatures to prevent damage to the materials. Contact with chemicals that may affect its durability and performance should be avoided. During transport, the product should be secured in such a way as to prevent displacement and damage. Use appropriate protective packaging.

Disposal

The appliance is subject to separate waste collection. Failure to comply with the obligation to hand over equipment to an authorised collection point may result in administrative sanctions. Electrical and electronic equipment compliant with Directive 2012/19/EU. The appliance, accessories and packaging must be disposed of correctly at the end of its use. Local ordinances or regulations regarding disposal must be observed. The product must be taken to a separate collection point, where it will be recovered and recycled.

Product life

The expected product life is up to 2 years if used in accordance with the information in the instructions for use. The product can be used after this period if it is undamaged and meets the specified requirements for safe use.

Warranty

- The product is covered by a 24-month warranty from the date of purchase, against any defects in material or workmanship.
- The warranty covers the free replacement and/or repair of parts found to be manufacturing defects.
- The warranty does not cover parts subject to normal wear and tear.
- The warranty is no longer valid if the product has been tampered with/modified, if the defect is the result of misuse or if the damage cannot be attributed to the manufacturer.

NOTE

Any serious incident involving the device must be reported to the manufacturer and to the competent authority of the Member State in which the user or patient is residing.

Technical specifications

Dimensions MATERIAL (length x width x height): Extension flaps: PUMP (length x width x height):	190cm×82cm×7cm 35cm (for both sides) 25cm×12cm×10cm
PRODUCT WEIGHT Mattress and pump	3 kg
BEARING CAPACITY MAXIMUM PATIENT WEIGHT:	135 kg
ELECTRICS OPERATING POWER SUPPLY: FUSES:	~ 220V 50Hz, 0.1A MAX. T1A 250V

ENVIRONMENTAL CONDITIONS OPERATING CONDITIONS Ambient temperature: Relative humidity:	10°C to 40°C 10% to 75%
STORAGE AND TRANSPORT CONDITIONS Ambient temperature Relative humidity:	-18°C to +43°C 10% to 95%

Product content: PVC (Polyvinyl chloride) - mattress; Plastic - pump casing.

Colour: Brown.

Explanation of symbols

	Conformity of the device with EU Regulation 2017/745 on medical devices		Medical device
	Follow instructions for use		Maximum user weight
	Batch/series number		Model
	Serial number		The product in question must not be disposed of as unsorted waste. WEEE product.
	Use before date (YYYYMMDD)		Permissible relative humidity range for unit operation
	Type B application part		Permissible environmental temperature range for operation
	Protect from moisture		Keep out of sunlight
	Date of production		Warning
	Unique identification code of the device		Medical device manufacturer

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ACTION S.A.