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Up to date user manual

Using the QR code, download the up to date user manual containing the latest, updated information regarding the use, maintenance and troubleshooting of the device. The manual is available at welltec.pro/r/manual/ock

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Dear Customer,

Thank you for your trust and for choosing the **Welltec** oxygen concentrator. We hope that it will provide you a clean oxygen for many years.

Table of contents

Safety notes and warnings	3-7
Safety when inhaling oxygen	8
Safety when nebulizing	9
Symbols	10
Contents of the package	11
Product introduction	12
First use	13-14
Usage	15-17
Remote controller	18
Nebulizer usage	19
Inhalation headset usage	20
Cleaning and maintenance	21-23
Specification	24
Troubleshooting	25-27
Error table	28
Pulse oximeter introduction	29
Pulse oximeter safety notes and warnings	30-31
Pulse oximeter information	32
Pulse oximeter usage	33
Pulse oximeter charging	34
Pulse oximeter cleaning and maintenance	34-35
Pulse oximeter specification	36
Pulse oximeter troubleshooting	37

Safety notes and warnings

To prevent loss of health or damage to property, read the following instructions.

Risk of health loss or death

Caution

Risk of health loss, hazardous substances

Prohibited

Should be done

To prevent oxygen concentrator's possible failures or power shut down, people in urgent need of oxygen and seriously ill patients must prepare other oxygen-supply devices for emergency use (etc: oxygen cylinders, oxygen bags).

The concentrator is suitable for oxygen supplementation and is not intended to be used for life support or for use on the unconscious.

Do not use the equipment near flammable materials as grease oil, detergent etc. In a certain pressure, oil, grease or grease substances which contacted with oxygen will self-ignition and intense combustion. These substances must be kept away from oxygen concentrator, piping, connectors, and all other oxygen devices. Do not use any lubricant other than the manufacturer's recommendation.

Oxygen is the gas that sustains combustion. **Absolutely no smoking near a working oxygen concentrator.** Keep the oxygen concentrator away from matches, burning cigarettes and other objects with high temperatures or open flames. Textiles and other normally non-combustible materials become flammable in oxygen-enriched air. Accessories should also not lay on such materials when connected to an operating device and not in use. Ignoring this warning could result in severe fire, property damage, personal injury, or death.

Personal and family use should be in accordance with the guidance of a physician when the concentration is higher than 93%.

Precautions and warnings

-  If any object or liquid enters the device, disconnect the power plug immediately and have it professionally inspected before using it again.
-  Do not place the oxygen concentrator in parallel or in series with other concentrators or oxygen treatment equipment.
-  Do not disassemble.
-  Do not block the air in/outlet under any circumstances.
-  Check the electrical connection regularly.
-  Turn off the device before disconnecting the power plug.
-  Only use the original plug without an extension cord.
-  After cleaning the filter elements, dry them completely before using them again.
-  Make sure the device is properly connected to the electrical outlet.
-  Do not expose the device to direct sunlight, wind or rain.
-  Keep a distance of 50 cm or more between the concentrator and walls, windows, furniture or other similar objects.
-  Make sure there are no humidifying devices in the same room or within 2 meters.
-  Do not place next to a source of heat or steam.
-  Do not use in explosive or corrosive atmosphere.
-  Do not allow the device to get wet.
-  Do not pierce or incinerate the device.
-  Protect the device from running water.
-  Do not spray the device with insecticides.
-  Do not touch with wet hands.
-  Do not clean with solvents.

Precautions and warnings

-  Before switching on, check that the supply voltage from the rating plate corresponds to the local mains voltage.
-  Place the device on a flat and stable surface to prevent water leakage.
-  The temperature near the air outlet will be higher when the concentrator is working, keep away from the air outlet to avoid burns.
-  Be careful when transporting or moving it and prevent tipping it over or tilting it at an angle greater than 5°.
-  Do not move the oxygen concentrator while it's working.
-  Do not step, sit or lie on the oxygen concentrator.
-  Wait 5 minutes before restarting the device.
-  Always grasp the plug when disconnecting the device.
-  The device should be used in a well-ventilated area. Do not use the device in confined spaces or in an environment with obstructed airflow.
-  Do not use in a smoky space or in a polluted air.
-  Do not lubricate the connections, nasal cannula or other accessories of the oxygen concentrator to avoid the risk of fire and burns.
-  Contact professional medical personnel if you feel discomfort or have had an accident while using the concentrator.
-  Geriatric, pediatric or other patients who are unable to communicate discomfort may require additional monitoring and/or an alert system to communicate the discomfort and/or urgent medical need to a responsible caregiver to avoid harm to health.

Precautions and warnings

Be careful to transport or convey, don't let it reversal or thwart, tilt angle not greater than 5°.

Geriatric, paediatrics or any other patient unable to communicate discomfort can require additional monitoring and or a distributed alarm system to convey the information about the discomfort and or the medical urgency to the responsible care giver to avoid harm.

Use only water-based lotions or salves that are oxygen-compatible before and during oxygen therapy. Never use petroleum or oil-based lotions or salves to avoid the risk of fire and burns.

Do not touch the cabinet or air outlet with your hands during the operation of the oxygen concentrator to avoid injury caused by overheat.

The temperature near the air outlet will be higher while the concentrator is operating, do not close to the air outlet to avoid being burned.

Do not move the oxygen concentrator when charged.

Personal and family use should be in accordance with the guidance of a physician when the concentration is higher than 93%.

Smoking during oxygen therapy is dangerous and is likely to result in facial burns or death. Do not allow smoking within the same room where the oxygen concentrator or any oxygen carrying accessories are located.

Do not lubricate fittings, connections, nasal cannula, or other accessories of the oxygen concentrator to avoid the risk of fire and burns.

Open flames during oxygen therapy are dangerous and is likely to result in fire or death. Do not allow open flames within 2 m of the oxygen concentrator or any oxygen carrying accessories.

The device should be placed on a flat and stable surface to prevent water leakage.

Do not tread, sit on or lie on the oxygen concentrator.

Precautions and warnings

Attention!

If any of the following events occur, turn off the device immediately:

- One or more switches do not work.
- A short circuit has occurred.
- Disturbing overheating of the power cord / plug.
- Burning smell or disturbing sounds and vibrations.
- Other malfunctions.

If the power cord is damaged or destroyed, it must be replaced by the manufacturer, service or another authorized entity.

The device can be used by people with limited physical, sensory and mental performance on condition of assurance appropriate supervision and instruction in the field of use. Children should not play or use the device.

Read the device's operating instructions and warnings carefully. To the extent permitted by applicable law, Welltec accepts no liability, express or implied, for any: loss of profits, use of the Product, functionality, contracts, transactions, revenues or anticipated savings, increased costs, or expenses, or no other indirect, consequential, or special loss or damage, in particular caused by using the device contrary to the operating instructions.

Attention!

The device is intended for use in rooms with a temperature from 5 to 40 °C, humidity below 75%, barometric pressure range from 860 hPa to 1060 hPa and altitude below 800 m.a.s.l. Operation of the device outside the specified temperatures, humidity, altitude, and atmospheric pressures may reduce the oxygen concentration.

Safety when inhaling oxygen

Absolutely no smoking when using the oxygen concentrator.

Do not move the oxygen concentrator while it's in use.

Follow physician's instructions if the concentrator is used to support treatment.

Not to be used by patients allergic to plastics or at risk of oxygen poisoning.

Clean or replace the sponge filter if the oxygen outlet is clogged. This will extend the life of the oxygen concentrator.

To maintain normal operation of the concentrator, replace the filter sponge with one provided or the other certified by the manufacturer.

The oxygen flow rate should not be too high. If the user or caregiver believes that the oxygen content is insufficient, contact the supplier or physician immediately and adjust the flow as directed by the physician. Patients with severe lung disease should consult a physician to determine the appropriate flow rate.

Excessive oxygen inhaling can cause some problems, e.g. oxygen poisoning (including CO₂ retention), irritating and dry cough, nausea, vomiting and headaches, damage to the nasal passage or nosebleeds. Please use the device with care and always consult your doctor before use to avoid personal injury.

The nebulization set and nasal cannulas included with the device are test samples. The capacity of the nebulization cup is 6-12 ml and the length of the nasal cannula is 1.6 m or 6 m. Use the nasal cannula personally, clean and disinfect it often. It is recommended to use a nebulization set and nasal cannulas produced by qualified manufacturers.

There should not be too much water in the humidifier bottle to prevent overflow. Keep the water level between the „MAX“ and „MIN“ marks and change the water frequently. Different bottles may affect the performance of the oxygen concentrator, so replace the bottle with the original one supplied by the manufacturer.

Safety when nebulizing

Before use, check that the nebulization set and the nebulization port on the oxygen concentrator are firmly connected. Check for leaks that could affect the nebulization of the liquid.

The water in the nebulization cup must be clean, free of precipitation and impurities, otherwise its inlet or outlet may be blocked.

The water added in nebulization cup should not exceed 12 ml each time.

To prevent cross infection, it is recommended that each patient use different nebulization set.

After the nebulization set is connected, the flow level of the oxygen concentrator should be adjusted to the minimum.

Use only the provided nebulization set or one certified by manufacturer, other nebulization sets may affect the nebulization effect.

Some symbols may appear on your device.

Symbol	Description	Symbol	Description
	General warning sign		No open flame: Fire, open ignition source and smoking prohibited
	Caution		No smoking
	Keep dry		Fragile-handle with care
	This side up		Increase or decrease (Knob)
	ON (power)		Class II equipment
	OFF (power)		Alternating current
	Type BF applied part		Date of manufacturing
	Refer to instruction manual/booklet		Manufacturing
	Consult instructions for use		Serial number
	Degree of protection provided by enclosure		Batch code
	Part Number		Fuse
	Max. Stack Quantity		WEEE symbol
	CE mark		Authorized representative in the European community

Contents of the package

Oxygen concentrator

Remote control

Pulse oximeter

Nasal cannula 1,6 m (2 pcs.)

Nasal cannula 6 m (2 pcs.)

2x Oxygen mask

Nebulization set (2 pcs.)

Inhalation headset (2 pcs.)

Splitter for connecting 2 people

Sponge filter (8 pcs.)

Fuse (4 pcs.)

Pulse oximeter usb cable

Outlet sponge filter (4 pcs.)

Instruction manual

Product introduction

- 1.** Port for connecting a humidifier bottle with an oxygen outlet
- 2.** Nasal cannula port
- 3.** Control panel
- 4.** Air outlet
- 5.** Handle
- 6.** Casters
- 7.** Humidifier bottle
- 8.** Nebulization outlet

- 9.** Air filter
- 10.** Fuse socket
- 11.** Power switch
- 12.** Air tube
- 13.** Nasal cannula port
- 14.** Oxygen outlet
- 15.** Humidifier bottle mounting
- 16.** Humidifier bottle holder
- 17.** Power cord

First use

NOTE. The device should stand upright for 4 to 6 hours before first use. This prevents damage to the compressor.

Before turning on, check if the concentrator has been connected to power supply.


1. Remove the lid of the humidification bottle the bottle with clean water (or distilled water). Make sure the water level is between the „MAX“ i „MIN“ marks on the bottle, then cover the top.

2. Mount the humidification bottle in the bottle holder of the concentrator's humidifier and secure it with a band.

3. Connect the humidifying bottle (A) to the oxygen outlet (D) with the connecting tube (B) and connect the other port of the humidifying bottle (C) to the nasal cannula or oxygen mask (as shown on the right).

First use

4. Set the power switch on the back of the case to the „I“ (ON) position. If power is connected, the power indicator will turn green.
5. Put on the nasal cannula, oxygen mask or inhalation headset.
6. Press the  button on the front panel to start the flow of oxygen.
7. Adjust the cannula to the correct position to facilitate oxygen inhalation for maximum comfort. The concentrator reaches 90% oxygen concentration in 12 minutes.
8. Turn off the concentrator and unplug it after use. Keep the nasal cannula and oxygen tubing clean for reuse.
9. If the concentrator will not be used for a long time, empty the water from the humidifier bottle. After cleaning and drying, keep the humidifier bottle.

Warning!


- During operation, the concentrator makes a „TSSSS“ sound every few seconds, this is a normal ventilation sound.
- Uptake time and oxygen flow adjustments should be as directed by your physician.
- If the air flow is imperceptible or low (less than 0.5 L/min), check that the tubes at the connectors are properly connected or that the tubes or accessories are not clogged, kinked or damaged. If, after making sure that all the equipment is in order, the air flow is still not felt or there is little, contact the seller of the device.

1. Power ON/OFF


Turn the „I/O“ power switch on the back of the concentrator to the „I“ (ON) position. The oxygen concentrator is in standby mode. Press the  (1) button on the control panel to start the device.

Within three seconds before restarting the concentrator, the system will automatically perform a test during which the device may beep and the LCD will remain off. The oxygen concentration can reach a stable state after about 12 minutes of concentrator operation.

2. 3. 5. Timer

Press the  or  button to program the time after which the device will turn off (possible settings are 1 - 10 hours), then press the button to confirm the setting. After starting the timer, the system proceeds to the countdown and the concentrator's LCD display (5) shows the remaining time. When the countdown reaches 0, the oxygen concentrator will enter standby mode.

4. Flow knob


Turn the knob on the control panel to adjust the flow. Left to increase the flow or right to decrease it. The higher the flow value, the greater the air flow and the lower the oxygen concentration. The recommended maximum flow rate of the OCK3 series oxygen concentrator is 3 l/min. The recommended maximum concentration is 90% - 96%. The recommended maximum flow rate of the OCK5 series oxygen concentrator is 5 l/min. The recommended maximum concentration is 90% - 96%.


6. Total run time

The total run time is the sum of the operating hours of the oxygen concentrator counted from the first use. 

7. 8. Flow and concentration

The display shows the current oxygen flow values. It can be changed by turning the knob (4). The higher the flow, the lower the oxygen concentration. The higher the oxygen concentration, the lower the flow. 

Function between oxygen concentration and nominal flow (nominal outlet pressure is 0 kPa):

Effect of altitude on oxygen concentration delivered by the device:

Altitude (m.a.s.l.)	≤800	800 - 1500	1500 - 3000	3000 - 4000
OCK3	≥90%	≥82%	≥75%	≥65%
OCK5	≥90%	≥80%	≥72%	≥63%

9. Power indicator

This indicator turns green when the oxygen concentrator is connected to a power source and the switch on the back of the housing is in the „I“ position. 

10. Air pressure / concentration / fan / temperature indicator


- If the air pressure in the system is higher than 260 kPa or lower than 20 kPa, this indicator will turn yellow and the corresponding error code will be displayed on the screen.
- When the oxygen concentration drops below 82%, the indicator will turn yellow.
- When the cooling fan is disconnected or damaged, the indicator will turn yellow and the corresponding error code will be displayed on the screen.
- When the system temperature exceeds the maximum allowable temperature limit (65°C), the oxygen generator stops operating. The indicator will turn yellow and the corresponding error code will appear on the screen.
- When the flow knob or oxygen path is blocked and flow adjustment is not possible, the indicator will turn yellow and the corresponding error code will appear on the screen.

A table listing the individual error codes is on page 26.

11. Power failure alarm indicator

If there is a power failure while the concentrator is operating, this indicator lights up red and the device sounds the alarm. Unplug the power cord and plug it back in when the power is back to normal.

12. Ionization function

A lit O_2^- icon means that the ionization function is turned active.

Remote controller

Overview:

1. Working indicator (the indicator lights green when using buttons of remote controller)
2. Timing 1 hour button
3. Turn on/off the device
4. Timing 2 hours button
5. Timing 3 hours button

Press the  button on the remote control to start the device. The oxygen concentrator must be plugged in and the power switch on the back of the case must be in the „I“ (ON) position. Press the  button again to turn off the device.

Press the corresponding timer button on the remote control to set the timer length.

If the remote control does not work, the battery may be exhausted. Periodically replace the battery with the same type with the correct orientation (positive and negative poles).

Warning!

- The concentrator should be in the remote range when using the controller.
- The maximum remote control distance is 5 m.
- Take the batteries out of the controller if unused for a long time.

Nebulizer usage

1. Turn the „I/O“ power switch on the back of the oxygen concentrator to the „I“ (ON) position.
2. Fill the nebulization cup with the medicinal preparation in the dose recommended by the physician. The water used for the preparation should always be pure.
3. Connect the cup with the medicinal preparation to the rest of the nebulization set.
4. Unscrew the cap on the nebulizer port on the concentrator. Connect the spray tube to the port.
5. Press the  button on the control panel to start the device. Then turn the airflow knob clockwise to reduce the flow to a minimum.
6. Put on the nebulization set and start therapy. Follow your doctor's instructions to use the proper nebulization time.
7. After nebulization is complete, press the  button on the control panel to turn off the device. Remove the nebulization set and turn off the power to the device by turning the power switch to the „O“ (OFF) position.


Nebulization performance

The nebulization outlet pressure is 60 kPa to 250 kPa, and the maximum air flow rate is 3 l/min. The spray rate of the nebulizer kit is ≥ 0.2 ml/min.

Warning!

- Oxygen inhaling and nebulization cannot be used at the same time.
- Adjust the flow rate to minimum value when using the nebulization function.

Inhalation headset usage

Overview:

- | | |
|---|--------------------------|
| 1. Snap hook | 4. Fast flow designation |
| 2.  flow selection marks | 5. Upper cover |
| 3. Slow flow designation | 6. Lower cover |

Inhalation headset usage

Thanks to the possibility of putting it on the head, this set allows for comfortable oxygen inhalation. Essential oils can also be added to the water cup for aromatherapy.

Adjust the flow rate of headset

The inhalation set has two levels of oxygen flow: slow and fast. To set the desired mode, rotate the upper cover (5) to align the  mark of the lower cover (6) with the desired level (2). A „click“ sound will be heard when the markers are correctly aligned.


Wearing the headset

When inhaling oxygen, wear the headset as shown in the illustration.

Water cup usage

Adding essential oil: turn the top cover (1) according to the arrow „OPEN  CLOSE“ (2) to open the cup, add the proper amount of essential oil, and then cover the top cover to the „CLOSE“ position. Install the water cup correctly and then you can turn on the concentrator.

Cleaning and maintenance

It is recommended that the operating time of the oxygen concentrator is not less than 30 minutes each time. Do not turn the oxygen concentrator on or off frequently. Unplug the concentrator at least 5 minutes after shutdown for the compressor protection.

Warning!

Before cleaning the device, disconnect it from the power supply.

Housing cleaning

With normal use, the housing of the device should be cleaned once a month.

The housing of the concentrator cannot be washed with water. Wipe its surface with a clean, soft cotton swab or sponge. Do not pour liquid into the slots of the casing. You can use a mild, neutral household detergent to wipe the housing.

Humidifier bottle maintenance and cleaning

Change the water in the bottle daily. Empty the bottle each time you finish using the device.

Use clean water to clean the humidifier bottle. If the bottle is dirty, you can clean it with a neutral detergent or a 1:10 solution of vinegar and hot water. After removing the dirt, wash the bottle with warm water, then refill with clean or distilled water before the next use.

Headset cleaning

Rinse or immerse the disassembled water cup in warm water for about 15 minutes. Add the right amount of vinegar to the water to clean it more thoroughly. Then clean the cup with diluted cleaning agent and rinse it with clean water.

Do not clean the set components with very hot water to prevent deformations caused by heat. After cleaning, all components must be thoroughly dried. Use a clean, soft cloth moistened with medical alcohol to wipe the air outlet.

Cleaning and maintenance

Air filter cleaning or replacement

Cleaning and replacing the air filter is important to protect and extend the life of the compressor and concentrator molecular sieve. It is recommended to clean the filter after approximately 100 hours of concentrator operation, but it should not be cleaned more than 5 times. If the filter is damaged or heavily worn, replace it with a new one.

Filter removal

The filter sponge is located on the back of the oxygen concentrator, as shown, pull the filter out and remove the sponge (1) from the filter.

Sponge cleaning

Clean the sponge with diluted detergent and wash with water. The sponge must be completely dry before it can be installed in the filter. Recommended cleaning agents are neutral household detergent or a 1:10 solution of vinegar and hot water.

Warning!

Do not use the oxygen concentrator without a sponge filter installed or with a wet or damp sponge. This will seriously damage the oxygen concentrator.

Cleaning and maintenance

Cleaning or replacing the outlet filter

It is recommended to clean the outlet filter after approximately 100 hours of operation of the concentrator. If the filter is damaged or heavily worn, replace it with a new one. The air filter is located in a drawer on the side of the concentrator.

Pull the filter (1) out of the device and clean it with diluted detergent and water.

The sponge must be completely dry before it can be installed. Recommended cleaning agents are a neutral household detergent or a 1:10 solution of vinegar and hot water.

Fuse replacement

If you suspect that there is an overload at the start of the device (i.e. when the mains works normally, and after starting the device there is an alarm - red flashing power failure indicator and a sequence of triple and double beeps), replace the fuse. Turn off the concentrator and unplug it, then use a tool to remove the screw (1) and replace the fuse with a new one (2) before restarting the concentrator.

Specification

Model	OCK3	OCK5
Factory model	KSW-3	KSW-5
Rated flow	3 l/min	5 l/min
Max oxygen concentration (V / V%)	93% ± 3%	
Operating mode	Continuous	
Control	• Control panel • Remote control	
Display	Yes	
Functions	• Timer • Nebulization • Ionization • Total running time • Error display	
Nebulization	≥0.2 l/min	
Nebulizer port pressure range	60 kPa - 250 kPa	
Nebulization outlet pressure	30 kPa - 80 kPa	
Equipment	pulse oximeter nasal cannula 1,6m nasal cannula 6m 2x oxygen mask 2x nebulization set 2x inhalation headset splitter for connecting 2 people 3x sponge filter 3x fuse	
Operating range - temperature	from 5 to 40°C	
Operating range - humidity	up to 75 %	
IP protection level	IP21	
Storage - temperature	from -40 to +55°C	
Storage - humidity	up to 93%	
Fuse	F5AL250V	
Noise level	53 dB(A)	54 dB(A)
Power	350 W	450 W
Power supply	220-240 V, 50 Hz	
Weight	18 kg	18 kg
Dimensions	410 x 265 x 530 mm	
Package size	655 x 472 x 322 mm	

Oxygen concentrator alarm and indicator system

The alarm system of this device is designed to monitor its operation and notify in the event of situations such as power failure or abnormal pressure. All alarms on this device are technical alarms.

Alarm	Cause	Indicator	Possible solutions
Power failure. Triple and double beeps sequence.	Mains power was disconnected during operation.	The power failure indicator flashes red.	Disconnect the device from the power supply. Check that the mains power supply is normal and the cable is not damaged. If the alarm still occurs after reconnecting, contact the seller.
Low concentration. Single beep.	The oxygen concentration is less than 82%.	The alarm indicator lights up yellow.	Contact the seller. For those who need oxygen urgently, reserve oxygen should be prepared.
Incorrect air pressure. Single beep.	The internal pressure of the concentrator is greater than 260 kPa or less than 20 kPa.	The alarm indicator lights up yellow and the error code "E05" or "E02" is displayed on the screen.	Turn off the concentrator and check and make sure that the air outlet and inlet are not blocked. Then restart your device. If the alarm is still active, contact the seller.
High temperature. Single beep.	The device system temperature is higher than the maximum allowed value.	The alarm indicator lights up yellow and the error code "E35" is displayed on the screen.	Turn off the concentrator and check and make sure that the air outlet and inlet are not blocked. Then restart your device. If the alarm is still active, contact the seller.

Alarm	Cause	Indicator	Possible solutions
Oxygen sensor communication failure. Single beep.	The oxygen concentrator is not receiving a signal from the oxygen sensor.	The alarm indicator lights up yellow and the error code "E31" is displayed on the screen.	Contact the seller. For those who need oxygen urgently, reserve oxygen should be prepared.
Low voltage. Single beep.	The power supply voltage of the device is 85% lower than the standard voltage.	The alarm indicator lights up yellow and the error code "E03" is displayed on the screen.	Disconnect the device from the power supply. Check that the mains power supply is normal and the cable is not damaged. If the alarm still occurs after reconnecting, contact the seller.

Please check the solutions below before reporting a problem.

Problem	Possible causes	Possible solutions
When turned on, the display, alarm system and oxygen concentrator do not work.	• Power plug not properly connected. • No mains power. • The fuse is broken.	• Check that the plug is pressed firmly into the socket. • Check the mains power. • Replace the fuse (see page 21).
The desired performance cannot be obtained.	• The nasal cannula is blocked or damaged. • The oxygen mask is blocked or damaged. • The humidifier bottle is blocked or damaged. • The oxygen tube is kinked.	• Remove the nasal cannula or mask, clean it and remove any kinks in the tube. • Replace the cannula or mask with a new one. • Remove the humidifier bottle, clean or replace it.

Troubleshooting

Please check the solutions below before reporting a problem.

Problem	Possible causes	Possible solutions
Nebulization is insufficient or not working properly.	• The nebulization set is installed incorrectly. • The nebulization port is blocked.	• Install the nebulization set correctly. • Replace nebulization set. If the problem persists, contact the seller.
When turned on, the display works, but the oxygen concentrator not function.	• Air compressor protection. • Blocked inlet or outlet. • The ambient temperature is below 5°C.	• Contact the seller if the concentrator does not work within 45 minutes of restarting. • Clear the filter. Check that the air inlet/outlet is not blocked. • Run the device at a higher ambient temperature.
The remote control doesn't work.	• The hub is out of range. • Battery run out.	• Reduce the distance to the hub (max. 5 m). • Replace the battery.

If the problem persists, contact the device seller.

Error table

Kod	Error cause
E02	During operation, the pressure drops below the lower pressure limit (20kPa).
E03	The power supply voltage of the device is 85% lower than the standard voltage.
E05	During operation, the pressure exceeds the upper pressure limit (260 kPa).
E31	Unable to receive data from the oxygen sensor.
E35	Compressor temperature control switch is not connected or exceeds the acceptable range.

If the problem persists, contact the device seller.

Warning!

It is strictly forbidden to open the concentrator housing for maintenance by non-professional personnel.

Pulse oximeter introduction

Overview:

1. Power On/Off switch:

- Turning on/off: Press the button to turn on the device. Press the button for 2 seconds to turn it off.
- Press this button while the machine is running to light up the display.

2. Bluetooth icon

- Visible: means that the oximeter is connected to a Bluetooth device.
- Not visible: means that no device is currently connected via Bluetooth.

3. Heart symbol

4. Pulse rate bar graph

5. PR value: Pulse rate per minute

6. Saturation value (SpO₂).

7. Battery capacity Indicator

-  - filling indicates battery level.
-  - the battery is discharged.

8. Micro USB connector: charging battery with specified charging cable connected to adaptor.

9. Rubber finger cover

Pulse oximeter safety notes and warnings

This device is not intended for treatment.

Do not disassemble.

This device is not intended for use on infants or neonates.

Do not use the device while using a defibrillator or other high-frequency equipment.

The presence of strong ambient light may cause inaccurate saturation (SpO₂) measurements.

Users who are allergic to rubber cannot use this product.

Explosion hazard: Do not use this device in an explosive or corrosive atmosphere.

In some circumstances, the device may interpret movement as a good quality heart rate. The patient's movements should be limited as much as possible.

Before using the pulse oximeter, carefully read this manual for all safety, operation, and specification information.

Measurement data displayed on the device are for information purposes only and cannot be directly used for diagnostic interpretation.

Do not use the device in situations where alarms are required. The pulse oximeter does not support alarms.

Do not connect to any external device or use for measurements while the device is charging.

Check and change the measurement site based on different circumstances, especially during long-term continuous monitoring. Check your finger every 2 hours. Change the measurement finger if unusual changes are seen. If the pulse oximeter is too tight because the measurement site is too large or becomes too large due to swelling, change finger or stop the measurement. Excessive pressure sustained for a long time can cause venous stasis, leading to interstitial edema and tissue ischemia.

Pulse oximeter precautions and warnings

-  Do not allow the device to get wet.
-  Protect from running water.
-  Do not clean the device with alcohol or solvents.
-  Do not use the device near sources of electromagnetic interference, such as mobile phones or radio transmitters.
-  Do not use the device if satisfactory results cannot be achieved.
-  Avoid placing the sensor on limbs with an arterial catheter or IV line.
-  Do not pierce or incinerate the device.
-  Do not use near an open flame.
-  Do not spray the device with insecticides.
-  Always handle the device carefully to avoid damage caused by drops, impacts, strong vibrations or other mechanical forces.
-  Check the device before use. Do not use it if there is significant damage.
-  Do not use the tester to access accurate saturation (SpO₂) data.

Pulse oximeter information

The pulse oximeter uses spectrophotometry by detecting blood redness and infrared light absorption to obtain oxygen saturation (SpO₂) and pulse rate (PR) values.

Oxygen saturation (blood saturation) is a term referring to the fraction of oxygen-saturated hemoglobin (HbO₂) relative to the total hemoglobin (Hb) in the blood, i.e. blood oxygen saturation, which is an important physiological parameter of the respiratory and circulatory systems. Many respiratory diseases can cause a decrease in oxygen saturation. The body's automatic regulatory dysfunction caused by anesthesia, major surgical injuries and injuries caused by some medical examinations can lead to decreased oxygen saturation, causing patients to experience dizziness, weakness, vomiting and other symptoms. Severe cases of low saturation can be life threatening. The value of saturation is of great importance in clinical medicine. Diagnosing the patient's blood oxygen saturation quickly will help the doctor find problems.

Intended use and contraindications

The pulse oximeter is designed to measure blood saturation and pulse rate of adults and children in home and outpatient settings.

Contraindications: the device cannot be used by users who are allergic to rubber.

Applications:

This device is suitable for the following users:

- people with vascular diseases, such as: coronary artery disease, hypertension, hyperlipidemia, cerebral thrombosis;
- people with respiratory diseases such as: asthma, tracheitis, chronic bronchitis, chronic cor pulmonale, chronic obstructive pulmonary disease;
- elderly over 60 years of age;
- people working more than 12 hours a day;
- people working in extreme environments or in environments that could lead to hypoxia;
- long-term alcoholics.

Pulse oximeter usage

1. Insert your finger fully under the pulse oximeter cover to obtain correct measurements.
2. Press the ON / OFF button to turn on the pulse oximeter.
3. The saturation (SpO₂) and pulse rate (PR) values will be displayed on the screen.

Notes:

- **The pulse oximeter cannot be used while charging.**
- There must be no nail polish or other cosmetics on the finger.
- Do not shake your finger, hand or body during measurement.
- Patches or other tapes around the site finger may affect the accuracy of the oxygen saturation and pulse rate measurements.
- The oximeter will turn off if no action is taken within 10 seconds of turning it on.

The following factors may influence the accuracy of measurements:

- Ambient light;
- Physical movement (patient or imposed motion);
- Diagnostic testing;
- Low perfusion;
- Electromagnetic interference, such as produced by MRI scanner;
- Electrosurgical units;
- Dysfunctional haemoglobin, such as carboxyhemoglobin (COHb) and methemoglobin (MetHb);
- The oximeter placed improperly;
- Drop of arterial blood flow to immeasurable level caused by shock, anemia, low temperature or vasoconstrictor;
- Dopamine, procaine, propranolol, lidocaine, bupivacaine and other drugs can cause more severe bias of SpO₂
- For anemic or toxic hypoxia, the SpO₂ value is only for reference, as some patients with severe anemia still show higher SpO₂ value.

Pulse oximeter charging

The oximeter uses a built-in non-removable rechargeable lithium battery. Symbol  will flash when battery capacity is low. Charge the battery according to the following steps.

1. Take out the cable and connect it to oximeter's Micro USB connector.
2. Plug the other end of cable to USB adaptor. Battery symbol appear and flashes on the screen while charging.

The performance of the rechargeable lithium battery may deteriorate over time. If the operating time of the battery is noticeably shorter than that stated in the specifications, contact the place where you bought the device.

Pulse oximeter cleaning and maintenance

Use only the approved substances and methods listed in this section to clean or disinfect the pulse oximeter.

We do not guarantee the effectiveness of the listed chemicals or methods as disinfectants. Consult your hygiene and epidemiology specialist for the most effective disinfection method

To avoid damaging the pulse oximeter, follow these rules:

- Always dilute substances according to instructions or use the lowest possible concentration.
- Do not immerse the pulse oximeter in liquid.
- Do not spill liquids on the pulse oximeter or accessories.
- Do not allow liquid to enter the oximeter housing.
- Never use abrasives (such as steel wool or silver polish) or erosive cleaners (such as acetone or acetone-based cleaners).

Pulse oximeter cleaning and maintenance

Recommended cleaning agents:

- Ethanol (75%)
- Isopropanol (70%)

Warning!

Turn off the pulse oximeter and disconnect the charging cable before cleaning.

Do not disinfect the pulse oximeter by heat, high pressure, gas vapors, or liquid immersion. Clean and disinfect the pulse oximeter as directed.

All servicing must be performed by an authorized service. The user or any other person should not disassemble or repair the pulse oximeter.

If liquid is spilled on the device or accessories, clean them immediately. If the pulse oximeter is not working properly, contact the seller.

Cleaning

- 1.** Make sure the pulse oximeter is turned off.
- 2.** Clean the exterior of the pulse oximeter with a soft cloth moistened with cleaning solution.
- 3.** If necessary, wipe off all the cleaning solution with a dry cloth after cleaning.
- 4.** Dry the pulse oximeter in a ventilated, cool place.

Disinfection

Disinfection may cause damage to the oximeter. We recommend that you disinfect the oximeter only when necessary. Clean the oximeter before disinfection. Recommended cleaning agent is Ethanol (75%) The recommended times of disinfection is three hundred.

Note: Do not use ethylene oxide or formaldehyde for disinfection.

Pulse oximeter specification

Model		F01
Operating mode		Continuous
Saturation (SpO₂%)	Measuring range	0% ~ 100%
	Measurement precision	70% ~ 100%, ±2%
	Measurement resolution	1%
Pulse rate (PR)	Measuring range	18 bpm ~ 250 bpm
	Measurement precision	25 bpm ~ 250 bpm, ±3 bpm
	Measurement resolution	1 bpm
	Refreshing rate	1 s
Control		• Control Panel
Display		Yes - LCD
Wavelengths		Red light: 660±3 nm Infrared light: 940±5 nm
Maximum power consumption of the converter		<30 mW
Working range - temperature		from 5 to 40°C
Working range - humidity		from 30% to 80%
Recommended max charged and discharged ambient temperature		Charged: 35°C, discharged: 45°C
Liquid protection level		IP22, non AP/APG type
Storage - temperature		from -20°C to +60°C
Storage - humidity		from 10% to 95%
Battery capacity		300 mAh
Charging voltage and current		5V ~ 160 mA
Charging connector		micro USB connector
Charge time		2.5 hours (charge to 90% of capacity)
Run time		20 hours (Fully-charged new battery, ambient temperature: 25°C)
Shutdown delay		max 5 minutes after the low battery prompt first occurs
Weight		< 50g
Dimensions		72 x 32 x 27 mm

Pulse oximeter troubleshooting

Before reporting a bug, check the solutions below.

Problem	Possible causes	Possible solutions
The oximeter cannot start.	- Lower battery capacity or no power, - Oximeter might be damaged.	- Charge the battery. - Check the power supply.
SpO ₂ or PR value cannot be shown normally.	- Finger might not be inserted properly, - High levels of ambient light, - User is in low perfusion or user's oxyhemoglobin is too low to be measured, - Oximeter might be damaged.	- Do not use the oximeter in the environment with high ambient light, - Go to a hospital for diagnosis.
SpO ₂ or PR value is unstable.	- Finger might not be inserted properly, - Finger is trembling or patient is moving.	- Insert finger properly. - Please keep still.
The display turns off suddenly.	No finger inserted after turning on the oximeter, or it automatically shut down caused by lowest battery capacity.	Press ON/OFF button to turn it on or charge the battery.

Manufactured for:

WELLTEC sp. z o.o.
Orzechowa 3, 72-010 Przęsocin, Poland
welltec.pro

Manufactured by:

Jiangsu Konsung Bio-Medical Science And Technology Co., Ltd.
NO.8, Shengchang West Road, Danyang Development Zone
Jiangsu Province, 212300, P.R. China

Shanghai International Holding Corp. GmbH (Europe)
Eiffestrasse 80, 20537 Hamburg, Germany

Warning! Do not dispose of the device in the household waste!

In accordance with the guidelines of the European Directive 2012/19 / EU on old, used electrical and electronic devices and its application in national law, unusable electrical equipment has been eliminated and electronic must be collected separately and recycled in an environmentally sound manner.

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