70*100mm

zenwire

Alcohol Breath Tester MR806

User Manual

Warnings

1. This product cannot be used to determine whether a person is fit to drive (e.g., drunk driving).
2. Results are for reference only. Even with the same alcohol concentration, the level of intoxication can vary depending on physical condition and individual differences.
3. Regardless of whether the breathalyzer is used, the manufacturer and the seller are not responsible for drunk driving or any legally prohibited behavior.
4. It takes about 20 minutes for alcohol to be absorbed from the digestive system into the blood. Therefore, measurements should be taken at least 20 minutes after drinking. Testing immediately after drinking may only reflect the alcohol present in the mouth, not the blood alcohol concentration.
5. For more accurate results, wait at least 3 minutes between repeated measurements.
6. Do not disassemble, modify, or wash the device. Store it in a place free of strong odors.
7. When the battery level is low, measurement results may have some error compared to actual values.
8. Keep the device out of direct sunlight, avoid excessive impact or vibration, and stay away from heat sources.
9. Do not use the device in environments with strong winds, smoke, or high humidity, as this may affect the accuracy of results.
10. Avoid eating, smoking, or using mouthwash within 15 minutes before testing, as residues may cause false readings.
11. This device is not a medical instrument and should not be used for health or treatment purposes.
12. Do not allow children to use the device unsupervised.
13. Clean the mouthpiece regularly and replace it if it becomes dirty or damaged to ensure accuracy.
14. Strong fragrances (perfume, air freshener, cleaning chemicals) near the sensor may interfere with measurement accuracy.
15. Always use high-quality batteries (3×AAA). Remove the batteries if the device will not be used for a long time to prevent leakage.
16. Do not immerse the device in water or expose it to rain.

Product Description

1	Blowing Nozzle
2	Indicator Light
3	Exhaust Port
4	LCD Display Screen
5	Power On/Off Button
6	Buzzer Sound Outlet
7	LED Flashlight

Usage Method

- 1. Battery Installation**
For first-time use, prepare 3 × 1.5V AAA dry batteries. Remove the battery cover, insert the batteries correctly according to the polarity markings, then replace the cover securely.
- 2. Power On**
Press the Power Button to turn on the device.
- 3. Preheating**
After startup, the device enters a preheating countdown (1–10 seconds). Please wait until preheating is completed.
- 4. Blow Detection**
After preheating, the device automatically enters the blow air detection countdown (10 seconds). Blow steadily into the nozzle when prompted. Wait for the test result to be displayed on the screen.
- 5. Set Display Unit**
After steps 2–4 are completed and the result appears on the screen (within 15 seconds), press the Power Button once to cycle through different units: %BAC (USA standard), g/L (European standard), mg/L BrAC (European/Japan standard)

- 6. Emergency LED Light**
Press the Power Button again while the device is off to turn on the emergency flashlight. Press again to turn it off.
- 7. Automatic Cleaning Mode**
If the time between two tests is less than 3 minutes, and residual alcohol is detected inside the chamber, the screen will flash and the device will automatically enter high-temperature cleaning mode. Once the display stops flashing and returns to normal, cleaning is complete and the device is ready for the next test.

Pay Attention

When the breathalyzer is first removed from the package or has not been used for a long period, gases and other impurities may accumulate on the sensor. In this case, press the Power Button several times to allow the device to preheat and clean itself before accurate use.

Measuring Method

- Preheating**
1. Press and hold the Power Button for 2 seconds until you hear a beep and the display turns on.
 2. The screen will show "Preheat" with a 10-second countdown.
 3. When the countdown reaches 0, preheating is complete.
- Testing**
1. After preheating, the device enters the "Blow" state with a 10-second countdown.
 2. Place your mouth on the nozzle and blow steadily for 4 seconds until you hear a beep. If blowing is too short or interrupted, the device will not measure properly.
 3. After a few seconds, the measurement result will be displayed. If no alcohol is detected, the result will show 0.

Measurement Results

- Green Light – Safe / Below Limit**
Blood Alcohol Concentration (BAC) is below the legal limit.
Conditions: (< 0.3 g/L, < 0.030% BAC, < 0.15 mg/L)
Indicator: Green light ON
- Yellow Light – Drinking Detected**
Indicates alcohol is present and driving may be impaired.
Conditions: (0.3–0.5 g/L, 0.030–0.080% BAC, 0.15–0.25 mg/L)
Indicator: Yellow light + Alarm Sound
- Red Light – Intoxicated**
Indicates intoxication / drunk driving level.
Conditions: (≥ 0.5 g/L, ≥ 0.080% BAC, 0.25 mg/L)
Indicator: Red light + Rapid Alarm Sound

Notification


1. Before first use, insert the batteries.
2. If the product is new or has not been used for a long time, the sensor must be cleaned. Follow these steps:
 - ① Press the Power Button
 - ② Wait without blowing into the nozzle
 - ③ The device will automatically power off when cleaning is complete. Repeat the above steps 3 times to ensure proper cleaning.
3. If the device will not be used for a long time, remove the batteries to prevent leakage.

Measurement Units

The breathalyzer supports different units commonly used in various countries: mg/100 mL → Commonly used in China, %BAC → Commonly used in the USA, g/L → Commonly used in Europe, mg/L BrAC → Commonly used in Japan
Example: 0.05% BAC means 0.05 g (50 mg) of alcohol per 100 mL of blood.

	BAC Alcohol concentration in the blood				BrAC Alcohol concentration in the exhaled breath	
	mg/100ml	%	‰	g/l	mg/l	PPM
1	10	0.01	0.1	0.1	0.05	25
2	20	0.02	0.2	0.2	0.10	50
3	30	0.03	0.3	0.3	0.15	75
4	40	0.04	0.4	0.4	0.20	100
5	50	0.05	0.5	0.5	0.25	125
6	60	0.06	0.6	0.6	0.30	150
7	70	0.07	0.7	0.7	0.35	175
8	80	0.08	0.8	0.8	0.40	200
9	90	0.09	0.9	0.9	0.45	225
10	100	0.10	1.0	1.0	0.50	250
11	110	0.11	1.1	1.1	0.55	275
12	120	0.12	1.2	1.2	0.60	300
13	130	0.13	1.3	1.3	0.65	325
14	140	0.14	1.4	1.4	0.70	350
15	150	0.15	1.5	1.5	0.75	375
16	160	0.16	1.6	1.6	0.80	400
17	170	0.17	1.7	1.7	0.85	425
18	180	0.18	1.8	1.8	0.90	450
19	190	0.19	1.9	1.9	0.95	475
20	200	0.20	2.0	2.0	1.00	500

E-mail: sklep@zenwire.eu
Manufacturer: Zenwire Sp.z.o.o. Kolobrzaska 4, Poznań 60-418, Made in PRC.

Business users in the European Union

If you need to dispose of electrical or electronic devices, please contact your nearest point of sale or supplier for additional information.

Disposal of waste outside the European Union

This symbol is only valid in the European Union. If you need to dispose of his product, please contact local authorities or the seller for information on the correct method of disposal.

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