# PL

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# EN

## Instructions for use

Sigvaris Aesthetic Medical Care VEST  
Compression Garments – Vest

### INTENDED USE:

A medical device with a compression and a stabilising effect, supporting the healing process, intended for use immediately after plastic, aesthetic and reconstructive surgeries and/or during the postoperative recovery period, as indicated by a physician.

**TARGET GROUP:** Device designed for persons showing the symptoms described in the indications.

### INDICATIONS:

- after a gynecostasia surgery:
- after oncological surgeries
- bilateral mastectomy
- other male breast surgeries

### CONTRAINDICATIONS:

- obstructive and restrictive diseases of the respiratory system, reduced respiratory volume conditions and reduced respiratory efficiency
- breast wound healing disorders, dehiscence of the surgical wound on the chest
- skin lesions with infection of the skin around the breast and chest area

### SAFETY PRECAUTIONS:

Do not use the product directly on wounds, ulcers, or damaged, afflicted or inflamed skin. Before using the product, treat inflammation and damaged or afflicted skin. In the case of postoperative use, a wound, an incision, or a skin puncture site should be covered with a dressing, as recommended by a physician or qualified medical personnel. The dressing on the wound (puncture, in the case of liposuction, or incision, in the case of surgery) should remain on it as protection until the wound heals or at least until secretion stops. Most often, an absorbent dressing or layer dressing is used on a wound, consisting of a non-stick layer applied directly on the skin, followed by an absorbent layer and, alternatively, a liquid-impermeable layer. A dressing on a wound should be changed, replaced or supplemented as recommended by a physician or qualified medical personnel. A dressing can be attached to the skin outside a wound to limit its displacement and relocation during movement and use. In the case of persistent or intensified pain, sensory disorders, an allergic reaction, itch, redness, skin maceration, irritation or any other adverse effect related to the use of the device, discontinue use and consult a doctor. The appropriate use of the product may require the advice of specialist medical staff. The place of application, and the location of the product, should be in accordance with the intended purpose. The product should not be applied on the neck or face (except for a face mask designed for this purpose); the product should not be tied around the neck, face or limbs. When applying a product with metal elements, such as hooks and zip fasteners, please make sure that these elements are not applied directly on the skin or pressed by the product, as it may result in point pressure, which leads to the formation of wounds, bedsores, or maceration. Direct contact of metal with the skin may also lead to the development of allergic contact dermatitis or skin allergy in sensitive individuals. The selection of the correct size determines the proper functioning of the product. If the product is too small in size or too tight, it may lead to circulatory disorders, especially to the impaired outflow of venous blood or impaired lymph movement and, as a result, to oedemas. In extreme cases, circular pressure may disturb the circulation by blocking the arterial inflow. Great caution is required when using the product immediately after surgery in patients under the influence of anaesthetic agents used for general anaesthesia, strong pain medication, and after using local or regional anaesthesia due to the lack of protective sensation. Great caution is required when using the product in patients with sensory disturbances, neuropathy and diabetes (diabetic polyneuropathy), polyneuropathies of other origins, nerve compression syndromes, arterial diseases (diabetic angiopathy, arteritis), certain skin diseases, e.g. vasculitis, connective tissue diseases, confirmed allergic contact dermatitis to a product component. Each use of the product should be preceded by verification of the signs of wear or ageing. Only an undamaged product guarantees optimal performance. Do not use the device continuously for 24 hours unless otherwise advised by your doctor. The use of ointments, oils, creams and similar substances may have a negative effect on the product's properties, reduce its durability and effectiveness and, in rare cases, may lead to the development of allergic contact dermatitis.

**INFORMATION ON INITIAL PREPARATION OF THE PRODUCT:** It is recommended to wash the product before the first use.

**CARE:** Wash in the washing bag; Before washing (if present), fasten the velcro and/or remove the straps.

**SIZE DETERMINATION:** Size chart is available on the packaging.

### METHOD OF PUTTING ON:

Undo the zipper and the hooks and clasps below. Put the product on: pass your arms through the holes and place the product so that it is around the chest with the hooks and clasps and the zipper at the front. Secure the hooks and clasps first, then it zip up. Pull the lower hem to stretch the product downward. Adjust the straps.

### DISPOSAL:

Throw the product away following the local rules for getting rid of waste.

**All serious incidents that occur when using the product shall be reported to the manufacturer and the competent supervisory authority.**

**CH** | SIGVARIS S.A., ul. Mazowiecka 11 lok. 49, 00-052 Warszawa, Poland  
**CH** | REP | SIGVARIS AG, Gröblistrasse 8, CH-9014 St.Gallen/Switzerland

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SIGVARIS, SIGVARIS GROUP and MOBILIS are registered trademarks of SIGVARIS AG, CH-9014 St.Gallen/Switzerland, in many countries worldwide. © 2025 Copyright by SIGVARIS AG, CH-9014 St.Gallen/Switzerland V3 / 07.2025
**CH** | SIGVARIS S.A., ul. Mazowiecka 11 lok. 49, 00-052 Warszawa, Pologne
**CH** | REP | SIGVARIS AG, Gröblistrasse 8, CH-9014 St.Gallen/Suisse