DRE Automatic

Analog Gas Manifold

Submittal Data Sheet

Project Information

Project Number	Approval



SPECIFICATIONS

The fully automatic manifold shall be a TMC series. No manual resetting of valves or levers shall be required. The unit shall switch from "Bank in Use" to "Reserve" bank without fluctuation in line delivery pressure. Simultaneously, the "Reserve in Use" alarm shall be triggered by the manifolds circuit board. The manifold shall continue to provide gas, in the event of a power failure, until both banks are depleted. After the switchover, the "Reserve" bank shall then become the "Bank in Use". The manifold shall be capable of being upgraded after installation; to be used with low or medium pressure portable bulk vessels, to upgrade to high flow line regulator(s), from single to dual line regulators and for use at higher or lower delivery pressures.

The control panel shall incorporate a set of LED's for each bank, green for "Bank in Use", amber for "Ready" and red for "Empty". Analog gauges are also provided so that line and both bank pressures may be observed.

All manifold regulators, piping and control switching equipment shall be cleaned for use with oxygen service and installed in a steel cabinet to provide protection and minimize tampering.

FEATURES

- > Five year parts and one year labor limited warranty*
- > Fully automatic no resetting of valves or levers
- > Input power 120 to 240 VAC, 50 to 60 Hz single point connection.
- > Easy to service piping design
- > Patented single solenoid pressure dierential changeover
- > 400 psi dierential rated solenoid
- May be converted from high pressure cylinder use to use with low or medium pressure liquid portable bulk vessels
- Includes 3/4" source or main line ball valve with copper tube extension
- > Dual line pressure regulators
- > Optional single point relief valve vertical kit
- > Double "Z" brackets for one man installation.
- > Cabinet weight 70 lbs



Maximum rated flow capacity of line regulators only, not the manifold cabinet, flowing to atmosphere. (Without restricting line pressure drop)

Gas Service		Standard Line Regulators	High Capacity Line Regulators
Oxygen or Medical Air	Delivery Pressure and Flow Option	1L	1H, 2H, 3H
		2,500 SCFH (1,180 l/min)	4,500 SCFH (2,120 l/min)
Nitrogen	Delivery Pressure and Flow Option		3H
		N/A	6,000 SCFH (2,830 l/min)

Maximum recommended flow due to the chill down nature of the gas.

Gas Service		Without Heaters	With Heaters
	Delivery Pressure and Flow Option	1L	1W, 1X
Nitrous Oxide or Carbon Dioxide		40 SCFH	500 SCFH
		(19 l/min)	(236 l/min)

Note


1W - Models incorporate 1L Line Regulators 1X - Models incorporate 1H Line Regulators

Manifold Cabinet Flow Capacity

			Manifold Line Regulator Delivery Pressure and Flow Option Average Flow Rate in SCFH (l/min)			
Static Delivery Pressure Setting PSI	Pressure Drop	Pressure Flowing psi				
			1L	1H	2H	3H
	3	50	195 (92 l/min)	640 (302 l/min)		
53	5	48	430 (203 l/min)	1,260 (595 l/min)		
	7	46	635 (300 l/min)	1,650 (779 l/min)		
	10	43	875 (413 l/min)	2,430 (1,147 l/min)		
	3	82			1,010 (477 l/min)	
85	5	80			1,610 (760 l/min)	
	7	78			2,670 (1,261 l/min)	
	10	75			3,120 (1,473 l/min)	
	10	165				1,230 (581 l/min)
175	20	155				2,535 (1,197 l/min)
1/3	30	145				4,140 (1,955 l/min)
	35	140				4,500 (2,125 l/min)

Flow rates shown were obtained using Nitrogen, flowing through the right primary regulator, which is considered the most restrictive flow path. (Worst case condition). Testing was performed with an average inlet pressure to the manifold cabinet at 1,425 psi.





Design Lengths	Total # of Cylinders	4	6	8	10	12	16	20
STAGGERED DESIGN (5" CENTERS)		4'-6"	5'-4"	6'-2"	7'-0"	7'-10"	9'-6"	11'-2"
OVERALL MANIFOLD LENGTH		(1.32m)	(1.63m)	(1.88m)	(2.13m)	(2.39m)	(2.90m)	(3.33m)
VERTICAL CROSSOVER (5" CENTERS)		3'-7"		4'-6"		5'-4"	8'-7"	10'-3"
OVERALL N	MANIFOLD LENGTH	(1.10m)	N/A	(1.32m)	N/A	(1.63m)	(2.62m)	(3.12m)