


User Instruction Sheet

Select the suitable glove size using the label. Before usage, inspect the gloves for any defects. If the gloves are defective, dispose them immediately. If in doubt do not use the gloves, get a new pair of gloves.

EN 420:2003+A1:2009

- Name and full address of manufacturer or his authorized representative.
Advansafety, Biyagama Export Processing Zone – A, Walgama, Malwana, Sri Lanka
- Name & address of the notified body responsible for both EU Type Examination and on-going conformity:
SATRA Technology Europe Ltd., Bracetown Business Park, Clonee, Dublin 15, D15 YN2P, Ireland, NB No.: 2777
- Glove designation (name or reference):
240-050BS - Latex, Powdered, 5g Disposable Gloves. (White, Green, Yellow, Pink, Purple, Blue)
- Information on the available size range Full Dipped Sizes: 7,8,9,10
- Reference to the relevant specific European standards
EN 420:2003+A1:2009
EN ISO 374-1:2016 + A1:2018
EN16523-1:2015
EN 374-4:2013
EN ISO 374-5:2016

Pictogram followed by the performance levels.



Certified performance level of the product as follows. Product is considered to be Category III of PPE hand protection and certified in accordance with PPE regulation (EU) 2016/425

Requirements for EN 374-1:2016 levels of performance

Levels of performance EN 374 - 1:2016

Level	Measured breakthrough time (minutes)
1	>10
2	>30
3	>60
4	>120
5	>240
6	>480

EN 374-1:2016

Method	Chemical Name	Performance level
EN 16523-1:2015	Sodium hydroxide 40% (K) 1310-73-2	1
EN 16523-1:2015	Hydrogen Peroxide 30% (P) 7722-84-1	1

EN 374-4 Degradation results

Sample description:		LRP/29.03.2019/374/4 and LRP/29.03.2019/420/01		
Challenge Chemical:	30% Sodium hydroxide (CAS: 7722-84-1)			
Test Temperature / OC:	(23 ± 1)			
Degradation / %:	Glove 1	Glove 2	Glove 3	
	11.8	22.0	4.2	
Mean degradation (DR) / %:	12.7			
Standard degradation (sDR) / %:	9.0			
UoM / ± %:	27.8			
Appearance of samples after testing:	Swollen and discoloured			

Sample description:		LRP/29.03.2019/374/4 and LRP/29.03.2019/420/01		
Challenge Chemical:	40% Sodium hydroxide (CAS: 1310-73-2)			
Test Temperature / OC:	(23 ± 1)			
Degradation / %:	Glove 1	Glove 2	Glove 3	
	-12.0	1.5	32.3	
Mean degradation (DR) / %:	-14.3			
Standard degradation (sDR) / %:	17.0			
UoM / ± %:	25.4			
Appearance of samples after testing:	No change			

- Degradation results indicate the change in puncture resistance of the gloves after the exposure to the challenge chemical

User Instruction Sheet

6. Basic explanation: The end user needs to know what the levels means under the pictogram.

Protection is claimed for the whole glove.

7. Product does not contain any known substances that may cause harm to the wearer's health and since this is a latex product there may be allergies in some individuals.

8. Instructions for decontamination: When reusing the gloves, it is recommended to dispose after use that day. Gloves may be cleaned and rinsed while being worn.

9. Disposal: Treat contaminated use gloves as bio-hazard and to be disposed professionally. Information can be obtained from the relevant waste disposal authorities. After carrying into contact with chemicals dispose the product in accordance with disposal regulation for the relevant chemicals.

10. Declaration of conformity can be viewed by visiting this link. <http://www.advansafety.online>

11. Warning:

a. Gloves not be worn when there is a risk of entanglement by moving parts of machines.

b. This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals.

c. The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400mm-where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical is used in a mixture.

d. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation.

e. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact may reduce. The actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves.

f. Before usage, inspect the gloves for any defect or imperfections.

g. EN ISO 374-5:2016 : The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

h. Explanation for EN ISO 374-5:2016 results: Protection against bacteria and fungi – Pass

Protection against viruses – Not tested

12. Storage & maintenance:

The gloves should be stored in the original packing at a dry and clean place. Please avoid exposing the glove to high temperature, humidity or direct sunlight light. Corrugated boxes with the inside poly bags should be used when transportation.

13. Shelf life of product:

Recommended shelf life of product is 5 years under controlled environment condition, product should store away from direct sunlight & away from humidity. Temperature of warehouse should not exceed more than 35 degree Celsius, product should store under proper packaging.

Declaration of conformity is available at www.advansafety.online

Manu. Date Xx/xx/xxxx
Exp. Date Xx/xx/xxxx