

System 98XT

The System 98XT IABP with CardioSync 2™ featuring R-Trac™

Datascope introduces the System 98XT with the autotiming supervisor, “R-Trac”, a highly advanced second generation IAB timing control system. It combines the sophistication of CardioSync 2 with rapid QRS identification and the fastest IAB inflation and deflation speeds available*.



For the Patient...

- R-Trac™ automatically matches the most appropriate deflation timing method to the patient's rhythm.
- CardioSync 2 Software and Fast Pneumatics provide highly accurate and reliable ventricular support, enhancing augmentation and LV unloading over a broad range of transient and sustained rhythm disturbances.
- CardioSync 2 improves tracking in ECG with automated management of ectopic patterns, including isolated PVC's, bigeminy, and couplets. Timing rapidly adapts to sudden changes in heart rate.
- CardioSync 2 delivers reliable tracking in pressure trigger during atrial fibrillation

For the Clinician...

- Automated controls and a significantly larger, brighter display enhance ease-of-use:
 - Rapid start-up function
 - Single ECG trigger mode
 - Advanced pressure trigger:
 - Beat-to-beat threshold adaptation
 - Auto reassessment of timing
 - Display and printer preferences menu
 - Keypad selectable signal sources
 - Advanced information display:
 - Auto R-wave deflation mode advisory
 - Auto scaling and positioning of pressure waveform
 - On-screen helium/battery indicators

* Inflation/deflation speed measured from command to 90% IAB inflation/deflation.



Datascope®

innovation is the best medicine

SYSTEM 98XT WITH CARDIOSYNC 2 SUMMARY TECHNICAL SPECIFICATIONS

EL Display: 8.3" (21.1 cm)W x 6.2" (15.8 cm)H; 160° viewing angle; Rotates 330°; Tilts 180°; Detachable; Laptop-like closure for storage and protection; Remote monitor mount (optional)

Preferences Menu: User may select display *sweep speed* (25 or 50 mm/sec), *brightness* (low, med., high); *balloon waveform* (on/off); *ECG inflation markers* (on/off); *AP waveform auto-scaling* (on/off); *flashing alarms* (on/off)

ECG ECG Trigger: Threshold dynamically adjusted by system for improved sensitivity and selectivity of the R-wave detection; Minimum = 120µV ± 20µV at normal gain; 40µV at max. gain

Pressure Trigger: *In automatic mode:* adjusted to 50% of the difference between peak systolic and end diastolic (avg. over multiple cycles); *In variable mode:* User adjustable between 7 and 30 mmHg ± 3 mmHg

Pacer A Trigger: R-wave detection (as above) except blanking is extended to 100 ms

Pacer V/A-V Trigger: V Pacer: fixed at rate up to 185 bpm (no demand pacing)

A-V Pacer: fixed at rate up to 125 bpm (no demand pacing) with A-V intervals between 80-224 ms

Internal Trigger: Variable mode: 40-120 bpm; Normal mode: 80 ± 1 bpm

Tall T-Wave Rejection: Rejects all T-Waves where Q-T interval is <300 ms and the amplitude is <70% of QRS input amplitude (ECG and Pacer A mode)

Pacer Rejection: Rejects all pulses of amplitude ± 2.0 mV to ± 700 mV and durations between 0.1 ms to 2.0 ms with:
1) No tail; 2) 100 ms time constant tail < 1 mV;
3) 25 ms time constant tail < 1 mV; 4) 4 ms time constant tail < 2 mV

ECG Leads: I, II, III, AVR, AVL, AVF, V (12 lead compatibility)
On-keypad switch between direct leads and external monitor

ECG Gain (default): 1 V output per 1 mV input ±5% (waveform automatically scaled to occupy ECG display window)

Gain (variable): 0.15 cm to 3.0 cm/mV ±20% (autoscaling disabled)

Frequency Response: 0.5-12 Hz (display); 0.5-135 Hz (Output to Monitor)

Defibrillator Protection: Discharge levels ≤360 J (trace returns to screen in 2 sec. max)

ESIS: Automatic suppression with internal ECG amplifier

AP Source Selection: On-keypad switch (transducer/external)

AP Input Sensitivity: 5.0µV/V/mm Hg

Excitation: +5 VDC ± 5% @ 130 mA (max)

Frequency Response: 0.5 to 12 Hz (-3 dB max atten.)

Timing Modes: Automatic (R-Trac ON or OFF); Manual

ECG Connector: 6-pin 5-lead AAMI standard

Pressure Connector: 6-pin AAMI standard

External Mon. Inputs: ECG: 1 V/mV (nominal); AP: 1 V/100 mm Hg (nominal)

ECG/BP Output: Phone Jack: ECG 1V / 1 mV, BP 1V/100 mm Hg

Mains Voltage/Frequency: 100-120 VAC ± 10% or 220-240 VAC ± 10% ; 50/60 Hz ± 3 Hz

Internal Battery: 24 VDC (nominal), 17.2 Amp-hour, 2.25 hrs. (min. @ 120 bpm)

Type: Maintenance free; Sealed lead-acid

System Compressor: Dual head diaphragm pump with brushless DC motor

IAB Helium (He*): Medical-grade; 0.5 L (2900 psi) or 0.69 L (2200 psi)

He+ Endurance (nominal): 1.8 mo. (.5 L); 2 mo. (0.69 L): continuous operation - 24 hrs./day

Condensate Removal: Fully automatic condensate removal and disposal

PC-IABP (optional): Software for remote clinical assistance and training

Service Diagnostics: Built-in software for system analysis and troubleshooting

Modem Data Rate: Up to 9,600 baud

Modem Certifications: Registered with FCC and accepted in ≥89 countries

Size on Cart: 43.3"H x 16.8"W x 22.3"D (109 cm x 42.7 cm x 56.6 cm)

Size off Cart: 26.9"H x 10.8"W x 20.5"D (68.3 cm x 27.4 cm x 52.1 cm)

Console Weight: 76.8 lbs. (34.8 Kg)

Monitor: 9.45 lbs. (4.3 Kg)

Hospital Cart Weight: 61.0 lbs. (27.7 Kg)

Internal Battery: 34 lbs. (15.4 Kg)

UTS Version: Rugged base which attaches to docking station (DS)

Lightweight DS: 17.8" (45.1 cm) x 20.0" (50.8 cm)
(Optional)

Basic DS: 21.5" (54.6 cm) x 21.5" (54.6 cm)

Ext Battery: 24 VDC (nominal), 17.2 Amp-hour, 1.8 hrs. (min. @ 120 bpm)
(Optional)

Aux. DC Input: 20.5-35.0 VDC

Shock and Vibration: ECRI AIII.3.1, 3.2.2, 3.3, 3.4 Mil Std 810E method 514.4, Cat. I; RTCA DO-160C, 1989, Section B, Curve N

Op. Temp.: 10°C - 40°C

Op. Humidity: 5 - 95% (R.H.) non-condensing

Op. Altitude: Up to 12,000 feet (3,657 m); automatic altitude correction for IAB pressure

Printer Type: Thermal array - 50 mm wide

Printer Menu: On-screen selection of waveforms; strip length; timed print; print on alarm; trend data; and alarm/trigger log

Waveforms (any 1 or 2): ECG (lead number & size labeled), arterial pressure, BPW

Event Markers: Trigger pulse (upper), inflation interval

Trend Graphs: ≤8 hrs of heart rate; peak systolic, end diastolic, mean, augmented diastolic

Additional Data: Space for patient I.D.; Date; Time; Alarm and Advisory messages; Timing mode; Trigger; HR, Frequency; Syst/Dias (in 1:1); Assisted and Unassisted Syst/Dias (in 1:2 or 1:3); Augmented Pressure; Mean

Doppler: 8 MHz non-directional probe

Saline Pole: Height adjustable

Simulator*: Emulates ventricular rhythms; Pacers; HR; A-fib


Remote Mount*: Mounts keypad/monitor on bypass pump

Storage Case*: 3-pocket design holds cables and spare items

Pump Cover*: Padded vinyl with pocket

* (optional)

- IEC601-1:1988/EN60601-1:1990
- EC Medical Device Directive 93/42/EEC
- IEC 601-2-27:1994 and IEC 601-2-34:1994
- UL 2601-1:1997
- CSA C22.2 No.601.1-M90



Metrics at 25°C , unless specified

Datascope Corp.
15 Law Drive
Fairfield, NJ 07004
800.777.4222 (tel.)
973.244.6100 (tel.)
800.258.8762 (fax)
973.244.6299 (fax)