



Snow / Siker Mirror / Whitehead Blades

American Profile

Laryngoscopes

With Fiber Optic Illumination

SNOW


A modification which includes the Miller spatula and the Wis-Foregger (Wisconsin modification) flange to reduce trauma during intubation, especially for use with the Sanders tube. Curving back from the distal end, it provides a convenient visual field.

SIKER MIRROR

For laryngoscopy of patients with anatomical variations which make intubation with conventional blades difficult or impossible. The basic portion of the blade makes an angle of 135° with a stainless steel mirror located at that angle on the flange facing the spatula. The distal portion of the blade is three inches long. Because the mirror inverts the reflected image, the operator should exercise caution until familiar with its use. A curved stylette is recommended for use with the Siker blade. Placing the blade in a warm water bath for ten minutes before use will avoid any mirror fogging problems.

WHITEHEAD MODIFICATION OF WISCONSIN

A modification of the Wisconsin blade with a reduced flange to increase visualization, facilitate intubation and reduce possibility of damage to the incisors.



CAT#	SIZE	DESCRIPTION	A	B
5-5249-03	3	Medium Adult	162mm	16mm



CAT#	SIZE	DESCRIPTION	A	B
5-5245-04	4	Large Adult	158mm	18mm



CAT#	SIZE	DESCRIPTION	A	B
5-5244-03	3	Medium Adult	162mm	15mm

All above blades are made from 303/304 surgical stainless steel and are provided with 4mm, replaceable fiber optic bundles. GreenLine® blades and handles are fully compliant with ISO standard 7376-3 and interchange perfectly with companies conforming to the "Green System". (i.e. Anesthesia Associates, AMS, Heine, Propper, Rusch and Welch Allyn)