Penlon Prima SP2 Anesthesia System

## The affordable, low maintenance anesthesia workstation.

The Prima SP2 Anesthesia Machine provides sophistication in a flexible, easy-to-use system, designed to meet all your anesthesia requirements.





## The durable, easy-to-use Penlon Prima SP2 is designed for you.

The Penlon Prima SP2 has an open architecture, is perfect for ultra-low flow, and is backed by a manufacturer that has never discontinued service for its anesthesia machines or ventilators.

- Balanced ultra-low flow anesthesia that provides increased patient safety with a mechanical anti-hypoxic device
- Modular construction with high levels of workstation integration
- Open architecture gives you ample space for monitors and accessories
- Standard with three gases (O<sub>2</sub> N<sub>2</sub>O Air)
- Available with a pull-out writing tablet









The fully integrated Penlon Prima SP2 features the AV-S touchscreen ventilator, the high performance A200SP absorber, and Sigma Delta vaporizer(s)

With the Penlon Prima SP2 you get high performance and a low cost of ownership.

#### **Feature-packed**

system includes a touchscreen ventilator with three spontanous breathing modes and ultra-low flow delivery Easy to use with intuitive features for a smooth, pain-free transition from your old

# Two-year manufacturer's parts warranty

protects your investment and reinforces peace of mind

# LOADED WITH ESSENTIAL FEATURES that can MAXIMIZE

**PERFORMANCE** 

# LOW COST OF OWNERSHIP HELPS REDUCE ANNUAL COSTS

# Save on annual service

with only two preventative maintenance visits per year

# Integrated spirometry

(not consumable) can save you as much as \$720 per year

#### Reduced software upgrade fees

help you save money throughout the life of your machine

# Get the features you need every day (and with the Prima SP2, you won't have to worry about a discontinuation letter looming around the corner).

The superior Prima SP is constructed of durable components (more aluminum, less flimsy plastic) and comes with standard features (like three gases) for which competitors charge extra. See below for how the Prima SP2 stacks up against competing machines.







Feature	Penlon Prima SP2	Fabius GS	Aestiva/5
Dimensions (h x d x w)	54" x 29" x 31.5" + abs.	51" x 32" x 35" + abs.	54" x 33" x 37" + abs.
Weight (lbs)	220 w/ 2 vaps.	260 w/ 2 vaps.	387 w/ 3 vaps.
Top Shelf (w x d)	28" x 19.5"	30" x 19"	26" x 16"
Work Surface (w x d)	25.6" x 11.9"	27" x 13"	18" x 12"
Writing Tablet	Optional	Yes	Optional
Drawers ( $h x d x w$ )	Up to 3 (7.5" x 20.5" x 15.4")	3	1 (4" x 10" x 15")
Flow Meter Type	Dual cascade flow tubes	LED and combined flow	Dual flow tubes
Gases (Standard)	$O_2 \cdot N_2 0 \cdot Air$	$O_2 \cdot N_2 0 \cdot Air$	$O_2 \cdot N_2 0$ (air is optional)
Cylinders (Standard)	Four	Optional	Optional
Electrical Outlets x 4	Standard	Optional	Optional
Anti-Hypoxic Device	Mechanical, 27 to 33%	S-ORC	Chain link – 25

## The easy-to-use AV-S Ventilator is designed for all patient types.



The AV-S anesthesia ventilator features a user-friendly touchscreen and comprehensive support modes that provide optimum therapy for all patient profiles. With the clearly-labeled, simple layout of the AV-S interface, you'll be able to make quick decisions on the fly, even in the most hectic operating environment. You can also easily configure its screen to display two waveforms, including a compliance loop not available on many competing ventilators.

5

Feature	AV-S	E-Vent	7900	7100
Screen Type/Size	Color Touchscreen/8.4"	EL Black and Orange/6.4"	EL Black and Orange/6.4"	LCD/5.5"
Volume	Yes	Yes	Yes	Yes
Pressure	Yes	Yes	Yes	Optional
Spontaneous Monitoring	Yes	Yes	Yes	Yes
Pressure Support	Yes	No	Optional	Optional
Control	Touchscreen and Com-Wheel	Com-Wheel only	Com-Wheel only	Com-Wheel only
Dual Waveform	Yes	No	No	No
Print From Screen Ability	Yes	No	No	No
Tidal Volume (ml)	20–1600	20–1400	20–1500	45–1500
Frequency/min.	4–100	4–60	4–100	4–65
I:E Ratio	3:1 to 1:8	4:1 to 1:4	2:1 to 1:8	2:1 to 1:6
Pressure v. Limits	10–80	15–70	12–100	12–100
Electronic PEEP cmH <sub>2</sub> 0	4–30	0–20	4–30	4–30
Fresh Gas Compensation	Yes	Decoupling	Yes	Yes
Gas Powered	Yes	No, Electric	Yes	Yes
Battery Back-up	45 min.	45 min.	45 min.	45 min.
Oxygen Monitor	Yes	Yes	Yes	Yes



## The AV-S features three advanced spontaneous breathing modes.




The AV-S ventilator provides three support modes that can be utilized as the patient attempts to breathe spontaneously. Patient recovery is accelerated by increased tidal volume and SpO<sub>2</sub>, and reduced EtCO<sub>2</sub>.

**SIMV** (Synchronized Intermittent Mandatory Ventilation) Provides a minimum level of volume, and allows spontaneous breaths and set mandatory breaths.


#### SIMV — Spontaneously Breathing Patient

Negative pressure<sup>†</sup> in the Trigger Window\* (generated by the patient's spontaneous breath) results in a synchronized mandatory breath at a preset volume and rate.



#### SIMV — Non-breathing Patient


If the patient makes no effort to breathe during a cycle, a mandatory breath, at the end of the Trigger Window,\* will still be delivered at the preset volume and rate.



**SMMV** (Synchronized Mandatory Minute Ventilation)
Provides a set level of minute volume ventilation and allows spontaneous breaths, combined with a synchronized mandatory breath, to achieve the set minute volume.


#### SMMV — Spontaneously Breathing Patient

Negative pressure<sup>†</sup> in the Trigger Window\* (generated by the patient's spontaneous breath) results in a synchronized mandatory breath, ensuring that the set minute volume is achieved.



#### **SMMV** — Non-breathing Patient

If the patient makes no effort to breathe during a cycle, a mandatory breath, at the end of the Trigger Window,\* will still be delivered at the preset volume and rate.




#### **PSV** (Pressure Supported Ventilation)

PSV assists each spontaneous breath with a preset pressure reducing the effort required to breathe. Negative pressure<sup>†</sup> (generated by the patient's spontaneous breath) results in synchronized pressure support.

### PSV is used to support spontaneously breathing patients only.

If the patient makes no attempt to breathe, the ventilator will not provide support and apnea alarm will be activated



<sup>†</sup>Negative relative to PEEP \*Trigger Window = 60% of Cycle Time



# The A200SP Absorber combines ease of use with advanced system integration.



The high performance A200SP absorber provides you with an efficient  $\mathrm{CO}_2$  absorption module and integrated, permanent spirometry sensors that don't have to be replaced year after year. It's autoclavable<sup>†</sup> and it can easily be configured for right-handed and left-handed operators. You can also utilize its quick-release canister to change absorbent while ventilating your patient.

Feature	A200SP	Cosy	Absorber	ABS
Heated Circuit	Optional	Optional	No	No
Large, Easy-to-Reach Controls	Yes	No	No	No
Integrated Ascending Bellows	Yes	No Bellows	Yes	Yes
Vent Bag/Vent Auto Switching	Yes	No	Yes	Yes
Integrated O <sub>2</sub> Sensor	Yes	Yes	Yes	Yes
Single Canister, Low Compliance	Yes	Yes	No	Yes
Quick Canister Release w/ Auto Seal	Yes	Yes	No	Yes
Loose Absorbent or Standard Pre-Pack	Yes	No	Yes	No
Autoclavable <sup>+</sup>	Yes	Yes	Yes	Yes
Robust Arm Bag	Standard	Optional	Optional	Optional
Pressure Gauge	Standard	Optional	Optional	Optional

†Excludes covers, manometer and oxygen sensor.



# The lightweight Sigma Delta Vaporizer delivers accurate concentrations at ultra-low flows.



Count on the Sigma Delta vaporizer to deliver accurate concentrations at varying temperatures and vapor concentrations, particularly at low flows. It is service-free for as many as 10 years, and at only 11 pounds, the lightweight Sigma Delta is much easier to handle than rival vaporizers.

Feature	Sigma Delta	Vapor 2000	Tec 5	Tec 7
Weight (lbs)	11	16	15	15
Back Bars	Selectatec, North American Drager, Cage	Selectatec, Cage	Selectatec, Cage Selectatec	
Major Service Intervals	10 Years	10 Years	3 Years	10 Years
In-House Service by User	Optional	No	No No	
3-Year Warranty	Yes	Yes	Yes	Yes
Ergonomic Push/Turn Control Dial	Yes	No	No	No



#### Technical Specifications

n: con a di con		AVCAR IN A L W CL			ADDOCT AL.	
Prima SP2 Anesthesia System		AV-S Multi-Mode Ventilator continued			A200SP Absorber	
Physical Size (H x W x D)	54" x 28.4" x 31.5"	SIMV, SMMV, PSV Trigger		0.2 to 0.4 L/min	<b>Physical</b> Size (H x W x D)	16.6" x 9.0" x 16.9"
312e (H X W X D)	(1370 x 720 x 800 mm)	Trigger Window	60		312e (FLX W X D)	(420 x 230 x 430 mm)
Weight	165.4 lb (75 kg)	Tidal Volume			Weight (empty)	(420 x 230 x 430 fffff) 33 lb (15 kg)
Top Shelf		Minute Volume			Absorbent Capacity	
10p 3Hell	(720 x 480 mm)	Inspiratory Time			Absorberit Capacity	2.5 lb (1.5 kg)
Work Surface	25.6" x 11.8"	Support Pressure			Sigma Delta Vaporizer	
WOIR Suitace	(650 x 300 mm)	Support riessure	3 10 20 CITII 1	20 (FEEF Teleficed)	Physical	
Drawers (up to three)	,	Alarms—optional, user set			Weight	Approx. 5 kg
Drawers (up to triree)	(190 x 520 x 390 mm)	Tidal Volume–Minimum		0 to 1600 ml	Canacity	Volume at MAX mark: 250 ml nominal
	(190 x 320 x 390 Hilli)			20 to 1600 ml	Сарасіту	Volume at MIN mark: 35 ±10 ml
Features		Minute Volume–Minimum			Note: After draining, approxima	
Number of Vaporizers	2	Minute Volume–Maximum			retained by the wick	itely 00 ±10 mil of liquid is
Gases		Apnea Adjustable Re			retained by the wick	
	Four	Low and High O <sub>2</sub> Concentratio			Operating Flow Range	0.2 to 15 liters / min
Oxygen Fail Safe		High Airway Pressure			Operating Temperature Range	
Anti-Hypoxic Device		riigir Aii way riessure	10 to t	50 CITII 1 <sub>2</sub> O adjustable	Operating remperature nange	(58 to 95°F)
Standards		Default Settings	Adult	Pediatric		(36 to 93 1)
Staridards	All relevant	Volume, Tidal Volume (Vt):	600 ml	150 ml	Size (H x W x D)	
AV-S Multi-Mode Ventilator		Volume, Rate (BPM):	10	150 1111	Cagemount Dimensions	8.6" x 5.2" x 6.2"
Physical		Volume, I:E Ratio:	1:2	1:2	cagemount Dimensions	(219 x 133 x 158 mm)
Size, Control Unit Only (H x W x D)	7 3" v 11 <i>4</i> " v 11 8"	Volume, Pmax:	38 cmH <sub>3</sub> O	38 cmH <sub>2</sub> O	Selectatec® Compatible w/	(217 × 133 × 130 11111)
Size, Control Offic Offiy (FTX W X D)	(185 x 290 x 300 mm)	voidifie, i friax.	30 CITII 1 <sub>2</sub> O	30 CITII 1 <sub>2</sub> O	Interlock Dimensions	9.5" x 4.7" x 7.5"
Size with Adult Bellows (H x W x D)		Pressure, Tidal Volume (Vt):	600 ml	150 ml	Interiock Dimensions	(242 x 120 x 190 mm)
Size With Addit bellows (FFX W X D)	(385 x 290 x 300 mm)	Pressure, Rate (BPM):	10	150 1111	Dräger Plug-In® Compatible D	
Weight (Control Unit)	16.8 lb (7.6 kg)	Pressure, I:E Ratio:	1:2	1:2	Drager riag in Compatible L	(242 x 100 x 190 mm)
Weight (w/Adult Bellows)	10.8 lb (0.0 kg)	Pressure, P-Target:	10 cmH <sub>3</sub> O	10 cmH <sub>2</sub> O		(242 × 100 × 130 IIIIII)
	8.4 inches (210 mm)	riessuie, riarget.	10 CITII 120	10 CITII 120	Sigma Alpha Vaporizer	
Bellows (Latex-Free)	20 to 1600 ml	SIMV, Tidal Volume (Vt):	600 ml	200 ml	Physical	
Bellows (pediatric option)		SIMV, Rate (BPM):	6	10	Weight	Approx. 6 kg
Power	90 to 264 VAC 47 to 63 Hz	SIMV, Nate (BI W).	2 sec	1 sec		ume at MAX mark: 330 ml working volume
Drive Gas		SIMV, Trigger:	0.4 L/min	0.4 L/min		lume at MIN mark: 70 ml refill indicator
Dive das	Oxygen or 7th	Silviv, mgger.	0.4 L/111111	0.4 [/111111	Note: After draining, approxima	
Functions		SMMV, Minute Volume (Vt):	3.6 L	2 L	retained by the wick	itely 60 ±10 mil of liquid is
Tidal Volume (Vt)	20 to 1600 ml	SMMV, Rate (BPM):	6	10	retained by the wick	
Rate (BPM)	4 to 100 bpm	SMMV, Inspiratory Time:	2 sec	1 sec	Operating Flow Range	0.5 to 10 liters / min
Sigh Mode		SMMV, Trigger:	0.4 L/min	0.4 L/min	Operating Temperature Range	
Sigit Woode	every 50 breaths	Sivilviv, mgger.	0.1 2/111111	0.1 2/111111	operating remperature names	(64 to 86°F)
I:E Ratio	1:0.3 to 1:8	PSV, Support Pressure:	10 cmH <sub>3</sub> O	10 cmH <sub>3</sub> O		(0110001)
Pressure Limit		PSV, Inspiratory Time:	2 sec	1 sec	Size (H x W x D)	
Fresh Gas Compensation		. 51,spiracory rime.	2 300	1 300	Cagemount Dimensions	8.6" x 5.2" x 6.2"
restrads compensation	Adjustment	Alarms—Automatic			cagemoune Dimensions	(219 x 133 x 158 mm)
Ventilation Modes Off, Stand	dby, Volume, Pressure Controlled,	Alarm Mute		30 seconds	Selectatec® Compatible w/	(21) X 133 X 130 11111)
vertellation violacs	Spontaneous, SIMV, SMMV, PSV	Low Drive Gas Pressure	Less	than 235 kPa (34 psi)	Interlock Dimensions	9.5" x 4.7" x 7.5"
(for us	se in anesthesia procedures only)	High Cont. Airway Press				(242 x 120 x 190 mm)
Pressure Control	10 to 70 cmH <sub>2</sub> O	Low Pressure		O (PEEP referenced)	Dräger Plug-In® Compatible D	· · · · · · · · · · · · · · · · · · ·
	tive Volume and Pressure Alarms,	Low Tidal Volume				(242 x 100 x 190 mm)
	t Function—Automatic switch to	Incorrect Rate or Ratio, Mains F				(2.27.1007.13011111)
	Mode if apnea alarm is triggered	co.recertace of riddo, mains i		backup		
	4 to 30 cmH <sub>2</sub> O	Low Battery		5 minutes use		
Oxygen Monitor		Vent Inop				
5.7, gc.11410111101	T der cell type	· cemop	IIICCI	a. o. bactery rundic		



1-800-477-2006 • www.dremed.com

9

