



YOUR OWN CLIMATE





ARALAB is a company specialized in designing, developing, manufacturing and servicing of **environmental simulation chambers** and **controlled environment rooms**.

Since 1985 we have been perfecting ways to create and control **temperature, humidity, light, air flow** and many other environmental conditions.



What drives us is simple: to develop solutions that go beyond the expectations of users so customers can enjoy the best equipment for their research and testing applications.



Aralab is ISO:9001 certified for its Quality Management System
Aralab is ISO:14001 certified for its Environmental Management System



TESTING





ENVIRONMENTAL AND TEMPERATURE TESTING CHAMBERS

ARALAB TESTING

Environmental and Temperature
Testing chambers

- Automotive components
- Electronics
- Aeronautics
- Coatings
- Building materials
- Calibration and Metrology
- Quality Control
- Pharmaceutical and Medical devices
- R&D
- HVAC-R
- Solar PV




KEY FEATURES

- Environmental conditions controlled with consistent precision through the years
- The most reliable technology in climate simulation
- Internal aerodynamic optimization to ensure uniformity of climatic conditions
- Highly resistant stainless-steel interior for maximum durability and easy cleaning
- Flexible interior with height adjustable and removable stainless-steel shelves
- Time saving features with easily configurable testing programs that can run, start and stop automatically
- Compliant with international standards and requirements EN, IEC, DIN, ISO, NP and UNE


SOME TESTING REFERENCES



TESTA CT & TT TESTING CHAMBERS



300 TO 2.000 LITERS “REACH-IN” ENVIRONMENTAL CHAMBERS FOR CLIMATIC AND/OR TEMPERATURE TESTING



TESTA CT & TT 300 TO 2000 LITRES TESTING CHAMBERS



'REACH-IN' ENVIRONMENTAL CHAMBERS FOR CLIMATIC AND/OR TEMPERATURE TESTING APPLICATIONS

'REACH-IN' TESTING CHAMBERS	
PERFORMANCE IN CLIMATIC TESTING RANGE (TEMPERATURE + HUMIDITY CONTROL)	
TEMPERATURE RANGE	10°C to 90°C
TEMPERATURE UNIFORMITY	± 0,1°C to ± 1,0°C ^(1b)
TEMPERATURE FLUCTUATION	± 0,1°C to ± 0,3°C ^(1b)
HUMIDITY RANGE	10% RH to 98% RH
HUMIDITY FLUCTUATION	± 0,5% RH to ± 2% RH
PERFORMANCE IN TEMPERATURE TESTING (TEMPERATURE CONTROL WITHOUT HUMIDITY)	
TEMPERATURE RANGES	-75°C and -50°C, up to 180 °C
TEMPERATURE UNIFORMITY	± 0,5°C to ± 1,5°C
TEMPERATURE FLUCTUATION	± 0,1°C to ± 0,5°C
TEMPERATURE RATE OF CHANGE HEATING	From 2,5°C to 10°C / minute
TEMPERATURE RATE OF CHANGE COOLING	From 2°C to 3,5°C / minute 5°C / minute in "5K models" 10°C/minute available in "10K models"
MODELS / DIMENSIONS AVAILABLE	
INTERNAL TEST VOLUMES	300, 500, 1.000 , 1.500 and 2.000 liters
OTHER TECHNICAL DATA	
NOISE LEVEL	55 to 68 dBA
ELECTRICAL CONNECTION	3/N/PE AC 400V ± 10% 50Hz

TESTA CT & TT 300 TO 2000 LITRES TESTING CHAMBERS



'REACH-IN' ENVIRONMENTAL CHAMBERS FOR CLIMATIC AND/OR TEMPERATURE TESTING APPLICATIONS



CONSTRUCTION

Interior: AISI 304 hermetical welded, vapor tight, stainless steel


Exterior: Zinc mild steel with epoxy coating finish (RAL 7035)

Insulation: Rock Wool

Interior light: Halogen lamp (available with optional window)

Door: Double silicone joints and anti-condensation heating frames (optional window)

Security, Communications and Mains: on the side of the chamber with high/low safety thermostat; mains switch; audible and visual alarms; RJ45 communications ports



ENVIRONMENTAL CONTROL

Heating: stainless steel electric heaters located in the air treatment tunnel

Cooling: scroll compressor group. -75°C chambers and models with cooling rate upgrades will have water cooled condenser as standard and air / water as option.

Humidification and Drying: Thermostatic bath with dew point control. For drying, an additional dry coil

Temperature control: Two sensors: one PT100-A located on the air return channel, and one PT100-A, movable sensor for flexible placing inside chamber

Humidity control: Psychrometric, Capacitive, or both.

TESTA 'WALK-IN' 8.000 TO 12.000L TESTING CHAMBERS



'WALK-IN' ENVIRONMENTAL CHAMBERS FOR CLIMATIC AND/OR TEMPERATURE TESTING APPLICATIONS



TESTA 8.000 TO 12.000 LITRES TESTING CHAMBERS



'WALK-IN' ENVIRONMENTAL CHAMBERS FOR CLIMATIC AND/OR TEMPERATURE TESTING APPLICATIONS

ARALAB TESTA 'WALK-IN' TESTING CHAMBERS	
PERFORMANCE IN CLIMATIC TESTING RANGE FOR TESTA CHAMBERS (TEMPERATURE + HUMIDITY CONTROL)	
TEMPERATURE RANGE	10°C to 95°C
TEMPERATURE UNIFORMITY	± 0,5°C to ± 1,0°C
TEMPERATURE FLUCTUATION	± 1°C
HUMIDITY RANGE	10% RH to 95% RH
HUMIDITY FLUCTUATION	± 1% RH to ± 3% RH
PERFORMANCE IN TEMPERATURE TESTING FOR TESTA CHAMBERS (TEMPERATURE CONTROL)	
TEMPERATURE RANGES	-60°C and -40°C, up to 150 °C
TEMPERATURE UNIFORMITY	± 0,5°C to ± 1,5°C
TEMPERATURE FLUCTUATION	± 0,1°C to ± 0,5°C
TEMPERATURE RATE OF CHANGE HEATING	From 1,5°C to 5°C / minute
TEMPERATURE RATE OF CHANGE COOLING	From 1°C to 7°C / minute
MODELS / DIMENSIONS AVAILABLE	
INTERNAL TEST VOLUMES	8.000, 12.000, or 18.000 liters
OTHER TECHNICAL DATA	
NOISE LEVEL	65 to 75 dBA
ELECTRICAL CONNECTION	3/N/PE AC 400V ± 10% 50Hz

TESTA 8.000 TO 12.000 LITRES TESTING CHAMBERS



'WALK-IN' ENVIRONMENTAL CHAMBERS FOR CLIMATIC AND/OR TEMPERATURE TESTING APPLICATIONS



CONSTRUCTION

Interior: AISI 304 hermetical welded, vapor tight, stainless steel

Exterior: Zinc mild steel with epoxy coating finish (RAL 7035)

Insulation: Rock Wool

Interior light: Halogen lamp (available with optional window)

Door: Double silicone joints and anti-condensation heating frames (optional window)

Security, Communications and Mains: on the side of the chamber with high/low safety thermostat; mains switch; audible and visual alarms; RS232 (or RJ45) communications ports



ENVIRONMENTAL CONTROL

Heating: stainless steel electric heaters located in the air treatment tunnel


Cooling: scroll compressor group. -40°C and -60°C chambers and models with cooling rate upgrades will have water cooled condenser as standard and air / water as option.

Humidification and Drying: Thermostatic bath with dew point control. For drying, an additional dry coil


Temperature control: Two sensors: one PT100-A located in air treatment tunnel and one PT100-A, movable sensor for flexible placing inside chamber

Humidity control: Psychrometric, Capacitive, or both.


TEMPERATURE SHOCK – FITOTERM 150 CTE



2 ZONES TEMPERATURE SHOCK TESTING



TEMPERATURE SHOCK – FITOTERM 150 CTE THERMAL SHOCK CHAMBER



TEMPERATURE TESTING CHAMBER FOR TEMPERATURE SHOCK TESTING

FITOTERM 150 CTE THERMAL SHOCK

SPECS

TEMPERATURE RANGES	-75°C to 200°C
TEMPERATURE IN HOT (UPPER) CHAMBER	Ambient to +200°C
TEMPERATURE IN COLD (LOWER) CHAMBER	-75° C to +60° C
TRANSFER TIME BETWEEN HOT AND COLD CHAMBER	≤ 5 seconds
TEMPERATURE FLUCTUATION (TIME)	≤ ± 0,5°C
TEMPERATURA PRECISION (SPACE)	≤ ± 1,0°C

MODELS/DIMENSIONS AVAILABLE

INTERNAL TEST VOLUME (NET)	125 liters
TEST SPACE DIMENSIONS	410 mm x 470 mm x 650 mm
ADMISSABLE WEIGHT LOAD	50 Kg

OTHER TECHNICAL DATA

NOISE LEVEL	60 to 70 dBA
ELECTRICAL CONNECTION	3/N/PE AC 400V ± 10% 50Hz



TEMPERATURE SHOCK - 150 CTE THERMAL SHOCK CHAMBER



TEMPERATURE TESTING CHAMBER FOR TEMPERATURE SHOCK TESTING



CONSTRUCTION

Interior: AISI 304 hermetical welded, vapor tight, stainless steel

Exterior: Zinc mild steel with epoxy coating finish (RAL 7035)

Insulation: Rock Wool

Door: Double silicone joints and anti-condensation heating frames

Observation Window: located at Hot (upper) chamber door. Multilayer anti condensation glass

Interior light: Halogen lamp

Security, Communications and Mains: on the side of the chamber with high/low safety thermostat; mains switch; audible and visual alarms; RS232 (or RJ45) communications ports



ENVIRONMENTAL CONTROL

Heating: stainless steel electric heaters located in the air treatment tunnel

Cooling: scroll compressor group. Water cooled condenser as standard.

Temperature control: Two sensors: one PT100-A located in air treatment tunnel and one PT100-A, movable sensor for flexible placing inside chamber

Cool-up Feature: accelerate cooling of the Hot (upper) chamber by fresh air renewal

TESTA 600 & 1200 PH-CI CURING CHAMBERS




REACH-IN CHAMBERS WITH TEMPERATURE AND HUMIDITY CONTROL



TECHNICAL DATA FOR FITOCLIMA 600 & 1200 PH-CI CHAMBERS

TEMPERATURE RANGE ^[1]	-5°C to 60
TEMPERATURE PRECISION	± 0,5°C
TEMPERATURE UNIFORMITY	± 1,0°C
HUMIDITY RANGE ^[1]	20 to 98% rH
HUMIDITY PRECISION	± 1% rH
HUMIDITY UNIFORMITY	± 2% rH
SHELVES (STANDARD CONFIGURATION)	FitoClima 600: 4 stainless steel reinforced shelves FitoClima 1.200: 8 stainless steel reinforced shelves
STANDARD WIRE SHELF SIZE	620mm x 520mm
STORAGE ^[2]	0.33m ² and 18 Kg weight load (per shelf)
STORAGE (STANDARD CONFIGURATION)	FitoClima 600 (w/ 4 shelves): 1,33 m ² FitoClima 1.200 (w/ 8 shelves): 2,67 m ²
AIRFLOW	0.2 m/s uniform across the shelves
INTERNAL VOLUMES	542 and 1.200 liters


COMMON 'TESTA REACH-IN' ACCESSORIES



Door with observation window



Air Dryer



Reinforced SS shelves (+100Kg)



Latex Gloves working ports



Safety stack light signal



Integrated Water Tank



Gas Sensors



Freeze-Thaw testing system

- (...)
- FitoLog Software pack
- Cooling / Heating rate upgrades
- Water or Air cooled condensers
- Safety door lock
- Additional entry ports
- ATEX and EUCAR safety systems

TESTA_E

New Line of Environmental and
Temperature Testing chambers

**BENCHTOP AND COMPACT
REACH-IN
TEMPERATURE SHOCK**



[Click image to see Highlights movie](#)

30 LITER BENCHTOP CHAMBER WITH TEMPERATURE CONTROL




TESTA_e MINI TT 30 E65

TEMPERATURE RANGE		-65°C to +150°C
TEMPERATURE FLUCTUATION		± 1°C
TEMPERATURE DEVIATION		± 2.0°C (temperature ≤ 100°C); ± 3.0°C (temperature > 100°C)
HUMIDITY RANGE		-
HUMIDITY DEVIATION		-
TEMPERATURE HEAT-UP RATE		-65°C to +150°C; ≤ 40min
TEMPERATURE PULL-DOWN RATE		+20°C to -60°C; ≤ 50min



TESTA_e TEMPERATURE SHOCK TS120



2 ZONES TEMPERATURE SHOCK CHAMBER





NORMAL INTERIOR VOLUME (L)		About 120
TEST SYSTEM		2 zones basket moving up-and-down driven by a motor and a screw
HIGH TEMP. ZONE	TEMPERATURE RANGE	+60°C - 220°C/ (20~85) °C
	TEMPERATURE HEAT-UP RATE	≥ 14°C/ min
LOW TEMP. ZONE	TEMPERATURE RANGE	-80°C - +70°C
	TEMPERATURE HEAT-UP RATE	≥ 2°C/ min
	TEMPERATURE PULL-DOWN RATE	≥ 5°C/ min





CLIMAPLUS AND FITOLOG®




PROPRIETARY CONTROLLERS
AND SOFTWARE




CLIMAPLUS HMI CONTROLLER



- Capability for creating 50 programs of 50 segments each
- Internal non volatile memory for storing test data
- Automatic restart of tests due to power failure, without losing data and restarting test where it was interrupted
- Real-time monitoring of all functions and control of equipment.
- Manage control settings via MODBUS/TCP
- Possibility of programming a delay of the beginning of test
- Monitoring and recording of all alarms
- Possibility of performing events by external commands
- Several outputs for connecting computers or other devices
- Alarms management
- Graphic representation of the tests and conditions
- Remote access through VNC server
- Possibility of running computer test programs and export them to the controller






FitoLog and FitoLogView: Records and displays in real time all data and details related to the set-points, running variables and equipment behaviour. It also retrieves information about the active components of the chamber, running processes, errors, alarms and allows the configuration of periodic or alarm triggered remote notifications (by email or SMS, depending on existing connections and accessories).

FitoProgram: This application simplifies the creation of programs and their integration on the chamber ClimaPlus controller. Easily create environments that simulate your testing environment and safely send the programs to any Aralab chamber.

Testing secured with alarms, notifications, fast diagnostics and prompt troubleshooting: With FitoLog it is possible to gather data from each of the chambers systems, which makes it a very useful tool to diagnose any necessary maintenance. This tool works as the “black box” of the equipment, giving Aralab technicians the necessary data to remotely carry out a fast and efficient diagnostic. All that is needed is a FitoLog file.





The screenshot shows the FitoLog software interface. On the left is a configuration sidebar for email settings. The main area displays a dashboard with several data panels for different sensors:

- Schneider -75**: Temperature -3,5 °C, Humidity 0,0 RH, Ventilation 100 %, DryBulb Temp -4,3 °C, Mobile Temp -4,2 °C, Cap Temp 0,2 °C, Psy Humidity 0,0 RH, Cap Humidity 93,9 RH.
- Schneider -45**: Temperature 25,0 °C, Humidity 59,7 RH, Ventilation 100 %, DryBulb Temp 25,0 °C, Mobile Temp 25,6 °C, Cap Temp 26,3 °C, Psy Humidity 59,7 RH, Cap Humidity 54,1 RH.
- 25000 HP 1070**: ALARM Fail Comm, Temperature °C, Humidity RH, Radiation x10µm/s, CO2 x10ppm, Soil Moist %vol, Light Bank 1-3 %, Ventilation %, Air Inlet %, Air Outlet %, Mobile °C, Temp 2 °C, Hum 2 RH, Events.
- 01200PH 1590**: ALARM Fail Comm, Temperature °C, Humidity RH, Radiation x10µm/s, CO2 x10ppm, Soil Moist %vol, Light Bank 1-3 %, Ventilation %, Mobile °C, Temp 2 °C, Hum 2 RH, Events.
- 2500 HP**: Temperature °C, Humidity RH, Radiation x10µm/s, CO2 x10ppm, Soil Moist %vol, Light Bank 1-3 %, Metal Flake %, Ventilation %, Mobile °C, Temp 2 °C, Hum 2 RH, Events.

Notify users when an alarm occurs or automatically send log files with a given periodicity to store data in a remote location.

TESTING



OTHER TESTING CHAMBERS



PSYCHROMETRIC TEST ROOMS FOR AIR CONDITIONERS, GAS BOILERS AND HEAT PUMPS

WALK-IN ROOMS FOR
TESTING ACCORDING TO ISO 14511

ISO 13253 | ISO 5151
DIN EN 14825 | DIN EN 16147
JIS B 8615-1
ASHRAE 41.1 | ASHRAE 41.2
ASHRAE 41.3 | ASHRAE 41.6



● ● ● ● INDOOR ROOM		
TEMPERATURE		+5°C to +45°C
RH RANGE		20% to 90% RH
TEMPERATURE UNIFORMITY		± 2,0°C
HUMIDITY UNIFORMITY		± 2 % RH
● ● ● ● OUTDOOR ROOM		
TEMPERATURE		-30°C to +60°C
RH RANGE		20% to 90% RH
TEMPERATURE UNIFORMITY		± 2,0°C
HUMIDITY UNIFORMITY		± 2 % RH

[Click here to see Highlights movie](#)



TESTACAL

ENVIRONMENTAL CHAMBER FOR METROLOGY AND CALIBRATIONS

300 Liter Reach-in chamber.

Highly accurate and uniform
temperature and humidity
control



Temperature uniformity

in Space @ low temp. point	°C	± 0,5
in Space @ 0°C	°C	± 0,50
in Space @ +25°C	°C	± 0,20
in Space @ +50°C	°C	± 0,30
in Space @ +75°C	°C	± 0,60
in Space @ high temp point	°C	± 1,5



BATTERY TESTING CHAMBER - EUCAR

REACH-IN CHAMBER
TESTA 500 BT



EUCAR HAZARD LEVEL PROTECTIONS

HAZARD LEVEL 0-3	HAZARD LEVEL 4	HAZARD LEVEL 5	HAZARD LEVEL 6
Status indicator	Mechanical Door Lock	Fire extinguishing	Permanent gas inertization with nitrogen or argon. O ₂ measurement
Emergency Stop	Reversible pressure release	Fire detection w/ CO Sensor	Gas measurement H ₂ , HF, O ₂ , CO ₂ , HC
Security Door lock	Especial entry port	ATEX Extration / renovation	Explosion Pressure relief






SOLAR AND UV RADIATION SIMULATION CHAMBER

ENVIRONMENTAL
AND RADIATION SIMULATION

TESTA 1.000 WITH RADIATION
SIMULATION MODULE







COLD BEND TESTING OF CABLES AND WIRES

REACH-IN CHAMBERS
TESTA 300 and 500 CB



TESTING



SOME EXAMPLES OF SPECIAL OR 'TURN-KEY' PROJECTS



COMBINED VIBRATION AND TEMPERATURE/ CLIMATIC TESTING

INTEGRATION WITH HORIZONTAL,
VERTICAL OR MULTIAXIAL SHAKERS






COMBINED VIBRATION AND TEMPERATURE/ CLIMATIC TESTING

INTEGRATION WITH HORIZONTAL,
VERTICAL OR MULTIAXIAL SHAKERS





**VIBRATION
AND CLIMATIC
TESTING
IN 'MULTI-POSITION'
PLATFORM, WITH
HORIZONTAL,
VERTICAL OR
MULTI-AXIAL**





SOLAR PANELS TESTING CHAMBERS

REACH-IN AND WALK-IN
CHAMBERS FOR SOLAR
PHOTOVOLTAIC TESTING
STANDARDS





CORROSION GAS TESTING

WALK-IN CHAMBERS FOR TESTING NOXIOUS GASES





TESTING OF REFRIGERATED DISPLAY CABINETS

ENVIRONMENTAL
TESTING ROOMS





ENVIRONMENTAL TESTING 'DRIVE-IN' CHAMBERS


AUTOMOTIVE
ENVIRONMENTAL SIMULATION
DRIVE-IN CHAMBERS





ENVIRONMENTAL TESTING 'DRIVE-IN' CHAMBERS

AUTOMOTIVE ENVIRONMENTAL SIMULATION DRIVE-IN CHAMBERS





COMBINED VIBRATION AND TEMPERATURE/ CLIMATIC TESTING

INTEGRATION WITH HORIZONTAL,
VERTICAL OR MULTIAXIAL SHAKERS





SUN RADIATION AND TEMPERATURE TESTING

200m³ DRIVE-IN CHAMBER





CUSTOM ENVIRONMENTAL TEST CHAMBERS

**TURN-KEY CHAMBERS WITH
BESPOKE PERFORMANCES AND
DIMENSIONS**



THANK YOU FOR YOUR TIME!

WWW.ARALAB.PT
ARALAB@ARALAB.PT

FIND OUT MORE:

