

PRODUCT CATALOGUE

light & medicine



UV-C
flow
germicidal lamps

UV-C
direct radiation
germicidal lamps

LED NGP
X-ray film
viewers

DiCO
digital images
viewing station

FOTOVITA
SAD phototherapy
lamps

ULTRA*Vi*OL

Table of contents:



Welcome
to our
product
catalogue.

About ULTRAVIOL.....	4 - 5
■ UV-C flow germicidal lamps	
ASEPTOR Basic series.....	6 - 9
GERMIPROTECT GP 4x55 series.....	10 - 11
NBVE series.....	12 - 15
■ UV-C direct radiation germicidal lamps	
NBV series.....	16 - 17
NBV IP65 series.....	18 - 19
NBV Multi-directional series.....	20 - 21
■ LED-NGP X-ray film viewers	
LED-NGP series.....	22 - 23
LED-NGP WS series.....	24 - 25
■ DiCO digital images viewing station	
DiCO series.....	26 - 29
■ FOTOVITA SAD phototherapy lamps	
FOTOVITA series.....	30 - 31

About ULTRAVIOL

history, mission & vision

Founders



Wiesław Pietras



Radosław Purgał

Successors



Piotr Pietras




Mateusz Purgał


ULTRAVIOL is an experienced manufacturer and supplier of medical equipment. The company was founded in 1993. ULTRAVIOL's appliances work in every prestigious clinic and hospital both in Poland and worldwide.

ULTRAVIOL company meets the highest requirements set for manufacturers of medical equipment. The enterprise obtained the certificate of compliance for quality management system ISO 13485 issued by TÜV Nord Poland sp. z o.o.

Germicidal UV-C lamps manufactured by ULTRAVIOL have been used to fight the world COVID-19 pandemic by thousands of medical units, beauty salons, hospitals and healthcare facilities.



Sales to over 70 countries



MEDICA 2019 - MESSE DÜSSELDORF



ARAB HEALTH 2023 - DUBAI

ASEPTOR Basic®

a new UV-C flow germicidal lamps for safe air disinfection

Our over 30 years of experience in designing and manufacturing UV-C germicidal lamps and awareness of customers' needs enabled us to create the new quality in safe and effective air disinfection. ASEPTOR is the latest innovation by ULTRAVIOL company.

UV-C flow germicidal lamps provide one of the most efficient methods of safe air disinfection (reduction of viruses, bacteria, moulds and fungi). Bulbs built in ASEPTOR emit UV-C irradiation (wavelength - 253.7 nm) known for its germicidal effectiveness. UV-C radiation that can kill every microorganism does not come out from disinfection chamber. That is why UV-C germicidal flow lamps are commonly used in occupied rooms and areas where air purity determines the quality and safety of services and work.




99.9 %
microorganisms
reduction


236 <20 dB
255 <30 dB
295 <37 dB
fan noise
level


236 - 90 m³
255 - 150 m³
295 - 200 m³
disinfected
cubature


9000 h
useful lifetime
of the bulbs



How does the single purpose UV-C flow germicidal lamp ASEPTOR Basic work?

1. Contaminated air is drawn by a fan – through a filter catching the dust and other contaminations into the disinfection chamber.
2. The UV-C radiation intensity and the time during which the air remains in the disinfection chamber are selected so that the air blown out from the lamp is practically free of microorganisms.
3. The forced flow of the air results in its smooth circulation and causes disinfection of the air in the whole room. This is a safe way to get rid of viruses, bacteria, molds and fungi from the air.

The original technical solutions applied in ASEPTOR are patented

Applications:



Healthcare



Public Facilities



Offices




Industry

Germicidal effectiveness of UV-C radiation

Microorganisms exposed to UV-C radiation are quickly deactivated. This phenomenon is called germicidal effect. As a result of numerous scientific researches it is scientifically proven statement, that the strongest biocidal effect occurs in the range of 250-270 nm radiation wave. The germicidal mechanism is based on the absorption of UV-C radiation energy by nucleic acids and proteins. This process triggers chemical reactions in cell nucleus and results in microorganisms deactivation (all the viruses, including SARS-CoV-2, bacteria, molds, fungi and many others).

UV-C radiation is a shortwave one, therefore it is also high-energy radiation. The energy of the photons absorbed by the nucleic acids causes disruption of DNA molecular bonds and creates pyrimidine dimers. UV-C radiation inactivates DNA and RNA of microorganisms.



ASEPTOR effectiveness research by Prof. Wacław Dąbrowski Institute of Agriculture and Food Biotechnology- State Research Institute

Modern & convenient solutions



Stability and easy smooth movement in every direction



Wall-mounted version ensures aesthetic installation



Fast & easy anti-dust filter exchange



Technical data

Lamp type	ASB 236	ASB 255	ASB 295
Supply voltage	230 V, 50 Hz	230 V, 50 Hz	230 V, 50 Hz
Power consumption	80 W	120 W	200 W
UV-C bulbs type (Philips / Osram)	2 x 36 W (PL-L TUV/HNS-L 2G11)	2 x 55 W (PL-L TUV/HNS-L 2G11)	2 x 95W (PL-L TUV/HNS-L 2G11)
Useful lifetime of the UV-C bulbs	9000 h	9000 h	9000 h
Fan capacity	80 m³/h	130 m³/h	165 m³/h
Air flow capability	35 m³/h	60 m³/h	80 m³/h
Cubage of disinfected room	90 m³	150 m³	200 m³
Effective area of the lamp	35 m²	60 m²	80 m²
Fan noise level	<20 dB	<30 dB	<37 dB
Protection against electric shock	I	I	I
Ingress Protection Code	IP 20	IP 20	IP 20
Class for the medical environment	B - home	B - home	B - home
Group according to PN-EN 55011 Clause 5	1	1	1
Compliance with PN-EN 60601-1	YES	YES	YES
Lamp body dimensions (L x W x H)	890 x 140 x 215 mm	1035 x 155 x 250 mm	1170 x 250 x 155 mm
Overall dimensions - wall-mounted version W	890 x 140 x 215 mm	1035 x 155 x 250 mm	1170 x 250 x 155 mm
Overall dimensions - mobile version M	600 x 600 x 1070 mm	600 x 600 x 1250 mm	600 x 600 x 1385 mm

Finishing of the lamp



Powder coating white colour



INOX version



RAL colours on special demand

Wall-mounted version

Mobile version



ASB 236 W C



ASB 255 W C INOX



ASB 295 W C



ASB 236 M C



ASB 255 M C INOX



ASB 295 M

Optional variants of lamp equipment



C - advanced working time LED counter in ASEPTOR Basic

working time counter (counter accuracy to 1 hour)

automatic brightness level - display brightness depending on the light intensity level in the room

quiet alert after 8800 h - visual signalling of the last 200 h of bulbs work

Ch-F check filter alert. Every 1500 hours of operation, the counter informs about the need to check the cleanliness of the filter and to replace it if it is visibly dirty.

acoustic and visual signalization of the UV-C bulbs exchange time (after 9000 h)

bulb error alarm - audio and visual signaling (bulb/Err)



RC - remote control for ASEPTOR Basic

The remote control (RC) is used for remote switching on/off the ASEPTOR Basic UV-C flow germicidal lamps. Available for the newly ordered lamps. Installation possible only by the manufacturer.

GERMIPROTECT GP 4x55®

UV-C flow germicidal lamps for large volume rooms

The UV-C flow germicidal lamp GERMIPROTECT is intended to be used in large - volume rooms such as: open space offices, hotels, cinemas, gyms, restaurants, kitchens, railway stations, waiting rooms and various crowded areas.

UV-C flow germicidal lamps enable air disinfection during personnel and customers presence. They are entirely safe for people. The UV-C irradiation emitted by UV-C bulbs built in GERMIPROTECT irreversibly destroys viruses, fungi, yeasts, moulds and other microorganisms present in the air. GERMIPROTECT creates a kind of barrier preventing from infections spreading and developing. The device has a low maintenance costs.



99.9 %
microorganisms
reduction

250 m³
disinfected
cubature

220 W
UV-C bulbs
power

9000 h
useful lifetime
of the bulbs



The original technical solutions applied in GERMIPROTECT are patented.

How does the single purpose UV-C flow germicidal lamp GERMIPROTECT work?

1. Contaminated air is drawn by a fan into the disinfection chamber.
2. The UV-C radiation intensity and the time during which the air remains in the disinfection chamber are selected so that the air blown out from the lamp is practically free of microorganisms.
3. The forced flow of the air results in its smooth circulation and causes disinfection of the air in the whole room. This is a safe way to get rid of viruses, bacteria, molds and fungi from the air.

Mounting types



GP 4x55 N
wall-mounted version



GP 4x55 S
ceiling-mounted version



GP 4x55 P
mobile version

Technical data

Lamp type	GERMIPROTECT GP 4x55
Supply voltage	230 V, 50 Hz
Power consumption	240 W
UV-C bulbs type (Philips / Osram)	4 x 55 W (PL-L TUV/HNS-L 2G11)
Useful lifetime of the UV-C bulbs	9000 h
Fan capacity	260 m³/h
Device capacity	100 m³/h
Cubage of disinfected room	250 m³
Effective area of the lamp	100 m²
Protection against electric shock	I
Ingress Protection Code	IP 20
Lamp body dimensions (L x W x H)	940 x 350 x 250 mm
Overall dimensions GP 4x55 N - wall-mounted (L x W x H)	940 x 292 x 350 mm
Overall dimensions GP 4x55 S - ceiling-mounted (L x W x H)	940 x 350 x 286 mm
Overall dimensions GP 4x55 P - on mobile stand (L x W x H)	940 x 350 x 900 mm

NBVE®

UV-C flow germicidal lamps in two variants

NBVE series of the UV-C germicidal lamps is the one of the most known product line by ULTRAVIOL. These devices are available in two types: NBVE single purpose UV-C flow germicidal lamps and NBVE dual purpose UV-C flow germicidal lamps.

NBVE single purpose UV-C flow germicidal lamps:

Single purpose UV-C flow germicidal lamps provide one of the most efficient methods of safe air disinfection (reduction of viruses, bacteria, molds and fungi). UV-C bulbs built in the device emit UV-C irradiation (wavelength - 253.7 nm). Its germicidal effectiveness means that UV-C can kill every microorganism.

NBVE dual purpose UV-C flow germicidal lamps:

Dual purpose UV-C flow germicidal lamps with external UV-C bulbs of direct action guarantee a full range of disinfection. They enable intensive disinfection of the air in the presence of people (UV-C flow chamber – function I) and direct disinfection of the whole room when the staff and patients stay outside (UV-C direct radiation tube – function II). Functions can be turned on/off separately or both at once.



99.9 %
microorganisms
reduction



90 m³
disinfected
cubature



RC
remote control
optional




WT
weekly timer
optional

How does the single purpose UV-C flow germicidal lamp NBVE work?

1. Contaminated air is drawn by a fan – through a filter catching the dust and other contaminations into the disinfection chamber.
2. The UV-C radiation intensity and the time during which the air remains in the disinfection chamber are selected so that the air blown out from the lamp is practically free of microorganisms.
3. The forced flow of the air results in its smooth circulation and causes disinfection of the air in the whole room. This is a safe way to get rid of viruses, bacteria, molds and fungi from the air.



Top view
NBVE dual
purpose
germicidal lamp
with 2 UV-C
bulbs




How does the dual purpose UV-C flow germicidal lamp NBVE work?

1. Contaminated air is drawn by a fan – through a filter catching the dust and other contaminations into the disinfection chamber.
2. The UV-C radiation intensity and the time during which the air remains in the disinfection chamber are selected so that the air blown out from the lamp is practically free of microorganisms.
3. The forced flow of the air results in its smooth circulation and causes disinfection of the air in the whole room. This is a safe way to get rid of viruses, bacteria, molds and fungi from the air.
4. Additional external bulbs enable direct irradiation of the surfaces. They can be used only during the absence of people.



Mounting types



Finishing of the lamp



Optional variants of lamp equipment



LW - digital counter with 4-field LED display and acoustic signalling



RC for NBVE - remote control

working time counter (counter accuracy to 1 hour)

LED display shows the current status of the hour counter. During the first hour of work minutes and seconds are displayed separated by a flashing point

Ch-F check filter alert (description p.9).

alert after 8950 h, intermittent acoustic signal informing about the approaching end of bulb(s) effective lifetime (50 hours of effective work left)

acoustic and visual signalization of the UV-C bulbs exchange time (after 9000 h)

The remote control RC NBVE is used for remote switching on/off the single purpose UV-C flow germicidal lamps. Available for newly ordered lamps. Installation possible only by the manufacturer.

Technical data

Lamp type	NBVE 60	NBVE 110	NBVE 60/30	NBVE 110/55	NBVE 60/60	NBVE 110/110
Supply voltage	230 V, 50 Hz					
Power consumption	85 W	115 W	115 W	145 W	145 W	185 W
UV-C bulbs type	2 x 30 W	2 x 55 W	2 x 30 W internal 1 x 30 W external	2 x 55 W internal 1 x 55 W external	2 x 30 W internal 2 x 30 W external	2 x 55 W internal 2 x 55 W external
Useful lifetime of the UV-C bulbs	9000 h					
Radiation intensity of the external UV-C tube at the distance of 1 m	-	-	100 µW/cm ²	150 µW/cm ²	100 µW/cm ²	150 µW/cm ²
Fan capacity	132 m ³ /h	199 m ³ /h	132 m ³ /h	199 m ³ /h	132 m ³ /h	199 m ³ /h
Cubage of disinfected room	25-50 m ³	45-90 m ³	25-50 m ³	45-90 m ³	25-50 m ³	45-90 m ³
Effective area of the lamp	10-20 m ²	18-36 m ²	10-20 m ²	18-36 m ²	10-20 m ²	18-36 m ²
Fan noise level	~32 dB	~36 dB	~32 dB	~36 dB	~32 dB	~36 dB
Protection against electric shock	I					
Ingress Protection Code	IP 20					
Lamp body dimensions (L x W x H)	1125 x 130 x 215 mm		1125 x 130 x 285 mm		1125 x 130 x 355 mm	
Overall dimensions wall-mounted version N	1190 x 145 x 215 mm		1190 x 145 x 285 mm		1190 x 145 x 355 mm	
Overall dimensions ceiling-mounted version S	1190 x 130 x 330 mm		1190 x 130 x 400 mm		1190 x 130 x 400 mm	
Overall dimensions mobile version P	600 x 600 x 1300 mm		600 x 600 x 1300 mm		600 x 600 x 1300 mm	



WT - weekly timer for NBV direct radiation and NBVE dual purpose UV-C flow germicidal lamps



RC for NBVE 2 - purpose - remote control

weekly timer (counter accuracy to 1 second)

programmer - enables programming 17 different configurations for ON and OFF time of the germicidal lamp

display - shows the day of the week and the hour with accuracy to 1 second

The remote control RC NBVE 2-purpose is used for remote switching on/off the dual purpose UV-C flow germicidal lamps. Available for the newly ordered lamps. Installation possible only by the manufacturer.

NBV®

Direct radiation UV-C germicidal lamps

Direct radiation germicidal lamps NBV provides one of the most efficient methods of supporting disinfection process (reducing the population of microorganisms). These devices equipped with UV-C bulbs emit radiation of wavelength 253.7 nm. This spectrum reveals the strongest biocidal characteristics and irreversibly deactivates bacteria, viruses, moulds, fungi and all other microorganisms. Due to their high efficiency germicidal lamps are used wherever high level of microbiological purity of surface and air is required. The devices can be used only during the absence of people in the room.



99.9 %
microorganisms
reduction

RC
remote control
optional

LW
LED digital
counter
optional

BT
bluetooth
switch app
optional p.21

How does the direct radiation UV-C germicidal lamp work?



1. Contaminated air and surface should be irradiated by UV-C direct radiation lamps - only during the absence of people.
2. Effective time required to sterilise the room depends on the power of the lamp, time of irradiation, microorganisms existing in the space and dimensions of the room.
3. The usage of the lamp results in reduction of viruses, bacteria, molds and fungi from the air and irradiated surface.

Mounting types



NBV 30 N
wall-mounted version



NBV 2x30 S
ceiling-mounted version



NBV 2x30 P
mobile version

Technical data

Lamp type	NBV 15	NBV 30	NBV 55	NBV 2x30	NBV 2x55
Supply voltage	230 V, 50 Hz				
Power consumption	18 W	33 W	60 W	66 W	115 W
UV-C bulbs type	15 W	30 W	55 W	2 x 30 W	2 x 55 W
Useful lifetime of the UV-C bulbs	9000 h				
Radiation intensity of the external UV-C tube at the distance of 1 m	0.9 W/m²	2.3 W/m²	2.9 W/m²	3.6 W/m²	3.6 W/m²
Effective area of the lamp	6-8 m²	12-15 m²	15-18 m²	18-22 m²	22-27 m²
Protection against electric shock	I				
Ingress Protection Code	IP 20				
Lamp body dimensions (L x W x H)	500 x 85 x 135 mm	960 x 85 x 135 mm	960 x 85 x 145 mm		

NBV IP65®

Direct radiation UV-C germicidal lamps for industrial use

Direct radiation germicidal lamps NBV IP65 are designed to prevent primary and secondary infections in the food, pharmaceutical and cosmetic industries, warehouses and everywhere the production process takes place. Direct radiation germicidal lamps applied in the rooms where airborne pathogenic microorganisms (pathogens) live, significantly reduce the probability of spreading the infection by the air. Raising the level of microbiological purity of the air and surface helps to destroy and reduce impact of existing outbreaks of pathogens.

In germicidal industrial lamps we apply antibacterial fluorescent UV-C bulbs cover - laminated with the special protective foil. It prevents glass from splashing in case of UV-C bulb breaks or damages. Shrink-wrapped fluorescent bulbs meet the EU and HACCP requirements. Adapted foil does not affect the efficiency of antibacterial fluorescent lamps and simultaneously protects UV-C bulbs against bruises. The ultraviolet radiation does not penetrate through the usual glass, plexiglass and similar materials, therefore anti-splashing foil is the only material that can be used to secure the UV-C bulb.



99.9 %
microorganisms
reduction

40 m²
disinfected
surface

RC
remote control
optional

BT
bluetooth
switch app
optional p.21

How does the direct radiation UV-C germicidal lamp for industrial use work?



1. Contaminated air and surface should be irradiated by UV-C direct radiation lamps - only during the absence of people.
2. Effective time required to sterilise the room depends on the power of the lamp, time of irradiation, microorganisms existing in the space and dimensions of the room.
3. The usage of the lamp results in reduction of viruses, bacteria, molds and fungi from the air and irradiated surface.

Mounting types



NBV 2x30 IP65 N
wall-mounted version / ceiling mounted version both hangers in standard

NBV 2x30 IP65 S

NBV 2x30 IP65 P
mobile version

Technical data

Lamp type	NBV 2x15 IP65	NBV 2x30 IP65	NBV 2x36 IP65	NBV 2x55 IP65	NBV 2x75 IP65
Supply voltage	230 V, 50 Hz				
Power consumption	40 W	65 W	75 W	115 W	160 W
UV-C bulbs type	2 x 15 W	2 x 30 W	2 x 36 W	2 x 55 W	2 x 75 W
Useful lifetime of the UV-C bulbs	9000 h				
Radiation intensity of the external UV-C tube at the distance of 1 m	1.0 W/m ²	2.1 W/m ²	2.8 W/m ²	3.6 W/m ²	6.8 W/m ²
Effective area of the lamp	spot	10-20 m ²	20-25 m ²	25-30 m ²	30-40 m ²
Protection against electric shock	I				
Ingress Protection Code	IP 65				
Lamp body dimensions (L x W x H))	520 x 220 x 160 mm	980 x 220 x 160 mm	1280 x 220 x 160 mm	980 x 220 x 160 mm	1280 x 220 x 170 mm

NBV Multi - directional®

Direct radiation UV-C germicidal lamps

Direct radiation germicidal lamps NBV provides one of the most efficient methods of supporting disinfection process (reducing the population of microorganisms). UV-C bulbs emit radiation of wavelength 253.7 nm. This radiation reveals the strongest biocidal characteristics and irreversibly deactivates bacteria, viruses, moulds, fungi and all other microorganisms. Due to their high efficiency germicidal lamps are used wherever high level of microbiological purity of surface and air is required.

The devices can be used only during the absence of people in the room. NBV Multi - directional lamp enables efficient disinfection of the air and surface in the whole room. The most important device feature is guarantee of 360° disinfection of the air and surface. NBV Multi - directional UV-C lamp is equipped with Bluetooth switch app and motion detector that provide complete safety of the device and protect against accidental undesirable switch on.





360 °
multi - directional
irradiation



BT
bluetooth
switch app
standard





MD
motion detector
standard



COUNTER
UV-C bulb work
time

How does the direct radiation UV-C germicidal lamp work?

1. Contaminated air and surface should be irradiated by UV-C direct radiation lamps - only during the absence of people.
2. Effective time required to sterilise the room depends on the power of the lamp, time of irradiation, microorganisms existing in the space and dimensions of the room.
3. The usage of the lamp results in reduction of viruses, bacteria, molds and fungi from the air and irradiated surface.




Technical data

Lamp type	NBV 8x36	NBV 8x75
Supply voltage	230 V, 50 Hz	
Power consumption	300 W	630 W
UV-C bulbs type	8 x 36 W	8 x 75 W
Useful lifetime of the UV-C bulbs	9000 h	
Radiation intensity of the external UV-C tube at the distance of 1 m	4.5 W/m²	9.2 W/m²
Effective area of the lamp	100 m²	160 m²
Protection against electric shock	I	
Ingress Protection Code	IP 20	
Lamp body dimensions (L x W x H)	500 x 500 x 1620 mm	

Complete safety with BT - Bluetooth Switch & MD - Motion Detector

The device is equipped with a programmable switch module controlled by the NBV App via mobile device with Android (Bluetooth).



setting the delay until the device is switched on (**time to leave the room - minimum 1 minute**)

setting the operating time of the device

acoustic signal informs about the approaching moment of switching on

working-time counter displays UV-C bulbs useful lifetime with notification of the need of replacing them (**visualization on the NBV App and an acoustic signal**)

motion detector turns off UV-C radiation regardless of the operation of the programmable BT switch module.

LED-NGP®

X-ray film viewers

LED-NGP X-ray film viewers are medical devices designed for the analysis of medical images on X-ray films. This is one of the basic methods of diagnosing. X-ray film viewers manufactured in LED technology have extremely efficient parameters such as: value of luminance intensity, uniformity of screen luminance, low operating cost - long lifetime, low energy consumption. The devices support the process of analyzing X-ray films by doctors. The products are made in the 1st class of protection against electric shock. They can be used in operating theaters, doctor's offices, X-ray laboratories, etc. LED-NGP X-ray film viewers are aesthetic devices with very slim casing. They can fit every modern space.




6000 cd/m²
high luminance

≥ 95 %
light
uniformity

35 mm
ultra slim

45 %
of the
standard power
consumption

Types of LED-NGP X-ray film viewers



Specific features of LED-NGP

Ultra slim
casing- 35 mm only

Step-less luminance
adjustment with
ON/OFF switch

*CLICK
Automatic-ON/OFF
when changing
the film- no glare

Technical data

X-ray film viewer type	LED-NGP-11	LED-NGP-21	LED-NGP-31	LED-NGP-41
Supply voltage	90 - 260 V, 50 - 60 Hz			
Power consumption	65 W	125 W	180 W	250 W
Luminance	600 - 6000 cd/m ²			
Uniformity	≥ 95 %			
Screen dimensions	36 x 43 cm	72 x 43 cm	108 x 43 cm	144 x 43 cm
Step-less luminance adjustment (each frame separately)	10 - 100 %			
Colour temperature of the lamps light	> 6500 K			
Protection against electric shock	I			
IP Code	IP 20			
Total mass	4.5 kg	8.0 kg	12.0 kg	16.0 kg
Overall dimensions (L x W x H) mm	430 x 35 x 520	795 x 35 x 520	1160 x 35 x 520	1520 x 35 x 520

LED-NGP WS®

Recessed X-ray film viewers

LED-NGP WS X-ray film viewers are medical devices designed for the analysis of medical images on X-ray films. This is one of the basic methods of diagnosing. X-ray film viewers manufactured in LED technology have extremely efficient parameters such as: value of luminance intensity, uniformity of screen luminance, low operating cost - long lifetime, low energy consumption. The devices support the process of analyzing X-ray films by doctors. The products are made in the 1st class of protection against electric shock. They are integral part of operating theaters, doctor's offices, X-ray laboratories, etc. WS marking means built - in version. The devices can be fit into the wall panels in reference to previously sent plans and drawings.



6000 cd/m²
high luminance

≥ 95 %
light
uniformity

50 mm
ultra slim
recess

45 %
of the standard
power
consumption

Types of LED-NGP WS X-ray film viewers

LED-NGP-11 WS

LED-NGP-21 WS

LED-NGP-31 WS

LED-NGP-41 WS

Installations examples



Technical data

X-ray film viewer type	LED-NGP-11 WS	LED-NGP-21 WS	LED-NGP-31 WS	LED-NGP-41 WS
Supply voltage	90 - 260 V, 50 - 60 Hz			
Power consumption	65 W	125 W	180 W	250 W
Luminance	600 - 6000 cd/m ²			
Uniformity	≥ 95 %			
Screen dimensions	36 x 43 cm	72 x 43 cm	108 x 43 cm	144 x 43 cm
Step-less luminance adjustment (each frame separately)	10 - 100 %			
Colour temperature of the lamps light	> 6500 K			
Protection against electric shock	I			
IP Code	IP 20			
Total mass	6.0 kg	10.0 kg	14.0 kg	18.0 kg
Recess size (mm)	435 x 50 x 545	795 x 50 x 545	1155 x 50 x 545	1515 x 50 x 545
Overall dimensions (L x W x H) mm	470 x 74 x 580	830 x 74 x 580	1190 x 74 x 580	1550 x 74 x 580

DiCO®

Digital images viewing station (PACS, RIS HIS) for operating room (operating theater)

DiCO viewing station is a complete medical computer for operating theater, also known as a digital X-ray films viewer. It is compatible with the PACS, RIS and HIS. DiCO viewing station works via LAN and can be part of the integrated operation room. The hermetic and washable housing distinguishes DiCO as a perfect viewing device to be used in the hospital.

The highest quality of the DiCO station is assured by the fact, that ULTRAVIOL compiles with the ISO 13485 norm for medical devices. Entire process of designing, manufacturing and assembling is strictly controlled and supervised. ULTRAVIOL offers various configurations, types of finish and provides complex technical consulting services. The combination of our know-how and experience resulted in sales to hundreds medical facilities in Poland and worldwide.







i3/i5/i7
intel
generation 12


EASILY
folding
keyboard


**SCHOTT
CONTURAN®**
anti-reflective
glass coating


DICOM
standard
of monitors

Examples of DiCO® station configuration

DICO 1M	1 monitor 21" 24" 27" 32"		PC computer system: - Intel mainboard; i3/i5/i7 processor - RAM 8GB/16GB/32GB - hard drive 240 GB SSD /500 GB HDD or larger on special request - LAN 10/100/1000 MBit socket protected against flooding - graphic card – professional, for high accuracy of the image reproduction - DVD+/-RW - WINDOWS 10® PRO 64 bit
DICO 2M	2 monitors 21" 24" 27"		
DICO 1M/1B	1 monitor 21" 24" 27" 32" 43" 49" 55" single LED X-ray film viewer		
DICO 1M	1 large monitor 43" - 70"		Additionally - Input socket of video signals: 1 x HDMI i 1 x Display Port 2 x USB 2.0 - output of the HDMI or DVI-D video signals
DICO 2M	1 monitor 43" or 49" 1 monitor 21" or 24"		
DICO 2M/1B	1 monitor 43" or 49" 1 monitor 21" lub 24" single LED X-ray film viewer		Additionally - input sockets of video signals: HDMI, DVI-D, VGA 2 x USB 2.0 - output of the HDMI or DVI-D video signals
R - additional monitor on the straight-line arm	1 monitor 19" - 24"		Shows the image „cloned” from the main monitor of the DiCO® staton (optionally - touch screen)
DICO 2M (60"/17")	1 monitor 21" 24" 27"		Visualization system allows simultaneous real-time display of the video signals from DVI-D, VGA, Composite, S-Video, Component inputs on the monitor with 8MP resolution. Advanced control system allows to select any configuration of the displayed images (e.g. 1, 2 or 4 images) and immediate switching them with touchscreen control monitor.



Łódź, Poland



Kraków, Poland

Details of DiCO® station



High-quality medical (reference) monitor EIZO with screen of diagonal from 21" to 70" and resolution min. 2 MP. Compatible with DICOM



Integrated CD/DVD drive for wall mounted and mobile version.



Top-class computer processor guarantees the highest performance and reliability



2x USB sockets protected against flooding



Integrated medical keyboard with aluminium housing, easy to disinfect, silicone with antibacterial coating and touchpad (can be supplied as a separate article)



Silicone medical mouse, Optional version with keyboard and a medical mouse tray.



CD/DVD drive for the recessed version



Reliable and safe mechanism of keyboard folding



High-grade X-ray film viewer made in LED technology, single or double frame, excellent parameters:
luminance 6000 cd / m²
superior light uniformity > 90%
adjustable light intensity 10-100%.
Optional LED X-ray film viewer




DiCO® station meets the requirements of Medical Device Regulation EU 2017/745 (MDR) and is registered in the URPLWiPB Polish competent authority for medical devices in Warsaw and in the European database EUDAMED.

DiCO® station complies with the requirements of the standards: EN 60601-1 (medical requirements safety) and EN 60601-1-2 (electromagnetic compatibility).

Antireflex

Protection of the monitor is made of special glass SCHOTT CONTURAN® with anti-reflective coating 8-times reducing undesirable glare.



Anti-reflective coating



Ordinary glass

FOTOVITA®

SAD phototherapy lamps

FOTOVITA is a medical device for treatment of Seasonal Affective Disorder symptoms. It was designed and tested in cooperation with the scientists from Ludwik Rydygier Collegium Medicum in Bydgoszcz Nicolaus Copernicus University in Toruń Poland.

FOTOVITA lamps have been helping and bringing relief to their users for over 15 years. The illuminance of white light emitted by FOTOVITA SAD lamps can be characterized as similar to the sunlight. It has a particular effect on the human brain, specifically on pineal gland, which is responsible for production of melatonin. People suffer from its overproduction during autumn-winter season, due to the lack of sunlight. The FOTOVITA light reaches the brain through the eyes. Its physical properties help stimulate the pineal gland what results in restoring the harmony of sleep-wake cycle (Circadian rhythm). The process is entirely safe and has been the subject of numerous scientific researches and dissertations all over the world.



CALM
sleep

BETTER
mood

FAST
treatment
effect

5
year
warranty

Technical data



01 FOTOVITA FV-10 S
small

Supply voltage	230 V, 50 Hz
Power consumption	80 W
Luminous flux density	3 200 lx at 50 cm 1 800 lx at 75 cm
Fluorescent lamp useful lifetime	10 000 h
Protection against electric shock	I
Dimensions	285 x 195 x 525 mm
Mass	2.8 kg



02 FOTOVITA FV-10 M
medium

Supply voltage	230 V, 50 Hz
Power consumption	120 W
Luminous flux density	5 000 lx at 50 cm 2 500 lx at 75 cm
Fluorescent lamp useful lifetime	10 000 h
Protection against electric shock	I
Dimensions	285 x 195 x 635 mm
Mass	3.1 kg



03 FOTOVITA FV-10 L
large

Supply voltage	230 V, 50 Hz
Power consumption	230 W
Luminous flux density	10 000 lx at 50 cm 5 000 lx at 75 cm
Fluorescent lamp useful lifetime	10 000 h
Protection against electric shock	I
Dimensions	345 x 280 x 570 mm
Mass	5.0 kg



ULTRA*Viol*

OVER 30 YEARS

ON THE MEDICAL MARKET



34 Stępowizna Str.
95-100 Zgierz, Poland



www.ultraviol.pl/en
office@ultraviol.pl



+48 42 717 77 45
+48 725 779 277