



We help the elderly, and anybody suffering cognitive impairment, fight against dementia, slow down the loss of mental abilities, keep their autonomy for longer, and delay dependency.

We facilitate the therapists and caregivers follow their evolution, and build upon the results of the practice




**PERSONALIZED  
COGNITIVE  
STIMULATION**

**IN-SITE AND  
AT-HOME  
PRACTICE**

**PROVEN  
EFFICACY  
AND RESULTS**

An ideal support for any institution or specialist practicing cognitive stimulation, in-site, at home, or remotely, to prevent and treat the loss of mental skills

# Developed countries are aging rapidly, and at the same time the population suffering from dementia is growing



Prevalence of dementia by age					
EDAD	65-70	70-75	75-80	80-85	>85
Prevalence	15%	26%	47%	63%	65%


Along with the growth of the elderly population, life expectancy also increases, exceeding 80 years.

The probability of developing dementia grows exponentially with age

And the bad news is that, at the present time, dementia cannot be cured. Those who suffer from it end up being totally dependent on others.


(\*) Main countries of Latin-America: México, Argentina, Chile, Colombia, Perú, Venezuela, Uruguay

## Current therapeutic strategy in the absence of a cure for dementia: prevent it and, when detected, slow down its progression.




**Doing cognitive stimulation exercises emerges as the most effective strategy  
to prevent the loss of mental abilities and, when it is detected,  
to slow down its evolution as much as possible.**

# Smartbrain: cognitive stimulation personalized to any need



- 1 Helps prevent and delay the onset of cognitive impairment
- 2 When an impairment is detected, it allows slowing down and mitigating its evolution, maintaining the autonomy and quality of life of the affected persons for a longer period of time.
- 3 Helps the specialists complete the diagnosis, and to follow closely the evolution of their patients over time, helping them to move away from dependence

Smartbrain features thousands of exercises to work out all the mental skills, to prevent and treat any case of cognitive impairment



More than 100.000 different exercises  
Multiple levels of difficulty, from the very basic to the really tough one  
Several languages: English, Spanish, Italian, Swedish, Hebrew, Catalan

**Smartbrain has been developed with a solid scientific foundation  
its therapeutic efficacy has been scientifically proven,  
and is backed by many testimonials from users**

**Medical studies published in peer-reviewed  
journals show that it slows down  
the process of cognitive impairment**

Alzheimer's  
Parkinson's  
Ageing  
AIDS treatment

**The users of Smartbrain tell us that they are  
effectively getting the desired results**

Chemobrain  
Ictus and brain damage  
Down syndrome and intellectual disability  
Multiple sclerosis

#### **Proven efficacy for treating cognitive impairment caused by Alzheimer's**

- "A randomized pilot study to assess the efficacy of an interactive, multimedia tool of cognitive stimulation in Alzheimer's disease" - L. Tarraga, M.Boada, O. Lopez, J. Becker et al. – Journal of Neurology, Neurosurgery and Psychiatry 2006;77

#### **Proven efficacy for treating cognitive impairment caused by Parkinson's**

- "A blind randomized controlled study of the efficacy of cognitive training in Parkinson's Disease" - A. Prats, H. Guerra, M. Crespo et al - Movement Disorders, Vol. 26, No. 7, 2011
- "Cognitive Training and Transcranial Direct Current Stimulation for Mild Cognitive Impairment in Parkinson's Disease: A Randomized Controlled Trial" - Blake J. Lawrence, Natalie Gasson, Andrew R. Johnson, Leon Booth, Andrea M. Loftus. – Hindawi Parkinson's Disease, Volume 2018, Article ID 4318475


#### **Proven efficacy for treating cognitive impairment caused by the therapy for AIDS**

- "A pilot study of the effects of internet-based cognitive stimulation on neuropsychological function in HIV disease" – J. Becker, M. Dew, H. Aizenstein, O. Lopez, L. Tárraga et al. Disability and Rehabilitation 2012


**The only cognitive stimulation program recommended by CEAFA  
Spanish confederation of associations helping people suffering Alzheimer's**



Smartbrain makes it easy for professionals to do cognitive stimulation exercises personalized to each situation and follow the treatment with as many people or groups as needed.



# Smartbrain's functionalities and tools make it easy to offer personalized cognitive stimulation services with several options: In-site, at-home, or telemedicine



Smartbrain presents professionals with a very simple and intuitive user interface, with a dashboard that informs them of the global situation of all their users.

The screenshot shows a dashboard titled 'Users of Juan Ramon Llorente'. The table lists 10 users with the following data:

Name/ alias	Last working session	Global scores	Trend	Notifications
Android TV	2024/02/14	67.53%	-	0
Asocioaciones jugadores veteranos	2024/09/26	22.21%	-	0
Castellsj	2024/02/14	79.41%	-	0
fenomeno	2024/09/30	84.40%	-	0
jrllorente9	2024/08/17	58.07%	-	0
juan	-	No data	-	0
Juan Ramon	2024/04/25	21.59%	↘	0
prueba AJFCB ingles	2023/12/06	43.84%	-	0
Prueba AJFCB italiano	2024/01/02	27.33%	-	0
Prueba Beatriz Moreno	2024/07/22	86.36%	↗	0

Annotations on the left side of the dashboard:

- Space to customize with corporate logo
- Users are identified with an alias and icon/photo

Annotations on the right side of the dashboard:

- Data search and download tools
- Quick editing of users' profiles

Annotations at the bottom of the dashboard:

- Date of last execution of the working plan
- Overall scores achieved to date
- Trend of overall scores
- Notifications of relevant events for each user

Smartbrain offers more than 100.000 exercises to work out all mental skills, and makes it simple to define users' stimulation plans quickly and easily

The image shows the Smartbrain software interface. On the left, a sidebar titled 'Activities' lists various cognitive exercises. On the right, a detailed view of an activity titled 'Pairing additions' is shown, including its description, objectives, and difficulty level. The interface is designed for customization, as indicated by a callout 'Space to customize with corporate logo' pointing to the top left. Another callout 'Activities for stimulating all cognitive skills' points to the sidebar. Three large blue arrows point from the text callouts to the corresponding features in the interface. A bottom callout 'Multiple difficulties: from the simplest to the diabolical' points to the difficulty level dropdown, and another callout 'All activities available in several languages' points to the language dropdown.

Space to customize with corporate logo

Activities for stimulating all cognitive skills

Multiple difficulties: from the simplest to the diabolical

All activities available in several languages

Several options for configuring the working plans

Detailed information on each activity

Smartbrain

Cuida tu memoria

Activities

Activities

10

Search: memory

Category Activity

Calculation	Pairing additions
Calculation	Solve multiplicaciones
Calculation	Solve subtractions
Calculation	Solve additions TV
Calculation	Solve additions
Language	Pairing geometric shapes and their names
Language	Pairing numbers and their names
Language	Pairing objects and their names AJFCB
Language	Pairing objects and their names
Language	Recognize hours named in the same order

Showing 11 to 20 of 115 records  
(filtered from 206 total entries)

1 2 3 4 5 ... 12 >

Activity data

Category Activity

Calculation	Pairing additions
-------------	-------------------

Description

The user must find the pairs of additions and their results that are hidden under the cards that appear on the screen

Objectives

Memory. Calculation. Visual memory. Immediate memory. Attention. Working memory. Numbers recognition.

Difficulty level N° of repetitions Language

1	1	Español
---	---	---------

Preview

# Professionals have multiple options for setting the parameters of the working plans, to customize them to the individual needs of each user

The image shows the Smartbrain platform interface for 'Settings: Juan Ramon'. The top navigation bar includes 'Users', 'Results', 'Notifications 0', 'Activities', 'Templates', 'Help', 'Messages', and a user profile 'Juan Ramon Llorente'. The main content area has tabs for 'Profile', 'Working plan' (which is selected), 'Global scores', 'Last executions', and 'Notifications'. The 'Working plan' tab contains several configuration options:

- Working days:** A list of days of the week with checkboxes for Monday, Tuesday, Wednesday, Saturday, and Sunday. A green arrow points to the 'Wednesday' checkbox.
- Minutes of work out per working session:** A dropdown menu set to 30. A green arrow points to the dropdown arrow.
- Maximum daily work out time:** A dropdown menu set to 60. A green arrow points to the dropdown arrow.
- Don't repeat audio instruction:** An unchecked checkbox.
- Allow skipping repetitions:** A checked checkbox.
- Test template:** A dropdown menu set to 'No test template'. A green arrow points to the dropdown arrow.
- Test frequency:** A dropdown menu set to 'Select frequency'. A green arrow points to the dropdown arrow.

Below these settings is a section for 'Notifications settings' with three checkboxes:

- The user has not done any working session in the [1] last weeks
- The user has not done more than [3] working sessions in the [1] last weeks
- The user has taken more than [6] minutes to complete one of the activities

A blue 'Change' button is located at the bottom of this section. The interface is branded with the 'Smartbrain' logo and the tagline 'Cuida tu memoria'.

Space to customize with corporate logo

Selection of working days

Configuration of notifications to be received


Setting of working time

Adaptation to user's particularities

Creation of templates to facilitate assignment


Option to set periodic exams

**For users, using Smartbrain is very simple and enjoyable:  
they simply answer the exercises that are automatically presented to them  
by pressing the touch screen or with the mouse**



**Users receive audio instructions on what to do,  
and can also see them in writing in case of deficiency or forgetfulness.  
Positive interjections are heard to indicate successes and failures**

**The scores of the executions are automatically recorded, and Smartbrain allows to follow them in full detail, with graphs and the option to download the data**



The image shows the 'Global scores - fenomeno' page of the Smartbrain application. The interface includes a navigation bar with 'Users', 'Results', 'Notifications' (with 0 notifications), 'Activities', 'Templates', 'Help', 'Messages', and a user profile for 'Juan Ramon Llorente'. The main content area has tabs for 'Profile', 'Working plan', 'Global scores' (which is selected), 'Last executions', and 'Notifications'. On the left, there are filters for 'Starting date' and 'Ending date' with search and refresh icons. Below these are sections for 'Overall results for the selected period' and 'Results of the last 10 sessions'. The 'Overall results' section shows bar charts for Attention (75.31%), Orientation (77.50%), Calculation (47.36%), Solve additions (English) (100.00%), Amounts and figures (Español) (0.00%), Order the figures (Español) (15.75%), and Language (96.70%). The 'Results of the last 10 sessions' section shows a bar chart for each session with its date, time, and percentage. A blue double-headed arrow is positioned between the two main sections. To the left of the interface, five blue ovals with arrows pointing to specific features are listed: 'Space to customize with corporate logo' (pointing to the logo area), 'Tools for searching and downloading scores' (pointing to the filters), 'Consultation of scores in each category' (pointing to the 'Overall results' section), 'Consultation of scores in each activity' (pointing to the 'Results of the last 10 sessions' section), and 'Consultation of global scores on all executions' (pointing to the double-headed arrow). To the right, a blue oval with an arrow pointing to the interface is labeled 'Periodic examinations for further control of evolution'.

Space to customize with corporate logo

Tools for searching and downloading scores

Consultation of scores in each category

Consultation of scores in each activity

Periodic examinations for further control of evolution

Consultation of global scores on all executions

**Either if done in-site or at-home, Smartbrain allows the offering of a cognitive stimulation practice that can be personalized to the individual need of each patient, userfriendly, easy to use and effective**



**The full interaction creates a continuous and feedback cycle, which allows the therapists detect any change, risk or need, and do long-term monitoring of the therapy**

**Smartbrain is already being successfully used by individuals and public and private entities in Spain, the United States, Italy and Latin America**

**Cloud operation, which allows the use of Smartbrain 24 hours a day, every day, from any place where Internet access is available, with any computer, tablet, Smart TV, smartphone, digital whiteboard, etc.**



## Smartbrain is also a high social impact project due to the significant socio-economic results that its dissemination allows to achieve

Helps to promote an active aging, with greater personal autonomy, and maintaining the quality of life of the elderly

For every year that cognitive abilities are maintained, the dependency rate is reduced by 20%

Family members and caregivers of affected persons maintain their ability to socialize and their quality of life

Both families and public health services can achieve significant savings in dependency care costs

