

Philips HeartStart XL Defibrillator/Monitor

Product information



The capabilities you need, the flexibility you demand

The Philips HeartStart XL defibrillator/monitor, at just under 14 pounds, is easily transported throughout the hospital to the site where cardiac care is needed. A combination manual device with automated external defibrillator capabilities, HeartStart XL enables the first caregiver on the scene, whether an ALS or BLS clinician, to deliver lifesaving defibrillation therapy.

In AED mode, voice prompts and text messages guide BLS users through the defibrillation process, while HeartStart XL continuously monitors and displays the patient's ECG. Upon the arrival of ALS personnel, HeartStart XL is easily switched from AED to manual mode, allowing operators to access the unit's advanced therapeutic features such as selectable energy (from 2 to 200 Joules), non-invasive pacing (optional), SpO₂ (optional) and synchronized cardioversion.

Easy to use, compact and Withstands hospital use

Features

- Manual and AED operation
- ECG monitoring with pads or electrodes
- Selectable limits and alarms
- Liquid crystal display
- 50 mm strip chart printer
- Data storage/event summary
- Mark events
- SMART Biphasic waveform for defibrillation

Options

- Adult paddles with paddle-to-patient contact indicators (PCI). Pediatric paddles underneath
- SpO₂ pulse oximetry with alarms
- Non-invasive pacing
- 5-Lead ECG



Compact and rugged, the HeartStart XL easily fits on a standard hospital stretcher and withstands the rigors of hospital use and patient transport.


SMART Biphasic Technology

Philips' patented SMART Biphasic waveform is clinically proven superior to high-energy, monophasic waveforms for efficacy and for minimizing post-defibrillation heart dysfunction.^{1,2,3}

- Impedance Compensation. Measures chest impedance and delivers a low-energy shock based on patient's unique physical requirements.
- Quick Charging. Charges to highest energy level, 200 Joules, in less than 3 seconds.
- Synchronized Cardioversion. Philips' SMART Biphasic waveform has undergone clinical testing, demonstrating its effectiveness for cardioversion of atrial fibrillation.⁴

Easy to Use

- 1-2-3 Operation. Defibrillation is intuitive for all users.
- AED Mode. Voice and text prompts guide users through the defibrillation process. 150 Joules non-escalating, pre-set energy level.
- Paddles (optional). Anterior/anterior adult paddles convert to pediatric by removing the outer contacts.
- Multi-function Defibrillator Pads. Adult and pediatric pads for defibrillation, ECG monitoring, pacing, and synchronized cardioversion.



Data Card for 2 hours
of event storage

rugged the rigors of and patient transport.

Data Management and Reporting

All patient data, including continuous ECG and events (i.e., shock and alarms), are stored in the unit's internal memory and on an optional data card. When used with Philips HeartStart Event Review data management software, users can edit, store, and print reports required for quality control and reporting.

Lightweight, Compact, Durable

- Grab and Go. Less than 14 pounds (6.35 kg).
- Compact. Easily fits on a standard hospital stretcher.
- Rugged. Withstands the rigors of hospital use and patient transport.

Service Solutions

- You Choose the Warranty. Five-year unit exchange (North America only), or one-year on-site service with extended warranty options, or 5-year BioMedical Warranty (U.S. only).
- Dependability. Part of the complete family of Philips resuscitation products used by healthcare providers for over 40 years.



XL Defibrillator/Monitor

A combination manual device with automated external defibrillator capabilities, HeartStart XL enables the first caregiver on the scene, whether an ACLS or BLS clinician, to deliver lifesaving defibrillation therapy.





Philips, Leading the Way

Our dedication to excellence in design, manufacturing and customer support makes us a trusted supplier of medical equipment, serving the healthcare community for more than 100 years. HeartStart XL is part of our cardiac resuscitation family of products, which includes ALS monitor/defibrillators and automated external defibrillators used in private and public settings. Each is tailored to the needs and skills of a particular type of user, extending the reach of care from the home to the hospital. To learn more about HeartStart Defibrillators and Philips Medical Systems, visit www.philips.com/heartstart or call 800-453-6860.

Philips Medical Systems is part of Royal Philips Electronics

On the web

www.philips.com/heartstart

Via email

medical@philips.com

By fax

+31 40 27 64 887

By postal service

Philips Medical Systems
3000 Minuteman Road
Andover, MA 01810-1085

Asia

Tel: +852 2821 5888

Europe, Middle East, Africa

Tel: +31 40 27 87246

Latin America

Tel: +1 954 628 1000

North America

Tel: +1 800 934 7372

© Koninklijke Philips Electronics N.V. 2005

All rights are reserved. Reproduction in whole or in part is prohibited without the prior written consent of the copyright holder.

Philips is a registered trademark of Koninklijke Philips Electronics N.V.

Philips Medical Systems North America Corporation reserves the right to make changes in specifications or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Printed in The Netherlands.

4522 962 08141/861 * NOV 2005

References

1. Bardy GH et al. "Multicenter Comparison of Truncated Biphasic Shocks and Standard Damped Sine Wave Monophasic Shocks for Transthoracic Ventricular Defibrillation". *Circulation* 1996;94:2507-2514.
2. Gliner BE, White RD. "Recurrence of Out-of-Hospital VF Following Low-Energy Biphasic and High-Energy Monophasic Defibrillation Shocks". *J American College of Cardiology (Abstract)* 1999;33:127A.
3. Reddy RK et al. "Transthoracic Ventricular Defibrillation Causes Fewer ECG ST-Segment Changes After Shock". *Ann Emerg Med* 1997;30:127-134.4.
4. Benditt, DG et al. "Biphasic Waveform Cardioversion as an Alternative to Internal Cardioversion for Atrial Fibrillation Refractory to Conventional Monophasic Waveform Transthoracic Shock". *Am J Cardiol*, December 15, 2001; 88(12):1426-1428.