User's Manual

RS PRO Portable Digital Microscope

Stock number: 196-4077


EN

IT

FR

DE

ES

User's Manual

RS PRO Portable Digital Microscope

Stock number: 196-4077

EN


Table of Contents

1.	Precaution.....	2
2.	Preparation.....	6
	● Standard Accessories.....	6
	● Part Names & Functions.....	7
3.	key features... ..	8
4.	Use the Microscope.....	9
	● Install Battery.....	9
	● Insert a Micro-SD Card.....	10
	● Display Icons & Meanings.....	11
	● Turn on the Microscope... ..	12
	● Turn off the Microscope.....	12
	● How to Capture... ..	13
	● Download to Computer.....	19
	● Preview on Computer.....	19
	● Charge Battery... ..	20
	● Output to TV.....	21
	● Mount Tripod.....	22
5.	Menu Introductions.....	23
	● Set Resolution... ..	23
	● Timer	24
	● Quality.....	25
	● Sharpness.....	27
	● Language.....	27
	● Power Save... ..	27
	● Date Stamp... ..	28
	● TV out... ..	28
	● Date/Time... ..	29
	● Format SD.....	29
	● System Reset... ..	30
	● Version... ..	30

6. Specifications	31
7. Computer System Minimum Requirement.....	31
8. Troubleshooting... ..	32

PC Software Introductions 33

Install the software..... 33

Start Microscope... .. 34

Software functions... .. 35

- Capture Photo... .. 35
- Capture Video... .. 35
- Time..... 36

Calibration and Measurement..... 37

- Calibration
- Measurement... .. 41
- Change Calibration..... 42
- Draw and text..... 45

Software Interface Instructions..... 32

- Main Interface Icons..... 46
- Main Menu
- Options... .. 48
- Capture... .. 49
- Measurement Interface Instructions..... 50

Welcome

Thank you for purchasing a Portable Digital Microscope. Please take a moment to read through this manual. Its contents will help you to get the best use of the smart microscope.

1. Precaution

Before using Portable Digital Microscope, please ensure that you read and understand the safety precautions described below. Always ensure that the Portable Digital Microscope is operated correctly.

- Never operate Portable Digital Microscope where operation of such devices is restricted. Improper use creates the risk of serious accident.
- Do not try to open the case of Portable Digital Microscope or attempt to modify Portable Digital Microscope in any way. Maintenance and repairs should be performed only by authorized service providers.
- Keep Portable Digital Microscope away from water and other liquids. Do not use Portable Digital Microscope with wet hands. Never use Portable Digital Microscope in the rain or snow. Moisture creates the danger of fire and electric shock.
- Keep Portable Digital Microscope and its accessories out of the reach of children and animals to prevent accidents or damage to the Portable Digital Microscope.

If you notice smoke or a strange odor coming from Portable Digital Microscope, turn the power off immediately. Take Portable Digital Microscope to the nearest authorized service center for repair. Never attempt to repair Portable Digital Microscope on your own.

- Portable Digital Microscope may become slightly warm during use. This does not indicate malfunction.
- If the exterior of Portable Digital Microscope needs cleaning, wipe it with a soft, dry cloth.


FCC Notice

- This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

CAUTION

To comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules, this device, when needed, must only be connected to computer equipment certified to comply with the Class B limits.

■ Standard accessories

Portable Digital
Microscope

Lithium Battery

AC adaptor

CD

USB cable


TV cable

Stand

Cleaning fabric


Carry bag


Manual

■ Parts Names & Functions

No.	Functions	No.	Functions
1	Capture	12	LEDs brightness dial
2	Zoom in/out; Screen icon display on/off	13	Battery compartment
3	TV-out	14	TFT display
4	USB port	15	UP Color mode switch
5	DC-in	16	Menu setup
6	Aux LEDs	17	OK
7	Lens	18	DOWN
8	Tripod mounting hold	19	LEFT / BACK
9	MicroSD card slot	20	Power on/off
10	Wrist grip hook	21	Speaker
11	Focus wheel	22	Microphone
23	Battery charging indicator		


3. Key Features

- 3 inches TFT panel
- 5M high quality image sensor (up to 12M by interpolation)
- 20x-200x-500x magnifications
- Photo & Video
- Timer
- Measurement (when working on computer)
- Support MicroSD card up to 16GB
- Bundle software PortableCapture Plus for Windows7/8/10, Mac 10.14 or later.
- Lithium battery (rechargeable & changeable)

4. Use the Microscope

4.1 Install battery

- 1). Open the battery door by sliding it outward.
- 2). Insert the lithium battery to the battery compartment according to the polarities shown in the compartment.

Note: Before using the microscope, you must charge its battery. Please refer to section 4.9 on how to charge the battery.

4.2 Insert a MicroSD Card (not included)


You must insert a MicroSD card (not supplied) to record and store your pictures and videos.

Inserting the MicroSD memory card:


- 1) Turn off the microscope
- 2) Insert the MicroSD card into the MicroSD card slot; gently press in until card is latched.

Note: Do not force the card into the slot; forcing may damage microscope and the MicroSD memory card. Align MicroSD card position (as marked) and try again. The MicroSD memory card must be formatted before capture.

See the following instruction.

4.3 Display Icons & Meanings


ICON	Meaning
	Still image capture
	Timer capture
00013/01618	Captured image qty/total qty that could be captured (based on the current resolution)
	MicroSD inserted
5.0M	Current resolution
	Battery electricity indicator

Note: The icon displays can be switched on/off by pressing Zoom Wheel(2)downward vertically.

4.4 Turn on the Microscope

- 1) Unfold the display panel, and then the microscope will power on directly.

- 2) If need to power on when panel unfolded, press POWER ON/OFF button.


4.5 Turn off the Microscope

- 1) Fold the display panel and then the microscope will turn off automatically
- 2) Press ON/OFF button
- 3) By Power Save. Please refer to Section 5.2.2 for Power Save setup.


4.6 How to Capture?

- 1) Turn on the microscope by simply unfolding its display panel to a proper position, normally at an angle of 90 degree with the viewing end as follows:


- 2) Adjust the LED brightness dial to get a proper brightness for the LEDs at the viewing end.


- 3) Hold the microscope main body by right hand; put

the viewing end closely to target object; adjust focus wheel to get clear image by left hand.


- 4) When a clear image got, hold the microscope steadily and press Capture button (1) one time.

5) Record Video

(1) Press LEFT/BACK (19) to enter video mode as follows. There is video icon on the left corner.

(2) Press Capture button (1) to start recording video. Press again to stop.

6) Play Back

(1) Press LEFT/BACK (19) to enter play back mode as follows. There is review icon on the left corner.


(2) Press UP/DOWN to select among captured pictures and video clips.

(3) Delete

At play back mode, press  to enter sub-menu as follows, then confirm delete.


(4) Protect

At play back mode, press  to enter sub-menu as follows, then confirm protect pictures from delete.

(5) Slide Show

At play back mode, press  to enter sub-menu as follows, select interval time, and then confirm by OK.


7) Color Mode


At preview status, press   to switch among Color/Positive/Negative.

8) **Zoom in:** push Zoom Wheel upwardly

9) **Zoom out:** push Zoom Wheel downwardly


Note: The zoom ranges from 1x – 2.5x

4.7 Download to Computer

- 1) Connect the microscope to computer by USB cable.
- 2) Choose Mass Storage and then you can view stored photos & videos.

4.8 Preview on Computer

- 1) Connect the microscope to computer by USB cable.
- 2) Choose PC Camera mode
- 3) Start PortableCapture software. Please refer to Part II on how to use the microscope by PortableCapture software.


Note:

- 1) When working on computer, the button functions are no longer valid, but most of functions are available by PortableCapture software.
- 2) Voice recording is not available in PC camera mode
- 3) Images and videos captured from PC camera mode are not stored in micro SD card, they store in the PC.
- 4) The microscope can be used as a webcam device for microscopic live video.

4.9 Charge Battery


Before using the microscope, you must charge its battery.

1. Plug the small end of the power adaptor into the DC-in port located on the bottom side of the unit. Then plug the other end of the power adaptor into a standard AC outlet. The battery indicator LED at the end corner of the unit will be in red to show that the battery is charging.

2. When the battery is finished charging, the battery indicator LED will be in off. Unplug the unit.
3. You can now use your microscope for up to four hours of continuous use. When the battery is low, there will be low-power indicator on the right down corner of the panel, front LED light will become weak and image quality becomes poor. Then you have to connect the power adaptor and plug it in to recharge the battery.


4.10 Output to TV

1. Insert the 1/8-inch plug on the end of the supplied video cable into the microscope's Video port.
2. Insert the RCA plug on the video cable into one of your television's yellow video input connectors. Refer to your television's manual for more information on locating these connectors.
3. Turn on the television. If necessary, switch your television to the appropriate input channel. Items that you magnify appear on the television screen.

4.11 Mount Tripod

1) Put the stand end with screw toward the mount hole on the microscope properly and rotate the end into the hold carefully.

Do not force the end into the hole; forcing may damage the microscope.

2) Adjust the height of the stand at a proper position.

5. Menu Introductions

5.1 Main Menu I

Press  to enter main menu I as follows:

Here you can setup Resolution/Timer/Quality/Sharpness.

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2) Press OK to confirm selection.

5.1.2 Timer

1) From the main menu I , choose Timer and then enter Timer setup:

Frames: input quantity of frames you want to take by UP/DOWN and LEFT/RIGHT.

The above setup means 32 frames.

Interval: input interval time between frame capture by UP/DOWN and LEFT/RIGHT.

The above setup means interval time is 15 seconds between each two captures.

- 3) Press OK to confirm selection and press BACK (19) back to preview status. Then press Capture (1) to start auto capture. When at auto capture status, there is

at the left corner of the screen.

- 4) When all frames captured, microscope will return to normal preview.
 5) If you need to stop auto capture during the capture process, press Capture one time.

5.1.3 Quality

From the main menu I , choose Quality and then select among Super Fine/Fine/Normal.

The three selections will result in different image compress ratio.


5.1.4 Sharpness

From the main menu I , choose Sharpness and then select among Super Fine/Fine/Normal.

5.2 Main Menu II

Press  to enter main menu I, and then press LEFT to enter Main Menu II as follows:

Here you can setup Language/Power Save/Date Stamp/TV Out/Date and Time/Format SD/System Reset/Version.

5.2.1 Language

From the main menu II, choose Language and then select a target language by UP/DOWN.

5.2.2 Power Save

1) From the main menu II, choose Power Save and then select among the following a target time for the microscope to power off without key operation.

5.2.3 Date Stamp

1) From the main menu II, choose Date Stamp to enter the following menu:

2) Select Date/Time, Date or Off.

Date/Time: both date and time will be marked on images taken.

Date: on date will be marked on images taken

Off: None of date and time info will be marked on images taken.


5.2.4 TV Out

From the main menu II, choose TV Out and then select, between NTSC and PAL, the correct mode compatible with your TV.

5.2.4 Date/time


From the main menu **II** , choose Date/Time to setup date and time.

- 1) Setup date and time by LEFT/RIGHT/UP/DOWN
- 2) Setup date mode among YY/MM/DD, DD/MM/YY, and MM/DD/YY.

5.2.5 Format SD

From the main menu **II** , choose Format SD to format the MicroSD card inserted.

5.2.6 System Reset

From the main menu **II**, choose Sys. Reset to reset the system back to default.

5.2.7 Version

From the main menu **II**, choose Version to view the current machine version.

6. Specifications

Image sensor:	5M high quality Image Sensor
Magnifications:	20x-200x-500x
Resolutions:	12M,9M,5M,3M,1.3M,VGA
Lens:	High quality microscopic lens
Display:	3inches TFT panel, 4:3 ratio
Light sources:	8 LEDs with adjustable brightness
Power sources:	Li-ion 3.7V / 1050mAh battery Working time: 4 hours; Charging time: 5 hours
Frame rate:	30fps
Adaptor:	Input:100-240V, 50/60Hz; Output: 5V, 1A
OSD languages:	English/Germany/French/Spanish/Italian/ Portugese/Japanese/Chinese
Bundle software	Win/Mac: PortableCapture Plus with measurement
Software languages:	English/Germany/French/Spanish
Product size:	130x103x29mm
Weight:	195g (with battery)

7. Computer System Minimum Requirement

Operating system	Windows 7/8/10, Mac 10.14 or later
CPU	Pentium 11 or above
RAM	512MB
Interface	USB port2.0

8. Troubleshooting

Problems	Solution
The unit does not turn on.	Charge the battery.
The screen is blank.	<p>Ensure that the power is on.</p> <p>Charge the battery.</p> <p>Verify that the item you want to magnify is positioned correctly under the lens.</p>
There are smudges or blemishes on the screen	Clean the screen and lens carefully with soft fabric
Microscope is connected to a television, but no image appears on the television screen.	<p>Ensure that the television is on.</p> <p>Switch to the appropriate video input channel on your television.</p> <p>Ensure that the video cable is securely connected to the Microscope Video-out port and plugged into the correct video input on your television.</p>
<p>The battery indicator LED is not in red, even though the power adaptor is connected and plugged into a wall socket.</p> <p>-or-</p> <p>The battery is not charging.</p>	<p>Plugged into the Microscope's DC-in port.</p> <p>Ensure that the outlet has power and is not malfunctioning. If the outlet is connected to a wall switch, make sure that the switch is on.</p>

■ PC Software Introduction

1. Install the Software on PC (Win/Mac)

Insert the driver CD into CD-ROM Drive and then open the CD.

Double click to install **PortableCapture Plus** and the installation wizard will guide you through the whole process.


Note:

The software **PortableCapture Plus** is also available at www.cameradownload.net.

2. Start Microscope

Connect your Microscope to your PC USB port, start the

software by clicking the  icon generated on the desktop after installation and you will see the following screen pop up.

When the Microscope is disconnected from your PC USB port, the following message will pop up.


**No Device detected, please
connect your Microscope
directly to your PC USB port.**

3. Software Functions

■ Capture Photo


- Click 'Photo' shortcut icon in the main menu bar.
- Click 'Capture > Photo'.

■ Capture Video

- Click 'Video' shortcut icon in the main menu bar.
- Click 'Capture > Video'.


■ Timer (Photo / Video)

Click 'Timer' shortcut icon in the main menu bar.

* Under Timer photo, you can set the start time, interval time and total picture number.

* Under Timer video, you can set the start time, interval time and total video number.


■ Calibration and Measurement

Please make calibration before measurement.

Calibration

Display Crossing from 'Options>Crossing'.

Microscope default at crossing “on”


It's used to check the dial on the calibration ruler.


1. Focus microscope on the attached calibration ruler


and then adjust the magnification till the picture is clear. Take a picture for the selected area when the dial is parallel to the crossing as below.

2. Double-click on the picture to enter into measurement interface as below.


3. Click Calibration icon , and then move mouse to the picture. Click on the beginning and ending points across a known value. Now a dialog will popup as below and you need to enter the exact value of the measurement sample into 'Actual dimension'. Then the software will automatically calibrate the magnification rate. Click 'OK' and the adjusted magnification rate will show in the upper right corner.


* Remark: The measured line should be paralleled.

4. Now you can use 'Any angle Line' under Measurement

icon  to measure the value again. If the value you measure by software is nearly the same as the exact value of the dial means the calibration finish.

Measurement

Double-click on the picture to enter into measurement interface. Under the 'Measurement' icon , you can choose any desired options to measure.

1. Any Angle Line Measurement

Simply click from one point and drag to the next point, and then release the click.

2. Continuous Line Measurement

The Continuous Line measurement allows you to measure any continuous desired length.

Simply click and drag from one point to other desired points, and then release the click.

3. Radius Circle Measurement

Simply click and extend out to the desired radius, the radius, circumference and area of the circle will show up accordingly.

4. Diameter Circle Measurement

Simply click and extend out to the desired diameter, the diameter, circumference and area of the circle will show up accordingly.

5. Three Points Angle Measurement


Simply define any three points angle and click, the angle value will then display.

Change Calibration


If you need to change the magnification rate and adjust focus to capture picture, please calibrate the magnification rate again.

1. Double click to open the picture, Click ' > Reset magnification > Reset picture magnification' as below. And then start calibration again (Calibration details please refer to page24).

2. Change the magnification rate and capture several pictures, click on the picture to enter into measurement

interface. And click  > Set magnification > Set picture magnification', now a dialog will show up as below; enter the magnification rate into the dialog and click 'OK'.


3. If you want to use the same magnification rate to

capture pictures, just click ' > Set magnification > Set capture magnification'. Now all pictures will use the same magnification rate and you can start measurement directly.


Draw and Text


You can draw or add a text on the picture and save for later analysis.

1. Click icon to choose any kinds of drawing under the pull-down manual as below.

2. Click icon and set the font for the text in the popup dialogue as below.

Main Interface Icons

Capture photo – It's used to capture photo one by one.

Start / Stop record video.


Timer for photo / video – Set the start time, interval time and total photo/video number.

Full screen mode,
Press ESC or double-click to exit.

Main Menu

File


Photos Directory: set the directory to store photos taken;

Videos Directory: set the directory to store video taken;

Options

Resolution: change preview size

Date/Time: set to display system date/time on preview window

Language: choose different OSD language

Crossing: display cross on preview window

Full-screen mode: enter into full screen. Click ESC or double click to exit

Capture


Photo: take a photo

Video: record a video clip

Measurement Interface Instructions

Save as

Redo

Copy to
clipboard

Draw

Copy file

Text note

Previous

Measurement

Next

Calibration

Undo

This symbol on the product or in the instructions means that your electrical and electronic equipment should be disposed at the end of its life separately from your household waste. There are separate collection systems for recycling in the EU.


For more information, please contact the local authority or your retailer where you purchased the product.

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