



SG Diagnostics

Committed to Quality of Health



SG Diagnostics COVID-19 Antigen Rapid Test Kit (for self-testing)

The COVID-19 Antigen Rapid Test Kit is an *in vitro* assay for the qualitative detection of antigens from SARS-CoV-2 in nasal or throat swab sample. It is intended to aid in the rapid diagnosis of a COVID-19 infection.

How it helps you



Rapid detection time of less than 15 minutes



97.57% sensitivity and 99.29% specificity



Able to detect asymptomatic patients



Able to detect new COVID-19 variants



Easy to store at room temperature

COVID-19 Antigen Rapid Test Kit is a **fast** and **accurate** test that can be done by anyone, anywhere.

- On-site testing, no lab needed.
- Hassle-free with all essential materials provided
- Gives quick result within 15 minutes

Suitable to use at :



Airports



Corporations



Events



Hotels



Homes



Schools



Scan to watch
Instructional Video


Steps for COVID-19 Antigen Test

1




Carefully insert the swab into the nasal or throat. Rotate the swab at least 5 times. Withdraw the swab carefully.

2




Insert the swab into the elution tube and roll the swab at least 10 times while pressing the swab tip against the bottom and side of the elution tube. Squeeze the elution tube as you remove the swab. Use the nozzle attached to the tube to cover the elution tube tightly.

3



10 - 15 minutes



Positive

Negative

Invalid

Remove the test cassette from the sealed pouch just prior to the testing and lay flat on a table. Add 4 drops of test sample into the sample well (●) by squeezing the elution tube. Wait for the colored band(s) to appear. The result should be read after 10 minutes. Do not interpret the result after 15 minutes.

Box Content

- Antigen test cassette
- Sample eluent
- Disposable swab
- Instructions for use

Specification

20 tests/kit, 5 tests/kit, 2 tests/kit

Specimen Type

Nasal or throat swab

Total Time

15 min

Storage Condition

4°C - 30°C, sealed