Instructions for Use

ADORA™

Round Microtextured Mammary Implants

904748 v1 - 11/2025

EUROSILICONE™

CONTENTS

EN	Round Microtextured Silicone Gel-Filled Mammary Implants.....	3
FR	Implants mammaires ronds micro-texturés remplis de gel de silicone	8
DE	Runde, mikrotexturierte, mit Silikongel gefüllte Brustimplantate.....	13
IT	Impianti mammari rotondi microtesturizzati riempiti in gel di silicone	18
ES	Implantes mamarios redondos microtexturizados rellenos de gel de silicona	23
PT	Implantes Mamários com Enchimento de Gel de Silicone Microtexturizados Redondos	29
DA	Runde mikroteksturerede silikonegelfyldte brystimplantater.....	34
NL	Ronde microgetextureerde siliconengelgevulde borstimplantaten	39
SK	Okrúhle mikrotextúrované prsné implantáty plnené silikónovým gélom.....	44
TR	Yuvarlak Mikro Dokulu Silikon Jel Dolgulu Meme İmplantları	49
PL	Okrągłe implanty piersi z mikroteksturą i wypełnieniem z żelu silikonowego.....	54

EN

Round Microtextured Silicone Gel-Filled Mammary Implants

Only to be implanted in an appropriate medical environment by appropriately trained medical doctors who are qualified or accredited in accordance with national law.

INTENDED PURPOSE/USE

The Round Microtextured Silicone Gel-Filled Mammary Implants is intended for cosmetic breast augmentation surgery or breast reconstruction surgery.

DESCRIPTION

The Round Microtextured Silicone Gel-Filled Mammary Implants is a sterile, single-use medical device that is surgically implanted to replace tissue lost by trauma, surgical intervention or to augment tissue for therapeutic purposes. The round silicone elastomer shell is produced with the outer surface micro textured to offer a disruptive contact surface for collagen interface and filled with a cohesive gel. Each gel-filled mammary implant is supplied in a cardboard box and packaged within a sealed, double primary pack. The packaging also contains removable labels indicating the characteristics of the implant (Lot number, Reference number), implant card, patient and user leaflets.

SUMMARY OF REPORTED PARAMETERS FOR SURFACE CLASSIFICATION

Surface finish classification	Roughness grade	Sa	Scx
SLO	Micro-textured (<50)	37 µm	0.22

LINK TO THE IFU ON THE MANUFACTURER'S WEBSITE

<https://www.gcaesthetics.com>.

MEDICAL INDICATIONS

The Round Microtextured Silicone Gel-Filled Mammary Implants are indicated for

- Surgical Augmentation and contour correction of congenital anomalies of the breast,
- Surgical Reconstruction of the breast following subcutaneous mastectomy and other suitable mastectomy procedures or trauma,
- Surgical correction of Combined breast and chest wall abnormalities,
- Replacement of breast implants for medical reasons.

COSMETIC INDICATIONS

The Round Microtextured Silicone Gel-Filled Mammary Implants are indicated for

- Cosmetic Augmentation Surgery,
- Replacement of breast implants for cosmetic reasons.

CONTRA-INDICATIONS

The use of these implants is contraindicated in patients who have one or more of the following conditions:

- FIBROCYSTIC disease,
- INSUFFICIENT TISSUE covering due to radiation damage on the chest wall, tight thoracic skin grafts or radical resection of the pectoralis major muscle,
- EXISTING LOCAL OR METASTATIC CARCINOMA of the breast,
- Suppressed/compromised immune system,
- A history of sensitivity to foreign materials or atopy,
- PHYSIOLOGICALLY/PSYCHOLOGICALLY unsuitable patient,
- ACTIVE INFECTION anywhere in the body,
- Recent history of breast abscess,
- A history of compromised wound healing or regional swelling,
- KELOID SCARS,
- Any other serious medical condition.
- These devices are not intended for persons who are less than 18 years old or pregnant or nursing women

INTENDED USER

The Round Microtextured Silicone Gel-Filled Mammary Implants is intended to be used by suitably qualified surgeons who have received appropriate training on how to use the device safely.

IMPORTANT: The surgeon is responsible for conducting the medical assessments regarding the patient's eligibility for implantation and for deciding on the most suitable surgical technique for both the patient and the type of implant chosen. The user should consider any specific risks that may be applicable to activities of the patient (for example profession, sports or other activities regularly performed by the patient) before the implantation. The manufacturer recommends the surgeons to perform a post implantation monitoring in order to identify any potential undesirable side-effects.

The recommendation for the patient has to consider any previous procedures, accidents, conditions, medications, or other simultaneous treatments of the consumer that may affect the procedure (for example skin, diseases, traumas, and auto-immune diseases).

TARGET POPULATION

The Round Microtextured Silicone Gel-Filled Mammary Implants are intended to be used in adult populations (excluding pregnant and nursing women), with no known contraindications.

INTENDED CLINICAL BENEFITS (FOR MEDICAL PURPOSE)

The clinical benefits of the Round Microtextured Silicone Gel-Filled Mammary Implants are as follows:

- Replace breast tissue that has been removed following mastectomy
- Replace breast tissue that has been removed due to trauma
- Enhance breast tissue affected by tuberous or congenital breast abnormality

EXPECTED EFFECTS (FOR NON-MEDICAL PURPOSE)

The expected effects of the Round Microtextured Silicone Gel-Filled Mammary Implants are as follows:

- Improve patient's quality of life (based on self-esteem, body image, sexual satisfaction)
- Improve patient's psychology in terms of levels of depression and anxiety

PERFORMANCE CHARACTERISTICS

The Round Microtextured Silicone Gel-Filled Mammary Implants are designed and manufactured according to the current EN ISO 14607 standard concerning mammary implants. The performance characteristics of the Round Microtextured Silicone Gel-Filled Mammary Implants in terms of:

- Shell integrity
- Implant resistance
- Compatibility between the filling gel and the shell
- Diffusion from the implant
- Volume

- Surface texture
- Biocompatibility

Had been validated and comply with the requirements of EN ISO 14607.

MATERIALS

Shell and patch: Silicone elastomer, Filling: Silicone gel

LINK TO THE SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The summary of Safety and effectiveness will be available on the Eudamed Website when it is live <https://ec.europa.eu/tools/eudamed/#/screen/home> and the company Website under section entitled "Additional Helpful Information" via the link (<https://www.gcaesthetics.com/reconstruction-breast-surgery>)

For access via Eudamed:

- Click on the link above for the Eudamed home page
- Click on the module called Certificates (Issued or Refused)
- Within the search criteria under Certificate number box insert MDR 740809 and press Search at the bottom of the page
- On the right-hand side of the page under View Detail click on the eye icon.
- On the left-hand side of the page click on SSCP
- Download the SSCP master document by clicking on the link

PATIENT INFORMATION

Any surgical procedure can have complications and risks. Mammary implant surgery is known to provide psychological satisfaction to patients, but like any surgical procedure it can have potential complications and risks. Mammary implantation is an elective procedure, and the patient should be well counselled on the risk/benefit relationship, by the surgeon. Each of the possible complications and warnings should be discussed with the patient prior to the decision to proceed with surgery.

Each patient should receive the patient information leaflet (provided together with the products) and/or the patient information booklet during their initial consultation, to allow time for the patient to read and understand the information on risks, follow-up recommendations and benefits associated with breast implants in order to make an informed decision about proceeding with the surgery. The patient information booklet also contains a 'Patient Informed Consent Form' which allows the patient to retain copies of their implant details such as lot number for their records. This document is available in PDF format from our website www.gcaesthetics.com.

LIFETIME

The patient should be informed that this device is not a lifetime device and may require removal or replacement, which could entail another surgical procedure. Patients must be informed that the implant must not be considered definitive due to possible individual and potential physiological reactions, the inherent to the design of silicone implants or the various implant procedures. The implant may be removed or replaced at any time, which requires revision surgery. The manufacturer guarantees no implant duration period for its implants.

From clinical PMCF study on equivalent device, the cumulative Kaplan-Meier risk of first occurrence of explantation for patients within 10 years post-implantation was 17.6% (95%CI = 16.9 – 18.3%).

MRI COMPATIBILITY

Although Round Microtextured Silicone Gel-Filled Mammary Implants has not been specifically tested for use in an MRI, please note that Round Microtextured Silicone Gel-Filled Mammary Implants are all manufactured from medical implant grade silicone materials which are compatible with MRI Scan. The implantable silicone grade material is the same as that used in other manufacturers' silicone implants where patients have undergone MRI Scans, and no compatibility issues have been reported to date.

STEROID USE

The patient should be informed to consult a physician before using steroid drugs in the implant area to avoid extrusion of the implant.

BREAST FEEDING

The patient should be informed of potential complications with breast feeding following the implantation. The surgery procedure may impact milk production during nursing. Specifically, a peri-areolar incision can significantly reduce the possibility of successful breastfeeding.

INTERFERENCE

During a mammography, the radiologist must be informed that there is an implant in the breast in order to adapt the mammography pressure accordingly. The surgeon or doctor must be notified about the implant in the event of any surgical intervention in the breast area.

The presence of breast implants could delay the detection of breast cancer by self-exam. For this reason, patients must be informed that they need to consult a surgeon and/or doctor for appropriate medical monitoring, as well as regular breast cancer screening. Patients must also be informed that they need to consult a doctor if they suspect complications particularly in the event of trauma or compression caused, for example, by violently massaging the breasts, a sporting activity or wearing seat belts.

IMPLANT CARD

After surgery, the surgeon should provide to the patient:

- the implant card which should be filled out according to the leaflet Instructions for Healthcare Professionals filling out the manufacturer International Implant Card.
- the patient leaflet which is available on the Information website for patients referenced on the card.

The patient should be informed that the implant card must be kept on their person at all the times and enables:

- the patient to identify the implanted devices and to get access to other information related to the implanted device (e.g. via EUDAMED, and other websites).
- the patients to identify themselves as persons requiring special care in relevant situations e.g. security checks.
- emergency clinical staff or first responder to be informed about special care/needs for relevant patients in case of emergency situations.

WARNINGS

ALTERATION OF IMPLANTS

No alteration should be made to the mammary implant before implantation. Alteration of implants voids all warranties, express or implied.

DAMAGE OF IMPLANTS - BREAKAGE AT SURGERY

Extreme care should be taken in the use and handling of implants to minimize the potential for breakage of shells. All implants have been produced by established manufacturing techniques and under strict quality control standards, but there may be some breakage of implants during handling or in surgery, both in initial and any subsequent surgeries.

Extreme care must be taken to avoid unintentional damage to the implant during the implantation or explantation surgery.

- Do not contact the implant with sharp surgical instruments or devices such as scalpels, forceps, hemostats, suture needles and hypodermic needles.
- Do not contact the implant with blunt surgical instruments such as clamps, retractors and dissectors.
- Do not contact the implant with cautery devices.
- Do not use excessive manipulation, force or stress.

Implants should be carefully inspected for structural integrity prior to use. Damaged products should not be implanted, do not attempt to repair damaged products.

An extra product should be available at the time of surgery in the event of damage to, or contamination of, the implant. Recommended procedures for testing, examination and handling of products should be meticulously followed to assure proper use of implants. Patients should be instructed to inform other treating physicians of the presence of implants to minimize the risk of damage.

INTERFERENCE WITH MAMMOGRAPHY

Standard positioning techniques have shown significant limitations when used for imaging augmented breasts. The implant can interfere with the detection of early breast cancer through mammography by obscuring some underlying breast tissue and/or by compressing underlying tissue which can 'hide' suspicious lesions in the breast. Patients should be instructed to request radiologists who are experienced with the most current radiological techniques and equipment for imaging breasts with implants, and to inform their radiologists of the presence, type, and placement of implants.

RUPTURE OF THE IMPLANT

Gel Rupture: Gel rupture can occur post-operatively due to damage during handling or surgery. Some ruptures may remain undetected unless surgery is performed for another reason (e.g., size exchange). Despite the cohesive properties of the gel material, extravasation from the surgical pocket can occur under pressure, potentially requiring additional surgery to retrieve the gel. Inflammation and formation of silicone granulomas have been reported.

Shell Rupture: Rupture of the implant shell can occur due to contracture, trauma, or excessive manipulation. If shell rupture is suspected, the implant should be removed. While improvements in implant shell and silicone gel design have minimized the risk of silicone bleeding through an intact shell, some silicone gel may still diffuse or "bleed" under certain conditions.

SINGLE USE

Implants are intended for SINGLE USE ONLY. DO NOT REUSE EXPLANTED PRODUCTS. DO NOT RESTERILIZE ANY PRODUCT.

Explanted products should not be reused because recleaning and re-sterilization procedures may not adequately remove biological residues, such as blood, tissue and other matter, which could retain resistant pathogens and may also affect the performance of the implant.

USE OF DRUGS

The manufacturer can neither predict, nor warrant, the safety of use of any drugs, including, but not limited to, aesthetic, steroid, antibiotic and vitamin solutions. If such use is contemplated, the appropriate drug manufacturer should be consulted.

COMPLICATIONS

Possible complications known to be associated with gel-filled mammary implants are as follows:

These include but are not limited to: Asymmetry; Calcium Deposits; Capsular Contracture; Compartment Syndrome; Delayed Wound Healing; Dissatisfaction with results ; Double capsular; Dysesthesia; Extrusion; External capsulotomy; Fibroadenoma; Galactorrhea; Gel bleed perspiration; Gel breakage; Granuloma; Hematoma and oedema; Inflammation/Erythema; Infection; Late seroma; Lymphadenopathy; Lymphorrhea; Malrotation of the implant; Nerve injury; Pain; Patient intolerance to any foreign implant; Phlebitis superficial; Ptosis; Reaction to medication; Rupture of the implant; Scar hyperpigmentation and hypertrophy; Silicone allergy; Silicone migration; Displacement of implant; Superficial wound; Upper limb lymphedema; Wrinkling... leading to explantation of mammary implant.

NB: Each patient shall receive the patient information leaflet (provided together with the products) and the patient information booklet during their initial consultation, to allow time for the patient to read and understand the information on risks, follow-up recommendations and benefits associated with breast implants in order to make an informed decision about proceeding with the surgery.

DEFINITIONS FOR THE SURGEON:

See below for definitions of specific complications explained for the surgeon:

ACCIDENTAL RUPTURE OF IMPLANTS IN SURGICAL POCKET

If this occurs, the gel mass may be penetrated with index finger of double gloved hand, whilst exerting pressure on the breast with other hand. The gel can be manipulated out of the cavity, within the hand, and the outer glove can then be pulled over the mass for disposal. Blot the surgical pocket with gauze sponges. Isopropyl alcohol will aid the removal of gel from instrumentation.

ASIA SYNDROME

Autoimmune/Inflammatory syndrome induced by adjuvants. Patients can present myalgia, myositis, muscle weakness, arthralgia and/or arthritis, chronic fatigue, neurological manifestations, cognitive impairment, pyrexia, dry mouth. Diagnosis is confirmed by tests and markers of immune response. Removal of inciting agent generally induces improvement.

In the case that your patient is diagnosed with ASIA syndrome, a report should be made to the manufacturer via the e-mail address: productsupport@gcaesthetics.com and an individual treatment should be initiated.

BREAST IMPLANT ASSOCIATED ANAPLASTIC LARGE CELL LYMPHOMA (BIA-ALCL) (and other lymphomas)

European safety information, the US FDA and current scientific literature have identified an association between breast implants and the development of a cancer of the immune system called **breast implant-associated anaplastic large cell lymphoma (BIA-ALCL)**. *This specific entity is included in the WHO 2016 classification under the terminology "ALCL – AIM". This cancer occurs more commonly in patients with textured breast implants rather than smooth implants and typically develops many years after breast implant surgery. At this time rates of BIA-ALCL are considered to be low.*

In most patients, BIA-ALCL is found within the fluid or scar tissue next to the implant. The main symptoms of BIA-ALCL are persistent swelling aforementioned symptoms, or if BIA-ALCL is suspected, it is advised that patients should be assessed to rule out peri-implant ALCL. When testing for BIA-ALCL, collect fresh seroma fluid and representative portions of the capsule and send for pathology tests to rule out ALCL. Diagnostic evaluation should include cytological evaluation of seroma fluid with Wright Giemsa stained smears and cell block immunohistochemistry testing for cluster of differentiation (CD) and Anaplastic Lymphoma Kinase (ALK) markers. Most cases of breast implant-associated ALCL are treated by removal of the implant and the capsule surrounding the implant and some cases have been treated by chemotherapy and radiation. Although treatment is usually successful BIA-ALCL can be fatal if left untreated, hence, early detection and treatment is key.

In the case that your patient is diagnosed with BIA-ALCL, a report should be made to the manufacturer via the e-mail address: productsupport@gcaesthetics.com and an individual treatment program should be initiated by a multidisciplinary team.

BIA-SCC

Breast implant-associated squamous cell carcinoma is a rare complication of breast implantation that can result in significant morbidity and mortality. BIA-SCC is a rare, epithelial-based tumour distinct from BIA-ALCL. Patients can present with unilateral breast seroma, fluid collection, pain, erythema, and capsular contracture (Baker Grade IV).

Initial treatment typically involved removal of the implant and the capsule surrounding the implant and some cases have been treated by chemotherapy with limited success. Aggressive first-stage surgery seems to have a better prognosis and outcome. BIA-SCC can be fatal, early detection and aggressive treatment is key.

In the case that your patient is diagnosed with BIA-SCC, a report should be made to the manufacturer via the e-mail address: productsupport@gcaesthetics.com and an individual treatment program should be initiated by a multidisciplinary team.

BREAST IMPLANT ILLNESS

A small proportion of women who have breast implants, for both aesthetic and reconstructive purposes, self-identify as having a number of symptoms they believe arise from the presence of their breast implants. Whilst not a medical diagnosis, they refer to their symptoms as 'Breast Implant Illness' (BII).

These multiple generic symptoms are varied, and all cases have been self-diagnosed and self-reported. They include, but are not limited to, flu- like symptoms such as extreme tiredness, brain fog, joint aches, immune-related symptoms, sleep disturbance, depression, hormonal issues, headaches, hair loss and chills.

There are a variety of other reasons these symptoms might be experienced, including background illnesses or hormonal changes. In addition, there are a number of scientific studies investigating similar symptoms experienced by women in the general population with and without breast implants. On average, around 50% of women who self-identify as having 'BII' feel that their symptoms improve after removal – sometimes temporarily and sometimes permanently. It therefore appears that removing breast implants does not necessarily improve symptoms in everyone. To date, there is no research demonstrating which symptoms may or may not improve with implant removal.

Currently, there are no tests that can confirm 'BII'. Research continues in this area, specifically around patients who have autoimmune disease or a predisposition to autoimmune disease. Breast implantation is contraindicated in women with a suppressed or compromised immune system – see CONTRAINDICATIONS section above. Medical advice should be sought if patients suspect they have BII. Their symptoms may not be related to the implants and other medical investigations should not be overlooked or ignored in determining the cause of their symptoms.

CARCINOMA DETECTION

Carcinoma is a type of cancer that originates in epithelial tissues and the implant can interfere with the early detection. It is essential for patients to have regular follow-up examinations to monitor for any signs or symptoms that may indicate the presence of carcinoma. In the case that your patient is diagnosed with carcinoma, a report should be made to the manufacturer via the e-mail address: productsupport@gcaesthetics.com and an individual treatment program should be initiated by a multidisciplinary team.

INFECTION

Pre-existing infection not resolved before implant placement increases the risk of periprosthetic infection. Do not expose the implant or filling accessories to contaminants, which increases the risk of infection.

Infection is an inherent risk following any type of invasive surgery. Infection around a breast implant may occur within days, weeks, or even years, after surgery. Signs of acute infection reported in association with breast implants include erythema, tenderness, fluid accumulation, pain and fever. Signs of subclinical infection may be difficult to detect. Postoperative infections should be treated aggressively according to standard medical practices to avoid more serious complications. Infection that is unresponsive to treatment or necrotizing infection may require implant removal. Capsular contracture may be related to infection in the area surrounding the implant.

LYMPHOMA DETECTION

Lymphoma is a type of cancer that originates in the lymphatic system and the implant can interfere with the early detection. It is essential for patients to have regular follow-up examinations to monitor for any signs or symptoms that may indicate the presence of Lymphoma. In the case that your patient is diagnosed with lymphoma, a report should be made to the manufacturer via the e-mail address: productsupport@gcaesthetics.com and an individual treatment program should be initiated by a multidisciplinary team.

INSTRUCTIONS FOR USE

STORAGE CONDITIONS

The Round Microtextured Silicone Gel-Filled Mammary Implants must be stored in normal conditions, must be protected from humidity and direct sunlight and must be stored with the arrow symbols pointing up. Where these conditions are followed, Eurosilicone products have a shelf life of five years.

SUPPLIED STERILE

Implants are supplied in sterile form (Ethylene Oxide Sterilization or Dry Heat Sterilization, corresponding symbol on product label), processed by validated strictly controlled sterilization cycles with ethylene oxide or dry heat. Sterility is verified in accordance with standards. Sterility of the implant is maintained only if the package is intact and undamaged. If the sterile packaging has been damaged or unintentionally opened before use DO NOT USE. Authorization must be requested to return the product.

SINGLE USE

Implants are intended for SINGLE USE ONLY. DO NOT REUSE EXPLANTED PRODUCTS. DO NOT RESTERILIZE ANY PRODUCT. Explanted products should not be reused because recleaning and re-sterilization procedures may not be adequate to remove biological residues, such as blood, tissue and other matter, which could retain resistant pathogens and may also affect the performance of the implant.

PACKAGING

Sterile product is supplied in a sealed, double primary package. Sterility is not guaranteed if the package has been damaged or opened. Tear-Off patient record labels are attached to the primary package. These labels are to be attached to the patients' records. The manufacturer recommends that the patient is provided with one of these labels to maintain the traceability of the implanted product.

TO OPEN PACKAGED STERILE PRODUCT

- Peel open outer package under clean, aseptic conditions, over sterile field, allowing sealed inner package to fall gently into the field.
- Attach patient record portion of the inner label to patient's chart.
- Using aseptic precautions, peel open the inner package.

PRODUCT EXAMINATION AND HANDLING

- Product should be visually examined for any evidence of particulate contamination, damage or leakage.
- The implant should be kept submerged in sterile water or normal saline prior to implantation to prevent contact with airborne and surgical field particulate contaminants.
- Do not immerse the implant in Betadine or Iodine containing solutions. If Betadine or Iodine containing solutions are used in the pocket, ensure that it is rinsed thoroughly so no residual solution remains in the pocket.
- Gel-filled mammary implants may contain air bubbles. This occurs normally during sterilization and has no effect on product integrity or performance.
- Use disposable, non-shedding wipes only to ensure there are no particles left on the surface of the implant.

CONTAMINATION OF IMPLANTS

Care must be taken to prevent surface contaminants such as talc, dust and skin oils from coming into contact with the implant. Products should be inspected for contamination prior to insertion. Contamination at the time of surgery increases the risk of periprosthetic infection, and possibly capsular contracture. Textured silicone surfaces are potentially more susceptible to contamination than smooth silicone surfaces. Extra care should be taken in the handling of textured surface implants. The manufacturer accepts no liability for products contaminated by other substances after the product leaves our possession

SURFACE CONTAMINANTS (TALCS, DUST, LINT, OILS) ON SURFACE OF IMPLANTS CAN CAUSE FOREIGN BODY REACTION, HANDLE WITH CARE WITH SURGICAL GLOVES (RINSED FREE OF TALC) WITH STRICT ASEPTIC TECHNIQUE. DO NOT IMPLANT CONTAMINATED PRODUCT.

Back-up implants must be readily available at the time of surgery for use in the event of contamination.

SURGICAL PROCEDURE

Correct surgical procedures and techniques are the responsibility of the medical profession. Each surgeon must evaluate the suitability of the procedure based upon current accepted techniques, individual judgement and experience. Proper size and shape of implants must be determined for the individual patient by the surgeon. An incision should be of appropriate length to accommodate the style, size, and profile of the implant. This will reduce the potential for creating excessive stress to the implant during insertion.

SPECIFIC GUIDELINES FOR THE USE OF MAMMARY IMPLANTS

1. Size selection should be decided by base dimensions, the patient's dimensions, and individual breast characteristics rather than preferred volume or projection.
2. The preformed shapes of the breast implants are firmer than other cohesive gel filled implants, and the surgical incisions should be longer than those required for compressible implants.
3. Where appropriate the new inframammary fold may be required to be approximately 1cm below the natural fold.
4. The recommended approach is via the inframammary fold. It is the responsibility of the surgeon to decide if axillary or trans areolar placements may be possible. We caution that insertion and orientation may be more difficult due to the increased firmness of the gel
5. Sub-muscular placement (with possible lower ligament release) is preferred where possible. As the gel is highly cohesive it has a more defined radius on the upper apex margin. Sub-muscular placement reduces the probability of this apex radius being palpable.
6. DO NOT use force during placement of breast implants. Over manipulation of breast implants may result in permanent distortion of the cohesive gel and loss of the implant's shape.
7. Maintenance of correct alignment and avoidance of rotation in the post-operative phase is aided by:
 - a) Achieving a surgical pocket closely aligned to the implant.
 - b) Achieving a "dry" surgical pocket.
 - c) Practicing the procedure of implant stabilization by using post-operative external immobilization.

MATERIALS DEVICE DESCRIPTION

These medical devices are manufactured from medical grade silicone materials which are suitable for long-term implantation. Master Files have been filed with the U.S FDA. The following three tables provide quantitative and qualitative information on the materials and substances which patients can be exposed based on chemical characterization of representative devices. The materials and device have been subjected to biocompatibility testing and evaluation, and risk assessments to demonstrate their biological safety. However, individual responses to chemicals may vary, and all reactions cannot be predicted.

Table 1. Breast Implant Device Materials

Device materials	Implant component
Silicone dispersions	Shell
Silicone dispersions	Closure Patch
Platinum catalyst	Shell and Gel fill
Silicone adhesive	Shell
Organotin catalyst	Shell (silicone adhesive)
Silicone gel	Gel fill

Table 2. Chemicals Released by Breast Implants

Compounds	Whole device (µg/g or ppm) ¹	Compounds	Whole device (µg/g or ppm)
Volatiles² - Chemicals that can be released by breast implants as a gas			
1,1-Dichloroethene	Not Detected	1,2-Dichloroethane	Not Detected
Cis-1,2-dichloroethene	Not Detected	Bromodichloromethane	Not Detected
1,1-Dichloroethane	Not Detected	Toluene	Not Detected
Trans-1,2-dichloroethene	Not Detected	1,1,2-Trichloroethane	Not Detected
Chloroform	Not Detected	Tetrachloroethene	Not Detected
1,1,1-Trichloroethane	Not Detected	1,3-Dichloropropane	Not Detected
Carbon tetrachloride	Not Detected	Dibromochloromethane	Not Detected
1,1-Dichloropropene	Not Detected	Chlorobenzene	Not Detected
Benzene	Not Detected	1,1,1,2-Tetrachloroethane	Not Detected
Ethanol, 2-(trimethylsilyl)-	2,41 max	1,3,5-Trimethylbenzene	Not Detected
Ethylbenzene	Not Detected	Sec-butylbenzene	Not Detected
m/p-Xylene	Not Detected	1,3-Dichlorobenzene	Not Detected
o-Xylene	Not Detected	1,4-Dichlorobenzene	Not Detected
Styrene	Not Detected	1,2-Dichlorobenzene	Not Detected
Bromoform	Not Detected	N-Butylbenzene	Not Detected
Isopropylbenzene	Not Detected	1,2-Dibromo-3-chloropropane	Not Detected
Bromobenzene	Not Detected	1,2,4-Trichlorobenzene	Not Detected
1,2,3-Trichloropropane	Not Detected	Hexachlorobutadiene	Not Detected
n-Propylbenzene	Not Detected	Naphthalene	Not Detected
2/4-Chlorotoluene	Not Detected	1,2,3-Trichlorobenzene	Not Detected
1,2,4-Trimethylbenzene	Not Detected	p-Isopropyltoluene	Not Detected
Tert-butyl benzene	Not Detected		
Total volatiles	2,41 ppm		
Extractables³ - Chemicals that can be released by breast implants following soaking in water and/or solvent			
Cyclic siloxanes (D4, D5, D6...)	Not Detected	Linear siloxanes (L3, L4, L5...)	Not Detected
Total extractables	Not Detected		
¹ The data reported in this table has been calculated cumulatively for shell and gel. ² Not Detected means that the level of the individual volatile was below the quantitation limit of the test method. The quantitation limits data is summarised below:			
Component	Individual quantitation limit	m & p-xylene and 2,4-Chlorotoluene	³ Extractable -Quantitation limit
Shell	0.239 ppm	0.465 ppm	4.645 ppm
Gel	0.3 ppm	0.590 ppm	6.0 ppm

Table 3. Heavy Metals Found in Breast Implants

Heavy metals	Concentration (µg/g or ppm)	Heavy metals	Concentration (µg/g or ppm)
Aluminium	Not Detected	Manganese	Not Detected
Antimony	Not Detected	Mercury	Not Detected
Arsenic	Not Detected	Molybdenum	Not Detected
Barium	Not Detected	Nickel	Not Detected
Beryllium	Not Detected	Platinum	Not Detected
Cadmium	Not Detected	Selenium	Not Detected
Chromium	Not Detected	Thallium	Not Detected
Cobalt	Not Detected	Tin	Not Detected
Copper	Not Detected	Vanadium	Not Detected
Lead	Not Detected	Zinc	0,043
Magnesium	Not Detected		
Not Detected means that the level of the individual element was below 0.002 ppm which is the assumed test method quantitation limit for the shell, and 0,0615 ppm the test method quantitation limit for the gel.			

DISPOSAL OF DEVICE

All mammary implants must be disposed of according to requirements related to medical wastes involving an infectious risk.

REPORTING OF SERIOUS INCIDENT:

Any serious incident (incident that directly or indirectly led, might have led or might lead to any of the following: (a) the death of a patient, user or other person, (b) the temporary or permanent serious deterioration of a patient's, user's or other person's state of health, (c) a serious public health threat) should be reported to the manufacturer via the e-mail address: productsupport@gcaesthetics.com and to the competent authority of the Member State in which the user and/or patient is established.

RETURN OF GOODS POLICY

Prior to the return of any product:

- a) Authorization must be requested in writing.
- b) Decontamination certificate must be supplied with any product returned which has been removed from the packaging.
- c) If it is necessary to puncture the shell of any sealed shell product to assist with safe decontamination sterilization procedures, the area of mechanical interference should be marked with indelible marker on the surface of the product and reference made on the decontamination certificate.
- d) Products File Notes and returned products must be sent to the manufacturer by your distributor only.

WARRANTY

The manufacturer warrants that reasonable care was used in the manufacture of these products and will replace any product that the manufacturer's investigation shows to be defective at time of shipment. Patient selection, surgical procedures, post-surgical treatment and stresses, and handling of the devices are totally the responsibility of the customer. The manufacturer has no control over the conditions of use and cannot warrant good effect or against ill effect following the use of the device and shall not be responsible for any incidental or consequential loss or damage or expenses directly or indirectly arising from the use of this device. All other warranties (whether implied by law or otherwise) are excluded to the fullest extent permitted by law.

CAUTION:

Federal (USA) law restricts this device to sale by or on the order of a physician.


HARMONISED STANDARD AND COMMON SPECIFICATIONS APPLIED

The following common specification and relevant harmonized standards are applied for all mammary implants: Common specification 2022/2346; EN ISO 13485; EN ISO 14971; EN ISO 62366-1; EN ISO 14630; EN ISO 14607; EN ISO 10993 series; EN ISO 15223-1; EN ISO 20417; EN ISO 11607; EN ISO 11135 or EN ISO 20857 (depending on sterilization method).

BIBLIOGRAPHIC REFERENCES

Literature references are available upon request from the manufacturer.

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Eurosilicone SAS