Instructions for Use

LUNA XT

Anatomical Silicone Gel-Filled Breast Implants

EUROSILICONE™

CONTENTS

EN	LUNA XT Micro textured Anatomical Silicone Gel-Filled Breast Implants	3
FR	LUNA XT Implants mammaires anatomiques micro texturés remplis de gel de silicone	9
DE	LUNA XT Mikrotexturierte anatomische mit Silikon gel gefüllte Mammaimplantate	17
IT	LUNA XT Impianti mammari anatomici microtesturizzati riempiti in gel di silicone	24
ES	Implantes mamarios anatómicos microtexturizados rellenos de gel de silicona LUNA XT	32
PT	Implantes Mamários Anatómicos com Enchimento de Gel de Silicone Microtexturado LUNA XT	39
DA	LUNA XT anatomiske silikongefyldte brystimplantater med mikrostruktur	46
NL	LUNA XT Microgetextureerde anatomische met siliconengel gevulde borstimplantaten	53
EL	Ανατομικά εμφυτεύματα στήθους γέλης σιλικόνης με μικροϋφή LUNA XT	60
LV	LUNA XT Anatomiski mikrotekstūras krūšu implantī ar silikongēla pildījumu	68
SK	LUNA XT – Anatomické prsné implantáty s mikrotexúrou vyphené silikónovým géloom	74
RO	Implanturi mamare anatomic cu gel de silicon microtexturate LUNA XT	81
TR	LUNA XT Mikro dokulu Anatomik Silikon Jel Dolgulu Meme İmplantları	88
RU	Анатомические грудные имплантаты LUNA XT с микротекстурой, заполненные силиконовым гелем	95
CS	Anatomické prsní implantáty vyplněné silikonovým gelem LUNA XT s mikrotexurou	103
ET	LUNA XT mikrotekstuuriga anatoomilised silikoonist geeliga täidetud rinnaimplantaadid	109
HU	LUNA XT mikroteksturált anatómiai, szilikongélel töltött mellimplantátumok	116
PL	Implanty anatomiczne piersi z mikroteksturą wypełnioną żelem silikonowym LUNA XT	123
SL	LUNA XT mikro teksturirani anatomici prsni vsadki, napolnjeni s silikonskim gelom	130
SV	LUNA XT mikrostrukturerade anatomiska silikongefyllda bröstimplantat	137L
	LUNA XT Anatomiski implantati za grudi s mikroteksturom punjeni silikonskim gelom	144
BG	Микроструктурни анатомични гръдни импланти със силиконов гел LUNA XT	151
HR	LUNA XT Anatomiski implantati za grudi s mikroteksturom punjeni silikonskim gelom	158
FI	LUNA XT -mikroteksturoidut anatomiset silikongeilitäytteiset rintaimplantit	165
NO	LUNA XT Anatomiske gelylte brystimplantater av silikon, med mikrotekstur	172
IS	LUNA XT dropalaga silikonhlaupfylltir brjóstapúðar með mikró-áferð	179

LUNA XT Micro textured Anatomical Silicone Gel-Filled Breast Implants

DESCRIPTION

The Micro textured Anatomical Breast Implant is a sterile, single-use medical device that is surgically implanted to replace tissue lost by trauma, surgical intervention or to augment tissue for therapeutic purposes.

The Micro textured Anatomical Breast Implant is intended for breast reconstruction surgery.

The anatomic silicone elastomer shell is produced with the outer surface micro textured to offer a disruptive contact surface for collagen interface and filled with a cohesive gel.

Each gel-filled mammary implant is supplied in a cardboard box and packaged within a sealed, double primary pack. The packaging also contains removable labels indicating the characteristics of the implant (Lot number, Reference number), implant card, patient and user leaflets.

INDICATIONS

The Micro textured Anatomical Breast Implant is indicated for

- Surgical Augmentation and contour correction of congenital anomalies of the breast,
- Surgical Reconstruction of the breast following subcutaneous mastectomy and other suitable mastectomy procedures or trauma,
- Surgical correction of Combined breast and chest wall abnormalities,
- Replacement of breast implants for medical reasons.

CONTRA-INDICATIONS

The use of these implants is contraindicated in patients who have one or more of the following conditions:

- FIBROCYSTIC disease,
- INSUFFICIENT TISSUE covering due to radiation damage on the chest wall, tight thoracic skin grafts or radical resection of the pectoralis major muscle,
- EXISTING LOCAL OR METASTATIC CARCINOMA of the breast,
- Suppressed/compromised immune system,
- A history of sensitivity to foreign materials or atopy,
- PHYSIOLOGICALLY/PSYCHOLOGICALLY unsuitable patient,
- ACTIVE INFECTION anywhere in the body,
- Recent history of breast abscess.
- A history of compromised wound healing or regional swelling,
- Any other serious medical condition.

INTENDED USER

The Micro textured Anatomical Breast Implant is intended to be used by suitably qualified surgeons.

IMPORTANT: The surgeon is responsible for conducting the medical assessments regarding the patient's eligibility for implantation and for deciding on the most suitable surgical technique for both the patient and the type of implant chosen.

TARGET POPULATION

The Micro textured Anatomical Breast Implant is intended to be used in adult populations (excluding pregnant and nursing women), with no known contraindications.

INTENDED CLINICAL BENEFITS

The clinical benefits of the Micro textured Anatomical Breast Implant are as follows:

- Replace breast tissue that has been removed following mastectomy
- Replace breast tissue that has been removed due to trauma
- Enhance breast tissue affected by tuberous or congenital breast abnormality
- 10-year Kaplan Meier complication estimate 16.5% Capsular Contracture, 3.8% Rupture and 17.6% Re-intervention rates
- Improve patient's quality of life (based on self-esteem, body image, sexual satisfaction)
- Improve patient's psychology in terms of levels of depression and anxiety

PERFORMANCE CHARACTERISTICS

The Micro textured Anatomical Breast Implants are designed and manufactured according to the current EN ISO 14607 standard concerning mammary implants. The performance characteristics of The Micro textured Anatomical Breast Implants in terms of:

- Shell integrity
- Implant resistance
- Compatibility between the filling gel and the shell
- Diffusion from the implant
- Volume
- Surface texture
- Biocompatibility

Had been validated and comply with the requirements of EN ISO 14607.

MATERIALS

Shell and patch: Silicone elastomer, Filling: Silicone gel

LINK TO THE SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The summary of Safety and effectiveness will be available on the Eudamed Website when it is live <https://ec.europa.eu/tools/eudamed/#/screen/home> and the company Website (<https://www.gcaesthetics.com/>)

PATIENT INFORMATION

Any surgical procedure can have complications and risks. Mammary implant surgery is known to provide psychological satisfaction to patients, but like any surgical procedure it can have potential complications and risks. Mammary implantation is an elective procedure, and the patient should be well counselled on the risk/benefit relationship, by the surgeon. Each of the possible complications and warnings should be discussed with the patient prior to the decision to proceed with surgery.

Each patient should receive the patient information booklet during their initial consultation, to allow time for the patient to read and understand the information on risks, follow-up recommendations and benefits associated with breast implants in order to make an informed decision about proceeding with the surgery. The patient information booklet also contains a 'Patient Informed Consent Form' which allows the patient to retain copies of their implant details such as lot number for their records. This document is available in PDF format from our website www.gcaesthetics.com/we-care/.

LIFETIME

The patient must be informed of the fact that implants have a limited lifetime and that inserting an implant may require removal or replacement, which could entail another surgical procedure. ***It has been demonstrated in published scientific literature that breast implants may last between 10-20 years. Therefore, Eurosilicone's mammary implants, when implanted and used under normal conditions, are expected to have a product lifetime similar to that stated in the aforementioned published literature.***

Clinical monitoring of Eurosilicone implants has enabled, to date, to establish implant durability of 82.4% after 10 years (risk of a first reoperation = 17.6%, using the Kaplan Meier method).

MRI COMPATIBILITY

Although Eurosilicone implants have not been specifically tested for use in an MRI, please note that Eurosilicone's Gel-filled Mammary Implants are all manufactured from medical implant grade silicone materials which are compatible with MRI Scan. The implantable silicone grade material is the same as that used in other Eurosilicone silicone implants where patients have undergone MRI Scans and no compatibility issues have been reported to date.

STEROID USE

The patient should be informed to consult a physician before using steroid drugs in the implant area to avoid extrusion of the implant.

INTERFERENCE

During a mammography, the radiologist must be informed that there is an implant in the breast in order to adapt the mammography pressure accordingly. The surgeon or doctor must be notified about the implant in the event of any surgical intervention in the breast area.

The presence of breast implants could delay the detection of breast cancer by self-exam. For this reason, patients must be informed that they need to consult a surgeon and/or doctor for appropriate medical monitoring, as well as regular breast cancer screening. Patients must also be informed that they need to consult a doctor if they suspect complications particularly in the event of trauma or compression caused, for example, by violently massaging the breasts, a sporting activity or wearing seat belts.

IMPLANT CARD

After surgery, the surgeon should provide to the patient:

the implant card which should be filled out according to the leaflet Instructions for Healthcare Professionals filling out the EUROSILICONE International Implant Card.

the patient leaflet which is available on the Information website for patients referenced on the card.

The patient should be informed that the implant card must be kept on their person at all the times and enables:

the patient to identify the implanted devices and to get access to other information related to the implanted device (e.g. via EUDAMED, and other websites).

the patients to identify themselves as persons requiring special care in relevant situations e.g. security checks.

emergency clinical staff or first responder to be informed about special care/needs for relevant patients in case of emergency situations.

WARNINGS

ALTERATION OF IMPLANTS

No alteration should be made to the mammary implant before implantation. Alteration of implants voids all warranties, express or implied.

DAMAGE OF IMPLANTS - BREAKAGE AT SURGERY

Extreme care should be taken in the use and handling of implants to minimize the potential for breakage of shells. All implants have been produced by established manufacturing techniques and under strict quality control standards, but there may be some breakage of implants during handling or in surgery, both in initial and any subsequent surgeries.

Extreme care must be taken to avoid unintentional damage to the implant during the implantation or explantation surgery.

Do not contact the implant with sharp surgical instruments or devices such as scalpels, forceps, hemostats, suture needles and hypodermic needles.

Do not contact the implant with blunt surgical instruments such as clamps, retractors and dissectors.

Do not contact the implant with cautery devices.

Do not use excessive manipulation, force or stress.

Implants should be carefully inspected for structural integrity prior to use. Damaged products should not be implanted, do not attempt to repair damaged products. An extra product should be available at the time of surgery in the event of damage to, or contamination of, the implant. Recommended procedures for testing, examination and handling of products should be meticulously followed to assure proper use of implants. Patients should be instructed to inform other treating physicians of the presence of implants to minimize the risk of damage.

INTERFERENCE WITH MAMMOGRAPHY

Standard positioning techniques have shown significant limitations when used for imaging augmented breasts. The implant can interfere with the detection of early breast cancer through mammography by obscuring some underlying breast tissue and/or by compressing overlying tissue which can 'hide' suspicious lesions in the breast. Patients should be instructed to request radiologists who are experienced with the most current radiological techniques and equipment for imaging breasts with implants, and to inform their radiologists of the presence, type, and placement of implants and if anatomical implants, they should be aware of orientation marks on the anatomical implant.

SINGLE USE

Implants are intended for SINGLE USE ONLY. DO NOT REUSE EXPLANTED PRODUCTS. DO NOT RESTERILIZE ANY PRODUCT.

Explanted products should not be reused because recleaning and re-sterilization procedures may not adequately remove biological residues, such as blood, tissue and other matter, which could retain resistant pathogens and may also affect the performance of the implant.

USE OF DRUGS

Eurosilicone can neither predict, nor warrant, the safety of intraluminal introduction of any drugs, including, but not limited to, aesthetic, steroid, antibiotic and vitamin solutions. If such use is contemplated, the appropriate drug manufacturer should be consulted.

COMPLICATIONS

Possible complications known to be associated with gel-filled mammary implants are as follows:

Other potential complications associated with all surgical procedures should be discussed with the patient. These include but are not limited to: infection (see below); hematoma (see below); serous fluid accumulation (see below); loss of sensation (see below); reaction to medication; nerve damage; patient intolerance to any foreign implant; and poor wound healing.

ACCIDENTAL RUPTURE OF IMPLANTS IN SURGICAL POCKET

If this occurs, the gel mass may be penetrated with index finger of double gloved hand, whilst exerting pressure on the breast with other hand. The gel can be manipulated out of the cavity, within the hand, and the outer glove can then be pulled over the mass for disposal. Blot the surgical pocket with gauze sponges. Isopropyl alcohol will aid the removal of gel from instrumentation.

ASYMMETRY

Asymmetry may result from improper initial placement, displacement, or failure to correct pre-existing asymmetry through variation in individual implant size. Pre-existing asymmetry may not be entirely correctable. Asymmetry may also be a symptom of capsular contracture, fluid accumulation, infection, post-operative breast dysplasia, unilateral discrepancy in muscle development, or deflation, requiring further investigation.

BREAST IMPLANT ASSOCIATED ANAPLASTIC LARGE CELL LYMPHOMA (BIA-ALCL)

European safety information, the US FDA and current scientific literature have identified an association between breast implants and the development of a cancer of the immune system called *breast implant-associated anaplastic large cell lymphoma (BIA-ALCL)*. *This specific entity is included in the WHO 2016 classification under the terminology "ALCL – ALM". This cancer occurs more commonly in patients with textured breast implants rather than smooth implants, and typically develops many years after breast implant surgery. At this time rates of BIA-ALCL are considered to be low.*

In most patients, BIA-ALCL is found within the fluid or scar tissue next to the implant. The main symptoms of BIA-ALCL are persistent swelling aforementioned symptoms, or if BIA-ALCL is suspected, it is advised that patients should be assessed to rule out peri-implant ALCL. When testing for BIA-ALCL, collect fresh seroma fluid and representative portions of the capsule and send for pathology tests to rule out ALCL. Diagnostic evaluation should include cytological evaluation of seroma fluid with Wright Giemsa stained smears and cell block immunohistochemistry testing for cluster of differentiation (CD) and Anaplastic Lymphoma Kinase (ALK) markers. Most cases of breast implant-associated ALCL are treated by removal of the implant and the capsule surrounding the implant and some cases have been treated by chemotherapy and radiation. Although treatment is usually successful BIA-ALCL can be fatal if left untreated, hence, early detection and treatment is key.

In the case that your patient is diagnosed with BIA-ALCL, a report should be made to Eurosilicone via the e-mail address: return@gcaesthetics.com and an individual treatment program should be initiated by a multidisciplinary team.

BREAST IMPLANT ILLNESS

A small proportion of women who have breast implants, for both aesthetic and reconstructive purposes, self-identify as having a number of symptoms they believe arise from the presence of their breast implants. Whilst not a medical diagnosis, they refer to their symptoms as 'Breast Implant Illness' (BII).

These multiple generic symptoms are varied, and all cases have been self-diagnosed and self-reported. They include, but are not limited to, flu-like symptoms such as extreme tiredness, brain fog, joint aches, immune-related symptoms, sleep disturbance, depression, hormonal issues, headaches, hair loss and chills.

There are a variety of other reasons these symptoms might be experienced, including background illnesses or hormonal changes. In addition, there are a number of scientific studies investigating similar symptoms experienced by women in the general population with and without breast implants. On average, around 50% of women who self-identify as having 'BII' feel that their symptoms improve after removal – sometimes temporarily and sometimes permanently. It therefore appears that removing breast implants does not necessarily improve symptoms in everyone. To date, there is no research demonstrating which symptoms may or may not improve with implant removal.

Currently, there are no tests that can confirm 'BII'. Research continues in this area, specifically around patients who have autoimmune disease or a predisposition to autoimmune disease. Breast implantation is contraindicated in women with a suppressed or compromised immune system – see CONTRAINDICATIONS section above. Medical advice should be sought if patients suspect they have BII. Their symptoms may not be related to the implants and other medical investigations should not be overlooked or ignored in determining the cause of their symptoms.

BREAST FEEDING

The implantation may affect the milk production during nursing, particularly, peri-areolar incision may reduce considerably the possibility of breast feeding.

CALCIUM DEPOSITS

Calcification commonly occurs in mature breast tissue with or without implantation. Calcification is also known to occur after implantation of a foreign body, although the etiology is unknown, and reported cases are rare. Microcalcification after implantation typically occurs on/or around the fibrous capsule in thin plaques or accumulations. Extensive microcalcifications may cause breast hardness and discomfort and may necessitate surgical intervention.

CAPSULAR CONTRACTURE

Fibrous capsular contracture is a common complication following breast implant surgery. While formation of a fibrous tissue capsule surrounding the implant is a normal physiological response to a foreign body, not all capsules contract. Contracture of the fibrous capsular tissue surrounding the implant may result in firmness, discomfort or pain in the breast, distortion of the breast, palpability of the implant, or displacement of the implant. The etiology of capsular contracture is unknown but is most likely multifactorial. Contracture develops to varying degrees, unilaterally or bilaterally, and may occur within weeks to years after surgery. Severe cases may require surgical intervention. Capsular contracture may recur after capsulotomy or capsulectomy. Lumps perceived as capsular scarring have delayed the diagnosis of palpable tumor. Questionable lumps must be promptly investigated.

DELAYED WOUND HEALING

Some patients experience delayed healing, and for others the incision site may not heal well. It may open from injury or infection. If the implant is exposed, further surgery will be required. Tissue necrosis is the development of dead tissue around the implant. It will delay wound healing, may cause wound infection and may require surgical correction and/or implant removal. Tissue necrosis has been reported following the use of steroid drugs, chemotherapy, radiation to breast tissue, and smoking, but in some cases, it may occur without any known cause.

DISSATISFACTION WITH RESULTS

The complications of incorrect size, misplaced scar location, hypertrophic scarring or unfavorable scarring and those listed below are usually related to surgical technique. Careful size selection, creation of an appropriate and adequate size surgical pocket, and use of current accepted surgical procedures are the surgeon's responsibility.

DISPLACEMENT OF IMPLANT AND MALROTATION

Implants may displace or rotate with accompanying discomfort and/or distortion in breast shape. Difficult placement techniques may increase the risk of displacement by reducing pocket size and placement accuracy. Displacement may require surgical intervention.

EXPLANATION OF MAMMARY IMPLANTS

Although there are no definite reports in the medical literature, it is anticipated that a large number of mammary implants will be surgically explanted or replaced. This is not unlike many other plastic surgery procedures which are commonly repeated to maintain patient satisfaction. Due to the wide variety of patients' physical responses to breast implant surgery, differences in surgical techniques and medical treatments, as well as potential complications, patients should be advised that these should not be considered lifetime implants and explant surgery may be indicated at any time. Eurosilicone makes no representations for the term of implantation of the device.

EXTERNAL CAPSULOTOMY

Treatment of capsular contracture by external manual compression may cause the shell to weaken or rupture. There have been reports of rupture and Eurosilicone recommends against the external capsulotomy procedure and is not responsible for the integrity of the implant should manual compression or other external stress techniques be used.

EXTRUSION

Unstable or compromised tissue covering and/or interruption of wound healing may result in exposure and extrusion of the implant. Causes or contributing factors may include infection, wound dehiscence, necrosis with or without infection, capsular contracture, closed capsulotomy, unsuitable skin flap, improper size and placement of implant, and/or tissue erosion associated with implant folds. The incidence of extrusion has been shown to increase when the prosthesis has been placed in injured areas: scarred, heavily irradiated or burned tissue or crushed bone areas; where severe surgical reduction of the area has been performed; and where steroids are used in the surgical pocket.

GRANULOMA

Granuloma formation is a common tissue response to the presence of foreign materials. It is possible for a tissue reaction to be caused by the presence of the implant leading to formation of silicone granuloma.

HEMATOMA/ EARLY ONSET SEROMA

Hematoma and serous fluid accumulation are complications associated with any type of invasive surgery. Postoperative hematoma and seroma may contribute to infection and/or capsular contracture. Postoperative hematoma and seroma may be minimized by meticulous attention to hemostasis during surgery, and also possibly by postoperative use of a closed drainage system. Persistent, excessive bleeding must be controlled before the device is implanted. Any postoperative evacuation of hematoma or seroma must be conducted with care to avoid contamination or damage to the implant.

INFECTION

Pre-existing infection not resolved before implant placement increases the risk of periprosthetic infection. Do not expose the implant or filling accessories to contaminants, which increases the risk of infection.

Infection is an inherent risk following any type of invasive surgery. Infection around a breast implant may occur within days, weeks, or even years, after surgery. Signs of acute infection reported in association with breast implants include erythema, tenderness, fluid accumulation, pain and fever. Signs of subclinical infection may be difficult to detect. Postoperative infections should be treated aggressively according to standard medical practices to avoid more serious complications. Infection that is unresponsive to treatment or necrotizing infection may require implant removal. Capsular contracture may be related to infection in the area surrounding the implant.

PAIN

Any unexplainable pain, not associated with any of the complications listed herein, must be promptly investigated.

PHLEBITIS SUPERFICIAL

In rare cases, damage to the vein walls may occur during implantation of the mammary prosthesis. Subsequently, this may result in phlebitis superficial, which is an inflammation of the veins caused by a blood clot below the surface of the skin. This complication is rarely serious and with proper care, can be resolved rapidly.

PTOSIS

Ptosis, also known as sagging of the breast, is a deformity which can occur pre- and post- breast surgery. Pre-breast surgery, due to the breast having no bony or cartilaginous support, time and gravity can affect the natural breast including glandular and skin stretching. Post- breast surgery, the size, weight and position of breast implants can have an effect on the degree of ptosis. If significant, this complication can often lead to re-intervention.

RUPTURE OF GEL IMPLANTS

Gel rupture can occur post-operatively from damage to the implant during handling or surgery. Rupture of the shell can also occur from contracture, trauma or excessive manipulation. Ruptures of unknown etiology have also been reported. Because of the cohesiveness of the gel material, some gel ruptures can go undetected unless surgery occurs for another reason (e.g. size exchange). Despite the cohesive properties of the gel, extravasation out of the surgical pocket can occur under pressure, in which case additional surgery may be required to retrieve the gel. Inflammation and formation of silicone granulomas have been reported. If shell rupture is suspected, the implant should be removed.

SENSATION

The risk of temporary or permanent dysesthesia exists following any invasive surgical procedure. Careful surgical technique can minimize, but not preclude, this risk. Dysesthesia of the nipple/areolar complex, and less frequently, the breast area in general, has been reported following implantation and may be temporary or permanent. The risk of neurological impairment increases with more extensive surgery. Surgical implantation of a mammary prosthesis may interfere with the ability to breast feed. However, it should be noted that previous breast reconstruction surgery, may be the initial cause of this interference.

SILICONE ALLERGY

Although silicone allergy is extremely rare, this risk cannot be completely ruled out.

UPPER LIMB LYMPHEDEMA

Lymphedema is a chronic condition characterized by swelling of the arm, hand, or breast, which can result from breast cancer treatment, due to the accumulation of lymphatic fluid in the interstitial tissues. This condition is known to significantly compromise quality of life due to body image changes, alterations in arm function, and increase the risk of other complications including infection and cellulitis.

Recent studies suggest immediate breast reconstruction may be associated with a reduced risk of post-mastectomy lymphedema. Treatment options include conservative management or surgical treatment using a super-microsurgical technique. However, no complete cure for lymphedema is known.

WRINKLING AND FOLDS

Thin or inadequate overlying tissue, patients with little or no subcutaneous fat, implants which are too large for the surgical pocket, or the anatomical structure of the patient, and subcutaneous placement may contribute to palpable, or visible, implant edge, implant rippling, creasing, wrinkles or folds. Folds may result in thinning and erosion of adjacent tissue and erosion of the implant. Signs of skin inflammation, such as tenderness and erythema, may indicate thinning or erosion and must be promptly investigated. Palpable wrinkling / and or folds may be confused with palpable tumor, and questionable cases must be promptly investigated.

INSTRUCTIONS FOR USE

STORAGE CONDITIONS

Eurosilicone products must be stored in normal conditions, must be protected from humidity and direct sunlight and must be store with the arrows pointing up. In these conditions, Eurosilicone products have a shelf life of five years.

SUPPLIED STERILE

Implants are supplied in sterile form (Ethylene Oxide Sterilization), processed by validated strictly controlled sterilization cycles with ethylene oxide. Sterility is verified in accordance with standards. Sterility of the implant is maintained only if the package is intact and undamaged. If the sterile packaging has been damaged or unintentionally opened before use DO NOT USE.

SINGLE USE

Implants are intended for SINGLE USE ONLY. DO NOT REUSE EXPLANTED PRODUCTS. DO NOT RESTERILIZE ANY PRODUCT. Explanted products should not be reused because recleaning and re-sterilization procedures may not adequately remove biological residues, such as blood, tissue and other matter, which could retain resistant pathogens and may also affect the performance of the implant.

PACKAGING

Sterile product is supplied in a sealed, double primary package. Sterility is not guaranteed if the package has been damaged or opened. Tear-Off patient record labels are attached to the primary package. These labels are to be attached to the patients' records.

TO OPEN PACKAGED STERILE PRODUCT

Peel open outer package under clean, aseptic conditions, over sterile field, allowing sealed inner package to fall gently into the field. Attach patient record portion of the inner label to patient's chart. Using aseptic precautions, peel open the inner package.

PRODUCT EXAMINATION AND HANDLING

Product should be visually examined for any evidence of particulate contamination, damage or leakage.

The implant should be kept submerged in sterile water or normal saline prior to implantation to prevent contact with airborne and surgical field particulate contaminants.

Do not immerse the implant in Betadine or Iodine containing solutions. If Betadine or Iodine containing solutions are used in the pocket, ensure that it is rinsed thoroughly so no residual solution remains in the pocket.

Gel-filled mammary implants may contain air bubbles. This occurs normally during sterilization and has no effect on product integrity or performance.

Use disposable, non-shedding wipes only to ensure there are no particles left on the surface of the implant.

CONTAMINATION OF IMPLANTS

Care must be taken to prevent surface contaminants such as talc, dust and skin oils from coming into contact with the implant. Products should be inspected for contamination prior to insertion. Contamination at the time of surgery increases the risk of periprosthetic infection, and possibly capsular contracture.

Textured silicone surfaces are potentially more susceptible to contamination than smooth silicone surfaces. Extra care should be taken in the handling of textured surface implants. Eurosilicone accepts no liability for products contaminated by other substances after the product leaves our possession.

SURFACE CONTAMINANTS (TALCS, DUST, LINT, OILS) ON SURFACE OF IMPLANTS CAN CAUSE FOREIGN BODY REACTION, HANDLE WITH CARE WITH SURGICAL GLOVES (RINSED FREE OF TALC) WITH STRICT ASEPTIC TECHNIQUE. DO NOT IMPLANT CONTAMINATED PRODUCT.

Back-up implants must be readily available at the time of surgery for use in the event of contamination.

SURGICAL PROCEDURE

Correct surgical procedures and techniques are the responsibility of the medical profession. Each surgeon must evaluate the suitability of the procedure based upon current accepted techniques, individual judgement and experience. Proper size and shape of implants must be determined for the individual patient by the surgeon. An incision should be of appropriate length to accommodate the style, size, and profile of the implant. This will reduce the potential for creating excessive stress to the implant during insertion.

SPECIFIC GUIDELINES FOR THE USE OF EUROSILICONE MAMMARY IMPLANTS

1. Size selection should be decided by base dimensions, the patient's dimensions, and individual breast characteristics rather than preferred volume or projection.
2. The preformed shapes of the Eurosilicone implants are firmer than other cohesive gel filled implants, and the surgical incisions should be longer than those required for compressible implants.
3. Where appropriate the new inframammary fold may be required to be approximately 1cm below the natural fold.
4. The recommended approach is via the inframammary fold. It is the responsibility of the surgeon to decide if axillary or trans areolar placements may be possible. We caution that insertion and orientation may be more difficult due to the increased firmness of the gel.
5. Sub-muscular placement (with possible lower ligament release) is preferred where possible. As the gel is highly cohesive it has a more defined radius on the upper apex margin. Sub-muscular placement reduces the probability of this apex radius being palpable.

6. DO NOT use force during placement of Eurosilicone implants. Over manipulation of Eurosilicone implants may result in permanent distortion of the cohesive gel and loss of the implant's shape.
7. Maintenance of correct alignment and avoidance of rotation in the post-operative phase is aided by:
 - a) Achieving a surgical pocket closely aligned to the implant.
 - b) Achieving a "dry" surgical pocket.
 - c) Practicing the procedure of implant stabilization by use of post-operative external immobilization.
8. All mammary implants have the potential in the immediate post-operative phase to rotate and the above guidelines will assist in avoidance of this with any anatomically shaped implants.
9. CORRECT ORIENTATION IS VITAL when using any anatomical design of Eurosilicone implant, due to the conical shape of these implants having an apex. Correct orientation during implantation is aided by orientation markers as below.

MIA HR1-3, MIA MR1-3, MIA LR1-3

Palpable orientation markers are placed on the superior and posterior surfaces of the implant. Please refer to the schematic diagram below to see the position of the markers. Schematic diagram showing orientation markers.

MATERIALS DEVICE DESCRIPTION

These medical devices are manufactured from medical grade silicone materials which are suitable for long term implantation. Master Files have been filed with the U.S FDA.

The following three tables provide quantitative and qualitative information on the materials and substances which patients can be exposed based on chemical characterization of representative devices. The materials and device have been subjected to biocompatibility testing and evaluation, and risk assessments to demonstrate their biological safety. However, individual responses to chemicals may vary, and all reactions cannot be predicted.

Table 1. Breast Implant Device Materials

Device materials	Implant component
Silicone dispersions	Shell
Silicone adhesive	Shell
Silicone gel	Gel fill
Silicone dispersions	Closure Patch
Silicone Elastomer	Orientation markers
Platinum catalyst	Shell and Gel fill
Organotin catalyst	Shell (silicone adhesive)

Table 2. Chemicals Released by Breast Implants

Compounds	Whole device (µg/g or ppm)	Compounds	Whole device (µg/g or ppm)
Volatiles¹ - Chemicals that can be released by breast implants as a gas			
Benzene	Not Detected	Ethylbenzene	Not Detected
Bromobenzene	Not Detected	Hexachlorobutadiene	Not Detected
Bromodichloromethane	Not Detected	Isopropylbenzene	Not Detected
Bromoform	Not Detected	p-Isopropyltoluene	Not Detected
Butylbenzene (N- / tert- / sec-)	Not Detected	N-propylbenzene	Not Detected
Carbon tetrachloride	Not Detected	Naphthalene	Not Detected
Chlorobenzene	Not Detected	Styrene	Not Detected
Chloroform	Not Detected	1,1,1,2-Tetrachloroethane	Not Detected
2,4-Chlorotoluene	Not Detected	Tetrachloroethene	Not Detected
Dibromochloromethane	Not Detected	Trichloroethene	Not Detected
1,2-Dibromo-3-chloropropane	Not Detected	Trichloroethane (1,1,1- / 1,1,2-)	Not Detected
Dichlorobenzene (1,2- / 1,3- / 1,4-)	Not Detected	1,2,3-Trichloropropane	Not Detected
Dichloroethane (1,1- / 1,2-)	Not Detected	Trimethylbenzenes (1,2,4- / 1,2,5-)	Not Detected
Dichloroethene (1,1- / cis-1,2- / trans-1,2-)	Not Detected	Toluene	Not Detected
Dichloropropane (1,2- / 1,3-)	Not Detected	Trichlorobenzenes (1,2,4- / 1,2,3-)	Not Detected
1,1-Dichloropropene	Not Detected	o-Xylene	Not Detected
Ethanol, 2-trimethylsilyl	2.91 max.	m- & p-Xylene	Not Detected
Total volatiles	2.91		
Extractables² - Chemicals that can be released by breast implants following soaking in water			
Cyclic siloxanes (D4, D5, D6...)	Not Detected	Linear siloxanes (L3, L4, L5...)	Not Detected
Total extractables	< 24.2		
¹ Not Detected means that the level of the individual volatile was below the quantitation limit of the test method. The quantitation limit was 0.86 µg/g for m- & p-Xylene and 2,4-Chlorotoluene, and 0.44 µg/g for all other compounds.			
² Not Detected means that the level of the individual extractable was below 24.2 µg/g, the quantitation limit of the test method.			

Table 3. Heavy Metals Found in Breast Implants

Heavy metals	Concentration (ppm)	Heavy metals	Concentration (ppm)
Arsenic	Not Detected	Chromium	Not Detected
Cadmium	Not Detected	Copper	Not Detected
Lead	Not Detected	Tin	Not Detected
Mercury	10 max.	Nickel	Not Detected
Vanadium	Not Detected	Platinum	Not Detected
Molybdenum	Not Detected	Zinc	0.25 max.
Selenium	Not Detected	Magnesium	0.41 max.
Cobalt	Not Detected	Manganese	Not Detected
Antimony	Not Detected	Thallium	Not Detected
Barium	Not Detected		

Not Detected means that the level of the individual element was below 0.25 ppm, the quantitation limit of the test method.

DISPOSAL OF DEVICE

All mammary implants must be eliminated according to requirements related to medical wastes involving an infectious risk.

REPORTING OF SERIOUS INCIDENT:

Any serious incident (incident that directly or indirectly led, might have led or might lead to any of the following: (a) the death of a patient, user or other person, (b) the temporary or permanent serious deterioration of a patient's, user's or other person's state of health, (c) a serious public health threat) should be reported to Eurosilicone via the e-mail address: retum@gcaesthetics.com and to the competent authority of the Member State in which the user and/or patient is established.

RETURN OF GOODS POLICY

Prior to the return of any product:

- Authorization must be requested in writing.
- Decontamination certificate must be supplied with any product returned which has been removed from the packaging.
- If it is necessary to puncture the shell of any sealed shell product to assist with safe decontamination sterilization procedures, the area of mechanical interference should be marked with indelible marker on the surface of the product and reference made on the decontamination certificate.
- Products File Notes and returned products must be sent to Eurosilicone S.A.S. by your distributor only.

WARRANTY

The manufacturer warrants that reasonable care was used in the manufacture of these products and will replace any product that the manufacturer's investigation shows to be defective at time of shipment. Patient selection, surgical procedures, post-surgical treatment and stresses, and handling of the devices are totally the responsibility of the customer. Eurosilicone has no control over the conditions of use and cannot warrant good effect or against ill effect following the use of the device and shall not be responsible for any incidental or consequential loss or damage or expenses directly or indirectly arising from the use of this device. All other warranties (whether implied by law or otherwise) are excluded to the fullest extent permitted by law.

CAUTION:

Federal (USA) law restricts this device to sale by or on the order of a physician.

BIBLIOGRAPHIC REFERENCES

Literature references are available upon request from Eurosilicone.

LABELLING SYMBOL:

Indicates the products dimensions

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- Einungis dreifingaraðili má senda skráðar upplýsingar um vöruna og þær vörur sem skila á til Eurosilicone S.A.S.
Eurosilicone SAS