SILICONE GEL-FILLED SMOOTH MAMMARY IMPLANTS

EN	FR	DE	IT	PT	ES	DA
NL	FI	EL	SV	CS	ET	HU
LV	LT	PL	SK	SL	NO	TR
BG						

EN

DESCRIPTION

Gel-filled mammary implants are designed for use in cosmetic augmentation mammoplasty or reconstruction of the breast. The round silicone elastomer shell is produced smooth and filled with a high cohesive gel. Each gel-filled mammary implant is supplied in a cardboard box and packaged within a sealed, double primary pack. The packaging also contains removable labels indicating the characteristics of the implant (Lot number, Reference number).

INDICATIONS

These gel-filled mammary implants are indicated for the following:

- Cosmetic Augmentation Surgery.
- Augmentation and contour correction of congenital anomalies of the breast.
- Reconstruction of the breast following subcutaneous mastectomy and other suitable mastectomy procedures or trauma.
- Combined breast and chest wall abnormalities.
- Replacement of implants for medical and cosmetic reasons.

CONTRA-INDICATIONS

The use of these implants is contraindicated in patients who have one or more of the following conditions:

- FIBROCYSTIC disease.
- INSUFFICIENT TISSUE covering due to radiation damage on the chest wall, tight thoracic skin grafts or radical resection of the pectoralis major muscle.
- EXISTING LOCAL OR METASTATIC CARCINOMA of the breast.
- Suppressed/compromised immune system.
- A history of sensitivity to foreign materials or or atopy.
- Previous UNSATISFACTORY HISTORY of augmentation/reconstruction.
- PHYSIOLOGICALLY/PSYCHOLOGICALLY unsuitable patient.
- ACTIVE INFECTION anywhere in the body.
- Recent history of breast abscess.
- A history of compromised wound healing.
- Any other serious medical condition which, in the opinion of the surgeon, would compromise the health of the patient.
- Since there is no data concerning the concomitant use of these devices with dermal fillers, such use of dermal fillers in the proximity of the device (e.g., in the décolletage skin area) is contraindicated. Physicians should also advise patients against future use of dermal fillers in proximity to these devices.
- These devices are not intended for the treatment of children or pregnant or nursing women.

N.B. It is the SURGEONS RESPONSIBILITY to make a medical assessment as to patient suitability for implantation, and to decide the surgical technique suitable to both the patient and the chosen implant type and design.

PATIENT INFORMATION

Any surgical procedure can have complications and risks. Mammary implant surgery is known to provide psychological satisfaction to patients, but like any surgical procedure it can have potential complications and risks. Mammary implantation is an elective procedure and the patient should be well counselled on the risk/benefit relationship, by the surgeon. Each of the possible complications and warnings should be discussed with the patient prior to the decision to proceed with surgery.

Each patient should receive the patient information booklet during their initial consultation, to allow time for the patient to read and understand the information on risks, follow-up recommendations and benefits associated with breast implants in order to make an informed decision about proceeding with the surgery. The patient information booklet also contains a 'Patient Informed Consent Form' which allows the patient to retain copies of their implant details such as lot number for their records. This document is available in PDF format from our website www.gcaesthetics.com/we-care/.

LIFETIME

Patients should be informed that mammary implants should not be considered as lifetime devices due to a variety of factors including potential differences in physiological reactions, implantation procedures, inherent nature of silicone implant designs and external mechanical influences. It is anticipated that the implant may have to be removed or replaced, requiring revision surgery during the patient's lifetime. However, it has been demonstrated in published literature that smooth breast implants may last between 10-20 years. Therefore these smooth mammary implants, when implanted as instructed and used under normal conditions, would be expected to have a product lifetime similar to that stated in published literature (which is estimated to be 10-20 years).

Clinical monitoring of GC Aesthetics' textured gel-filled mammary implants has enabled, to date, to establish implant durability of 82.4% after 10 years (re-operations took place for 17.6% of implants originally placed). In addition, GC Aesthetics will perform a long-term clinical study on its smooth mammary implants to gather short and long term safety and performance data along with survivorship data over 10 years.

MRI COMPATIBILITY

N.B. Although these implants have not been specifically tested for use in an MRI. They are all manufactured from medical implant grade silicone materials which are compatible with MRI Scan. The implantable silicone grade material is the same as that used in other silicone implants where patients have undergone MRI Scans and no compatibility issues have been reported to date.

STEROID USE

The patient should be informed to consult a physician before using steroid drugs in the implant area to avoid extrusion of the implant.

WARNINGS

ALTERATION OF IMPLANTS

No alteration should be made to the mammary implant before implantation. Alteration of implants voids all warranties, express or implied.

DAMAGE OF IMPLANTS - BREAKAGE AT SURGERY

Extreme care should be taken in the use and handling of implants to minimise the potential for breakage of shells. All implants have been produced by established manufacturing techniques and under strict quality control standards, but there may be some breakage of implants during handling or in surgery, both in initial and any subsequent surgeries.

Extreme care must be taken to avoid unintentional damage to the implant during the implantation or explantation surgery.

- Do not contact the implant with sharp surgical instruments or devices such as scalpels, forceps, hemostats, suture needles and hypodermic needles.
- Do not contact the implant with blunt surgical instruments such as clamps, retractors and dissectors.
- Do not contact the implant with cautery devices.
- Do not use excessive manipulation, force or stress.

Implants should be carefully inspected for structural integrity prior to use. Damaged products should not be implanted, do not attempt to repair damaged products. An extra product should be available at the time of surgery in the event of damage to, or contamination of, the implant. Recommended procedures for testing, examination and handling of products should be meticulously followed to assure proper use of implants. Patients should be instructed to inform other treating physicians of the presence of implants to minimise the risk of damage.

INTERFERENCE WITH MAMMOGRAPHY

Standard positioning techniques have shown significant limitations when used for imaging augmented breasts. The implant can interfere with the detection of early breast cancer through mammography by obscuring some underlying breast tissue and/or by compressing overlying tissue which can 'hide' suspicious lesions in the breast. Patients should be instructed to request radiologists who are experienced with the most current radiological techniques and equipment for imaging breasts with implants, and to inform their radiologists of the presence, type, and placement of implants and if anatomical implants, they should be aware of orientation marks on the anatomical implant.

SINGLE USE

Implants are intended for SINGLE USE ONLY. DO NOT REUSE EXPLANTED PRODUCTS. DO NOT RESTERILIZE ANY PRODUCT. Explanted products should not be reused because recleaning and resterilization procedures may not adequately remove biological residues, such as blood, tissue and other matter, which could retain resistant pathogens.

USE OF DRUGS

The manufacturer can neither predict, nor warrant, the safety of intraluminal introduction of any drugs, including, but not limited to, anaesthetic, steroid, antibiotic and vitamin solutions. If such use is contemplated, the appropriate drug manufacturer should be consulted.

COMPLICATIONS

Possible complications known to be associated with gel-filled mammary implants are as follows:

Other potential complications associated with all surgical procedures should be discussed with the patient. These include but are not limited to: infection (see below); haematoma (see below); serous fluid accumulation (see below); loss of sensation and nerve damage (see below); reaction to medication; patient intolerance to any foreign implant; and poor wound healing.

ACCIDENTAL RUPTURE OF IMPLANTS IN SURGICAL POCKET

If this occurs, the gel mass may be penetrated with index finger of double gloved hand, whilst exerting pressure on the breast with other hand. The gel can be manipulated out of the cavity, within the hand, and the outer glove can then be pulled over the mass for disposal. Blot the surgical pocket with gauze sponges. Isopropyl alcohol will aid the removal of gel from instrumentation.

ASYMMETRY

Asymmetry may result from improper initial placement, displacement, or failure to correct pre-existing asymmetry through variation in individual implant size. Pre-existing asymmetry may not be entirely correctable. Asymmetry may also be a symptom of capsular contracture, fluid accumulation, infection, postoperative breast dysplasia, unilateral discrepancy in muscle development, or deflation, requiring further investigation.

BREAST IMPLANT-ASSOCIATED ANAPLASTIC LARGE CELL LYMPHOMA (BIA-ALCL)

European safety information, the FDA and current scientific literature have identified an association between breast implants and the rare development of anaplastic large cell lymphoma (ALCL), a type of non-Hodgkin's lymphoma. Specifically, women with textured breast implants may have a low but higher risk of developing ALCL in an area adjacent to the implant. This specific entity is included in the WHO 2016 classification under the terminology "ALCL – AIM". These safety reports advise that an extremely small number of cases have been reported in either the fluid or scar tissue adjacent to the implant.

The main symptoms of BIA-ALCL in women with breast implants reported are, persistent swelling or pain, late onset persistent peri-implant seroma¹, capsular contracture or a mass in the vicinity of the breast implant. If women present with any of the aforementioned symptoms, or if BIA-ALCL is suspected, it is advised that patients should be assessed to rule out peri-implant ALCL. As per the recommended guidelines, when testing for BIA-ALCL, collect fresh seroma fluid and representative portions of the capsule and send for pathology tests to rule out ALCL. Diagnostic evaluation should include cytological evaluation of seroma fluid with Wright Giemsa stained smears and cell block immunohistochemistry testing for cluster differentiation (CD) and Anaplastic Lymphoma Kinase (ALK) markers. Most cases of breast implant-associated ALCL are treated by removal of the implant and the capsule surrounding the implant and some cases have been treated by chemotherapy and radiotherapy.

In the case that your patient is diagnosed with BIA-ALCL, a report should be made to the manufacturer via the e-mail address: ukvqilance@gcaesthetics.com and an individual treatment programme should be initiated by a multidisciplinary team.

'BREAST IMPLANT ILLNESS'

A small proportion of women who have breast implants, for both aesthetic and reconstructive purposes, self-identify as having a number of symptoms they believe arise from the presence of their breast implants. Whilst not a medical diagnosis, they refer to their symptoms as 'Breast Implant Illness' (BII).

These multiple generic symptoms are very varied and all cases have been self-diagnosed and self-reported. They include, but are not limited to, flu-like symptoms such as extreme tiredness, brain fog, joint aches, immune-related symptoms, sleep disturbance, depression, hormonal issues, headaches, hair loss and chills.

There are a variety of other reasons these symptoms might be found including background illnesses or hormonal changes. In addition, there are a number of scientific studies investigating similar symptoms experienced by women in the general population both with and without breast implants. On average, around 50% of women who self-identify as having 'BII' feel that their symptoms improve after implant removal – sometimes temporarily and sometimes permanent. It therefore appears that removing breast implants does not necessarily improve symptoms in everyone. To date, there is no research demonstrating which symptoms may or may not improve with implant removal.

Currently, there are no tests that can confirm 'BII'. Research continues in this area, specifically around patients who have autoimmune disease or a predisposition to autoimmune disease. Breast implantation is contraindicated in women with a suppressed or compromised immune system - see CONTRA-INDICATIONS section above. Medical advice should be sought if patients suspect they have BII. Their symptoms may not be related to the implants and other medical investigations should not be overlooked or ignored in determining the cause of their symptoms.

CALCIUM DEPOSITS

Calcification commonly occurs in mature breast tissue with or without implantation. Calcification is also known to occur after implantation of a foreign body, although the etiology is unknown, and reported cases are rare. Microcalcification after implantation typically occurs on or around the fibrous capsule in thin plaques or accumulations. Extensive Microcalcification may cause breast hardness and discomfort, and may necessitate surgical intervention.

CAPSULAR CONTRACTURE

Fibrous capsular contracture is a common complication following breast implant surgery. While formation of a fibrous tissue capsule surrounding the implant is a normal physiological response to a foreign body, not all capsules contract. Contracture of the fibrous capsular tissue surrounding the implant may result in firmness, discomfort or pain in the breast, distortion of the breast, palpability of the implant, or displacement of the implant. The etiology of capsular contracture is unknown, but is most likely multifactorial. Contracture develops to varying degrees, unilaterally or bilaterally, and may occur within weeks to years after surgery. Severe cases may require surgical intervention. Capsular contracture may recur after capsulotomy or capsulectomy. Lumps perceived as capsular scarring have delayed the diagnosis of palpable tumour. Questionable lumps must be promptly investigated.

DELAYED WOUND HEALING

Some patients experience delayed healing, and for others the incision site may not heal well. It may open from injury or infection. If the implant is exposed, further surgery will be required. Tissue necrosis is the development of dead tissue around the implant. It will delay wound healing, may cause wound infection and may require surgical correction and/or implant removal. Tissue necrosis has been reported following the use of steroid drugs, chemotherapy, radiation to breast tissue, and smoking, but in some cases it may occur without any known cause.

DISSATISFACTION WITH RESULTS

The complications of incorrect size, misplaced scar location, hypertrophic scarring are usually related to surgical technique. Careful size selection, creation of an appropriate and adequate size surgical pocket, and use of current accepted surgical procedures are the surgeons' responsibility.

DISPLACEMENT OF IMPLANT

Implants may displace with accompanying discomfort and/or distortion in breast shape. Difficult placement techniques may increase the risk of displacement by reducing pocket size and placement accuracy. Displacement may require surgical intervention.

EXPLANTATION OF MAMMARY IMPLANTS

Although there are no definite reports in the medical literature, it is anticipated that a large number of mammary implants will be surgically explanted or replaced. This is not unlike many other plastic surgery procedures which are commonly repeated to maintain patient satisfaction. Due to the wide variety of patients' physical responses to breast implant surgery, differences in surgical techniques and medical treatments, as well as potential complications, patients should be advised that these should not be considered lifetime implants and explant surgery may be indicated at any time.

EXTERNAL CAPSULOTOMY

Treatment of capsular contracture by external manual compression may cause the shell to weaken or rupture. There have been reports of rupture and the manufacturer recommends against the external capsulotomy procedure and is not responsible for the integrity of the implant should manual compression or other external stress techniques be used.

EXTRUSION

Unstable or compromised tissue covering and/or interruption of wound healing may result in exposure and extrusion of the implant. Causes or contributing factors may include infection, wound dehiscence, necrosis with or without infection, capsular contracture, closed capsulotomy, unsuitable skin flap, improper size and placement of implant, and/or tissue erosion associated with implant folds. The incidence of extrusion has been shown to increase when the prosthesis has been placed in injured areas: scarred, heavily irradiated or burned tissue or crushed bone areas; where severe surgical reduction of the area has been performed; and where steroids are used in the surgical pocket.

GRANULOMA

Granuloma formation is a common tissue response to the presence of foreign materials. It is possible for a tissue reaction to be caused by the presence of the implant leading to formation of silicone granuloma.

HAEMATOMA/ EARLY ONSET SEROMA

Haematoma and serous fluid accumulation are complications associated with any type of invasive surgery. Postoperative haematoma and early onset seroma² may contribute to infection and/or capsular contracture. Postoperative haematoma and early onset seroma may be minimised by meticulous attention to haemostasis during surgery, and also possibly by postoperative use of a closed drainage system. Persistent, excessive bleeding must be controlled before the device is implanted. Any postoperative evacuation of haematoma or early seroma must be conducted with care to avoid contamination or damage to the implant.

INFECTION

Pre-existing infection not resolved before implant placement increases the risk of periprosthetic infection. Do not expose the implant or filling accessories to contaminants, which increases the risk of infection.

Infection is an inherent risk following any type of invasive surgery. Infection around a breast implant may occur within days, weeks, or even years, after surgery. Signs of acute infection reported in association with breast implants include erythema, tenderness, fluid accumulation, pain and fever. Signs of subclinical infection may be difficult to detect. Postoperative infections should be treated aggressively according to standard medical practices to avoid more serious complications. Infection that is unresponsive to treatment or necrotizing infection may require implant removal. Capsular contracture may be related to infection in the area surrounding the implant.

PAIN

Any unexplainable pain, not associated with any of the complications listed herein, must be promptly investigated.

¹ Late seroma is defined as a predominant serous accumulation of periprosthetic liquid developing ≥ 1 year post-implantation.

² Early onset seroma is defined as periprosthetic fluid accumulation developing < 1 year post-implantation.

PHLEBITIS SUPERFICIAL

In rare cases, damage to the vein walls may occur during implantation of the mammary prosthesis. Subsequently, this may result in phlebitis superficial, which is inflammation of the veins caused by a blood clot just below the surface of the skin. This complication is rarely serious and with proper care, usually resolves rapidly.

PTOSIS

Ptosis, also known as sagging of the breast, is a deformity which can occur pre- and post-breast surgery. Pre-breast surgery, due to the breast having no bony or cartilaginous support, time and gravity can affect the natural breast inducing glandular and skin stretching. Post-surgery, the size, weight and position of breast implants can have an effect on the degree of ptosis. If serious, this complication can often lead to re-intervention.

RUPTURE OF GEL IMPLANTS

Gel rupture can occur post-operatively from damage to the implant during handling or surgery. Rupture of the shell can also occur from contracture, trauma or excessive manipulation. Ruptures of unknown etiology have also been reported. Because of the cohesiveness of the gel material, some gel ruptures can go undetected unless surgery occurs for another reason (e.g. size exchange). Despite the cohesive properties of the gel, extravasation out of the surgical pocket can occur under pressure, in which case additional surgery may be required to retrieve the gel. Inflammation and formation of silicone granulomas have been reported. If shell rupture is suspected, the implant should be removed.

SENSATION

The risk of temporary or permanent dysesthesia exists following any invasive surgical procedure. Careful surgical technique can minimise, but not preclude, this risk. Dysesthesia of the nipple/areolar complex, and less frequently, the breast area in general, has been reported following implantation and may be temporary or permanent. The risk of neurological impairment increases with more extensive surgery. Surgical implantation of a mammary prosthesis may interfere with the ability to breast feed. However, it should be noted that previous breast reconstruction surgery, may be the initial cause of this interference.

UPPER LIMB LYMPHEDEMA

Lymphedema is a chronic condition characterized by swelling of the arm, hand, or breast which can result from breast cancer treatment due to the accumulation of lymphatic fluid in the interstitial tissues. This condition is known to significantly compromise quality of life due to body image changes, alterations in arm function, and increase the risk of other complications including infection and cellulitis.

Recent studies suggest immediate breast reconstruction may be associated with a reduced risk of post-mastectomy lymphedema. Treatment options include conservative management or surgical treatment using a super-microsurgical technique. However, no complete cure for lymphedema is known.

WRINKLING AND FOLDS

Thin or inadequate overlying tissue, patients with little or no subcutaneous fat, implants which are too large for the surgical pocket or the anatomical structure of the patient, and subcutaneous placement may contribute to palpable, or visible, wrinkles and folds. Folds may result in thinning and erosion of adjacent tissue and erosion of the implant. Signs of skin inflammation, such as tenderness and erythema, may indicate thinning or erosion and must be promptly investigated. Palpable wrinkling / and or folds may be confused with palpable tumour, and questionable cases must be promptly investigated.

INSTRUCTIONS FOR USE

STORAGE AND TRANSIT CONDITIONS

There is no special Storage or Transit conditions associated with Gel-filled Mammary Implants.

SUPPLIED STERILE

Implants are supplied in sterile form using dry heat sterilisation, processed by validated strictly controlled sterilisation cycles. Sterility is verified in accordance with standards. Sterility of the implant is maintained only if the package is intact and undamaged. If the package has been opened or damaged the product must not be used.

PACKAGING

Sterile product is supplied in a sealed, double primary package. Sterility is not guaranteed if the package has been damaged or opened. Tear-Off patient record labels are attached to the primary package. It is recommended that these labels are attached to the patients' records, as applicable, and also that a label is given to the patient to retain for their own information.

TO OPEN PACKAGED STERILE PRODUCT

1. Peel tab on corner of outer tray to open the lid of the outer primary package under clean, aseptic conditions.
2. Remove the inner primary package from the outer primary package. This can be done by either of the following methods:
 - a) The inner tray can be removed using the recessed area.
 - b) Invert the outer primary package over the sterile field, allowing the sealed inner primary package to gently fall into the field.
3. Peel tab on corner of inner tray to open the lid of the inner primary package.

PRODUCT EXAMINATION AND HANDLING

1. Product should be visually examined for any evidence of particulate contamination, damage or leakage.
2. The implant should be kept submerged in sterile water or normal saline prior to implantation to prevent contact with airborne and surgical field particulate contaminants.
3. **Do not** immerse the implant in Betadine or Iodine containing solutions. If Betadine or Iodine containing solutions are used in the pocket, ensure that it is rinsed thoroughly so no residual solution remains in the pocket.
4. Gel-filled mammary implants may contain air bubbles. This occurs normally during sterilisation and has no effect on product integrity or performance.
5. Use disposable, non-shedding wipes only to ensure there are no particles left on the surface of the implant.

CONTAMINATION OF IMPLANTS

Care must be taken to prevent surface contaminants such as talc, dust and skin oils from coming into contact with the implant. Products should be inspected for contamination prior to insertion. Contamination at the time of surgery increases the risk of periprosthetic infection, and possibly capsular contracture.

SURFACE CONTAMINANTS (TALCS, DUST, LINT, OILS) ON SURFACE OF IMPLANTS CAN CAUSE FOREIGN BODY REACTION, HANDLE WITH CARE WITH SURGICAL GLOVES (RINSED FREE OF TALC) WITH STRICT ASEPTIC TECHNIQUE. DO NOT IMPLANT CONTAMINATED PRODUCT.

Back-up implants must be readily available at the time of surgery for use in the event of contamination.

SURGICAL PROCEDURE

Correct surgical procedures and techniques are the responsibility of the medical professional. Each surgeon must evaluate the suitability of the procedure based upon current accepted techniques, individual judgement and experience. Proper size and shape of implants must be determined for the individual patient by the surgeon. An incision should be of appropriate length to accommodate the style, size, and profile of the implant. This will reduce the potential for creating excessive stress to the implant during insertion.

N.B. size of the mammary implant may be presented in terms of mass (grams) or volume (millilitres) on the product labelling. Due to the nature of the silicone gel-fill material mass can be considered equivalent to volume.

DISPOSAL OF DEVICE

There are no special disposal instructions for this device. Normal hospital protocol should be followed to dispose of this device.

RETURN OF GOODS POLICY

Prior to the return of any product:

- a) Authorisation must be requested in writing.
- b) Decontamination certificate must be supplied with any product returned which has been in contact with body fluid.
- c) If it is necessary to puncture the shell of any sealed shell product to assist with safe decontamination sterilisation procedures, the area of mechanical interference should be marked with indelible marker on the surface of the product and reference made on the decontamination certificate.

WARRANTY

The manufacturer warrants that reasonable care was used in the manufacture of these products, and will replace any product that the manufacturer's investigation shows to be defective at time of shipment. Patient selection, surgical procedures, post-surgical treatment and stresses, and handling of the devices are totally the responsibility of the customer. The manufacturer has no control over the conditions of use and cannot warrant good effect or against ill effect following the use of the device, and shall not be responsible for any incidental or consequential loss or damage or expenses directly or indirectly arising from the use of this device. All other warranties (whether implied by law or otherwise) are excluded to the fullest extent permitted by law.

CAUTION

Federal (USA) law restricts this device to sale by or on the order of a physician.

BIBLIOGRAPHIC REFERENCES

Literature references are available upon request.

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NAGOR LTD
GC Aesthetics (Management) Limited