

[Company logo]	DNA Quantitation - Qubit: dsDNA High Sensitivity		Page 1 of 4
	Qubit dsDNA High Sensitivity Quant Kit Lot No.		Version #
Operator Name		Signature	Date

1. Formulate Qubit Working Solution Formulation						
NOTE:	The Qubit dsDNA HS Assay is a useful tool for quantifying cfDNA samples after cfDNA Extraction , whereas the Qubit dsDNA BR Assay is more appropriate for quantifying pools after Barcoding . Please see the Qubit: dsDNA Broad Range Quantitation bench worksheet for the latter procedure.					
NOTE:	Qubit dsDNA HS Assay reagents are very sensitive to light--keep in the dark whenever possible. Qubit Reagent and Buffer are stored at Room Temperature in the dark. Qubit Standards are stored at 4°C in the dark.					
___ 1.	Warm Qubit HS Standards to room temperature (15°C to 25°C) protected from light for 30 minutes.					
	Start time:		End time:			
___ 2.	Vortex cfDNA sample tubes for 5 seconds and centrifuge for 5 - 10 seconds.					
___ 3.	Vortex Qubit dsDNA HS Buffer and Qubit dsDNA HS Standards 1 and 2.					
___ 4.	Label two Qubit assay tubes for standards "S1" and "S2."					
___ 5.	Label each Qubit sample tube <Sample ID>					
___ 6.	Prepare Qubit Working Solution in a LoBind 1.5 mL microcentrifuge tube according to the following formulation:					
	Reagent	µL/sample	x	Total # of samples (incl. 2 standards)	+ 20% Overage	Total µL
	Qubit dsDNA HS Buffer	195-199			1.2	
	Qubit dsDNA HS Reagent	1-5	x		=	
___ 7.	Vortex Qubit Working Solution for 5 seconds and centrifuge for 5 - 10 seconds.					

2. Prepare Qubit Standards and Samples						
___ 1.	Dispense 190 µL of Qubit Working Solution into each Qubit standard tube.					
___ 2.	Dispense 199 µL of Qubit Working Solution to each Qubit sample tube.					
___ 3.	Dispense 10 µL of Qubit dsDNA HS Standard 1 into the S1 tube and 10 µL of Standard 2 into the S2 tube.					
___ 4.	Dispense 1 µL of each sample to its respective Qubit sample tube.					
NOTE:	Ensure that bubbles are not introduced to any of the standards or samples.					
NOTE:	Hold the Qubit assay tubes from the top, as temperature fluctuations will affect the fluorescence measurements.					
___ 5.	Vortex for 5 seconds, then centrifuge for 5 to 10 seconds.					
___ 6.	Confirm that each Qubit standard/sample tube is free of bubbles and particulate, or if present, remove them.					
___ 7.	Incubate at room temperature (15°C to 25°C) in the dark for 5 minutes.					
	Start time:		End time:			
NOTE:	Do not store Qubit assay samples or standards for more than 3 hours.					


[Company logo]	DNA Quantitation - Qubit: dsDNA High Sensitivity	Page 3 of 4
	Qubit dsDNA High Sensitivity Quant Kit Lot No.	Version #
Analyst Name	Signature	Date

Qubit 2.0 Screen Examples

- After inserting a sample tube and clicking “Read Next Sample”, below is what the Qubit screen will look like. This is the Raw Qubit value (circled in orange)
- To calculate the ng/ μ L stock concentration click “Calculate Stock Conc.” (circled in green)



- After clicking “Calculate Stock Conc.”, the screen will look as below
- Adjust “Volume of Sample Used” (circled in red) to 1 μ L. 1 μ L should line up the red dashes on the sides.
- Select ng/ μ L as the units (circled in green)
- Screen will now display your ng/ μ L concentration (circled in yellow)



[Company logo]	DNA Quantitation - Qubit: dsDNA High Sensitivity	Page 4 of 4
	Qubit dsDNA High Sensitivity Quant Kit Lot No.	Version #
Analyst Name	Signature	Date

Qubit 3.0 Screen Examples

- Before inserting a sample, below is what the screen will look like
- Select 1 μL as Sample Volume (green box)
- Select $\text{ng}/\mu\text{L}$ as Output Sample Units (yellow box)
- Insert Sample tube and Select "Read tube"



- Qubit concentration in $\text{ng}/\mu\text{L}$ will be shown in the bottom circle (red box)

