

GC Application Note

Quantification of nitrosamines in water by automated PAL SPME-Arrow and GC/MS



FOOD SAFETY





Quantification of nitrosamines in water by automated PAL SPME-Arrow and GC/MS


Peter Egli, Beat Schilling, BGB Analytik AG, Adliswil, Switzerland
Guenter Boehm, CTC Analytics AG, Zwingen, Switzerland

Short summary:

- A method for the quantitative analysis of a number of different N-nitrosamines from aqueous samples has been worked out.
- Carboxen/PDMS SPME fibers and Carbon WR SPME Arrow were compared. Carbon WR SPME Arrow gave a 3-5 x higher extraction yield/detector signal.
- The limits of detection (= S/N > 3) range from 5 ng/L for N-nitroso-dipropylamin to 50 n/L for N-nitroso-pyrolindin-amine.
- Repeatability @ 1 µg/L ranges from 4-12%.

Introduction:

Most nitroso compounds are suspected carcinogens in various animals including human beings. Human exposure to carcinogenic N-nitrosamines may result directly from ingestion or inhalation of preformed compounds from the environment. They are formed by the reaction of amines, especially secondary and tertiary amines or amino group containing compounds with nitrite (present e.g. as food preservative). N-nitroso-dimethylamine (NDMA) and N-nitroso-diethylamine (NDEA) are the most common nitrosamines found in food materials.



Experimental:

1. Immersion Extraction

Sample	10ml water in 20ml HS vial with 4g KCl
Fiber	Carbon WR SPME Arrow (PAL Sytem No: ARR11-C-WR-120-20-P1) Carboxen/PDMS SPME Fiber (part of Supelco kit, No: 57284-U)
Pre conditioning	0:30 min
Pre incubation time	1:00 min
Incubation temp	50°C
Agitation speed	500 rpm
Needle penetration	22 mm
Fiber penetration	30 mm
Extraction time	30 min
Desorption time/temp	2:00 min/ 260°C

2. Chemicals

Nitrosamines calibration standard (Restek No: 31898)
N-nitroso-dimethylamine
N-nitroso-methylethylamine
N-nitroso-diethylamine
N-nitroso-dipropylamine
N-nitroso-dibutylamine
N-nitroso-pyrolidone
N-nitroso-piperidine
Diluted in Methyl tert-butyl ether (Fluka No: 20249)

3. GC/MS

GC	Varian 3400
MS	Varian Saturn ion trap
Column	60 m; 0.25 mm ; 25 µm BGB Wax 20M
Carriergas	Hydrogen 10.0 psi
Temp. program	50°C for 1min, then 10°C/min to 240°C
Injector	260°C
Mass Range	42 bis 250 m/z

With the PAL RTC and RSI the entire SPME process is fully automated. This guarantees process safety and high reproducibility.



Overlaid Chromatogram Plots

Plot 1: c:\... \desktop\spme\nitrosamine\iniam in 22 .sm s Ions: 74+88+102+102+130+84+114+100 all

Plot 2: c:\... \desktop\spme\nitrosamine\iniam in 26 .sm s Ions: 74+88+102+102+130+84+114+100 all




Fig. 1: Chromatograms of 1 µL of a 1 µg/L standard solution and blank direct injection





Fig. 2: Comparison of relative detector (MS) response for different nitrosamines @ 1 µg/L between Carbon WR SPME Arrow and Carboxen/PDMS SPME

	N-nitroso-dimethylamine	N-nitroso-methylethylamine	N-nitroso-diethylamine	N-nitroso-dipropylamine	N-nitroso-dibutylamine	N-nitroso-piperidine	N-nitroso-pyrrolidine
Average (n=5)	153621	1105607	430555	566316	218648	218468	264031
% RSD	11.9	8.8	5.9	4.2	6.8	14.1	10.4

Tab. 1: Reproducibility for the quantification of different nitrosamines @ 1 µg/L with Carbon WR SPME Arrow.



N-nitroso-diethylamine @ 50 ng/L, S/N 24 N-nitroso-dipropylamine @ 50 ng/L, S/N 34 N-nitroso-pyrrolidinamine @ 50 ng/L, S/N 3 (=LOD)

Fig. 3: Chromatograms for selected N-nitrosamines, @ 1000 ng/L and 50 ng/L with Carbon WR SPME Arrow.






Fig. 4, 5, 6: Calibration curves for N-nitroso-diethylamine (top), N-nitroso-dipropylamine (middle) and N-nitroso-pyrrolidineamine (bottom)

Conclusions:

- A method for the quantitative analysis of a number of different N-nitrosoamines from aqueous samples has been worked out.
- Based on previous publications (ref.1) Carboxen/PDMS SPME fibers have been selected as sorbent material
- Carboxen/PDMS SPME fibers and Carbon WR SPME Arrow were compared. Carbon WR SPME Arrow gave a 3-5 x higher extraction yield/detector signal.
- The limits of detection (= S/N > 3) range from 5 ng/L for N-nitroso-dipropylamine to 50 ng/L for N-nitroso-pyrrolidineamine.
- Repeatability @ 1 µg/L ranges from 4-12% @ 1 µg/L.
- Time for one sample was 30 min extraction + 30 min GC runtime. With overlapped extraction (e.g. with PAL Sample Control <http://www.palsystem.com/index.php?id=243>, or many other chromatographic data systems) the total runtime for 7 samples was approximately 240 min.

References:

- [1] Grebel JE, Young CC, Suffet IH, J. Chrom. A, 1117 (2006) 11-18

Legal Statements

CTC Analytics AG reserves the right to make improvements and/or changes to the product(s) described in this document at any time without prior notice.

CTC Analytics AG makes no warranty of any kind pertaining to this product, including but not limited to implied warranties of merchantability and suitability for a particular purpose.

Under no circumstances shall CTC Analytics AG be held liable for any coincidental damage or damages arising as a consequence of or from the use of this document.

© 2016 CTC Analytics AG. All rights reserved. Neither this publication nor any part hereof may be copied, photocopied, reproduced, translated, distributed or reduced to electronic medium or machine readable form without the prior written permission from CTC Analytics AG, except as permitted under copyright laws.

CTC Analytics AG acknowledges all trade names and trademarks used as the property of their respective owners.

PAL is a registered trademark of CTC Analytics AG | Switzerland

Imprint

Date of print: 02.2016

CTC Analytics AG
Industriestrasse 20
CH-4222 Zwingen
Switzerland
T +41 61 765 81 00
Contact: info@ctc.ch