

## An Innovative Approach to Coffee Characterization, Using Size Exclusion Chromatography and High Resolution Gas Chromatography

*A combination of high performance SEC and GC were applied to the analysis of six coffee samples. Data obtained from the two techniques were statistically evaluated by principal components analysis (PCA). The results confirmed the potential of these techniques for coffee analysis.*

### Key Words:

- coffee quality
- size exclusion chromatography
- principal components analysis
- headspace extraction

Coffee beans are one of the largest traded commodities at international levels — important to producers and consumers. Coffee bean quality traditionally has been judged by specially trained individuals who taste, smell, and visually examine the beans (sensory analysis).

Today's international business trade requirements are stricter. However, methods for detecting trading fraud and chemical standards for characterizing coffee beans have yet to be developed. Efficient methods for characterizing coffee beans are becoming a necessity.

New techniques for characterizing coffee beans combine GC or reversed phase LC with sensory analysis. High resolution GC, using headspace extraction, is the preferred method for analyzing components related to aroma. This method provides a better profile of volatile chemical compounds when compared to distillation techniques, which tend to form inhibitor compounds that interfere with some analytical protocols. HPLC techniques define flavor components rather than aroma components.

Investigators at the Universidade Federal do Rio de Janeiro, Instituto de Química, Brazil, characterized coffee beans by using size exclusion chromatography (SEC), headspace-GC, and principal component analysis (PCA). Six samples of green coffee beans from various countries of origin were first identified as high or low quality by sensory analysis — beans A, B, and C were high quality; and beans D, E, and F were low quality. These beans were then analyzed by the following methods.


### Size Exclusion Chromatography

To extract a sample for SEC analysis, a peeled bean was granulated and passed through a 0.75mm mesh filter, then mixed with distilled water at 80°C for 15 minutes. A 0.5g sample was withdrawn, cooled to ambient temperature, filtered into a flask, and taken to 50mL volume with distilled water. The samples were separated on a TSKgel G3000SW column, with UV detection and differential refraction index detection (Figure A).

SEC excludes, or fractions, compounds according to their molecular weight — molecules with larger molecular weights are eluted from the column before those with smaller molecular weights. The TSKgel G3000SW column was calibrated using protein standards, which split the chromatogram into several regions. The chromatograms in Figure A present a distinct profile showing low and high molecular weight regions.

**Figure A. Low-Quality Coffee Bean by SEC**

Column: **TSKgel G3000SW, 30cm x 7.5mm ID, 10µm particles**  
 Cat. No.: **805789**  
 Mobile Phase: 0.05%  $N_3Na$  in twice-distilled water  
 Flow Rate: 0.5mL/min  
 Det.: UV, 272nm, 0.16 AUFS  
 differential refraction index, 64 scale factor of 40  
 Inj.: 10µL of prepared sample from coffee bean D



The regions up to 21 minutes consist of larger molecules with molecular weights greater than  $1.2 \times 10^4$ . The regions from 21 to 40 minutes consist of molecules with low molecular weights, less than  $1.2 \times 10^4$ . The latter regions include very important substances — total chlorogenic acid, trigonelline and caffeine — shown respectively in Figure A as peaks 1, 2, and 3 (1).

The compounds shown as peaks A, B, and C contribute 90% of the varieties' quality attributes among the samples — the peak areas are the principal parameters representing the variability of the analyses. Therefore, these small peaks represent the final quality characteristics of coffee preparations.

In addition to the conventional size exclusion mechanism, the presence of some hydrophobic interactions may explain the better resolution of the low molecular weight regions, relative to the higher molecular weight regions.

### Headspace-Gas Chromatography

Headspace-GC analysis was performed according to the method described by De Maria et al. (2). A coffee bean sample of 0.5g was roasted at 280°C for 14 minutes in a vacuum-sealed test tube. 4mL of the resulting gaseous phase was injected into the GC. These samples were separated on a bonded phase SUPELCOWAX™ 10 column. Figure B shows the large quantity of peaks resulting from this procedure, which demonstrates the complexity of the samples. However, the number of extraneous peaks decreases the ease in identifying those that are important quality characteristics.

### Principal Components Analysis

Results from the chromatographic analyses were processed by the statistical method PCA. In PCA the following methods are used to assess discrimination stability: PCA, modified PCA, and discriminate function analysis with cross validation. System integrity and reproducibility are maintained as long as the original sample conditions are kept constant.


Highly correlated variables are analyzed using nonlinear multiple regression analysis, partial regression analysis, and least square techniques. Sensory test data were correlated to the analytical techniques used, for each green coffee sample. Variation of a complex system of variables can be expressed by a simple, graphical representation using PCA.

Results were analyzed statistically using the peak area but, due to the complexity of the headspace-GC chromatograms, a relationship between the techniques and the chromatograms was needed to correlate the data. Variation coefficients were calculated from characteristic peaks on ten repetitive runs per sample. Peaks having a variation coefficient of less than 15% were selected and subjected to PCA, individually and as a group.

To examine sample dispersion, a scatterplot graph was made with the two principal values corresponding to the main variation peaks (Figure C). From these graphs we can tell that the dispersion was high. These results were seen during previous attempts, and dispersion was higher when each peak from each chromatographic technique was analyzed individually. We suggest that further grouping can remain according to sample characteristics, as long as testing continues and sample size increases.

**Figure B. Coffee Bean by Headspace-GC**

Column: SUPELCOWAX 10, 30cm x 0.25mm ID, 15µm film  
Cat. No.: 24287  
Oven: 40°C (10 min) to 190°C (15 min) at 3°C/min  
Carrier: hydrogen, 40cm/sec  
Det.: FID, 280°C  
Inj.: 4mL headspace of coffee bean



797-0570

**Figure C. Scatterplot Graph Showing Sample Dispersion for Each Coffee Bean**




Table 1 offers a relationship of the percent variance obtained through PCA for each peak area — SEC compared to traditional reversed phase chromatography. SEC performed as well as reversed phase chromatography. In fact, SEC presents certain advantages — SEC is capable of simultaneous analysis of the important peaks, whereas reversed phase HPLC requires different chromatographic conditions for each compound (1). Additionally, SEC uses water as the mobile phase, which is attractive from a cost containment perspective.

Further studies are necessary to provide detailed identification of components. Additional complementary analytical techniques, such as mass spectrometry and GC or LC, may help obtain a detailed characterization profile of most aroma and bouquet compounds for coffee.

**Table 1. Average Results for CGA, Trigonelline, and Caffeine**

Samples	CGA		Methods Trigonelline		Caffeine	
	SEC	RP	SEC	RP	SEC	RP
Mean	8.0	9.0	1.1	1.1	1.1	1.2
Std. Dev.	0.32	0.31	0.05	0.04	0.03	0.04
Coef. Var. (%)	4.0	3.4	4.5	3.6	2.7	3.3
B	6.7	7.7	1.2	1.2	1.1	1.3
C	7.4	8.1	1.0	1.1	1.1	1.1
D	7.2	8.0	1.6	1.6	0.7	0.7
E	6.8	7.1	1.5	1.4	1.1	1.3
F	6.8	7.0	1.2	1.2	1.2	1.3

SEC — size exclusion chromatography; RP — reversed phase liquid chromatography, a different technique was used for each component

A — results are means for 12 replicates originating the respective standard deviation and coefficient of variation; B-F — results are means for duplicate determinations.

**Ordering Information:**

Description	Cat. No.
<b>TSKgel G3000SW Gel Filtration Column</b> 30cm x 7.5mm ID, 10µm particles	<b>805789</b>
<b>SUPELCOWAX 10 Capillary Column</b> 30m x 0.25mm ID, 0.15µm film	<b>24287</b>

**References**

1. De Maria, C. A. B., L. C. Trugo, and R. F. A. Morreira, *Food Chemistry*, **52**: 447 (1995).
2. De Maria, C. A. B., L. C. Trugo, R. F. A. Morreira, and C. C. Werneck, *Food Chemistry*, **49**: 141 (1994).

References not available from Supelco.

Information in this bulletin used with permission from Química Nova, **20**: 1:5-8, (1997) Professor Luiz C. Trugo, Universidade Federal do Rio de Janeiro, Instituto de Química, Brazil, and Sociedade Brasileira de Química, Professors Ricardo F.A. Moreira and Carlos A.B. De Maria.

SUPELCOWAX is a trademark of Sigma-Aldrich Co.

**For more information on these analyses and PCA, please contact Professor Luiz C. Trugo directly — email: cabm@acd.ufrj.br.**

## Free Technical Literature From Supelco

Supelco produces a wide variety of technical literature, all of which is free. You may request this literature several different ways:

- Use the business reply card located at the back of our catalog.
- Call our Ordering and Customer Service department (800-247-6628 or 814-359-3441).
- Use our ChromFax™ service.

## How to Use ChromFax

ChromFax is Supelco's automated digital data system, which allows you to instantly obtain technical literature any time of the day or night. Supelco technical bulletins, application notes, lab hints, and product specification sheets can be sent to your Fax machine.

1. To use ChromFax you must dial from a touch-tone phone. If you do not have a touch-tone phone, call our Ordering and Customer Service department and we will Fax the publication to you.
2. Have the phone number of your Fax machine ready, and dial **800-652-8069 or 814-359-5748**. Remember you must have the *complete* Fax number (011 + country/city code + Fax number) when you are calling from outside the US or Canada, otherwise the document cannot be transmitted.
3. When ChromFax answers, follow the simple recorded instructions. Select the **ChromFax number** of the document you want. Be certain you order documents by their ChromFax numbers and NOT their literature numbers.
4. Enter your Fax number and hang up. The requested information will be sent to you within minutes. To report a problem with the ChromFax system contact our Technical Service department.

## Technical Literature Available on ChromFax

The listing of documents available on ChromFax is updated every month. We encourage you to request a list of available ChromFax documents each time you use the ChromFax system.

Contact our Technical Service Department  
(phone 800-359-3041 or 814-359-3041, FAX 800-359-3044  
or 814-359-5468) for expert answers to your questions.

BULLETIN 919

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website ([www.sigma-aldrich.com](http://www.sigma-aldrich.com)), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA · Sigma-Aldrich de Argentina, S.A. · Buenos Aires 1119 AUSTRALIA · Sigma-Aldrich Pty. Ltd. · Castle Hill NSW 2154 AUSTRIA · Sigma-Aldrich Handels GmbH · A-1110 Wien  
BELGIUM · Sigma-Aldrich N.V./S.A. · B-2880 Bornem BRAZIL · Sigma-Aldrich Quimica Brasil Ltda. · 01239-010 São Paulo, SP CANADA · Sigma-Aldrich Canada, Ltd. · 2149 Winston Park Dr., Oakville, ON L6H 6J8  
CZECH REPUBLIC · Sigma-Aldrich s.r.o. · 186 00 Praha 8 DENMARK · Sigma-Aldrich Denmark A/S · DK-2665 Vallensbaek Strand FINLAND · Sigma-Aldrich Finland/YA-Kemia Oy · FIN-00700 Helsinki  
FRANCE · Sigma-Aldrich Chimie · 38297 Saint-Quentin-Fallavier Cedex GERMANY · Sigma-Aldrich Chemie GmbH · D-82041 Deisenhofen GREECE · Sigma-Aldrich (o.m.) Ltd. · Ilioupoli 16346, Athens  
HUNGARY · Sigma-Aldrich Kft. · H-1067 Budapest INDIA · Sigma-Aldrich Co. · Bangalore 560 048 IRELAND · Sigma-Aldrich Ireland Ltd. · Dublin 24 ISRAEL · Sigma Israel Chemicals Ltd. · Rehovot 76100  
ITALY · Sigma-Aldrich s.r.l. · 20151 Milano JAPAN · Sigma-Aldrich Japan K.K. · Chuo-ku, Tokyo 103 KOREA · Sigma-Aldrich Korea · Seoul MALAYSIA · Sigma-Aldrich (M) Sdn. Bhd. · Selangor  
MEXICO · Sigma-Aldrich Química S.A. de C.V. · 50200 Toluca NETHERLANDS · Sigma-Aldrich Chemie BV · 3330 AA Zwijndrecht NORWAY · Sigma-Aldrich Norway · Torshov · N-0401 Oslo  
POLAND · Sigma-Aldrich Sp. z o.o. · 61-663 Poznań PORTUGAL · Sigma-Aldrich Quimica, S.A. · Sintra 2710 RUSSIA · Sigma-Aldrich Russia · Moscow 103062 SINGAPORE · Sigma-Aldrich Pte. Ltd.  
SOUTH AFRICA · Sigma-Aldrich (pty) Ltd. · Jet Park 1459 SPAIN · Sigma-Aldrich Quimica, S.A. · 28100 Alcobendas, Madrid SWEDEN · Sigma-Aldrich Sweden AB · 135 70 Stockholm  
SWITZERLAND · Supelco · CH-9471 Buchs UNITED KINGDOM · Sigma-Aldrich Company Ltd. · Poole, Dorset BH12 4QH  
UNITED STATES · Supelco · Supelco Park · Bellefonte, PA 16823-0048 · Phone 800-247-6628 or 814-359-3441 · Fax 800-447-3044 or 814-359-3044 · email: [supelco@sial.com](mailto:supelco@sial.com)

H

Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.

BMB