

MASSPREP OLIGONUCLEOTIDE SEPARATION TECHNOLOGY STANDARD

I. INTRODUCTION

The pre-packaged MassPREP™ Oligonucleotide Separation Technology (OST) Standard is designed for verification of HPLC/UPLC® instrument and column performance for analysis of synthetic oligonucleotides. Approximately equimolar amounts of 15, 20, 25, 30 and 35 nucleotide (nt) long oligodeoxythymidines are lyophilized and packaged in 1.5 ml LC vials. These vials are vacuum-sealed in foil pouches to reduce degradation that can occur by excessive exposure to light and air. Approximately 1 nmole of each oligonucleotide is present in the vial.

II. RECOMMENDED USAGE

Waters MassPREP OST standard is intended for testing the column and HPLC/UPLC systems separation performance. For selection of Waters XBridge™ OST C18, 2.5 µm HPLC or ACQUITY UPLC® OST C18, 1.7 µm UPLC column, see document 720002008EN on waters.com/library. Both columns and LC system performance can be tested using this dedicated MassPREP OST standard. The most common problems in oligonucleotide analysis that can be detected using OST standard are listed in the Table 1.

CONTENTS

I. INTRODUCTION

II. RECOMMENDED USAGE

III. PREPARATION PROCEDURE

IV. STORAGE AND STABILITY

V. ORDERING INFORMATION

Table 1: System Troubleshooting using MassPREP OST Standard

Chromatogram Appearance	Potential problem
Peaks elute outside of the expected retention time window	Incorrectly prepared mobile phase, aged mobile phase, incorrect column temperature setting; excessive gradient delay (method transfer between different LC systems).
Tailing peaks	Poor tubing connections (especially between column and detector); column bed deterioration, column contamination.
Inconsistent peak width, split peaks	Inadequate gradient mixing, incompatible sample solvent or weak wash (purge solvent for Alliance® HT).
Lost resolution	Aged or contaminated column. (Has column backpressure changed?)

III. PREPARATION PROCEDURE

The following procedure is provided as a general guideline for MassPREP OST standard reconstitution. The described method only serves as a starting point. Depending on the specific application, one may consider using other solvents and/or dilutions. Recommended chemicals to prepare sample diluent and LC mobile phase are following: Acetic Acid (2M):Triethylamine (2M), Fluka, P/N 09748; 1,1,1,3,3,3-hexafluoro-2-propanol (HFIP), Fluka, P/N 52512; triethylamine (TEA), Sigma, P/N 417283.

- Prepare 100 mM of triethylammonium acetate (TEAA) by diluting a 2M:2M stock solution of Acetic Acid:Triethylamine 20 fold in deionized or HPLC grade water.
- Add 0.5 mL of 100 mM TEAA in MassPREP OST standard vial. The final concentration is ~2 pmole/μL for each oligonucleotide. Vortex the vial briefly to dissolve and homogenize the sample.
- For MassPREP OST standard analysis in UPLC mode, prepare the following mobile phases:
 - Mobile phase A is 15 mM TEA, 400 mM HFIP in water. Add 8.31 mL (13.44 g) of HFIP into 191.3 g of water, add 416 μL of TEA; the final volume is 200 mL, the buffer pH is ~7.9.
 - Mobile phase B is 50% A, 50% MeOH. Prepare TEA-HFIP buffer as described above, and add 158.2 g of MeOH (200 mL).
 - The LC solvents are highly volatile; the mobile phase should be used only 24 hours only before making a fresh one. Keep the mobile phase containers sealed to minimize the buffer evaporation.
- Prime the ACQUITY UPLC system and connect a 2.1 x 50 mm ACQUITY UPLC OST C18, 1.7 μm column (P/N 186003949). Setup flow rate at 0.2 mL/min, and column temperature at 60 °C. Equilibrate with initial mobile phase conditions (38% B) for ~20 minutes.
- Inject 10 μL of the MassPREP OST standard. Run the gradient from 38 to 50% B in 12 minutes. The example of chromatogram is shown in Figure 1.
- The MassPREP OST sample can also be used to troubleshoot XBridge OST C18, 2.5 μm columns configured to a HPLC system. The resolution is not expected to match the performance obtained with a ACQUITY UPLC OST C18, 1.7 μm column and UPLC System.
- The alternative conditions for HPLC (UPLC) oligonucleotide analysis use TEAA ion-pairing system:
 - Mobile phase A is 100 mM TEAA (measure 190 g of water, add 10 mL of 2M:2M Acetic Acid:Triethylamine Stock Solution. The pH adjustment is not necessary.
 - Mobile phase B: 20% acetonitrile, 80% TEAA. Prepare 200 mL of TEAA as above, and add 50 mL (39.3 g) of acetonitrile.
 - The most useful analytical column for HPLC is a 2.1 x 50 mm, XBridge OST C18, 2.5 μm column (P/N 186003952). Flow rate is set to 0.2 mL/min, column temperature to 60 °C. Gradient is 40-60% B in 25 minutes. Approximately 20-40 μL of the MassPREP OST sample is typically injected on column (80 pmole per peak). Example of chromatogram is not shown.




Figure 1: ACQUITY UPLC analysis of MassPREP OST standard on an ACQUITY UPLC OST C18, 1.7 μm column. The main components are labeled. Small peaks eluting between labeled oligonucleotides are N-1, N-2, etc. failure sequences generated during the oligonucleotide syntheses. The ACQUITY UPLC system is equipped with 50 μL standard mixer and PDA detector (260 nm).

IV. STORAGE AND STABILITY

It is recommended that the MassPREP OST Standard remain in the foil pouch (in the dark) at room temperature until use. Once reconstituted in appropriate solvents, the MassPREP OST Standard solution is stable for one week when stored at 4 °C. Stability can be extended by freezing samples at -20 °C.

V. ORDERING INFORMATION

Description	Particle Size	Pore Size	Dimension	Part No.
ACQUITY UPLC OST C18 *	1.7 μm	135 Å	2.1 x 50 mm	186003949
ACQUITY UPLC OST C18 *	1.7 μm	135 Å	2.1 x 100 mm	186003950
Custom ACQUITY UPLC OST C18 *	--	--	--	186003951
XBridge OST C18	2.5 μm	135 Å	2.1 x 50 mm	186003952
XBridge OST C18	2.5 μm	135 Å	4.6 x 50 mm	186003953
XBridge OST C18	2.5 μm	135 Å	10 x 50 mm	186003954
Custom XBridge OST C18	--	--	--	186003955
MassPREP OST Standard			1 Vial	186004135

* For use on Waters ACQUITY UPLC Systems

† MSDS information: Waters MassPREP OST Standard contains a non hazardous peptide in powder form. It is for laboratory use only; solely for purposes of scientific experimentation, analysis, research or development. It should be used only as directed and in accordance with good laboratory practice. It is non hazardous as defined by OSHA's Hazard Communication Standard, 29 CFR 1910.1200(c) and Appendix A and B, and are thus exempt from the requirements of the Hazard Communication Standard, including the requirement to supply MSDSs.

Sales Offices:

Austria and European Export
(Central South Eastern Europe,
CIS and Middle East) 431 877 18 07

Australia 2 9933 1777

Belgium 32 2 726 1000

Brazil 55 11 5094 3788

Canada 800 252 4752

China 8621 6495 6999

CIS/Russia +7 495 3367000

Czech Republic 42 02 617 11384

Denmark 45 46 59 8080

Finland +358 9 5659 6288

France (33) 1 30 48 72 00

Germany 49 6196 400600

Hong Kong 852 29 64 1800

Hungary 36 1 350 5086

India and India Subcontinent
91 80 2 837 1900

Ireland 353 1 448 1500

Italy 39 02 274 211

Japan (81) 3 3471 7191

Korea (82) 2 820 2700

Mexico 5255 5200 1860

The Netherlands +31 (0)76-50 87 200

Norway 47 63 84 60 50

Poland (48) 22 833 4400

Puerto Rico 787 747 8445

Singapore 65 6273 1221

Spain 34 93 600 93 00

Sweden 46 8 555 11500

Switzerland 41 62 889 2030

Taiwan 886 2 2543 1898

United Kingdom 44 208 238 6100



© 2007 Waters Corporation. Waters, The Science of What's Possible, MassPREP, XBridge, ACQUITY, Alliance and UPLC are trademarks of Waters Corporation.

October 2007 715001677 Rev A VW-PDF

Waters Corporation
34 Maple Street
Milford, MA 01757 U.S.A.
T: 1 508 478 2000
F: 1 508 872 1990
www.waters.com