

# Addressing the World Helium Shortage For Gas Chromatography

Helium Conservation and  
Converting GC Methods to  
Nitrogen and Hydrogen  
Carrier Gas

James D. McCurry, Ph.D.  
Senior Scientist  
Agilent Technologies

# Presentation Outline

## Carrier Gas Decision Tree

- Decision making guide to fit your carrier gas requirements

## Helium Conservation


- Smarter helium use with new hardware/software tools
- No need to revalidate existing GC methods

## Migrating Existing Helium GC Methods to H<sub>2</sub> and N<sub>2</sub>

- Best practices for obtaining the same results and minimizing method revalidation

# Carrier Gas Decision Tree

Continue using helium, but in a smarter way



Method revalidation not required!


# Reducing Helium Use With Conservation

## New 7890B Helium Conservation Module

- Automatically switches carrier gas supply to N<sub>2</sub> Standby during idle time
- Integrates into the new 7890B Sleep and Wake function
- Combined with Helium Gas Saver to **GREATLY** reduce helium consumption
- Better alternative to just “shutting off the GC”
  - No system contamination with ambient air exposure
  - Faster re-start of heated zones

# Helium Conservation Module


Seamlessly integrated onto 7890 GC hardware and software



- Built on 5<sup>th</sup> generation EPC
- Fully controlled by Agilent data systems
- Purge channel prevents cross contamination of gases
- Precise pressure control between tank and GC
- Switch between gases within 15-30 min for most detectors


# How Does It Work?

## Helium Savings Mode (Nitrogen Carrier, or Sleep Mode)




# How Does It Work?

## Normal Operation Mode (Helium Carrier or Wake Mode)



# How It Works: Configuring Sleep/Wake Operation


## Simple, Straight Forward Setup






# Performance: No Change in Chromatography After N<sub>2</sub> Carrier Sleep Method. GC/FID Analysis


Day 1 - Original He carrier gas run




Day 2 – First He carrier gas run after overnight N<sub>2</sub> sleep method



Day 3 – First He carrier gas run after overnight N<sub>2</sub> sleep method



# Performance: Pass MS Tune **Within 15min** After Switching From N<sub>2</sub> To He As Carrier. GC/MSD



Time (min)	Counts of Nitrogen Ion			
	5 mL/min He	Relative to Saturation	2 mL/min He	Relative to Saturation
3	1735168	20.69%	8388096	100.00%
4	1033280	12.32%	4959232	59.12%
5	590080	7.03%	1618944	19.30%
6	354112	4.22%	722944	8.62%
7	228480	2.72%	333696	3.98%
10	56984	0.68%	102576	1.22%
15	9052	0.11%	17080	0.20%

# Helium Savings Calculator – Single GC Channel

## Extend helium supply and lower cost using conservation techniques



Agilent Technologies

**Method:** ASTM D4815 - Ethanol in Gasoline

**Column:** PDMS 30m x 0.53mm x 2.65um

### GC Flow Conditions

He Carrier Flow (mL/min):	8
He Split flow (mL/min):	70
Gas Saver Flow (mL/min):	20
Gas Saver On (min):	3
Run Time(min.):	20
Gas Volume in Cylinder (L):	8000
Runs per Day:	30
He Cylinder Cost (\$):	300
N2 Cylinder Cost (\$):	60

Parameter	No Conservation	Helium Conservation
Daily He Usage (L)	112	21
<b>He Cylinder Life (days)</b>	<b>71</b>	<b>376</b>
Daily N <sub>2</sub> Usage (L)	0	24
N <sub>2</sub> Cylinder Life (days)	0	340
Yearly He Cost (\$)	\$1,537	\$292
Yearly N2 Cost (\$)	\$0	\$64
<b>Yearly Total Gas Cost (\$)</b>	<b>\$1,537</b>	<b>\$356</b>


## Example

- ASTM Method D4815
  - Widely used to measure ethanol in gasoline
  - Helium cylinder last 2 months under normal operation
- Helium Conservation
  - Helium cylinder life extended to 12 months
  - 4x yearly gas costs per year



# Carrier Gas Decision Tree

## Migrating GC methods to nitrogen and hydrogen



# Safety Considerations for Hydrogen Migration

## GC, GC/MS: Both offer H<sub>2</sub> enabled features

- Agilent H<sub>2</sub> safety letter and safety manuals available
- GC levels of safety design
  - Safety Shutdown
  - Flow Limiting Frit
  - Oven ON/OFF Sequence
- Newer version 6890, 7890 GC and 5773, 5975 MSD offer greater safety than older versions of GC and GC/MSD
- Plumbing Considerations
  - Use chromatographic quality stainless steel tubing
  - Do not use old tubing (H<sub>2</sub> is known as scrubbing agent)
  - Especially don't use old copper tubing (brittleness is a safety concern)

# Source of Hydrogen

## H<sub>2</sub> Generator – preferred

- Very clean H<sub>2</sub>, >99.9999% available
- More consistent purity
- Built-in safety considerations
- Make sure to buy a good one with a low spec for water and oxygen
- Parker's H-MD are used in LFS and SCS



## H<sub>2</sub> Cylinder

- Consider gas clean filter
- Possible to add safety device



# Use N<sub>2</sub> as carrier gas

## Many helium GC methods suited to nitrogen conversion

- Readily available and less expensive gas
- No safety concerns
- Suitable for simple routine analysis (with sufficient resolution)
- More inert than H<sub>2</sub>, especially with PLOT/Micropacked columns
  - Some compounds catalytically reduced in H<sub>2</sub>
- Most helium GC methods have too much resolution
  - lower column efficiency with N<sub>2</sub> won't affect results
- 2-D GC ideally suited to nitrogen
  - column combinations designed to solve specific resolution problems

## Potential issues

- Reduced chromatographic resolution at high flows
- Not suitable for GC/MSD and certain GC detector applications



# Helium Carrier Gas Alternatives

## Important Theoretical Considerations Relating To Peak Efficiency


Sharp, narrow peaks in a chromatogram is an indication of a high efficiency GC column.

- Remember that efficiency is represented mathematically by the symbol " $N$ " called *Theoretical Plates*, and that the larger  $N$  is, the better the resolving power of the column (i.e., higher resolution).
- Resolution is described mathematically by the symbol  $R_s$  and its numeric value tells how well two adjacent peaks are separated from each other.

$$R_s = \frac{\sqrt{N}}{4} \left( \frac{k}{k+1} \right) \left( \frac{\alpha-1}{\alpha} \right)$$

A resolution value of 1.5 tells us that two peaks are baseline separated. The greater (higher) the  $R_s$  value, the more separation that has been achieved.

### Calculating Efficiency



We would like to know the actual time the component spends in the stationary phase.

$$t'_R = t_R - t_m \quad n = \frac{(t'_R)^2}{5.545 W_h^2}$$

$t'_R$  = Corrected Retention Time.

Let's relate "n" to the length of the column.


$$\text{Plates per meter (N)} = \frac{n}{L} \quad \text{or}$$

$$\text{Height equivalent to a theoretical plate (HETP)} = \frac{L}{n}$$

n = effective theoretical plates.

Thus, the more efficient the column, the bigger the "N" the smaller the "HETP".

### Efficiency & Carrier Gas Linear Velocity



Efficiency is a function of the carrier gas linear velocity or flow rate.

The minimum of the curve represents the smallest HETP (or largest plates per meter) and thus the best efficiency. "A" term is not present for capillary columns.

- Plot of HETP versus linear velocity is known as the **Van Deemter** plot.
- The linear velocity value at the minimum of the curve is the optimum value for achieving the best efficiency.



# Helium Carrier Gas Alternatives

## Let's Make This Easy


Sharp, narrow peaks in a chromatogram is an indication of a high efficiency GC column.

- Remember that efficiency is represented mathematically by the symbol "N" called *Theoretical Plates*, and that the larger N is, the better the resolving power of the column (i.e., higher resolution).
- Resolution is described mathematically by the symbol  $R_s$  and its numeric value tells how well two adjacent peaks are separated from each other.

$$R_s = \frac{\sqrt{N}}{4} \left( \frac{k}{k+1} \right) \left( \frac{\alpha-1}{\alpha} \right)$$

A resolution value of 1.5 tells us that two peaks are baseline separated. The greater (higher) the  $R_s$  value, the more separation that has been achieved.

### Calculating Efficiency



Let's relate "n" to the length of the column.

Plates per meter (N) =  $\frac{n}{L}$  or

Height equivalent to a theoretical plate (HETP) =  $\frac{L}{n}$

Thus, the more efficient the column, the bigger the "N" the smaller the "HETP".


We would like to know the actual time the component spends in the stationary phase.

$$t'_R = t_R - t_m \quad n = \frac{5.545}{\left( \frac{t'_R}{W_h} \right)^2}$$

$t'_R$  = Corrected Retention Time.

n = effective theoretical plates.

### Efficiency & Carrier Gas Linear Velocity



Efficiency is a function of the carrier gas linear velocity or flow rate.

The minimum of the curve represents the smallest HETP (or largest plates per meter) and thus the best efficiency. "A" term is not present for capillary columns.

- Plot of HETP versus linear velocity is known as the Van Deemter plot.
- The linear velocity value at the minimum of the curve is the optimum value for achieving the best efficiency.



# Helium Carrier Gas Alternatives

## Let's Make This Easy

- **Goal: change carrier gas while keeping other method conditions the same**
  - use the same column
  - use the same oven program
  - adjust column flow or holdup time to:
    - maintain same peak elution order
    - maintain same peak retention times (or as close as possible)
- Easier method revalidation using this approach
  - minimal changes to timed integration events
  - minimal changes to peak identification table
- For N<sub>2</sub>, test resolution of key components
  - adjust GC conditions (temp, flow) if needed
- Use tools built into the Agilent Chemstation to guide us through the process

# The Tyranny of Van Deemter

## Why Nitrogen Gets a Bad Rap for Capillary GC





- $N_2$  actually provides the best efficiency, but at a slower speed
- Most helium GC methods have too much resolution
  - Lower  $N_2$  efficiency at “typical” helium flows can still provide good enough resolution
- Most GC methods now use constant flow
  - Efficiency losses with temp programming are not as severe

# Many Helium GC Have Excess Resolution


## EN14103 – GC Analysis of FAME content in biodiesel

HP-INNOWax, 30m x 0.25mm ID x 0.25  $\mu$ m



# Helium Carrier Gas Alternative

## Test Case 1: ASTM D6584 for Free and Total Glycerin in Biodiesel



COC Inlet:	Oven Track Mode
Pre-column:	Ultimetal 2m x 0.53mm ID
Column:	Ultimetal DB5HT, 15m x 0.32mm ID x 0.1 df
Column Flow:	Helium at 3.0 mL/min (50 deg C)
Column Pressure:	7.63 psi constant pressure mode
Initial Column Temp:	50 °C for 1 min.
Oven Ramp 1:	15 °C/min to 180 °C
Oven Ramp 2:	7 °C/min to 230 °C
Oven Ramp 3:	30 °C/min to 380 °C, hold 10 min.
Detector:	FID with 25 mL/min N <sub>2</sub> makeup

# Set the Control Mode: Flow or Holdup Time

Try the same flow or holdup time of the original Helium method

The screenshot displays the 'Setup Method' dialog box for an Agilent 7890A. The top section features a graph with 'Run Time, min' on the x-axis (0 to 30) and two y-axes: temperature in °C (0 to 400) and flow rate in mL/min (0 to 3). Three data series are plotted: 'Oven: °C\*' (red solid line), 'Column 1 He: psi' (black solid line), and 'Column 1 He: mL/min' (green dashed line). Below the graph is a toolbar with icons for ALS, Valves, Inlets, Columns, Oven, Detectors, Events, Signals, Configuration, Counters, and Readiness. A table on the left lists method selections:

#	Selection
1	Agilent 123-BD11: 425 °C: 15 m x 320 µm x 0.1 µm Main Segment Heated By Oven inSeg Heated By Oven: 2 m x 530 µm x 0 µm In: Back COC Inlet He Out: Front Detector FID
2	Agilent 123-456: 400 °C: 30 m x 320 µm x 0.5 µm In: Front SS Inlet He Out: Back Detector FID

The 'Control Mode' section is checked 'On'. The 'Setpoint' table is as follows:

Parameter	Setpoint
Flow	3 mL/min
Pressure	7.6296 psi
Average Velocity	52.99 cm/sec
Holdup Time	0.47179 min

Below the setpoint table, a dropdown menu is set to 'Constant Pressure'. A 'Post Run' field is set to '7.6 psi'. A 'Column #1 Configuration' section contains 'Change Column...' and 'Calibrate Column...' buttons. At the bottom are 'OK', 'Apply', 'Cancel', and 'Help' buttons.

# Configure Inlet for Carrier Gas in Chemstation

Setup Method

Agilent 7890A Agilent 7890A Sample Prep Program

Graph: Run Time, min (0 to 30). Y-axis: °C (0 to 400), psi (0 to 8), mL/min (0 to 3). Legend: Oven: °C\* (red solid line), Column 1 He: psi (black solid line), Column 1 He: mL/min (green dashed line).

Configuration Options:

- ALS
- Valves
- Inlets
- Columns
- Oven
- Detectors
- Events
- Signals
- Configuration
- Counters
- Readiness

Miscellaneous | Columns | Modules | ALS

**Front Inlet**  
SS Inlet: He

**Back Inlet**  
COC Inlet: He

**Front Detector**  
FID  
Makeup: N2  
Set Lit Offset with GC Keyboard.

**Back Detector**  
FID  
Makeup: N2  
Set Lit Offset with GC Keyboard.

**Aux EPC 1,2,3**  
Aux EPC 1: N2  
Aux EPC 2: H2  
Aux EPC 3: He

**PCM A**  
PCM A-1: He  
PCM A-2: He  
Channel B Control Mode: Forward Pressure

Buttons: OK, Apply, Cancel, Help

**Select H<sub>2</sub> or N<sub>2</sub>**

# New Windows 7 Method Translation Calculator

## Another useful tool for carrier gas calculations

The screenshot displays the Agilent Method Translator software interface. The window title is "Agilent Technologies Method Translator". The interface is divided into several sections:

- Speed gain:** Set to 1.0000. Radio buttons for "Speed gain", "Translate", and "Best Efficiency" are visible.
- Last file imported:** A text field and icons for file operations.
- Original Method Parameters:** Gas: He. Parameters include Length (30 m), Inner Diameter (320 µm), Film Thickness (0.25 µm), Phase Ratio (320), Inlet Pressure (7.0569 psi), Outlet Flow (1.3158 mL/min), Average Velocity (24.342 cm/s), Outlet Pressure (14.696 psi), Holdup Time (2.0541 min), and Outlet Velocity (30.468 cm/s).
- Calculated Method Parameters:** Gas: N2. Parameters include Length (30 m), Inner Diameter (320 µm), Film Thickness (0.25 µm), Phase Ratio (320), Inlet Pressure (6.4601 psi), Outlet Flow (1.2921 mL/min), Average Velocity (24.342 cm/s), Outlet Pressure (14.696 psi), Holdup Time (2.0541 min), and Outlet Velocity (29.919 cm/s).
- Pressure Units:** Set to PSI.
- Pressure Capacity:** Original Column Capacity: 2.48; Translated Column Capacity: 2.48.
- Temperature Ramps:** Radio buttons for "Isothermal" and "Ramps". A table shows the ramp schedule for both methods.

#	Ramp Rate (°C/min)	Final Temp (°C)	Final Time (min)
Init		60	1
1	5.0000	200	1

#	Ramp Rate (°C/min)	Final Temp (°C)	Final Time (min)
Init		60	1
1	5.0000	200	1


Total Run Time: 30.00 min

- Flexible tool helps convert existing helium methods to alternative carrier
- Built into the New OpenLAB CDS Software
- Can also run as Windows 7 program
- Download from the Agilent Helium Update Page:


[www.agilent.com/chem/heliumupdate](http://www.agilent.com/chem/heliumupdate)



# Wider Retention Time Variation Using the Same Flow as the Original Helium Method



# Same Holdup Time ( $T_r$ ) Gives Consistent Retention Times Compared to Original Helium Method



# ASTM D6584 - Quantitative Results For Alternative Carrier Gas

Carrier gas has no effect on reported results

	Weight Percent		
	Helium	Hydrogen	Nitrogen
<b>Glycerin</b>	<b>0.015</b>	<b>0.014</b>	<b>0.013</b>
<b>Monoglycerides</b>	<b>0.226</b>	<b>0.216</b>	<b>0.223</b>
<b>Total Glycerin</b>	<b>0.097</b>	<b>0.095</b>	<b>0.098</b>

# Test Case 2: Analysis of Oxygenates in Gasoline Using 2-D Gas Chromatography

## ASTM Method D4815 – Oxygenated Additives

- Ethers and alcohols from 0.1 wt% to 15 wt%
- Usually only one or two additives in a sample

Preliminary separation removes light hydrocarbons from sample


- Polar TCEP micro-packed columns retains ethers and alcohols
- Back flush TCEP column to non-polar capillary column (HP-1) to complete analysis

# Configuration and Operation for D4815


Nitrogen carrier gas using original ASTM GC flows conditions

## D4815 Method

<b>carrier gas</b>	<b>nitrogen</b>
Inlet	Split/Splitless
inlet temperature	200 Deg C
<b>TCEP column flow</b>	<b>5 mL/min</b>
<b>split vent flow</b>	<b>70 mL/min</b>
split ratio	15:1
<b>HP-1 column flow</b>	<b>3 mL/min</b>
FID temperature	250 deg C
oven temperature	60 deg C isothermal
Run time	16 min



# Analysis of MtBE and Ethanol in Gasoline using N<sub>2</sub> Carrier Gas



# ASTM Precision Specifications

## D4815 Precision Measures


Compound	Mass %	Repeatability		Reproducibility	
		Spec	Observed	Spec	Observed
Ethanol	0.99	0.06	0.01	0.23	0.01
Ethanol	6.63	0.19	0.03	0.68	0.04
MtBE	2.10	0.08	0.01	0.20	0.01
MtBE	11.29	0.19	0.05	0.61	0.08

## Accuracy Evaluation

Sample	MtBE mass %	
	known	found
SRM2294 #1	10.97	10.61
SRM2294 #2	10.97	10.60
AccuStd Check	12.00	11.81

# Migration to H<sub>2</sub>:

## -- Specific Considerations for GC/MS






# H<sub>2</sub>: Chromatographic Method Migration

## Be aware:

- Consider flow limitation due to MSD pumping capacity
- Ensure >35 cm/sec flow rate (see Van Deemter Curves)
- Keeping similar peak elution order
- Consider column sample capacity

## *H<sub>2</sub> isn't a inert gas*

- Consider full inert flow path
- Use lowest temp possible
- Avoid methylene chloride, carbon disulfide as solvent



# H<sub>2</sub>: Additional HW modification on GC/MS

## Check Magnetic (5975 only)

- Ensure SN is printed on it
- Call CE if not



## Use H2 optimized draw out lens



- PN: G2589-20045



# H<sub>2</sub> for GC/MS: Initial Setup

## When start: customer should expect

- High background that looks like hydrocarbons ☹️
- Reduced signal to noise (worse MDL) ☹️
- Significant tailing for many compounds ☹️





# H<sub>2</sub> for GC/MS: Initial Setup

## Overnight system clean up:


- After setup, purging and pump down
- Set the source to max temp for your source
- Reduce the EMV to 800V
- Leave the FILAMENT ON overnight.
- System will be cleaned in the morning
- Make a run with matrix to check
- Ready to go



# H<sub>2</sub> for GC/MS: Analytical Result Expectations

- Sensitivity reduction: 2 – 5 times
- Trace conc. “reactive” compounds (e.g. flavor) may lost
- Possible spectrum infidelity
  - Observed ion ratio change for some compounds
  - Extra/missing ions, but often time not qualifier ions

## Lindane



	NIST Library Match	
	Helium	Hydrogen
Dichlorvos	924	867
Mevinphos	925	852
Ethalfuralin	927	897
Trifluralin	897	856
Atrazine	914	901
BHC-gamma (Lindane, gamma HCH)	925	858
Chlorpyrifos-methyl	922	897
Heptachlor	926	851
Malathion	916	836
Chlorpyrifos	901	890
DDE, p,p'-	934	903
Dieldrin	926	903
Hexazinone	901	832
Propargite	851	821
Leptophos	884	847
Mirex	903	881
Fenarimol	905	872
Coumaphos	864	804
Etofenprox (Ethofenprox)	906	842



# Summary: Helium Conservation Benefits

- **Seamless integration**

  - No need to revalidate existing GC methods

  - Fully integrated with 7890B and CDS (OpenLAB, Mustang, Mass Hunter)

  - Carrier gas ID and set points are a part of the method for compliance and transfer

  - Easily implemented using new Agilent Sleep/Wake functions

- **Greater reliability**

  - Based on proven 5<sup>th</sup> generation AUX EPC

  - 7890 provides warning if set points are not reached

  - For hydrogen users, nitrogen substitution when not running GC

- **Greater performance**

  - Purge channel prevents cross contamination of gases

  - Delivers more stable gas pressure control from the tank regulator to the inlet EPC module

  - Acts as an intermediate pressure regulator from the tank to inlet EPC to ensure greater analytical precision

# Summary – Migration To H<sub>2</sub> and N<sub>2</sub>

If you still need a helium alternative:

- For resolution critical methods, H<sub>2</sub> offer the best alternative
  - Agilent GC and GC/MS systems have many built-in safety features
- For many GC applications, N<sub>2</sub> offers a cheap, easy alternative without any safety worries
  - Many existing helium methods have too much resolution
  - N<sub>2</sub> can be used without changing any of the existing GC conditions
    - keep the holdup time the same as the original method
  - 2-D methods have high resolution built-in, so N<sub>2</sub> is ideally suited as a carrier gas
    - Valve-based or Deans switch, not GCxGC
- For more information on Helium Carrier Gas  
[www.agilent.com/chem/heliumupdate](http://www.agilent.com/chem/heliumupdate)