

# Updating Old Methods: Is the Gain Worth the Pain?

Becoming a Better Chromatographer  
HPLC educational webinar

Mark Powell  
Applications Engineer  
Columns and Supplies Technical Support



# Why Would You Consider Updating? What Is the Goal?

- Do I need to update?
- Reduce analysis time
  - Productivity (free up time for other tasks)
  - More analyses (free up instrument time)
  - More samples (higher throughput)
- Improve resolution
- Improve sensitivity
- Increase column life
- Save solvent
  - Cost
  - Disposal




# What Can Be Changed? Summary of Allowable Adjustments: USP General Chapter <621>


Parameters for System Suitability		
	Isocratic	Gradient
Particle Size	L/dp: -25% to +50% or N: -25% to +50%	No Changes allowed
Column Length		
Column Inner Diameter	Flexible, w/ constant linear velocity	No Changes allowed
Flow Rate	Based on dp: $F_2 = F_1 \times [(dc_2^2 \times dp_1) / (dc_1^2 \times dp_2)]$ Additional adjustments: $\pm 50\%$ , provided N decreases $\leq 20\%$	No Changes allowed
Injection volume	May be adjusted, as far as is consistent with precision and detection limits	May be adjusted, as far as is consistent with precision and detection limits
Column Temperature	$\pm 10^\circ\text{C}$	$\pm 10^\circ\text{C}$
Mobile phase pH	$\pm 0.2$ units	$\pm 0.2$ units
Salt Concentration	within $\pm 10\%$ if the permitted pH variation is met	within $\pm 10\%$ if the permitted pH variation is met
Ratio of Components in Mobile Phase	Minor component ( $\leq 50\%$ ): $\pm 30\%$ relative, but cannot exceed $\pm 10\%$ absolute; may only adjust 1 minor component in ternary mixtures	No Changes allowed *  * Not specified in <621>, assume no changes are allowed
Wavelength of UV-Visible Detector	No changes allowed	No changes allowed

# What Are the Column Options?

- Smaller particle size
  - Higher efficiency -> shorter column -> faster method
  - Increase resolution
  - Better sensitivity
  - Consider pressure limit of instrument
- Smaller diameter
  - Solvent savings
  - Depends on instrument configuration and plumbing
- Bonded phase
  - Match USP designation
  - More robust column life
  - Consider a different bonded phase?




# Agilent InfinityLab Poroshell 120 Particle Sizes




InfinityLab Poroshell 120  
**1.9 µm**

Highest UHPLC  
performance



InfinityLab Poroshell 120  
**2.7 µm**


UHPLC performance at  
lower pressure




InfinityLab Poroshell 120  
**4 µm**

Improved HPLC  
performance

# Particle Size: When to Use What Size



# Decreasing Particle Size Increases Efficiency



## Columns:

Agilent InfinityLab Poroshell 120 EC-C18, 2.1 x 50 mm, 1.9 um

Agilent InfinityLab Poroshell 120 EC-C18, 2.1 x 50 mm, 2.7 um

Agilent InfinityLab Poroshell 120 EC-C18, 2.1 x 50 mm, 4 um

**Mobile phase A:** 0.2% formic acid in water

**Mobile phase B:** Acetonitrile

**Gradient:** 5-16% B in 7 min

**Flow rate:** 0.5 mL/min

**Detection:** 240 nm @ 80 Hz

**Sample:** 1 μL of 0.06 mg/mL each of gallic acid, gallic acid gallate, gallo catechin, gallo catechin gallate, epigallocatechin, epigallocatechin gallate, catechin, catechin gallate, caffeine, caffeine gallate, epicatechin, epicatechin gallate, epigallocatechin gallate, gallo catechin gallate, epicatechin gallate, catechin gallate


- Higher N improves resolution as particle size is decreased

$$N \propto \frac{L}{d_p}$$


Particle	Pressure	R <sub>smin</sub>
1.9 μm	226 bar	2.2
2.7 μm	131 bar	1.3
4 μm	53 bar	0.7

# Agilent InfinityLab Poroshell 120 phases


Best all around	Best for low pH mobile phases	Best for high and mid pH mobile phases	Best for alternative selectivity	Best for more polar compounds	HILIC for polar compounds	Chiral phases
Poroshell 120 <b>EC-C18 (L1)</b> 1.9 µm, 2.7 µm, 4 µm	Poroshell 120 <b>SB-C18 (L1)</b> 2.7 µm	Poroshell 120 <b>HPH-C18 (L1)</b> 1.9 µm, 2.7 µm, 4 µm	Poroshell 120 <b>Phenyl-Hexyl (L11)</b> 1.9 µm, 2.7 µm, 4 µm	Poroshell 120 <b>SB-Aq (L96)</b> 2.7 µm	Poroshell 120 <b>HILIC (L43)</b> 1.9 µm, 2.7 µm, 4 µm	Poroshell 120 <b>Chiral-CF</b> 2.7 µm
Poroshell 120 <b>EC-C8 (L7)</b> 1.9 µm, 2.7 µm, 4 µm	Poroshell 120 <b>SB-C8 (L7)</b> 2.7 µm	Poroshell 120 <b>HPH-C8 (L7)</b> 2.7 µm, 4 µm	Poroshell 120 <b>Bonus-RP (L60)</b> 1.9 µm, 2.7 µm, 4 µm	Poroshell 120 <b>EC-CN (L10)</b> 2.7 µm	Poroshell 120 <b>HILIC-Z</b> 2.7 µm	Poroshell 120 <b>Chiral-CD</b> 2.7 µm
			Poroshell 120 <b>PFP (L43)</b> 2.7 µm		Poroshell 120 <b>HILIC-OH5</b> 2.7 µm	Poroshell 120 <b>Chiral-V</b> 2.7 µm
						Poroshell 120 <b>Chiral-T</b> 2.7 µm



4 µm



2.7 µm



1.9 µm



# Benefits of Transferring to Agilent InfinityLab 120 Poroshell


- Superficially porous particles
  - Can run at high speed without sacrificing resolution
  - Existing instruments can turn over more samples per hour
  - 30-50% increase in efficiency near the backpressure of a 3.5-5.0  $\mu\text{m}$  particle
  - Boosts the sample throughput on low-pressure instruments
- Short column lengths (or narrow column diameters)
  - Cut solvent flow by 50-80%
  - Reducing solvent consumption
  - Cuts purchase and disposal costs



# Van Deemter Equation

$$h = L/N$$

$$h = A + B/u + C \cdot u$$



- **A term**: eddy diffusion and flow distribution
  - Particle size & packing quality important
  - Narrow particle size distribution
- **B term**: longitudinal diffusion
  - Diffusion in the mobile phase
- **C term**: mass transfer
  - shorter diffusion paths
  - better with superficially porous particles
  - more effect on large molecules
- **u**: linear velocity
  - velocity of mobile phase through column
  - $u = L/t_0$  in cm/sec

Lower  $h$  (reduced plate height) = higher efficiency


# Van Deemter Equation

## Eddy diffusion


$$W_{eddy} \sim \lambda d_p$$

$\lambda$ : Quality of column packing


Differences in diffusion paths due to:



Different paths



Poor column packing



Broad particle size distribution


# Van Deemter Equation

## Axial or longitudinal diffusion

Increase in peak width due to self-diffusion of the analyte

At low flow the analyte remains in the mobile phase for a long time


- High increase in peak width
- Increased height of a theoretical plate



# Van Deemter Equation


Resistance to mass transfer

$$W_C \sim d_p^2$$




Totally porous particle (TPP)

Different diffusion paths



Superficially porous particle (SPP)

# Van Deemter Curves






Optimize your flow rate for Agilent InfinityLab Poroshell 120:

- For 2.1 mm id, we suggest 0.42 mL/min
- For 3.0 mm id, we suggest 0.85 mL/min
- For 4.6 mm id, we suggest 2 mL/min

Overlay of van Deemter plots: the optimal flow rate for Agilent InfinityLab Poroshell 120 is faster than for 5 or 3.5 μm columns

# Faster Ibuprofen Analysis on Agilent InfinityLab Poroshell 120 EC-C8 4.6 x 50 mm, 2.7 $\mu\text{m}$



“Transfer and Optimization of HPLC Methods to Superficially Porous UHPLC” Pittcon 2012, W. Long, A Brooks Oral Presentation 420-3.

# Transfer of Isocratic USP Methods (Naproxen)




Parameters for System Suitability		
	Isocratic	Gradient
Particle Size	L/dp: -25% to +50% or N: -25% to +50%	No Changes allowed
Column Length		

L (mm)	dp (µm)	L/dp	%	N	%	<621> compliant
150	5	30,000	100%	10,639	100%	Yes
150	4	37,500	125%	19,054	179%	Yes
100	4	25,000	83%	13,186	124%	Yes
100	2.7	37,037	123%	21,046	198%	Yes
50	2.7	18,519	62%	11,281	106%	Yes



# Scaling USP Naproxen Method from a 5 μm TPP to 4 μm SPP

System Suitability Method Requirement:  $N > 4000$ ,  $R_s > 11.5$






Mobile Phase: 50:49:1 MeCN:H<sub>2</sub>O Acetic Acid

Flow Rate: 1.2 mL/min

Peak 1. Naproxen 2. Butyrophenone

# Scaling USP Naproxen Method from a 5 $\mu\text{m}$ TPP to 2.7 $\mu\text{m}$ SPP

System Suitability Method Requirement:  $N > 4000$ ,  $R_s > 11.5$



Mobile phase: 50:49:1 MeCN:H<sub>2</sub>O Acetic Acid  
Flow rate: 1.2 mL/min, **All Pressures < 300 bar**  
Peak 1. Naproxen 2. Butyrophenone

# Retention Factor - Gradients

$$k^* = \frac{t_G \cdot F}{S \cdot \Delta\Phi \cdot V_m}$$


With gradient separations, the retention factor is influenced by

- $F$  = flow rate
- $t_G$  = gradient time (minutes)
- $\Delta\Phi$  = change in volume fraction of B mobile phase
- $V_m$  = column volume
- $S$  = constant (4 - 6 for small molecules, 10 - 1000 for peptides and proteins)

To keep the retention factor constant, changes in the denominator need to be offset by proportional changes in the numerator, and vice versa.

# Retention Factor - Gradients

$$k^* = \frac{t_G \cdot F}{S \cdot \Delta\Phi \cdot V_m}$$



Increasing the solvent strength  
= Increasing the % organic in the  
mobile phase

Linear solvent strength gradient  
= % per min is a constant

$$\Delta\Phi = 80\%$$

$$t_G = 40 \text{ min}$$

$$\frac{\Delta\Phi}{\Delta t_G} = 2\%/min$$

# Maintaining $k^*$

Keep relative peak position and shorten analysis

## Any decrease in

- Column length
- Column volume (i.d.)
- $\Delta\%B$  (same column)



## Can be offset by a proportional

- Decrease in  $t_G$  or  $F$
- Increase in  $\Delta\%B$
- Decrease in  $t_G$  or  $F$
- Increase in  $\Delta\%B$
- Decrease in  $t_G$  or  $F$



$$k^* \propto \frac{t_G \cdot F}{S \cdot \Delta\Phi \cdot V_m}$$


# Reduce Analysis Time, Keep Gradient Steepness the Same

Sample: 1. Aldicarb sulfoxide, 2. Oxamyl, 3. Methomyl, 4. Aldicarb sulfone, 5. Carbofuran-3-hydroxy, 6. Aldicarb, 7. Propoxur, 8. Carbofuran, 9. Carbaryl, 10. Methiocarb, 11. ISTD (BDMC)

Column: Agilent ZORBAX Eclipse Plus-C18  
4.6 x 150 mm, 5  $\mu$ m

Gradient Time: 20 min  
Flow Rate: 1.0 mL/min


Analysis  
Time: 14 min




Column: Agilent InfinityLab Poroshell  
120 EC-C18 4.6 x 50 mm, 2.7  $\mu$ m

Gradient Time: 6.7 min  
Flow Rate: 1.0 mL/min

Analysis  
Time: 5 min



# Gradient Method Transfer– 4.6 x 150mm, 5 μm Agilent ZORBAX SB-C18




Flow Rate 1.0 ml/min  
 Injection Volume 15uL  
 Temperature 30° C  
 Wavelength 246nm  
 Sample rate 2.5 Hz

Time (min)	% Acetonitrile
0	50
10	90
13.5	90
13.6	50
15	50

# Maintaining Peak Position & Resolution

Have shortened column & gradient time – need to do so by the SAME factor  
1/3 column length – 1/3 gradient time  
ex: RRHT column – 4.6 x **50 mm**, 1.8 $\mu$ m, SB-C18



Flow Rate 1.0 ml/min  
Injection Volume 5 $\mu$ L  
Temperature 30° C  
Wavelength 246nm  
Sample rate **13.74 Hz**

Time (min)	% Acetonitrile
0	50
<b>3.33</b>	90
<b>4.5</b>	90
<b>4.53</b>	50
<b>5</b>	50



# Gradient Transfer: 4.6 x 250 mm, 5 µm to 4.6 x 100 mm, 2.7µm


Mobile phase:

A: 0.1% formic acid in water

B: 0.1% formic acid in ACN

Time	%B
0	8
33	33
34	33


Column: **Agilent Eclipse Plus C1860**  
**4.6 x 250mm, 5 µm**  
Flow rate: 1 mL/min



Sulfadiazine,  
Sulfathiazole  
Sulfapyridine  
Sulfamerazine,  
Sulfamethazine,  
Sulfamethazole,  
Sulfamethoxypyridazine,  
Sulfachloropyridazine  
Sulfamethoxazole,  
Sulfadimethoxine

Time	%B
0	8
12	33
13.2	33

Column: **4.6 x 100mm**  
**Agilent InfinityLab Poroshell**  
**120 EC-C18, 2.7 µm**  
Flow rate: 1 mL/min




- Faster analysis
- More sensitive (smaller particle sharpens peak)
- Both run under 400 bar

- Column length decreased from 250 to 100 mm ( $250 / 100 = 0.4$ )
- Gradient time points decreased proportionally:  $33 * 0.4 = 13.2$  (adjusted to ~12 minutes)

# Resolution

## Baseline separations for rugged methods





If we consider peaks of equal height:

- 1 - minimum for a measureable separation
- 0.6 - required to discern a valley between two equal-height peaks
- 1.5 - considered to be a baseline separation
- 1.7 or greater - desirable for rugged methods

# Notoginseng Analysis: 4 $\mu\text{m}$ Agilent Poroshell and 5 $\mu\text{m}$ Totally Porous

Same dimensions, Same conditions, higher resolution



Mobile phase:


A) water

B) acetonitrile

Gradient for 4.6 x 150 mm:


Time (min)	%A
0	81
12	81
60	64
61	10
65	10
66	81
70	81

# Notoginseng Analysis: Speed Optimization



Gradient for 4.6 x 150 mm:

Time (min)	%A
0	81
12	81
60	64
61	10
65	10
66	81
70	81





Gradient for 4.6 x 100 mm:

Time (min)	%A
0	81
8	81
40	64
40.5	10
43.5	10
44	81
47	81

- Column length decreased from 150 to 100 mm ( $150 / 100 = 0.67$ )
- Gradient time points decreased proportionally
- $12 * 0.67 = 8$  minutes
- $60 * 0.67 = 40$  minutes

# Notoginseng Analysis: Speed Optimization

- Flow rate increased from 1 mL/min to 1.5 mL/min ( $1.5 / 1.0 = 0.67$ )
- Gradient time points decreased proportionally
- $45 * 0.67 = 30$  minutes





Gradient for 4.6 x 100 mm:

Time (min)	%A
0	81
8	81
40	64
40.5	10
43.5	10
44	81
47	81

# Achieving Better Sensitivity

Totally porous 3.5  $\mu\text{m}$  to 2.7  $\mu\text{m}$  superficially porous





Same method can be used with either column due to similar selectivity

2.7  $\mu\text{m}$  columns:

- Taller, narrower peaks
- >40% more sensitivity, as noted by the S/N of ibuprofen
- Conditional peak capacity >20% higher than the 3.5  $\mu\text{m}$  column
- 15 compounds, 2 min

# Amino Acid Analysis on Agilent InfinityLab Poroshell 120 HPH-C18



Column: Agilent InfinityLab Poroshell HPH-C18 or Agilent ZORBAX Eclipse Plus C18

Column Temperature: 40 °C





Mobile Phase A: 10 mM Na<sub>2</sub>HPO<sub>4</sub>:  
10 mM Na<sub>2</sub>B<sub>4</sub>O<sub>7</sub>, pH 8.2: 5 mM NaN<sub>3</sub>

Mobile Phase B: Acetonitrile:  
Methanol: Water (45:45:10, v: v: v)

Injection Diluent: (0.25 mL H<sub>3</sub>PO<sub>4</sub> +  
100 mL H<sub>2</sub>O)

- Poroshell HPH uses a special coating process
- Particle resists attack at high pH
- Extends column lifetime in alkaline mobile phase

# Amino Acid Analysis on HPH-C18 in pH 8.3 Phosphate Buffer






# Competitor High pH column vs Agilent InfinityLab Poroshell 120 HPH-C18: Stress Test


A: 10 mM ammonium bicarbonate, pH 10

B: acetonitrile


5-95%B over 5 minutes




Other column




Agilent InfinityLab  
Poroshell  
120 HPH-C18  
2.1 x 50 mm, 2.7  $\mu$ m



**Still maintains resolution after 2000 injections!**




# Solvent Use




4.6 mm

1.5 mL/min



3.0 mm

0.62 mL/min




2.1 mm

0.21 mL/min

Amino Acid Analysis on Agilent InfinityLab Poroshell 120, 100 mm, 2.7 $\mu$ m HPH-C18

# System Considerations




**Dwell Volume:** from formation of gradient to top of column

-minimize for faster equilibration and more efficient gradient formation

**Extra Column Volume** from injection to detector (flow cell) outside of the column

-minimize to reduce band broadening, for sharper peaks and better resolution

# Paraquat and Diquat on ZORBAX Eclipse XDB-C8 Ion Pairing Conditions



Conditions:

Column: ZORBAX Eclipse XDB-C18  
4.6 x 150 mm, 5 mm

Mobile Phase: 70% water + heptanesulfonic acid  
phosphoric acid to pH = 3.0 :30% methanol  
Temperature: 30°C

Flow Rate: 1 mL/min

Detection: UV 257 nm (Paraquat)  
UV 308 nm (Diquat)

Sample: 50 ng on column  
1. Paraquat  
2. Diquat

## Ion Pairing methods

- Advantage: Uses standard system and C18 columns
- Disadvantage: Permanently contaminates system, lengthy method development, restricted to positive or negative mode MS

# Agilent InfinityLab Poroshell 120 HILIC-Z Approach

- **No loss in throughput:** InfinityLab Poroshell HILIC columns operates as quickly and reliably as reverse phase
- **Easy to implement:** Uses the same solvents and buffers as reverse phase
- Positive or Negative ESI


## InfinityLab Poroshell HILIC-Z 2.1 x 100 2.7 $\mu$ m

Mobile phase A: 20 mM ammonium formate, pH 3

Mobile phase B: MP A in 90% acetonitrile


Min	%B
0	80
1	80
5	73
6	80
8	80

ESI Positive



# Updating Chiral Separations

- 4 chiral stationary phases on 2.7µm Agilent Poroshell particles
- Fast, efficient analysis using superficially porous particles (<5 min!)
- High efficiency provides superior peak shape and resolution compared to totally porous chiral phases
- Use common reversed-phase solvents for maximum method flexibility - even LCMS



Column: 4.6x100mm, 2.7µm  
Mobile phase: 100/0.2/0.05: Methanol/Acetic Acid/Ammonium Hydroxide  
Flow Rate: 1.0 mL/min Detection: UV 230 nm

*Method guidance available:  
Chiral Application Notebook  
[5991-8450EN](#)*

## Best for Chiral

InfinityLab Poroshell  
**Chiral-V**  
2.7 µm

InfinityLab Poroshell  
**Chiral-T**  
2.7 µm


InfinityLab Poroshell  
**Chiral-CD**  
2.7 µm

InfinityLab Poroshell  
**Chiral-CF**  
2.7 µm

# Fast, High Efficiency Chiral Separations

## Traditional Chiral Separation— totally porous particle


Chirobiotic V2 (250 x 4.6 mm, 5 µm)



1. D-(+)-Amphetamine, 2. L-(-)-Amphetamine, 3. D-(+)-Methamphetamine  
4. L-(-)-Methamphetamine 100/0.1/0.02 MeOH/HOAc/NH<sub>4</sub>OH with a  
1.0 mL/min flow rate at room temperature and UV at 220 nm

## Agilent InfinityLab Poroshell 120 Chiral Separation— superficially porous particle

InfinityLab Poroshell 120 Chiral-V (100 x 4.6 mm, 2.7 µm)



1. D-(+)-Amphetamine, 2. L-(-)-Amphetamine, 3. D-(+)-Methamphetamine  
4. L-(-)-Methamphetamine 100/0.1/0.02 MeOH/HOAc/NH<sub>4</sub>OH with a  
1.0 mL/min flow rate at room temperature and UV at 220 nm

# Tips for Method Transfer

Adjusting flow for different column diameters

$$Flow_{column\ 1} \times \left( \frac{diameter_{column\ 2}}{diameter_{column\ 1}} \right)^2 = Flow_{column\ 2}$$

$$1\ mL/min \times \left( \frac{2.1\ mm}{4.6\ mm} \right)^2 = 0.21\ mL/min$$



# Tips for Method Transfer

## Adjusting injection volume

$$V_m = \pi \times r^2 \times L \times \sim 0.6$$

$$Injection_{column\ 1} \times \left( \frac{Volume_{column\ 2}}{Volume_{column\ 1}} \right) = Injection_{column\ 2}$$

Original 4.6 x 250 mm: ~2.5 mL

Transferred to 2.1 x 100 mm: ~0.21 mL

$$30\ \mu L \times \left( \frac{0.21\ mL}{2.5\ mL} \right) = 2.5\ \mu L$$

# Agilent LC Column Navigator

<http://navigator.chem.agilent.com/>




**LC COLUMN  
and SAMPLE PREP  
NAVIGATOR**




Let us help you find the best Agilent column and sample prep products for your application.

**BEGIN ▶**



**WHAT CAN WE HELP YOU WITH?**  
Choose from these options to narrow your search.

- Option 1: I already have a column or sample prep method.** ▶
- Option 2: Column recommendation based on method parameters. ▶
- Option 3: I would like to see compound-specific references. ▶
- Option 4: USP method.** ▶



Request your Poroshell Selection Poster

# Contact Agilent Chemistries and Supplies Technical Support



1-800-227-9770 Option 3, Option 3:

Option 1 for GC and GC/MS columns and supplies

Option 2 for LC and LC/MS columns and supplies

Option 3 for Sample preparation, filtration and QuEChERS

Option 4 for Spectroscopy supplies

**Available in the USA & Canada 8-5 all time zones**

[gc-column-support@agilent.com](mailto:gc-column-support@agilent.com)

[lc-column-support@agilent.com](mailto:lc-column-support@agilent.com)

[spp-support@agilent.com](mailto:spp-support@agilent.com)

[spectro-supplies-support@agilent.com](mailto:spectro-supplies-support@agilent.com)