



Mestrelab Research

Chrom Reaction Optimization 1.1

STARTING GUIDE



Document Number
P/N 312 R1.1



COPYRIGHT

©2023 MESTRELAB RESEARCH S.L.

All rights reserved. No parts of this work may be reproduced in any form or by any means - graphic, electronic, or mechanical, including photocopying, recording, taping, or information storage and retrieval systems - without the written permission of the publisher.

Products that are referred to in this document may be either trademarks and/or registered trademarks of the respective owners. The publisher and the author make no claim to these trademarks.

While every precaution has been taken in the preparation of this document, the publisher and the author assume no responsibility for errors or omissions, or for damages resulting from the use of information contained in this document or from the use of programs and source code that may accompany it. In no event shall the publisher and the author be liable for any loss of profit or any other commercial damage caused or alleged to have been caused directly or indirectly by this document.

Table of Contents

TABLE OF CONTENTS	2
1. INSTALLATION	3
2. PRE-REQUISITE (WHEN ANALYZING POSITIVE AND NEGATIVE MASS FUNCTIONS)	5
3. THE WORKFLOW	6
3.1. INPUT	6
3.2. PROCESSING.....	10
3.3. PLUGINS.....	10
3.3.1. <i>The Input tab</i>	11
3.3.2. <i>The Analysis tab</i>	16
3.3.3. <i>The Quality Controls tab</i>	19
3.3.4. <i>The Output tab</i>	20
3.4. DESIGN	27
3.5. OUTPUT	27
4. THE OUTPUT FOLDER	28
4.1. THE HTML FILE	29
4.2. THE CSV REPORT	30
4.3. OTHER OUTPUT.....	30
5. THE MNOVA GEARS RESULTS VIEWER	30
5.1. THE MGEARS VIEWER	30
5.1.1. <i>The Well plate view</i>	32
5.1.2. <i>The Results section</i>	33
5.2. THE CHROM REACTION OPTIMIZATION VIEWER	36

A chemist sometimes needs to optimize a reaction. This typically involves running this reaction numerous times whilst systematically varying physical conditions or reagents in order to find the best combination in which the product yield is maximized, or the level of some undesired by-product is minimized.

Parallel reaction technology and high throughput LC-MS analysis can allow the chemist to run many such combinations; however, this approach can also easily lead to new bottlenecks, in particular, the analysis of the resulting large body of data.


Chrom Reaction Optimization is an Mgears plugin that helps solve this problem by quickly processing large chromatography, GC/LC-MS datasets and generating automatic reports for further evaluation and decision making. Chrom Reaction Optimization tracks a defined set of chemicals across a large number of samples and that reports on their levels.

This document is to help you get started with Chrom Reaction optimization.


1. Installation

Before installing Chrom Reaction Optimization, make sure you already have Mnova MSChrom (minimum version: 15.0) and Mgears (minimum version: 2.5) installed and running with valid licenses.


Go to **Files>Advanced Plug-ins>Available**. Tick the Chrom Reaction Optimization plugin, then press **Install**.




Another option is to drag and drop the Chrom Reaction Optimization installer into the Mnova interface. The following dialog will open. Click on **Install**.



Restart Mnova.




Chrom Reaction Optimization must now be installed. You can check your license status by going to **Files>Help>License Manager**. A green check must appear in the plugin's status column.




You are now good to go.

2. Pre-requisite (when analyzing positive and negative mass functions)

When you wish to analyze positive and negative mass functions together, it is crucial to ensure that your Mnova Mass preferences are appropriately configured to detect both positive and negative ionization chromatograms. To do so, go to **Files>Preferences>Mass** (1), navigate to the **Setup** tab (2), locate and click on the **Add** button  to incorporate your mass functions (3). Select the **Positive** function (4) then confirm your choice by pressing **OK** (5).




Repeat the operation and choose now the **Negative** function. Both “+ MS TIC” and “- MS TIC” should now appear in the setup.



Click **OK** to save this configuration. You can proceed now with setting up your reaction optimization analysis.

3. The workflow


Launch Mgears from the Mnova **Automation** ribbon. The dialog with the usual six tabs will open.



3.1. Input


Chrom Reaction Optimization (Chrom RO) analysis follows the general Mnova Gears workflow based on its five-step setup, the first step being the choice of input data. Mgears can read data from your **Disk**, **Database**, or from **Real-Time** acquisition. In this guide, we will work with data from disk directories (*please refer to the [Mnova Gears manual](#) for more details about other input types*).

Click on the button and select the data files to upload as your **Main Directory**.




Use the **Automatic Detection** mode to detect your input files.


Experiment selection is recommended if your data folder contains different types of data file in order to restrict detection to GC/LC/MS and avoid analysis of other undesired file types. In this case, you must select the **GC/LC/MS** experiment type as shown below.




Important! When using **Automatic Detection**, make sure your data file extension(s) is included in the **Mask Manager** so that Mgears can recognize the appropriate files correctly. Otherwise, you can manually add your standard file extension or custom file formats and press **OK**.




Using a **Mapping file** is convenient (but optional) when input files are listed in a **.txt** or **.csv** document along with other metadata. Mgears will map the information found in the file provided with the information found in the selected input directory (on Disk, DB, or RT folder). To add and use a **Mapping file**, click on the **Filter** icon at the bottom left-hand side of the Mgears dialog.



Click on **...** and upload your **.txt** or **.csv** document with the data and parameters you want to work with.




Now, click on to set up the parameters to read the mapping file. The **Advanced Filtering Options** dialog will open.




In the **Advanced Filtering Options** dialog, you must set the column with the data file names for matching. You can either indicate the number of the column directly by clicking on this box , or open the assistant to visualize the CSV and select the corresponding column, as shown in the image below. Click on **OK** to save the selection.

Note. In the **CSV column Detection** dialog, only the first 10 rows of the filtering file will be loaded.



Important! Remember that in order to be correctly read by Mgears, the CSV separator (comma, semicolon, tab, or vertical line) should be configured in **Settings** tab (in this case, the CSV is using a comma as a separator, therefore we select the comma in the **Settings**).



Other filtering options can be configured for the mapping file. You can, for example, set the **Starting Row** at which matching should start, decide whether to **Read Until** a specific row, or **Allow Gaps** in the results.

The **well plate position** (A1, A3...) can also be matched.

Read Column with the Position in the Well Plate
Column: 5

Read Column with ID of Compound
Column: 1

Read Column with SMILES Compound
Column: 13

OK Cancel

Once you've completed configuring your filtering options, click on **OK**.

You can now perform an automatic inspection of the selected directory by clicking on to check if your input files have been found and filtered correctly.

Mnova Gears

Input Processing Plugins Design Output Settings

Using Automatic Detection Experiments Selection

Main Directory: C:/Users/Usuario/Desktop/Datasets/By Plugin/Chrom RO/Dataset/Input

Optional Directory: ...

Experiments Found: 10

List of detected experiments:
*Result filtered using file: C:/Users/Usuario/Desktop/Datasets/By Plugin/Chrom RO/Dataset/Mapping file.csv(Using Column: 1 Delimiter: ,)

Experiment: REACTION-01
MS: C:/Users/Usuario/Desktop/Datasets/By Plugin/Chrom RO/Dataset/Input/Reaction-01.raw/_FUNC001.DAT
Position: A1

Experiment: REACTION-02
MS: C:/Users/Usuario/Desktop/Datasets/By Plugin/Chrom RO/Dataset/Input/Reaction-02.raw/_FUNC001.DAT
Position: A2

Experiment: REACTION-03
MS: C:/Users/Usuario/Desktop/Datasets/By Plugin/Chrom RO/Dataset/Input/Reaction-03.raw/_FUNC001.DAT
Position: A3

Experiment: REACTION-04
MS: C:/Users/Usuario/Desktop/Datasets/By Plugin/Chrom RO/Dataset/Input/Reaction-04.raw/_FUNC001.DAT
Position: A4

Experiment: REACTION-05
MS: C:/Users/Usuario/Desktop/Datasets/By Plugin/Chrom RO/Dataset/Input/Reaction-05.raw/_FUNC001.DAT
Position: A5

Experiment: REACTION-06
MS: C:/Users/Usuario/Desktop/Datasets/By Plugin/Chrom RO/Dataset/Input/Reaction-06.raw/_FUNC001.DAT
Position: A6

Experiment: REACTION-07
MS: C:/Users/Usuario/Desktop/Datasets/By Plugin/Chrom RO/Dataset/Input/Reaction-07.raw/_FUNC001.DAT
Position: A7

Experiment: REACTION-08
MS: C:/Users/Usuario/Desktop/Datasets/By Plugin/Chrom RO/Dataset/Input/Reaction-08.raw/_FUNC001.DAT
Position: A8

Experiment: REACTION-09
MS: C:/Users/Usuario/Desktop/Datasets/By Plugin/Chrom RO/Dataset/Input/Reaction-09.raw/_FUNC001.DAT
Position: A9

Experiment: REACTION-10
MS: C:/Users/Usuario/Desktop/Datasets/By Plugin/Chrom RO/Dataset/Input/Reaction-10.raw/_FUNC001.DAT
Position: A10


OK

Input folder

- Reaction-01.raw
- Reaction-02.raw
- Reaction-03.raw
- Reaction-04.raw
- Reaction-05.raw
- Reaction-06.raw
- Reaction-07.raw
- Reaction-08.raw
- Reaction-09.raw
- Reaction-10.raw

3.2. Processing


In the **Processing** tab, you can optionally upload a script to apply customized processing options.




In this analysis, we will not be setting up any **Processing** options.

3.3. Plugins

In the **Plugins** section, select and add the Chrom Reaction Optimization plugin.



Then, click on **Chrom Reaction Optimization Plugin Settings** to configure your analysis and report parameters.




A dialog with four tabs should appear.


3.3.1. The Input tab

In the **Input** tab, you must declare the chemical reaction components that need to be monitored. For each component, you must provide at least one of the following: a **Mol.file**, **SMILES** code, **molecular formula**, **MZ**, or a **Retention Time (RT)**. This information will be used to identify and assign chromatogram peaks to the components, so the choice of the information provided in this regard is critical to an accurate assignment.


Click on to add a reaction component.



Note. The first time you click on the **Add** button you will get a message informing you that changes in the reaction components may affect the analysis formulae used in the evaluation and the output configuration. Make sure to review those and amend them according to your changes before launching the analysis.




A new line is added to the table. Fill in the **Name** and **Type** of the compound you wish to monitor. The compound can be a reaction start material, product, by-product, or an internal standard that will serve for relative quantification. When a compound is known but undesired or uninteresting to the analysis, it is possible to set it as “ignored”, and in this way the peaks from that compound will be excluded from the analysis.




Note. Duplicate component names are not allowed within this system. Whenever you attempt to input the same name twice manually, a warning message will promptly appear to alert you. Likewise, if you happen to be using older settings that contain duplicated component names, you will receive a warning message when loading these settings into Mgears. This message will inform you that the duplicated name has been replaced by another to ensure the uniqueness of component names within the system.

When adding a new component




When loading old settings



Complete the table row by entering the information that will be used for peak assignment and quantification.

To add a Mol.file, simply double-click on the cell. A dialog will open, allowing you to browse your files and select the Mol.file. Once you've selected your file, click **OK**, and the file path will appear in the cell.

To delete a previously added file, you can click on **X**, then press **OK**.



To add a **SMILES** code, **molecular formula**, **MZ**, or a **Retention Time (RT)**, simply type your value inside the corresponding cell.

When providing a **Mol.File**, **SMILES**, or **Molecular Formula**, the component is matched using the Molecule Match feature. When providing an **MZ**, the component is matched using the main peak of the auxiliary EIC for that **MZ**.

In cases where a particularly stringent assignment method is needed, you can provide both an **RT range** and an **MZ** or **molecular formula** for a component. In this scenario, the **RT** information takes precedence, and matching **MZ** peaks must fall within the specified **RT range** to be assigned. This approach helps eliminate false-positive **MZ** matches.

Detection method considerations:

For detection and peak assignment **by mass** (i.e., when a structure, a SMILES code, a molecular formula or an **MZ** value is provided for matching), you will need to:

- Make sure that the component is ionizing properly. If ionization is poor, it is recommended to set detection by **RT**.
- Use the **MZ** of the strongest ionization. In some cases, the strongest ionization is a fragmentation peak (e.g., for Boc-protected compounds, $M\text{-Boc} = M\text{-}100$).
- Check that the selected **MZ** for a component is “exclusive” (i.e., that no other compound is likely to give the same fragmentation peak).

Otherwise, you can define a **RT** or an **RT range** for detection and peak assignment.

A **Color** can be attributed to each component. This color will be used in the results file to highlight the peaks in the spectra, and to designate the compound in the graphics showing the composition of the samples (e.g., Pie and Bar Charts).

Input Analysis Quality Controls Output									
Reaction Components									
	Name	Type	Mol. File	Smiles	Mol. Formula	MZ	RT	Color	External Calibration
1	SM	Start Material				261.38			


When requiring an absolute quantification, it is possible to add an external calibration output file or a calibration value to be taken into account in calculations. To do so, double click on the cell.

Input Analysis Quality Controls Output									
Reaction Components									
	Name	Type	Mol. File	Smiles	Mol. Formula	MZ	RT	Color	External Calibration
1	SM	Start Material				261.38			<input type="text"/>

Double-click


A new **Calibration** dialog will open. Here, you can:

- Select and add a [Chrom Calibration](#) output file (a JSON file): calibration information will be imported and used in the current analysis.




CSV Input/Output Controls									
	Mol. Formula	MZ	RT	Color	External Calibration				
		261.38			C:/Users/Usuario/Desktop/Results/MSCalibration/Compound_148_calibration.json				

- Manually enter a calibration value. This value is the gradient of a Peak Area vs. Concentration plot. When using this option, the calculation of the compound concentration is achieved by dividing the Peak Area by the gradient.



To remove a previously added calibration file or value, check the **None** option, then **OK**. You can also check the **Use for all components** option if you want to delete all calibration information at once.



Complete your components table by repeating all the previous steps.

You can always delete unwanted components from your table by selecting the row and clicking on , or you can otherwise edit the various columns in the table as needed.

Reaction Components									
Name	Type	Mol. File	Smiles	Mol. Formula	MZ	RT	Color	External Calibration	
1 SM	Start Material				261.38		Blue		
2 P1	Product				296.2		Green		
3 P2	By Product				255.04		Red		
4 P3	Ignored					1.2870	Grey		

Minus icon highlighted in the bottom left corner of the table.

In the CSV section, the path to the CSV mapping file will be displayed if one has been uploaded in the [Mgears Input](#) tab. If no file has been uploaded, this section says “No CSV Input”.


3.3.2. The Analysis tab

3.3.2.1. The LC/MS Evaluation tab


In the LC/MS Evaluation tab, the initial step involves configuring the chromatogram extraction preferences.

To get started, click on the **Add** button , which will open a dialog box containing the **Chromatogram Parameters**.

Select the analytical trace type (TIC, DAD, DAD range, PDA, or other relevant trace type) for chromatogram extraction. Then, decide whether you want to perform chromatogram alignment automatically or manually. We strongly recommend opting for manual alignment. In the case of manual alignment, you will introduce a specific value, which will be applied to the chromatogram to align it with TIC. (Please refer to this [article](#) to see how you can calculate the time-shift required to align chromatograms manually with Mnova).



In this dialog, it is also possible to **Apply Trace Baseline Correction**. The **Filter** and **Smooth Factor** can be set manually by (1) typing the value, (2) using the cursor, or automatically by (3) autodetection.



Once you have adjusted your chromatogram parameters to your satisfaction, remember to click **OK** to save your settings.

In the main LC/MS Evaluation tab, you can configure some other parameters related to the chromatograms. For instance:

- If you wish to modify the peak detection parameters, you can access the **Peak Detection Options** by pressing this  button.
- If you are using a mass to assign a component, it will be necessary to set the **Tolerance** for the auxiliary EIC extractions generated by the automation.

EIC Tolerance:

- If you are using a molecular formula to assign a component, you can access and modify the **MolMatch settings** by pressing this  button.
- If you are using an RT or RT range that are defined prior to chromatogram alignment, check the **RT alignment correction** option to make sure those values are corrected.


RT alignment correction

- When re-analyzing a dataset, you can check the **Force RT Range as Integrated area** option to use the RT range as the total reaction component peak width. It's important to note that this option is accessible exclusively when initiating a re-analysis from the Mgears result viewer as explained in the [Results section](#). Enabling this option ensures consistent quantitation of a component across all samples on the plate, as it uses the same peak range for measurement in every case.

Force RT Range as Integrated Area

3.3.2.2. The Unknowns tab


Chrom RO allows you to evaluate and report unknown peaks found in your spectra. If you wish to do this, you must enable the **Evaluate unknowns** option and set the peak **Area threshold**. The unknown peaks will be reported in the Mgears viewer, and the reports generated by the analysis.




Additionally, you have the option to **Report unknowns in last CSV columns**. When this option is enabled, the RT and areas of the unknown peaks in each sample will be added to the last columns of the CSV output.

3.3.2.3. The Formulas tab

Within the Chrom RO tab, you have the ability to craft your own custom output values by creating formulas based on the default outputs generated by the plugin. To get started with this process, click on the pencil icon  to enter a new formula. This action will open a new dialog box.




Inside the dialog box, click the **Add** button . This will insert a new line into the table. Edit the parameter's name as needed. Then, enter your custom formula in the designated space and decide how many decimal places you want to display in the resulting value.




The formulas you create can only be generated using a limited set of parameters, including component areas (and concentrations, if calibration data is available), as well as the total area and total concentration. To see the available parameters for each analysis, click this icon .

Please note that the total area or concentration does not include the Internal Standards, as they are not part of the reaction. However, you can use the internal standard area or concentration for calculations.



Tip. Clicking on an item in the dropdown menu will copy the parameter name to the clipboard, allowing you to conveniently paste it into the formula editor cell.

After you have written your formula, you can test its validity by clicking the checkmark button. If the formula is not valid, a warning message will appear, prompting you to make the necessary adjustments.




3.3.3. The Quality Controls tab

In the **Quality Controls** tab, you can activate a series of quality tests to identify various anomalous situations that may require manual correction. Failed controls will trigger a flag for the sample and report the components or peaks that failed the test, facilitating a review-by-exception approach. There are five options available which can be configured:

- **Peak Not Assigned:** If enabled, a flag will be raised when peak assignment for a component fails.
- **Assignment Overlap:** If enabled, a flag will be obtained when two or more components are assigned to the same peak.
- **RT Outlier:** If enabled, a threshold must be defined (a factor of the Standard Deviation; xSTD). A flag is obtained when the distance [RT-mean RT] exceeds the accepted threshold.

- **Multiple Peak in RT range:** This control is only available when an RT range is defined for at least one component. If enabled, a flag is obtained when multiple peaks are detected in the specified RT range.
- **RT/MS Mismatch:** This control is only available when an RT value or range and a mass (SMILES, structure, molecular formula, or mz provided) are provided for at least one component. If enabled, a flag is obtained when the matched mz peak does not elute at the define RT or RT range.

Enable the option(s) of interest, then edit the flag text and click on the color selection button to choose the flag color. Click on **OK** to save your preferences.




An additional control on the overall well plate is also available, the **Accepted RT variation**. If enabled, a threshold must be defined, and a flag will be obtained when the spread (max-min) of the RT of a component across all samples exceeds the accepted threshold.

Finalize your plugin settings setup by pressing **OK** and move to the next step.

3.3.4. The Output tab

Within the **Output** tab, you have the flexibility to customize the CSV output configuration to align it with your specific requirements. To begin, you can choose whether to name this file in accordance with the results folder you will configure in the Mgears **Output** tab or opt for the default naming convention, which is "results.csv."




Next, simply click on **Configure Output** to access the dialog box where you can select which parameters to include in your output. This action will open a new window, where each tab corresponds to a distinct source of information. These sources encompass the outcomes of the automation process, data extracted from the input CSV file (if it is part of your workflow), and parameters derived from the raw datafile's metadata, typically found in the **Parameters table** within Mnova. For all three data sources, you will have the flexibility to export the values to the output CSV file, the Mgears Viewer, and the HTML report.

Note. If you have added a CSV mapping file in the **Input** tab, Mgears will automatically detect this file and will allow you to configure the output CSV based on its structure, as is the case in this analysis.

3.3.4.1. The Automation Variables tab

By default, this table will include:

- The **Name** of the sample and the **Well** plate position information, if provided, during the mapping file configuration (in the **Mgears Input** tab).
- For each component (other than Internal Standard or Ignored components) declared in the **Chrom RO Input** tab:
 - **RT**
 - **Area**
 - **%Area**
 - **Relative Area to the Internal Standard** (If an Internal Standard is present)
 - **Concentration** (If there is a calibration for the component)
 - **%Concentration** (If there is a calibration for the component)
 - **Relative Concentration to the Internal Standard** (If a calibration is applied to both the Internal Standard and component)



	Parameter	Output Column	Output Header	Print to HTML
1	Name	1	Name	<input checked="" type="checkbox"/>
2	Well	0	Well Position	<input type="checkbox"/>
3	SM (RT)	5	SM (RT)	<input checked="" type="checkbox"/>
4	SM (Area)	6	SM (Area)	<input type="checkbox"/>
5	SM (% Area)	7	SM (% Area)	<input checked="" type="checkbox"/>
6	SM (Rel. Area)	0	SM (Rel. Area)	<input type="checkbox"/>
7	SM (Conc.)	0	SM (Conc.)	<input type="checkbox"/>
8	SM (% Conc.)	0	SM (% Conc.)	<input type="checkbox"/>
9	SM (Rel. Conc.)	0	SM (Rel. Conc.)	<input type="checkbox"/>
10	P1 (RT)	0	P1 (RT)	<input checked="" type="checkbox"/>
11	P1 (Area)	0	P1 (Area)	<input type="checkbox"/>

Test Output CSV: 

OK Cancel

Annotations:

- Always
- If provided during the mapping file configuration
- Always
- When an internal standard is declared
- If absolute quantification is required
- If calibration is applied to both Internal Standard and component

- For a component declared as Internal Standard (If present):
 - **RT**
 - **Area**
 - **Relative Area** (Its value is expected to be 1)

- **Concentration** (If there is a calibration for the component)
- **Relative Concentration** (Its value is expected to be 1)

IS (RT)	0	IS (RT)	<input type="checkbox"/>	Always when an internal standard is declared
IS (Area)	0	IS (Area)	<input type="checkbox"/>	
IS (Rel. Area)	0	IS (Rel. Area)	<input type="checkbox"/>	
IS (Conc.)	0	IS (Conc.)	<input type="checkbox"/>	If absolute quantification is required
IS (Rel. Conc.)	0	IS (Rel. Conc.)	<input type="checkbox"/>	

- If an Internal Standard is present, the **Sum of relative areas** of all components is also included.

Sum (Rel. Area)	0	Sum (Rel. Area)	<input type="checkbox"/>
-----------------	---	-----------------	--------------------------

- If the option to **Evaluate Unknowns** is enabled, the **Unknowns Area**, total **percentage of Unknown Peaks area** and a text describing the unknown **Peaks** [list of RT(%Area, mz)] are also included.

Unknowns (Area)	0	Unknowns (Area)	<input type="checkbox"/>
Unknowns (% Area)	0	Unknowns (% Area)	<input type="checkbox"/>
Unknowns (Peaks)	0	Unknowns (Peaks)	<input type="checkbox"/>

For each **Parameter**, it is possible to configure an **Output Column** and **Output Header** and decide whether to include them in the HTML report (**Print to HTML**).

To choose the **Output column**, you can either:

- Select the cell and type the value:

Automation Variables			
Parameter	Output Column	Output Header	
1 Name	0	Name	<input type="checkbox"/>
2 Well	0	Well Position	<input type="checkbox"/>
3 SM (RT)	0	SM (RT)	<input type="checkbox"/>
4 SM (Area)	0		
5 SM (% Area)	0	SM (% Area)	<input type="checkbox"/>

Click and type value

B. Double-click on the cell to open the **CSV Column Detection** assistant and select the desired column:

The first screenshot shows the 'CSV Input/Output Configuration' dialog with a table where the 'Output Column' for 'SM (RT)' is 0. A red circle '1' and the text 'Double-click' point to this cell.

The second screenshot shows the 'CSV Column Detection' dialog. It asks to 'Click or Select column with values for: Output Column for SM (RT)'. A table lists columns 1-5. Column 2 is highlighted in red, with a red circle '2' over it. A red circle '3' is over the 'OK' button.

The third screenshot shows the 'CSV Input/Output Configuration' dialog again, but now the 'Output Column' for 'SM (RT)' is 2. A red circle '4' is over this cell.


C. Click on this button to automatically generate output column values for all the parameters in the table:

The first screenshot shows the 'CSV Input/Output Configuration' dialog with a 'Print to HTML' button highlighted by a red circle '1'.

The second screenshot shows the 'CSV Input/Output Configuration' dialog with the 'Print to HTML' button highlighted by a red circle '2'. The table below shows the resulting configuration with auto-generated output column values.


Parameter	Output Column	Output Header	Print to HTML
1 Name	1	Name	<input type="checkbox"/>
2 Well	2	Well Position	<input type="checkbox"/>
3 SM (RT)	3	SM (RT)	<input type="checkbox"/>
4 SM (Area)	4	SM (Area)	<input type="checkbox"/>
5 SM (% Area)	5	SM (% Area)	<input type="checkbox"/>
6 P1 (RT)	6	P1 (RT)	<input type="checkbox"/>
7 P1 (Area)	7	P1 (Area)	<input type="checkbox"/>
8 P1 (% Area)	8	P1 (% Area)	<input type="checkbox"/>

If you already have assigned certain columns to other parameters in the other tabs (**CSV File** and/or **Mnova Parameters Table**), a message will appear to inform you that the automatic assignment will start right after the last assigned column, as seen in the image below.




To clear **Output Column** values, you can press this button .

To edit the CSV **Output Header**, you must select the cell and type the header of your choice.




To include the parameters to the HTML report, tick the corresponding box in the final column.

You can quickly select/unselect all parameters using the buttons / , respectively.



Note. The **Automation Variables** tab will be updated whenever a component is added or modified in the [Reaction Components table](#) of the plugin.

In this analysis, we will use the following configuration for the **Automation Variables**.




3.3.4.2. The CSV File tab

Here you can manually add other parameters included in your input mapping file to the output CSV. Click on the add button, . A new row will be added to the table.

Click and edit the parameter name, **Input column** and **Header**, and **Output column** and **Header**, etc.

Enable the **Print in Viewer/HTML** if you wish to include the parameter to the viewer and the HTML report. The selected parameter will be added in the **Metadata** tab of the [Results section](#) in **Mgears Viewer**.

Repeat the workflow to add all the parameters you wish to include.



You can always select a row and click on to remove an entry.

3.3.4.3. The Mnova Parameters Table

To include a parameter from the **Mnova parameters table** in the final report, a datafile must be open in Mnova (1). Click on the add button **+** (2). The program will collect the Mnova parameters from the active spectrum and display them in a new dialog. Select the parameters you want to export (3) and click on **OK** to save your choices (4). The selected parameters will appear in the table as shown below (5).

Add a value for the **Output column**, edit the **Output Header**, and enable printing the results in the Mgears viewer and HTML report if needed.

1 An Mnova datafile must be open

2 +

3 Acquired Date

4 OK

5 Add a value


Parameters retrieved from the active document

Parameter	Output Column	Output Header	Print to Viewer/HTML
1 Acquired Date	8	Acquisition date	<input type="checkbox"/>

Parameter	Selected
1 Data Path	<input type="checkbox"/>
2 Acquired Date	<input checked="" type="checkbox"/>
3 Sample Name	<input type="checkbox"/>
4 Sample ID	<input type="checkbox"/>
5 Task Code	<input type="checkbox"/>
6 Job Code	<input type="checkbox"/>
7 Title	<input type="checkbox"/>
8 Description	<input type="checkbox"/>
9 Comment	<input type="checkbox"/>
10 User	<input type="checkbox"/>
11 Operator	<input type="checkbox"/>
12 Submitter	<input type="checkbox"/>

You can always select a parameter and press **-** to remove it from your CSV output.

Now that you have completed the **CSV Input/Output configuration**, you can click on at the bottom left-hand side of the configuration dialog to view a simulation of the output CSV. An overview of the CSV output will be displayed. The columns that will contain automation data are colored grey.




3.4. Design


In the **Design** tab, you can use Mnova layout templates to produce final, customizable reports or even select your own scripts.

3.5. Output


Choose a directory to which Mgears can save the analysis results.



Optionally, enable the **Add Nickname to the Results Folder** and type the nickname of your choice.




Enable the Mnova document creation to be able to review the resulting spectra in Mnova.



You can also choose to create a PDF with the results or save the results to a database (*please refer to the [Mnova Gears manual](#) for more details on the output configuration options*).

Now that everything is set, click on  Run to start the analysis.




4. The output folder

The results folder is saved under the directory previously specified and contains all the output generated by the current evaluation.


Analysis date and time_Nickname

2022-04-25T12.29.22_LastAnalysis




4.1. The HTML file

In the HTML report, an overview of the results is made available. Each sample is reported in a row, along with the predefined parameters (as per the [output settings](#)).




The HTML format allows you to change the display of the columns and adapt them to your preferences. Click on **Columns** and select/unselect the parameters you wish to display/hide.

You can also **Copy**, **Print**, or save these results into **CSV** or **PDF** formats.



4.2. The CSV report

The output folder named “ChromRO” contains the CSV file with all the output and input information as predefined in the [output settings](#).



2022-04-25T12.29.22_LastAnalysis

Name

- ChromRO
- css
- data
- documents

overallResults

results

	CSV Input				Automation variables			Mnova parameters		Unknowns	
	A	B	C	D	E	F	G	H	I	J	
1	Name	Temperature	Catalyst	Solvent	SM (RT)	SM (Area)	SM (% Area)	Acquisition date	Unknown 1 (RT)	Unknown 1 (Area)	
2	REACTION-01	0C	NaBH3(CN)	EtOH	0.957	2468128377	28.99	2021-06-14T12:12:17.000			
3	REACTION-02	20C	NaBH3(CN)	EtOH	0.96	2452735496	27.44	2021-06-14T12:16:17.000			
4	REACTION-03	0C	NaBH3(CN)	THF	0.957	2404991750	45.29	2021-06-14T12:20:02.000			
5	REACTION-04	20C	NaBH3(CN)	THF	0.956	2475581697	40.8	2021-06-14T12:23:49.000			
6	REACTION-05	0C	NaBH(OAc)3	EtOH	0.961	2498589554	67.37	2021-06-14T12:27:33.000			
7	REACTION-06	20C	NaBH(OAc)3	THF	0.964	2464925096	45.18	2021-06-14T12:31:16.000			
8	REACTION-07	20C	Ti(OiPr)4 - NaBH4	THF	0.956	2495984154	38.58	2021-06-14T12:34:59.000			
9	REACTION-08	40C	Ti(OiPr)4 - NaBH4	THF	0.96	2546771623	24.56	2021-06-14T14:04:56.000	0.477366122	6352361384	
10	REACTION-09	20C	Ti(OiPr)4 - NaBH4	THF - EtOH	0.964	2499719297	30.39	2021-06-14T14:08:40.000			
11	REACTION-10	40C	Ti(OiPr)4 - NaBH4	THF - EtOH	0.959	2465778118	57.6	2021-06-14T12:46:10.000			

4.3. Other output

- A “documents” directory, containing the output Mnova files (unless Mgears is configured to save Mnova files in another location).
- A log file of the execution.
- A copy of the settings used in the current evaluation.
- A resume file of the steps followed in the execution.
- A CSS folder, a data folder, a JS folder, and an images folder.

5. The Mnova Gears Results Viewer

Chrom RO has two different – but interconnected – result viewers, the Mgears Viewer, which displays overall results for the sample/well, and the Chrom Reaction Optimization Viewer, which displays individual components RT statistics for each reaction dataset.

5.1. The Mgears Viewer

Open the Mgears Viewer, then click on  to open your analysis results.




When the experiment is open, the **Mgears Viewer** shows the individual results datasets, a well plate view, and the associated numerical results and metadata.

Results file name

Datasets list

Well plate display

Results



The screenshot displays the Mgears Viewer interface. At the top, the title bar reads "Mgears Viewer - 2023-10-03T14:30:03_FormulaTest". Below the title bar is a toolbar with various icons and an "Analyze Again" button. The main area is divided into three sections:

- Datasets list:** A table with columns "#", "Title", and "Document". It lists three datasets: REACTION-01, REACTION-02, and REACTION-03.
- Well plate display:** A grid showing numerical values for a well plate. The first row is labeled "A" and the first column is labeled "1".
- Results:** A detailed view for "REACTION-01" showing a table with columns "Name", "Type", "RT", "Area", and "% Area".

#	Title	Document
1	REACTION-01	REACTION-01.mnova C:/Users/Usuario/Desktop/Results/Chrom React
2	REACTION-02	REACTION-02.mnova C:/Users/Usuario/Desktop/Results/Chrom React
3	REACTION-03	REACTION-03.mnova C:/Users/Usuario/Desktop/Results/Chrom React


	1	2	3	4	5	6	7	8	9	10	11	12
A	29.20	27.71	45.43	41.10	67.31	45.31	39.02	24.80	30.74	58.26		

Name	Type	RT	Area	% Area
1 SM	Start Material	0.960	2643937169	29.20
2 P1	Product	0.576	2851935555	31.49
3 P2	By Product	0.485	3434576553	37.93
4 P3	Ignored	-	-	-

Click on a specific dataset to view the corresponding results and spectrum in Mnova.

5.1.1. The Well plate view

The **Well Plate** view is highly adaptable, allowing you to present results in a manner that suits your preferences. You can choose to showcase numerical data, such as the percentage area of various reaction components or the outcomes of your custom formulas. Alternatively, you can opt for graphical representations, including bar graphs and pie charts, to visualize the results.




Some examples of results presentations are given below:


SM (% Area)

	1	2	3	4	5	6	7	8	9	10
A	29.20	27.71	45.43	41.10	67.31	45.31	39.02	24.80	30.74	58.26


Bar Chart Areas



SM (% Area) (Bubbles - Enhanced graphics*)




Pie Chart Areas



SM (% Area) (Heatmap - With enhanced graphics*)

	1	2	3	4	5	6	7	8	9	10
A	29.20	27.71	45.43	41.10	67.31	45.31	39.02	24.80	30.74	58.26

Controls Chart



* The Enhanced graphics can be configured in the Mgears Viewer Settings



5.1.2. The Results section

The **Results** section includes different tabs with information relating to:

- Reaction components **Results** with the Matched peak's RT, Area, and %Area.
- **Unknowns** results, if [unknowns' analysis](#) is enabled in the plugin settings.
- **Formulas** results, if [custom formulas](#) are introduced to the analysis.
- Extracted **Metadata**, if [output parameters](#) are set to appear in the Mgears viewer.
- **Controls** results, if [quality tests](#) are enabled.
- **RT Statistics** for each component peak within the well, encompassing the RT and its deviation from the mean RT calculated across all samples as well as a score to evaluate the RT outliers (defined as $[RT - \text{Mean}]/\text{STD}$).

Within the **Results** section, you will also find a set of tools designed to facilitate the correction and refinement of your analysis. For example, if a peak assignment has encountered an error or is incorrect, you have the option to effortlessly remove the assigned peak and proceed to manually reassign this component to another peak within your spectrum or spectra. You can also open the **Chrom Reaction Optimization** dialog box to edit some settings and quickly relaunch the analysis on the same samples.

5.1.2.1. Removing peak assignment

To address an incorrectly assigned component, simply click on the corresponding component row (1), and then press the button (2). Mgears will recalculate and update the results for the selected well (3).

Result in Well Plate: Pie Chart Areas


Results Unknowns Formulas Metadata Controls RT Statistics

Name	Type	RT	Area	% Area
SM	Start Material	1.002	45712667	0.50
P1	Product	0.576	2851935555	31.49
P2	By Product	0.485	3434576553	37.93

M. ? X
Recalculating...
Cancel


Deleting this peak assignment to SM

All the Mgears tables and spectra are automatically updated, as can be seen below:




5.1.2.2. Re-assigning a peak

To manually assign a peak to a component, select the component row (1), then press the button (2). A dialog with a list of detected chromatographic peaks appears. Select the peak you wish to assign (3) and press **OK** (4). Mgears will recalculate and update the results for the selected well (5).



All the Mgears tables and spectra are automatically updated, as can be seen below:




5.1.2.3. Editing analysis settings


To edit analysis settings and quickly relaunch the analysis on the same samples, press on the **Brick** icon  in the **Results** section. This will open the **Chrom Reaction Optimization Settings** dialog box from which you can add a new component or modify the RT, color, and external calibration of the added components in the **Input** tab. In the **LC/MS Evaluation** tab, only the **Force RT Range as Integrated Area** option can be modified. In addition, you can modify the Unknowns window, add, delete, or modify Formulas, edit the Controls and modify the Output.

If you enable the **Force RT Range as Integrated Area** option the peak is integrated in exactly the provided RT range.

When you have finished, click on **OK** to save your changes.




Now click on **Analyze again** in the Mgears viewer and select **Analyze Again All Results**. Mgears will recalculate and update the results for the whole well plate.



The new results will be automatically saved to your output folder.

5.2. The Chrom Reaction Optimization Viewer

You can launch the **Chrom Reaction Optimization Viewer** directly from the Mgears Viewer interface by clicking the  button. This viewer includes a comprehensive **Peaks Statistics** table featuring essential data such as peak RT mean, standard deviation, variation range, and the outcome of the Control test, if enabled.



The screenshot shows the 'Chrom Reaction Optimization Viewer' window. It contains a table titled 'Peaks Statistics' with the following data:

	Component	RT mean	RT STD	RT Range	RT Variation Control Passed
1	SM	0.961	9.14e-4	0.003	✓
2	P1	0.571	2.51e-3	0.007	✓
3	P2	0.480	2.12e-3	0.007	✓
4	P3	0.000	0.00e+0	NaN	✓

For more details on Mnova Gears' options, please refer to the [Mnova Gears manual](#).