

# Streamlining Cannabis Testing Using GC-MS and GCxGC-MS

David E. Alonso and Joe Binkley | LECO Corporation, St. Joseph, Michigan 49085




## Background

- The composition of cannabis is very important in determining its utility, potency, and medicinal effectiveness
- Cannabis is difficult to analyze because of its chemical diversity, but more importantly the wide concentration ranges of its constituents
- Analysis often requires:
  - Tedious sample preparation
  - Specialized (targeted) instrumental analysis

## Objectives

- Reduce sample manipulation & shorten preparation times
- Explore the possibility of streamlining cannabis screening
- Maximize the information obtained per analysis
- Implement the use of GC-TOFMS and GCxGC-TOFMS to effectively profile cannabis samples
- Use automated data processing software to quickly and confidently identify compounds
- Ultimately, the question to be answered... *Is the product safe?*



## Samples



## Instrument Parameters

Gas Chromatograph	Agilent 7890, LECO Dual Stage Quad Jet Modulator & L-PAL 3 Autosampler
Injection	Variable (250 °C)
Carrier Gas	He @ 1.0 mL/min, Constant Flow
Columns (1 <sup>st</sup> Dimension)	Rxi-5 MS, 30 m x 0.25 mm i.d. x 0.25 µm (Restek, Bellefonte, PA, USA)
Columns (2 <sup>nd</sup> Dimension)	Rxi-17 SIL MS 0.6 m x 0.25 mm i.d. x 0.25 µm (Restek, Bellefonte, PA, USA)
Temperature Program	40 °C (5 min), ramped 10 °C/min to 300 °C (2 min) Secondary oven maintained +5 °C relative to primary oven
Modulation	2s with temperature maintained +15 °C relative to secondary oven
Mass Spectrometer	LECO Pegasus® BT 4D
Ion Source Temperature	250 °C
Ionization Mode	EI
Mass Range (m/z)	45-600
Acquisition Rate	10 spectra/s (1D), 200 spectra/s (2D)

## Sample Preparation Methods

Name	Formula	R.T. (s)	Mass Δ (Da)	Similarity
Ethylbenzene	C <sub>8</sub> H <sub>10</sub>	323	0.00	946
o-Xylene	C <sub>8</sub> H <sub>10</sub>	331	0.00	939
2-Heptanone	C <sub>7</sub> H <sub>14</sub> O	349	0.00	813
p-Xylene	C <sub>8</sub> H <sub>10</sub>	354	0.00	956
Heptanal	C <sub>7</sub> H <sub>14</sub> O	359	0.00	926
Benzene, (1-methylethyl)-	C <sub>8</sub> H <sub>10</sub>	421	0.00	824
Benzene, 1,2,3-trimethyl-	C <sub>9</sub> H <sub>12</sub>	453	0.00	809
Octanal	C <sub>8</sub> H <sub>16</sub> O	459	0.02	937
Heptanoic acid	C <sub>7</sub> H <sub>14</sub> O <sub>2</sub>	520	-0.01	873
3-Methylcoumaran	C <sub>9</sub> H <sub>10</sub> O	582	0.04	872
Octanoic acid	C <sub>8</sub> H <sub>16</sub> O <sub>2</sub>	610	-0.01	911
2-Decanone	C <sub>10</sub> H <sub>20</sub> O	636	0.01	812
Dec-(5E)-en-1-ol	C <sub>10</sub> H <sub>20</sub> O	692	0.03	878
Nonanoic acid	C <sub>9</sub> H <sub>18</sub> O <sub>2</sub>	695	0.00	839
Hexanoic acid, hexyl ester	C <sub>12</sub> H <sub>24</sub> O <sub>2</sub>	797	-0.04	821
4-Chloro-4'-methoxybiphenyl	C <sub>17</sub> H <sub>13</sub> ClO	869	0.12	779
n-Hexadecanoic acid	C <sub>16</sub> H <sub>32</sub> O <sub>2</sub>	1185	-0.01	923
Linoleic acid	C <sub>18</sub> H <sub>34</sub> O <sub>2</sub>	1284	0.01	888
Stearic acid	C <sub>18</sub> H <sub>36</sub> O <sub>2</sub>	1297	0.01	824
γ-Tocopherol	C <sub>29</sub> H <sub>50</sub> O <sub>2</sub>	1727	0.04	828
dl-α-Tocopherol	C <sub>29</sub> H <sub>50</sub> O <sub>2</sub>	1774	0.04	866
Phytosterol	C <sub>27</sub> H <sub>48</sub> O	1923	0.04	893
α-Amyrin	C <sub>30</sub> H <sub>50</sub> O	2010	0.03	887

Figure 1: List of miscellaneous compounds in cannabis with deconvoluted (top) and NIST library (bottom) spectra shown for 4-Chloro-4'-methoxybiphenyl.

## GC-TOFMS Results: Terpenes

Name	Formula	R.T. (s)	M.D. (Da)	Similarity	Name	Formula	R.T. (s)	M.D. (Da)	Similarity	Name	Formula	R.T. (s)	M.D. (Da)	Similarity
3-Thujene	C <sub>10</sub> H <sub>16</sub>	386	0.00	896	Perillal	C <sub>10</sub> H <sub>16</sub> O	654	0.00	773	cis-Calamenene	C <sub>10</sub> H <sub>16</sub>	913	0.00	861
2-Pinene	C <sub>10</sub> H <sub>16</sub>	394	0.00	967	Verbenone	C <sub>10</sub> H <sub>16</sub> O	660	0.00	896	(E)-α-Bisabolene	C <sub>10</sub> H <sub>16</sub>	921	0.00	946
Camphene	C <sub>10</sub> H <sub>16</sub>	410	0.00	952	Geraniol	C <sub>10</sub> H <sub>16</sub> O	690	0.07	814	α-Calacorene	C <sub>10</sub> H <sub>16</sub>	929	0.00	878
4(10)-Thujene	C <sub>10</sub> H <sub>16</sub>	433	0.00	925	3-Terpinolene	C <sub>10</sub> H <sub>16</sub>	772	0.00	805	Caryophyllene oxide	C <sub>10</sub> H <sub>16</sub> O	961	0.00	936
β-Myrcene	C <sub>10</sub> H <sub>16</sub>	448	0.00	954	α-Cubebene	C <sub>10</sub> H <sub>16</sub>	776	0.00	909	Guaiol	C <sub>10</sub> H <sub>16</sub> O	966	0.00	935
α-Phellandrene	C <sub>10</sub> H <sub>16</sub>	463	0.00	865	Ylangene	C <sub>10</sub> H <sub>16</sub>	796	0.00	926	(-)-Globulol	C <sub>10</sub> H <sub>16</sub> O	975	0.00	918
β-Caradiene	C <sub>10</sub> H <sub>16</sub>	470	0.00	925	Capene	C <sub>10</sub> H <sub>16</sub>	800	0.00	926	Humulene epoxide 2	C <sub>10</sub> H <sub>16</sub> O	979	0.00	868
D-Limonene	C <sub>10</sub> H <sub>16</sub>	487	0.00	866	β-Guaiane	C <sub>10</sub> H <sub>16</sub>	805	0.00	780	Selin-6-en-4e-ol	C <sub>10</sub> H <sub>16</sub> O	985	0.01	814
1,8-Cineole	C <sub>10</sub> H <sub>16</sub> O	491	0.00	879	β-Yalagene	C <sub>10</sub> H <sub>16</sub>	808	0.00	749	(-)-10-epi-Eudesmol	C <sub>10</sub> H <sub>16</sub> O	987	0.00	928
β-Cidmene	C <sub>10</sub> H <sub>16</sub>	503	0.00	970	β-Cubebene	C <sub>10</sub> H <sub>16</sub>	810	0.00	798	Aristolene	C <sub>10</sub> H <sub>16</sub>	1002	0.00	756
γ-Terpinene	C <sub>10</sub> H <sub>16</sub>	516	0.00	835	Selina-5,11-diene	C <sub>10</sub> H <sub>16</sub>	817	0.00	884	(-)-α-Panasinsen	C <sub>10</sub> H <sub>16</sub>	1012	-0.24	767
3,8-p-Menthadiene	C <sub>10</sub> H <sub>16</sub>	545	0.00	758	Cadinene	C <sub>10</sub> H <sub>16</sub>	828	-0.02	893	α-Bisabolol	C <sub>10</sub> H <sub>16</sub> O	1022	0.01	922
L-Fenchone	C <sub>10</sub> H <sub>16</sub> O	547	0.00	795	Santalene	C <sub>10</sub> H <sub>16</sub>	832	0.00	866	Ledene	C <sub>10</sub> H <sub>16</sub>	1029	0.00	855
Linalool	C <sub>10</sub> H <sub>16</sub> O	552	0.01	941	Isocaryophyllene	C <sub>10</sub> H <sub>16</sub>	837	0.00	960	Junipercamphor	C <sub>10</sub> H <sub>16</sub> O	1038	0.01	853
1,3,8-p-Menthatriene	C <sub>10</sub> H <sub>16</sub>	567	0.00	814	α-Guaiane	C <sub>10</sub> H <sub>16</sub>	848	0.00	825	Aromadendrane-4,10-diol	C <sub>10</sub> H <sub>16</sub> O <sub>2</sub>	1101	0.01	824
Fenchol	C <sub>10</sub> H <sub>16</sub> O	570	0.00	964	Aromadendrene	C <sub>10</sub> H <sub>16</sub>	852	0.00	930	Isaaromadendrene epoxide	C <sub>10</sub> H <sub>16</sub> O	1127	0.00	789
Pinanone	C <sub>10</sub> H <sub>16</sub> O	579	0.00	885	(E)-β-Farnesene	C <sub>10</sub> H <sub>16</sub>	854	0.02	884	Spathulanol	C <sub>10</sub> H <sub>16</sub> O	1144	0.00	751
cis-p-Mentha-2,8-dien-1-ol	C <sub>10</sub> H <sub>16</sub> O	588	0.06	797	Humulene	C <sub>10</sub> H <sub>16</sub>	864	0.00	955	Corymbolone	C <sub>10</sub> H <sub>16</sub> O <sub>2</sub>	1169	0.00	794
Alloaromadendrene	C <sub>10</sub> H <sub>16</sub>	595	0.00	889	Alloaromadendrene	C <sub>10</sub> H <sub>16</sub>	869	0.00	959	m-Camphorene	C <sub>10</sub> H <sub>16</sub>	1188	0.01	851
(-)-2-Bornanone	C <sub>10</sub> H <sub>16</sub> O	601	0.01	850	γ-Muurolene	C <sub>10</sub> H <sub>16</sub>	878	0.00	855	p-Camphorene	C <sub>10</sub> H <sub>16</sub>	1208	0.01	818
Pinocaradiene	C <sub>10</sub> H <sub>16</sub> O	617	0.00	813	α-Selinene	C <sub>10</sub> H <sub>16</sub>	895	0.00	943	Phytol	C <sub>10</sub> H <sub>16</sub> O	1274	0.02	951
endo-Bornol	C <sub>10</sub> H <sub>16</sub> O	620	0.00	963	β-Bisabolene	C <sub>10</sub> H <sub>16</sub>	897	0.00	966	Neophyladiene	C <sub>10</sub> H <sub>16</sub>	2120	-0.02	887
Terpinen-4-ol	C <sub>10</sub> H <sub>16</sub> O	629	0.00	895	γ-Cadinene	C <sub>10</sub> H <sub>16</sub>	907	0.00	932					
α-Terpineol	C <sub>10</sub> H <sub>16</sub> O	640	0.01	943	δ-Cadinene	C <sub>10</sub> H <sub>16</sub>	912	0.00	892					




Figure 2: Deconvoluted (top) and NIST library spectra (bottom) for B-Ocimene and m-Camphorene.

## GC-TOFMS: Cannabinoids

Name	Formula	R.T. (s)	Mass Delta (Da)	Similarity
Cannabidiol	C <sub>21</sub> H <sub>32</sub> O <sub>2</sub>	1342	0.01	882
Cannabidiol	C <sub>21</sub> H <sub>32</sub> O <sub>2</sub>	1445	0.07	905
Delta 8-THC	C <sub>21</sub> H <sub>32</sub> O <sub>2</sub>	1483	0.01	837
Delta 9-THC	C <sub>21</sub> H <sub>32</sub> O <sub>2</sub>	1491	-0.01	928
Cannabigerol	C <sub>21</sub> H <sub>32</sub> O <sub>2</sub>	1510	0.01	900
Cannabinol	C <sub>21</sub> H <sub>32</sub> O <sub>2</sub>	1522	0.01	902
Cannabichromene	C <sub>21</sub> H <sub>32</sub> O <sub>2</sub>	1578	N/A	873





Figure 3: List of representative cannabinoids with deconvoluted (top) and NIST library (bottom) spectra shown for Delta 9-THC.

## GCxGC-TOFMS Contour Plot Showing Separation of Several Compound Classes



## GCxGC-TOFMS - More Compounds Identified

Name	Formula	R.T. (s)	Similarity
Decanal	C <sub>10</sub> H <sub>20</sub> O	Not Found	
Myrtenal	C <sub>10</sub> H <sub>16</sub> O	648	745
Selina-3,7(11)-diene	C <sub>15</sub> H <sub>24</sub>	929	805
α-Calacorene	C <sub>15</sub> H <sub>20</sub>	929	878

Name	R.T. (s)	Similarity
Decanal	648 s, 1.262 s	723
Myrtenal	648 s, 1.445 s	890
Selina-3,7(11)-diene	928 s, 1.323 s	869
α-Calacorene	928 s, 1.414 s	881

## Summary

- A variety of sample introduction methods can be used to prepare cannabis products for analysis
- Comprehensive profiling of cannabis is an attractive alternative to targeted, panel-based methods
- GCxGC-TOFMS facilitates fast and confident cannabis profiling.
- Compound identification: Similarity searches, mass Δ calculations, and RI filtering