

# The Chromatography Checklist

Rita Steed  
Alexander Ucci  
May 23, 2018



# Checklist

- Is your sample ready?
  - ✓ Solubility: sample solvent and starting mobile phase conditions
  - ✓ Filters
  - ✓ Sample Clean-Up
- Supplies
  - ✓ In-line filters, safety caps, quick connect fittings
- Instrument
  - ✓ Maintenance
  - ✓ Role
  - ✓ Housekeeping
- Method
  - ✓ Established or new
  - ✓ Conditions
- Column
  - ✓ Choice
  - ✓ Documentation
  - ✓ Guards
- Final checklist



## Sample Clean-Up

Filtration, Solid Phase Extraction, QUECHERS, and more!




# Why Perform Sample Clean-Up?

- To acquire desired sensitivity/selectivity
- To reduce contamination/carryover issues
- Use of sensitive and expensive instruments: Protect your investment!!!




Salt build-up in LC-MS ion source from unextracted salts

Pesticides in Avocado without SP



Pesticides in Avocado with SP



Curtain plate after injection of 25 samples with extractions from raisins without cleanup

# Which Sample Clean-Up technique is right for YOU?

## Solid Phase Extraction (SPE)

Multi-step approach for highest level of sample cleanup

## QuEChERS (dSPE)


Sample cleanup by extraction of bulk interferences

## Captiva EMR-Lipid (PPT and lipid removal)

Removes precipitated proteins by in-well protein precipitation and also removes lipids

## Filtration

Simple and fast removal of particulates



# Captiva Filtration and it's Benefits



Filtration is basic sample preparation method for all kinds of samples


Physically removes particulates from the sample

Prevents blocking of capillaries, frits, and column inlet (especially for UHPLC columns with 1.8 and 2.7  $\mu\text{m}$  particle sizes)

Results in less downtime of the instrument for repairs

Results in less wear and tear on the critical moving parts of the injection valves

96-well plate formats available



Unfiltered, centrifuged, and filtered plasma extracts  
Zorbax RRHD Eclipse Plus C18, 2.1 x 50 mm, 1.8  $\mu\text{m}$  column, PN 959757-902

Captiva Syringe Filters Guide 5991-1230EN

[Syringe Filter Selection Tool](#)



# Captiva EMR-Lipid





- One of Agilent's newest products with a 2 in 1 benefit of removing proteins and lipids
- Simple pass-through format
- Solvent-retention frit in 1 mL cartridge/96-well plate format for in well protein precipitation (*in situ*)
  - Unique cartridge/well construction minimizes clogging – and **ensures protein and lipid removal** (no cloudy samples)
- 3 and 6 mL cartridge format for larger samples
  - Do not contain solvent retain frit which allow for gravity flow
  - Protein precipitation performed offline (QUECHERS, etc.)
- Unique cartridge/well construction minimizes clogging – and **ensures protein and lipid removal** (no cloudy samples)
- High analyte recoveries
- Effective use will reduce ion suppression, increase analyte sensitivity, and detection, and extend the lifetime of your analytical column





# Enhanced Matrix Removal: EMR-Lipid

EMR-Lipid sorbent technology effectively traps lipids through two mechanisms:

- ✓ **Size exclusion** – Unbranched hydrocarbon chains (lipids) enter the sorbent; bulky analytes do not
- ✓ **Sorbent chemistry** – Lipid chains that enter the sorbent are trapped by hydrophobic interactions




## Effective phospholipid removal




# Protein Precipitation vs. Captiva EMR-Lipid RSD and Peak Area

## Captiva EMR-Lipid



## Protein Precipitation



Lipids cause reproducibility problems resulting in high RSD values

Using Captiva EMR-Lipid → low RSD values and higher peak areas

Higher peak area due to less ion suppression → can lead to lower detection limits

# QuEChERS

Screening of pesticide residues in fruit and vegetables

- Developed to make sample cleanup of food faster, simpler, less expensive, and greener

Now used with other matrices and compound classes as well

QuEChERS: Quick Easy Cheap Effective Rugged Safe

Commercially available kits allow for ease of use and convenience leading to increased throughput

Consists of two steps, and thus **2 kits**:

Step 1: Liquid Extraction



Step 2: Dispersive SPE / Interference Removal



# Bond Elut Solid Phase Extraction

The Trusted Name in Solid Phase Extraction

## SPE Benefits

- Selectivity ranges support a range of sample prep goals, from targeted matrix removal to target analyte concentration
- Easily automatable for increased throughput and reproducibility
- **Removes the widest set of matrix interferences for cleaner samples**
- Depth and breadth of published applications makes choosing, optimizing and implementing SPE methods easier
- Allows for lower detection limits and longer instrument uptime from cleaner extracts

## Bond Elut SPE Products Deliver

- **Widest selection of formats to accommodate different sample types** and sizes and to fit into lab workflows
- Over 50 unique polymeric, silica, and specialty sorbents for a wide range of application requirements
- Reliable, consistent performance to ensure best performance for your lab

## SPE Cartridge Options



## 96-well plate SPE



# Productivity Benefits with Sample Preparation

**More Matrix Removal = Less Matrix Entering System = Time and Cost Savings!**


- ✓ Less matrix build-up
  - Less interferences
  - Improved S/N
  - Better reproducibility
- ✓ Better chromatography
  - Less time spent on data analysis/manual integration
  - Less time spent on re-runs/recalibrations
- ✓ Less maintenance
  - Less instrument down-time
  - Saves \$\$ on consumables/services
- ✓ Less troubleshooting
  - “Is it my column or my MS”?
  - Less instrument down-time



# LC Supplies



# InfinityLab Solvent Bottles and Inlet Filters



- Inlet filters- Not a replacement for good mobile phase hygiene
- Glass solvent inlet filter (20  $\mu\text{m}$ ), 5041-2168
- Stainless steel solvent inlet filter (recommended for LC-MS), 01018-60028
- Agilent solvent filters show uniform particle size and superior pre-size homogeneity
- Agilent solvent inlet filters are packed in ultraclean antistatic bags with inner metallic coating, while other vendors use normal plastic bags.
- LC/MS analysis shows potential **contamination** through slip agent **when using 3<sup>rd</sup> party filters.**
- Optional bottle tags and bottle labels with Agilent's InfinityLab solvent bottles


# InfinityLab Stay Safe Caps for Solvent Bottles



The venting valve plays an important role in maintaining your workspace, keeping it safe and clean. It blocks harmful vapors and assures the required ventilation while solvent is being pulled to the LC instrument.

The time strip has a lifetime of 6 months and will alert you when to replace the venting valve

We also sell caps for waste container  
We sell a variety of different caps with different amount of ports



# Pump



## Items to have on hand:

- PTFE frits
- Pump piston seals



p/n: 01018-22707

## Typical Schedule for Pump

Item	Typical Schedule	Comments
Solvent inlet filter	Replace every 6 - 12 months	
PTFE frits in purge valve + gold seal	Every 12 months	
Piston seals	Every 12 months	When changing the seal, check the piston for scratches--replace if scratched
Inlet valve cartridge, outlet ball valve	Every 24 months	

Adjust according to your samples, conditions, and performance goals

# Detector – Which One & Why



## ➤ UV/DAD

- Popular, simple to use, reliable, sensitive
- Sample must have UV absorbance

## ➤ MS

- Sensitive
- Sample must be ionizable

## ➤ RI

- Refractive Index; difference between analyte and mobile phase
- Need strict temperature control

## ➤ ELSD

- Independent of a compound's absorbance, fluorescence, or electro-activity
- Enables detection of semi-volatile and thermally sensitive compounds

## ➤ FLD

- More selective and can be more sensitive
- Compounds must fluoresce; Compounds often derivatized

## ➤ ECD

- Very sensitive
- Can produce severe noise

# Detector Care

## MS Detectors



Flush the nebulizer	Daily after use to flush the tubing, valves, and nebulizer
Replace the nebulizer needle	When plugged
Clean the spray chamber	Daily or when carryover is suspected
Check the rough pump fluid level	Check weekly for color and level; replace every six months



Edwards pump  
Oil 6040-0834



MS40+ pump  
Oil 6040-1361



# Autosampler and Column Compartment



## Maintenance points on Autosampler

- Needle
- Loop capillary
- Needle seat
- Injection valve rotor seal
- Metering device seal



## In-line Filters

- In-line filters can help extend the life of your column
- Traps particulates that can plug column frits
- Not intended to be a replacement for good sample cleanup

## Guard Columns


- Protects analytical column from sample contaminants that are strongly or permanently retained




# Making Great Connections

Connection problems can lead to:

- *Downtime*
- *High cost of operation*
  
- Poor chromatography results
  - Broad or tailing peaks
  - Loss of resolution
- Expensive maintenance cost
  - Overtightening
  - Column damage
  - Leaks, added troubleshooting



# Simplifying column and fitting connections with InfinityLab Quick Connect and Quick Turn Fittings



Quick Connect fitting



Quick Turn fitting

- The unique spring-loaded design applies a constant force to eliminate dead volume. The spring-loaded design constantly pushes the tubing against the receiving port, delivering a reproducible connection with no dead volume for consistent chromatographic performance.
- Finger tight to 1300 bar: seals with a simple turn of the lever
- For instrument connections that are tight, you can use the InfinityLab QuickTurn fitting
  - Finger tight to 600 bar, wrench tight to 1300 bar
- Both fittings contain stem length that is adjustable through the spring which makes the fitting compatible with all types of LC columns



**Quick Connect fitting requires capillary with spring and a PEEK component**



**Quick Turn fitting needs capillary with long socket due to its internal spring action**

# Agilent A-Line Vials

**Maximum inertness:** The inert performance of Agilent A-Line vials, results in reduced analyte peak variability, so you can have the utmost confidence in your results.

**Consistent performance:** Vial-to-vial lot-to-lot Agilent A-Line vials demonstrate consistent performance, so you spend less time troubleshooting and rerunning samples.

**Certification of analysis:** Agilent A-Line vials come with a certificate of analysis, so you can be sure that they will perform even in the most demanding environments.

**Designed to fit a range of caps:** Agilent A-Line vials can be used with your existing 2 mL autosampler caps, for easier inventory management.

**Fewer septa issues:** Agilent septa are continually being improved to limit leaching, coring, sticking, push-through, hardness, and adsorption/absorption.



# Checklist

- Is your sample ready?
  - ✓ Solubility: sample solvent and starting mobile phase conditions
  - ✓ Filters
  - ✓ Sample Clean-Up
- Supplies
  - ✓ In-line filters, safety caps, quick connect fittings
- Instrument
  - ✓ Maintenance
  - ✓ Role
  - ✓ Housekeeping
- Method
  - ✓ Established or new
  - ✓ Conditions
- Column
  - ✓ Choice
  - ✓ Documentation
  - ✓ Guards
- Final checklist



# LC Instrument



- Maintenance
- Role
- Housekeeping



# Is your System Maintenance Up To Date?

## Typical Schedule\*

### PUMPS

Item	Typical Schedule	Comments
Solvent inlet filter	Replace every 6 - 12 months	
PTFE frits in purge valve + gold seal	Every 12 months	
Piston seals	Every 12 months	When changing the seal, check the piston for scratches--replace if scratched
Inlet valve cartridge, outlet ball valve	Every 24 months	

### AUTOSAMPLER

Item	Typical Schedule	Comments
Needle and needle seat	Every 12 months	
Rotor seal	Every 12 months	
Metering device seal	Every 24 months	

### COLUMN COMPARTMENT

Item	Typical Schedule	Comments
Column switching valve rotor seal	Every 12 months	
Column fittings	Every 5 to 10 column changes	A-line fittings last a lot longer than traditional fittings

### DETECTORS

Item	Typical Schedule	Comments
Lamps	Every 2000 hours	Watch for a noisy baseline
Flow cell	Check cleanliness every 6 months	Low light intensity could be caused by a dirty flow cell

\*Adjust according to your samples, conditions, and performance goals

**Check out** - A Gram of Prevention:  
Simple Tips for Maintaining Instrument  
Performance

<https://agilentseminar.webex.com/ec3300/eventcenter/recording/recordAction.do?siteurl=agilentseminar&theAction=poprecord&recordID=8057207&internalRecordTicket=4832534b00000004c45844572945e7b28ba731ae9f5854e008a76b4598b318dcb3f39f8efb42e39d>



<https://www.youtube.com/watch?v=vFU VHssMnx4>



<https://www.youtube.com/watch?v=iTiIOMH51Uc&index=11&list=PLThrdl2ragolmT3J-W5r8ailvJN94DJMR>

# Role – Pump

## Performance Characteristics

### Common to isocratic & gradient

- Flow accuracy
- Pressure pulsation

### Gradient Pump Only

- Delay volume in low and high pressure mixing
  - Determine your dwell volume (see appendix)
- Composition accuracy & precision

## *Influence on...*

- *RT & peak area precision*
- *Baseline noise*
  
- *Gradient shape and precision*
  
- *RT & peak area precision*



# Role – Injector

## Performance Characteristics

- Injection volume precision
- Wide linearity
- Minimum carryover


## *Influence on...*

- *Precision of peak area/height*
- *Accuracy of peak area/height (when using different injection volumes)*
- *Precision of peak/area height*



Peak Tailing – Injector Seal Failure

# Role - Thermostatted Column Compartment

Performance Characteristics	<i>Influence on...</i>	
<ul style="list-style-type: none"><li>➤ Temperature accuracy</li><li>➤ Temperature precision</li></ul>	<ul style="list-style-type: none"><li>• <i>Elution order</i></li><li>• <i>Peak identification</i></li><li>• <i>Elution order</i></li><li>• <i>RT precision</i></li><li>• <i>Peak identification</i></li></ul>	 <p>The figure displays two chromatograms side-by-side. The top chromatogram is labeled '40° C' and shows a series of peaks. Two peaks at the end of the run are labeled 'A' and 'B' in blue, with 'A' eluting before 'B'. The bottom chromatogram is labeled '65° C' and shows the same series of peaks. The two peaks at the end are labeled 'B' and 'A' in blue, with 'B' eluting before 'A'. The x-axis for both is 'Time in Minutes' with major ticks at 0, 5, and 10.</p>

# Role - HPLC UV/Vis Detectors

Performance Characteristics	<i>Influence on...</i>
<u>Variable Wavelength &amp; Diode Array Detector</u>	
<ul style="list-style-type: none"><li>➤ Low noise, wander, and drift</li><li>➤ Wide linear range</li><li>➤ Wavelength accuracy &amp; precision</li></ul>	<ul style="list-style-type: none"><li>• <i>Detection limit, quantitation limit</i></li><li>• <i>Quantitation at low &amp; high concentration</i></li><li>• <i>Accuracy &amp; precision of peak areas/heights</i></li></ul>
<u>Diode Array Detector Only</u>	
<ul style="list-style-type: none"><li>➤ Spectral resolution</li><li>➤ Spectral sensitivity</li></ul>	<ul style="list-style-type: none"><li>• <i>Accuracy of spectra, peak id by spectra</i></li><li>• <i>Accuracy of spectra, peak id by spectra at low concentrations</i></li></ul>

# Housekeeping

## Daily/System Start-up


- Mobile phase – how fresh is yours?
- Purge pumps ~5 minutes
- Condition pump for ~15 minutes

## Weekly

- Seal wash solvent
- Buffer flush
- Visual inspection of solvent filters
- Purge with composition of your application
- Condition with composition of your application



# Growth in H<sub>2</sub>O over 1 week in a clear bottle



# Avoid growth of algae



## Algae

The presence of algae in aqueous media can cause a variety of problems in your HPLC. Buffers, light and a pH between 4 and 8 will accelerate their growth.

For more information refer to your pump reference manual.

# Method & Column



- Method
  - ✓ Established or New
  - ✓ Conditions
- Column
  - ✓ Choice
  - ✓ Documentation
  - ✓ Guards
  - ✓ Equilibration



# Method Conditions

## Mobile Phase

- ✓ HPLC or MS grade solvents
- ✓ Buffer – right choice, column, LC, LC/MS, filtered?
- ✓ Preparation procedure
- ✓ Fresh mobile phase
- ✓ Bottles covered? No paraffin
  - Volatile components
- ✓ Label bottles, content + date
- ✓ Amber bottle for aqueous
- ✓ Make sure system flushed before introducing new mobile phase
- ✓ pH

## Temperature

## Pressure

## Standards, test mix

# Buffer Options

<b>Non-volatile:</b>		pK <sub>a</sub>	Buffer Range
Phosphate	$\text{H}_3\text{PO}_4 \rightleftharpoons \text{H}_2\text{PO}_4^-$	pK <sub>1</sub> = 2.1	1.1 – 3.1
	$\text{H}_2\text{PO}_4^- \rightleftharpoons \text{HPO}_4^{2-}$	pK <sub>2</sub> = 7.2	6.2 – 8.2
	$\text{HPO}_4^{2-} \rightleftharpoons \text{PO}_4^{3-}$	pK <sub>3</sub> = 12.3	11.3 – 13.3
Citrate	$\begin{array}{c} \text{CH}_2\text{COOH} \\   \\ \text{HOCCOOH} \\   \\ \text{CH}_2\text{COOH} \end{array}$	pK <sub>1</sub> = 3.1	2.1 – 4.1
		pK <sub>2</sub> = 4.7	3.7 – 5.7
		pK <sub>3</sub> = 5.4	4.4 – 6.4
Borate	$\text{H}_3\text{BO}_3$	pK <sub>1</sub> = 9.2	8.2 – 10.2
<b>Volatile:</b>			
Trifluoroacetate	$\text{F}_3\text{CCOOH}$	pK <sub>1</sub> = 0.5	xx – 1.5
Formate	$\text{HCOOH}$	pK <sub>1</sub> = 3.8	2.8 – 4.8
Acetate	$\text{CH}_3\text{COOH}$	pK <sub>1</sub> = 4.8	3.8 – 5.8
Ammonium	$\text{NH}_4^+$	pK <sub>1</sub> = 9.2	8.2 – 10.2

# Mobile Phase Preparation

➤ Small changes in mobile phase strength can have a large affect on retention

- ✓ HPLC grade or better
- ✓ Buffer prep procedure
  - Be consistent
    - Document process

Volume % of solvents can depend on preparation

Specified volume ACN added to a 1 L volumetric and made to volume with H<sub>2</sub>O

≠


Specified volume H<sub>2</sub>O added to a 1 L volumetric and made to volume with ACN

≠

500 ml H<sub>2</sub>O added to 500 ml ACN

- ✓ Degree of contraction is affected by the relative quantities of each
- ✓ Temperature

# Mobile Phase Preparation Effect on Chromatography




HPLC System: Agilent 1100 with quaternary pump  
Column: ZORBAX Eclipse XDB-C8 Rapid-Resolution (3.5 $\mu$ m), 4.6 x 50 mm  
Agilent Part No. 935967-906  
Mobile Phases: Dial-a-Mix= A: water B: MeOH, pump 50% B  
Vol/Vol=250 mL water + 250 mL MeOH, pump 100%  
Mix-to-Mark = 250 mL MeOH, fill to 500 mL with water, pump 100%  
Premixed (w/w) = 200 g MeOH + 200 g water, pump 100%  
Detection: UV 254 nm  
Flow: 1 mL/ min.  
Temperature: ambient

1. Uracil
2. Butylparaben
3. Naphthalene
4. Dipropylphthalate
5. Acenaphthene

- Method used to prepare MP can significantly affect the elution pattern
- **Be consistent**
  - w/w is more accurate than v/v

Effect of Mobile Phase Preparation on Chromatography,  
Pub. No. 5988-6476EN

# Changes in Volatile Buffer Concentration Can Cause Shifts in Retention Time and Peak Shape



# Column

- Choice
  - ✓ Conditions, flow rate, pressure, pH
- Performance report
- Datasheet or column guide
- Benchmark with your system
- Suitable fittings
- Equilibrated
- In-line filters and/or guards
- Sample
- Store properly when done

# Column Specifications

## InfinityLab Poroshell 120 column specifications

InfinityLab Poroshell 120	Pore Size	Temperature Limit	pH Range	Endcapped	Carbon Load	Surface Area	USP Designation
EC-C18	120 Å	60 °C	2.0–8.0	Yes	10%	130 m <sup>2</sup> /g	L1
EC-C8	120 Å	60 °C	2.0–8.0	Yes	5%	130 m <sup>2</sup> /g	L7
SB-C18	120 Å	90 °C	1.0–8.0	No	9%	130 m <sup>2</sup> /g	L1
SB-C8	120 Å	80 °C	1.0–8.0	No	5.5%	130 m <sup>2</sup> /g	L7
HPH-C18	100 Å	60 °C	3.0–11.0	Yes	Proprietary	95 m <sup>2</sup> /g	L1
HPH-C8	100 Å	60 °C	3.0–11.0	Yes	Proprietary	95 m <sup>2</sup> /g	L7
Bonus-RP	120 Å	60 °C	2.0–8.0	Yes	9.5%	130 m <sup>2</sup> /g	L60
PFP	120 Å	60 °C	2.0–8.0	Yes	5.1%	130 m <sup>2</sup> /g	L43
Phenyl-Hexyl	120 Å	60 °C	2.0–8.0	Yes	9%	130 m <sup>2</sup> /g	L11
SB-Aq	120 Å	80 °C	1.0–8.0	No	Proprietary	130 m <sup>2</sup> /g	L96
EC-CN	120 Å	60 °C	2.0–8.0	Yes	3.5%	130 m <sup>2</sup> /g	L10
HILIC-Z	100 Å	80 °C	2.0–12.0	No	Proprietary	95 m <sup>2</sup> /g	L114
HILIC	120 Å	60 °C	0.0–8.0	No	NA	130 m <sup>2</sup> /g	L3
HILIC-OH5	120 Å	45 °C	1.0–7.0	Proprietary	Proprietary	130 m <sup>2</sup> /g	L86
Chiral-V	120 Å	45 °C	2.5–7.0	Proprietary	Proprietary	130 m <sup>2</sup> /g	L88
Chiral-T	120 Å	45 °C	2.5–7.0	Proprietary	Proprietary	130 m <sup>2</sup> /g	L63
Chiral-CD	120 Å	45 °C	3.0–7.0	Proprietary	Proprietary	130 m <sup>2</sup> /g	L45
Chiral-CF	120 Å	45 °C	3.0–7.0	Proprietary	Proprietary	130 m <sup>2</sup> /g	NA

# Column Documentation

## Performance Report

SERIAL NUMBER: USDAZ01333

PART NUMBER: 959758-902

COLUMN TYPE: ZORBAX RRHD Eclipse Plus C18 2.1 x 100 mm, 1.8 µm


PACKING LOT #: B09089

### TEST CONDITIONS

MOBILE PHASE = 60% Acetonitrile / 40% Water  
COLUMN PRESSURE = 517.2 Bar  
COLUMN FLOW = 0.50 ml / min  
LINEAR VELOCITY = 0.436 cm / sec  
TEMPERATURE = AMBIENT (Nominally 23 °C)  
INJECTION VOLUME = 1 µl

### QUALITY CONTROL PERFORMANCE RESULTS FOR NAPHTHALENE

TEST VALUES	SPECIFICATIONS
THEORETICAL PLATES = 22337	MIN = 21000
SELECTIVITY = 1.90	RANGE = 1.82 - 1.92
USP TAILING FACTOR = 1.08 (@ 5% Peak Height)	RANGE = 0.98 - 1.20
k' = 4.58	



Sample components with concentrations diluted in mobile phase in the following elution order.

Peak #	Conc (ug/ml)	Sample Component
1	10	Uracil
2	400	Phenol
3	50	4-Chloro Nitrobenzene
4	80	Naphthalene

Manufacturing test chromatogram is done on a modified LC system to minimize ECV and will differ from a typical lab HPLC

- Don't expect to get the exact same result as the performance report
- Test column performance on your instrument to have as a reference

# Column Documentation

## Column Guide

This booklet provides general information for all ZORBAX, Poroshell, Pursuit, and Polaris reversed-phase columns.


For additional detailed information about your specific phase or family, see: [agilent.com/chem/columnchoices](http://agilent.com/chem/columnchoices)

### Getting Started

A QC Column Performance Report, including a test chromatogram, is enclosed with every Agilent column. The QC test system has been modified from a standard system to minimize system dead volume, so it may vary from the system used in your lab. This allows a better evaluation of the column and assures a more consistent product. A properly configured LC system will generate similar results to the chromatogram on your QC Performance Report.

Modern columns are robust and are designed to operate for long periods under normal chromatographic conditions. You can maximize column performance by running it within specifications. Always review the specifications before putting in place a final method.

## Data Sheet




### Agilent ZORBAX Eclipse Plus Phenyl-Hexyl Rapid Resolution Threaded Column

Data Sheet

#### General Description


Eclipse Plus Phenyl-Hexyl columns are designed for superior peak shape with basic compounds, and deliver high efficiency and excellent peak shape with all sample types. Eclipse Plus Phenyl-Hexyl is especially useful for the separation of acidic, basic, and other highly polar compounds by reverse-phase liquid chromatography. Eclipse Plus Phenyl-Hexyl packing is made by first chemically bonding a dense monolayer of dimethylphenylhexylsilane stationary phase to a specially prepared, improved ultra-high purity (>99.995% SiO<sub>2</sub>) ZORBAX Rx-SIL porous silica support. This special silica support (Type B) is designed to reduce or eliminate strong adsorption of basic and highly polar compounds. The bonded-phase packing is then doubly endcapped using proprietary reagents and procedures to obtain maximum deactivation of the silica surface. Eclipse Plus Phenyl-Hexyl columns can be used for acidic and neutral samples, but are especially suited for separating basic compounds that produce poor peak shapes on other columns. These columns can be used for a wide range of applications and over a pH range of 2 to 8, accommodating most popular mobile phases.



# Column Documentation - Benchmark

Benchmark new column on your system

1. Standard mix; test mix (5188-6529, 01080-68704; QC reference material;
2. Criteria like retention time, peak area, peak tailing, resolution, response, system pressure, etc.
3. Theoretical plates
  - Monitor column over time
  - Troubleshoot



## Chromatographic conditions

Sample: RRLC Checkout sample  
(p/n 5188-6529)  
Column: Agilent Poroshell 120  
EC C18, 3 mm x 50 mm,  
2.7 µm  
Mobile phase: A = Water  
B = Acetonitrile  
Gradient: 0 min 20% B  
8 min 80% B  
Flow rate: 1.2 mL/min  
Stop time: 8 min  
Post time: 4 min  
Injection volume: 1 µL  
Column temperature: 30 °C  
DAD: 245/10 nm  
Ref 400/100 nm  
Flow cell: 10 mm  
Peak width: <0.025 min (10 Hz)

# Initial Column and System Equilibration\*

In an appropriate vessel, test highest % organic/buffer ratio to verify that buffer will not precipitate. With stirring, add organic to buffer first, not vice versa.

Equilibrate column with, in order:

- 100% organic modifier (if brand new)
- mobile phase minus buffer
- buffered mobile phase containing highest % organic modifier (gradient high end)
- buffered mobile phase containing lowest % organic modifier (gradient low end).

Inject standard or sample several times until RTs stable, or for gradient methods, precede former with 1 or 2 blank gradients.

\*Or follow instructions in your column user guide

# Do You Need a Guard Column?

## If you do, is it the correct one?



The **ZORBAX High Performance Guard Cartridge** components assemble quickly and easily to provide a high efficiency, low dead volume guard column that seals, with hand tightening, **up to 340 bar or 200 bar with a PEEK fitting.**

For use with columns having a 5um, 3um or 3.5um packing and **400 bar pressure limit**



**Agilent Fast Guards** (3/pk), stainless steel **UHPLC** guards packed with **1.8um Zorbax** or **Poroshell 120** materials.

- Single replacement guard column (no cartridge)
- Rated to **600 bar – 1300 bar** to match column

# Checklist

## Supplies

---

- Critical supplies on hand

## Sample

---

- Is it ready for chromatography?

## Instrument

---

- Maintenance up to date

## Method conditions

---

## Column

---

- Right choice for your sample and conditions

## Final checklist

---

- Shutdown (see appendix)
  - Short term
  - Long Term

# Contact Agilent Chemistries and Supplies Technical Support



1-800-227-9770 Option 3, Option 3:

Option 1 for GC/GCMS Columns and Supplies

Option 2 for LC/LCMS Columns and Supplies

Option 3 for Sample Preparation, Filtration and QuEChERS

Option 4 for Spectroscopy Supplies



[gc-column-support@Agilent.com](mailto:gc-column-support@Agilent.com)

[lc-column-support@agilent.com](mailto:lc-column-support@agilent.com)

[spp-support@agilent.com](mailto:spp-support@agilent.com)

[spectro-supplies-support@agilent.com](mailto:spectro-supplies-support@agilent.com)

# Resources for Support

- Agilent University <http://www.agilent.com/crosslab/university>
- Tech support <http://www.agilent.com/chem/techsupport>
- Resource page <http://www.agilent.com/chem/agilentresources>
  - Quick Reference Guides
  - Catalogs, Column User guides
  - Online Selection Tools, How-to Videos
- InfinityLab Supplies Catalog ([5991-8031EN](http://www.agilent.com/chem/infinitylab))
- Your local FSE and Specialists
- Youtube – [Agilent Channel](https://www.youtube.com/agilent)
- Agilent Service Contracts



# APPENDIX

# Determining the Dwell Volume of Your System

**Replace column with short piece of HPLC stainless steel tubing**

**Prepare mobile phase components**

**A. water - UV-transparent**

**B. water with 0.2% acetone - UV-absorbing**


**Monitor at 265 nm**

**Adjust attenuation such that both 100% A and 100% B are on scale**

**Run gradient profile 0 - 100% B/10 min at 1.0 mL/min**

**Record**

# Measuring Dwell Volume



Intersection identifies dwell time (t<sub>D</sub>)


$$V_D = t_D \times F$$
$$V_D = \text{Dwell Volume}$$

# Measuring Dwell Volume

If using gradient conditions - report dwell volume ( $V_D$ )  
 $V_D$  varies from instrument to instrument

## Dwell Volume Impact

A chromatogram generated on one instrument ( $V_{D1}$ ) can have a very different profile if generated on another instrument ( $V_{D2}$ )



High Pressure Mixing:  $V_D$  = mixing chamber + connecting tubing + injector

Low Pressure Mixing:  $V_D$  = the above + pump heads + associated plumbing

001011P1.PPT

# Correcting for Dwell Volume

1. Measure the Dwell Volume of your HPLC System

$$V_D = 1.0 \text{ mL}$$

2. Draw Effective Gradient Profile at First Flow Rate  
Calculate the time delay (imposed isocratic hold)  
caused by dwell volume


$$V_D = t_D \cdot F \quad 1.0 \text{ mL} = t_D \cdot 1.0 \text{ mL / min}$$

where  $F = 1.0 \text{ mL / min}$  for 4.6 x 150 mm column

$$V_D = 1.0 \text{ mL}$$

$$t_D = F/V_D \quad t_D = 1.0 \text{ mL / min} / 1.0 \text{ mL}$$

$$t_D = 1.0 \text{ min}$$



000884P2.PPT

# Correcting for Dwell Volume

If  $V_{D1} > V_{D2}$

Compensate for longer  $V_{D1}$  by adding an isocratic hold to  $V_{D2}$ , such that  
Hold +  $V_{D2} = V_{D1}$

If  $V_{D1} < V_{D2}$

Delay injection, such that  $V_{D2} - \text{delay} = V_{D1}$

( very difficult to accomplish in practice )

# Shutdown State and Instrument Flushing

## Shutdown State

Next day use—using same buffers

- Pump mobile phase very slowly (for example, 0.01 – 0.1 mL/min).

When flushing column or for longer term column storage

- Flush with 20/80 organic/water, then 80/20 organic/water or 100% organic.

## Instrument flushing

- ✓ Replace column with capillary tubing. Leave disconnected from detector.
- ✓ Flush pumps with water, then connect capillary tubing to detector.
- ✓ Inject water 2-3 times at maximum injection volume setting.
- ✓ Flush all pumps with 100% organic for long term storage

➤ Check your instrument manual for manufacturer's guidance

# Buffer Preparation – General Guidance

1. Dissolve salt in organic-free water in 1- or 2-L beaker. Use appropriate volume to leave room for pH adjustment solution. Equilibrate solution to room temperature for maximum accuracy.
2. Calibrate pH meter. Use 2-level calibration and bracket desired pH. Use appropriate audit solution to monitor statistical control (for example, potassium hydrogen tartrate, saturated solution, pH = 3.56).
3. Adjust salt solution to desired pH. Minimize amount of time electrode spends in buffer solution (contamination). Avoid overshoot and readjustment (ionic strength differences can arise).
4. Transfer pH-adjusted buffer solution quantitatively to volumetric flask, dilute to volume, and mix.
5. Filter through 0.45  $\mu\text{m}$  filter. Discard first 50 – 100 mL filtrate. Rinse solvent reservoir with small volume of filtrate and discard. Fill reservoir with remaining filtrate or prepare premix with organic modifier.
  - Agilent Solvent Filtration Kit, 250-mL reservoir, 1000-mL flask, p/n 3150-0577
  - Nylon filter membranes, 47 mm, 0.45  $\mu\text{m}$  pore size, p/n 9301-0895