

# Extraction of a Full Suite of Semivolatile Compounds from Drinking Water Using Automated Solid Phase Extraction

Melissa Lever, Biotage.

Key Words: US EPA 525.2, SPE, solid phase extraction, drinking water.

## Introduction

Drinking water is an important part of environmental exposure, especially for small children. Countries around the world have put regulations in place to monitor drinking water quality for a wide range of hazardous compounds. Methods such as SL 392- 2007 in China, the EN methods in Europe and US methods such as method 525.2 cover a large suite of analytes of concern, extracted using solid phase extraction (SPE) disks and using GC/MS for detection.

US EPA method 525.2 may use solid phase extraction (SPE) to extract the analytes of interest from water samples. It includes a variety of quality control measures to ensure the method is under control throughout the analysis.<sup>1</sup> The Biotage® Horizon 5000 (previously known as SPE-DEX 5000) was used in this study to extract the EPA Method 525.2 analytes from six prepared water samples as described in section 9.3, Initial Demonstration of Laboratory Accuracy and Precision. The test involves measuring 4-7 samples of reagent water spiked at approximately the mid-point of the calibration curve, 2-5 µg/L of the full suite of analytes. For each analyte and surrogate, the spike recovery, expressed as a percentage of the true value, should be 70–130% and the relative standard deviation (RSD) should be <30% to meet method criteria. All six samples were extracted using the same procedure and calibration parameters.

### » Biotage

- » Biotage® Horizon 5000 (previously known as SPE-DEX 5000)
- » DryVap® Automated In-line drying and concentration system
- » Drydisk® Separation membranes
- » Atlantic® C18 High Capacity Disk



The Biotage® Horizon 5000 was used for extraction of the analytes from the water samples. The Biotage® Horizon 5000 is an automated system that conditions the solid phase extraction disk, loads the sample through the disk, rinses the sample bottle and elutes the sample all without user intervention. The 47 mm diskholder was used and high capacity C18 disks were chosen. High capacity disks were used because some of the more water soluble compounds are better retained and suffer less breakthrough with this disk. Ethyl acetate and methylene chloride were used for elution.



Biotage® Horizon 5000 (previously known as SPE-DEX 5000) Automated Extractor.



The method run on the Biotage® Horizon 5000 is shown in Table 1. The method information, run time, sample identification and other information are stored in a file when the sample is run. This can be printed in a report or exported to the laboratory LIMS for archiving.

## Method Summary

1. Obtain six 1 liter samples of drinking water.
2. Add dechlorinating agent to each 1 liter sample.
3. Acidify each 1 liter water sample to pH <2 using concentrated HCl.
4. Add surrogate and internal standard compounds into samples.
5. Start extraction method shown in Table 1 and collect extract (≈32 mL).
6. Add extract to the DryDisk® holder and start automated drying and concentration process on the DryVap system (the DryVap® system automatically dries and concentrates extract to 0.9 mL).
7. Quantitatively bring extract volume to 1.0 mL using methylene chloride once the extracts are evaporated to less than 1 mL.
8. Add external standard into the 1 mL extract.
9. Transfer the extract to a 2.0 mL GC vial.
10. Analyze by GC/MS.

**Table 1.** Extraction Method.

Step	Solvent	Solvent Volume (mL)	Purge Time (s)	Pump Rate (#)	Sat. Time (s)	Soak Time (s)	Drain Time (s)	
1. Condition SPE Disk	Methylene Chloride	15	60	2	1	20	30	
2. Condition SPE Disk	Ethyl Acetate	11	60	2	1	20	30	
3. Condition SPE Disk	Methanol	11	60	2	1	60	2	
4. Condition SPE Disk	Reagent Water	9	30	2	1	5	5	
5. Condition SPE Disk	Reagent Water	9	60	2	1	30	0	
Step	Sample Flow Rate (#)		Done Loading Sample Delay (s)					
6. Load Sample	2		45					
Step	Dry Time (s)		Pump Rate (#)		N <sub>2</sub> Blanket			
7. Air Dry Disk Timer	60		6		Off			
Step	Solvent	Solvent Volume (mL)	Purge Time (s)	Pump Rate (#)	N <sub>2</sub> Blanket	Sat. Time (s)	Soak Time (s)	Elute Time (s)
8. Elute Sample Container	Ethyl Acetate	8	60	2	Off	1	30	45
9. Elute Sample Container	Methylene Chloride	8	15	2	Off	1	30	45
10. Elute Sample Container	Methylene Chloride	8	15	2	Off	1	30	45
11. Elute Sample Container	Methylene Chloride	8	15	6	Off	2	30	60

## Gas Chromatography Mass Spectrometry System

Column ZB Semi-volatiles, 30 m x 0.35 mm i.d., 0.25 µm film thickness (Phenomenex)

Flow Rate 9 psig helium ramped up with the oven temperature to maintain a constant flow

Temperature Ramp

Temperature (°C)	Rate (°C/min)	Hold (min)
60	0	2.00
270	20	0.00
320	6	3.00

The GC/MS used was a 6890 GC with a 5973 MSD (Agilent).

Total Run Time: 23.83 minutes

Injection Method: 1.0 µL injected, Temperature 280°C, Pulsed splitless

- Inlet pulse pressure 25.0 psi for 1.00 min
- Purge flow to split vent 50 mL/min for @2.00 min

## Results and Discussions

Six replicate laboratory fortified blanks (LFBs) were extracted as described in US EPA Method 525.2, following the procedure in the method summary in this note. Drinking water was spiked with standards and surrogates at a concentration of 5 µg/L.

The results are shown in Table 2 for each of the six replicate samples.

**Table 2.** Recoveries of spiked analytes in drinking water.

Compound	Sample 1 (% Rec)	Sample 2 (% Rec)	Sample 3 (% Rec)	Sample 4 (% Rec)	Sample 5 (% Rec)	Sample 6 (% Rec)	AVG	SD
Acenaphthene d10	70.0	73.0	77.6	81.0	75.0	76.8	75.6	3.82
Phenanthrene d10	75.6	84.2	88.0	95.4	84.4	85.8	85.6	6.41
Chrysene d12	77.0	83.8	90.2	98.2	84.0	86.0	86.5	7.14
2,4-Dinitrotoluene	113	112	109	110	112	109	111	1.60
2,6-Dinitrotoluene	113	112	108	109	113	107	110	2.75
2-Nitro-m-xylene	95.6	95.6	90.0	87.0	98.0	94.6	93.5	4.12
4,4'-DDD	94.4	93.8	94.0	90.4	94.0	93.6	93.4	1.48
4,4'-DDE	91.4	89.8	92.8	87.2	90.4	89.6	90.2	1.88
4,4-DDT	94.4	93.8	94.0	90.4	94.0	93.6	93.4	1.48
a-BHC	99.2	95.2	97.8	93.8	98.4	96.6	96.8	2.05
Acenaphthene	110	109	101	105	109	106	107	3.28
Acenaphthylene	95.6	95.6	93.8	90.2	98.0	97.0	95.0	2.76

Compound	Sample 1 (% Rec)	Sample 2 (% Rec)	Sample 3 (% Rec)	Sample 4 (% Rec)	Sample 5 (% Rec)	Sample 6 (% Rec)	AVG	SD
Acetochlor	108	109	104	103	107	108	106	2.24
a-Chlordane	92.8	91.4	93.4	88.8	93.0	92.0	91.9	1.68
Alachlor	99.8	97.0	99.0	95.2	98.8	98.4	98.0	1.67
Aldrin	85.0	81.6	84.6	81.2	82.4	83.6	83.1	1.58
Ametryn	98.6	95.2	96.8	93.8	98.0	96.4	96.5	1.77
Atrazine	104	98.0	100	96.6	98	100	99.4	2.47
b-BHC	99.0	96.8	98.6	96.0	99.6	97.0	97.8	1.43
Benz(a)anthracene	89.2	87.6	89.0	85.4	88.6	88.8	88.1	1.44
Benzo(a)pyrene	73.8	74.0	76.2	70.4	74.8	75.4	74.1	2.02
Benzo(b)fluoranthene	95.0	93.8	95.4	89.0	92.0	92.0	92.9	2.38
Benzo(ghi)perylene	95.0	93.0	95.2	91.4	93.0	93.2	93.5	1.42
Benzo(k)fluoranthene	91.0	88.8	91.8	89.0	92.2	92.4	90.9	1.60
Bis(2-ethylhexyl)adipate	95.4	92.6	93.4	90.8	94.4	94.2	93.5	1.61
Bis(2-ethylhexyl)phthalate	96.0	94.0	95.0	92.4	94.8	95.4	94.6	1.26
Bromacil	105	103	104	100	104	102	103	1.86
Butaclor	98.0	98.4	97.6	92.6	100	97.4	97.4	2.54
Butyl benzyl phthalate	100	99.0	98.6	93.4	101	99.0	98.5	2.65
Butylate	103	101	101	97.0	104	101	101	2.44
Caffeine	82.4	88.8	77.2	78.0	84.8	79.0	81.7	4.51
Chlorobenzilate	101	101	99.2	95.4	104	102	100	2.85
Chloroneb	112	110	114	111	113	109	111	1.81
Chlorothalonil	107	105	106	101	107	105	105	2.00
Chlorpropham	125	122	119	122	122	121	122	1.90
Chlorpyrifos	101	95.0	101	96.2	98.6	98.2	98.3	2.41
Chrysene	92.8	90.8	92.8	89.2	91.8	91.6	91.5	1.36
cis-Permethrin	97.8	96.4	96.4	94.4	95.8	96.0	96.1	1.10
Cyanazine	104	99.4	101	98.2	101	100	101	1.94
Cycloate	113	111	110	110	113	110	111	1.46
Dacthal	99.0	96.4	98.6	95.2	99.6	97.8	97.8	1.68
d-BHC	99.6	97.6	101	96.0	100.2	97.4	98.6	1.83
Diazinon	92.8	89.8	90.2	88.0	91.6	91.0	90.6	1.65
Dibenz(ah)anthracene	92.2	91.6	92.2	92.8	91.6	95.6	92.7	1.51
Dichlorvos	119	115	110	112	119	115	115	3.57
Dieldrin	95.6	95.2	95.0	90.2	95.8	93.6	94.2	2.12
Diethyl phthalate	115	114	113	112	114	111	113	1.52
Dimethoate	81.0	90.0	72.8	77.2	82.2	80.0	80.5	5.72
Dimethyl phthalate	109	107	107	106	109	105	107	1.41

Compound	Sample 1 (% Rec)	Sample 2 (% Rec)	Sample 3 (% Rec)	Sample 4 (% Rec)	Sample 5 (% Rec)	Sample 6 (% Rec)	AVG	SD
Di-n-butyl phthalate	101	98.8	100	95.0	99.6	98.2	98.9	2.18
Di-n-octyl phthalate	99.6	96.0	98.6	94.8	97.6	97.4	97.3	1.73
Diphenamid	101	98.0	98.8	95.4	100	97.6	98.5	2.01
Disulfoton	96.4	92.8	96.2	91.8	97.4	97.4	95.3	2.42
Disulfoton sulfone	102	103	99.6	97.4	105	102	102	2.68
Endosulfan I	97.8	92.8	97.0	93.4	97.0	93.4	95.2	2.26
Endosulfan II	96.6	97.6	98.8	92.2	98.6	95.2	96.5	2.49
Endosulfan Sulfate	99.2	99.8	98.8	95.4	101	100	99.1	2.04
Endrin	115	111	110	105	114	111	111	3.47
Endrin Aldehyde	88.8	88.4	84.4	81.6	89.4	84.8	86.2	3.10
Endrin Ketone	97.6	97.6	98.8	91.8	101	98.2	97.5	3.11
EPTC	107	104	99.4	99.4	107	104	103	3.35
Ethoprop	119	119	112	116	118	115	117	2.66
Etridiazole	107	107	105	103	110	106	106	2.41
Fenamiphos	116	119	108	109	119	117	115	5.09
Fenarimol	104	104	96.6	96.8	102	103	101	3.31
Fluoranthene	98.2	95.2	97.8	93.8	96.6	95.2	96.1	1.70
Fluorene	104	104	102	100	105	103	103	1.78
Fluridone	110	113	103	104	112	110	109	3.97
g-Chlordane	90.6	89.4	90.8	88.2	91.4	90.2	90.1	1.14
Heptachlor	96.4	91.0	95.6	92.6	94.2	94.8	94.1	1.99
Heptachlor epoxide A	96.8	92.2	95.2	91.0	93.8	92.4	93.6	2.14
Heptachlor epoxide B	96.2	95.2	98.0	92.6	96.8	94.0	95.5	1.96
Hexachlorobenzene	92.0	87.2	91.0	87.2	89.8	91.0	89.7	2.06
Hexazinone	100	100	98.0	94.6	101	100	99.1	2.47
Indeno(1,2,3-cd)pyrene	94.0	94.4	95.2	91.0	92.6	93.4	93.4	1.48
Isophorone	109	101	96.6	97.4	107	100	102	5.00
Lindane (g-BHC)	99.4	95.6	99.8	95.2	99.4	96.4	97.6	2.12
Malathion	117	116	108	110	115	117	114	3.91
Merphos	93.4	90.6	104	101	108	108	101	7.39
Methoxychlor	100	99.8	98.6	94.2	99.2	97.8	98.3	2.18
Methyl paraoxon	98.8	97.2	97	93.8	94.6	93.4	95.8	2.17
Metolachlor	102	99.4	100	97	102	99.4	100	1.98
Metribuzin	100	96.4	96.2	94.6	99.4	97	97.3	2.16
Mevinphos	127	124	116	122	126	122	123	3.99
MGK-264-A	96.0	94.6	95.8	92.6	97.2	96.0	95.4	1.59
MGK-264-B	96.0	94.6	95.8	92.6	97.2	96.0	95.4	1.59

Compound	Sample 1 (% Rec)	Sample 2 (% Rec)	Sample 3 (% Rec)	Sample 4 (% Rec)	Sample 5 (% Rec)	Sample 6 (% Rec)	AVG	SD
Molinate	111	109	108	107	111	107	109	1.78
Naphthalene	89.6	89.8	86.8	80.2	92.0	90.4	88.1	4.24
Napropamide	102	101	98.6	94.2	103	100	99.8	3.09
Norflurazon	99.6	100	96.8	93.8	102	99.6	98.7	2.89
Pebulate	108	106	104	102	109	105	105	2.64
Pentachlorophenol	106	104	105	102	108	107	105	2.08
Perylene-d12	79.2	78.4	81.4	76.6	79.4	80.6	79.3	1.69
Phenanthrene	97.4	94.2	96.0	91.4	97.4	95.8	95.4	2.28
Prometon	79.2	76.6	79.8	73.4	75.2	77.6	77.0	2.42
Prometryn	99.4	95.6	97.4	94.4	97.8	97.0	96.9	1.75
Pronamide	100	96.6	98.8	95.4	99.4	98.0	98.0	1.75
Propachlor	115	112	112	112	113	110	112	1.54
Propazine	103	100	100	98	101	100	101	1.77
Pyrene	96.2	95.2	95.8	90.8	96.4	95.4	95.0	2.09
Pyrene-d10	95.8	94.6	95.6	90.2	96.8	93.8	94.5	2.33
Simazine	98.8	95.8	98.0	95.0	97.4	94.0	96.5	1.86
Simetryn	99.8	97.0	102	99.6	98.4	102	99.9	2.07
Stirafos	112	111	105	104	114	114	110	4.36
Tebuthiuron	120	121	110	118	120	118	118	3.83
Terbacil	119	117	112	111	118	116	115	3.42
Terbufos	122	114	121	110	113	106	115	6.21
Terbutylazine	101	98.4	98.0	96.0	99.4	98.0	98.5	1.79
Terbutryn	101	98.4	100	96.2	100.6	99.6	99.3	1.73
Terphenyl-d14	130	119	111	102	119	116	116	9.26
Thiobencarb	101	99.6	95.4	96.8	98.6	99.8	98.6	2.13
trans-Nonachlor	92.6	92.0	92.6	89.4	92.8	91.2	91.8	1.30
trans-Permethrin	94.4	92.0	94.2	90.0	93.2	92.4	92.7	1.63
Triademefon	102	99.4	97.2	92.6	98.6	97.0	97.8	3.18
Tricyclazole	99.6	98.4	88.8	89.6	97.6	95.2	94.9	4.63
Trifluralin	110	108	107	107	108	107	108	1.20
Triphenylphosphate	101	101	99.0	95.0	102	101	99.8	2.57
Vernolate	106	105	101	100	109	104	104	3.22

The recoveries are within the 70—130% recovery as specified in the method quality control criteria in all cases for the more than 100 compounds measured. The precision was excellent and the standard deviations ranged for 2-3% for most of the analytes.

## Conclusion

The Atlantic® high-capacity C18 disks provided excellent recovery of the large suite of compounds extracted in water. The compounds included in the method had excellent performance and an average recovery of 98.8% was achieved with a 5 µg/L spike. The spike recovery criterion of 70–130% was achieved in all cases. The Biotage® Horizon 5000 system provided uniform performance and a hands-off approach to the extraction step. The reproducibility of the six runs was excellent and the average of the standard deviation values was 2.57%. The Atlantic high-capacity C18 disk allowed even more water soluble compounds, such as caffeine, to be successfully retained with good recovery, demonstrating the utility of solid phase extraction. In combination with the Biotage® Horizon, the samples were reliably extracted with excellent precision.

## References

1. Method 525.2, Determination of Organic Compounds in Drinking Water by Liquid-Solid Extraction and Capillary Column Gas Chromatography/Mass Spectrometry, US EPA [https://www.epa.gov/sites/production/files/2015-10/documents/method\\_525-2\\_rev-2\\_1995.pdf](https://www.epa.gov/sites/production/files/2015-10/documents/method_525-2_rev-2_1995.pdf), accessed April 12, 2017.

### EUROPE

Main Office: +46 18 565900  
Toll Free: +800 18 565710  
Fax: +46 18 591922  
Order Tel: +46 18 565710  
Order Fax: +46 18 565705  
order@biotage.com  
Support Tel: +46 18 56 59 11  
Support Fax: + 46 18 56 57 11  
eu-1-pointsupport@biotage.com

### NORTH & LATIN AMERICA

Main Office: +1 704 654 4900  
Toll Free: +1 800 446 4752  
Fax: +1 704 654 4917  
Order Tel: +1 704 654 4900  
Order Fax: +1 434 296 8217  
ordermailbox@biotage.com  
Support Tel: +1 800 446 4752  
Outside US: +1 704 654 4900  
us-1-pointsupport@biotage.com

### JAPAN

Tel: +81 3 5627 3123  
Fax: +81 3 5627 3121  
jp\_order@biotage.com  
jp-1-pointsupport@biotage.com

### CHINA

Tel: +86 21 68162810  
Fax: +86 21 68162829  
cn\_order@biotage.com  
cn-1-pointsupport@biotage.com

### KOREA

Tel: + 82 31 706 8500  
Fax: + 82 31 706 8510  
korea\_info@biotage.com  
kr-1-pointsupport@biotage.com

### INDIA

Tel: +91 22 4005 3712  
india@biotage.com

Distributors in other regions  
are listed on [www.biotage.com](http://www.biotage.com)

### Literature Number: AN114-HOR.V.1

© 2019 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at [www.biotage.com/legal](http://www.biotage.com/legal). Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.