

Mass spectrometry

High precision copper and zinc isotopes with Neoma MS/MS MC-ICP-MS



Authors

G. Craig, M. Pfeifer,
J. Roberts, C. Bouman

Thermo Fisher Scientific,
Bremen, Germany

Keywords

MS/MS MC-ICP-MS, Pre-cell mass filter, Copper, Zinc, Mass bias, Neoma

Goal

To demonstrate the ability of Thermo Scientific™ Neoma™ MS/MS MC-ICP-MS to measure copper and zinc isotopes at high precision.

Introduction

Copper and zinc are important metals that exhibit distinct isotopic fractionations as a result of various physical, chemical and biological processes. As a result, copper and zinc isotopes have been widely used in a range of applications, varying from cosmochemistry to medical metallomics. In this application note, we discuss the various application spaces for copper and zinc isotopes and show the capabilities of the Neoma MS/MS MC-ICP-MS for measuring these isotope systems at high precision.

Copper and zinc are important elements for studying metal cycling in the lithosphere, atmosphere, biosphere and hydrosphere. Cu is a siderophile and highly chalcophilic element with moderate volatility. It is the most refractory of the chalcophile elements. As a multivalent element in nature, Cu isotopes respond to changes in redox. Zn is a slightly siderophile and moderately volatile element. It is somewhat lithophile in nature because it often occurs in silicate, forming strong bonds with oxygen. Zinc is not multivalent in nature and thus changes in redox will not fractionate Zn isotopes in a geological context. Within the geosciences discipline, Cu and Zn isotopes have been used to trace planetary formation^{1,2}, core-mantle differentiation³, mineralization and chemical weathering^{4,5}, and biogeochemical cycling^{6,7}.

Copper and zinc are also useful tools for studying metal dispersion in the environment. Their isotope ratios have been used as a tracer of anthropogenic Cu and Zn contamination in surface and aquatic environments⁸⁻¹⁰. Trace metals are released by anthropogenic activities (e.g. agriculture, industry, mining, metallurgy, energy production and transport) and dispersed among the different natural compartments such as soil, water and atmosphere. To develop correct environmental management strategies, either for appropriately targeted remediation actions or to maintain emissions at sustainable levels, it is imperative to identify their sources and understand their transport and fate in the environment.

Copper isotope ratios in groundwater have also been used as a mineral exploration vector. Relatively large surface and groundwater $\delta^{65}\text{Cu}$ dispersion haloes can exist around mineral deposits^{11,12}.

Cu and Zn are promising elements in the biomedical field of isotope metallomics because of their small number of functional roles in biology and also because their turnover rate in the body is relatively short. Copper plays a major role in oxidizing iron and controlling electron fluxes, while Zn is a cofactor of hundreds of important enzymes. Cu and Zn isotope ratios have been used to identify multiple cancers and other diseases¹³⁻¹⁶.

Methods

Experimental set-up

Here we present 3 analytical runs of Cu and Zn solutions (Table 1). Samples are measured via standard-sample bracketing relative to itself.

The cup configuration for all analytical runs is shown in Table 2.

Table 1. Experimental set-up of three analytical runs for Cu and Zn isotope analysis

Analytical run	1	2	3
Solution	100 ppb Zn and 100 ppb Cu in 3% HNO ₃	100 ppb Zn and 100 ppb Cu in 3% HNO ₃	1 ppm Zn and 1 ppm Cu in 3% HNO ₃
Plasma	Wet	Dry	Wet
Sample introduction system	SIS spray chamber for Neoma MC-ICP-MS	ESI® Apex Omega™ Q desolvating nebulizer system	SIS spray chamber for Neoma MC-ICP-MS
Nebulizer	ESI 100 µL/min microFAST MC™ syringe injection autosampler	ESI 100 µL/min microFAST MC syringe injection autosampler	ESI PFA 100 µL min
Cones	X skimmer Ni cone, Jet sampler Ni cone	X skimmer Ni cone, Jet sampler Ni cone	X skimmer Ni cone, Jet sampler Ni cone
Resolution	Low	Medium	Low
Method	75 cycles of 4 second integration time (5 min total measurement time)	45 cycles of 4 second integration time (3 minute total measurement time)	50 cycles of 8 second integration time (7 min total measurement time)

Table 2. Cup configuration for Cu and Zn isotope analysis

L5	L4	L3	L2	L1	C	H1	H2	H3	H4	H5
⁶² Ni	⁶³ Cu	⁶⁴ Zn	⁶⁵ Cu	¹³¹ Xe ⁺⁺	⁶⁶ Zn	⁶⁷ Zn	¹³⁵ Ba ⁺⁺	⁶⁸ Zn	¹³⁸ Ba ⁺⁺	⁷⁰ Zn
10 ¹¹ Ω	10 ¹¹ Ω	10 ¹¹ Ω	10 ¹¹ Ω	10 ¹¹ Ω	10 ¹¹ Ω	10 ¹¹ Ω				

The instrument settings are shown in Table 3. Note that these experiments could equally have been carried out with the base Neoma MC-ICP-MS (rather than Neoma MS/MS MC-ICP-MS). In all runs, the pre-cell mass filter was in full ion transmission mode and there was no gas within the collision/reaction cell.

Table 3. Instrument settings for 3 analytical runs described in Table 1

Analytical run	1	2	3
RF power [W]	1,200	1,200	1,200
Cool gas [L min ⁻¹]	14.0	14.0	14.0
Aux gas [L min ⁻¹]	0.798	0.798 (0.8)	0.798
Neb gas [L min ⁻¹]	1.057	1.062	1.059
Desolvator settings	–	Ar – 4.6 L/min, N ₂ – 5 mL/min	–
Wien filter electric field [V]	206.5	213	213
Wien magnetic field [%]	30	30	30
Wien center mass	⁶⁶ Zn	⁶⁶ Zn	⁶⁶ Zn
Pre-filter slit position [%]	95	80	80
CCT entry [V]	-50.5	-50.4	-50.4
CCT bias [V]	-0.1	-0.1	-0.1
RF amplitude [%]	96.8	96.8	96.8
CCT exit 1 [V]	-180	-180	-180
CCT exit 2 [V]	-73.1	-73.1	-73.1
CCT gas	None	None	None

Interference corrections were applied automatically with Thermo Scientific™ Qtegra™ Intelligent Scientific Data Solution (ISDS) Software (Table 4). All isotope ratios were mass bias corrected using ⁶⁶Zn/⁶⁴Zn.

Table 4. Interference corrections applied to each isotope of interest

Isotope of interest	Interferent	Measured interference-free isotope used for correction
⁶⁴ Zn	¹²⁸ Xe ²⁺	¹³¹ Xe ²⁺
	⁶⁴ Ni ⁺	⁶² Ni ⁺
⁶⁵ Cu	¹³⁰ Xe ²⁺	¹³¹ Xe ²⁺
	¹³⁰ Ba ²⁺	¹³⁸ Ba ²⁺
⁶⁶ Zn	¹³² Xe ²⁺	¹³¹ Xe ²⁺
⁶⁸ Zn	¹³⁶ Xe ²⁺	¹³¹ Xe ²⁺

Results

Sensitivity

The absolute sensitivity of the Neoma MS/MS MC-ICP-MS is shown in Table 5. In low resolution and wet plasma, the sensitivity for 118–125 V/ppm and 46–48 V/ppm for copper and zinc respectively. For copper, this is almost an order of magnitude higher than sensitivities expressed for other available MC-ICP-MS in the market (ca. 13–34 V/ppm for Cu^{17,18}). For zinc, the sensitivity of Neoma MS/MS MC-ICP-MS is a factor of 3 better than other available MC-ICP-MS at low resolution in wet plasma (16 V/ppm)¹⁹.

In dry plasma and medium resolution, the Neoma MS/MS MC-ICP-MS out-performs other available MC-ICP-MS in the market. The sensitivity of Neoma MS/MS MC-ICP-MS for copper in medium resolution was approximately double the sensitivity achieved with other MC-ICP-MS in low resolution (90–100 V/ppm)²⁰. For zinc, the sensitivity of Neoma MS/MS MC-ICP-MS in medium resolution was 197 V/ppm, more than an order of magnitude better than other MC-ICP-MS instruments also measuring in medium resolution (17 V/ppm²¹), and more than double what can be achieved in low resolution (70–80 V/ppm²²).

Table 5. Sensitivity of each analytical run for copper and zinc

Analytical run	1 (Low resolution, wet plasma)	2 (Medium resolution, dry plasma)	3 (Low resolution, wet plasma)
Cu sensitivity (V/ppm)	125	208	118
Zn sensitivity (V/ppm)	48	197	46

Mass bias

Mass bias of Neoma MS/MS MC-ICP-MS follows the exponential law (Figure 1). All samples plot on a straight line in log-log space and the gradient of the slopes match the predicted value based on exponential mass bias law.




Figure 1. Log-log plots of interference-corrected samples. The expected slopes based on exponential law are (A) 0.9113, (B) 1.4244, and (C) 1.6958.

Cu isotope reproducibility

Replicate measurements of a copper/zinc solution are shown in Figure 2. The performance of Neoma MS/MS MC-ICP-MS is the same under wet plasma in low resolution and dry plasma in medium resolution – this is because the experiments have been set-up such that they produce the same number of total counts across the measurement. Internal precision is similar to external reproducibility in all experiments. At higher concentration (1 ppm Cu/Zn), the external reproducibility was 11.5 ppm (2SD) and 3.45 ppm (2SE).

Zn isotope reproducibility

A similar picture is observed for zinc isotopes, with the external reproducibility of the 1,000 ppb (Analytical run 3) of 17.5 ppm and 22.7 ppm (2SD) for $\delta^{67}\text{Zn}$ and $\delta^{68}\text{Zn}$, respectively (Figure 3). The standard errors are 4.9 ppm and 6.3 ppm (2SE) for $\delta^{67}\text{Zn}$ and $\delta^{68}\text{Zn}$, respectively.




Figure 2. Replicate measurements of copper/zinc solution. $\delta^{65}\text{Cu}$ has been mass bias corrected using $^{66}\text{Zn}/^{64}\text{Zn}$. (Blue) Analytical run 1 – wet plasma, low resolution, 100 ppb Cu/Zn, (Red) Analytical run 2 – dry plasma, medium resolution, 100 ppb Cu/Zn, (Gray) Analytical run 3 – wet plasma, low resolution, 1,000 ppb Cu/Zn.




Figure 3. Replicate measurements of copper/zinc solution. $\delta^{67}\text{Zn}$ and $\delta^{68}\text{Zn}$ have been mass bias corrected using $^{66}\text{Zn}/^{64}\text{Zn}$. (Blue) Analytical run 1 – wet plasma, low resolution, 100 ppb Cu/Zn, (Red) Analytical run 2 – dry plasma, medium resolution, 100 ppb Cu/Zn, (Gray) Analytical run 3 – wet plasma, low resolution, 1,000 ppb Zn.




Figure 3 continued. Replicate measurements of copper/zinc solution. $\delta^{67}\text{Zn}$ and $\delta^{68}\text{Zn}$ have been mass bias corrected using $^{66}\text{Zn}/^{64}\text{Zn}$. (Blue) Analytical run 1 – wet plasma, low resolution, 100 ppb Cu/Zn, (Red) Analytical run 2 – dry plasma, medium resolution, 100 ppb Cu/Zn, (Gray) Analytical run 3 – wet plasma, low resolution, 1,000 ppb Zn.

Conclusion

Neoma MS/MS MC-ICP-MS in full ion transmission mode can be used to measure copper and zinc isotopes at ppm-level precision. Instrumental sensitivity is more than double that achievable by other MC-ICP-MS instruments available. Mass bias is predictable, following the exponential law. This allows copper and zinc isotopes to be accurately corrected for mass bias, producing accurate as well as high precision isotope data.

References

- Herzog, G. F. *et al.* Isotopic and elemental abundances of copper and zinc in lunar samples, Zagami, Pele's hairs, and a terrestrial basalt. *Geochim Cosmochim Acta* **73**, 5884–5904 (2009).
- Moynier, F. *et al.* Isotopic composition of zinc, copper, and iron in lunar samples. *Geochim Cosmochim Acta* **70**, 6103–6117 (2006).
- Savage, P. S. *et al.* Copper isotope evidence for large-scale sulphide fractionation during Earth's differentiation. *Geochem Perspect Lett* **1**, 53–64 (2015).
- Mathur, R. *et al.* Modern and paleofluid pathways revealed by Cu isotope compositions in surface waters and ores of the Pebble porphyry Cu-Au-Mo deposit, Alaska. *Economic Geology* **108**, 529–541 (2013).
- Asael, D. *et al.* Copper isotope fractionation in sedimentary copper mineralization (Timna Valley, Israel). *Chem Geol* **243**, 238–254 (2007).
- Conway, T. M. *et al.* The biogeochemical cycling of zinc and zinc isotopes in the North Atlantic Ocean. *Global Biogeochem Cycles* **28**, 1111–1128 (2014).
- Takano, S. *et al.* Isotopic constraints on biogeochemical cycling of copper in the ocean. *Nature Communications* 2014 5:1 **5**, 1–7 (2014).
- Araújo, D. F. *et al.* Assessment of the metal contamination evolution in the Loire estuary using Cu and Zn stable isotopes and geochemical data in sediments. *Mar Pollut Bull* **143**, 12–23 (2019).
- Viers, J. *et al.* Are Cu isotopes a useful tool to trace metal sources and processes in acid mine drainage (AMD) context? *Chemosphere* **193**, 1071–1079 (2018).
- Gelly, R. *et al.* Lead, zinc, and copper redistributions in soils along a deposition gradient from emissions of a Pb-Ag smelter decommissioned 100 years ago. *Science of The Total Environment* **665**, 502–512 (2019).
- Mathur, R. *et al.* Use of Cu isotopes to distinguish primary and secondary Cu mineralization in the Cañariaco Norte porphyry copper deposit, Northern Peru. *Miner Depos* **47**, 755–762 (2012).
- Kidder, J. A. *et al.* Hydrogeochemical mineral exploration in deeply weathered terrains: An example from Mumbwa, Zambia. *Science of The Total Environment* **810**, 151215 (2022).
- Mahan, B. *et al.* Isotope metallomics approaches for medical research. *Cellular and Molecular Life Sciences* 2020 77:17 **77**, 3293–3309 (2020).
- Albarede, F. *et al.* Medical applications of Cu, Zn, and S isotope effects. *Metallomics* **8**, 1056–1070 (2016).
- Hastuti, A. A. M. B. *et al.* High-precision isotopic analysis of serum and whole blood Cu, Fe and Zn to assess possible homeostasis alterations due to bariatric surgery. *Anal Bioanal Chem* **412**, 727–738 (2020).
- Costas-Rodríguez, M. *et al.* Isotopic analysis of Cu in blood serum by multi-collector ICP-mass spectrometry: a new approach for the diagnosis and prognosis of liver cirrhosis? *Metallomics* **7**, 491–498 (2015).
- Zhang, Y. *et al.* Copper Isotope Ratio Measurements of Cu-Dominated Minerals Without Column Chromatography Using MC-ICP-MS. *Front Chem* **8**, 609 (2020).
- Wang, J. *et al.* High-precision copper isotopic analysis using a Nu Sapphire MC-ICP-MS. *J Anal At Spectrom* **37**, 2589–2598 (2022).
- Sauzéat, L. *et al.* Inter-comparison of stable iron, copper and zinc isotopic compositions in six reference materials of biological origin. *Talanta* **221**, 121576 (2021).
- Larner, F. *et al.* A new separation procedure for Cu prior to stable isotope analysis by MC-ICP-MS. *J Anal At Spectrom* **26**, 1627–1632 (2011).
- Moore, R. E. T. *et al.* High Precision Zinc Stable Isotope Measurement of Certified Biological Reference Materials Using the Double Spike Technique and Multiple Collector-ICP-MS. *Anal Bioanal Chem* **409**, 2941–2950 (2017).
- Arnold, T. *et al.* Measurement of zinc stable isotope ratios in biogeochemical matrices by double-spike MC-ICPMS and determination of the isotope ratio pool available for plants from soil. *Anal Bioanal Chem* **398**, 3115 (2010).

Appendix

A1. Analytical run 1

	MBC		MBC		MBC	
	$\delta^{65}\text{Cu}$	2SE	$\delta^{67}\text{Zn}$	2SE	$\delta^{68}\text{Zn}$	2SE
	0.016	0.016	-0.048	0.036	0.018	0.029
	0.013	0.021	0.019	0.040	0.015	0.022
	0.027	0.018	0.037	0.036	0.000	0.027
	-0.006	0.017	-0.044	0.045	-0.039	0.021
	0.012	0.017	-0.023	0.043	0.106	0.033
	-0.012	0.017	0.012	0.044	-0.008	0.026
	-0.019	0.020	0.035	0.039	-0.077	0.031
	-0.002	0.022	-0.002	0.035	-0.045	0.031
	0.009	0.021	0.040	0.043	-0.042	0.023
	0.016	0.019	0.005	0.038	-0.048	0.028
	$\delta^{65}\text{Cu}$	2SE	$\delta^{67}\text{Zn}$	2SE	$\delta^{68}\text{Zn}$	2SE
Mean	0.005	0.019	0.003	0.040	-0.012	0.027
2SD	0.029	-	0.065	-	0.103	-

A2. Analytical run 2

	MBC		MBC		MBC	
	$\delta^{65}\text{Cu}$	2SE	$\delta^{67}\text{Zn}$	2SE	$\delta^{68}\text{Zn}$	2SE
	0.010	0.057	0.019	0.049	0.118	0.057
	0.009	0.034	0.007	0.037	-0.068	0.047
	-0.004	0.040	0.012	0.047	-0.010	0.052
	0.009	0.037	0.008	0.044	-0.021	0.044
	-0.018	0.043	0.038	0.053	0.162	0.057
	0.018	0.036	0.008	0.051	0.013	0.057
	0.013	0.047	-0.033	0.047	-0.057	0.055
	-0.001	0.033	0.006	0.046	0.089	0.052
	0.004	0.038	-0.001	0.051	-0.016	0.053
	-0.005	0.036	-0.006	0.040	0.000	0.047
	$\delta^{65}\text{Cu}$	2SE	$\delta^{67}\text{Zn}$	2SE	$\delta^{68}\text{Zn}$	2SE
Mean	0.004	0.040	0.006	0.047	0.021	0.052
2SD	0.021	-	0.036	-	0.153	-

A3. Analytical run 3

	MBC		MBC		MBC	
	$\delta^{65}\text{Cu}$	2SE	$\delta^{67}\text{Zn}$	2SE	$\delta^{68}\text{Zn}$	2SE
	-0.0077	0.0075	-0.0048	0.0119	0.0015	0.0122
	0.0063	0.0064	0.0171	0.0102	0.0069	0.0099
	-0.0038	0.0071	0.0122	0.0119	0.0139	0.0114
	0.0007	0.0079	0.0129	0.0104	0.0199	0.0100
	-0.0015	0.0089	0.0197	0.0130	0.0051	0.0095
	0.0109	0.0085	0.0086	0.0086	-0.0072	0.0084
	0.0086	0.0077	0.0113	0.0094	-0.0046	0.0100
	-0.0006	0.0062	-0.0060	0.0097	-0.0235	0.0085
	-0.0039	0.0065	0.0036	0.0098	0.0134	0.0094
	-0.0027	0.0094	0.0159	0.0112	0.0103	0.0119
	-0.0001	0.0087	0.0012	0.0098	-0.0041	0.0092
	-	-	0.0041	0.0129	-0.0044	0.0111
	-	-	0.0206	0.0086	0.0049	0.0104
	$\delta^{65}\text{Cu}$	2SE	$\delta^{67}\text{Zn}$	2SE	$\delta^{68}\text{Zn}$	2SE
Mean	0.0006	0.0077	0.0090	0.0106	0.0025	0.0101
2SD	0.0115	-	0.0175	-	0.0227	-

 Learn more at thermofisher.com/ms-ms

General Laboratory Equipment - Not For Diagnostic Procedures. © 2023 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified ESI and Apex Omega are trademarks of Elemental Scientific Inc. This information is presented as an example of the capabilities of Thermo Fisher Scientific products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representatives for details. **AN002360-EN 0723C**