

APPLICATIONS

A Fast Isothermal GC Analysis of Spark Ignition Fuel Oxygenate Additives Using an Alternative Phase

Tim Nelson
Phenomenex, Inc., 411 Madrid Avenue., Torrance, CA 90501 USA



Tim Nelson
Global Marketing Manager
Fuels and Chemicals

Tim is an avid outdoorsman who loves to hike and ski. His most recent exploration is tall ship sailing in our local Pacific Ocean. Tim loves history and everything about the stars and space.

Introduction

Oxygenates are chemical additives used to enhance the quality of spark ignition fuels and to boost the performance of engines. They are effective modifiers to improve acceleration of the engine which can also lower the vehicle's fuel mileage and help keep engines clean for longer while performing more efficiently. Oxygenates help gasoline burn efficiently during combustion to provide optimum torque to the vehicle. Another benefit to an efficient combustion burn is that it lowers harmful carbon monoxide levels and other harmful emissions from vehicles.

Below is a list of some of the common oxygenates found in spark ignition fuels:

Alcohols

- Methanol (MeOH)
- Ethanol (EtOH)
- Isopropyl alcohol (IPA)
- n-Butanol (BuOH)
- Gasoline grade tert-butanol (GTBA)

Ethers

- Methyl tert-butyl ether (MTBE)
- Tert-Amyl methyl ether (TAME)
- Tert-Hexyl methyl ether (THEME)
- Ethyl tert-butyl ether (ETBE)
- Tert-Amyl ethyl ether (TAEE)
- Diisopropyl ether (DIPE)


Historically MTBE had been used as the primary oxygenate additive in gasoline (petrol), from approximately 1974 as a replacement for lead up to the year 2000 when itself became an environmental concern. In place of MTBE, ethanol and ETBE are now commonly used as oxygenate additives in spark ignition fuels, though some regions of the world still use MTBE. Methanol is also being used in some proposed new low emission formulations.

A 100% dimethylpolysiloxane phase GC column is typically used for the analysis of oxygenates however this technical work describes an alternative selectivity, specifically a Zebron ZB-5ms 5% phenyl-arylene phase.

GC Method Parameters

Column: Zebron™ ZB-5ms
Phase: 5 % Phenyl-Arylene 95 % Dimethylpolysiloxane
Dimensions: 30 meter x 0.25 mm x 0.50 µm
Part No.: 7HG-G010-17
Injection: Split 100:1 @ 200 °C, 0.2 µL
Carrier Gas: Helium @ 1.4 mL/min (Constant Flow)
Oven Program: 30 °C Isothermal
Detector: Flame Ionization (FID) @ 250 °C
Sample: 1. Methanol
 2. Ethanol
 3. 2-Propanol
 4. tert-Butanol
 5. n-Propanol
 6. Methyl tert-butyl ether (MTBE)
 7. Isopropyl ether
 8. sec-Butanol
 9. tert-Butyl ethyl ether (ETBE)
 10. Methylcyclopentane
 11. Isobutanol
 12. tert-Pentanol
 13. 1,2-Dimethoxyethane
 14. Benzene
 15. n-Butanol
 16. tert-Amyl methyl ether

Figure 1.
Separation of 16 Common Spark Ignition Fuel Oxygenates



App ID 15434

Results and Discussion

The Zebron™ ZB-5ms has a special Engineered Self Cross-linking™ (ESC) stationary phase that provides low bleed when used to separate oxygenates. The column's low polar 5 % phenyl-arylene phase provided great separation for all the test components, with a fast sub-7 minute analysis time @ 30 °C isothermal run.

Alcohols in general are more polar than ethers, so the slightly polar Zebron™ ZB-5ms helped separate the more polar elements from the nonpolar ones, for example there was good separation of tert-pentanol from the methylcyclopentane, which can coelute when run on a 100% dimethylpolysiloxane phase.

Conclusion

The Zebron ZB-5ms provided optimal separation of critical pairs with symmetric peaks and a fast run time of under 7 minutes at an isothermal temperature program. The Zebron ZB-5ms with the 5% phenyl-arylene phase is a great alternative selectivity to provide increased separation of polar and aromatic elements from the less polar components.

Ordering Information

Zebron ZB-5ms GC Columns			
ID(mm)	df(μm)	Temp. Limits °C	Part No.
10-Meter			
0.10	0.10	-60 to 325/350	7CB-G010-02
0.18	0.18	-60 to 325/350	7CD-G010-08
12-Meter			
0.20	0.33	-60 to 325/350	7DE-G010-14
15-Meter			
0.25	0.25	-60 to 325/350	7EG-G010-11
20-Meter			
0.18	0.18	-60 to 325/350	7FD-G010-08
0.18	0.32	-60 to 325/350	7FD-G010-51
0.18	0.36	-60 to 325/350	7FD-G010-53
25-Meter			
0.20	0.33	-60 to 325/350	7GE-G010-14
30-Meter			
0.25	0.25	-60 to 325/350	7HG-G010-11
0.25	0.50	-60 to 325/350	7HG-G010-17
0.25	1.00	-60 to 325/350	7HG-G010-22
0.32	0.25	-60 to 325/350	7HM-G010-11
0.32	0.50	-60 to 325/350	7HM-G010-17
0.32	1.00	-60 to 325/350	7HM-G010-22
60-Meter			
0.25	0.10	-60 to 325/350	7KG-G010-02
0.25	0.25	-60 to 325/350	7KG-G010-11
0.32	0.25	-60 to 325/350	7KM-G010-11

Note: If you need a 5 in. cage, simply add a (-B) after the part number, e.g., 7HG-G010-11-B. Some exceptions may apply. Agilent 6850 and some SRI and process GC systems use only 5 in. cages.

Find Your Liner Online!

Easily search by part numbers, applications, injection mode, or system manufacturer for quick selection in under 1 minute!

www.phenomenex.com/FindLiner



Inlet Consumables Are Available Online!

Need inlet seals, septa, or syringes? Explore hundreds of available parts online at:

www.phenomenex.com/InletGC



BE-HAPPY™
guarantee

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms & Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Zebron and Engineered Self Cross-linking (ESC) are trademarks of Phenomenex.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2019 Phenomenex, Inc. All rights reserved.