

Capillary GC Analyses of Chlorinated Pesticides in Dry Pet Food

Based on federal and state regulations for identifying and quantifying low levels of pesticides in food and environmental samples, we selected three capillary columns to screen for chlorinated pesticides. A nonpolar and two low/intermediate polarity phases were chosen to evaluate differences in component elution order and retention times. Low level screening analyses were performed effectively by using split/splitless injection and an electron capture detector (ECD). Each column separated the 18 pesticides in approximately 35 minutes. Example chromatograms are shown.

Key Words:

• pesticides • aldrin • endrin • methoxychlor • pet food

Federal and state regulations require that pesticides in food and environmental samples be identified and quantified at low levels. Based on these regulations, we selected three capillary GC columns to screen for pesticides at low levels. We chose a nonpolar phase, PTE™-5, and two low/intermediate polarity phases, SPB™-608 and SPB-1701, to illustrate the differences in component elution order and retention times for chlorinated pesticides.

Eighteen chlorinated pesticides (Table 1) were spiked into, then extracted from, dry pet food purchased at a local grocery store. Extracts were prepared by weighing out 50 grams of pet food, manually crushing it, and adding 100mL of acetonitrile (1). *The untreated extracts contained no pesticides.*

Table 1. Chlorinated Pesticide Standards Mixture
(TCL Pesticide Mix, Cat. No. 4-8913, 2000µg/mL each component in toluene:hexane, 50:50)

Aldrin
α-BHC
β-BHC
γ-BHC
δ-BHC
4,4'-DDD
4,4'-DDE
4,4'-DDT
Dieldrin
Endosulfan I
Endosulfan II
Endosulfan sulfate
Endrin
Endrin aldehyde
Endrin ketone
Heptachlor
Heptachlor epoxide
Methoxychlor

Samples of the spiked and unspiked extracts were injected onto each capillary column under the conditions listed in Figure A. Figure A shows chromatograms of the extracted pesticides from each column. The low/intermediate polarity SPB-608 and SPB-1701 columns selectively eluted the analytes, based on dipole-dipole and hydrogen bonding interactions between the solute and the stationary phase. Each column separated all of the analytes. The nonpolar PTE-5 column eluted the analytes by boiling point. One pair of analytes (endosulfan sulfate and 4,4'-DDT) coeluted from the nonpolar column. Analysis time for each column was 31-36 minutes. Elution order differences among the three columns are evident in Figure A. In a two-column analysis, these differences can be used to help confirm the identities of the analytes.

In general, recovery of the chlorinated pesticides was poor. The extraction procedure apparently was unable to remove these analytes from the fat in the pet food.* Table 2 lists the recovery values for the pesticides, determined using the SPB-608 column. Values ranged from 5% to 122%. The limits of detection for the standards ranged from 0.1 to 10.0 ppm. These values should improve with a more effective extraction procedure.

Based on these evaluations, we determined that the three stationary phases, PTE-5, SPB-608, and SPB-1701, exhibited differences in retention times, resolution, and elution order for 18 common chlorinated pesticides. Using these columns, screening analyses for low levels of these pesticides can be performed effectively, in approximately 35 minutes, with split/splitless injection and an electron capture detector (ECD). Pesticide recovery can be affected by the fat content of the sample, thus development of an effective extraction procedure is important.




Table 2. Recovery of Chlorinated Pesticides from Dry Pet Food (SPB-608 Column)

Pesticide	Recovery (%)
Aldrin	122
α-BHC	17
β-BHC	14
γ-BHC	5
δ-BHC	84
4,4'-DDD	23
4,4'-DDE	91
4,4'-DDT	61
Dieldrin	51
Endosulfan I	65
Endosulfan II	19
Endosulfan sulfate	—
Endrin	60
Endrin aldehyde	—
Endrin ketone	—
Heptachlor	6
Heptachlor epoxide	41
Methoxychlor	—

Figure A. Chlorinated Pesticides from Dry Pet Food

Stationary Phases: PTE-5, SPB-608, SPB-1701
 Column Dimensions: 30m x 0.25mm ID, 0.25µm phase film
 Catalog Nos.: 24135-U (PTE-5), 24103-U (SPB-608), 24113 (SPB-1701)
 Oven: 100°C to 280°C at 6°C/min
 Carrier: helium, 40cm/sec.
 Det: ECD, 250°C
 Sample: 1µL of 0.1µg/mL extract, split/splitless 45 sec, 250°C

- | | |
|-------------|--|
| 1. Aldrin | 10. Endosulfan I |
| 2. α-BHC | 11. Endosulfan II |
| 3. β-BHC | 12. Endosulfan sulfate |
| 4. γ-BHC | 13. Endrin |
| 5. δ-BHC | 14. Endrin aldehyde |
| 6. 4,4'-DDD | 15. Endrin ketone |
| 7. 4,4'-DDE | 16. Heptachlor |
| 8. 4,4'-DDT | 17. Heptachlor epoxide |
| 9. Dieldrin | 18. Methoxychlor
(not observed with PTE-5 column) |



795-0848, 0849, 0850

Ordering Information:

Description	Cat. No.
Fused Silica Capillary Columns	
all 30m x 0.25mm ID, 0.25µm phase film	
PTE-5	24135-U
SPB-608	24103-U
SPB-1701	24113
Chlorinated Pesticides Mixture (TCL Pesticides Mix)	
2000µg/mL each analyte in Figure A in toluene:hexane, 50:50	48913

Reference

1. Kaphalia, B.S. *Assoc. Official Analytical Chemists* **73** (4) 1990.

* For a suitable extraction procedure, refer to AOAC Methods, 16th edition (Order from AOAC International, 481 North Frederick Avenue, Suite 500, Gaithersburg, Maryland 20877-2504 USA. Tel.: +1-301-924-7077; FAX: +1-301-924-7089.)

PTE and SPB are trademarks of Supelco, Inc.

Fused silica columns manufactured under HP US Pat. No. 4,293,415.

Contact our Technical Service Department
 (phone 800-359-3041 or 814-359-3041, FAX 814-359-5468)
 for expert answers to your questions.

NOTE 91

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website (www.sigma-aldrich.com), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA · Sigma-Aldrich de Argentina, S.A. · Buenos Aires 1119 AUSTRALIA · Sigma-Aldrich Pty. Ltd. · Castle Hill NSW 2154 AUSTRIA · Sigma-Aldrich Handels GmbH · A-1110 Wien
 BELGIUM · Sigma-Aldrich N.V./S.A. · B-2880 Bornem BRAZIL · Sigma-Aldrich Quimica Brasil Ltda. · 01239-010 São Paulo, SP CANADA · Sigma-Aldrich Canada, Ltd. · 2149 Winston Park Dr., Oakville, ON L6H 6J8
 CZECH REPUBLIC · Sigma-Aldrich s.r.o. · 186 00 Praha 8 DENMARK · Sigma-Aldrich Denmark A/S · DK-2665 Vallensbaek Strand FINLAND · Sigma-Aldrich Finland/YA-Kemia Oy · FIN-00700 Helsinki
 FRANCE · Sigma-Aldrich Chimie · 38297 Saint-Quentin-Fallavier Cedex GERMANY · Sigma-Aldrich Chemie GmbH · D-82041 Deisenhofen GREECE · Sigma-Aldrich (o.m.) Ltd. · Ilioupoli 16346, Athens
 HUNGARY · Sigma-Aldrich Kft. · H-1067 Budapest INDIA · Sigma-Aldrich Co. · Bangalore 560 048 IRELAND · Sigma-Aldrich Ireland Ltd. · Dublin 24 ISRAEL · Sigma Israel Chemicals Ltd. · Rehovot 76100
 ITALY · Sigma-Aldrich s.r.l. · 20151 Milano JAPAN · Sigma-Aldrich Japan K.K. · Chuo-ku, Tokyo 103 KOREA · Sigma-Aldrich Korea · Seoul MALAYSIA · Sigma-Aldrich (M) Sdn. Bhd. · Selangor
 MEXICO · Sigma-Aldrich Química S.A. de C.V. · 50200 Toluca NETHERLANDS · Sigma-Aldrich Chemie BV · 3330 AA Zwijndrecht NORWAY · Sigma-Aldrich Norway · Torshov · N-0401 Oslo
 POLAND · Sigma-Aldrich Sp. z o.o. · 61-663 Poznań PORTUGAL · Sigma-Aldrich Quimica, S.A. · Sintra 2710 RUSSIA · Sigma-Aldrich Russia · Moscow 103062 SINGAPORE · Sigma-Aldrich Pte. Ltd.
 SOUTH AFRICA · Sigma-Aldrich (pty) Ltd. · Jet Park 1459 SPAIN · Sigma-Aldrich Quimica, S.A. · 28100 Alcobendas, Madrid SWEDEN · Sigma-Aldrich Sweden AB · 135 70 Stockholm
 SWITZERLAND · Supelco · CH-9471 Buchs UNITED KINGDOM · Sigma-Aldrich Company Ltd. · Poole, Dorset BH12 4QH
 UNITED STATES · Supelco · Supelco Park · Bellefonte, PA 16823-0048 · Phone 800-247-6628 or 814-359-3441 · Fax 800-447-3044 or 814-359-3044 · email: supelco@sial.com

H

Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.