

APPLICATIONS

Detailed Hydrocarbon Analysis of Refinery Naphtha Streams for Determination of Compound Group-types

Dr. Ramkumar Dhandapani and Tim Nelson
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA



Tim Nelson
Global Marketing Manager
Fuels and Chemicals

Tim is an avid outdoorsman who loves to hike and ski. His most recent exploration is tall ship sailing in our local Pacific Ocean. Tim loves history and everything about the stars and space.

Introduction

Petroleum fuels are required to meet federal and state environmental regulations to mitigate their impact on the environment. In the United States, the refineries use test methods that are approved by the EPA, and in other countries there are similar national regulations to ensure that the required performance criteria are met. There are different regulated test methods to measure the physical characteristics of the fuel as well as its specific chemical components. Depending on global location, spark ignition fuel is called gasoline, petrol, or benzin, and it is distilled from crude oil as part of the naphtha stream (C1-C12). Straight-run product is blended with other processed streams which makes it critical to have a detailed understanding of the hydrocarbon components for a consistent product.

There are hundreds of chemical components in spark ignition fuels, and so to make the testing and regulation manageable they are measured in categories by chemical group-type. These category group-types are based on molecular similarities that give them common performance behaviors. These tests are referred to as Detailed Hydrocarbon Analysis (DHA) which includes several test methods that measure different hydrocarbon subset or group-type categories. DHA identifies the important individual components in spark ignition fuels and other refinery light hydrocarbon process streams. Common refinery naphthas are virgin, alkylate, reformate, FCC and Coker, and finished spark ignition fuels. This technique is sometimes also referred to as PONA, PIONA, or PIANO analysis since specific test methods are measuring the group-type paraffins, isoparaffins, olefins, naphthenes, aromatics, and oxygenates. The DHA analysis helps optimize the production processes in addition to meeting the regulatory requirements.

Many of the DHA test methods utilize gas chromatography to separate and measure the group-type components. This technical work describes the use of a Zebtron™ ZB-DHA-PONA GC column which contains a 100 % Dimethylpolysiloxane phase that is specially designed to provide excellent peak symmetry and resolution of hydrocarbon components for a true boiling point separation of complex mixtures in spark ignition fuels.


GC Conditions for All Analyses

Column: Zebtron™ ZB-DHA-PONA
Phase: 100 % Dimethylpolysiloxane
Dimensions: 100 meter x 0.25 mm x 0.50 µm
Part No.: [7MG-G042-17](#)
Injection: Split 40:1 @ 300 °C, 0.2 µL
Recommended Liner: Zebtron^{PLUS} Straight Z-Liner™
Carrier Gas: Helium @ 1.55 mL/min (Constant Flow)
Oven Program: 35 °C for 14 min, to 60 °C at 1.1 °C/min. for 19 min, to 280 °C at 2 °C/min. for 5 min
Detection: Flame Ionization (FID) @ 320 °C

Table 1. Zebtron ZB-DHA-PONA Column Dimensional Options for Detailed Hydrocarbon Analysis (ASTM and alternative methods)

Part Number	Dimensions	Stationary Phase	Benefits
7JE-G042-17	50 meter x 0.2 mm x 0.50 µm	100 % Dimethylpolysiloxane	50 meter column provides shorter run time while the 0.2 µm tighter ID provides higher efficiency.
7MG-G042-17	100 meter x 0.25 mm x 0.50 µm	100 % Dimethylpolysiloxane	100 meter length provides high plate count/efficiency.
7QG-G042-22	150 meter x 0.25 mm x 1.0 µm	100 % Dimethylpolysiloxane	150 meter length with a 1.0 µm thicker film provides better separation of lower boiling fractions and maintains high efficiency from the column length.
7AG-G042-22	5 meter x 0.25 mm x 1.0 µm (tuning column)	5 % Phenyl, 95 % Polydimethylsiloxane	Optional tuning column provides phenyl selectivity in addition to true boiling point separation. This helps in resolving certain aromatics from alkanes and alkenes.


Figure 1.
Separation of Paraffins on a Zebron™ ZB-DHA-PONA GC Column



- Sample:**
1. Pentane
 2. Hexane
 3. Heptane
 4. Octane
 5. Nonane
 6. Decane
 7. Undecane
 8. Dodecane
 9. Tridecane
 10. Tetradecane
 11. Pentadecane

App ID 25306


Figure 2.
Separation of Isoparaffins on a Zebron ZB-DHA-PONA GC Column



- Sample:**
- | | |
|----------------------------|-------------------------|
| 1. Isopentane | 19. 4-Methylheptane |
| 2. 2,3-Dimethylbutane | 20. 3-Methylheptane |
| 3. 2-Methylpentane | 21. 3-Ethylhexane |
| 4. 3-Methylpentane | 22. 2,5-Dimethylheptane |
| 5. 2,2-Dimethylpentane | 23. 3,3-Dimethylheptane |
| 6. 2,4-Dimethylpentane | 24. 2,3-Dimethylheptane |
| 7. 2,2,3-Trimethylbutane | 25. 3,5-Dimethylheptane |
| 8. 3,3-Dimethylpentane | 26. 3,4-Dimethylheptane |
| 9. 2-Methylhexane | 27. 2-Methyloctane |
| 10. 2,3-Dimethylpentane | 28. 3,3-Diethylpentane |
| 11. 3-Methylhexane | 29. 3-Methyloctane |
| 12. 3-Ethylpentane | 30. 2,2-Dimethyloctane |
| 13. 2,2-Dimethylhexane | 31. 3,3-Dimethyloctane |
| 14. 2,5-Dimethylhexane | 32. 2,3-Dimethyloctane |
| 15. 2,2,3-Trimethylpentane | 33. 3-Ethylheptane |
| 16. 2,4-Dimethylhexane | 34. 2-Methylnonane |
| 17. 2,3-Dimethylhexane | 35. 3-Methylnonane |
| 18. 2-Methylheptane | |

App ID 25307


Figure 3.
Separation of Aromatic Compounds on a Zebron ZB-DHA-PONA GC Column



- Sample:**
- | | |
|---------------------------------|-----------------------------------|
| 1. Benzene | 20. 1-Methyl-3-n-propylbenzene |
| 2. Toluene | 21. 1-Methyl-4-n-propylbenzene |
| 3. Ethylbenzene | 22. n-Butylbenzene |
| 4. m-Xylene | 23. 1-Methyl-2-n-propylbenzene |
| 5. p-Xylene | 24. 1,2-Diethylbenzene |
| 6. o-Xylene | 25. 1,4-Dimethyl-2-ethylbenzene |
| 7. Isopropylbenzene | 26. 1,3-Dimethyl-5-ethylbenzene |
| 8. n-Propylbenzene | 27. 1,2-Dimethyl-4-ethylbenzene |
| 9. 1-Methyl-3-ethylbenzene | 28. 1,3-Dimethyl-2-ethylbenzene |
| 10. 1-Methyl-4-ethylbenzene | 29. 1,2-Dimethyl-3-ethylbenzene |
| 11. 1,3,5-Trimethylbenzene | 30. 1,2,4,5-Tetramethylbenzene |
| 12. 1-Methyl-2-ethylbenzene | 31. 2-Methylbutylbenzene |
| 13. 1,2,4-Trimethylbenzene | 32. trans-1-Butyl-2-methylbenzene |
| 14. tert-Butylbenzene | 33. n-Pentylbenzene |
| 15. Isobutylbenzene | 34. t-1-Butyl,3,5-dimethylbenzene |
| 16. sec-Butylbenzene | 35. t-1-butyl-ethylbenzene |
| 17. 1-Methyl-3-isopropylbenzene | 36. 1,3,5-Triethylbenzene |
| 18. 1-Methyl-4-isopropylbenzene | 37. 1,2,4-Triethylbenzene |
| 19. 1-Methyl-2-isopropylbenzene | 38. n-Hexylbenzene |


App ID 25309

Figure 4.
Separation of Naphthenes on a Zebron ZB-DHA-PONA GC Column



- Sample:**
1. Cyclopentane
 2. Methylcyclopentane
 3. Cyclohexane
 4. 1,1-Dimethylcyclopentane
 5. cis-1,3-Dimethylcyclopentane
 6. trans-1,3-Dimethylcyclopentane
 7. trans-1,2-Dimethylcyclopentane
 8. Methylcyclohexane
 9. Ethylcyclopentane
 10. ctc-1,2,4-Trimethylcyclopentane
 11. ctc-1,2,3-Trimethylcyclopentane
 12. cct-1,2,4-Trimethylcyclopentane
 13. trans-1,4-Dimethylcyclohexane
 14. 1-Ethyl-1-methylcyclopentane
 15. trans-1,2-Dimethylcyclohexane
 16. ccc-1,2,3-Trimethylcyclopentane
 17. Isopropylcyclopentane
 18. cis-1,2-Dimethylcyclohexane
 19. n-Propylcyclopentane
 20. ccc-1,3,5-Trimethylcyclohexane
 21. 1,1,4-Trimethylcyclohexane
 22. ctt-1,2,4-Trimethylcyclohexane
 23. ctc-1,2,4-Trimethylcyclohexane
 24. 1,1,2-Trimethylcyclohexane
 25. Isobutylcyclopentane
 26. Isopropylcyclohexane
 27. n-Butylcyclopentane
 28. Isobutylcyclohexane
 29. t-1-Methyl-2-propylcyclohexane
 30. t-1-Methyl-2-(4MP)cyclopentane

Figure 5.
Separation of Olefins on a Zebron ZB-DHA-PONA GC Column



- Sample:**
1. 3-Methyl-1-butene
 2. 1-Pentene
 3. 2-Methyl-1-butene
 4. 2-Methyl-1,3-butadiene
 5. trans-2-Pentene
 6. cis-2-Pentene
 7. 4-Methylpentene-1
 8. 1-Hexene
 9. 2-Methylpentene-2
 10. trans-2-Hexene
 11. cis-2-Hexene
 12. 1-Heptene
 13. trans-3-Heptene
 14. cis-3-Heptene
 15. trans-2-Heptene
 16. cis-2-Heptene
 17. 1-Octene
 18. trans-2-Octene
 19. cis-2-Octene
 20. 1-Nonene
 21. trans-3-Nonene
 22. cis-3-Nonene
 23. trans-2-Nonene
 24. cis-2-Nonene
 25. 1-Decene

APPLICATIONS

Results and Discussion

Detailed hydrocarbon analysis involves analysis of a matrix that can contain a variety of compounds including paraffins, isoparaffins, olefins, naphthenes, and related compounds which makes the separation complex. Presented in this tech note is the analysis of PIONA compounds by Zebtron ZB-DHA-PONA column which is specially designed for the separation of detailed hydrocarbons. The polyimide layer provides temperature stability and flexibility. The middle ultra clean silica layer provides great peak shape for polar and nonpolar compounds. The inner layer of a 100 % dimethylpolysiloxane that delivers great separation of the paraffins, isoparaffins, olefins, naphthalenes, and aromatics (**Figure 6**).

The Zebtron ZB-DHA-PONA column has a special Engineered Self Cross-linkingTM (ESC) stationary phase that provides low bleed and exceptional column life. The dimensions used in this test (100 meter length x 0.25 mm ID x 0.50 μ m film thickness) are important to provide the ideal separation and resolution. The smaller 0.25mm diameter along with the column length provides an efficient mass transfer of the complex analytes for great separation and peak shapes. The column's low polar 100 % Dimethylpolysiloxane phase provided great separation for all of the DHA test components (**Figure 1-5**). In Table 1. there is a list of additional column dimensions that can be used for other test methods which will provide optimal separation of detailed hydrocarbons in the refinery naphtha process streams.

Conclusion

The Zebtron ZB-DHA-PONA GC column provides the optimal separation of DHA critical pairs with symmetric peaks. The Zebtron ZB-DHA-PONA GC column with the 100 % Dimethylpolysiloxane phase is a great selectivity to provide separation of paraffins, isoparaffins, aromatics, naphthalene, and olefins.

Ordering Information

ZebtronTM ZB-DHA-PONA GC Columns

ID (mm)	df (μ m)	Temp. Limits °C	Part No.
5-Meter			
0.25	1.00	-60 to 340/360	7AG-G042-22
50-Meter			
0.20	0.50	-60 to 360/370	7GE-G042-17
100-Meter			
0.25	0.50	-60 to 360/370	7MG-G042-17
150-Meter			
0.25	1.00	-60 to 340/360	7QG-G042-22



Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy


Terms and Conditions
Subject to Phenomenex Standard Terms & Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks
Zebtron, Z-Liner, and Engineered Self Crosslinking (ESC) are trademarks of Phenomenex.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2019 Phenomenex, Inc. All rights reserved.

Figure 6.
Benefits of ZB-DHA-PONA GC Column



Australia
t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria
t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium
t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada
t: +1 (800) 543-3681
info@phenomenex.com

China
t: +86 400-606-8099
cninfo@phenomenex.com

Denmark
t: +45 4824 8048
nordicinfo@phenomenex.com

Finland
t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France
t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany
t: +49 (0)6021-58830-0
anfrage@phenomenex.com

India
t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Ireland
t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy
t: +39 051 6327511
italiainfo@phenomenex.com

Luxembourg
t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico
t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands
t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand
t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway
t: +47 810 02 005
nordicinfo@phenomenex.com

Portugal
t: +351 221 450 488
ptinfo@phenomenex.com

Singapore
t: +65 800-852-3944
sginfo@phenomenex.com

Spain
t: +34 91-413-8613
espinfo@phenomenex.com

Sweden
t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland
t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan
t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

United Kingdom
t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA
t: +1 (310) 212-0555
info@phenomenex.com

All other countries
Corporate Office USA 
t: +1 (310) 212-0555
info@phenomenex.com



 ...breaking with traditionSM

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com