



Halogenated hydrocarbons

Analysis of chlorinated volatiles and C₁-C₆ hydrocarbons

Application Note

Environmental

Authors

Agilent Technologies, Inc.

Introduction

Gas chromatography with an Agilent CP-SilicaPLOT column separates nine chlorinated volatiles and hydrocarbons in 14 minutes.

Chloromethane (methyl chloride) and vinylchloride are separated well from butane and can be quantified accurately because of the high retention, good separation and peak shape that the CP-SilicaPLOT column produces.



Agilent Technologies


Conditions

Technique : GC-capillary
Column : Agilent CP-SilicaPLOT, 30 m x 0.32 mm, fused silica
PLOT (df = 4 μ m) (Part no. CP8567)
Temperature : 30 °C (2 min) \rightarrow 250 °C, 10 °C/min
Carrier Gas : He, 210 kPa (2.1 bar, 30 psi)
Injector : Split, 1:100
T = 200 °C
Detector : FID
T = 300 °C
Sample Size : 1.0 mL
Concentration Range : 100 ppm in N₂

Courtesy : Jim Luong, The Dow Chemical Company, Canada

Peak identification

1. methane
2. ethane
3. propane
4. butane
5. chloromethane (methyl chloride)
6. vinylchloride
7. pentane
8. chloroethane (ethyl chloride)
9. hexane



www.agilent.com/chem

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A01289



Agilent Technologies