

Solvents

Analysis of water in acetone

Application Note

Materials Testing & Research

Authors

Agilent Technologies, Inc.

Introduction

Gas chromatography with an Agilent PorapLOT Q-HT column identifies water in an acetone sample in six minutes.



Agilent Technologies


Conditions

Technique : GC-wide-bore
Column : Agilent PoraPLOT Q-HT, 0.53 mm x 25 m fused silica
PLOT PoraPLOT Q-HT (df = 20 µm)
(Part no. CP7559)
Temperature : 150 °C
Carrier Gas : He, 150 kPa (1.5 bar, 21.7 psi)
Injector : Split
Detector : µ-TCD
Concentration : 4000 ppm in acetone

Peak identification

1. water
2. acetone

Water elutes very quickly from the hydrophobic porous polymers.
Water can be analyzed directly without using titration techniques.



www.agilent.com/chem

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A01301



Agilent Technologies