

# Temperature Stability and Cooling Performance of Double-Shot Pyrolyzer (PY-2020D)

In analysis with pyrolysis gas chromatography (Py-GC), because the pyrolysis temperature of a sample greatly influences the pyrogram obtained, the temperature control of the pyrolysis furnace is one of the most important performances. Fig. 1 shows a temperature profile of pyrolysis furnace of Double-Shot Pyrolyzer. When the set temperature of pyrolysis furnace was 600°C, the temperature fluctuation was 600+/-0.5°C, while the set temperature was 40°C, it was 40+/-1°C.

Also, in Evolved Gas Analysis in which heating and cooling are repeated, the cooling rate of the pyrolysis furnace greatly influences the efficiency in continuous analysis. Double-Shot Pyrolyzer employs a forced cooling of the pyrolysis furnace with cooling gas (nitrogen or air), and can achieve a rapid cooling from 600°C down to 40°C within 20 min.

1) Double-Shot Pyrolyzer® Technical Note, PYT-004




Fig. 1. Structure of Double-Shot Pyrolyzer




Fig. 2 Temperature Profile of Pyrolysis Furnace

**Keywords :** Heating Furnace, Cooling Rate, Temperature Control

**Products used :** Multi-functional pyrolyzer

**Applications :** General Polymer Analysis

**Related technical notes :**

Please forward your inquiries via our web page or send us a fax message.

**R&D and manufactured by :**  
**Frontier Laboratories Ltd.**

Phone: (81)24-935-5100 Fax: (81)24-935-5102  
<http://www.frontier-lab.com/>