



**Axel Semrau®**

## CHRONECT LC-GC **CHRONECT Workstation PAH**



**Product information**

---

## PAH analytics according to the imPAHct method

Polycyclic aromatic hydrocarbons (PAHs) are among the most frequently analyzed contaminants in diverse matrices. The reason for this is the partly very high carcinogenicity of the substance group and their broad distribution. PAHs are formed during partial combustion and pyrolysis processes. They enter the environment through fumes, bound to dust particles or mineral oil residues. PAHs contaminate foodstuff either from the environment (e.g. dust particles on broad-leaf vegetables) or by roasting and smoking processes.

European guidelines set limiting values which vary depending on the foodstuff. For the marker substance benzo[a]pyrene the value lies between 1–6 µg/kg while the total limit value for the PAH4 (this term summarizes benzo[a]pyrene, benzo[a]anthracene, chrysene and benzo[b]fluoranthene) is 1–35 µg/kg.

Due to the required sensitivity and complex matrices, PAH determination is a great challenge for analytics. As an analytical method either HPLC with fluorescence detection or GC-MS are used. Sample preparation is often extensive and time-consuming.

Recent requirements to laboratories include a reduction in sample handling time (Turn-Around-Time) and an increase in the sample throughput. Still, no losses within the analytical parameters, e.g. sensitivity, precision and accuracy are desired. Consequently, highly automated analytical procedures are necessary.

The LC-GC technology by Axcel Semrau® is established for the routine analytics of foodstuff and therefore provides an optimal basis for such an application system. Axcel Semrau® has developed a LC-LC-GC-MS/MS solution and a method for PAH analytics, the imPAHct method: **innovative multidimensional PAH clean-up technology**.

To measure PAHs, a simple yet fast extraction of the sample is performed. The extract is quickly manually treated and then transferred to the LC-GC system by Axcel Semrau®. Here, the extract is purified by a two-step LC separation (LC-LC Clean-Up) and subsequently analyzed via GC-MS. Edible oils may be measured straightforwardly.

Applying this intelligent matrix management, it is possible to determine the PAHs which are prescribed by EFSA fully automated. Depending on the detector, detection limits are reached which are even lower than the required European values for infant formula by a factor of 100.


The LC-GC solutions by Axcel Semrau® are pre-installed in application laboratories, tested and delivered to the customer ready-to-use. Thus, the quickest possible continuation of the routine measurement is assured.

### Example chromatograms


The LC-GC coupling provides two chromatograms at the same time:

- the signal of the UV detector from the HPLC purification
- the signal from the GC-MS

The illustrations on the following pages show examples of these chromatograms.



**Figure 1:** LC-LC Clean-Up: HPLC chromatogram of the first purification (upper chromatogram) and HPLC chromatogram of the two step HPLC purification (bottom chromatogram).



**Figure 2:** LC-LC-GC-MS/MS measurement of cocoa butter which was doted to 0.04 µg/kg with PAH.

## Specifications

- LOD for all EFSA PAHs between 0.01 – 0.02 µg/kg
- LOQ for all EFSA PAHs < 0.06 µg/kg
- sensitivity values apply if using a Bruker EVOQ GC-TQ
- repeatability < 10 % for entire method
- duration of analysis: 45 min

## Benefits of the imPAHct method

- LOD 1/100 of the limit value for infant formula (example substance benzo[a]pyrene)
- no consumables necessary due to LC-LC Clean-Up (e.g. SPE cartridges)
- longer maintenance intervals due to Clean-Up technology
- fast results
- high protection of GC due to
- high degree of automation
- excellent reproducibility
- applicable to cosmetics and all food-stuff
- investment safety
- expandable to other applications

## System components

The application system consists of the following components:

- CHRONECT Robotic consisting of:
  - CTC PAL3 with an 85 cm axis
  - Software CHRONOS
- Agilent 1260 Infinity II HPLC pump with UV detector and degasser
- Bruker EVOQ GC-TQ
- CHRONECT LC-GC interface for the LC-GC coupling
- data system with control and evaluation software
- accessories and consumables
- Factory Acceptance Test
- instruction and commissioning
- Site Acceptance Test
- training, support

**imPAHct and the CHRONECT Workstation PAH are developments by Axel Semrau®.**

Subject to technical changes

### **Axel Semrau GmbH & Co. KG**

Stefansbecke 42  
45549 Sprockhövel  
Tel.: + 49 2339 / 12090  
Fax: + 49 2339 / 6030

[www.axel-semrau.de](http://www.axel-semrau.de)  
[info@axel-semrau.de](mailto:info@axel-semrau.de)