

# Trap Selection, Fast GC, and Troubleshooting Strategy for Purge & Trap GC Volatiles

Updated: 18-Jun-2014



## Agenda

Overview

Basic P&T Operation

Adsorbents

Trap Selection

Fast GC for Volatiles

Troubleshooting Strategy

Summary / Resources

[sigma-aldrich.com/analytical](http://sigma-aldrich.com/analytical)

## Overview

Why this Topic?


- Purge & trap (P&T) in combination with gas chromatography (GC) is a widely used analytic technique for the analyses of **volatiles**, particularly in the environmental laboratory industry
- Many analysts performing this work may not fully understand, nor appreciate, the **underlying physical/chemical aspects**



# Overview

## Goals

- Gain a basic understanding of what adsorbents and traps do
- Learn how to **select the proper trap** for their applications
- Learn the Principles of **Fast GC** for volatiles
- Learn how to employ a **sound troubleshooting strategy** for use when things go wrong



## Basic P&T Operation

### Steps Prior to the Start of GC Analysis

- **Purge**: helium is bubbled through the sample, carrying volatile analytes to the trap
- **Dry-purge (optional)**: dry helium (bypassing the sample) is passed across the trap, removing any residual moisture from the system
- **Desorb pre-heat**: the trap is heated with no helium flow, causing analytes to be thermally desorbed from the adsorbent(s)



The purge parameters are typically method-specified.

## Basic P&T Operation

### Steps During GC Analysis

- **Desorb:** helium flow is added to the trap in the opposite direction as the purge mode, resulting in analytes being swept onto the head of the GC column
- **Bake:** the combination of helium flow, heat, and time are used to remove any non-desorbed compounds from the trap, venting them away from the GC, leaving the trap ready for the next sample



The desorb parameters are typically method-specified.

## Adsorbents

Commonly Used in Purge Traps

- **Carbon** adsorbents (Carbopack, Carbotrap, Carboxen, Carbosieve, activated charcoal) trap analytes based on analyte size/shape
- **Tenax** (a porous polymer) traps analytes predominantly based on analyte size/shape
- **Silica gel** traps analytes predominantly based on analyte polarity




Each adsorbent has specific advantages and drawbacks.

# Adsorbents

## Particles and Pores

- Pore composition determines **adsorption and desorption characteristics**
- Specialty carbons can be engineered with:
  - More or less of any pore type to serve a **specific purpose**
  - Through pores or closed pores, which influences **microporous strength**
  - Tapered pores (from macro- to meso- to micro-) which increases **thermodynamic efficiency**



Pore: any cavity present on a solid surface with a depth:width ratio of ~10:1.

Macropore diameter >500 Å; Mesopore diameter 20 – 500 Å; Micropore diameter <20 Å

# Adsorbents

## Carbon



### Advantages

- High desorption temperature (250 °C)
  - Rapid transfer to the GC, **less band broadening** and **improved chromatography**, especially the gases
  - **Better removal of large compounds** stuck in pores during the bake step
- Hydrophobic nature: **less water** trapped and pushed over to the GC
- Can be dry-purged: to **remove residual moisture** from the system
- Engineered pore composition: **application specific**

### Drawbacks

- The stronger carbon adsorbents (Carboxens and Carbosieves) have a great **affinity for methanol**
  - **Fills pores**, leaving none available for gases (blow off back of trap)
  - May **displace gases** (get adsorbed, get displaced, blow off back of trap)
- Adsorbed methanol will be desorbed during the desorption step, leaving the trap ready for the next run (however, the gases will be biased low for that run)
- Solution is to keep total amount of methanol in samples, standards and blanks <50 µL (<20 µL is desirable)



# Adsorbents

## Tenax



### Advantages

- High desorption temperature (250 °C)
- The ability to trap **analytes with a wide range** of sizes/shapes (as small as pentane; larger than trichlorobenzene)
- **Well-studied** (breakthrough volume data for numerous analytes)
- **Referenced** in many methods and numerous literature and **easily obtainable**

### Drawbacks

- **Not suitable for the gases** and other small molecules
- Oxidation of the polymer framework can produce **high background** levels as residual monomers left in the adsorbent are released
  - Detected as aromatic ring compounds
  - Voids can be created, leading to **channeling**

Tenax is a porous polymer based on 2,6-diphenylene-oxide.

# Adsorbents

## Silica Gel



### Advantages

- Great **affinity for polar analytes** (ketones) that are very water-soluble and may not trap well on other adsorbents
- **Referenced** in many methods and numerous literature
- **Easily obtainable**

### Drawbacks

- **Lower recommended desorption temperature** (180 °C)
- Its great **affinity for moisture**: more water is trapped and subsequently pushed over to the GC
- **Cannot be dry-purged** without also removing polar analytes

Unlike carbon and Tenax, silica gel traps analytes predominantly based on analyte polarity instead of analyte size/shape.

## Trap Selection

Adsorbent Bed Order (regardless of adsorbents)

- Sample comes into contact with the **weakest adsorbent first** and the **strongest adsorbent last**
  - Larger compounds trap on the first bed while the smaller compounds pass through the first bed(s) and are trapped on later bed(s)
  - If a **large compound** gets onto a **strong adsorbent**, it may become **irreversibly retained**
- Traps work similar to a sieve system



Traps can also be designed without a strong adsorbent bed (small compounds blow off the back of the trap during the purge mode).

# Trap Selection

Which Trap is 'Best'?

- **No absolute 'best' trap** for all situations
- It really depends on
  - The **method** being followed (the list of analytes)
  - Whether it is desirable to trap **small analytes** (the gases)
  - The **hardware** being used
    - Is moisture management available?
    - Is it being used or by-passed?
    - If being used, is it functioning properly?



# Trap Selection

## Handling Moisture

- **Water is a problem** for volatile analyses
  - Chromatographically **interferes with the gases**
  - **Damages quads** in GC-MS systems
- **Dry-Purge**
  - 1-2 minute will remove water to an acceptable level (without loss of adsorbed analytes)
  - 4-5 minute dry-purge will remove virtually all water (with some loss of adsorbed analytes)
- **Moisture Management**
  - Tubing section which grabs moisture as it passes
  - Must be maintained for proper operation



Note: Traps that contain a silica gel bed cannot be dry-purged.

Note: Do not use a dry-purge with a properly operating moisture management.

# Trap Selection

Most Methods (dry-purge; no moisture management)

## “K” Trap


- Supelco's most popular trap
- Able to adsorb/desorb small analytes (gases) to less volatile analytes (1,3,5-trichlorobenzene)

## “I” Trap

- Has a fourth bed in addition to the three beds in the “K” trap
- The fourth bed is Carboxen 1000
  - Weaker material than any bed in a “K” trap
  - For more efficient adsorption/desorption of larger analytes

### Adsorbent Key

Tenax TA	Green
Silica Gel 15	Purple
Activated Charcoal	Blue
3% SP-2100*/ Chromosorb W AW	Yellow
Carbopack B	Red
Carbopack C	Dark Green
Carbosieve S-III <sup>®</sup>	Cyan
Carboxen 1000	Gold
Carboxen 1001	Pink



# Trap Selection

Most Methods (moisture management; no dry-purge)

## “C” Trap or “#10” Trap


- Great traps as long as moisture is being managed by the concentrator
- Both traps have similarly arranged beds using a combination of different materials
  - Tenax polymer bed for analytes larger in size than pentane (BTEX and halocarbons)
  - Silica gel bed for polar analytes (ketones)
  - Carbon bed for small analytes (gases)
    - “C” trap has activated charcoal
    - “#10” trap has a stronger carbon

## “K” Trap or “I” Trap

- Will also work

### Adsorbent Key

Tenax TA	
Silica Gel 15	
Activated Charcoal	
3% SP-2100*/ Chromosorb W AW	
Carbopack B	
Carbopack C	
Carbosieve S-III <sup>®</sup>	
Carboxen 1000	
Carboxen 1001	



# Trap Selection

## Medium-Level Soil Methods

100 µL methanol extract added to 5 mL water

### “E” Trap or “D” Trap


- This amount of methanol is not a problem
- Drawbacks are lower (180 °C) desorption temperature (peaks not as sharp) and increased moisture transferred to the GC
- Try “D” trap if increased moisture is an issue (no silica gel bed)

### “K” Trap or “I” Trap

- Will also work

#### Adsorbent Key

Tenax TA	Light Green
Silica Gel 15	Purple
Activated Charcoal	Blue
3% SP-2100*/Chromosorb W AW	Yellow
Carbopack B	Red
Carbopack C	Dark Green
Carbosieve S-III <sup>®</sup>	Cyan
Carboxen 1000	Gold
Carboxen 1001	Pink





# Trap Selection

## BTEX-Only Methods

### “J” Trap


- Designed specifically for BTEX analysis
- Will not adsorb small compounds (methanol, gases, MTBE), which allows for a shorter analytical run

### “K” Trap or “I” Trap

- Will also work

#### Adsorbent Key

Tenax TA	
Silica Gel 15	
Activated Charcoal	
3% SP-2100*/ Chromosorb W AW	
Carbopack B	
Carbopack C	
Carbosieve S-III <sup>®</sup>	
Carboxen 1000	
Carboxen 1001	



# Trap Selection

BTEX+MTBE and/or GRO methods










## “M” Trap or “L” Trap


- Developed for underground storage tank (UST) and leaking underground fuel tank (LUFT) testing
- Methanol and the gases will not be trapped
  - Use “M” trap for BTEX, MTBE and analytes as small as n-pentane
  - Use “L” trap for BTEX and MTBE

## “K” Trap or “I” Trap

- Will also work

### Adsorbent Key

Tenax TA	
Silica Gel 15	
Activated Charcoal	
3% SP-2100*/Chromosorb W AW	
Carbopack B	
Carbopack C	
Carbosieve S-III™	
Carboxen 1000	
Carboxen 1001	



## Fast GC for Volatiles

### The Principles of Fast GC

- Decrease analysis time by using:
  1. Shorter column
  2. Quicker oven temperature ramp rate
  3. Higher carrier gas linear velocity

*But these changes also decrease resolution!*


- Offset the decrease in resolution by also using:
  4. Narrow I.D. column
  5. Hydrogen carrier gas
  6. Low film thickness



Fast GC is manipulating a number of column and instrument parameters to provide faster analysis times while maintaining resolution.

## Fast GC for Volatiles


Conventional GC (30 m x 0.25 mm I.D., 1.4  $\mu\text{m}$ ), 30:1 Split



17 minute run; poor responses for the gases.

## Fast GC for Volatiles


Fast GC (20 m x 0.18 mm I.D., 1.0  $\mu\text{m}$ ), 30:1 Split



9.5 minute run; still poor responses for the gases.

## Fast GC for Volatiles

Fast GC (20 m x 0.18 mm I.D., 1.0  $\mu\text{m}$ ), 100:1 Split



9 minute run; much improved responses and peak shapes for the gases.

## Fast GC for Volatiles

### Modified Column Parameters and Instrument Conditions

	Conventional GC	Fast GC (30:1 split)	Fast GC (100:1 split)
desorption process	40 mL/min at 210 °C for 2 min	40 mL/min at 210 °C for 2 min	<b>150 mL/min</b> at 210 °C for 2 min
column	SPB-624, 30 m x 0.25 mm I.D., 1.4 µm (24255)	SPB-624, <b>20 m x 0.18 mm I.D., 1.0 µm (28662-U)</b>	SPB-624, 20 m x 0.18 mm I.D., 1.0 µm (28662-U)
oven	40 °C (2 min), 7 °C/min to 135 °C, 30 °C/min to 230 °C (3 min)	<b>40 °C (1 min), 11 °C/min to 125 °C, 35 °C/min to 230 °C (2 min)</b>	40 °C (1 min), 11 °C/min to 125 °C, 35 °C/min to 230 °C (2 min)
carrier gas	helium, 1.1 mL/min	helium, <b>1.2 mL/min</b>	helium, <b>1.5 mL/min</b>
injection	30:1 split	30:1 split	<b>100:1 split</b>

Common mistake is not increasing the split when switching to a narrow I.D. column.  
Sensitivity maintained because very sharp/narrow peaks are produced (greater signal-to-noise).

## Fast GC for Volatiles

### Identical Column Parameters and Instrument Conditions

- sample/matrix: each analyte at 50 ppb in 5 mL water
- purge trap: VOCARB 3000 “K” (24940-U)
- purge: 40 mL/min at 25 °C for 11 min
- dry purge: 2 min
- desorption pre-heat: 205 °C
- bake.: 260 °C for 10 min
- transfer line/valve temp.: 110 °C
- inj. temp.: 150 °C
- detector: MS (interface at 200 °C), 35-400 m/z
- liner: 0.75 mm I.D., direct (SPME) type, straight design



## Fast GC for Volatiles

### Peak IDs

- |                                     |                                      |                                      |
|-------------------------------------|--------------------------------------|--------------------------------------|
| 1. Chloromethane                    | 14. 1,2-Dichloroethane-d4<br>(surr.) | 27. Tetrachloroethene                |
| 2. Vinyl Chloride                   | 15. Benzene                          | 28. Dibromochloromethane             |
| 3. Bromomethane                     | 16. 1,2-Dichloroethane               | 29. Chlorobenzene-d5 (I.S.)          |
| 4. Chloroethane                     | 17. Fluorobenzene (I.S.)             | 30. Chlorobenzene                    |
| 5. Trichlorofluoromethane           | 18. Trichloroethene                  | 31. Ethylbenzene                     |
| 6. 1,1-Dichloroethene               | 19. 1,2-Dichloropropane              | 32. Bromoform                        |
| 7. Methylene chloride               | 20. Bromodichloromethane             | 33. 4-Bromofluorobenzene<br>(surr.)  |
| 8. trans-1,2-Dichloroethene         | 21. 2-Chloroethyl vinyl ether        | 34. 1,1,2,2-Tetrachloroethane        |
| 9. 1,1-Dichloroethane               | 22. cis-1,3-Dichloropropene          | 35. 1,3-Dichlorobenzene              |
| 10. Chloroform                      | 23. Toluene-d8 (surr.)               | 36. 1,4-Dichlorobenzene-d4<br>(I.S.) |
| 11. Dibromofluoromethane<br>(surr.) | 24. Toluene                          | 37. 1,4-Dichlorobenzene              |
| 12. 1,1,1-Trichloroethane           | 25. trans-1,3-Dichloropropene        | 38. 1,2-Dichlorobenzene              |
| 13. Carbon tetrachloride            | 26. 1,1,2-Trichloroethane            |                                      |

## Fast GC for Volatiles

Taking Advantage of a <10 Minute GC Run

- **Ballistic Dry-Purge and Bake**

- EST: Encon Evolution
- Tekmar: Stratum PTC
- Tekmar: Velocity XPT



- **Two Concentrators per GC**

- EST: Centurion WS
- EST: PT<sup>2</sup>
- OI: PT Express



Example equipment that allows P&T to catch up with improvements in GC speed.

# Troubleshooting Strategy

## Best Practices

### Documentation


- Use maintenance log books
- May save weeks of troubleshooting

### Make one change at a time

- To uncover root cause
- Multiple changes may offset each other


### Keep a 'good' trap

- Remove and store the trap as a reference for when issues occur at a later time
- Replace the caps, place in original shipping container, label properly, and protect from vibration



# Troubleshooting Strategy

Overview / Flow



Supelco Technical Service strategy successfully used to help many customers over many years.

## Troubleshooting Strategy

### Potential Cause: Moisture Issue

- As samples are purged, moisture is carried into the concentrator
- If **active sites or cold spots** develop, **moisture will accumulate**, causing
  - Low/no response of the more water-soluble analytes
  - Decreased I.S. response throughout the tune window
- **Opening the system** allowed the moisture to escape
  - When the non-working trap was removed
  - Also why the non-working trap now seems to work



Perform preventative maintenance on moisture management  
or find/eliminate active site and/or cold spot.

## Troubleshooting Strategy

Potential Cause: Trap Issue

Try an **elevated bake-out**

- Make sure that a PTFE ferrule is NOT installed on the trap (use a VESPEL or graphite ferrule able to withstand 400 °C or higher)
- Bake trap with 40 mL/min helium at 375 °C for 4 hours
- Cool the trap down and discard the ferrule
- Replace with a new ferrule (whatever normally used)
- Evaluate performance (if it works, great; if not, it is a 'bad' trap)



Note: This procedure is only for use with a Supelco carbon trap (H, I, J, K, L, or M). These traps can easily handle this temperature. The reason that Supelco does not list this on the data sheet is that many customers use PTFE ferrules (maximum temperature is 300 °C) and we do not want to be responsible for melted PTFE stuck in instruments. Customers must verify that their instrument can handle this temperature.

# Troubleshooting Strategy

Potential Cause: Incorrect Parameter

- Autosampler, concentrator, GC inlet, column dimension, and detector combinations are **quite large**
- Impossible to write up what each individual parameter should be (depends on hardware being used)
- Recommend contacting a **Technical Service** (have **complete list of all operating parameters** handy)
  - Concentrator manufacturers (CDS, EST, O.I., Tekmar)
  - Trap manufacturers (CDS, EST, O.I., Supelco, Tekmar)
  - System manufacturers (Agilent, PerkinElmer, Shimadzu, Thermo, Varian)
  - Third-party (MDL [[www.ecsmdl.com](http://www.ecsmdl.com)])



## Troubleshooting Strategy

Potential Cause: System Issue

Make **50 and 200 ng/ $\mu$ L standards (with I.S./surr.) in methanol** and perform the following **five runs**

1. Inject 0.5  $\mu$ L (200 ng/ $\mu$ L) into the GC (if accessible)
2. Inject 2  $\mu$ L (50 ng/ $\mu$ L) into the top of the trap; step to desorb pre-heat
3. Inject 2  $\mu$ L (50 ng/ $\mu$ L) into the top of the trap; dry-purge 5 minutes then step to desorb pre-heat
4. Inject 2  $\mu$ L (50 ng/ $\mu$ L) into 5 mL water; set-up manually on concentrator (do not use autosampler)
5. Inject 2  $\mu$ L (50 ng/ $\mu$ L) into 5 mL water; set-up manually on the autosampler



A drastic loss of response indicates where the issue resides.  
Focus further troubleshooting in this area.



# Troubleshooting Strategy

## Troubleshooting Worksheet (T596001)

- Should work for any combination of:
  - Trap type
  - Sample matrix
  - Concentrator make/model
  - Moisture reduction strategy
  - GC conditions
- Fill out and use for phone, fax, or email correspondence with any of the Technical Services

**Purge & Trap Troubleshooting Worksheet**

Name: \_\_\_\_\_ Date: \_\_\_\_\_  
Company: \_\_\_\_\_ Tel: \_\_\_\_\_  
Address: \_\_\_\_\_ Fax: \_\_\_\_\_  
City, State, Zip: \_\_\_\_\_ E-mail: \_\_\_\_\_

-----Trap-----

Trap Type: \_\_\_\_\_ Catalog Number: \_\_\_\_\_ Lot: \_\_\_\_\_  
Was the Previous Trap OK: \_\_\_\_\_ Previous Trap Lot: \_\_\_\_\_  
Has This Application Ever Worked on This Instrument: \_\_\_\_\_

-----Sample-----

Method(s): \_\_\_\_\_ Sample Matrix: \_\_\_\_\_  
Purge Volume: \_\_\_\_\_ Temp: \_\_\_\_\_  
ICAL Range: \_\_\_\_\_  
Total  $\mu$ L Methanol (from STD, IS, SS) in Each Standard: \_\_\_\_\_  
Total  $\mu$ L Methanol (from IS, SS) in Each Sample: \_\_\_\_\_  
If Water Samples, Preserved ( Yes / No ) If Preserved, Using: \_\_\_\_\_

-----Concentrator-----

Autosampler Make / Model: \_\_\_\_\_ Concentrator Make / Model: \_\_\_\_\_  
If Water Samples, Purge Glassware ( Fritted / Fritless )  
Initial Trap Conditioning Temp: \_\_\_\_\_ Time: \_\_\_\_\_  
Purge Flow: \_\_\_\_\_ Temp: \_\_\_\_\_ Time: \_\_\_\_\_  
Dry-Purge Flow: \_\_\_\_\_ Temp: \_\_\_\_\_ Time: \_\_\_\_\_  
Desorb Pre-Heat Temp: \_\_\_\_\_  
Desorb Flow: \_\_\_\_\_ Temp: \_\_\_\_\_ Time: \_\_\_\_\_  
Bake Flow: \_\_\_\_\_ Temp: \_\_\_\_\_ Time: \_\_\_\_\_  
Transfer Line Temp: \_\_\_\_\_ Valve Temp: \_\_\_\_\_  
Moisture Reduction (i.e. MCM, MCS, etc.) Used ( Yes / No )  
Temp During Purge: \_\_\_\_\_ Desorb: \_\_\_\_\_ Bake: \_\_\_\_\_

-----GC-----

GC Make / Model: \_\_\_\_\_  
Injector Configuration ( Split / Direct ) Temp: \_\_\_\_\_ EPC Used ( Yes / No )  
Column Flow: \_\_\_\_\_ Split Flow: \_\_\_\_\_  
Liner Type: \_\_\_\_\_ Dimensions: \_\_\_\_\_  
Column Phase: \_\_\_\_\_ Dimensions: \_\_\_\_\_  
Oven Temperature Program: \_\_\_\_\_  
Detector Type / Make / Model: \_\_\_\_\_ Temp: \_\_\_\_\_  
If MS, Scan Range: \_\_\_\_\_ MS Pump Type ( Turbo / Diffusion )

-----Problem-----

How Do You Know the Trap is Defective: \_\_\_\_\_



## Summary / Resources

### Summary

- Select the proper trap for the application
- “K” trap will work for most applications
- Employ Fast GC to increase throughput
- Use a sound troubleshooting strategy

Web Site ([sigma-aldrich.com/environmental](http://sigma-aldrich.com/environmental))

- Sections for organics, metals, wet chemistry, sample collection, chromatograms, and newsletters
- From sample collection through sample preparation to analysis, we have you covered



Supelco site  
Bellefonte, Pennsylvania (USA)



Fluka site  
Buchs (Switzerland)



**Thank You!**

