

# Measuring Ethanol in Kombucha Tea by Headspace Solid Phase Microextraction (HS-SPME); a Walk Through Method Development and Validation

Katherine K. Stenerson, Olga I. Shimelis, Teresa Marsala  
MilliporeSigma, Bellefonte PA

## Introduction

Recently, AOAC adopted a method for first action status using headspace solid phase microextraction (HS-SPME) for analysis of ethanol in Kombucha tea.

Described here are some of the parameters that were optimized in development of an HS-SPME method to meet standard method performance requirements (SMPR) as designated by the working group within AOAC.

## Method

Per the SMPR, the method had to meet the following requirements:

- Analytical range of 0.1 – 2.0% alcohol by volume (ABV)
- LOQ  $\leq$  0.05% ABV
- Accuracy over analytical range: 97-102%
- Method repeatability:  $\leq$ 4% RSD
- Method reproducibility:  $\leq$ 6% RSD

In developing the final method (Tables 1 & 2), the following parameters were optimized:

- Fiber chemistry** – PDMS chosen over Carboxen/PDMS due to greater capacity - Figure 1.
- Equilibration/Extraction temperature & time** – optimized for best reproducibility, highest sensitivity and lowest method cycle time.
- Internal standard** – Ethanol-d6 chosen for best accuracy. Addition to samples was done with the diluent for best reproducibility. This I.S. then required use of GC-MS for analysis.
- Sample diluent** – salt addition and buffer increased response and reproducibility (Table 3). Buffer at pH 7 also normalized matrix between sample types.
- Calibration** - high sample dilution allowed for calibration from water (Figure 2).

### Table 1. Optimized HS-SPME Method

Sample/matrix: 400  $\mu$ L tea sample + 3.6 mL salt/buffer soln. (.05M Na<sub>2</sub>HPO<sub>4</sub> at pH 7 w/25% NaCl) containing I.S. in 10 mL headspace vial

Internal standard: ethanol-d6; at .08% ABV in .05M salt/buffer solution (chilled prior to use with samples)

SPME fiber: 100  $\mu$ m PDMS

Incubation: 7 min, 40 °C

Extraction: 2 min, headspace, 40 °C

Desorption: 3 min, 250 °C, split 10:1

Fiber post-bake: 5 min, 260 °C

### Table 2. GC-MS conditions

Column: Supelcowax-10, 30 m x 0.25 mm I.D., 0.50 $\mu$ m  
Oven: 40°C (5 min), 8°C/min to 70°C, 20°C/min to 250°C (5 min)  
MSD interface: 250°C  
Scan range: m/z 25-300  
Carrier: helium, 1 mL/min constant flow  
Injection: SPME fiber, split 10:1  
Liner: 0.75 mm I.D. SPME

Table 3. Effect of salt and buffer on response and reproducibility

dilution medium:	ethanol-d6 response (area counts)	
	salt/buffer	water
<b>Matrix</b>		
<b>water</b>	24958	19623
<b>raspberry lemon Kombucha</b>	24879	17289
<b>mango Kombucha</b>	26986	22006
<b>cherry Kombucha</b>	23684	17677
<b>cranberry Kombucha</b>	25785	14527
<b>average response</b>	<b>25258</b>	<b>18224</b>
<b>%RSD</b>	<b>5%</b>	<b>15%</b>

Figure 1. Comparison of ethanol uptake from water; Carboxen vs. 100  $\mu$ m PDMS fiber





Figure 2. Ethanol calibration; Kombucha vs. water: both diluted 10:1 into salt/buffer diluent.



- From 0.1 – 2% ABV, accuracy was 98-100%, and repeatability 1-3% RSD. (Table 4)
- LOD/LOQ was determined at a spike concentration of 0.1% ABV:
  - LOD = .01% ABV
  - LOQ = .03% ABV
  - Accuracy at the LOD and LOQ were 83% and 95%, respectively, for n=5 spiked replicates.
- Method accuracy was 98-100% using certified reference materials of ethanol in water and low alcohol beer. Repeatability was 2-5% RSD (Table 5).
- Method repeatability determined by triplicate analyses of kombucha samples done on two separate days was from 1-3% RSD. The ethanol levels in these samples ranged from 0.39 to 1.87% ABV (Table 6).

## Validation Results

Table 4. Method accuracy and repeatability in spiked Kombucha samples from 0.1 to 2.0% ABV

alcohol spiking level (%ABV)	amt. of alcohol measured in unspiked Kombucha (%ABV)	avg. %ABV measured	avg. % ABV less unspiked	repeatability %RSD (n=5)	% accuracy
<b>0.10</b>	.011	0.11	.098	3*	98
<b>0.50</b>	.011	0.51	0.50	2	100
<b>1.00</b>	.011	1.00	0.99	1	99
<b>2.00</b>	.026	1.99	1.96	1	98

\*n=7 replicates (used for LOD, LOQ calcs.)

Table 5. Method accuracy and repeatability using certified reference materials analyzed on multiple days using different instruments and analysts

sample	certified conc. (%ABV <sup>a</sup> )	# of replicates	# instruments/analysts	range <sup>b</sup> (%ABV)	avg. amt. measured (%ABV)	avg. % accuracy	reproducibility %RSD <sup>c</sup>
<b>low alcohol beer</b>	0.51	10	2/2	0.48-0.51	0.50	98	2
<b>alcohol in water, 80 mg/dL</b>	0.10	4	2/1	.095-0.11	.099	98	5
<b>alcohol in water, 200 mg/dL</b>	0.25	3	2/2	0.24-0.26	0.25	99	3
<b>alcohol in water, 400 mg/dL</b>	0.51	7	2/2	0.49-0.52	0.51	100	3

Table 6. Method repeatability from different Kombucha tea varieties analyzed in 2 sets of 3 replicates over 2 days.

Sample #	Producer #	Flavor #	"Best by" date	Testing dates	Avg. %ABV	% RSD
1	1	1	10/20/2017	11/29, 12/1/2017	1.87	3
2	1	2	12/25/2017	11/29, 12/1/2017	1.36	3
3	1	3	12/1/2017	11/29, 12/1/2017	0.81	1
4	2	1	1/20/2018	11/29, 12/1/2017	0.66	2
5	3	1	2/28/2018	11/29, 12/1/2017	1.66	1
6	3	2	6/12/2018	11/29, 12/1/2017	1.54	3
7	4	1	12/19/2017	11/30, 12/2/2017	0.39	3
8	5	1	2/6/2018	11/30, 12/2/2017	1.06	2
9*	6	1	NA	11/30, 12/2/2017	0.49	2

\*homemade Kombucha provided by a colleague

## Summary

In relation to the SMPR:

- Analytical range of 0.1 – 2.0% ABV:** demonstrated with  $\leq$ 3% RSD in Kombucha matrix spiked over designated range.
- LOQ  $\leq$  0.05% ABV:** demonstrated with LOQ of .03% ABV (verified with spikes at this level)
- Accuracy over analytical range 97-102%:** demonstrated with range of 98-100% from 0.1 – 2.0 %ABV
- Method repeatability of  $\leq$ 4% RSD:** demonstrated with repeatability of  $\leq$ 3% for multiple replicates analyzed over 2 days.
- Method reproducibility of  $\leq$ 6% RSD:** demonstrated using certified reference materials with % RSDs of  $\leq$ 5% for multiple analyses done on 2 instruments by 2 different analysts.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada

SigmaAldrich.com