

APPLICATIONS

Optimized Analysis of 1,4-Dioxane in Tap Water by GC/MS Using a Zebron™ ZB-624 GC Column

Matthew Trass, Tim Anderson, and Kristen Parnell
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA



Timothy Anderson
GC Product Manager
Tim was raised in Texas where it was completely too hot, then moved to Pennsylvania and Ohio where it was entirely too cold. He finally settled on California where the weather is just right.



Growing concern over the negative health effects associated with dioxane have led to increased global regulation and testing of this compound. This work presents an optimized method for testing of 1,4-dioxane in drinking water by EPA Method 522 using a Zebron ZB-624 GC column.

Introduction

EPA Method 522 is used to determine 1,4-dioxane (dioxane) levels in drinking water. Dioxane is a highly water-soluble solvent which has a multitude of industrial uses. Therefore, there is a high likelihood that dioxane can migrate into groundwater systems. In fact, dioxane has been detected in groundwater supplies in numerous locations throughout the United States. Regulatory agencies are especially concerned because dioxane does not readily bind to soils and is resistant to naturally occurring biodegradation processes.

Because of its serious environmental implications, a rapid and robust analytical testing method for dioxane is necessary. EPA Method 522 requires a 624 type column for analysis by GC/MS following SPE extraction of either a 100 mL or 500 mL water sample. Presented here is an optimized EPA 522 GC/MS protocol on a Zebron ZB-624 GC column. The method's effectiveness is demonstrated on a 100 mL tap water sample.

Materials and Methods

Tap water samples were collected and prepared using the EPA 522 extraction procedure as shown in **Table 1**, and demonstration of method efficiency was performed using a secondary option for preparation as shown in **Table 2**. Following sample preparation, samples were analyzed by GC/MS as shown in **Table 3**.

Table 1.
Sample preparation protocol following EPA 522 extraction (Option 1).

Option 1	
Step	Protocol
1	Open the tap and allow the system to flush until the water temperature has stabilized. Collect samples from the flowing system into a vessel containing sodium sulfite (chlorine reducer (50 mg/L)).
2	Add sodium bisulfate (1 g/L) and verify that the pH is < 4.
3	Fortify samples with 1,4-dioxane (analyte) and 1,4-dioxane-d8 (surrogate standard) at the appropriate concentration levels. (Samples are now ready for solid phase extraction (SPE)).

Table 2.
Sample preparation protocol using SPE extraction and a vacuum manifold.

Option 2	
Step	Protocol
1	Condition Supelco® Supelclean™ ENVI-Carb™ Plus cartridges (54812-U) with 1 mL of dichloromethane (DCM), 2 mL of methanol, 2 mL of methanol and 3 mL of water.
2	Using a reservoir, load 100 mL tap water sample and elute at approximately 10 mL/min.
3	Dry cartridge using high vacuum for 10 minutes.
4	Fill each cartridge with DCM and soak for 1 minute.
5	Continue to add DCM and elute dropwise until 1.5 mL of eluent has been collected.
6	Bring sample volume to 2 mL with DCM and add THF-d8 internal standard at the appropriate concentration.
7	Remove residual water by passing samples through Strata® Silica anhydrous sodium sulfate cartridges (Part No.: 8B-S124-JCH). (Samples are now ready for GC/MS analysis.)

Table 3.
GC/MS conditions for analysis of tap water samples.

Column:	Zebron ZB-624		
Dimensions:	30 meter x 0.25 mm x 1.4 µm		
Part No.:	7HG-G005-27		
Guard Column:	5 m Z-Guard™ (7AG-G000-00-GZ0)		
Injection:	Splitless @ 200 °C (hold 0.5 min), 2.0 µL		
Recommended Liner:	Zebron PLUS Single Taper Z-Liner™, 4 mm ID		
Liner Part No.:	AG2-0A13-05 (for Agilent® systems)		
Carrier Gas:	Helium @ 1.5 mL/min (constant flow)		
Oven Program:	30 °C for 1 min to 90 °C @ 20 °C/min to 250 @ 50 °C/min		
Detector:	MSD @ 250 °C		
Sample:	Peak	Analyte	SIM Ions Window (min)
	1	THF-d8	46*, 78, 80 4.00 – 4.75
	2	1,4-dioxane-d8	96*, 62, 64 4.75 – 7.8
	3	1,4-dioxane	88*, 58 4.75 – 7.8




Results and Discussion

As shown in **Figures 1** and **2**, the Zebtron™ ZB-624 GC column provided excellent peak shape and resolution for the EPA Method 522 analytes. The GC/MS method was modified from the original EPA 522 GC parameters to reduce analysis time. Specifically, the carrier gas flow rate was increased to 1.5 mL/min and the oven ramps were increased from 7 °C/min to 20 °C/min and 20 °C/min to 50 °C/min. These method modifications result in a faster analysis time of less than 8 minutes without sacrificing method performance. The neat seven-point calibration curve generated ranged from 5 to 2000 ng/mL and showed

good linearity for 1,4-dioxane ($R^2 = 0.9994$) and 1,4-dioxane-d8 ($R^2 = 0.9995$).


To demonstrate method efficiency, a tap water sample fortified with analyte and surrogate standards at low (3 µg/L) and high (25 µg/L) levels was processed. In this instance, extraction protocol Option 2 (shown in **Table 2**), which requires a 100 mL sample extracted with a 400 mg activated carbon SPE cartridge, was chosen. The sample extracts were quantified by comparison to a standard calibration curve (**Figures 3** and **4**). The results, listed in **Table 4**, show high recovery values that fall within the 70 – 130 % range specified by the EPA method.

Figure 1. Extracted ion chromatogram for 1000 ng/mL calibration curve, THF-d8 (internal standard) at 500 ng/mL.




App ID 24123

Figure 2. Extracted ion chromatogram for tap water extract, 25 ng/mL.



App ID 24125

Figure 3.
1,4-dioxane (analyte) calibration curve, 5 to 2000 ng/mL.



Conclusion

The Zebron™ ZB-624 GC column provided excellent peak shape and resolution for the EPA Method 522 GCMS analysis. The GC/MS method was modified from the original EPA 522 GC parameters, resulting in a faster analysis time and excellent linearity.

The extraction data collected shows that the modified GC method is effective for the analysis of real-world samples and will benefit a laboratory wishing to improve sample throughput without sacrificing method performance.

Figure 4.
1,4-dioxane-d8 (surrogate standard) calibration curve, 5 to 2000 ng/mL.




Table 4.
100 mL tap water LFSM mean recovery values.

Compound	µg/L	Mean Recovery (%) (n=3)	RSD (%) (n=3)
1,4-dioxane	3	96	4.2
1,4-dioxane-d8	3	109	5.9
1,4-dioxane	25	103	7.8
1,4-dioxane-d8	25	105	8.8



APPLICATIONS


Ordering Information

Zebron™ ZB-624 GC Columns

Length (m)	ID (mm)	df (µm)	Temp. Limits (°C)	Part No.
20	0.18	1.00	-20 to 260	7FD-G005-22
30	0.25	1.40	-20 to 260	7HG-G005-27
30	0.32	1.80	-20 to 260	7HM-G005-31
30	0.53	3.00	-20 to 260	7HK-G005-36
60	0.25	1.40	-20 to 260	7KG-G005-27
60	0.32	1.80	-20 to 260	7KM-G005-31
60	0.53	3.00	-20 to 260	7KK-G005-36
75	0.53	3.00	-20 to 260	7LK-G005-36
105	0.53	3.00	-20 to 260	7NK-G005-36

Note: If you need a 5 in. cage, simply add a (-B) after the part number, e.g., 7HG-G005-27-B. Some exceptions may apply. Agilent 6850 and some SRI and process GC systems use only 5 in. cages.

Zebron PLUS GC Inlet Liners

Description	Dimensions ID x L (mm)	Part No.	Unit.
For Agilent® or Thermo Scientific® GC Systems			
Single Taper Z-Liner™	 4 x 78.5	AG2-0A13-05	5/pk
		AG2-0A13-25	25/pk

Find more Zebron GC columns and
accessories at
www.phenomenex.com/GC



If Phenomenex products in this technical note do not provide at least an equivalent separation as compared to other products of the same phase and dimensions, return the product with comparative data within 45 days for a FULL REFUND.

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Zebron, Guardian, Z-Liner, and Z-Guard are trademarks and Strata is a registered trademark of Phenomenex. Supelclean and ENVI-Carb are trademarks and Supelco is a registered trademark of Sigma-Aldrich Co., LLC. Agilent is a registered trademark of Agilent Technologies, Inc. Thermo Scientific is a registered trademark of Thermo Fisher Scientific, Inc.

Disclaimer

Phenomenex is in no way affiliated with Sigma-Aldrich Co., LLC., Agilent Technologies, Inc., or Thermo Fisher Scientific, Inc.

© 2017 Phenomenex, Inc. All rights reserved.

Australia
t: +61 (0)2-9428-6444
f: +61 (0)2-9428-6445
auserinfo@phenomenex.com

Austria
t: +43 (0)1-319-1301
f: +43 (0)1-319-1300
anfrage@phenomenex.com

Belgium
t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
f: +31 (0)30-2383749
beinfo@phenomenex.com

Canada
t: +1 (800) 543-3681
f: +1 (310) 328-7768
info@phenomenex.com

China
t: +86 400-606-8099
f: +86 (0)22 2532-1033
phen@agela.com

Denmark
t: +45 4824 8048
f: +45 4810 6265
nordicinfo@phenomenex.com

Finland
t: +358 (0)9 4789 0063
f: +45 4810 6265
nordicinfo@phenomenex.com

France
t: +33 (0)1 30 09 21 10
f: +33 (0)1 30 09 21 11
franceinfo@phenomenex.com

Germany
t: +49 (0)6021-58830-0
f: +49 (0)6021-58830-11
anfrage@phenomenex.com

India
t: +91 (0)40-3012 2400
f: +91 (0)40-3012 2411
indiainfo@phenomenex.com

Ireland
t: +353 (0)1 247 5405
f: +44 1625-501796
eirinfo@phenomenex.com

Italy
t: +39 051 6327511
f: +39 051 6327555
italiainfo@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

Luxembourg
t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

Mexico
t: 01-800-844-5226
f: 001-310-328-7768
tecnicomx@phenomenex.com

The Netherlands
t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

New Zealand
t: +64 (0)9-4780951
f: +64 (0)9-4780952
nzinfo@phenomenex.com

Norway
t: +47 810 02 005
f: +45 4810 6265
nordicinfo@phenomenex.com

Puerto Rico
t: +1 (800) 541-HPLC
f: +1 (310) 328-7768
info@phenomenex.com

Spain
t: +34 91-413-8613
f: +34 91-413-2290
espinfo@phenomenex.com

Sweden
t: +46 (0)8 611 6950
f: +45 4810 6265
nordicinfo@phenomenex.com

United Kingdom
t: +44 (0)1625-501367
f: +44 (0)1625-501796
ukinfo@phenomenex.com

USA
t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com

**All other countries
Corporate Office USA** 
t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com