

## USEPA Procedures for Wastewater Analyses by Packed Column GC and HPLC

*This bulletin outlines the United States Environmental Protection Agency's procedures for analyzing 113 organic priority pollutants in wastewater. The chromatographic column or column packing listed by the EPA, and its Supelco equivalent, are described for each gas chromatographic-mass spectrometric, gas chromatographic, or high performance liquid chromatographic method. Confirmational columns and traps for concentrating analytes from samples are described for methods in which they are specified. The analyses are illustrated with chromatograms.*

### Key Words

- wastewater ● priority pollutants ● volatile pollutants
- acidic pollutants ● base/neutral pollutants ● pesticides
- PCBs

The United States Environmental Protection Agency (USEPA) has amended its priority pollutant protocol (1) by revising the gas chromatographic-mass spectrometric methods (GC-MS) and developing alternative gas or high performance liquid chromatographic methods (GC or HPLC) for organic water pollutant analyses. These alternative methods of analysis enable each regulated industry to choose the approach that suits its individual needs or capabilities. The EPA considers these methods, described in the *Federal Register* (2), to be the only validated analytical approach to monitoring priority pollutants in wastewater. Industries are required to use the methods both in evaluating their wastewater and in applying for National Pollution Discharge Elimination System (NPDES) permits (3). Deviations from these methods require EPA approval in advance. General information about the procedures follows, preceding illustrations and brief outlines of the analyses.

### GC-MS Analyses for Organic Pollutants

For the revised GC-MS analyses, the organic pollutants are classified as volatiles or as nonvolatile acids (phenols), base-neutrals, and pesticides and PCBs, as they were in the 1977 protocol (1). The chromatographic columns described in the original protocol, however, have been replaced with columns developed through advanced column technology. The newer packings are more thoroughly deactivated than the packings they replace, and use of the columns listed in the revised methods improves the detection of priority pollutants. The GC-MS analyses for the organic compounds are incorporated into two methods, Method 624 for volatiles and Method 625 for nonvolatiles.

### GC or LC Organic Pollutant Analyses

The GC or HPLC methods (Methods 601-613) simplify chromatographic analyses by separating 113 organic pollutants into 13 classes and describing an analysis for each class. Each GC or HPLC

analysis specifies the use of a chromatographic column that will resolve all of the pollutants in one class and a detection system that is specific for their identification (2). Certain methods also suggest using a confirmational column to substantiate the identity of the pollutants.

Supelco GC packings were used by the EPA to evaluate these methods for priority pollutant analyses. Supelco packings are listed by the EPA in many of the methods, and for most of the other methods Supelco packings are generically equivalent alternatives to the listed packings. Supelco packings and columns were used to obtain all of the chromatograms in this bulletin.

### Concentration of Organic Pollutants

Organic pollutants normally are present at ppb or lower levels in wastewater samples, and they must be concentrated prior to chromatographic analysis. The volatile organics are purged from the water sample and are trapped on solid adsorbents (GC-MS Method 624 or GC Method 601, 602, or 603). The recommended purge-trap method has been described by Bellar and Lichtenberg (4). After the organics are adsorbed, they are rapidly heat desorbed from the trap onto the chromatographic column. The EPA recommends a different adsorbent trap for each of the four methods of analysis for volatile pollutants. Each trap provides the best trapping efficiency for the volatile compounds listed in the corresponding method.

The nonvolatile organics are concentrated by solvent extraction of the water sample followed by concentration of the extract to 1 mL in a Kuderna-Danish concentrator. Subsequently, the extract is dried on a sodium sulfate column. When water samples are extremely dirty, specific cleanup procedures are used with each extraction method. Details of the extraction techniques and cleanup procedures are given in the *Federal Register* (2).


### Outlines of EPA Methods

Outlines of the GC-MS methods for organic pollutant analyses are presented first, followed by outlines of the GC and HPLC methods. Each method outline lists the class of organic compounds, the concentration technique, the column, the confirmational column (if any), the detector, and any qualifications or potential problems that have been indicated by the EPA. For methods in which the listed packing is not a Supelco packing, the Supelco equivalent is given in parentheses. We have used two-meter glass columns for those methods that call for a nominal six-foot glass column. Two-meter glass columns give more consistent results, through tighter control of the column length and internal diameter, and thus of the volume of chromatographic packing.

These methods are discussed in detail in the *Federal Register* (2). The chromatographic packings and columns used in these analyses are listed following the method outlines.

**GC-MS Method 624 – Volatile Pollutants**  
**Figure A. Volatiles by EPA Method 624**

Column: **60/80 Carboxpack™ B/1% SP™-1000, 8' x 1/8" OD SS**  
 Cat. No.: **11815 (15g packing)**  
 Oven: **45°C (3 min) to 220°C at 8°C/min, hold 15 min**  
 Carrier: **helium, 40mL/min**  
 Det.: **FID, 250°C**  
 Inj.: **1µL dodecane containing 200-500ng/µL each analyte, 200°C**



**Sample Concentration**

Purge and trap, using a 30.5cm x 0.105" ID x 0.125" OD SS trap containing 15cm Tenax®, 8cm Silica Gel 15, and 1cm 3% SP-2100.

**GC Column**

6' x 1/8" OD SS containing 60/80 Carboxpack B/1% SP-1000. Figure A shows the analysis on an 8' column also used in Method 601. A 2" x 1/8" OD SS precolumn containing 3% SP-1000 on 60/80 Chromosorb® W AW is recommended.

**Confirmational Column**

8' x 1/8" OD SS containing 60/80 Carboxpack C/0.2% CARBOWAX® 1500. A 2" x 1/8" OD SS precolumn containing 3% CARBOWAX 1500 on 60/80 Chromosorb W AW is recommended.

**Detector**

Mass spectrometer. For Figure A an FID was used to simulate the mass spectrometer analysis.


**Qualifications**

Acrolein and acrylonitrile are volatile compounds, but purging efficiencies for them are low and erratic. For these compounds, the EPA recommends direct aqueous injection or the modified purge-trap procedure described in EPA Method 603.

**GC-MS Method 625 – Nonvolatile Pollutants**  
**I. Base Neutrals**

**Figure B. Base-Neutrals by EPA Method 625**

Column: **3% SP-2250 on 100/120 SUPELCOPORT™, 2m x 2mm ID glass**  
 Cat. No.: **11756 (20g packing)**  
 Oven: **50°C (4 min) to 270°C at 8°C/min**  
 Carrier: **nitrogen, 30mL/min**  
 Det.: **FID**  
 Inj.: **1µL methylene chloride containing 0.125mg/mL or 0.063mg/mL (\*) each analyte**



**Sample Concentration**

Solvent extraction with methylene chloride followed by concentration to 1mL.

**GC Column**

2m x 2mm ID glass containing 3% SP-2250 on 100/120 SUPELCOPORT.

**Confirmational Column**

None listed.

**Detector**


Mass spectrometer. For Figure B an FID was used to simulate the mass spectrometer analysis.

**Qualifications**

Analysts must be able to chromatograph 100ng of benzidine to demonstrate the inertness of the GC-MS system before they use this procedure on wastewater samples. The EPA has set the detection limit for most of the base-neutrals at 20ng/µL. For 6- and 7-ring polycyclic aromatics, however, the limit is 50ng/µL.

## II. Acidics (Phenols)

**Figure C. Phenols by EPA Method 625**



### Sample Concentration

Solvent extraction with methylene chloride followed by concentration to 1mL.

### GC Column

2m x 2mm ID glass containing 1% SP-1240-DA\* on 100/120 SUPELCOPORT.

### Confirmational Column

None listed.

### Detector

Mass spectrometer. For Figure C an FID was used to simulate the mass spectrometer analysis.

### Qualifications

Analysts must chromatograph 50ng of pentachlorophenol to evaluate the inertness of their GC-MS system before they use this method to analyze wastewater samples for phenols.


\*DA – Deactivated for acidic compounds.

**In Method 625, capillary columns may be used as an alternative to packed columns if they meet the QC criteria listed in Sections 8.2, 12, and 13.1 (4). Two available protocols for capillary methodology are 2 & 3 below:**

1. *Method of Organic Chemical Analysis of Municipal and Industrial Wastewater* US EPA EMSL, Cincinnati, OH 45268, EPA-600/4-83-057, July 1982.
2. *Testing Methods for Evaluating Solid Waste* US EPA Office of Solid Waste and Emergency Response, Washington, D.C., 20460, SW-846, July, 1982, 2nd Edition.
3. A.D. Sauter, et al., *J. High Resolution Chromatography and Chromatography Communications* 4: 366 (1981).

## III. Pesticides and PCBs

**Figure D. Pesticides by EPA Method 625**



### Sample Concentration

Solvent extraction with methylene chloride followed by concentration to 1mL. The specific procedure is described in EPA Method 608.

### GC Column

2m x 2mm ID glass containing 3% SP-2250 on 100/120 SUPELCOPORT.

### Confirmational Column

None listed.

### Detector

Mass spectrometer (see **Qualifications**). For Figure D an FID was used to simulate the mass spectrometer analysis.

### Qualifications

The sample is analyzed on the GC column with an electron capture detector (ECD) as described in GC Method 608. A GC-MS analysis is performed only if the GC/ECD analysis indicates the presence of pesticides or PCBs.

Pesticides with multiple peaks (PCBs, toxaphene, and chlordane) can be separated from the other pesticides by sample cleanup with Florisil® column chromatography (5).

Analysts must be able to chromatograph 100ng of aldrin before they use the GC-MS procedure on wastewater samples.

### Quantification of Priority Pollutants by GC-MS

Internal standards may be used in any of the GC-MS analyses, if they are adequately resolved from the priority pollutants in the sample and from matrix interferences. Two groups of initial standards recommended by the EPA are deuterated and fluorinated compounds of various aromatics, phenols, amines, and nitroaromatics. A mass spectrum of either 4-bromofluorobiphenyl or decafluorotriphenylphosphine must be checked daily to meet performance criteria set forth in Method 625.

If external standards are used to calibrate, individual standards of each compound must be analyzed. A three-point calibration must be performed, with one point of the calibration curve near the limit of detection. The calibration must be verified daily with at least one of the standard solutions used in the initial calibration. If significant variation ( $\pm 10\%$ ) occurs, a new calibration should be conducted.

## GC Method 601 – Purgeable Halocarbons

### Sample Concentration

Purge and trap, using a 30.5cm x 0.105" ID x 0.125" OD SS trap containing 7.7cm activated charcoal, 7.7cm Silica Gel 15, 7.7cm Tenax, and 1cm 3% SP-2100.

### GC Column

8' x 1/8" OD SS containing 60/80 Carboxpack B/1% SP-1000.

### Confirmational Column

8' x 1/8" OD SS containing n-octane on 80/100 Porasil® C.

### Detector

Hall® electroconductivity.

### Qualifications

Electron capture detectors are not recommended for this analysis because they have limited sensitivity for mono- and di-substituted halocarbons.

## GC Method 602 – Purgeable Aromatics

### Sample Concentration

Purge and trap, using a 30.5cm x 0.105" ID x 0.125" OD SS trap containing 1cm 3% SP-2100 and 23cm Tenax.

### GC Column

6' x 1/8" OD SS containing 5% SP-1200/1.75% Bentone® 34 on 100/120 SUPELCOPORT.

### Confirmational Column

8' x 1/8" OD SS containing 5% 1,2,3-tris(2-cyano-ethoxy) propane on 60/80 Chromosorb W AW.

### Detector

Photoionization (10.2 ev lamp). For Figure F an FID was used to simulate the photoionization detector.

### Qualifications

FID detection can be used if interfering organics can be resolved from the purgeable aromatics.

## GC Method 603 – Acrolein and Acrylonitrile

### Sample Concentration

Purge and trap, using a 30.5cm x 0.105" ID x 0.125" OD SS trap containing 1cm 3% SP-2100 and 23cm Tenax.

### GC Column

6' x 1/8" OD SS containing Durapak®/Carbowax 400 on 100/120 Porasil C.\*

### Confirmational Column

6' x 1/8" OD SS containing 80/100 Chromosorb 101.

### Detector


FID.

\*NOTE: 100/120 Porasil C is no longer available from the manufacturer.

## Figure E. Purgeable Halocarbons by EPA Method 601


Column: **60/80 Carboxpack B/1% SP-1000, 8' x 1/8" OD SS**  
Cat. No.: **11815** (15g packing)  
Oven: 45°C (3 min) to 220°C at 8°C/min, hold 15 min  
Carrier: nitrogen, 40mL/min  
Det.: Hall, 85°C  
Inj.: 1µL dodecane containing 200-500ng/µL each analyte, 200°C

- |   |  |
|---|--|
| 1. Chloromethane                          | 15. Bromodichloromethane   |
| 2. Bromomethane                           | 16. 1,2-Dichloropropane & trans-1,3-Dichloropropane                        |
| 3. Vinyl chloride & Dichlorofluoromethane | 17. Trichloroethylene  |
| 4. Chloroethane                           | 18. cis-1,3-Dichloropropene, 1,1,2-Trichloroethane, & Chlorodibromomethane |
| 5. Methylene chloride                     | 19. 1-Chloro-2-bromopropane (IS)   |
| 6. Trichlorofluoromethane                 | 20. 2-Chloroethyl vinyl ether  |
| 7. 1,1-Dichloroethylene                   | 21. Bromoform  |
| 8. Bromochloromethane (IS)                | 22. 1,1,2,2-Tetrachloroethylene & 1,1,2,2-Tetrachloroethane                |
| 9. 1,1-Dichloroethane                     | 23. 1,4-Dichlorobutane (IS)  |
| 10. trans-1,2-Dichloroethylene            | 24. Chlorobenzene  |
| 11. Chloroform                            |  |
| 12. 1,2-Dichloroethane                    |  |
| 13. 1,1,1-Trichloroethane                 |  |
| 14. Carbon tetrachloride                  |  |



## Figure F. Purgeable Aromatics by EPA Method 602


Column: **5% SP-1200/1.75% Bentone 34 on 100/120 SUPELCOPORT, 6' x 1/8" SS**  
Cat. No.: **12134** (20g packing)  
Oven: 50°C (2 min) to 90°C at 6°C/min, hold  
Carrier: helium, 36mL/min  
Det.: FID  
Inj.: 1µL methanol containing 0.2mg/mL each analyte



## Figure G. Acrolein and Acrylonitrile by EPA Method 603

Column: **80/100 Chromosorb 101, 6' x 1/8" OD SS**  
Cat. No.: **20214** (50g packing)  
Oven: 100°C (5 min) to 140°C at 10°C/min, hold  
Carrier: helium, 45mL/min  
Det.: FID  
Inj.: 2µL water


Acrolein, 3mg/mL  
Acrylonitrile, 1mg/mL



## GC Method 604 – Phenols

### Figure H. Phenols by EPA Method 604

Column: 1% SP-1240-DA\* on 100/120 SUPELCOPORT,  
2m x 2mm ID glass  
Cat. No.: 11832 (20g packing)  
Oven: 70°C (2 min) to 200°C at 8°C/min, hold  
Carrier: helium, 30mL/min  
Det.: FID, 250°C  
Inj.: 0.5µL standard mix (Cat. No. 48859)



#### Sample Concentration

Solvent extraction with methylene chloride followed by concentration to 1mL.

#### GC Column (Free Phenols)

2m x 2mm ID glass containing 1% SP-1240-DA\* on 100/120 SUPELCOPORT. The analysis is similar to the GC-MS analysis in Figure C.

#### Confirmational Column

None listed.

#### Detector

FID.

#### Qualifications

When sample cleanup is required, the extracted phenols are derivatized with pentafluorobenzylbromide (PFB) and the extract is cleaned by silica gel column chromatography. PFB phenols are separated on a column different from that used to separate free phenols, and are detected by electron capture.

#### GC Column (Derivatized Phenols)

1.8m x 2mm ID glass containing 5% OV®-17 on 80/100 Chromosorb W AW-DMCS. (Supelco Equivalent: 5% SP-2250 on 80/100 SUPELCOPORT.)

#### Detector


Electron capture detector.

\*DA – Deactivated for acidic compounds.

## HPLC Method 605 – Benzidine and 3,3'-Dichlorobenzidine

### Figure I. Benzidines by EPA Method 605

Column: SUPELICOSIL™ LC-1, 5cm x 4.6mm, 5µm particles  
Cat. No.: 58237  
Mobile Phase: 0.1M acetate buffer:acetonitrile, 60:40  
Flow Rate: 0.8mL/min  
Det.: UV, 254nm  
Inj.: 10µL mobile phase containing 0.5µg each analyte



#### Sample Concentration

Solvent extraction with chloroform followed by concentration to 1mL.

#### HPLC Column

25cm x 4.6mm ID SS containing LiChrosorb® RP-2. (Supelco Column: 25cm x 4.6mm ID SS containing SUPELICOSIL LC-1).

#### Confirmational Column

None listed.

#### Detector

Electrochemical (0.8 volts). For Figure I, UV detection at 254nm was used (see **Qualifications**).

#### Qualifications

A UV detector may be used if there are no matrix interferences.

The EPA procedure suggests using a 25cm x 4.6mm ID column, but a more than adequate separation is readily achieved with a 5cm x 4.6mm ID or a 15cm x 4.6mm ID SUPELICOSIL LC-1 column.

## GC Method 606 – Phthalate Esters

### Sample Concentration

Solvent extraction with methylene chloride followed by concentration to 1 mL.

### GC Column

2m x 4mm ID glass containing 1.5% SP-2250/1.95% SP-2401 on 100/120 SUPELCOPORT. Two temperatures, 180°C and 220°C, are used.

### Confirmational Column

2m x 4mm ID glass containing 3% OV-1 on 100/120 SUPELCOPORT. (Supelco Equivalent: 3% SP-2100 on 100/120 SUPELCOPORT.)


### Detector

Electron capture detector.

### Qualifications

If the extract must be cleaned, the EPA recommends Florisil or silica gel column chromatography.

Figure J. Phthalates by EPA Method 606



## GC Method 607 – Nitrosamines

### Sample Concentration

Solvent extraction with methylene chloride followed by concentration to 1 mL.

### GC Column

2m x 4mm ID glass containing 10% CARBOWAX 20M/2% KOH on 80/100 Chromosorb W AW. Two temperatures, 110°C and 220°C, are used.

### Confirmational Column

2m x 4mm ID glass containing 10% SP-2250 on 100/120 SUPELCOPORT.


### Detector

Nitrogen-phosphorous detector. For Figure K an FID was used in place of a nitrogen-phosphorous detector.

### Qualifications

N-Nitrosodiphenylamine decomposes in the GC injection port and is detected as diphenylamine. Therefore, diphenylamine in the sample must be removed by Florisil column chromatography or it will interfere with the analysis.

Figure K. Nitrosamines by EPA Method 607



## GC Method 608 – Organochlorine Pesticides and PCBs

### Sample Concentration

Solvent extraction with methylene chloride followed by concentration to 1 mL.

### GC Column

2m x 4mm ID glass containing 1.5% SP-2250/1.95% SP-2401 on 100/120 SUPELCOPORT.

### Confirmational Column

2m x 2mm ID glass containing 3% OV-1 on 100/120 SUPELCOPORT. (Supelco Equivalent: 3% SP-2100 on 100/120 SUPELCOPORT.)


### Detector

Electron capture detector.

### Qualifications

The extract must be cleaned by Florisil column chromatography to eliminate possible phthalate interference (5).


Figure L. Pesticides by EPA Method 608





**GC Method 609 – Nitroaromatics and Isophorone**  
**Figure M. Nitroaromatics and Isophorone by EPA Method 609**

Column: 1.5% SP-2250/1.95% SP-2401 on 100/120 SUPELCOPORT, 2m x 4mm ID glass  
 Cat. No.: 11947 (20g packing)  
 Carrier: nitrogen, 44mL/min  
 Det.: FID (nitrobenzene, isophorone) or ECD (dinitrotoluenes)  
 Inj.: 2µL, 0.2mg/mL each analyte (nitrobenzene, isophorone) or 0.5µL, 2ng/µL each analyte (dinitrotoluenes)



**Sample Concentration**

Solvent extraction with methylene chloride followed by concentration to 1mL.

**GC Column**

4' x 4m ID glass containing 1.95% QF-1/1.5% OV-17 on 80/100 Gas-Chrom® Q. (Supelco Equivalent: 2m x 4mm ID glass containing 1.5% SP-2250/1.95% SP-2401 on 100/120 SUPELCOPORT.) Two temperatures, 85°C and 145°C, are used.

**Confirmational Column**

10' x 4mm ID glass containing 3% OV-101 on 80/100 Gas-Chrom Q. (Supelco Equivalent: 3% SP-2100 on 100/120 SUPELCOPORT.)

**Detector**

FID for isophorone and nitrobenzene; electron capture detector for nitrotoluenes.

**Qualifications**


The Supelco equivalent column is two meters long, rather than four feet long, to permit its application to Methods 606, 608, and 613.

**GC or LC Method 610 – Polynuclear Aromatic Hydrocarbons (PAHs)**


**Figure N. PAHs by HPLC by EPA Method 610**

Columns: SUPELCO SILLIC-PAH  
 Cat. No.: 58229 (top) or 59133  
 Mobile Phase: water (A):acetonitrile (B) gradient top – 40% B (5 min) to 100% B over 25 min, hold bottom – 60% B (0.3 min) to 100% B over 2.7 min, hold 1 min  
 Flow Rate: 1mL/min (top) or 3.0mL/min  
 Det.: UV, 254nm  
 Inj.: 2µL methanol/methylene chloride containing 0.4µg each analyte

25cm x 4.6mm column, 5µm particles




5cm x 4.6mm column, 3µm particles



**Figure O. PAHs by GC by EPA Method 610**

Column: 3% SP-2250 on 100/120 SUPELCOPORT, 2m x 2mm ID glass  
 Cat. No.: 11744 (20g packing)  
 Oven: 100°C (4 min) to 280°C at 8°C/min, hold  
 Carrier: nitrogen, 40mL/min  
 Det.: FID  
 Inj.: 1µL methanol/methylene chloride containing 0.2µg each analyte



## Method 610 (contd.)

### Sample Concentration

Solvent extraction with methylene chloride followed by concentration to 1 mL.

### HPLC Column

25cm x 4.6mm ID SS containing HC-ODS Sil-X reversed phase. (Supelco Equivalent: 25cm x 4.6mm ID SS containing SUPELCOSIL LC-PAH.)

### Confirmational Column

None listed.

### Detector

Fluorescence.

### GC Column

2m x 2mm ID glass containing 3% OV-17 on 100/120 Chromosorb W AW-DMCS. (Supelco Equivalent: 3% SP-2250 on 100/120 SUPELCOPORT.)

### Confirmational Column

None listed.

### Detector

FID.


### Qualifications

The HPLC procedure is recommended for complete resolution of the PAHs; the GC procedure cannot adequately resolve four pairs of these compounds. The EPA method suggests using a 25cm x 4.6mm ID HPLC column containing 5µm packing particles for this analysis, but shorter SUPELCOSIL LC-PAH columns (15cm x 4.6mm ID, 5µm particles or 5cm x 4.6mm ID, 3µm particles) can greatly reduce the analysis time, as Figure N shows.

## GC Method 611 – Haloethers

### Figure P. Haloethers by EPA Method 611

Column: 3% SP-1000 on 100/120 SUPELCOPORT, 2m x 2mm ID glass  
Cat. No.: 11746 (20g packing)  
Oven: 60°C (2 min) to 230°C at 8°C/min, hold  
Carrier: helium, 40mL/min  
Det.: FID  
Inj.: 1µL methanol containing 0.2mg/mL each analyte



795-0190

## Sample Concentration

Solvent extraction with methylene chloride followed by concentration to 1 mL.

### GC Column

2m x 2mm ID glass containing 3% SP-1000 on 100/120 SUPELCOPORT.

### Confirmational Column

2m x 2mm ID glass containing 60/80 Tenax.

### Detector

Hall electroconductivity. An FID was used to obtain Figure P.


### Qualifications

The sample should be cleaned on Florisil if interferences are suspected to exist in the sample.

## GC Method 612 – Chlorinated Hydrocarbons

### Figure Q. Chlorinated Hydrocarbons by EPA Method 612

Column: 1.5% OV-1/1.5% OV-225 on 80/100 SUPELCOPORT, 2m x 2mm ID glass  
Cat. No.: custom-prepared  
Carrier: argon:methane, 95:5, 30mL/min  
Det.: ECD  
Inj.: 2µL hexane containing 0.1-10ng/µL each analyte



795-0191

## Sample Concentration

Solvent extraction with methylene chloride followed by concentration to 1 mL.

### GC Column

2m x 2mm ID glass containing 1.5% OV-1/2.4% OV-225 on 80/100 SUPELCOPORT (replaces 1.5% OV-1/1.5% OV-225 packing used for Figure Q.) Two temperatures, 75°C and 160°C, are used.

### Confirmational Column

None listed.

### Detector

Electron capture detector.

### Qualifications


Sample cleanup, if required, is conducted on a Florisil column.



## GC Method 613—2,3,7,8-Tetrachlorobenzo-p-dioxin (TCDD)

### Figure R. 2,3,7,8-TCDD and Other TCDD Isomers by EPA Method 613

Column: SP-2331 capillary column, 60m x 0.32mm ID, 0.2µm film  
 Cat. No.: 24105  
 Oven: 200°C (1 min) to 250°C at 3°C/min  
 Carrier: helium, 30cm/sec  
 Det.: ECD  
 Inj.: 0.2µL n-dodecane containing 0.2ng each isomer, splitless injection, 250°C



713-1065

#### Sample Concentration

Solvent extraction with methylene chloride followed by concentration to 1mL. The sample is cleaned on a silica gel column, then on an alumina column.

#### GC Column

60m x 0.32mm ID capillary column containing SP-2331.

#### Detector

Mass spectrometer.

#### Qualifications

A separate method was established for TCDD analysis because of the extreme toxicity of this compound. This method should be reviewed thoroughly before it is used.

A 60m x 0.32mm ID capillary column containing SP-2331 resolves 2,3,7,8-TCDD from other TCDD isomers. For additional information on this analysis, request Application Note 113.

#### Quantification of Priority Pollutants by GC or LC

The EPA also has recommended calibration procedures for the GC or LC analyses of priority pollutants. Prior to any calibration, the chromatographer must demonstrate that the system is operational by injecting a known mixture of the components to be analyzed.

A blank containing all reagents used is carried through the extraction and cleanup procedures and is analyzed to ensure the absence of interferences. Concentrations of calibration standards are chosen to bracket the expected concentration of the compound in the sample. These standards establish the sensitivity limit of the detector and the linear range of the analytical system for each component.

Once a graph of the standard concentrations vs. response is established, it is checked daily to ensure that no significant variation has occurred in the standard curve. Many of the methods recommend that surrogate standards be included in all standards, samples, and blanks for evaluation of the recovery and precision of the sample workup and analysis. If doubt exists about the identity of a peak in the chromatographic analysis, mass spectrometry should be used for confirmation.

Internal standards, surrogates, and calibration standards required in the GC-MS or the GC or LC methods are listed in the chemical standards section of the Supelco catalog.


#### Use a Capillary Column to Analyze Difficult or Complex Samples

In US EPA Method 625, capillary columns may be used as an alternative to packed columns if they meet the QC criteria listed in Sections 8.2, 12, and 13.1 (see page 3). 0.25mm ID, 0.32mm ID, and 0.53mm ID SPB™-5 capillary columns provide similar resolution of acid and base/neutral pollutants for Method 625 (Figure S). In addition, you can use 0.53mm ID columns in packed column chromatographs. In a GC-MS system with a high capacity pump or jet separator, a 0.53mm ID column combines large sample capacity (up to 2000ng/component) with rapid analysis and good resolution.

### Figure S. Acid and Base/Neutral Pollutants

Column: PTE™-5 capillary column, 30m x 0.25mm ID, 0.25µm film  
 Cat. No.: 24135-U  
 Oven: 35°C (4 min) to 300°C at 10°C/min, hold 10 min  
 Carrier: helium, 40cm/sec (set at 250°C)  
 Det.: MS (mass range: 35-450 m/z; 1.18 sec/scan)


N-Nitrosodimethylamine	p-Chloro-m-cresol	4. Phenanthrene-d <sub>10</sub>
Phenol	Hexachlorocyclopentadiene	Phenanthrene
2-Chlorophenol	2,4,6-Trichlorophenol	Anthracene
bis(2-Chloroethyl) ether	2-Chloronaphthalene	Di-n-butyl phthalate
1,3-Dichlorobenzene	Acenaphthylene	Fluoranthene
1,4-Dichlorobenzene-d <sub>4</sub>	Dimethylphthalate	Benzidine
1,4-Dichlorobenzene	2,6-Dinitrotoluene	Pyrene
1,2-Dichlorobenzene	3. Acenaphthene-d <sub>10</sub>	Benzyl butyl phthalate
bis(2-Chloroisopropyl) ether	Acenaphthene	Benzo(a)anthracene
Hexachloroethane	2,4-Dinitrophenol	3,3'-Dichlorobenzidine
N-Nitroso-di-n-propylamine	4-Nitrophenol	5. Chrysene-d <sub>12</sub>
Nitrobenzene	2,4-Dinitrotoluene	Chrysene
Isophorone	Diethyl phthalate	bis(2-Ethylhexyl)phthalate
2-Nitrophenol	Fluorene	Di-n-octyl phthalate
2,4-Dimethylphenol	4-Chlorophenylphenyl ether	Benzo(b)fluoranthene
bis(2-Chloroethoxy)methane	4,6-Dinitro-o-cresol	Benzo(k)fluoranthene
2,4-Dichlorophenol	N-Nitrosodiphenylamine	Benzo(a)pyrene
1,2,4-Trichlorobenzene	1,2-Diphenylhydrazine	Indeno(1,2,3-cd)pyrene
2. Naphthalene-d <sub>8</sub>	4-Bromophenylphenyl ether	Dibenzo(a,h)anthracene
Naphthalene	Hexachlorobenzene	Benzo(ghi)perylene
Hexachlorobutadiene	Pentachlorophenol	



797-0665

0.25mm ID and 0.53mm ID SPB-608 capillary columns offer excellent resolving power and inertness for separating the pesticides listed in US EPA Method 608 (Figure T). In addition to routine testing for efficiency and inertness, each SPB-608 column is tested to ensure minimal breakdown of DDT and endrin, as required by the EPA method. Furthermore, you can use a 0.53mm ID column in instruments designed for packed columns.

**Figure T. Chlorinated Pesticides**






\*A halogenated gas such as chloromethane or chlorotrifluoromethane is used to determine linear velocity by ECD.

**Wide Bore Capillary Columns Are Compatible with Purge and Trap Devices**

Our 0.53mm ID and 0.75mm ID VOCOL™ columns are specifically manufactured and tested to rapidly separate the volatile pollutants listed in EPA Methods 601, 602, and 624 for wastewater (Figure U) and 502.2 and 524.2 for drinking water. Because these columns have high optimum carrier gas flow rates, they accept samples directly from purge and trap devices. VOCOL columns are compatible with less sensitive detectors and concentration-dependent hyphenated techniques, such as GC-FTIR. The 0.53mm ID columns have a capacity of about 2000ng/sample component; the 0.75mm ID columns have a capacity of about 10,000ng/component.

**Figure U. Volatile Priority Pollutants**



## Ordering Information:

### Packings and Columns for USEPA Wastewater Analyses

For chemical standards for these methods, refer to the current Supelco catalog.

EPA Method No.	Class of Pollutants	Packings and Columns	Cat. No.
<b>Gas or Liquid Chromatography</b>			
601	Purgeable Halocarbons	60/80 Carbo-pack B/1% SP-1000, 15g	11815
		<i>Packing for Confirmation Column:</i> n-Octane on 80/100 Porasil C, 75cc	11733-U
		<i>Capillary GC Columns:</i> VOCOL, 60m x 0.75mm ID glass VOCOL, 30m x 0.53mm ID fused silica	NA* 25320-U
602	Purgeable Aromatics	5% SP-1200/1.75% Bentone 34 on 100/120 SUPELCOPORT, 20g	12134
		<i>Packing for Confirmation Column:</i> 5% 1,2,3-TCEP on 60/80 Chromosorb W AW, 20g	11765-U
		<i>Capillary GC Columns:</i> VOCOL, 60m x 0.75mm ID glass VOCOL, 30m x 0.53mm ID fused silica	NA* 25320-U
603	Acrolein & Acrylonitrile	80/100 Porapak QS, 24g	20343
		<i>Packing for Confirmation Column:</i> 80/100 Chromosorb 101, 50g	20214
604	Phenols	<i>Free Phenols:</i> 1% SP-1240-DA on 100/120 SUPELCOPORT, 20g	11832
		<i>Derivatives:</i> 5% SP-2250 on 80/100 SUPELCOPORT, 20g	11737
		<i>Capillary GC Column:</i> SPB-5, 15m x 0.53mm ID fused silica	25304
605	Benzidines	<i>HPLC Column:</i> SUPELCOSIL LC-1, 25cm x 4.6mm	58296
		<i>Capillary GC Column:</i> SPB-5, 15m x 0.53mm ID fused silica	25304
606	Phthalates	1.5% SP-2250/1.95% SP-2401 on 100/120 SUPELCOPORT, 20g	11947
		<i>Packing for Confirmation Column:</i> 3% SP-2100 on 100/120 SUPELCOPORT, 20g	11738
		<i>Capillary GC Column:</i> SPB-5, 15m x 0.53mm ID fused silica	25304
607	Nitrosamines	10% CARBOWAX 20M/2% KOH on 80/100 Chromosorb W AW, 20g	11739
		<i>Capillary GC Column:</i> SPB-5, 15m x 0.53mm ID fused silica	25304
608	Organochlorine Pesticides & PCBs	1.5% SP-2250/1.95% SP-2401 on 100/120 SUPELCOPORT, 20g	11947
		<i>Packing for Confirmation Column:</i> 3% SP-2100 on 100/120 SUPELCOPORT, 20g	11738
		<i>Capillary GC Columns:</i> SPB-608, 30m x 0.25mm ID fused silica SPB-608, 15m x 0.53mm ID fused silica	24103-U 25310-U
609	Nitroaromatics & Isophorone	1.5% SP-2250/1.95% SP-2401 on 100/120 SUPELCOPORT, 20g	11947
		<i>Packing for Confirmation Column:</i> 3% SP-2100 on 100/120 SUPELCOPORT, 20g	11738
		<i>Capillary GC Column:</i> SPB-5, 15m x 0.53mm ID fused silica	25304
610	PAHs	3% SP-2250 on 100/120 SUPELCOPORT, 20g	11744
		<i>HPLC Column:</i> SUPELCOSIL LC-PAH, 25cm x 4.6mm	58229
		<i>Capillary GC Column:</i> SPB-5, 30m x 0.53mm ID fused silica	25305-U
611	Haloethers	3% SP-1000 on 100/120 SUPELCOPORT, 20g	11746
		<i>Capillary GC Column:</i> SPB-5, 15m x 0.53mm ID fused silica	25304

## Packings and Columns for USEPA Wastewater Analyses (contd.)

EPA Method No.	Class of Pollutants	Packings and Columns	Cat. No.
612	Chlorinated Hydrocarbons	1.5% OV-1/2.4% OV-225 on 80/100 SUPELCOPORT <i>Capillary GC Column:</i> SPB-5, 15m x 0.53mm ID fused silica	custom  25304
613	2,3,7,8-TCDD	<i>Capillary GC Columns:</i> SP-2331, 60m x 0.25mm ID fused silica SP-2331, 60m x 0.32mm ID fused silica	24104-U 24105
<b>GC/MS Methods</b>			
624	Purgeable Halocarbons	60/80 Carbowax B/1% SP-1000, 15g <i>Packing for Confirmation Column:</i> 60/80 Carbowax C/0.2% CARBOWAX 1500, 15g <i>Packing for Precolumn:</i> 3% SP-1000 on 60/80 Chromosorb W AW, 20g <i>Capillary GC Columns:</i> VOCOL, 60m x 0.75mm ID glass VOCOL, 30m x 0.53mm ID fused silica	11815  11826  11741  NA* 25320-U
625	Acids (Phenols), Base-Neutrals, Organochlorine Pesticides, PCBs	3% SP-2250 on 100/120 SUPELCOPORT, 20g <i>Acids:</i> 1% SP-1240-DA on 100/120 SUPELCOPORT, 20g <i>Acids, Base-Neutrals, Organochlorine Pesticides, PCBs:</i> PTE-5, 30m x 0.25mm ID fused silica (0.25µm phase film)	11756  11832  24135-U

\*No longer available; we recommend a 0.53mm ID fused silica column for this application.

**For chemical standards for these methods, refer to the current Supelco catalog.**

### References

1. *Sampling and Analysis Procedures for Screening of Industrial Effluents for Priority Pollutants* USEPA Effluent Guidelines Division, Washington (April, 1977).
2. *Federal Register* Vol. 44, No. 233, Dec. 3, 1979.
3. *Federal Register* Vol. 45, No. 98, May 19, 1980.
4. Bellar, T.A., and J.J. Lichtenberg, *J. Am. Water Works Assoc.*, **66**: 739-744 (1974).
5. *Manual of Analytical Methods for the Analysis of Pesticides in Human and Environmental Samples* EPA-600/8-80-038, USEPA, Health Effects Research Laboratory, Research Triangle Park, NC (1980).

References not available from Supelco.

### Trademarks

Carbowax, PTE, SP, SPB, SUPELCOPORT, SUPELCOSIL, VOCOL – Sigma-Aldrich Co.  
Bentone – National Lead Company, Baroid Sales Division  
CARBOWAX – Union Carbide Corporation  
Chromosorb – Celite Corp.  
Durapak, Porapak, Porasil – Waters Associates, Inc.  
Florisil – U.S. Silica Co.  
Gas-Chrom – Applied Science Laboratories, Inc.  
LiChrosorb – EM Laboratories, Inc.  
OV – Ohio Valley Specialty Chemical Company  
Tenax – Enka Research Institute Arnhem

Fused silica columns manufactured under HP US patent no. 4,293,415.

BULLETIN 775

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website ([www.sigma-aldrich.com](http://www.sigma-aldrich.com)), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA • Sigma-Aldrich de Argentina, S.A. • Buenos Aires 1119 AUSTRALIA • Sigma-Aldrich Pty. Ltd. • Castle Hill NSW 2154 AUSTRIA • Sigma-Aldrich Handels GmbH • A-1110 Wien  
BELGIUM • Sigma-Aldrich N.V./S.A. • B-2880 Bornem BRAZIL • Sigma-Aldrich Quimica Brasil Ltda. • 01239-010 São Paulo, SP CANADA • Sigma-Aldrich Canada, Ltd. • 2149 Winston Park Dr., Oakville, ON L6H 6J8  
CZECH REPUBLIC • Sigma-Aldrich s.r.o. • 186 00 Praha 8 DENMARK • Sigma-Aldrich Denmark A/S • DK-2665 Vallensbaek Strand FINLAND • Sigma-Aldrich Finland/YA-Kemia Oy • FIN-00700 Helsinki  
FRANCE • Sigma-Aldrich Chimie • 38297 Saint-Quentin-Fallavier Cedex GERMANY • Sigma-Aldrich Chemie GmbH • D-82041 Deisenhofen GREECE • Sigma-Aldrich (o.m.) Ltd. • Ilioupoli 16346, Athens  
HUNGARY • Sigma-Aldrich Kft. • H-1067 Budapest INDIA • Sigma-Aldrich Co. • Bangalore 560 048 IRELAND • Sigma-Aldrich Ireland Ltd. • Dublin 24 ISRAEL • Sigma Israel Chemicals Ltd. • Rehovot 76100  
ITALY • Sigma-Aldrich s.r.l. • 20151 Milano JAPAN • Sigma-Aldrich Japan K.K. • Chuo-ku, Tokyo 103 KOREA • Sigma-Aldrich Korea • Seoul MALAYSIA • Sigma-Aldrich (M) Sdn. Bhd. • Selangor  
MEXICO • Sigma-Aldrich Química S.A. de C.V. • 50200 Toluca NETHERLANDS • Sigma-Aldrich Chemie BV • 3330 AA Zwijndrecht NORWAY • Sigma-Aldrich Norway • Torshov • N-0401 Oslo  
POLAND • Sigma-Aldrich Sp. z o.o. • 61-663 Poznań PORTUGAL • Sigma-Aldrich Quimica, S.A. • Sintra 2710 RUSSIA • Sigma-Aldrich Russia • Moscow 103062 SINGAPORE • Sigma-Aldrich Pte. Ltd.  
SOUTH AFRICA • Sigma-Aldrich (pty) Ltd. • Jet Park 1459 SPAIN • Sigma-Aldrich Quimica, S.A. • 28100 Alcobendas, Madrid SWEDEN • Sigma-Aldrich Sweden AB • 135 70 Stockholm  
SWITZERLAND • Supelco • CH-9471 Buchs UNITED KINGDOM • Sigma-Aldrich Company Ltd. • Poole, Dorset BH12 4QH  
UNITED STATES • Supelco • Supelco Park • Bellefonte, PA 16823-0048 • Phone 800-247-6628 or 814-359-3441 • Fax 800-447-3044 • email: [supelco@sial.com](mailto:supelco@sial.com)

H

Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.

BNV