

IDENTIFICATION AND QUANTIFICATION OF HOST CELL PROTEIN IMPURITIES IN HIGH PURITY MONOCLONAL ANTIBODIES DOWN TO 1 PPM: AN INTER-LABORATORY STUDY

Catalin Doneanu¹, Scott J. Berger¹, Malcolm Anderson², Brad Williams³, Matt Lauber¹, Asish Chakraborty¹ and Weibin Chen¹
¹Waters Corporation, Milford, MA USA; ²Wilmslow, UK; ³Beverly, MA USA

To learn more about LC-MS based approaches for HCP analysis, please visit: www.waters.com/HCP

INTRODUCTION

- Host cell protein impurities (HCPs) in biopharmaceuticals raise immunogenicity, stability, and drug efficacy concerns, even when present at ppm levels.
- Regulatory guidances mandate that downstream processes are optimized to clear HCPs to the lowest levels that can be practically achieved, and that sensitive assays are in place to demonstrate clearance.
- The ELISA and immunoassay based HCP detection methods widely adopted today typically report total HCP burden, and often lack broad specificity required for comprehensive HCP monitoring.
- LC-MS based approaches for identification and quantification of HCPs in biopharmaceuticals have now been realized [Refs. 1-5, 7] to enable comprehensive HCP discovery and monitoring of individual HCP protein impurities.
- In this study, a novel analytical approach (2DLC-HDMS^E), that combines high pH/low pH reversed phase (RP/RP) separations with orthogonal ion mobility (IM) gas phase separations, is presented as a means to maximize analytical peak capacity and achieve robust low ppm HCP analysis.
- Our results indicate that LC-MS assays are now able to attain comparable sensitivity to traditional immunological based HCP assays, while offering discrete HCP identification and monitoring.

METHODS

Sample preparation

Three murine derived mAbs: a NIST reference mAb standard; Remicade, a commercial Infliximab (Innovator); and Inflectra (its approved biosimilar) biotherapeutics were analyzed. Samples (2.5 mg mAb) were denatured with RapiGest™ SF (60 °C, 15 min), DTT reduced (60 °C, 1h), alkylated with IAM (RT, 30 min) and digested with a mixture of Lys-C and porcine trypsin (Promega) overnight. Four protein digest standards were spiked post-digestion such that the amounts loaded on-column using a 250 µl injection were: 5,000 fmol ADH (yeast alcohol dehydrogenase), 1,000 fmol PHO (rabbit phosphorylase b), 250 fmol BSA (bovine serum albumin) and 50 fmol ENL (yeast enolase).

2D-LC configuration

ACQUITY™ UPLC® M-Class system with 2D technology was used for reversed-phase/reversed-phase (RP/RP) 2D-LC separations of the antibody digest.

First Dimension (1D) pH 10: 1.0 x 50 mm XBridge C₁₈ column (5 µm particles), 10 µL/min flow. Mobile phase: 20 mM ammonium formate in water (Solvent A) and ACN (Solvent B).

Online dilution (1:10) of the eluent from 1D prior to the trap column.

Trap column: 0.3 x 50 mm packed with 5-µm Symmetry C₁₈ particles.

Second Dimension (2D) pH 2.4: 0.3 x 150 mm analytical column CSH C₁₈ 1.7 µm, kept at 60 °C and operated at 10 µL/min.

Fractions were eluted in ten steps of increasing organic strength, with each step diluted online to 1:10 ratio with 0.1% TFA in water (pH=2.1) before trapping. Low pH separations in the second chromatographic dimension used a 40 min gradient from 3 to 40% acetonitrile (0.1% formic acid).

MS and data processing

Data-independent acquisition (MS^E) were performed on a SYNAPT G2-S mass spectrometer:

- Acquisition cycle time was 0.5 sec, m/z range: 100-1990.
- Fixed CE at 5 V for MS scans; CE ramp 20-45 eV for MS^E scans.
- For IMS (HDMS^E) enabled acquisition, a fixed wave velocity (650 m/s) and fixed wave height (40 V) applied.
- Ion mobility specific CE were applied during HDMS^E acquisition [6]
- Data Processing: ProteinLynx Global Server (PLGS) 3.0.2




Figure 1. Fluidic configuration for two-dimensional high pH/ low pH RP/RP chromatography, employing pre-trap on-line dilution.

SYNAPT G2-S HDMS system coupled to an ACQUITY UPLC M-Class with 2D technology.

ANALYTICAL DEVELOPMENT

Reproducibility of the 2D LCMS^E HCP Discovery Methodology




Figure 2. LC separation reproducibility of the 2DLC chromatographic system. Extracted mass chromatograms were generated for the monoisotopic peaks of two low-abundance HCP peptides. The peptides reproducibly elute in the same fraction at consistent retention times over the course of the triplicate HCP discovery experiment.

Ion Mobility (HDMS^E) acquisition simplifies MS^E spectral complexity of HCP peptide fragments of an HCP Peptide Precursor Ion (m/z 877.92, 2+ Charge)




Figure 3. Enabling Ion Mobility during 2DLC-MS^E acquisition (i.e. 2DLC-HDMS^E) separates peptide MS ions in the gas phase, and simplifies the resulting MS^E fragmentation spectra. The complex MS^E fragmentation spectra (TOP) seen without IMS separation, are greatly simplified by enabling IMS (MIDDLE), enabling clear ion detections. Confirmatory targeted MS/MS (2 m/z isolation window) independently validates fragment ion assignments. Example: NSFL1 Cofactor p47 peptide SYQDPSNAQFLESIR from NIST mAb.

HCP ANALYSIS OF NIST MAB REFERENCE STANDARD

Detecting low level (sub 20 PPM) HCPs in the NIST mAb Reference Standard from multiple peptides

No	Average HCP conc. (ppm)	HCP Identification	Sequence Coverage (%)	MW (kDa)	Peptide Sequence	Monoisotopic peak (charge)	RT (min)	Drift time (2D) (min)	Fraction no.
1	18	Catenin-interactor 1	7.3	100	ELVHAGDQDRLR SLNHYHFDHSDAGDINR YVGGSDYDSEGR	538.91 (+2) 475.01 (+3) 485.31 (+2)	11.7 10.2 11.8	3.3 1.9 3.5	1
2	14	Low affinity IgG gamma Fc region receptor	32.7	36.7	SGQKAKPTK TLDGKPTVYDQPK EMGELTLPFVYGR RFTFKK HPEALDELTDFQDQRI ATFDVYDGLDQYDGR EDYVLYLSEGLNHPNDSYDQYDGR	579.79 (+2) 578.47 (+3) 433.17 (+2) 313.26 (+2) 492.36 (+1) 482.36 (+2) 322.29 (+2)	12.2 6.2 12.2 12.7 11.7 11.7 11.7	1.2 3.5 5.1 3.2 4.1 4.1 6.6	2
3	13	Stress induced phosphoprotein 1	9.0	42.4	LARFVQDLER ALLAGDQDQDQYDQAK ALDLDQDQKADQDGR	384.39 (+2) 485.31 (+2) 433.17 (+2)	12.2 12.7 11.7	4.3 3.2 3.9	3
4	12	Glucose-6-phosphate isomerase	15.6	62.7	ELVHAGDQDRLR HFNALSTKAK AMVHALR VWFYDNDQTHAK TSLGSDYDSEGR	551.26 (+2) 394.82 (+2) 429.78 (+2) 529.42 (+2) 492.36 (+2)	12.7 12.7 7.3 12.7 11.7	3.5 3.6 2.7 3.4 3.9	4
5	7	NSFL1 cofactor p47	28.2	40.7	EVAVYGTEDR EGAVPEESYVNGR MALKERDQGR LAAGDQDQDQYDQAK SYQDPSNAQFLESIR SPVLEKQDQK DLSHQDQEELEFEGGR	678.44 (+3) 824.89 (+2) 695.80 (+2) 538.91 (+2) 877.92 (+2) 438.81 (+2) 365.29 (+2)	16.1 10.2 10.9 11.1 15.5 15.1 12.6	4.5 4.1 3.2 3.0 3.3 3.5 4.4	5
6	7	Beta-2-microglobulin	22.7	13.8	TYWQDQK VMSQDQK TRQDQYR	541.24 (+2) 345.24 (+2) 346.29 (+2)	16.9 17.9 17.9	3.2 3.2 3.3	6
7	5	Serine/threonine-rich signaling factor 7	13.7	30.8	TPVYDQDQYDQAK NPPQGFVYDPEQDQDQYDGR VLDRTYDQDQK	346.31 (+2) 791.37 (+1) 514.44 (+1)	15.3 26.2 12.9	3.3 4.4 3.9	7
8	5	Transketolase	8.3	67.4	TSPNPNVYDQDQYDQAK NMAQDQDQYDQAK VLDRTYDQDQK	827.31 (+1) 492.36 (+2) 461.99 (+2)	11.3 12.9 12.9	3.1 3.5 3.5	8
9	4	Adenylyl kinase 2	11.3	35.5	MEGFNAR AVGLQDQDQK	481.25 (+2) 492.36 (+2)	17.7 15.5	2.9 3.0	9
10	3	Prostaglandin reductase 1	8.2	15.4	TPQDQDQYDQAK LTPQDQDQYDQAK	433.17 (+2) 355.24 (+2)	17.2 17.2	3.4 3.4	10
11	2	Flavin reductase (NADPH)	12.6	22.3	YVGGSDYDSEGR IAIFGATGR	453.26 (+2) 453.26 (+2)	11.1 11.1	3.4 3.4	11
12	1	Penicillinase 5	25.7	21.8	YVGGSDYDSEGR THQDQDQYDQAK ALNWFQDQDQYDQAK	453.26 (+2) 788.41 (+1) 790.08 (+1)	11.1 12.2 29.1	3.4 4.5 4.0	12

Table I. List of peptides identified during 2D-LC/HDMS^E analysis of NIST mAb identifying HCPs at levels below 20 ppm. Detailed information is provided for each peptide regarding the observed monoisotopic m/z, charge state, ion mobility drift time, peptide LC retention time from the 2nd-dimension separation, as well as peptide elution fraction from high pH 1st dimension separation.

Consistency of Hi3 Quantification of NIST mAb Reference Standard HCPs across three laboratories

No	Accession crt	Protein Description	Identified without IMS	Average MW (kDa)	Measured HCP Concentration (ppm)
1	P00489	Glycogen phosphorylase rabbit (PHO)-1000 fmoles	Y	97.1	155 155 155
2	P00330	Alcohol dehydrogenase yeast (ADH)-5000 fmoles	Y	36.7	106 125 164
3	P05664	Fructose biphosphate aldolase A isoform	Y	39.3	113 126 109
4	P02769	Fructose biphosphate aldolase C isoform	Y	39.4	106 93 91
5	P02769	Bovine serum albumin (BSA)-250 fmoles	Y	66.3	35 28 20
6	Q99K99	Clathrin interactor	N	68.5	22 19 7
7	P08101	Low affinity IgG gamma Fc region receptor	Y	36.7	10 28 5
8	Q08064	Stress induced phosphoprotein 1	N	62.5	10 16 13
9	P06745	Glucose-6-phosphate isomerase	N	62.7	7 15 13
10	P01887	Beta-2-microglobulin	N	13.8	5 14 1 7
11	Q8C244	NSFL1 cofactor p47	N	40.7	7 10 5 7
12	P00924	Enolase 1 yeast (ENL)-50 fmoles	N	46.6	7 7 3 6
13	Q8L197	Serine/arginine-rich signaling factor 7	N	30.6	6 7 3 5
14	P40142	Transketolase	N	67.6	5 4 6 5
15	Q9W776	Adenylyl kinase 2	N	26.5	4 6 2 4
16	Q91Y89	Prostaglandin reductase 1	N	35.5	4 3 3 3
17	Q02302	Flavin reductase (NADPH)	N	22.2	1 3 1 2
18	P99029	Penicillinase 5	N	21.8	1 2 1 1

Table II. Fourteen NIST mAb HCPs were identified in common by all three laboratories performing 2DLC/HDMS^E acquisition. Individual HCP concentrations (range of 1–200 ppm) were calculated by the Hi3 methodology [Ref. 8] from three replicate injections, calibrated by the PHO protein standard, spiked in the NIST mAb post-digestion. Eleven HCPs (highlighted in yellow) required ion mobility based HDMS^E analysis for identification.



No	Accession crt	Protein Description	Sequence Coverage (%)	Average MW (kDa)	Amount on column fmoles	Concentration ng/mL	Concentration ppm	RSD (%)
1	P00489	Glycogen phosphorylase rabbit (PHO) - 1000 fmoles	41.2	97.1	1000	97	1942	92
2	P00330	Alcohol dehydrogenase yeast (ADH) - 5000 fmoles	44.6	36.7	1981	73	1454	69
3	P02769	Bovine serum albumin (BSA) - 250 fmoles	21.9	66.3	296	20	392	19
4	Q8GUQ1	Epidermal growth factor-like protein 8	8.2	33.3	792	26	527	25
5	Q8CBE3	WD repeat-containing protein 37	4.8	55.1	288	16	317	15
6	P00924	Enolase 1 yeast (ENL) - 50 fmoles	15.7	46.6	102	5	95	5
		Total				ng/mL	845	
		Total				ppm	40	
		mAb				ppm	99.99%	

No	Accession crt	Protein Description	Sequence Coverage (%)	Average MW (kDa)	Amount on column fmoles	Concentration ng/mL	Concentration ppm	RSD (%)
1	P00489	Glycogen phosphorylase rabbit (PHO) - 1000 fmoles	35.7	97.1	1000	97	3884	388
2	P00330	Alcohol dehydrogenase yeast (ADH) - 5000 fmoles	40.7	36.7	2073	76	3043	304
3	Q8GUQ1	Epidermal growth factor-like protein 8	8.2	33.3	687	23	915	92
4	P02769	Bovine serum albumin (BSA) - 250 fmoles	19.1	66.3	294	19	780	78
5	Q8CBE3	WD repeat-containing protein 37	4.8	55.1	142	8	313	31
6	P00924	Enolase 1 yeast (ENL) - 50 fmoles	18.3	46.6	58	3	108	11
		Total				ng/mL	1228	
		Total				ppm	123	
		mAb				ppm	99.99%	

Table III. Comparison of HCPs identified and quantified in Remicade and Inflectra. The same 2 HCPs were detected in common for both samples: epidermal growth-factor like protein 8 and WD repeat containing protein 37. Higher levels of the HCPs were detected for both proteins in the biosimilar version, for a total HCP burden roughly three times higher than the innovator.

References

1. Doneanu et al. *mAbs*, 2012, 4, 24-44.
2. Schenauer, *Anal Biochem*, 2012, 428, 150-157.
3. Protein Downstream Processing, *Methods in Molecular Biology* 1129, 2014, Editor NE Labrou, chapter 25, 341-350.
4. Zhang et al. *mAbs*, 2014, 6, 659-670.
5. Thompson et al. *Rapid Commun. Mass Spec*, 2014, 28, 855-860.
6. Distler et al. *Nat Meth*, 2014, 11, 167-170.
7. Doneanu, et al. *Anal Chem*, 2015, 87, 10283-10291. (Contains elements of the NIST mAb work in this poster)
8. Silva et al. *Mol Cell Proteom*, 2006, 144-127.

CONCLUSIONS

- An improved generic 2D-LC/HDMS^E assay was developed for quantitative HCP discovery analysis in biopharmaceuticals that incorporates ion mobility separations with 2D-LC/MS^E to extend assay sensitivity by an order of magnitude, allowing for the identification and quantification of HCPs at single digit ppm levels.
- A multi-laboratory trial of this methodology has established the quantitative robustness of the HCP discovery workflow.
- This methodology has been able to establish qualitative similarity of the HCP profile for an innovator biotherapeutic mAb and its approved biosimilar, but demonstrating quantitative differences in HCP levels.