


Solutions Offered by Shimadzu Corporation

Data Integrity Compliance in the Analytical Laboratory



● Strengthen Data Integrity with LabSolutions™ DB/CS

Data integrity refers to maintaining data completeness, consistency, accuracy, and reliability throughout the product life cycle. LabSolutions DB/CS software offers powerful support for ensuring compliance with data integrity requirements by comprehensively managing analytical instrument information.

● Strengthen Data Integrity for Various Instruments

LabSolutions DB/CS software enables stronger data integrity for all types of instruments in analytical laboratories, such as LC, GC, LC-MS, GC-MS, UV, FTIR, spectrofluorophotometers, protein sequencers, ICP-MS, AA, EDX, thermal analyzers, particle size analyzers, AG, TOC systems, Rigaku X-ray diffractometers and analytical balances.

● Configure Systems Optimized for the Given Laboratory Size

LabSolutions DB/CS systems can be optimized for a wide range of laboratory sizes and operating styles, from small standalone systems (with an instrument and computer) to networked or cloud-based systems with many analytical or measuring instruments.


Managing Data Securely

Systems can provide highly reliable data, including preventing data falsification by using a LabSolutions DB/CS database to centrally manage the data obtained from various analytical instruments in the laboratory.

- Instruments, users, and analytical data are managed in a database.
- Projects can be created in the database based on operating practices.

Commonizing User Management and User Privilege Settings for All Instruments


A wide variety of instruments other than LC and GC systems, such as LC-MS, GC-MS, UV, and FTIR systems, are also used in analytical laboratories. Managing the users and user privileges for all those instruments to ensure they are used appropriately is extremely tedious and time-consuming. LabSolutions DB/CS can reduce the time and trouble required for user management by commonizing the user management system for all instruments and specifying all user privilege settings.



Managing User Operations (Processes) Appropriately

Data integrity requires not only restricting access to laboratory instruments, but also assigning specific user privileges to laboratory personnel.

To ensure user privileges are specified appropriately, visualize the workflow by creating a process flow chart and then use LabSolutions DB/CS to register users and specify user privilege settings for implementing that workflow.




Example of User Privilege Groups and Workflow Recommended by LabSolutions DB/CS

Automating Report Creation


Spreadsheet software is often used, such as for quantitative calculations of analytical data, but that can prevent obtaining proper results due to errors or alterations, such as when transcribing analytical results, creating or storing files, or printing or storing reports.

The multi-data report function can create reports for the selected analytical results using the selected report template and then automatically create and save the report in the database. That eliminates the need to transcribe data, reduces the risk of errors, and can prevent data alterations or falsification. It also means reports with a combination of LC, GC, MS, FTIR, UV, analytical balance, and other instrument data can be created more efficiently using three easy steps.



Reducing Potential Risks in the Data Acquisition/Analysis Workflow - LabSolutions DB/CS Report Set

It is impossible to determine from simply looking at an acquired chromatogram or other analytical data whether the data was acquired using correct test parameter settings (such as instrument or data processing parameters) or rule out the possibility of the test parameters being falsified. Furthermore, log reviews are required to guarantee that electronic records do not include any fraudulent or suspicious data. By using LabSolutions DB/CS report sets, analytical data can be electronically managed and linked to metadata (results from processes that involve human intervention, such as specifying parameter settings or data analysis) and to the entire series of corresponding log data. In other words, metadata, log data, and analytical data can all be reviewed within a single report. In addition, by specifying all log events that should be reviewed in advance, all log events that need to be reviewed can be automatically extracted to enable efficiently checking for fraud or suspicious content in log data. Additionally, when review records are signed, they can be linked to data as a log review report. That means the entire series of process steps from data acquisition to data log reviews can be managed electronically.



Note: LabSolutions DB/CS report set functionality is compatible with all Shimadzu LC, GC, LC-MS, GC-MS, UV, FTIR, RF, SALD, TA, AG, TOC-L systems and Rigaku X-ray diffractometers.
 Note: Contact Shimadzu separately about the specifications and functionality of specific instruments related to strengthening data integrity.

For more information about Shimadzu solutions for strengthening data integrity, go to the following website.

<https://www.shimadzu.com/an/industries/small-molecule-pharmaceutical/data-integrity/index.html>

LabSolutions is a trademark of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.



Shimadzu Corporation
www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.
 Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".
 Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".
 Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.