

GLOBAL ANALYTICAL SOLUTIONS



GERSTEL

AppNote 2/1992

Applications of Mass Flow Controlled Multi Column Switching in On-Line Capillary GC/MS

Andreas Hoffmann, Ralf Bremer, Friedhelm Rogies
*Gerstel GmbH & Co.KG, Eberhard-Gerstel-Platz 1,
D-45473 Mülheim an der Ruhr, Germany*

Jacques A. Rijks^{1,2,3}, Jacek Staniewski^{2,3}

¹*Foundation Education and Development of Chromatography,
P.O. Box 4990, NL-5604 CD Eindhoven, The Netherlands*

²*Poznan Technical University, Institute of Chemical Technology and
Engineering, Pl. M. Skłodowskiej-Curie 2, 60-965 Poznan, Poland*

³*Eindhoven University of Technology, Laboratory of Instrumental
Analysis, P.O. Box 513, NL-5600 MB Eindhoven, The Netherlands*

KEYWORDS

PTV-injection, multi column switching, cryotrapping,
capillary GC-MS, industrial applications

ABSTRACT

The separation and analysis of low concentrations of organic compounds in complex sample matrices, such as petroleum products, waste and drinking water, food, beverages and pharmaceutical products is a rather complex analytical problem. Current methods are hampered by insufficient resolution obtained by single capillary columns even if they have rather high plate numbers.

In this paper the potential of a combination of programmed temperature sample introduction and mass flow controlled multi column dual oven capillary gas chromatography and on-line mass spectrometry will be discussed and illustrated.

The effect of cold trapping in between the columns for components with a moderate volatility will be demonstrated for different applications dealing with the determination of trace impurities in various main products such as gasoline, aromatics, steroids and aniline.

INTRODUCTION

Due to ever increasing demands for resolution, sample volume and detection limits in capillary GC nowadays, multi column switching is becoming increasingly important. Essentially these problems are related to the complexity of the samples and the required compatibility of sample size, input band width, column properties, detector specifications and the actual concentrations of the components of interest in the sample.

To realize adequate analysis a number of clean-up and/or enrichment steps, or solvent elimination for diluted samples prior to separation, are needed in order to achieve the required separation efficiency and detection limits. The identification by an on-line coupled MS-system is another breakpoint in the analysis of such complex mixtures, particularly if they elute in low concentrations close together with a major component. Multi column switching is, therefore, an attractive approach to solve these problems, particularly if high and low concentrations have to be determined simultaneously.

The aim of this paper is to demonstrate the possibilities of a combination of temperature programmed sample introduction, multi column dual oven capillary GC, coupled on-line with mass spectrometry in order to enhance identification. The effect of cryotrapping at the inlet of the second column will be emphasized.

EXPERIMENTAL

Instrumentation. The system consists of a combination of a manual or automatic injection (HP 7673, Hewlett-Packard, Avondale, USA), a temperature programmable cold injection system (CIS-3, Gerstel GmbH, Mülheim an der Ruhr, Germany), a multi column switching system (MCS A or MCS P, Gerstel GmbH, Mülheim an der Ruhr, Germany), 2 Ovens (HP 5890 series II, Hewlett-Packard, Avondale, USA) with a cryotrap system (CTS-1, Gerstel GmbH, Mülheim an der Ruhr, Germany) in between and a mass selective detector (HP 5971, Hewlett-Packard, Avondale, USA) as illustrated in **Figure 1**.




Figure 1. Schematic diagram of the applied system consisting of a manual or automated injector, a temperature programmable cold injection system with septumless sampling head (2), a GC (3) configured with monitor-FID, column switching device (4) and MCSA or MCSP-pneumatic, connected via a heated transferline with included cryotrap (5) to a second GC (6) with a mass selective detector (7).

The pneumatics of the multi column switching systems used in this study are presented in **Figure 2a** and **b**.




Figure 2a. Schematic design of a dual column MCS P system. Solvent flush prior to the first column NV1 and MFC2, both open.




Figure 2b. Schematic design of the multi column MCS A system. Solvent flush via splitline, MFC3 open.

The switching pneumatics of both systems are mass flow controlled. The carrier gas flow of the MCS P system is mass flow controlled, the carrier gas flow of the MCS A system is pressure controlled.

RESULTS

Multi column switching is the most powerful approach for analytical and preparative separation of components which can hardly or not at all be separated in single column systems. This is illustrated below for 5 different applications all dealing with the analysis of trace components, which elute very close to or together with a major component or the determination of trace impurities in highly pure products.

Example 2. Multi column analysis of m-,p-xylene in ethylbenzene.




Figure 4a. Precolumn chromatogram, 1 μ l, split ratio x:30, marked compound flushed.




Figure 4b. Main column chromatogram, without CTS.

Analytical conditions.

System:

MCS A, 2 columns, 2 ovens, with and without CTS

Columns:

Pre-column in GC-1	25 m Ultra-1 (Hewlett-Packard)
	$d_i = 0.32$ mm $d_f = 1.02$ μ m.
Main column in GC 2	30 m CP-Wax (Chrompack)
	$d_i = 0.32$ mm $d_f = 1.2$ μ m.

Pneumatics:

Carriergas	He	$p_i = 130$ kPa	split x:40
Control flow		$p_c = 50$ kPa	10 ml/min
FID	H ₂ , 30 ml/min	Air, 300 ml/min	
	N ₂ , 30 ml/min		

Temperatures:

CIS	60°C;	$\nearrow 260^\circ\text{C}$;	12°C/s.
Oven 1	70°C.		
Oven 2	50°C;	$\nearrow 180^\circ\text{C}$;	3°C/min.
CTS	200°C;	$\searrow -150^\circ\text{C}$;	12°C/s;
		$\nearrow 200^\circ\text{C}$;	12°C/s.

Detectors:

Monitor detector	FID
Main detector	FID

Compounds:

1. Benzene
2. Toluene
3. Ethylbenzene
4. m-Xylene
5. p-Xylene
6. Styrene




Figure 4c. Main column chromatogram, with CTS.

Example 4. Multi column analysis of impurities in aniline.




Figure 6a. Precolumn chromatogram, 1 μ l, split ratio x:10, marked compound flushed.




Figure 6b. Main column chromatogram (TIC), without CTS.

Analytical conditions.

System:

MCS A, 2 columns, 2 ovens, with and without CTS

Columns:

Pre-column in GC-1	25 m OV-17 (home made)
	$d_i = 0.32$ mm $d_f = 1.0$ μ m.
Main column in GC 2	50 m HP-1(Hewlett-Packard)
	$d_i = 0.32$ mm $d_f = 1.05$ μ m.

Pneumatics:

Carriergas	He	$p_i = 60$ kPa	split x:15
Control flow		$p_c = 45$ kPa	10 ml/min
FID	H ₂ , 30 ml/min	Air, 300 ml/min	
	N ₂ , 30 ml/min		

Temperatures:

CIS	60°C;	∇ 260°C;	12°C/s.
Oven 1	50°C;	∇ 200°C;	10°C/min.
Oven 2	40°C;	∇ 260°C;	10°C/min.
CTS	220°C;	∇ -150°C;	12°C/s;
		∇ 220°C;	12°C/s.

Detectors:

Monitor detector	FID
Main detector	MSD Scan 10-300 amu

Compounds:

1. Benzene
2. Cyclohexane
3. Cyclohexylamine
4. Cyclohexanol
5. Phenol
6. Aniline
7. Toluidine
8. Nitrobenzene
9. Dicyclohexylamine





Figure 6c. Main column chromatogram (TIC), with CTS.



GERSTEL GmbH & Co. KG

Eberhard-Gerstel-Platz 1
45473 Mülheim an der Ruhr
Germany

- +49 (0) 208 - 7 65 03-0
- +49 (0) 208 - 7 65 03 33
- gerstel@gerstel.com
- www.gerstel.com

GERSTEL Worldwide

GERSTEL, Inc.

701 Digital Drive, Suite J
Linthicum, MD 21090
USA

- +1 (410) 247 5885
- +1 (410) 247 5887
- sales@gerstelus.com
- www.gerstelus.com

GERSTEL AG

Wassergrabe 27
CH-6210 Sursee
Switzerland

- +41 (41) 9 21 97 23
- +41 (41) 9 21 97 25
- swiss@ch.gerstel.com
- www.gerstel.ch

GERSTEL K.K.

1-3-1 Nakane, Meguro-ku
Tokyo 152-0031
SMBC Toritsu-dai Ekimae Bldg 4F
Japan

- +81 3 5731 5321
- +81 3 5731 5322
- info@gerstel.co.jp
- www.gerstel.co.jp

GERSTEL LLP

10 Science Park Road
#02-18 The Alpha
Singapore 117684

- +65 6779 0933
- +65 6779 0938
- SEA@gerstel.com
- www.gerstel.com

GERSTEL (Shanghai) Co. Ltd

Room 206, 2F, Bldg.56
No.1000, Jinhai Road,
Pudong District

- Shanghai 201206
- +86 21 50 93 30 57
- china@gerstel.com
- www.gerstel.cn

GERSTEL Brasil

Av. Pascoal da Rocha Falcão, 367
04785-000 São Paulo - SP Brasil

- +55 (11)5665-8931
- +55 (11)5666-9084
- gerstel-brasil@gerstel.com
- www.gerstel.com.br

Information, descriptions and specifications in this Publication are subject to change without notice. GERSTEL, GRAPHPACK and TWISTER are registered trademarks of GERSTEL GmbH & Co. KG.

© Copyright by GERSTEL GmbH & Co. KG



Awarded for the active pursuit of environmental sustainability