

Table of contents

Summary

Overview:  
Dietary supplements and  
botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of  
supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

LC that takes your  
productivity to  
new heights

## HPLC - UV-Vis absorption and charged aerosol detection

Dietary supplements and botanical natural products applications notebook  
Complex samples, powerful chromatographic analysis



## Table of contents

## Summary

Overview:  
Dietary supplements and  
botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of  
supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

Reliable analysis of dietary supplements and [botanical](#) natural products is vital in ensuring the consistent quality and safety of raw materials and end products. However, the nature of these materials can make their analysis particularly demanding.

[Botanical](#) natural products, for example, are often complex mixtures. Their exact compositions may differ under the influence of variables such as genetics, cultivation conditions, adulteration, storage conditions, and processing methods.

Chromatographic analysis of both dietary supplements and [botanical](#) natural products presents specific challenges with respect to:

1. **Sample preparation:** the wide variety of sample matrices can affect the efficiency and reproducibility of extraction.
2. **Separation:** an exceedingly wide range of chemical structures will require different chromatographic approaches.
3. **Detection:** analytes lacking a [chromophore](#) are not detected by absorption detectors.
4. **Quantitation:** often there are no reference standards available.

This notebook offers a summary of applications and helpful information as to how high-performance liquid chromatography (HPLC) combined with charged aerosol detection, absorption detection or both can be used to meet the needs of scientists working with dietary supplements and [botanical](#) natural products.



## Table of contents

Summary

Overview:  
Dietary supplements and  
botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of  
supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

## What is meant by dietary supplements and botanical natural products?

### Dietary supplements

The [US FDA defines dietary supplements](#), in part, as “products taken by mouth that contain a dietary ingredient. Dietary ingredients include vitamins, minerals, amino acids, and [herbs](#) or [botanicals](#), as well as other substances that can be used to supplement the diet. Dietary supplements come in many forms, including tablets, capsules, powders, energy bars and liquids.”

### Natural products

The term natural product can be used in many ways. Some applications are specific, based on the [FDA's](#) definition of a chemical compound or substance produced by a living organism (bacteria, fungi, marine organisms or plants), that usually has a pharmacological or biological activity, and is of use in pharmaceutical drug discovery and drug design. Other definitions are broader, and natural product may be used as an umbrella term to describe complex plant-based preparations.



## Table of contents

Summary

Overview:  
Dietary supplements and  
botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of  
supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary



### Botanical natural products

[Botanical](#) natural products are referred to as supplements, [herbal](#) medicines, [herbal](#) drugs, herbs, [botanical](#) preparations, nutraceuticals, phytomedicines and [botanical](#) medicines. [Botanical](#) natural products are available to consumers in a variety of forms, including capsules containing raw or extracted material, extracts, teas (typically raw plant material that is extracted in hot water prior to use), tinctures (ethanolic extracts), and [traditional formulations](#), such as powders used in [traditional Chinese medicine](#) or [Ayurvedic](#) practices.

### Botanical drug products

The [FDA defines botanical drug products](#) as "...intended for use in the diagnosis, cure, mitigation, treatment or prevention of disease in humans..." and consisting of "...vegetable materials, which may include plant materials, algae, macroscopic fungi, or combinations thereof." [Botanicals](#) include drugs that are FDA-approved under the [botanical](#) drug pathway, while the category of plant-derived drugs includes both [botanicals](#) and other approved drugs that contain a mixture of natural plant-derived and synthetic or semi-synthetic substances.




Table of contents

- Summary
- Overview: Dietary supplements and botanical natural products
- Overview: Global market**
- Measurement and analysis
- Instrumentation
  - Sample preparation
  - Separation
  - Detection
- Authentication of supplements
- Application examples
  - Substances A-C
  - Substances D-G
  - Substances H-M
  - Substances N-T
  - Substances U-Z
- Literature
- Glossary

The increasing demand for dietary supplements is being driven by the growing health concerns of consumers, especially aging populations. Dietary supplements are one key segment of the nutraceuticals market, alongside functional foods and beverages. With respect to natural substances, a long history of use and a more favorable regulatory climate are key drivers in the market growth of [botanical](#) natural products.

The use of innovative technologies in production and enhanced quality monitoring by regulating bodies enable the growth of these markets and the delivery of safe drugs.



**Outer circle:** Dietary supplement market distribution by regions worth \$100.8 Billions in total (est. for 2023)  
**Inner circle:** Dietary supplement market distribution by ingredients (est. for 2023)

[From bcc Research – Nutraceuticals: Global Markets to 2023. Report Code FOD013G](#)

Estimated growth of plant-derived and [botanical](#) drugs markets. Compound annual growth rate (CAGR) shown for both.

[From bcc Research – Botanical and Plant-derived Drugs: Global Market, Report Code BIO022H](#)



Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature


Glossary

This Applications Notebook focuses on high performance liquid chromatography (HPLC) using two different HPLC/UHPLC detectors – UV-Vis Absorbance and Charged Aerosol Detection (CAD) – for the measurement of a wide variety of compounds in different sample matrices. There are differences in the performance of these detector technologies (shown in the table below), but they are also complementary. When used serially (UV-Vis then CAD) they can markedly extend the range of analytes measured in a sample.

	UV-Vis Absorbance	CAD
Type of Detector	<a href="#">Optical - selective*</a>	<a href="#">Nebulizer/particle charging - universal</a>
What Compounds are Measured?	Analytes must have a <a href="#">chromophore</a> , but can be volatile or non-volatile	Analytes must be non- or semi-volatile but do not need to contain a <a href="#">chromophore</a>
Detector Response	<a href="#">Dependent upon chromophore structure and absorption wavelength selected*</a>	<a href="#">Uniform for all non-volatile compounds, independent of chemical structure</a>
Calibration	Requires external standard	<a href="#">Single calibrant can be used for quantitation when individual standards are not available</a>

\*Low-wavelength UV detection is used for compounds that possess a weak [chromophore](#). While this approach increases the range of analytes detected, it lacks sensitivity and results in more complex chromatograms with a greater chance for analyte coelution and interference.

# Measurement and analysis



**CAD is a Universal Detector.**






# Using CAD for quantitation when standards are not available

- Table of contents
- Summary
- Overview: Dietary supplements and botanical natural products
- Overview: Global market
- Measurement and analysis
- Instrumentation
  - Sample preparation
  - Separation
  - Detection
- Authentication of supplements
- Application examples
  - Substances A-C
  - Substances D-G
  - Substances H-M
  - Substances N-T
  - Substances U-Z
- Literature
- Glossary

A critical issue for the quantitation of analytes in many [botanical](#) natural products is the lack of authentic standards.

With UV-Vis absorbance, the analytes can be volatile or non-volatile, but the response of an analyte depends upon its [chromophore](#) structure and the absorption wavelength used for detection. One compound may show a strong response at a specific absorption wavelength, while others may respond only weakly or not at all. Using a single calibrant to quantify all analytes in a sample may, therefore, lead to erroneous quantitation.

In contrast, the CAD response for non-volatiles is independent of the analyte's chemical structure. Consequently, a single calibrant can be used for quantitation of all non-volatile analytes in the sample without the need for multiple authentic standards. For semi-volatile analytes smart-chosen mobile phase additives could support the single calibrant quantification in many cases.



**CAD shows uniform response.**



Table of contents

Summary

Overview: Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T


Substances U-Z

Literature


Glossary

CAD response is, however, dependent upon the organic content of the mobile phase. This means the uniformity of response will be affected when gradients are employed. In order to use single calibrant quantitation, changes in the detector response during gradient elution must be minimized. This is readily achieved using a secondary make up flow post column. Here the inverse gradient ensures that the detector “sees” a constant mobile phase composition, thereby ensuring uniform response. Such a set-up is easily implemented using the [Thermo Scientific™ Vanquish™ Duo UHPLC System](#) for Inverse Gradient LC Workflows.

CAD response for ionizable compounds may be affected (to some extent) by the composition of the mobile phase, as a result of [salt formation](#). In general, this is a minor consideration when measuring the analyte composition of many dietary supplements and [botanical](#) natural products.



**Response of the CAD is dependent upon mobile phase composition during a gradient (red bars). An inverse gradient overcomes the effect of the gradient (gray bars) on detector response.**



**An inverse gradient delivered by the left pump ensures that the composition of the mobile phase entering the detector remains constant thereby maintaining response uniformity**

TN72806: Charged aerosol detection – factors affecting uniform response

TN73449: Why use Charged Aerosol Detection with Inverse Gradient?

SP73026-EN 0819C: Achieving standard free quantitation: Thermo Scientific Charged Aerosol Detectors





Table of contents

Summary

Overview: Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

# Collective power of chromatography

Instruments configured to meet your exact measurement needs


## Sample Preparation



## Separation



## Detection



## Table of contents

Summary

Overview:  
Dietary supplements and  
botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of  
supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

## Sample preparation

Sample preparation is a crucial step in any analytical workflow and the wide variety of sample matrices seen across dietary supplements and [botanical](#) natural products can make it especially challenging. However, modern sample preparation systems are designed to improve the speed, efficiency and reproducibility of this step, making it easy to extract even the most complex samples with confidence.

#### Automate sample preparation with an Accelerated Solvent Extractor system (ASE)

- Automates sample preparation for solid and semi-solid samples, using solvents at elevated temperatures.
- Operates above the boiling point of extraction solvents by using elevated pressure, thereby permitting fast liquid extractions at high temperature.
- Walk-away system that performs extraction and clean up on 24 samples unattended.
- Well-established and proven technique that is superior to Soxhlet and is approved for numerous regulatory methods.
- Ideal for the extraction of complex plant samples and dietary supplements prior to HPLC analysis.



[Thermo Scientific™ Dionex™ ASE™ 150 and ASE™ 350 Accelerated Solvent Extractor](#)

For more information about chromatography sample preparation and consumables click [here](#)



Table of contents


- Summary
- Overview: Dietary supplements and botanical natural products
- Overview: Global market
- Measurement and analysis
- Instrumentation
  - Sample preparation
  - Separation
  - Detection
- Authentication of supplements
- Application examples
  - Substances A-C
  - Substances D-G
  - Substances H-M
  - Substances N-T
  - Substances U-Z
- Literature
- Glossary

**Choosing between general purpose and specialty HPLC columns**

The choice of column used at the separation stage can significantly impact the quality of results obtained, so is an important element to consider.

Thermo Fisher Scientific offers columns with a wide range of bonded phases, particle sizes, particle morphologies and column dimensions. These comprehensive options are designed to meet any application need to deliver maximum analyte resolution. Some things to consider when choosing the most appropriate LC column for your requirements are:

- Fully porous particles, found in [Thermo Scientific™ Hypersil™ GOLD columns](#) and [Thermo Scientific™ Acclaim™ columns](#), allow for large injection volumes and high resolution of complex sample matrices.
- Superficially porous (solid core) particles, found in [Thermo Scientific™ Accucore™ columns](#), provide greater signal-to-noise characteristics with smaller injection volumes and/or when performing UHPLC separations without elevated system backpressure.
- A column's hydrophobic properties govern the separation of most analytes. While C18 is the most common phase selected in method development, consider hydrophobic retention trends with the carbon load on the column: C1 < C4 < C8 < C18 / C30.



Analyte	Suggested column(s)	Benefit
Lipids	<ul style="list-style-type: none"> <li>• <a href="#">Thermo Scientific™ Acclaim™ C30 column</a></li> <li>• <a href="#">Thermo Scientific™ Accucore™ C30 column</a></li> </ul>	Greater selectivity with shorter run times compared to a traditional C18
Large molecular weight (>4kDa)	<ul style="list-style-type: none"> <li>• <a href="#">Thermo Scientific™ Acclaim 300 C18 column</a></li> <li>• <a href="#">Thermo Scientific™ Accucore 150 C18 column</a></li> </ul>	Wider pore silica column minimizes sample carryover
Moderately polar analytes (including basic compounds)	<ul style="list-style-type: none"> <li>• <a href="#">Thermo Scientific™ Acclaim™ PA2 column</a></li> <li>• <a href="#">Thermo Scientific™ Accucore™ Polar Premium column</a></li> </ul>	Achieves greater retention by selecting a polar embedded column
Aromatic compounds	<ul style="list-style-type: none"> <li>• <a href="#">Thermo Scientific™ Accucore™ Biphenyl column</a></li> <li>• <a href="#">Thermo Scientific™ Hypersil GOLD™ PFP column</a></li> </ul>	Enhances steric selectivity and resolution with phenyl-based columns
Neutral and charged analytes	<ul style="list-style-type: none"> <li>• <a href="#">Thermo Scientific™ Acclaim™ Trinity P1 column</a></li> <li>• <a href="#">Thermo Scientific™ Acclaim™ Trinity P2 column</a></li> </ul>	Provides controlled ion exchange and reversed-phase (RP)/hydrophilic interaction liquid chromatography ( HILIC) properties - the Trinity line works well for samples with a mix of ion exchange and hydrophobic/hydrophilic properties

For details of the complete range of HPLC/UHPLC columns click [here](#)



- Table of contents
- Summary
- Overview: Dietary supplements and botanical natural products
- Overview: Global market
- Measurement and analysis
- Instrumentation
  - Sample preparation
  - Separation**
  - Detection
- Authentication of supplements
- Application examples
  - Substances A-C
  - Substances D-G
  - Substances H-M
  - Substances N-T
  - Substances U-Z
- Literature
- Glossary

The complexity and variety of potential analytes makes the analysis of dietary supplements and [botanical](#) natural products especially challenging.

Thermo Scientific Vanquish HPLC and UHPLC systems can solve even the toughest analytical challenges—they are dependable and consistently deliver high caliber results you can depend on.



**Thermo Scientific™ Vanquish™ Core HPLC system**  
 Absolute dependability to enable worry-free applications

**Thermo Scientific™ Vanquish™ Flex UHPLC system**  
 Complete flexibility for method development or fast and reliable UHPLC

**Thermo Scientific™ Vanquish™ Horizon UHPLC system**  
 Unrivalled performance and throughput for applications requiring high-end UHPLC

**Thermo Scientific™ Vanquish™ Duo HPLC and UHPLC systems**  
 increased productivity with two completely independent flow paths for higher throughput and improved sample characterization

[Find out more about Thermo Scientific Vanquish HPLC and UHPLC systems](#)



- Table of contents
- Summary
- Overview: Dietary supplements and botanical natural products
- Overview: Global market
- Measurement and analysis
- Instrumentation
  - Sample preparation
  - Separation
  - Detection**
- Authentication of supplements
- Application examples
  - Substances A-C
  - Substances D-G
  - Substances H-M
  - Substances N-T
  - Substances U-Z
- Literature
- Glossary

Choosing the most appropriate detection technology for your application is crucial in revealing all the components of interest in your sample. The Vanquish HPLC and UHPLC platforms offer a wide range of detection capabilities that can be easily integrated and combined to fit your analytical needs.

Although the focus here is on UV-Vis Absorbance and Charged Aerosol Detectors, a complete range of detectors is available for the analysis of dietary supplements and [botanical](#) natural products, including [mass spectrometry](#) solutions. These enable you to choose and configure the exact set-up for your specific requirements.

**Industry-leading diode array detection**



[Thermo Scientific™ Vanquish™ Diode Array Detectors](#)

**Cost-effective reliable UV-Vis detection**




[Thermo Scientific™ Vanquish™ Variable Wavelength Detectors](#)

**Unmatched all-round performance of CAD**



[Thermo Scientific™ Vanquish™ Charged Aerosol Detectors](#)

**Easily integrated mass detection**



[Thermo Scientific™ ISQ™ EC Single Quadrupole Mass Spectrometer](#) or [Thermo Scientific™ ISQ™ EM Single Quadrupole Mass Spectrometer](#)



Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

**State-of-the-art instrumentation is essential for successful sample analysis using charged aerosol detection**

Thermo Scientific HPLC and UHPLC systems combined with **Thermo Scientific Charged Aerosol Detectors (CAD)**, advanced column technologies and proven analytical methods deliver the precise automation and advanced data handling that helps you to:

- Characterize many classes of compounds
- Analyze compounds in a broad range of samples types
- Profile or quantify analytes

**Thermo Scientific Vanquish Charged Aerosol Detectors and Thermo Scientific™ Corona™ Veo™ Charged Aerosol Detectors provide:**

- Simple, intuitive operation
- Wide linear and dynamic range
- Sub-nanogram sensitivity
- Method flexibility covering micro-flow HPLC and UHPLC applications with a single nebulizer
- Adjustable evaporation temperature to optimize signal-to-noise ratio



**Vanquish Charged Aerosol Detector**



**Corona Veo Charged Aerosol Detector**



- Table of contents
- Summary
- Overview: Dietary supplements and botanical natural products
- Overview: Global market
- Measurement and analysis
- Instrumentation
  - Sample preparation
  - Separation
  - Detection**
- Authentication of supplements
- Application examples
  - Substances A-C
  - Substances D-G
  - Substances H-M
  - Substances N-T
  - Substances U-Z
- Literature
- Glossary



**Variable wavelength detectors offer sensitive, accurate UV-Vis detection**

Ensure high sensitivity and accuracy during UV-Vis detection of your analytes with Thermo Scientific variable wavelength detectors (VWDs) for HPLC. The optical bench design of this detector type minimizes baseline drift and noise.

**Diode array detectors combine flexibility and performance in UV-Vis detection**

Thermo Scientific Vanquish Diode Array Detectors (DADs), supporting 3D-field detection, and Multiple Wavelength Detectors (MWDs) are designed for the highest performance and flexibility. A wide spectral range and several programmable detector optimization parameters meet challenging method development requirements and application demands.

Thermo Scientific UV-Vis detectors offer several advantages, including excellent noise, drift, and linearity performance to support a wide operating range and trace compound detection. These absorbance detectors are easy to operate and streamline sample analysis by covering both the UV and Vis ranges of light.

Product Specification Sheet: [Vanquish Variable Wavelength Detectors](#)  
Product Specification Sheet: [Vanquish Diode Array Detectors and Multiple Wavelength Detectors](#)



## Table of contents

Summary

Overview:  
Dietary supplements and  
botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of  
supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T


Substances U-Z

Literature

Glossary

Safety, authenticity and sustainability are key issues when using plant-derived materials. Unintentional contamination or economically-motivated deliberate adulteration within the supply chain continue to be major problems affecting commercially available products. The need for authentication applies whether the sample is the source plant or a plant material such as a root, shoot or leaf.

As well as the analytical challenges already highlighted, authentication may be complicated by the use of poorly characterized material in previous studies. To address this, HPLC-based methods are used to evaluate [botanical metabolomics](#) approaches. Presented here is an untargeted HPLC metabolomic method, which uses [principal component](#) analysis to distinguish between authentic samples and an adulterant, based on differences in their metabolite patterns.



**Principal component analysis score plot of hoodia samples showing the ability to differentiate between authentic dried plant, capsules and tablets, and an adulterant *Opuntia leptocaulis* samples. See [Poster Note 70540](#) for greater detail.**

PN70540: Profiling Hoodia Extracts by HPLC with Charged Aerosol Detection, Electrochemical Array Detection, and Principal Component Analysis

For more background on the use of CAD to determine product authentication and adulteration see:

AN7317: Determination of olive oil purity based on triacylglycerols profiling by UHPLC-CAD and Principal Component Analysis





In the following pages, the analysis of a selection of dietary supplements and [botanical](#) natural products is described.

Table of contents
Summary
Overview: Dietary supplements and botanical natural products
Overview: Global market
Measurement and analysis
Instrumentation
Sample preparation
Separation
Detection
Authentication of supplements
Application examples
Substances A-C
Substances D-G
Substances H-M
Substances N-T
Substances U-Z
Literature
Glossary



## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

Acetylcarnitine is produced naturally by the body and is also taken as a dietary supplement. It is broken down in the blood by plasma esterases to carnitine, which functions as part of the system that transports fatty acids into the mitochondria for energy metabolism.

Lipoic acid is an essential cofactor for several enzyme complexes involved in aerobic metabolism. As well as occurring naturally in many foods, lipoic acid is also formed in the body and may additionally be taken as a supplement.


Both acetylcarnitine and lipoic acid are highly polar compounds and they are not well-retained on RP-HPLC columns. Presented here is an HPLC-UV method using an Acclaim Trinity P1 column for the simultaneous measurement of acetylcarnitine and lipoic acid in a dietary supplement. The cation-exchange, anion-exchange and reversed-phase retention mechanisms of the Acclaim Trinity P1 make it an ideal column for this application.

Column: Thermo Scientific Acclaim Trinity P1, 3  $\mu$ m, 3.0 x 50 mm


Mobile phase A: 3.00 g Monobasic sodium phosphate (25 mmol), 22 mg Tetrasodium pyrophosphate decahydrate (0.5 mmol) + 270  $\mu$ L 85% Phosphoric acid (4 mmol) + 196 g Acetonitrile + 750 g DI Water

Mobile Phase B: 196 g Acetonitrile + 750 g DI Water


Detector: Absorbance, UV 210 nm



Acetylcarnitine



Lipoic acid

**HPLC-UV analysis of dietary supplement.**1 – Acetylcarnitine;  
2 – Lipoic acid.

## Table of contents

Summary

Overview:  
Dietary supplements and  
botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of  
supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

The genus *Aristolochia* includes several plants known as birthwort or Dutchman's pipe. One species, *Aristolochia clematitis*, was used as a medicinal plant by the ancient Egyptians, Greeks and Romans, to treat asthma, hiccups, spasms, pains, and expulsion of afterbirth. In [traditional Chinese medicine](#) *Aristolochia* species are used to treat arthritis and edema. *Aristolochia*-containing supplements are also promoted to help with weight loss.

*Aristolochia* species are now considered to be toxic. There are reports of people suffering nephritis and rapid kidney failure after consuming certain weight loss supplements from China. The toxins responsible are most likely the aristolochic acids (AA I and AA II). These are potent carcinogens that can cause liver and urothelial cancer. Consequently, the [FDA has issued warnings](#) against consumption of AA-containing supplements.


Presented here is an HPLC-CAD method for the determination of AA1 and AA2 in supplements.

HPLC Column: Thermo Scientific Accucore C18,  
2.6  $\mu\text{m}$ , 2.1 x 150 mm


Mobile Phase A: 0.1 % Formic acid in 10 mM  
Ammonium formate

Mobile Phase B: 0.1 % Formic acid in Acetonitrile


Detector: Charged Aerosol



**HPLC-CAD analysis of a methanolic extract of *Aristolochia fangchi* root (black) and standards (red).**



**Aristolochic acid I (AAI)**



**Aristolochic acid II (AAII)**

Further information on [Thermo Scientific AppsLab Library of Analytical Applications](#)



## Table of contents

## Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

## Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

Sweet wormwood (*Artemisia annua*), a plant used in [traditional Chinese medicine](#), is the source of the drug artemisinin. Artemisinin and its semisynthetic derivatives are used in the treatment of malaria and parasitic worm infections.


Artemisinin does not contain a [chromophore](#). Low-wavelength UV absorbance detection lacks sensitivity and is prone to chromatographic interferences. Artemisinin is an ideal candidate for analysis by HPLC-CAD.

Presented here is an HPLC-CAD method that can be used for the analysis of artemisinin and related compounds.


HPLC Column: C18, 3  $\mu$ m, 4.6 x 75 mm

Mobile Phase: 60% Acetonitrile, pH 3 with Trifluoroacetic acid


Detector: Charged Aerosol




**Analysis of artemisinin and related compounds (3.3  $\mu$ g each on column).**



**Artemisinic acid**



**Artemisinin**



**Arteannuin B**

For HPLC-UV methods see:

[Separation of artemether and its impurities using reversed-phase HPLC-UV](#)

[Rapid analysis of artesunate and dihydroartemisinin using a Thermo Scientific Accucore RP-MS HPLC column](#)

[Fast analysis of artesunate and dihydroartemisinin using a Thermo Scientific Synchronis C18 HPLC column](#)



## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary


*Withania somnifera*, a member of the nightshade family, is commonly called Indian winter cherry, poison gooseberry or Indian ginseng. Known as ashwagandha, it is used in [Ayurvedic](#) medicine to treat stress, improve concentration and increase energy.

Ashwagandha is taken as a strengthening tonic, sexual tonic, memory enhancer, stress reliever and to cool the body. Purported active compounds include the steroidal lactones (e.g., withanolides and withaferins), [alkaloids](#) (anaferine, anahygrine, cuseohygrine, and isopelletierine) and [saponins](#). The withanolides are structurally similar to the ginsenosides of *Panax ginseng*.

Ashwagandha extracts can be analyzed by HPLC-UV detection, but this approach is limited as many compounds possess only weak [chromophores](#). HPLC-CAD provides a more comprehensive analysis.

As shown here, HPLC-UV absorbance detection fails to measure many compounds found in an ashwagandha extract.

Column: Thermo Scientific Accucore C8, 2.6  $\mu\text{m}$ , 4.6 x 150 mm  
 Mobile Phase A: DI Water  
 Mobile Phase B: Acetonitrile, Optima LCMS  
 Detector: Charged Aerosol and Absorbance, UV 230 nm



**Comparison between CAD and UV absorbance detection for HPLC separation of ashwagandha extract. The shaded area shows that many components are missed by the HPLC-UV approach.**

## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature


Glossary

*Bacopa monnieri*, commonly called brahmi, water hyssop, herb of grace or Indian pennywort, is used in [Ayurvedic](#) medicine to enhance mind power (the [Medhya effect](#)) and improve all aspects of mental functioning, including memory and comprehension. Among the purported bioactive compounds found in *Bacopa monnieri* are dammarane-type [triterpenoid saponins](#), known as bacosides.


Although extracts can be analyzed by HPLC-UV detection, this approach is limited because many compounds possess weak [chromophores](#). HPLC-CAD enables a more comprehensive and sensitive analysis.

Presented here is a comparison of HPLC with either CAD or UV detection of a *Bacopa monnieri* extract.

Column: Thermo Scientific Accucore C18, 2.6  $\mu\text{m}$ , 2.1 x 150 mm  
 Mobile Phase A: DI Water  
 Mobile Phase B: Acetonitrile, Optima LCMS  
 Detector: Charged Aerosol and Absorbance, UV 210 nm



**Comparison of HPLC with either CAD or UV detection of a *Bacopa monnieri* extract. Note UV detection misses early eluting analytes.**





## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature


Glossary

*Actaea racemosa* (synonym *Cimicifuga racemosa*), also called black cohosh, black bugbane, black snakeroot or fairy candle, is a species of flowering plant native to eastern North America. Black cohosh was used by Native Americans to treat gynecological and other disorders, while European settlers used it to treat snakebite, inflamed lungs and pain from childbirth. Today, it is taken as a supplement to address gynecological issues.


The active ingredients are believed to be the [triterpene glycosides](#), including actein, cimifugoside and related compounds. Extracts can be analyzed by HPLC-UV detection, but this approach is limited as the [triterpene glycosides](#) contain weak [chromophores](#). A more comprehensive and sensitive analysis is obtained using HPLC-CAD.

Presented here is a method for HPLC-CAD analysis of a black cohosh extract.

Column: Thermo Scientific Accucore C18, 2.6  $\mu\text{m}$ , 2.1 x 150 mm  
 Mobile Phase A: DI Water  
 Mobile Phase B: Acetonitrile, Optima LCMS  
 Detector: Charged Aerosol



HPLC-CAD analysis of black cohosh extract.



CAN113: Determination of Triterpene Glycosides in *Cimicifuga racemosa* (Black Cohosh) by HPLC-CAD

PN70543: Novel, Universal Approach for the Measurement of Natural Products in a Variety of Botanicals and Supplement, Part 2

PN70153: Can High Peak Capacity and Universal Detection Solve the Challenges in LC Characterization of Botanicals and Natural Products





## Table of contents

Summary

Overview:  
Dietary supplements and  
botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of  
supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature



Glossary

*Boswellia* is a genus of tree known for its fragrant resin. A well-known extract is frankincense, which comes from the resin of *Boswellia sacra*. More generally, *Boswellia* extracts have been used for centuries in [Asian and African folk medicine](#) to treat chronic inflammatory illnesses and various other health conditions.


The active ingredients are believed to be the pentacyclic [triterpene](#) boswellic acids. These compounds contain weak [chromophores](#) so are often measured using HPLC with low-wavelength UV detection.

Presented here is method for analyzing boswellic acid and related compounds in a *Boswellia* extract using HPLC-UV.

Column: C18, 5  $\mu$ m, 4.6  $\times$  150 mm  
 Mobile phase A: Acetonitrile/ DI Water/ Phosphoric acid (85%)  
 80:20:0.1 v/v/v  
 Mobile phase B: Acetonitrile  
 Detector: Absorbance, UV 210 or 250 nm


**Determination of boswellic acids in boswellia extracts, UV at 210 nm.**

1. 11-keto- $\beta$ -Boswellic acid
2. 3-acetyl-11-keto- $\beta$ -Boswellic acid
3.  $\alpha$ -Boswellic acid
4.  $\beta$ -Boswellic acid
5. 3-acetyl- $\alpha$ -Boswellic acid
6. 3-acetyl- $\beta$ -Boswellic acid


**Determination of boswellic acids in boswellia extracts, UV at 250 nm.**

1. 11-keto- $\alpha$ -Boswellic acid
2. 3-acetyl-11-keto- $\alpha$ -Boswellic acid







## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

Calcium and magnesium are essential nutrients and are commonly taken as dietary supplements. The type of calcium or magnesium salt present can affect the rate of absorption, so a variety of counterions may be included in the formulation. The analysis of mineral supplements provides an interesting analytical challenge in that both the anions and cations are functional ingredients that require determination. In the past, separate assays for anions and cations were used. However, advances in multi-mode chromatography now permit resolution of both anions and cations in one run using a single HPLC column.


Presented here is an HPLC-CAD method using an [Acclaim Trinity P2 column](#) for separation and detection of various cations and anions in a mineral supplement.

Column: Thermo Scientific Acclaim Trinity P2,  
3  $\mu$ m, 3 x 50 mm

Mobile phase A: DI Water

Mobile phase B: 100 mM Ammonium formate, pH 3.65

Detector: Charged Aerosol

**Analysis of a mineral supplement.**

1. Aspartate
2. Citrate
3. Unknown
4. Magnesium
5. Calcium

AppsLab 656: Simple gradient method for the analysis of calcium and magnesium in a dietary mineral supplement using HPLC-CAD  
AN20871: Separation of Calcium, Magnesium and Counterions in a Dietary Supplement Using Multi-mode Liquid Chromatography with Charged Aerosol Detection



## Table of contents

Summary

Overview:  
Dietary supplements and  
botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of  
supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

The edible cactus *Caralluma fimbriata* has been used in India for thousands of years as part of the diet and as an appetite suppressant, typically during long hunts and in times of famine. Today, *Caralluma* supplements are taken to suppress hunger and to enhance endurance.

The purported active components contained in *Caralluma* supplements are the [oxypregnane glycosides](#). These compounds contain weak [chromophores](#) so are often measured by HPLC with low-wavelength UV detection. A more comprehensive and sensitive analysis is obtained using HPLC-CAD.


Presented here is a method for the analysis of [oxypregnane glycosides](#) in *Caralluma* using HPLC with both CAD and UV detection.

Column: Thermo Scientific Accucore Vanquish C8, 2.6  $\mu\text{m}$ , 4.6  $\times$  150 mm


Mobile phase A: DI Water

Mobile phase B: Acetonitrile, Optima LCMS

Detector: Charged Aerosol and Absorbance, UV 210 nm



**Analysis of [oxypregnane glycosides](#) in a *Caralluma* supplement extract by HPLC with both CAD and UV detection. Note that many analytes are not detected by HPLC-UV detection and that this approach suffers from baseline perturbation throughout the gradient.**



**Pregnane [glycoside](#) from *Caralluma fimbriata***



## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of  
supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

Carnitine functions as part of the system that transports fatty acids into the mitochondria for energy metabolism. Carnitine is produced by the liver, but hepatic synthesis may not always meet the body's needs, so it is sometimes taken as a dietary supplement. Carnitine is highly polar and not well retained on RP-HPLC columns. The Acclaim Trinity P1 column provides cation-exchange, anion-exchange and reversed-phase retention mechanisms, and is found to be ideal for carnitine analysis.


Presented here is an HPLC-UV method using an Acclaim Trinity P1 column for the measurement of carnitine tartrate in a dietary supplement.

Column: Thermo Scientific Acclaim Trinity P1, 3  $\mu\text{m}$ , 3.0 x 50 mm


Mobile phase A: 3.00 g Monobasic sodium phosphate (25 mmol), 22 mg Tetrasodium pyrophosphate decahydrate (0.5 mmol) + 270  $\mu\text{L}$  85% Phosphoric acid (4 mmol) + 196 g Acetonitrile + 750 g DI Water

Mobile Phase B: 196 g Acetonitrile + 750 g DI Water


Detector: Absorbance, UV 210 nm



Carnitine



Tartaric acid

**HPLC-UV analysis of dietary supplement.**

1. Carnitine
2. Tartrate



## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary








*Lonicera japonica*, known as Japanese honeysuckle and golden-and-silver honeysuckle, is a twining vine native to eastern Asia. In [traditional Chinese medicine](#) *Lonicera japonica* is called ren dong téng or jin yin hua. The dried [rattan](#) (caulis ionicerae) leaves and flowers ([flos ionicerae](#)) are used to treat fever, headache, cough, thirst, sore throat and epidermal diseases.

Presented here is an HPLC-UV method that can simultaneously measure the six main active components (caffeic acid, chlorogenic acid, galuteolin, loganin, rutin and sweroside) of caulis ionicerae, in under 10 minutes. It is a viable alternative to the two HPLC methods required by The [Pharmacopeia](#) of the People's Republic of China for the regulation of this herbal material.

Column: Thermo Scientific Acclaim Phenyl-1, 3  $\mu$ m, 4.6 x 150 mm

Mobile phase: Acetonitrile, 0.4% Formic acid

Detector: Absorbance, UV 236 nm



**Analysis of a A) caulis ionicerae sample and B) mixture of standards by HPLC-UV.**

1. Loganin
2. Sweroside
3. Chlorogenic acid
4. Caffeic acid
5. Rutin
6. Galuteolin



## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature



Glossary

Choline is an essential nutrient involved in cell membrane integrity, cell membrane signaling, lipid transport and metabolism, neurotransmission, fetal brain development, and modulation of gene expression. Although choline can be produced endogenously in the liver, this synthesis is insufficient for human needs. Choline is naturally present in some foods, especially animal-based products (meat, poultry, fish, dairy products and eggs), and in cruciferous vegetables, certain beans, nuts, seeds and whole grains. However, most people in the United States consume less than the adequate intake for choline. Choline deficiency can lead to muscle damage, liver damage and nonalcoholic fatty liver disease. Choline supplements, including choline bitartrate, choline chloride, choline citrate, phosphatidylcholine and lecithin, are used to address this problem.

Quality testing is needed to ensure the safety and effectiveness of choline dietary supplements. Choline and its salts possess weak [chromophores](#), and so are not suitable candidates for detection by UV absorbance.

Presented here is a mixed mode (zwitterionic operated in HILIC mode) chromatographic method using CAD to measure choline and its salts. The method is compatible with LC-MS, which was used to identify the impurity O-(2-hydroxyethyl)choline, in choline supplements.

Column: Zwitterionic HILIC  
 Mobile phase A: Ammonium acetate  
 Mobile Phase B: Acetonitrile  
 Detector: Charged Aerosol



**HPLC-CAD measurement of sodium, chloride and sample impurities in a 2 mg/mL choline citrate sample. The sample impurity (RT~6 min) was studied using a flow split to Thermo Scientific Vanquish ISQ EM Single Quadrupole Mass Spectrometer.**

## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of  
supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature


Glossary

Cyclotides are bioactive peptides that typically contain 28-37 amino acid residues, forming a cyclotide peptide with a structure that is locked in place by three disulfide bonds. More than 100 cyclotides have so far been identified in plants. Cyclotides contain a conserved core of amino acids and a series of hypervariable loops. This suggests they may play an important role in nature, for example protecting the plant from predators and pathogens. Cyclotides are reported to have a wide range of biological activities including insecticidal, antimicrobial, antitumor, anti-HIV, protease inhibitory, hormone-like, cytotoxic and uterotonic properties. Recent research suggests that cyclotides are also novel candidates as scaffolds for peptide drug delivery.


Cyclotides possess strong [chromophores](#) so are usually measured using HPLC-UV detection.

Presented here is an HPLC-UV method for the measurement of the cyclotide cycloviolacin O2 in a plant extract.

Column:	Thermo Scientific Acclaim C18, 2.2 $\mu$ m, 3 x 75 mm
Mobile Phase A:	50 mM Ammonium formate, pH 4.4,
Mobile Phase B:	65% acetonitrile, 15% DI Water, 20% 100mM Ammonium formate pH 4.4
Detectors:	Absorbance, UV 280 nm



**HPLC-UV analysis of A) cycloviolacin O2 standard (100  $\mu$ g/mL) and B) a *Viola odorata* plant extract**



## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of  
supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z


Literature

Glossary





*Echinacea* is a group of herbaceous flowering plants in the daisy family, commonly called the coneflower. Three species of *Echinacea* are used as herbal remedies, including *Echinacea angustifolia*, *Echinacea pallida* and *Echinacea purpurea*. *Echinacea angustifolia* was widely used by the North American Plains Indians for its supposed medicinal qualities. Today, as an [herbal supplement](#), echinacea is consumed in teas, liquid extracts, as a dried herb, or capsules or pills to reduce many of the symptoms of the common cold, flu and some other illnesses, infections and conditions. However, [studies to date](#) have not reported any benefit of echinacea on common cold prevention or duration.

Presented here is an HPLC-UV absorbance method for measuring key phenolic compounds, derivatives of caffeic acid, in *echinacea* extracts.

Column: C18, 3  $\mu\text{m}$ , 4.6 x 150 mm  
 Mobile Phase A: Acetonitrile/DI Water (90:10 v/v)  
 Mobile Phase B: Acetonitrile/DI Water/Phosphoric acid (27:75:0.1 v/v/v)  
 Column Temp.: 25°C  
 Flow rate: 1.5 mL/min  
 Inj. volume: 25  $\mu\text{L}$   
 Gradient: 100 % A to 100 % B 20 min; hold 100% B 2 min  
 Detector: Absorbance, UV 330 nm

**Determination of phenolic compounds in Echinacea extracts.**

1. Cattaric acid
2. Chlorogenic acid
3. Echinacoside
4. Cichonic acid

**Cichoric acid****Chlorogenic acid****Echinacoside****Cattaric acid**



## Falcarinol, falcarindiol and falcarindiol-3-acetate

### Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature


Glossary

Some plants in the *Apiaceae* (or *Umbelliferae*) family (e.g., carrots, parsnips, parsley and celery) produce a group of bioactive C17-polyacetylene compounds, the polyacetylenic oxylipins, in response to pathogens. Three such compounds - falcarinol, falcarindiol and falcarindiol-3-acetate - are natural pesticides and highly toxic towards bacteria and fungi. They also exhibit a diverse range of biological activities in mammals, both beneficial (e.g., their cytotoxicity is proposed to reduce the risk of developing cancer) and detrimental (e.g., occupational allergic contact dermatitis).

Falcarinol and related compounds contain weak [chromophores](#) so are usually measured using low-wavelength UV absorbance (205 nm), but sample chromatograms tend to be complex due to the presence of many other compounds absorbing at this wavelength. HPLC-CAD overcomes limitations of HPLC-UV and is much more sensitive.



Presented here is an HPLC-CAD method for the analysis falcarinol and related compounds showing results from parsnip skin and parsnip core extracts.

Column: Thermo Scientific Acclaim RSLC 120 C18, 2.2  $\mu\text{m}$ , 2.1 x 150 mm  
 Mobile phase A: Methanol: DI Water: Acetic acid (500 : 500 : 4)  
 Mobile phase B: Acetone: Methanol: Tetrahydrofuran: Acetic acid (500 : 375 : 125 : 4)  
 Detector: Charged Aerosol



**Falcarinol**

PN-2023: Simple and Direct Analysis of Falcarinol and other Polyacetylenic Oxylipins in Carrots by Reverse Phase HPLC and Charged Aerosol Detection



**HPLC-CAD analysis of A) parsnip skin extract and B) parsnip core extract.**

1. Falcarindiol
2. Falcarindiol-3-acetate
3. Falcarinol







## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

Anthocyanins and their [aglycones](#), the anthocyanidins, are flavonoid plant pigments. Anthocyanins are particularly abundant in both the European bilberry (*Vaccinium myrtillus* L.) and North American blueberry (*Vaccinium corymbosum*). However, fresh bilberries contain four times more anthocyanin than fresh blueberries. Bilberry extracts are widely used in nutritional supplements and pharmaceuticals for improving visual acuity and treating circulatory disorders. [Chemical and pharmacological studies have identified anthocyanins as the main components responsible for the therapeutic effect of the extracts that are used in these supplements.](#)

Presented here is HPLC method using absorbance detection for the analysis of anthocyanins and anthocyanidins.


HPLC Column: Thermo Scientific Acclaim 120 C18, 2.2  $\mu\text{m}$ , 2.1 x 150 mm

Mobile Phase A: 10 % Formic acid

Mobile Phase B: 10 % Formic acid, 22.5% Methanol, 22.5% Acetonitrile


Detector: Absorbance, Vis 520 nm

Anthocyanidin	R3	R5	R6	R7	R3'	R4'	R5'
Pelargonidin	OH	OH	H	OH	H	OH	H
Cyanidin	OH	OH	H	OH	OH	OH	H
Delphinidin	OH	OH	H	OH	OH	OH	OH
Petunidin	OH	OH	H	OH	OCH <sub>3</sub>	OH	OH
Peonidin	OH	OH	H	OH	OCH <sub>3</sub>	OH	OH
Malvidin	OH	OH	H	OH	OCH <sub>3</sub>	OH	OCH <sub>3</sub>



**Separation of a fifteen anthocyanin standard and five anthocyanidins on the solid core Acclaim RSLC C18 column.**

- |           |                |              |               |                                    |
|-----------|----------------|--------------|---------------|------------------------------------|
| 1. Dp3Gal | 6. Pet3Gal     | 11. Pet3Ara  | 16. Mal3Glu   | Abbreviations are defined in AN281 |
| 2. Dp3Glu | 7. Cy3Ara      | 12. Peo3Glu  | 17. Mal3Ara   |                                    |
| 3. Cy3Gal | 8. Delphinidin | 13. Mal3Gal  | 18. Petunidin |                                    |
| 4. Dp3Ara | 9. Pet3Glu     | 14. Peo3Ara  | 19. Peonidin  |                                    |
| 5. Cy3Glu | 10. Peo3Gal    | 15. Cyanidin | 20. Malvidin  |                                    |



AN1042: Rapid Separation of Anthocyanins in Cranberry and Bilberry Extracts Using a Core-Shell Particle Column

AN281: Rapid and Sensitive Determination of Anthocyanins in Bilberries Using UHPLC

AN264: Fast Determination of Anthocyanins in Pomegranate Juice





## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T


Substances U-Z

Literature

Glossary

HPLC columns with fully porous 2.2µm and solid core-porous shell particles are two of the newer developments in separation science. Both the fully porous Thermo Scientific Accucore C18 HPLC Column and solid core Acclaim Rapid Separation LC (RSLC) 120, C18 Analytical Column can be used for high-speed, high-resolution separations without excessive backpressure.

Presented here is a simple and rapid HPLC-UV method for the simultaneous determination of anthocyanins and anthocyanidins in natural products.



**Separation of 15 anthocyanins and 5 anthocyanidins standards using a fully porous Accucore C18 column using conditions shown to left.**


1. Dp3Gal	6. Pet3Gal	11. Pet3Ara	16. Mal3Glu	Abbreviations are defined in AN281.
2. Dp3Glu	7. Cy3Ara	12. Peo3Glu	17. Mal3Ara	
3. Cy3Gal	8. Delphinidin	13. Mal3Gal	18. Petunidin	
4. Dp3Ara	9. Pet3Glu	14. Peo3Ara	19. Peonidin	
5. Cy3Glu	10. Peo3Gal	15. Cyanidin	20. Malvidin	

HPLC Column: Thermo Scientific Accucore C18, 2.6 µm, 2.1 × 150 mm

Mobile Phase A: 10 % Formic acid

Mobile Phase B: 10 % Formic acid, 22.5% Methanol, 22.5% Acetonitrile

Detector: Absorbance, Vis 520 nm



**Separation of anthocyanins in cranberry extract using a fully porous Accucore C18 column. See AN1042 for modifications to flow rate, temperature and gradient conditions.**

1. Dp3Gal	4. Cy3Gal	7. Peo3Gal	10. Mal3Gal	Abbreviations are defined in AN281.
2. Dp3Glu	5. Cy3Glu	8. Peo3Glu		
3. Unknown	6. Cy3Ara	9. Peo3Ara		

AN1042: Rapid Separation of Anthocyanins in Cranberry and Bilberry Extracts Using a Core-Shell Particle Column

AN281: Rapid and Sensitive Determination of Anthocyanins in Bilberries Using UHPLC.

AN264: Fast Determination of Anthocyanins in Pomegranate Juice



## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z


Literature

Glossary

Catechins are [flavonoids](#) found primarily in green tea, and in smaller amounts, in grapes, black tea, chocolate and wine. Catechins are potent antioxidants *in vitro* and some suggest they provide protection against certain diseases, such as cardiovascular disease and cancer. In North America, the consumption of green tea products increased due to reported health benefits. However, commercially available teas show a high variability in catechin content, so simple and rapid methods are needed to evaluate product quality.

Presented here is a study that evaluates a Thermo Scientific Accucore C18 High-Performance LC (HPLC) column to rapidly (<6 min) determine catechins in three different types of tea.


Column: Thermo Scientific Accucore C18, 2.6  $\mu\text{m}$ , 2.1  $\times$  150 mm  
 Mobile phase A: 2.5 % aq Acetonitrile  
 Mobile phase B: 0.1 % THF in Acetonitrile  
 Detector: Absorbance, UV 280 nm




**Analysis of catechins and caffeine in tea samples.**

- |                     |                |                             |
|---------------------|----------------|-----------------------------|
| 1. Gallic acid      | 4. Catechin    | 7. Epigallocatechin gallate |
| 2. Galocatechin     | 5. Caffeine    | 8. Galocatechin gallate     |
| 3. Epigallocatechin | 6. Epicatechin | 9. Epicatechin gallate      |


**Structures of catechins from *Camellia sinensis***




**Epicatechin**




**Epicatechin gallate**




**Epigallocatechin**



**Epigallocatechin gallate**



**Galocatechin**



**Galocatechin gallate**

AB150: [Rapid separation of catechins in tea using core-shell columns](#)

AN275: [Sensitive determination of catechins in tea by HPLC](#)

AN20536: [Analysis of Catechins Using an Accucore XL C8 4  \$\mu\text{m}\$  HPLC Column](#)

AN20583: [Determination of Catechins and Phenolic Acids in Red Wine by Solid Phase Extraction and HPLC](#)





## Flavonoids - Flos chrysanthemi indicum

### Table of contents

#### Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

*Flos chrysanthemi indicum*, *Dendranthema indicum* L. (D. *indicum*), is a common medicinal plant known in China as wild chrysanthemum (ye ju hua). The Chinese [Pharmacopoeia](#) (CP) Edition 2005 regulates its use as a [traditional](#) Chinese medicine (TCM). It is used with the belief that it improves eyesight and cures fever, swelling, erysipelas (a bacterial infection of the skin), sore throat, and headache.

Unfortunately, another plant, *Dendranthema lavandulaefolium* (Fish) Mak, in the same genus as *D. indicum* and growing in the same types of environments, appears similar to *Flos chrysanthemi indicum*. It is not approved for use as a TCM and often mislabeled as *Flos chrysanthemi indicum* when the plant is harvested. HPLC methods described in CP 2005 ([see AN207](#)) are inadequate and inaccurate and are not suitable for use in quality control of *Flos chrysanthemi indicum*.


Presented here is an HPLC-UV method for fingerprinting *Flos chrysanthemi indicum* based on measurement of discriminating peaks, including chlorogenic acid and [flavonoids](#) (luteolin-7-o-glucoside, linarin, luteolin and apigenin). Accelerated Solvent Extraction (ASE<sup>®</sup>) was used for isolation of the target components from the herbal samples.

Column: Thermo Scientific Acclaim C18,  
5 µm, 4.6 x 250 mm

Mobile Phase A: Acetonitrile


Mobile Phase B: 0.1 % Acetic acid

Detectors: Absorbance, UV 326 nm




**Chromatogram of a *Flos chrysanthemi indicum* sample (purchased from Tongrentang, Beijing).**


1. Chlorogenic acid
3. Luteolin-7-O-glucoside
6. Apigenin
11. Linarin
12. Luteolin-1




**Chlorogenic acid**




**Luteolin-7-O-glucoside**



**Linarin**



**Luteolin**



**Apigenin**



## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

Hesperidin, a [flavanone glycoside](#), is the predominant [flavonoid](#) in orange peel and other citrus fruits. Hesperidin is an antioxidant suggested to enhance the action of vitamin C to lower cholesterol levels. It is reported to have pharmacological action as an anti-inflammatory, antihistaminic and antiviral agent.

The [Pharmacopeia](#) of the People's Republic of China (PPRC) 2010 recommends its extraction from fruits with a Soxhlet extraction method using ligarine and methanol. This method is both time- and solvent-consuming, requiring  $\geq 5$  hours and  $>200$  mL of ligarine and methanol for each sample. The PPRC 2010 also recommends determination with a 12-minute RP-HPLC method ([See AB142](#)).


Presented here is a more efficient and cost-effective HPLC-UV absorbance method to determine hesperidin extracted from orange peel and other citrus fruits using Accelerated Solvent Extraction (ASE®).

Column: Thermo Scientific Accucore C18,  
2.6  $\mu\text{m}$ , 4.6  $\times$  150 mm


Mobile Phase A: DI Water

Mobile Phase B: Acetonitrile, Optima LCMS

Detector: Absorbance, UV 210 nm



**HPLC-UV chromatograms of A) mobile phase; B) hesperidin standard (2  $\mu\text{g/mL}$ ); C) orange peel sample; D) lemon peel sample (50-fold dilution).**



AB142: [Rapid determination of hesperidin in orange peel using accelerated solvent extraction and UHPLC](#)





## Table of contents

## Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z


Literature

Glossary


Procyanidins are members of the proanthocyanidin class of [flavonoids](#). They are polyphenols consisting of polymerized subunits of epicatechin or catechin and are structurally highly diverse because of the many possible combinations of subunits, type of bonding and branching. Procyanidins, after lignans, are the second most common class of natural phenolic substances found in nature. They are abundant in many foods, with apples and cocoa being most prominent in the Western diet. Procyanidins are purported to have many health benefits, including anti-inflammatory, hypoglycemic, insulin activation, antioxidant, hypocholesterolemic and anti-allergic properties. Particularly important may be a connection between procyanidin consumption and the lowering of risk of cardiovascular disease. To correlate dietary intake of procyanidins with an impact on disease prevention and amelioration, there is a need to develop new, as well as improved, analytical methodologies for pharmacological studies and the standardization of foods and dietary supplements.

Presented here is an HPLC-CAD method for the determination of individual procyanidins in various sample matrices, including extracts of crab apple, cocoa and cinnamon.

Column: Thermo Scientific Acclaim 120 C18, 3  $\mu$ m, 3  $\times$  150 mm,  
 Mobile Phase A: 0.05% Formic acid, 3% Acetonitrile, 0.2% Tetrahydrofuran  
 Mobile Phase B: 0.05% Formic acid, 50% Acetonitrile, 10% Tetrahydrofuran  
 Mobile Phase C: 90% Methanol  
 Detector: Charged Aerosol



**HPLC-CAD analysis of A-type and B-type procyanidins in plant extracts.**  
 CAT – catechin; EPI – epicatechin.





## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

Giant knotweed rhizome, the dried rhizome and root of *Polygonum cuspidatum* Sieb. et Zucc. is a common medicinal plant in China. [Chinese Pharmacopeia](#) Edition 2005 regulates its use as an [herbal medicine](#). It is used to treat angiocardopathy, skin inflammations and liver diseases, reduce fever, and relieve arthritis.

Purported active components, include [anthraquinones](#) (e.g., anthraglycoside A, anthraglycoside B, emodin, physcion, rhein, and chrysophanol) and stilbenes (e.g., resveratrol and polydatin).


Presented here is a method that uses an ASE® 200 Accelerated Solvent Extractor for efficient and reproducible sample preparation, and an HPLC-DAD method for simultaneous measurement of all key analytes.

HPLC Column: Thermo Scientific Acclaim 120 C18,  
5 µm, 4.6 x 250 mm


Mobile Phase A: Acetonitrile

Mobile Phase B: 20 mM Ammonium acetate


Detector: Absorbance, Vis 254 nm



	R
Resveratrol	H
Polydatin	Glucose




	R <sub>1</sub>	R <sub>2</sub>
Anthraglycoside A	Glucose	Me
Anthraglycoside B	Glucose	H
Emodin	H	H
Physcion	H	Me




Overlay of chromatograms of A) water, B) the mixed standard, and C-F) samples.

- |                      |                 |
|----------------------|-----------------|
| 1. Polydatin         | 5. Rhein        |
| 2. Resveratrol       | 6. Emodin       |
| 3. Anthraglycoside B | 7. Chrysophenol |
| 4. Anthraglycoside A | 8. Physcion     |



Chrysophanol



Rhein

AN232: Determination of Anthraquinones and Stilbenes in Giant Knotweed Rhizome by HPLC with UV detection.

## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature


Glossary

*Ginkgo biloba* (known as ginkgo) is an ancient species and the only member of the Ginkgophyta. Its leaf is used in [traditional medicine](#), while ginkgo nuts are used as food. [Extracts of \*Ginkgo biloba\* leaf sold as a dietary supplement](#) are marketed as being beneficial for cognitive function (improving memory and enhancing concentration), and in treating a number of other health issues, including high blood pressure, peripheral arterial disease, macular degeneration, tinnitus and altitude sickness. However, there is no conclusive evidence that ginkgo is helpful for any health condition.


The purported biologically active compounds present in ginkgo extracts are the sesquiterpenoid bilobalide and numerous diterpenoid ginkgolides. These compounds contain weak [chromophores](#) so measurement by HPLC with low-wavelength UV is limited, lacking sensitivity and generating complex chromatograms.

Presented here is a method for more comprehensive and sensitive analysis using HPLC-CAD.


Column: C18, 5  $\mu\text{m}$ , 4.6 x 250 mm  
 Mobile phase A: 5 % Acetonitrile in 0.1 % Trifluoroacetic acid  
 Mobile phase B: 70 % Acetonitrile in 0.1 % Trifluoroacetic acid  
 Detector: Charged Aerosol



The use of HPLC-CAD to profile a ginkgo supplement extract.



Ginkgolide



Bilobalide

Apps Lab 2366: AN1048: Novel, Universal Approach for the Measurement of Natural Products in a Variety of Botanicals and Supplements

USP Monograph: Content of flavonol glycosides in ginkgo tablets using a C18 HPLC column







## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of  
supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

The name ginseng applies to various herbs that have been used in [traditional medicines](#) for many centuries. Asian ginseng (*Panax ginseng*), also known as red ginseng or Korean ginseng, is one of several types of true ginseng. Another is American ginseng, *Panax quinquefolius*, which is traditionally used by Native Americans. However, Siberian ginseng (*Eleutherococcus senticosus*) is not a true ginseng. Ginseng is used to improve memory, fatigue, menopause symptoms and insulin response in people with mild diabetes. Ginseng-infused teas and energy drinks are consumed daily in China and neighboring countries as a tonic for vitality. The dried root is also used in foods.

The main bioactive compounds contained in *Panax ginseng* are the ginsenoside [triterpene saponins](#), including the protopanaxatriols (Rg1, Re and Rf) and protopanaxadiols (RB1, Rc, Rb2 and Rd). *Panax quinquefolius* contains the same suite of ginsenosides except for Rf. Gradient RP-HPLC with low-wavelength UV detection is typically used to measure ginsenosides as they lack strong [chromophores](#). However, this approach typically results in strongly sloping baselines that complicate peak integration, and interferences from minor components that have stronger UV [chromophores](#) than the ginsenosides. HPLC-CAD does not suffer these drawbacks.

Presented here us an HPLC-CAD method that can be used for the routine measurement of ginsenosides in ginseng samples.


Column: Thermo Scientific Acclaim RSLC PA, 2.2  $\mu$ m, 2.1 x 100 mm  
 Mobile phase A: DI Water  
 Mobile phase B: Acetonitrile  
 Detector: Charged Aerosol

AN1048: [Novel, Universal Approach for the Measurement of Natural Products in a Variety of Botanicals and Supplements](#)

CAN112: [Determination of ginsenosides in Panax ginseng by HPLC-CAD](#)


PN70153: [Can High Peak Capacity and Universal Detection Solve the Challenges in LC Characterization of Botanicals and Natural Products](#)

AN192: [Rapid analysis of ginseng using accelerated solvent extraction and HPLC](#)



**HPLC-CAD analysis of A) ginsenoside standards and B) ginseng extract.**

1. Rg1	4. Rh1	7. Rc	10. RD	13. Rh2
2. Re	5. Rg2	8. Rb2	11. SRg3	14. PPD
3. Rf	6. Rb1	9. Rb3	12. PPT	



## Table of contents

Summary

Overview:  
Dietary supplements and  
botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of  
supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature


Glossary

*Centella asiatica* (commonly called gotu kola) is a small herbaceous annual plant native to India, Sri Lanka, other parts of Asia, Northern Australia and the Western Pacific. It is employed as a medicinal herb in [Ayurvedic](#) medicine and traditional [Chinese medicine](#) to treat a wide variety of conditions. These include its use for improving memory and blood flow, as an agent for wound-healing, an anti-epileptic and an antidepressant, and for topical application in skin conditions such as ulcers, wounds and eczema. The chemical compounds of interest in gotu kola are usually considered to be the ursane- and oleanane-type [triterpenes](#) and the [triterpene glycosides](#).


Gradient RP-HPLC with low-wavelength UV detection is typically used to measure total [triterpenes](#) as they lack strong [chromophores](#). However, this approach often results in strongly sloping baselines that complicate peak integration, and interferences from minor components that have stronger UV [chromophores](#) than the [triterpenes](#).

Presented here is an HPLC-CAD method that does not suffer these drawbacks and can be used for the routine measurement of total [triterpenes](#) in gotu kola samples.

Column: Fused-Core C18, 2.7  $\mu\text{m}$ ,  
3.0  $\times$  100 mm,  
Mobile Phase A: 0.1 % Formic acid in DI Water  
Mobile Phase B: Acetonitrile  
Detector: Charged Aerosol




**Analysis of a gotu kola extract by HPLC-CAD.**



#### Ursane-type [triterpenes](#)

Asiatic acid: R = H, R<sub>1</sub> = OH  
 Madecassic acid: R = OH, R<sub>1</sub> = OH  
 Asiatikoside: R = H; R<sub>1</sub> = O-Glu(6-1)-Glu(4-1)-Rham  
 Madecassoside: R = OH, R<sub>1</sub> = O-Glu(6-1)-Glu(4-1)-Rham



#### Oleanane-type [triterpenes](#)

Terminolic acid: R = OH, R<sub>1</sub> = OH  
 Asiatikoside B: R = OH, R<sub>1</sub> = O-Glu(6-1)-Glu(4-1)-Rham

AN1048: [Novel, Universal Approach for the Measurement of Natural Products in a Variety of Botanicals and Supplements](#)

CAN111: [Determination of Triterpenes in Centella asiatica \(Gotu Kola\) by HPLC-CAD](#)





## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of  
supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

Hoodia, scientific name *Hoodia gordonii*, is a flowering succulent that grows in the Kalahari Desert in Africa. Historically, the San Bushmen used hoodia to suppress hunger during long hunts. Although dietary supplements are used as an appetite suppressant for weight loss, published clinical evidence that hoodia aids weight loss is currently lacking.

The alleged biological effects of hoodia may be due to several [oxypregnane steroidal glycosides](#), the [hoodigosides](#), that are abundant in *Hoodia gordonii*. One, P57, may be responsible for its putative appetite-suppressant effect. Measuring hoodigoside content is essential in determining whether commercial products actually contain *Hoodia gordonii*. As these compounds contain weak [chromophores](#), measurement by HPLC with low-wavelength UV is limited, lacks sensitivity and generates complex chromatograms.


Presented here is a sensitive HPLC-CAD method capable of profiling numerous hoodigosides in supplements and plant extracts.

Column: Thermo Scientific™ Accucore™ Vanquish™ C18, 1.5 μm, 2.1 ×100 mm


Mobile phase A: DI Water

Mobile phase B: Acetonitrile, Optima LCMS

Detector: Charged Aerosol



**Analysis of hoodigosides C, E, F, L, O, S and P57 standards by UHPLC-CAD.**



**P57 (aglycone = hoodigogenin)**



## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary


Since ancient times, hops (*Humulus lupulus* L.) have been a main ingredient of beer, used not only to provide bitterness, but also for their purported human health benefits (antioxidant, anti-inflammatory, immunosuppressive and chemopreventive properties). Hop supplements are used as a natural alternative to traditional hormone therapy to relieve menopausal symptoms in older women.

The chemistry of hops during beer production is complex. Isohumulones (iso- $\alpha$ -acids), formed from humulones ( $\alpha$ -acids), are the essential bitter constituents of hop resin. The poorly water-soluble  $\alpha$ -acids are isomerized to the water-soluble iso- $\alpha$ -acids during wort-boiling. Their antimicrobial effect leads to a sterile beverage, their tensioactive character stabilizes the foam, and they have a major influence on the general flavor, smell and smoothness of beer. The three major iso- $\alpha$ -acid variants, differing in their acyl side chain, include the iso-n-humulones, iso-cohumulones and iso-adhumulones. Due to the stereochemistry of iso- $\alpha$ -acids, they all occur as cis- and trans-isomers. The lifetimes of cis- and trans-isomers significantly differ from each other. Degradation products of iso- $\alpha$ -acids sensitively impact the important beer attributes mentioned above and avoidance of less stable iso- $\alpha$ -acid variants is beneficial.


Monitoring the isomerization progress, as well as the general content of iso- $\alpha$ -acids in beer during and after the brewing process, is mandatory in order to control important beer properties.

Presented here is an application that uses online solid phase extraction (SPE) with UHPLC-UV absorbance detection, to rapidly profile all key analytes in untreated beer samples.

Column: Thermo Scientific Hypersil GOLD C18, 1.9  $\mu$ m, 2.1  $\times$  100 mm  
 Mobile Phase A: 1% aq Formic acid containing 100 mg/L EDTA  
 Mobile Phase B: Acetonitrile, Optima LCMS  
 Detector: Absorbance, UV 270 nm



**Online SPE, UHPLC-UV analysis of a German Pilsner beer sample.**



AB153: Savor the flavor – robust iso- $\alpha$ -acids assaying in beer within ten minutes

AB155: Monitor the brewing process with LC – transformation of hop  $\alpha$ -acids into beer iso- $\alpha$ -acids

AB156: The everlasting paradigm – keep beer tradition or prevent beer from a skunky off-flavor?





## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of  
supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T


Substances U-Z

Literature

Glossary



Mangosteen (*Garcinia mangostana* L) is a tropical fruit indigenous to Southeast Asia. It is used broadly in [Ayurvedic](#) medicine to treat abdominal pain, cholera, diarrhea, dysentery, inflammation, wound infection, suppuration and chronic ulcers. Recently, mangosteen pericarp was proposed as an [adjunctive therapy for bipolar disorder and schizophrenia](#). Such therapeutic benefits have been mostly attributed to a unique family of compounds, the [xanthenes](#), that are most abundant in the pericarp of the fruit.

Presented here is a UHPLC-CAD method capable of analyzing several [xanthenes](#) in mangosteen pericarp. Conventional extraction methods, such as Soxhlet, are time consuming, labor-intensive and often lack reproducibility. To overcome these issues [Accelerated Solvent Extraction](#) (ASE<sup>®</sup>), an automated extraction technique that rapidly performs solvent extractions using high temperature and pressure, was used.



UHPLC-CAD analysis of mangosteen xanthenes. A) Standards and B) ASE extraction of mangosteen pericarp powder.

Column: Thermo Scientific Acclaim 120 C18, 2.2  $\mu\text{m}$ , 2.1  $\times$ 100 mm,  
Mobile Phase A: DI Water  
Mobile Phase B: Acetonitrile  
Detector: Charged Aerosol

 $\alpha$ -Mangostin

3-Isomangostin

PN70991: Fast Analysis of Selected Xanthenes in Mangosteen Pericarp Using Accelerated Solvent Extraction and Ultra High Performance Liquid Chromatography.

AB172: The Vanquish Platform: Major Improvement in Throughput and Resolution of Xanthenes in Mangosteen Pericarp



## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z


Literature

Glossary

The milk thistle plant, *Silybum marianum*, is native to Mediterranean Europe and has long been used as an herbal remedy to promote liver health and treat liver disorders, such as hepatitis and cirrhosis, and gallbladder problems. Recent research suggests that it may be beneficial for diabetes and dyspepsia.

The purported active ingredient in milk thistle is silymarin, which is particularly abundant in the plant's seeds. Silymarin consists of a mixture of different chemical species including a [flavonoid](#) (taxifolin) and several flavonolignans (silybin A, silybin B, isosilybin A, isosilybin B, silychristin, isosilychristin and silydianin).

Presented here is a high-resolution UHPLC-CAD method that measures not only the major silymarin components, but also numerous minor constituents that may be missed by typical HPLC-UV approaches.




High resolution analysis of milk thistle supplement extract by UHPLC-CAD.

Column: Thermo Scientific Acclaim Vanquish C18, 2.2  $\mu\text{m}$ , 2.1  $\times$ 250 mm


Mobile phase A: 0.1 % Formic acid in DI Water

Mobile phase B: 0.1 % Formic acid in Methanol


Detector: Charged Aerosol



Silybin



Isosilybin



Silychristin

AN1048: [Novel, universal approach for the measurement of natural products in a variety of botanicals and supplements](#)

PN70153: [Can High Peak Capacity and Universal Detection Solve the Challenges in LC Characterization of Botanicals and Natural Products](#)





## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary


Luo han kuo fruit (*Siraitia grosvenorii* Swingle), also called monk fruit, has long been used in [traditional Asian medicine](#) to treat coughs, sore throats and fatigue. Cucurbitane-type and other [triterpenes](#) isolated from the fruit were investigated for numerous potential health benefits, such as possible anti-cancer and antihyperglycemic effects. Many of these compounds are intensely sweet and, therefore, are also used as sugar substitutes and flavor enhancers. Extracts of Luo han kuo fruit used as sweeteners are acknowledged as Generally Recognized as Safe (GRAS) based on a [GRAS submission](#) to the FDA in January of 2010.

Presented here is a method that separates analytes on an Acclaim Trinity P1 column using HILIC conditions. CAD rather than UV absorbance was used to quantify these compounds as they lack a strong [chromophore](#).

Column: Thermo Scientific Acclaim Trinity P1,  
3  $\mu$ m, 2.1  $\times$  10 mm


Mobile phase: 81:19 Acetonitrile/  
10 mM Ammonium formate, pH = 3

Detector: Charged Aerosol



### HPLC-CAD analysis of steviol glycoside and mogroside V standards standards.

- |                   |                   |                |
|-------------------|-------------------|----------------|
| 1. Dulcoside A    | 4. Rebaudioside A | 7. Mogroside V |
| 2. Stevioside     | 5. Steviolbioside |                |
| 3. Rebaudioside C | 6. Rebaudioside B |                |



## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

A phenol is a compound that contains a hydroxyl group attached to an aromatic hydrocarbon group. Many phenolic phytochemicals are found in plants and classified as either **simple phenols** or **polyphenols**, based on the number of phenol hydroxyl groups and aromatic rings in the molecule. Representative **simple phenols**, include: carvacrol (a phenol), salicylic acid (a phenolic acid), resveratrol (a stilbenoid), caffeic acid (a hydroxycinnamic acid), tyrosol (a phenylethanoid), and **flavonoids**. **Polyphenols** include condensed tannins (e.g., proanthocyanidins); hydrolyzable tannins (e.g., ellagitannins) and phlorotannins. However, terminology and classification can be complex with the term polyphenol erroneously used to refer to all compounds with more than one phenol group.

Phenols have strong [chromophores](#) so are typically measured using HPLC with UV absorbance detection. Presented here is a general gradient HPLC-UV method for the measurement of multiple phenols and other compounds in a wide variety of [botanicals](#).


Column: Thermo Scientific Acclaim 120 C18, 3  $\mu\text{m}$ , 3 x 150 mm

Mobile phase A: 20 mM Monobasic sodium phosphate, 3 % Acetonitrile, 0.2 % Tetrahydrofuran, pH 3.35

Mobile phase B: 20 mM Monobasic sodium phosphate, 50 % Acetonitrile, 10 % Tetrahydrofuran, pH 3.45

Mobile Phase C: 90 % Methanol

Detector: Absorbance, UV 275 nm

**HPLC-UV analysis of multiple phenols.**

- |                               |                       |                       |
|-------------------------------|-----------------------|-----------------------|
| 1. Gallic acid                | 16. Syringaldehyde    | 31. Trans-resveratrol |
| 2. 4-Hydroxybenzyl alcohol    | 17. Umbelliferone     | 32. Luteolin          |
| 3. p-Aminobenzoic acid        | 18. p-Coumaric acid   | 33. Cis-resveratrol   |
| 4. 3,4-Dihydroxybenzoic acid  | 19. Salicylic acid    | 34. Quercetin         |
| 5. Gentisic acid              | 20. Sinapic acid      | 35. Kaempferol        |
| 6. 2-Hydroxybenzyl alcohol    | 21. Ferulic acid      | 36. Isorhamnetin      |
| 7. 4-Hydroxybenzoic acid      | 22. Ellagic acid      | 37. Eugenol           |
| 8. Chlorogenic acid           | 23. Rutin             | 38. Isoxanthohumol    |
| 9. p-Hydroxyphenylacetic acid | 24. Ethyl-vanillin    | 39. Cavaicol          |
| 10. Catechin                  | 25. 4-Hydroxycoumarin | 40. Thymol            |
| 11. Vanillic acid             | 26. Hesperidin        | 41. Carnisole         |
| 12. 4-Hydroxybenzaldehyde     | 27. Naringin          | 42. Xanthohumol       |
| 13. Syringic acid             | 28. Rosemarinic acid  | 43. Carnosic acid     |
| 14. Caffeic acid              | 29. Fisetin           |                       |
| 15. Vanillin                  | 30. Myricetin         |                       |

AN1063: Targeted Analyses of Secondary Metabolites in Herbs, Spices, and Beverages Using a Novel Spectro-Electro Array Platform

PN70019: The Spectro-Electro Array: A Novel Platform for the Measurement of Secondary Metabolites in Botanicals, Supplements, Foods and Beverages - Part 1: Theory and Concepts

AN1077: Determination of phenolic compounds in apple orchard soil







# Phytoestrogens – isoflavones and their aglycones

## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

A [phytoestrogen](#), also called a dietary estrogen, is a [xenoestrogen](#) obtained from a number of dietary plants. [Phytoestrogens](#) are a diverse group of naturally occurring nonsteroidal phytochemicals that, because of their structural similarity with estradiol (17- $\beta$ -estradiol), cause mild estrogenic and/or antiestrogenic effects. [Phytoestrogens](#) include [flavonoids](#), (e.g., the isoflavone daidzein), [coumestans](#) (e.g., coumestrol), prenylflavonoids and [mammalian lignans](#) (enterodiol and enterolactone). Presented here is a rapid HPLC-CAD method for routine determination of isoflavones and [aglycones](#).



### Isoflavone Analysis

Column: C18, 3  $\mu$ m, 4.6  $\times$  35 mm  
Mobile phase : 15% acetonitrile in 0.1% Acetic acid  
Detector: Charged Aerosol


### Isoflavone [Aglycone](#) Analysis

Column: C18, 3  $\mu$ m, 4.6  $\times$  35 mm  
Mobile phase : 30% acetonitrile in 0.1% Acetic acid  
Detector: Charged Aerosol

	R <sub>1</sub>	R <sub>2</sub>	R <sub>3</sub>
Daidzein	H	H	H
Daidzin	H	H	Glucose
Genistein	OH	H	H
Genistin	OH	H	Glucose
Glycitein	H	OCH <sub>3</sub>	H
Glycitin	H	OCH <sub>3</sub>	Glucose



HPLC-CAD analysis of soy isoflavones.



HPLC-CAD analysis of soy isoflavones and [aglycones](#).

LPN2930: Determination of the composition of natural products by HPLC with charged aerosol detection

For an HPLC-UV method see: An improved separation of isoflavones in red clover using a Thermo Scientific Acclaim 120 C18 HPLC column



## Table of contents

Summary

Overview:  
Dietary supplements and  
botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of  
supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary


Foods with the highest relative [phytoestrogen](#) content include nuts and oilseeds, soy products, cereals and breads, legumes, vegetables, and fruits. Their consumption may offer a range of health benefits with positive impacts on the cardiovascular, metabolic and central nervous systems, reduction in the risk of cancer, and by improving post-menopausal symptoms. However, [phytoestrogens](#) may also act as endocrine disruptors, adversely affecting health. Based on currently available evidence, it is not clear whether any potential health benefits of [phytoestrogens](#) outweigh their risks.

Presented here is an HPLC-CAD method for analysis of [coumestans](#) and [mammalian lignans](#).


Isoflavone [Aglycone](#) AnalysisColumn: C18, 3  $\mu\text{m}$ , 4.6  $\times$  35 mm

Mobile phase : 30% acetonitrile in 0.1% Acetic acid


Detector: Charged Aerosol




Coumestrol



Enterodiol



Enterolactone

HPLC-CAD analysis of [coumestans](#) and [mammalian lignans](#).

LPN2930: Determination of the composition of natural products by HPLC with charged aerosol detection

For a UHPLC-UV-CAD-MS approach see: [Profiling and quantitating the constituents of red clover extracts using UHPLC/UV-CAD/HRMS: A component of the safety assessment process](#)





## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of  
supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z


Literature

Glossary

Phytosterols (sterols and stanols) are steroid alcohols that are found in plants. They are key structural components of plant cell membranes, assuming the role that cholesterol plays in mammalian cells. In foods, phytosterols exist in the free form, as esters with fatty acids, and as glycolipids. Good food sources for phytosterols include unrefined vegetable oils, whole grains, nuts, seeds and legumes. Phytosterols are purported to have health benefits, such as lowering cholesterol and positively impacting cardiovascular disease. Consequently, foods and beverages supplemented with phytosterols are available in many countries. However, consumption of phytosterol-enriched foods may have undesirable effects too, such as reduced plasma carotenoid concentrations.


Phytosterols lack a good [chromophore](#), so UV detection below 210 nm is typically used for their analysis. Gas chromatographic approaches are sometimes used, but require extensive sample preparation (hydrolysis and derivatization).

Presented here is a direct, simple and sensitive HPLC-CAD approach for the measurement of several phytosterols in red palm oil extracts.



HPLC-CAD analysis of standards (blue) and red palm oil extract (red).

Column: Porous Shell C8, 2.7  $\mu\text{m}$ , 4.6  $\times$  150 mm,  
Mobile Phase A: Methanol/Water/Acetic acid (750 : 250 : 4)  
Mobile Phase B: Acetone/Methanol/Tetrahydrofuran/  
Acetic acid (500 : 375: 125 : 4)  
Detector: Charged Aerosol





## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of  
supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

Pomegranate (*Punica granatum*) fruit has been used since ancient times for medicinal purposes. Today, consumption of pomegranate juice is purported to have several health benefits (e.g., positively impacting arthritis, cancer and cardiovascular disease), possibly due to its high antioxidant content. Ellagitannins (hydrolysable tannins) are a family of water-soluble bioactive [polyphenols](#) particularly abundant in pomegranates. Amongst them, the punicalagins are reported to be responsible for more than half the potent antioxidant activity of the juice.


Presented here is an HPLC-UV absorbance method for the determination of punicalagin A and B.

Column: Thermo Scientific Acclaim Polar Advantage PA2,  
3  $\mu$ m, 3  $\times$  150 mm,


Mobile Phase A: 1% Formic acid

Mobile Phase B: Acetonitrile

Detector: Absorbance, UV 260 nm



**Analysis of punicalagin standards.**





## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

Pyrethrins are [terpenoid](#) esters derived from the flowers of *Chrysanthemum cinerariifolium*. The pyrethrin family includes six similar compounds that contain a cyclopropane core (see below). Pyrethrins both repel and kill insects by delaying the closure of voltage-gated sodium ion channels in their nerve cells. The insecticidal and insect repellent properties of these compounds have been known for millennia and chrysanthemum species have long been cultivated for this purpose. Interest in using pyrethrin-based insecticides is growing because of their low toxicity to humans (allowing home use) and favorable, fast biodegradability. However, they are also toxic to bees, with fatal doses as low as 0.02 µg, thus requiring very cautious application. Increasing pyrethrin use in agricultural and consumer products means there is a need for improved analytical techniques, both to assure product quality and to monitor the fate of pyrethrins in the environment.

Pyrethrins contain weak [chromophores](#), and are often measured using insensitive low-wavelength UV detection. Furthermore, published methods typically require excessively long run times to fully resolve closely eluting components.


Presented here is a comprehensive UHPLC-UV-CAD method that enables the resolution and detection of more compounds in pyrethrum oil in less time than previously possible.

Column: Thermo Scientific Acclaim Vanquish C18, 1.5 µm, 2.1 × 100 mm

Mobile phase A: DI Water


Mobile phase B: Acetonitrile, Optima LCMS

Detectors: Charged Aerosol and Absorbance, UV 220 nm



UHPLC-CAD analysis of pyrethrum oil extract.

	R <sub>1</sub>	R <sub>2</sub>
Cinerin I	CH <sub>3</sub>	CH <sub>3</sub>
Cinerin II	CO <sub>2</sub> CH <sub>3</sub>	CH <sub>3</sub>
Jasmolin I	CH <sub>3</sub>	CH <sub>2</sub> CH <sub>3</sub>
Jasmolin II	CO <sub>2</sub> CH <sub>3</sub>	CH <sub>2</sub> CH <sub>3</sub>
Pyrethrin I	CH <sub>3</sub>	CH=CH <sub>2</sub>
Pyrethrin II	CO <sub>2</sub> CH <sub>3</sub>	CH=CH <sub>2</sub>





## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of  
supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

*Schisandra chinensis* (Turcz.) Baill, (commonly known as Chinese magnolia-vine or schisandra) produces a fruit called magnolia berry or five-flavor-fruit (wù wèi zi). It is used both as a food and in [traditional Chinese medicine](#), where it is considered one of the 50 fundamental herbs. Schisandra has many proposed health benefits. It is promoted as an "adaptogen" for increasing resistance to disease and stress, boosting energy, and improving physical performance. It is also used for treating liver disease (hepatitis) and protecting the liver from poisons.

The major active components found in Schisandra are the lignanoids - schizandrin, schizandrin A and schizandrin B. The [Pharmacopeia of the People's Republic of China 2010](#) regulates its quality control using a UHPLC method for determining levels of these three analytes ([see AB139](#)).


Presented here is an efficient UHPLC method that can be used to measure schizandrin, schizandrin A and schizandrin B, for the quality control of Hupan tablets, a [traditional Chinese medicine](#) used for treating nonalcoholic fatty liver disease.

Column: Thermo Scientific Accucore C18, 2.6  $\mu\text{m}$ , 4.6  $\times$  150 mm


Mobile Phase A: DI Water

Mobile Phase B: Acetonitrile, Optima LCMS


Detector: Absorbance, UV 210 nm




**HPLC-UV analysis of schizandrin, schizandrin A and schizandrin B. A - standards; B - "Hupan" capsule extract.**



Schizandrin



Schizandrin A



Schizandrin B



## Table of contents

Summary

Overview:  
Dietary supplements and  
botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of  
supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary


Semen raphani, the seed of *Raphanus sativus* L., is a [Chinese medicinal](#) plant commonly used for treatment of dyspepsia and hypertension. Its anti-hypertension properties are attributed to the presence of sinapine thiocyanate. The [Pharmacopeia](#) of the People's Republic of China 2010 monitors the quality control of semen raphani with a RP-HPLC method for the determination of sinapine thiocyanate. The method specifies a stationary phase with phenyl groups bonded to silica ([see AB126](#)).

Presented here is an HPLC-UV method for the determination of sinapine thiocyanate in semen raphani seeds using an Acclaim Phenyl-1 column and UV absorbance detection.

Column: Thermo Scientific Acclaim Phenyl-1,  
3.0  $\mu\text{m}$ , 4.6  $\times$  150 mm

Mobile phase: Acetonitrile/3% Acetic acid, 10/90 (v/v)

Detector: Absorbance, 326 nm



**Separation of sinapine thiocyanate on the Acclaim Phenyl-1 column following the Chinese Pharmacopeial method. A) mobile phase, B) sinapine thiocyanate standard (10  $\mu\text{g}/\text{mL}$ ), C) semen raphani sample, and D) semen raphani sample spiked with sinapine thiocyanate (5  $\mu\text{g}/\text{mL}$ ).**

## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary


Stevia plant and extracts from stevia leaves have long been used as sweeteners in Asia and Latin America. Since the FDA approved the steviol [glycoside](#) rebaudioside A, purified from *Stevia rebaudiana* (Bertoni), as [Generally Recognized as Safe](#) for use as a sugar substitute, stevia products have become common table-top and beverage sweeteners.

Stevia plants contain many [terpene glycosides](#) that have different flavor profiles with both sweet and unpleasant bitter flavors. Two steviol [glycosides](#), stevioside and rebaudioside A, are largely responsible for the desired sweet flavor of the leaves, with rebaudioside A preferred for sweeteners.


Steviol [glycoside](#) determination is challenging for many reasons. The [glycoside](#) structures are quite similar, differing only in small changes in glycosylation which makes chromatographic separation difficult. Furthermore, they do not absorb strongly in the UV spectrum, and typical detection wavelengths for steviol [glycosides](#), such as 210 nm, are nonspecific and lack sensitivity. CAD can be used together with UV detection to improve steviol [glycoside](#) quantification. In addition, CAD has the advantage of measuring additional components in the sample that are not UV absorbing.

Presented here is a HPLC-CAD method for measurement of steviol [glycoside](#) and related compounds in [botanicals](#) and stevia products.

Column: Thermo Scientific Acclaim Trinity P1, 3  $\mu\text{m}$ , 2.1  $\times$  150 mm  
 Mobile Phase: 88:12 (v/v) Acetonitrile:10 mM Ammonium formate, pH 3  
 Detection: Charged Aerosol



Analysis of stevia [glycosides](#) by HPLC-CAD.



AN70278: [Analysis of Commercially Available Products Containing Stevia](#)

AN293: [Steviol Glycoside Determination by HPLC with Charged Aerosol and UV Detections Using the Acclaim Trinity P1 Column](#)

AN1048: [Novel, Universal Approach for the Measurement of Natural Products in a Variety of Botanicals and Supplements](#)





## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z


Literature

Glossary

St. John's Wort (*Hypericum perforatum*) is used in [folk medicine](#) to treat a variety of conditions, including depression, kidney and lung ailments, and insomnia, and to aid wound healing. Currently, St. John's Wort is most often used as a dietary supplement to treat mild depression. [However, studies have shown that it is not consistently effective.](#) Unfortunately, not only can it cause photosensitivity it may also interact with a numerous medications, causing major health issues.


Many biologically active compounds have been isolated from the extracts of St. John's wort, but the most potent appear to be the phloroglucinol-derivatives hyperforin and adhyperforin, and the naphthodianthrone hypericin and pseudohypericin.

Presented here is an HPLC-UV method capable of profiling key components in St. John's Wort supplements.



HPLC-UV analysis of a St. John's Wort extract.

Column:	Thermo Scientific Acclaim 120 C18, 3 $\mu$ m, 3 x 150 mm
Mobile phase A:	20 mM Monobasic sodium phosphate, 3 % Acetonitrile, 0.2 % Tetrahydrofuran, pH 3.35
Mobile phase B:	20 mM Monobasic sodium phosphate, 50 % Acetonitrile, 10 % Tetrahydrofuran, pH 3.45
Mobile Phase C:	90 % Methanol
Detector:	Absorbance, UV 254 nm



AN 1063: Targeted Analyses of Secondary Metabolites in Herbs, Spices, and Beverages Using a Novel Spectro-Electro Array Platform

AN335: Accelerated Solvent Extraction (ASE) of Active Ingredients From Natural Products

AN346: Totally Automated Sample Preparation Using Accelerated Solvent Extraction (ASE) Coupled with Gilson ASPEC: The Determination of Dianthrone in St. John's Wort



## Table of contents

Summary

Overview:  
Dietary supplements and  
botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of  
supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

Paclitaxel, first isolated from the Pacific yew (*Taxus brevifolia*), was approved as [Taxo<sup>®</sup>](#) by the National Cancer Institute (NCI) in 1992 for the treatment of ovarian cancer. An analysis of paclitaxel and related compounds, including cephalomannine (related compound A), 10-deacetyl-7-epipaclitaxel (related compound B) and 7-epipaclitaxel (related compound C) by RP-HPLC was published by both the United States [Pharmacopeia 2009](#) and Chinese [Pharmacopeia 2010](#) (see [AB119](#)). These methods each required longer than 70 minutes.



Presented here is a UHPLC-UV method using an Acclaim RSLC C18. This method is capable of resolving all key analytes in under 6 minutes and was developed to address the demand for a more rapid analysis.

Column: Thermo Scientific Acclaim 120 C18, 2.2  $\mu\text{m}$ , 2.1 x 100 mm

Mobile phase A: Water

Mobile phase B: Acetonitrile/Methanol 60:40 v/v

Detector: Absorbance, UV 227 nm



Overlay of chromatograms of A) mixture of paclitaxel and related compounds standards (5  $\mu\text{g}/\text{mL}$  for each) and B) paclitaxel injection sample.



# Tomatine and tomatidine

## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

Tomatine (or lycopersicin) is a steroidal [glycoalkaloid](#) particularly abundant in the stems and leaves of tomato plants where it acts as a natural fungicide and insecticide. Tomatine and its [aglycone](#) tomatidine are reported to have [health benefits](#), including possessing antibiotic and anticarcinogenic properties, and having a positive impact on cardiovascular disease.

Tomatine and tomatidine possess weak [chromophores](#) making analysis by HPLC-UV difficult. Presented here is an HPLC-CAD method that overcomes such analytical limitations.

HPLC Column: Thermo Scientific Accucore C8, 2.6  $\mu\text{m}$ , 2.1  $\times$ 100 mm

Mobile Phase A: 0.2 % Formic acid in Deionized water

Mobile Phase B: 0.2 % Formic acid in Acetonitrile


Column Temp.: 35  $^{\circ}\text{C}$

Flow Rate: 1.0 mL/min


Injection Volume: 0.2  $\mu\text{L}$

Gradient: 5-95% B in 10 min

Detectors: Charged Aerosol and Absorbance, UV 227 nm



**HPLC-CAD can be used to measure tomatine and tomatidine along with a number of impurities in commercially available standards. (Each standard - 0.25 mg/mL in methanol containing 0.2% formic acid). Compounds cannot be measured by UV absorbance detection.**



For an LC-MS approach see: [Plant Metabolomics: Tomato Metabolite Profiling and Identification Employing High-Resolution LC-MS Strategies](#)





## Turmeric – curcuminoids

## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of  
supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

Turmeric is a popular spice and dye made from the powdered dried root of *Curcuma longa*. It is widely used as a culinary additive to impart a distinctive yellow-orange color to Pakistani, Indian and Thai cuisines. It has been used in Asia for thousands of years, and is a major part of [Ayurvedic](#), Siddha medicine, [traditional Chinese medicine](#) and Unani medicine. In [Ayurvedic](#) medicine, it is called haridra and used as an anti-inflammatory, and in the treatment of arthritis, cancer, gastric ulcers, neurodegenerative diseases and allergies. Its purported bioactive components are the brilliant yellow pigments curcumin, desmethoxycurcumin and bis-desmethoxycurcumin, together with other minor curcuminoids.


HPLC analysis with a C18 column is typically used for measurement of curcuminoids. While C18 may be satisfactory for some applications, the selectivity of polar-embedded stationary phases, such as the Accucore Polar Premium phase provides a superior solution, capable of completely resolving the curcuminoids in under 3 minutes.

Presented here is a rapid HPLC-UV method for the measurement of three curcuminoids in spice extracts.

Column: Thermo Scientific Accucore Polar Premium, 2.6  $\mu\text{m}$ , 3.0 x 100 mm




Mobile phase: Methanol: 10 mM Phosphoric acid (80:20 v/v)

Detector: Absorbance, Vis 428 nm



**Analysis of ethanolic turmeric extract analyzed using HPLC with Vis absorbance.**

1. Curcumin
2. Desmethoxycurcumin
3. Bis-desmethoxycurcumin



AN20853: Separation of Curcuminoids from Turmeric – Comparison of Polar Embedded and C18 Solid HPLC Core Columns

PN70677: The Quantitative Analysis of Curcuminoids in Food and Food Additives Using Rapid HPLC With Electrochemical, UV, or Fluorescence Detection

AppsLab: Rapid analysis of pigments in turmeric on a Thermo Scientific Acclaim PolarAdvantage II (PA2) column





## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z


Literature

Glossary



Ursoleic acid is a pentacyclic [triterpenoid](#) found in the epicuticular waxes and peels of fruits and herbs, such as rosemary and thyme, but is particularly abundant in apple peel. There is a growing interest in ursolic acid because of its [purported health benefits](#), which include antioxidant, anticarcinogenic and anti-inflammatory effects.

Ursoleic acid and related compounds do not have a strong [chromophore](#), so HPLC-UV methods lack sensitivity. Presented here is a sensitive HPLC-CAD method using a C30 column that easily overcomes the poor resolution between ursolic acid and the structurally similar triterpenoid oleanolic acid that is often seen with other HPLC methods.

HPLC Column: Thermo Scientific Acclaim C30, 5  $\mu$ m, 4.6  $\times$  250 mm  
 Mobile Phase A: 1% aq Ammonium acetate  
 Mobile Phase B: Acetonitrile/Methanol 3:1 v/v  
 Detectors: Charged Aerosol



Measurement of four triterpenoids by HPLC-CAD. Black – standards (500 ng each on column); Blue – apple peel extract.



Ursoleic acid	$R_1 = \text{CH}_3$	$R_2 = \text{H}$	$R_3 = \text{CO}_2\text{H}$
Oleanolic acid	$R_1 = \text{H}$	$R_2 = \text{CH}_3$	$R_3 = \text{CO}_2\text{H}$
Glycyrrhetic acid	$R_1 = \text{H}$	$R_2 = \text{CO}_2\text{H}$	$R_3 = \text{H}$ Ketone at (*)



# Vitamins – simultaneous determination of water and fat soluble

## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature


Glossary

Vitamins are essential nutrients found in various natural foods and food supplements. Vitamins can be classified as water-soluble vitamins (WSV) or fat-soluble vitamins (FSV), based on their hydrophobicity. RP-HPLC is widely used to determine vitamins in food, supplements and beverages. The dramatically different hydrophobicity of WSV and FSV makes simultaneous liquid chromatography analysis with the same method difficult, so typically, separate methods are used for each class. To address this, a method for the separation of FSV and WSV was developed with two columns operated sequentially on one system. However, this required a complex hardware setup, and even more complex chromatography data system programming. Additionally, since the columns were run sequentially, sample throughput was low.

The work presented [here](#) shows an improved method capable of the simultaneous analysis of WSV and FSV in the same sample. The workflow is based on a novel Vanquish Flex Duo system for Dual LC. The system enables the independent and simultaneous use of two different columns and methods. Compared to the previous solution, this approach is remarkably simple to implement, allows the use of better optimized methods, and increases throughput thanks to the simultaneous use of two columns with two methods, and faster analysis cycles.


AN72592: [Simultaneous Determination of Water- and Fat-Soluble Vitamins in Tablets and Energy Drinks by Using a Novel Vanquish Flex Duo System for Dual LC](#)

AN72877: [Accelerated Method Development for the Separation of Water-soluble Vitamins by RP HPLC with UV Detection using an Empirical Approach to Predict Separation](#)



**Separation of water-soluble vitamins in a supplement tablet. Blue trace. 210 nm Black trace. 270 nm.**

- |                  |                     |                |
|------------------|---------------------|----------------|
| 1. Ascorbic acid | 4. Nicotinamide     | 7. Riboflavin. |
| 2. Thiamine      | 5. Pantothenic acid |                |
| 3. Pyridoxin     | 6. Folic acid       |                |



**Comparison of the standard mixture (black) with 100 µg/mL and the vitamin tablet (blue) at a wavelength of 280 nm.**

- |                    |                                 |                         |
|--------------------|---------------------------------|-------------------------|
| 1. Retinol         | 5. Menaquinone                  | 9. $\alpha$ -Tocopherol |
| 2. Retinyl acetate | 6. $\delta$ -Tocopherol         | 10. Phyllochinone       |
| 3. Ergocalciferol  | 7. $\alpha$ -Tocopheryl acetate |                         |
| 4. Cholecalciferol | 8. $\gamma$ -Tocopherol         |                         |



## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z


Literature

Glossary


Vitamin D, a fat-soluble vitamin, occurs naturally in a few foods (e.g., fatty fish and fish oils), is added to others and can also be consumed as a dietary supplement. It is also produced endogenously when skin is exposed to sunlight, triggering synthesis. Vitamin D is biologically inert and must undergo biochemical modification in the body for activation. The first reaction occurs in the liver, converting vitamin D to calcidiol (25-hydroxyvitamin D). The second reaction occurs primarily in the kidney and forms the physiologically active Vitamin D3 – calcitriol (1,25-dihydroxyvitamin D). Vitamin D is important for maintaining calcium levels and promoting bone health, and has other roles in the body, including modulation of cell growth, neuromuscular and immune function, and reduction of inflammation.

Presented here is an HPLC method using HILIC separation and UV absorbance detection to measure Vitamin D3 in an oil-based supplement.

HPLC Column: Thermo Scientific Acclaim HILIC-10, 5  $\mu$ m, 4.6  $\times$  150 mm  
 Mobile Phase: Heptane/Acetonitrile/Propan-2-ol/Acetic acid (98.4:1:0.5:0.1 v/v/v/v)  
 Detectors: Absorbance, UV 265 nm

**Analysis of a vitamin D3 in an oil-based supplement.**

1. Excipient
2. Excipient
3. Excipient
4. Vitamin D3



## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

*Zanthoxylum nitidum* (Roxb.) DC (ZN) (liang mian zhen) is an important [traditional Chinese medicine](#). The [Pharmacopeia](#) of the People's Republic of China (PPRC) 2010 regulates this dried root as an [herbal medicine](#) (see [AN 1008](#)). *ZN var. fastuosum* (ZNF) is another plant in the same genus as ZN. Although ZNF is not regulated in the PPRC, its dried root is still used in [Chinese folk medicine](#) because some of its reported medical benefits are the same as those reported for ZN.

The major active components of ZN and ZNF are [alkaloids](#). Nitidine and toddalolactone are the specific active components of ZN and ZNF, respectively. The PPRC 2010 method regulates the quality control of ZN using a HPLC method for determination of nitidine chloride but uses a thin-layer chromatography method for detection of chelerythrine chloride and toddalolactone. The presence of toddalolactone in ZN is not permitted.


The HPLC-UV method presented here is an efficient and comprehensive method for the quality control analysis of ZN through quantification of the main active components: nitidine chloride, chelerythrine chloride and toddalolactone.

Column: Thermo Scientific Acclaim LC PA, 3  $\mu$ m, 2.1  $\times$  150 mm

Mobile phase A: 25 mM Ammonium acetate (pH 4.5)


Mobile Phase B: Acetonitrile

Detector: Absorbance, UV 273 nm




#### Analysis of a *Zanthoxylum nitidum* (Roxb.) DC sample.


- A. Blank  
 B. Standards (20  $\mu$ g/mL each) and  
 C. *Zanthoxylum Nitidum* (Roxb.) sample  
 1. Toddalolactone  
 2. Nitidine chloride  
 3. Chelerythrine chloride.



**Nitidine chloride**



**Chelerythrine chloride**



**Toddalolactone**



# Thermo Scientific references

	Title	Authors	Publication
Table of contents			
Summary	<a href="#">Why use Charged Aerosol Detection with Inverse Gradient?</a>	Grosse, S.; Muellner, T.; Lovejoy, K.; Acworth, I.; Gamache, P.	TN73449, 2020
Overview: Dietary supplements and botanical natural products	<a href="#">Achieving standard free quantitation: Thermo Scientific Charged Aerosol Detectors</a>	Anon	SP73026, 2019
Overview: Global market	<a href="#">Determination of olive oil purity based on triacylglycerols profiling by UHPLC-CAD and Principal Component Analysis</a>	Green, H.; Li, Q.; De Pra, M.; Acworth, I.; Wang, S.	AN73174, 2019
Measurement and analysis	<a href="#">Simultaneous determination of water- and fat-soluble vitamins in tablets and energy drinks by using a novel Vanquish Flex Duo system for Dual LC</a>	Grosse, S.; De Pra, M.; Steiner, F.	AN72592, 2018
Instrumentation	<a href="#">Accelerated Method Development for the Separation of Water-soluble Vitamins by RP-HPLC with UV Detection using an Empirical Approach to Predict Separation</a>	Grosse, S.; Park, S.; De Pra, M.; Steiner, F.	AN72877, 2018
Sample preparation	<a href="#">Charged Aerosol Detection – Factors Affecting Uniform Response</a>	Menz, M.; Eggart, B.; Lovejoy, K.; Acworth, I.; Gamache, P.; Steiner, F.	TN72806, 2018
Separation	<a href="#">Totally Automated Sample Preparation Using Accelerated Solvent Extraction (ASE) Coupled with Gilson ASPEC: The Determination of Dianthrone in St. John's Wort</a>	Anon	AppsLab, 2017
Detection	<a href="#">Extraction of Herbal Marker Compounds Using Accelerated Solvent Extraction Compared to Traditional Pharmacopoeia Protocols</a>	Anon	AN362, 2017
Authentication of supplements	<a href="#">Quantitative Analysis of Toosendanin in the Fruit of <i>Melia toosendan</i> Sieb. Et Zucc (Meliaceae) by High Performance Liquid Chromatography Coupled with Charged Aerosol Detection</a>	Anon	AppsLab, 2017
Application examples	<a href="#">A Non-derivative Method for the Quantitative Analysis of Isosteroidal Alkaloids from Fritillaria by High Performance Liquid Chromatography Combined with Charged Aerosol Detect</a>	Anon	AppsLab, 2017
Substances A-C	<a href="#">Rapid Analysis of Natural Sweeteners Found in Food and Beverages Using an Advanced UHPLC System</a>	Hillbeck, D.	AppsLab, 2017
Substances D-G	<a href="#">Sensitive HILIC UHPLC-UV determination of steviol glycoside natural sweeteners</a>	Lamb, A.; Jones, J.	AppsLab, 2017
Substances H-M	<a href="#">Traditional Chinese Medicine HPLC Interactive Applications Notebook</a>	Acworth, I.	AN71120, 2016
Substances N-T	<a href="#">Novel, Universal Approach for the Measurement of Natural Products in a Variety of Botanicals and Supplements</a>	Acworth, I.; Bailey, B.; Plante, M.; Crafts, C.; Thomas, D.; Roman, M.	AN1048, 2016
Substances U-Z	<a href="#">Profiling Hoodia Extracts by HPLC with Charged Aerosol Detection, Electrochemical Array Detection, and Principal Component Analysis</a>	Acworth, I.; Bailey, B.; Plante, M.; Zhang, Q.; Thomas, D.	PN70540, 2016
Literature	<a href="#">Fast Determination of Anthocyanins in Pomegranate Juice</a>	Anon	AN264, 2016
Glossary	<a href="#">Sensitive Determination of Catechins in Tea by HPLC</a>	Anon	AN275, 2016
	<a href="#">Separation of Sinapine Thiocyanate in Semen Raphani Using an Acclaim Phenyl-1 Column</a>	Anon	AB126, 2016
	<a href="#">Determination of A-Type and B-Type Procyanidins in Apple, Cocoa and Cinnamon Extracts</a>	Glinski, J.; Thomas, D.; Wong, A.; Glinski, V.; Acworth, I.	PN71527, 2016
	<a href="#">Savor the Flavor – Robust Iso-a-acids Assaying in Beer Within Ten Minutes</a>	Heidorn, M.	AB153, 2016
	<a href="#">Monitor the Brewing Process with LC – Transformation of Hop <math>\alpha</math>-acids into Beer iso-<math>\alpha</math>-acids</a>	Heidorn, M.	AB155, 2016
	<a href="#">The Everlasting Paradigm – Keep Beer Tradition or Prevent Beer From a Skunky off-flavor?</a>	Heidorn, M.	AB156, 2016
	<a href="#">Mogroside V Determination by HPLC with Charged Aerosol and UV Detections</a>	Hurum, D.; Rohrer, J.	AU184, 2016

# Thermo Scientific references

	Title	Authors	Publication
Table of contents			
Summary	<a href="#">Steviol Glycoside Determination by HPLC with Charged Aerosol and UV Detections Using the Acclaim Trinity P1 Column</a>	Hurum, D.; Rohrer, J.	AN293, 2016
Overview: Dietary supplements and botanical natural products	<a href="#">Determination of Nitidine Chloride, Toddalolactone, and Chelerythrine Chloride by HPLC</a>	Jing, C.; Qun, X.; Rohrer, J.	AN1008, 2016
Overview: Global market	<a href="#">Natural Products Extraction Using Accelerated Solvent Extraction. White Paper</a>	Kettle, A.	WP71206, 2016
Measurement and analysis	<a href="#">Rapid Separation of Anthocyanins in Cranberry and Bilberry Extracts Using a Core-Shell Particle Column</a>	Perati, P.; De Borba B.; Rohrer, J.	AN1042, 2016
Instrumentation	<a href="#">Rapid Separation of Catechins in Tea Using Core-shell Columns</a>	Perati, P.; De Borba B.; Rohrer, J.	AB150, 2016
Sample preparation	<a href="#">Improved Universal Approach to Measure Natural Products in a Variety of Botanical and Supplements</a>	Plante, M.; Acworth, I.; Bailey, B.; Zhang, Q.; Thomas, D.	PN70543, 2016
Separation	<a href="#">HPLC Analysis of Six Active Components of Caulis Lonicerae Using a Phenyl-1 Column</a>	Qun, X.; Xiongfeng, H.; Rohrer, J.	AN299, 2016
Detection	<a href="#">Separation of Schizandrin, Schizandrin A, and Schizandrin B in a Tablet Sample</a>	Qun, X.; Xiongfeng, H.; Rohrer, J.	AB139, 2016
Authentication of supplements	<a href="#">Determination of Pyrethrins in Pyrethrum Oil Extracts by UHPLC with Charged Aerosol Detection</a>	Thomas, D.; Glinski, J.; Wong, A.; Acworth, I., Mohindra, D.	PN21431, 2016
Application examples	<a href="#">Separation of Curcuminoids from Turmeric – Comparison of Polar Embedded and C18 Solid HPLC Core Columns</a>	Tracy, M.	AN20853, 2016
Substances A-C	<a href="#">Separation of Calcium, Magnesium and Counterions in a Dietary Supplement Using Counterions in a Dietary Supplement Using Multi-mode Liquid Chromatography with Charged Aerosol Detection</a>	Tracy, M.; Liu, X.	AN20871, 2016
Substances D-G	<a href="#">The Spectro-Electro Array: A Novel Platform for the Measurement of Secondary Metabolites in Botanicals, Supplements, Foods and Beverages - Part 1: Theory and Concepts</a>	Ullucci, P.; Acworth, I.; Bailey, B.; Crafts, C.; Plante, M.	PN70019, 2016
Substances H-M	<a href="#">Targeted Analyses of Secondary Metabolites in Herbs, Spices, and Beverages Using a Novel Spectro-Electro Array Platform</a>	Ullucci, P.; Acworth, I.; Crafts, C.; Bailey, B.; Plante, M.	AN1063, 2016
Substances N-T	<a href="#">Rapid Determination of Hesperidin in Orange Peel Using Accelerated Solvent Extraction and UHPLC</a>	Xiongfeng, H.; Qun, X.; Jinshui, C.; Rohrer, J.	AB142, 2016
Substances U-Z	<a href="#">Cystine, an Essential Determinant of Protein Tertiary Structure, Is Also a Target for Electrochemical Manipulation</a>	Zhang, Q.; Bailey, B.; Plante, M.; Acworth, I.	PN71020, 2016
Literature	<a href="#">Fast Analysis of Selected Xanthones in Mangosteen Pericarp Using Accelerated Solvent Extraction and Ultra High Performance Liquid Chromatography</a>	Zhang, Q.; Bailey, B.; Plante, M.; Acworth, I.	PN70991, 2016
Glossary	<a href="#">The Quantitative Analysis of Curcuminoids in Food and Food Additives Using Rapid HPLC With Electrochemical, UV, or Fluorescence Detection</a>	Zhang, Q.; Thomas, D.; Acworth, I.	PN70677, 2016
	<a href="#">The Vanquish Platform: Major Improvement in Throughput and Resolution of Xanthones in Mangosteen Pericarp</a>	Zhang, X.; Bailey, B.; Plante, M.; Acworth, I.	AB172, 2016
	<a href="#">Improved analysis of resveratrol and related substances in dietary supplements using a Thermo Scientific Acclaim 120 C18 HPLC column</a>	Anon	AppsLab, 2015
	<a href="#">Analysis of alkaloids in bitter orange extract by HPLC-UV</a>	Anon	AppsLab, 2015
	<a href="#">Comparison of six stationary phases for the separation of catechins</a>	Anon	AppsLab, 2015
	<a href="#">Chromatography for Foods and Beverages: Vitamin and Antioxidant Applications Notebook</a>	Acworth, I.	AI71478, 2015
	<a href="#">USP Monograph: Content of flavonol glycosides in ginkgo tablets using C18 HPLC column</a>	Anon	AppsLab, 2015
	<a href="#">Determination of Hoodigosides by UHPLC-CAD on Vanquish</a>	Thomas, D.	AppsLab, 2015

# Thermo Scientific references

	Title	Authors	Publication
Table of contents	<a href="#">Aristolochic acids</a>	Anon	AppsLab, 2014
Summary	<a href="#">Separation of artemether and its impurities using reversed-phase HPLC-UV</a>	Anon	AppsLab, 2014
Overview: Dietary supplements and botanical natural products	<a href="#">An improved separation of isoflavones in red clover using a Thermo Scientific Acclaim 120 C18 HPLC column</a>	Anon	AppsLab, 2014
	<a href="#">Rapid analysis of pigments in turmeric on a Thermo Scientific Acclaim PolarAdvantage II (PA2) column</a>	Anon	AppsLab, 2014
	<a href="#">Simple gradient method for the analysis of calcium and magnesium in a dietary mineral supplement using HPLC-CAD</a>	Anon	AppsLab, 2014
Overview: Global market	<a href="#">Rapid analysis of vitamin D3 in supplements using the Thermo Scientific Acclaim HILIC-10 HPLC column</a>	Anon	AppsLab, 2014
Measurement and analysis	<a href="#">Rapid determination of Carnitine in a nutritional supplement using a Thermo Scientific Acclaim Trinity P1 HPLC column</a>	Anon	AppsLab, 2014
Instrumentation	<a href="#">Rapid determination of Acetylcarnitine and Lipoic Acid in a nutritional supplement using a Thermo Scientific Acclaim Trinity P1 HPLC column</a>	Anon	AppsLab, 2014
	Sample preparation	<a href="#">An improved analysis of alkaloids in goldenseal root on a Thermo Scientific Acclaim 120 C18 HPLC column with Thermo Scientific Dionex ASE Accelerated Solvent Extractor</a>	Anon
Separation	<a href="#">HPLC-UV method for the rapid determination of alkaloids</a>	Chander, P.	AppsLab, 2014
Detection	<a href="#">Determination of Phenolic Compounds in Apple Orchard Soil</a>	Jinshui, C.; Qun, X.; Lina, L.; Rohrer, J.	AN1077, 2014
Authentication of supplements	<a href="#">Separation of curcuminoids from turmeric - comparison of polar embedded and C18 solid core HPLC columns</a>	Anon	AppsLab, 2013
Application examples	<a href="#">Novel, Universal Approach for the Measurement of Natural Products in a Variety of Botanicals and Supplements. Part 2</a>	Bousova, K.; Acworth, I.; Bailey, B.; Plante, M.; Zhang, Q.; Thomas, D.; Roman, M.	PN70543, 2013
	Substances A-C	<a href="#">Rapid analysis of artesunate and dihydroartemisinin using a Thermo Scientific Accucore RP-MS HPLC column</a>	Anon
Substances D-G	<a href="#">Fast analysis of artesunate and dihydroartemisinin using a Thermo Scientific Synchronis C18 HPLC column</a>	Anon	AppsLab, 2012
Substances H-M	<a href="#">Fast analysis of arbutin using a Thermo Scientific Hypersil GOLD aQ HPLC column</a>	Anon	AppsLab, 2012
Substances N-T	<a href="#">Improved analysis of glycosides using a Thermo Scientific Hypersil GOLD aQ HPLC column</a>	Anon	AppsLab, 2012
Substances U-Z	<a href="#">Can High Peak Capacity and Universal Detection Solve the Challenges in LC Characterization of Botanicals and Natural Products</a>	Bauder, R.; Steriner, F.; Heidorn, M.; Martin, M.; McCleod, F.	PN70153, 2012
Literature	<a href="#">Analysis of Commercially Available Products Containing Stevia</a>	Crafts, C.; Bailey, B.; Plante, M.; Acworth, I.	AN70278, 2012
Glossary	<a href="#">Determination of Catechins and Phenolic Acids in Red Wine by Solid Phase Extraction and HPLC</a>	Dolci, M.	AN20583, 2012
	<a href="#">Analysis of Catechins Using an Accucore XL C8 4 µm HPLC Column</a>	Khan, A.	AN20536, 2012
	<a href="#">Sensitive HPLC Method for Triterpenoid Analysis Using Charged Aerosol Detection with Improved Resolution</a>	Plante, M.; Bailey, B.; Crafts, C.; Acworth, I.	PN70037, 2012

# Thermo Scientific references

Table of contents	Title	Authors	Publication
Summary	<a href="#">Phytosterols by HPLC with the Thermo Scientific Corona ultra Charged Aerosol Detection</a>	Plante, M.; Crafts, C.; Bailey, B.; Gamache, P.; Waraska, J.; Acworth, I.	AN1041, 2012
Overview: Dietary supplements and botanical natural products	<a href="#">Determination of the Composition of Natural Products by HPLC with Charged Aerosol Detection</a>	Acworth, I.; Bailey, B.; Gamache, P.; Waraska, J.	LPN2930, 2011
Overview: Global market	<a href="#">Simple and Direct Analysis of Falcarinol and other Polyacetylenic Oxylipins in Carrots by Reverse Phase HPLC and Charged Aerosol Detection</a>	Acworth, I.; Plante, M.; Bailey, B.; Craft, C.; Waraska, J.	LPN2923-01, 2011
Measurement and analysis	<a href="#">Rapid and Sensitive Determination of Anthocyanins in Bilberries Using UHPLC</a>	Anon	AN281, 2011
Instrumentation	<a href="#">Accelerated Solvent Extraction (ASE) of Active Ingredients From Natural Products</a>	Anon	AN335, 2011
Sample preparation	<a href="#">Extraction of Herbal Marker Compounds Using Accelerated Solvent Extraction Compared to Traditional Pharmacopoeia Protocols</a>	Anon	AN362, 2011
Separation	<a href="#">Determination of Triterpene Glycosides in <i>Cimicifuga racemosa</i> (Black Cohosh) by HPLC-CAD</a>	Roman, M.	CAN113, 2011
Detection	<a href="#">Determination of ginsenosides in <i>Panax ginseng</i> by HPLC-CAD</a>	Roman, M.	CAN112, 2011
Authentication of supplements	<a href="#">Determination of Triterpenes in <i>Centella asiatica</i> (Gotu Kola) by HPLC-CAD</a>	Roman, M.	CAN111, 2011
Application examples	<a href="#">Rapid Separation of Paclitaxel and Related Compounds in Paclitaxel Injection</a>	Anon	AB119, 2010
Substances A-C	<a href="#">Ultrafast determination of chlorophyll on a Thermo Scientific Acclaim PolarAdvantage II (PA2) RSLC column</a>	Anon	AppsLab, 2010
Substances D-G	<a href="#">Chromatographic fingerprinting of <i>Chrysanthema indicum</i> using HPLC</a>	Anon	AN207, 2009
Substances H-M	<a href="#">Determination of Anthraquinones and Stilbenes in Giant Knotweed Rhizome by HPLC with UV detection</a>	Anon	AN232, 2009
Substances N-T	<a href="#">Determination of the Punicalagins Found in Pomegranate by High Performance Liquid Chromatography</a>	Baugh, S.; Revell, J.; Eastman, K.	CAN106, 2009
Substances U-Z	<a href="#">Rapid Analysis of Ginseng Using Accelerated Solvent Extraction and HPLC</a>	Anon	AN192, 2007
Literature			
Glossary			

# Peer reviewed journal references

Table of contents	Title	Authors	Publication
Summary	<a href="#">Solvent and temperature effects of accelerated solvent extraction (ASE) coupled with ultra-high pressure liquid chromatography (UHPLC-DAD) technique for determination of thymoquinone in commercial food samples of black seeds (<i>Nigella sativa</i>)</a>	Ahmad, R.; Ahmad, N.; Shehzad, A.	Food Chem. <b>2020</b> , 309, 125740.
Overview: Dietary supplements and botanical natural products	<a href="#">Uncovering the antioxidant characteristics of black tea by coupling <i>in vitro</i> free radical scavenging assay with UHPLC-MS analysis</a>	Chen, N.; Han, B.; Fan, X.; Cai, F.; Ren, F.; Xu, M.; Zhong, J.; Zhang, Y.; Ren, D.; Yi, L.	J. Chromatogr. B <b>2020</b> , 1145, 122092.
Overview: Global market	<a href="#">Rapid Simultaneous Determination of Pentacyclic Triterpenoids by Mixed-Mode Liquid Chromatography-Tandem Mass Spectrometry</a>	Falev, D. I.; Kosyakov, D. S.; Ul'yanovskii, N. V.; Ovchinnikov, D. V.	J. Chromatogr. A <b>2020</b> , 1609, 460458.
Measurement and analysis	<a href="#">Determination of seven oligosaccharides and sucrose in <i>Pseudostellaria heterophylla</i> by pressurized liquid extraction and ultra-high performance liquid chromatography with charged aerosol detector and tandem mass spectrometry</a>	Hua, D.; Han, B.; Chen, C.; Chen, N.; Zhu, B.; Zhao, J.; Lia, S.	J. Chromatogr. A <b>2020</b> , 1609, 460441.
Instrumentation	<a href="#">Reinvestigation of <i>Herniaria glabra</i> L. saponins and their biological activity</a>	Kozachok, S.; Pecio, L.; Orhan, I. E.; Deniz, F. S. S.; Marchyshyn, S.; Oleszeka, W.	Phytochem. <b>2020</b> , 169, 112162.
Sample preparation	<a href="#">Quality assessment of <i>Moringa</i> seed shells based on fingerprinting using HPLC-DAD</a>	Li, X-F.; Shi, H-M.; Xu, M.; Meng, L.	Acta Chromatogr. <b>2020</b> , 32, 28-33.
Separation	<a href="#">Integration of micro-fractionation, high-performance liquid chromatography-ultraviolet detector-charged aerosol detector-mass spectrometry analysis and cellular dynamic mass redistribution assay to accelerate alkaloid drug discovery</a>	Wang, R.; Liu, Y.; Zhou, H.; Chen, Y.; Wang, J.; Zhang, X.; Yu, R.; Liang, X.	J. Chromatogr. A <b>2020</b> , 1616, 460779.
Detection	<a href="#">Preparation and Certification of a New Salvianolic Acid A Reference Material for Food and Drug Research</a>	Yang, D.; Su, B.; Bi, Y.; Zhang, L.; Zhang, B.; Song, J.; Lu, Y.; Du, G.	Nat. Prods. Bioprospect. <b>2020</b> , 10, 67-75.
Authentication of supplements	<a href="#">Development and validation of a novel high performance liquid chromatography-coupled with Corona charged aerosol detector method for quantification of glucosamine in dietary supplements</a>	Asthana, C.; Peterson, G. M.; Shastri, M. D.; Patel, R. P.	PLoS ONE, <b>2019</b> , 14, e0216039.
Application examples	<a href="#">Comprehensively qualitative and quantitative analysis of ginsenosides in <i>Panax notoginseng</i> leaves by online two-dimensional liquid chromatography coupled to hybrid linear ion trap Orbitrap mass spectrometry with deeply optimized dilution and modulation system</a>	Cao, J-L.; Ma, L-J.; Wang, S-P.; Deng, Y.; Wang, Y-T.; Li, P.; Wan, J-B.	Anal. Chim. Acta <b>2019</b> , 1079, 237-251.
Substances A-C	<a href="#">Nutritional composition, mineral content, antioxidant activity and quantitative estimation of water soluble vitamins and phenolics by RP-HPLC in some lesser used wild edible plants</a>	Datta, S.; Sinha, B. K.; Bhattacharjee, S.; Seal, T.	Helyon <b>2019</b> , 5, E01432.
Substances D-G	<a href="#">Characterization of Silver Fir Wood Decay Classes Using Sugar Metabolites Detected with Ion Chromatography</a>	Di Lella, S.; Tognetti, R.; La Porta, N.; Lombardi, F.; Nardin, T.; Larcher, R.	J. Wood Chem. Tech. <b>2019</b> , 39, 90-110.
Substances H-M	<a href="#">Isolation of secondary metabolites from <i>Geranium molle</i> L. with anticancer potential</a>	Graça, V. C.; Calheta, R.C.; Nunes, F. M.; Berthet, J.; Ferreira, I.; Santos, P. F.	Ind. Crops Prods. <b>2019</b> , 142, 111859.
Substances N-T	<a href="#">Determination of fatty acids in the seeds of <i>Lepidium apetalum</i> Willdenow, <i>Descurainia sophia</i> (L.) Webb ex Prantl, and <i>Draba nemorosa</i> L. by ultra-high-performance liquid chromatography equipped with a charged aerosol detector</a>	Kim, H. S.; Moon, B. C.; Sungyu Yang, Song, J-H.; Chun, J. M.; Kwon, B-I.; Yeong Lee, A.	J. Liq. Chromatogr. Rel. Technol. <b>2019</b> , 42, 128-136.
Substances U-Z	<a href="#">Determination of glucosinolates in broccoli-based dietary supplements by cyclodextrin-mediated capillary zone electrophoresis</a>	Lechtenberg, M.; Hensel, A.	J. Food Comp. Anal. <b>2019</b> , 78, 138-149.
Literature			
Glossary			

# Peer reviewed journal references

	Title	Authors	Publication
Table of contents			
Summary	<a href="#">Phytochemical analysis of <i>Brasolia</i>, <i>Eleanthus</i>, and <i>Sobralia</i>. Three genera of orchids with antibacterial potential against <i>Staphylococcus aureus</i></a>	Rykcaczevska, M.; Krauze-Baranowskab, M.; Żuchowski, J.; Krychowiak-Maśnickad, M.; Fikowicz-Krośkod, J.; Królickad, A.	Phytochem Letts. <b>2019</b> , 30, 245-253
Overview: Dietary supplements and botanical natural products	<a href="#">Impact of Nitrogen Fertilizer Levels on Metabolite Profiling of the <i>Lycium barbarum</i> L. Fruit</a>	Shi, Z.; Wei, F.; Wan, R.; Li, Y.; Wang, Y.; An, Y.; Qin, K.; Dai, G.; Cao, Y.; Feng, J.	Mol. <b>2019</b> , 24, 3879.
Overview: Global market	<a href="#">Rapid securing of reference substances from <i>Peucedanum japonicum</i> Thunberg by recycling preparative high-performance liquid chromatography</a>	Won, H. J.; Lee, S. M.; Kim, D-Y.; Kwon, O-K.; Park, M. H.; Kim, J-H.; Ryu, H. W.; Oha, S-R.	J. Chromatogr. B <b>2019</b> , 1133, 121835.
Measurement and analysis	<a href="#">Fast and non-derivative method based on high-performance liquid chromatography-charged aerosol detection for the determination of fatty acids from <i>Agastache rugosa</i> (Fisch. et Mey.) O. Ktze. seeds</a>	Yang, R.; Wu, Z.; Pu, Y.; Zhang, T.; Wang, B.	Nat. Prod. Res. <b>2019</b> , 33, 1969-1974.
Instrumentation			
Sample preparation	<a href="#">Integrated liquid chromatography-mass spectrometry and nuclear magnetic resonance spectra for the comprehensive characterization of various components in the Shuxuening injection</a>	Yu, Y.; Li, J.; Guo, L.; Di, C.; Qin, X.; Li, Z.	J. Chromatogr. A <b>2019</b> , 1599, 125-135.
Separation	<a href="#">A new cinnamide derivative and two new <math>\beta</math>-carboline alkaloids from the stems of <i>Picrasma quassioides</i></a>	Zhang, J.; Wang, C-X.; Song, W-J.; Li, S.; Fan, C-L.; Chen, G-D.; Hu, D.; Yao, X-S.; Gao, H.	Fitoterapia <b>2019</b> , 139, 104375.
Detection	<a href="#">Quality assessment of <i>Astragali radix</i> from different production areas by simultaneous determination of thirteen major compounds using tandem UV/charged aerosol detector</a>	Zhanga, C-E.; Lianga, L-J.; Yua, X-H.; Wua, H.; Tub, P-F.; Maa, Z-H.; Zhaoa, K-J.	J. Pharm. Biomed. Anal. <b>2019</b> , 165, 233-241.
Authentication of supplements	<a href="#">Comparative quality of the forms of decoction pieces evaluated by multidimensional chemical analysis and chemometrics: <i>Poria cocos</i>, a pilot study</a>	Zhu, L-X.; Xu, J.; Wu, Y.; Su, L-F.; Lam, K. Y. C.; Qi, E. R.; Dong, X-P.; Chen, H-B.; Liu, Y-D.	J. Food Drug Anal. <b>2019</b> , 27, 766-777.
Application examples			
Substances A-C	<a href="#">A multi-detector chromatographic approach for characterization and quantitation of botanical constituents to enable in silico safety assessments</a>	Baker, T. R.; Regg, B. T.	Anal. Bioanal. Chem. <b>2018</b> , 410, 5143-5154.
Substances D-G	<a href="#">Preparation and identification of oligosaccharides in lotus seeds and determination of their distribution in different parts of lotus</a>	Chen, L.; Hu, D.; Liang, X.; Zhao, J.	Electrophor. <b>2018</b> , 29, 2020-2028.
Substances H-M	<a href="#">Laser microdissection hyphenated with high performance gel permeation chromatography-charged aerosol detector and ultra performance liquid chromatography-triple quadrupole mass spectrometry for histochemical analysis of polysaccharides in herbal medicine: Ginseng, a case study</a>	Chen, Q. L.; Chen, Y. J.; Zhou, S. S.; Yip, K. M.; Xu, J.; Chen, H. B.; Zhao, Z. Z.	Int. J. Biol. Macromol. <b>2018</b> , 107, 332-342.
Substances N-T			
Substances U-Z	<a href="#">Development and application of bio-sample quantification to evaluate stability and pharmacokinetics of inulin-type fructo-oligosaccharides from <i>Morinda officinalis</i></a>	Chia, L.; Chena, L.; Zhangb, J.; Zhaoa, J.; Lia, S.; Zheng, Y.	J. Pharm. Biomed. Anal. <b>2018</b> , 156, 125-132.
Literature	<a href="#">Chemical characterization of a variety of cold-pressed gourmet oils available on the Brazilian market</a>	Cicero, N.; Albergamo, A.; Salvo, A.; Bua, G. D.; Bartolomeo, G.; Mangano, V.; Rotondo, A.; Stefano, V. D.; Di Bella, G.; Dugo, G.	Food Res. Int. <b>2018</b> , 109, 517-525.
Glossary	<a href="#">Effect of solvent extraction system on the antioxidant activities of three invasive alien species and quantification of phenolic by compounds by HPLC</a>	Datta, S.; Sinha, B. K.; Seal, T.	J. Pharmacog. Phytochem. <b>2018</b> , 7, 3963-3970.
	<a href="#">Chemical characterization and bioactive properties of decoctions and hydroethanolic extracts of <i>Thymus carnosus</i> Boiss.</a>	Martins-Gomes, C.; Taghouti, M.; Schäfer, J.; Bunzel, M.; Silva, A. M.; Nunes, F. M.	J. Funct. Foods <b>2018</b> , 43, 154-164.

# Peer reviewed journal references

	Title	Authors	Publication
Table of contents			
Summary	<a href="#">Application of high performance liquid chromatography with charged aerosol detection (LC-CAD) for unified quantification of synthetic cannabinoids in herbal blends and comparison with quantitative NMR results</a>	Popławska, M.; Błażewicz, A.; Kamiński, K.; Bednarek, E.; Fijatek, Z.; Kozerski, L.	Forensic Toxicol. <b>2018</b> , 36, 122-140.
Overview: Dietary supplements and botanical natural products	<a href="#">A rapid high-performance liquid chromatography method for the simultaneous estimation of water-soluble vitamin in ten wild edible plants consumed by the tribal people of North-eastern Region in India</a>	Seal, T.; Chaudhuri, K.; Pillai, B.	Phcog. Mag. <b>2018</b> , 14, 72-77.
Overview: Global market	<a href="#">Analysis of vitexin in aqueous extracts and commercial products of Andean <i>Passiflora</i> species by UHPLC-DAD</a>	Sepúlveda, P.; Costa, G. M.; Aragón, D. M.; Ramos, F.; Castellanos, L.	J. App. Pharm. Sci. <b>2018</b> , 8, 81-86.
Measurement and analysis	<a href="#">Identification of cis/trans isomers of menaquinone-7 in food as exemplified by dietary supplements</a>	Szterk, A.; Zmysłowski, A.; Bus, K.	Food Chem. <b>2018</b> , 243, 403-409.
Instrumentation	<a href="#">Techniques for the analysis of pentacyclic triterpenoids in medicinal plants</a>	Xu, C.; Wang, B.; Pu, Y.; Tao, J.; Zhang, T.	J. Sep. Sci. <b>2018</b> , 41, 6-9.
Sample preparation	<a href="#">Determination of the Marker Diarylheptanoid Phytoestrogens in <i>Curcuma comosa</i> Rhizomes and Selected Herbal Medicinal Products by HPLC-DAD</a>	Yingngam, B.; Brantner, A.; Jinarat, D.; Kaewamatawong, R.; Rungsevijitprapa, W.; Suksamrarn, A.; Piyachaturawat, P.; Chokchaisiri, R.	Chem. Pharm. Bull. (Tokyo) <b>2018</b> , 66, 65-70.
Separation			
Detection	<a href="#">A new exploration of <i>Dregea volubilis</i> flowers: Focusing on antioxidant and antidiabetic properties</a>	Das, B.; De, A.; Das, M.; Das, S.; Samanta, A.	S. African J. Bot. <b>2017</b> , 109, 16-24.
Authentication of supplements	<a href="#">A new approach to the rapid separation of isomeric compounds in a <i>Silybum marianum</i> extract using UHPLC core-shell column with F5 stationary phase</a>	Fibigr, J.; Šatinský, D.; Solich, P.	J. Pharma. Biomed. Anal. <b>2017</b> , 134, 203-213.
Application examples	<a href="#">A UHPLC method for the rapid separation and quantification of anthocyanins in acai berry and dry blueberry extracts</a>	Fibigr, J.; Šatinský, D.; Solich, P.	J. Pharma. Biomed. Anal. <b>2017</b> , 143, 204-213.
Substances A-C	<a href="#">A UHPLC method for the rapid separation and quantification of phytosterols using tandem UV/Charged aerosol detection - A comparison of both detection techniques</a>	Fibigr, J.; Šatinský, D.; Solich, P.	J. Pharm. Biomed. Anal. <b>2017</b> , 140, 274-280.
Substances D-G	<a href="#">Qualitative and quantitative characterization of two licorice root species (<i>Glycyrrhiza glabra</i> L. and <i>Glycyrrhiza uralensis</i> Fisch.) by HPTLC, validated by HPLC and DNA sequencing</a>	Frommenwiler, D. A.; Maire-Widmer, V.; Upton, R.; Nichols, J.; Heubl, G.; Reich, E.	J. Planar Chromatogr. <b>2017</b> , 30, 467-473.
Substances H-M			
Substances N-T	<a href="#">Comprehensive quantitative analysis of 32 chemical ingredients of a Chinese patented drug sanhuang tablet</a>	Fung, H. Y.; Lang, Y.; Ho, H. M.; Wong, T. L.; Ma, D. L.; Leung, C. H.; Han, Q. B.	Molecules <b>2017</b> , 22, E111.
Substances U-Z	<a href="#">Comprehensive characterization and identification of antioxidants in <i>Folium Artemisiae argyi</i> using high-resolution tandem mass spectrometry</a>	Han, B.; Xin, Z.; Ma, S.; Liu, W.; Zhang, B.; Ran, L.; Yi, L.; Ren, D.	J. Chromatogr. A <b>2017</b> , 1063, 84-92.
Literature	<a href="#">Investigations on the Constituents of SagaPro Tablets, a Food Supplement Manufactured From <i>Angelica archangelica</i> Leaf</a>	Kowal, N. M.; Eyjolfsson, R.; Olafsdottir, E. S.	Pharmazie. <b>2017</b> , 72, 3-4.
Glossary	<a href="#">Application of Charged Aerosol Detection in Traditional Herbal Medicines</a>	Liang, L.; Jiang, Y.; Tu, P.	In Charged Aerosol Detection for Liquid Chromatography and Related Separation Techniques; Gamache, P. H., Ed.; Wiley: New York, <b>2017</b> ; p 341.
	<a href="#">Phenolic profiles of Lauraceae plant species endemic to Laurisilva forest: A chemotaxonomic survey</a>	Llorent-Martínez, E. J.; Spínola, V.; Castilho, P. C.	Ind. Crops Prods. <b>2017</b> , 107, 1-12.

# Peer reviewed journal references

	Title	Authors	Publication
Table of contents			
Summary	<a href="#">Nutritional Evaluation Of Fresh And Dried Goji Berries Cultivated In Italy</a>	Niro, S.; Fratianni, A.; Panfili, G.; Falasca, L.; Cinquanta, L.; Alam, R.	It. J. Food Sci. <b>2017</b> , 29, Epub.
Overview: Dietary supplements and botanical natural products	<a href="#">The Estimation of the Traditionally Used Yarrow (<i>Achillea millefolium</i> L. Asteraceae) Oil Extracts With Anti-Inflammatory Potential in Topical Application</a>	Tadić, V.; Arsić, I.; Zvezdanović, J.; Zugić, A.; Cvetković, D.; Pavkov, S.	J. Ethnopharmacol. <b>2017</b> , 199, 138-148.
Overview: Global market	<a href="#">Optimization of ultrasonic-assisted extraction of fatty acids in seeds of <i>Brucea javanica</i> (L.) Merr. from different sources and simultaneous analysis using high-performance liquid chromatography with Charged Aerosol Detection</a>	Wu, Li, L.; Li, N.; Zhang, T.; Pu, Y.; Zhang, X.; Zhang, Y.; Wang, B.	Molecules <b>2017</b> , 22, E931.
Measurement and analysis	<a href="#">A modification on the vector cosine algorithm of Similarity Analysis for improved discriminative capacity and its application to the quality control of <i>Magnoliae Flos</i></a>	Yang, G.; Zhao, X.; Fan, G.	J. Chromatogr. A <b>2017</b> , 1518, 34-45.
Instrumentation	<a href="#">Simultaneous fingerprint, quantitative analysis and anti-oxidative based screening of components in <i>Rhizoma Smilacis Glabrae</i> using liquid chromatography coupled with Charged Aerosol and coulometric array detection</a>	Yang, G.; Zhao, X.; Wen, J.; Zhou, T.; Fan, G.	J. Chromatogr. B <b>2017</b> , 1049-1050, 41-50.
Sample preparation	<a href="#">Determination of three saponins in rhizoma and fibrous root of <i>Trillium tschonoskii</i> and <i>Trillium kamschaticum</i></a>	Yang, Y. J.; Sun, X. G.; Yang, J.; Li, Q.; Zhang, J.; Zhao, Y.; Ma, B. P.; Guo, B. L.	Zhongguo Zhong Yao Za Zhi <b>2017</b> , 42, 1146-1151.
Separation	<a href="#">Rapid separation and simultaneous quantitative determination of 13 constituents in <i>Psoraleae fructus</i> by a single marker using high-performance liquid chromatography with diode array detection</a>	Zhang, Y.; Chen, Z.; Xu, X.; Zhou, Q.; Liu, X.; Liao, L.; Zhang, Z.; Wang, Z.	J. Sep. Sci. <b>2017</b> , 40, 4191-4202.
Detection	<a href="#">Production of Schisandrin A and Schisandrin B from Callus and Suspension Cell Cultures of <i>Schisandra chinensis</i></a>	Zhou, Y. Q.; Li, T-C.; Cheng, Y-P.	Nat. Prod. Commun. <b>2017</b> , 12, Epub.
Authentication of supplements	<a href="#">Phenolic profile and antioxidant activity of <i>Coleostephus myconis</i> (L.) Rchb.f.: An underexploited and highly disseminated species</a>	Bessada, S. M. F.; Barreira, J. C. M.; Barros, L.; Ferreira, I. C. F. R.; Oliveira, M. B.	Ind. Crops Prods. <b>2016</b> , 89, 45-51.
Application examples	<a href="#">Qualitative and Quantitative Characterization of Phenolic and Diterpenoid Constituents in Danshen (<i>Salvia miltiorrhiza</i>) by Comprehensive Two-Dimensional Liquid Chromatography Coupled With Hybrid Linear Ion Trap Orbitrap Mass</a>	Cao, J-L.; Wei, J-C.; Hu, Y-J.; He, C-W.; Chen, M-W.; Wan, J-B.; Li, P.	J. Chromatogr. A <b>2016</b> , 1427, 79-89.
Substances A-C	<a href="#">LC-PDA-ESI-MSn analysis of phenolic and iridoid compounds from <i>Globularia</i> spp</a>	Friščić, M.; Bucar, F.; Pilepić, K. H.	J. Mass Spec. <b>2016</b> , 51, 1211-1236.
Substances D-G	<a href="#">Transformation of <i>Panax notoginseng</i> saponins by steaming and <i>Trichoderma longibrachiatum</i></a>	Ge, F.; Huang, Z.; Yu, H.; Wang, Y.; Liu, D.	J. Biotech. Biotechnology Equip. <b>2016</b> , 30, 165-172.
Substances H-M	<a href="#">Simultaneous Ultra Performance Liquid Chromatography Determination and Antioxidant Activity of Linarin, Luteolin, Chlorogenic Acid and Apigenin in Different Parts of Compositae Species</a>	Hwang, S. H.; Paek, J. H.; Lim, S. S.	Molecules <b>2016</b> , 21, 1609.
Substances N-T	<a href="#">Chemometrics applied to quality control and metabolomics for traditional Chinese medicines</a>	Liu, S.; Liang, Y-Z.; Liu, H-T.	J. Chromatogr. B, <b>2016</b> , 1015-1016, 82-91.
Substances U-Z	<a href="#">A non-derivative method for the quantitative analysis of isosteroidal alkaloids from <i>Fritillaria</i> by high performance liquid chromatography combined with charged aerosol detection</a>	Long, Z.; Guo, Z. M.; Acworth, I. N.; Liu, X. D.; Jin, Y.; Liu, X.G.; Liu, L.; Liang, L. N.	Talanta <b>2016</b> , 151, 239-244.
Literature	<a href="#">Verifying the botanical authenticity of commercial tannins through sugars and simple phenols profiles</a>	Malacarne, M.; Nardin, T.; Bertoldi, D.; Nicolini, G.; Larcher, R.	Food Chem. <b>2016</b> , 206, 274-283.
Glossary			



# Peer reviewed journal references

Table of contents	Title	Authors	Publication
Summary	<a href="#">Structure-based prediction of CAD response factors of dammarane-type tetracyclic triterpenoid saponins and its application to the analysis of saponin contents in raw and processed <i>Panax notoginseng</i></a>	Peng, M.; Zhang, T.; Ding, Y.; Yi, Y. X.; Yang, Y. J.; Le, J.	RSC Adv. <b>2016</b> , 6, 36987-37005.
Overview: Dietary supplements and botanical natural products	<a href="#">Quantitation of phenylpropanoids and iridoids in insulin-sensitising extracts of <i>Leonurus sibiricus</i> L. (Lamiaceae)</a>	Pitschmann.; A; Zehl, M.; Heiss, E.; Purevsuren, S.; Urban, E.; Dirsch, V. M.; Glasl, S.	Phytochem. Anal. <b>2016</b> , 27, 23-31.
Overview: Global market	<a href="#">HPLC determination of phenolic acids, flavonoids and ascorbic acid in four different solvent extracts of <i>Zanthoxylum acanthopodium</i>, a wild edible plant of Meghalaya state of India</a>	Seal, T.	Int. J. Pharm. Pharma. Sci. <b>2016</b> , 8, 103-109.
Measurement and analysis	<a href="#">Identification And Quantification Flavonoids In Two Wild Edible Plants, <i>Viburnum foetidum</i> And <i>Perilla ocimoides</i> Of North-Eastern Region In India, Using High Performance Liquid Chromatography With Diode Array Detection</a>	Seal, T.; Chaudhuri, K.	J. Med Plant Studies <b>2016</b> , 4, 79-85.
Instrumentation	<a href="#">Identification and Quantification flavonoids in three wild edible plants, <i>Houttuynia cordata</i>, <i>Solanum gilo</i> and <i>Solanum kurzii</i> of North-Eastern region in India, using High Performance Liquid Chromatography with Diode Array Detection</a>	Seal. T.; Chaudhuri, K.; Pillai, B.	J. Chem. Pharm. Res. <b>2016</b> , 8, 859-867.
Sample preparation	<a href="#">Identification and Quantification of phenolic acids by HPLC, in two wild edible plants viz. <i>Solanum gilo</i> and <i>Solanum kurzii</i> collected from North-Eastern region in India</a>	Seal. T.; Pillai, B.; Chaudhuri, K.	J. Chem. Biol. Phys. Sci. <b>2016</b> , 6, 1108-1121.
Separation	<a href="#">Phenolic Compounds in Chilean Mistletoe (<i>Quintral, Tristerix tetrandus</i>) Analyzed by UHPLC-Q/Orbitrap/MS/MS and Its Antioxidant Properties</a>	Simirgiotis, M. J.; Quispe, C.; Areche, C.; Sepúlveda, B.	Molecules <b>2016</b> , 21, 245-251.
Detection	<a href="#">Quantitative determination of 15 bioactive triterpenoid saponins in different parts of <i>Acanthopanax henryi</i> by HPLC with charged aerosol detection and confirmation by LC-ESI-TOF-MS</a>	Zhang, X. D.; Li, Z.; Liu, G. Z.; Wang, X.; Kwon, O. K.; Lee, H. K.; Whang, W. K.; Liu, X. Q.	J. Sep. Sci. <b>2016</b> , 39, 2252-2262.
Authentication of supplements	<a href="#">A Validated HPLC Method for Simultaneous Determination of Caffeoyl Phenylethanoid Glucosides and Flavone 8-C-glycosides in <i>Haberlea rhodopensis</i></a>	Zheleva-Dimitrova, D.; Nedialkov, P.; Giresser, U.	Nat. Prod. Commun. <b>2016</b> , 11, 791-792.
Application examples	<a href="#">Reversed-phase-liquid Chromatography Method for Separation and Quantification of Gallic Acid From Hydroalcoholic Extracts of <i>Qualea grandiflora</i> and <i>Qualea parviflora</i></a>	de Mesquita, M. L.; Leão, W. F.; Ferreira, M. R. A.; de Paula, J. E.; Espindola, L. S.; Soares, L. A. L.	Pharmacogn. Mag. <b>2015</b> , 11, S316-321.
Substances A-C	<a href="#">Ground green coffee beans as a functional food supplement – Preliminary study</a>	Dziki, D.; Gawlik-Dziki, U.; Pecio, L.; Różyło, R.; Świeca, M.; Krzykowski, A.; Rudy, S.	LWT - Food Sci. Tech. <b>2015</b> , 63, 691-699.
Substances D-G	<a href="#">Development of a reliable extraction and quantification method for glucosinolates in <i>Moringa oleifera</i></a>	Förster, N.; Ulrichs, C.; Schreiner, M.; Müller, C. T.; Mewis, I.	Food Chem. <b>2015</b> , 166, 456-464.
Substances H-M	<a href="#">Quantitative and qualitative investigations of pharmacopoeial plant material <i>Polygoni avicularis herba</i> by UHPLC-CAD and UHPLC-ESI-MS methods</a>	Granica, S.	Phytochem. Anal. <b>2015</b> , 26, 374-382.
Substances N-T	<a href="#">Qualitative and quantitative analyses of secondary metabolites in aerial and subaerial of <i>Scorzonera hispanica</i> (black salsify)</a>	Granica, S.; Lohwasser, U.; Jöhner, K.; Zidorn, C.	Food Chem. <b>2015</b> , 173, 321-331.
Substances U-Z	<a href="#">Variation of bioactive compounds content of 14 oriental strawberry cultivars</a>	Kim, S. K.; Kim, D. S.; Kim, D. Y.; Chun, C.	Food Chem. <b>2015</b> , 184, 196-202.
Literature	<a href="#">Three new compounds from the bark of <i>Antiaris toxicaria</i></a>	Li, X-S.; Zhu, J-J.; Zhao, H.; Li, S-L.; Hao, X-J.; Yao, X-S.; Tang, J-S.	Phytochem. Lett. <b>2015</b> , 13, 182-186.
Glossary	<a href="#">Decoding glycome of <i>Astragalus membranaceus</i> based on pressurized liquid extraction, microwave-assisted hydrolysis and chromatographic analysis</a>	Lv, G. P.; Hu, D. J.; Cheong, K. L.; Li, Z. Y.; Qing, X. M.; Zhao, J.; Li, S. P.	J. Chromatogr. A <b>2015</b> , 1409, 19-29.

# Peer reviewed journal references

	Title	Authors	Publication
Table of contents			
Summary	<a href="#">A New HPLC Method for Analysis of Natural Monacolin K in Red Yeast Rice Pharmaceutical Preparations</a>	Omm, S.	J. Pharmacogn. Nat. Prod. <b>2015</b> , 1, 1000106.
Overview: Dietary supplements and botanical natural products	<a href="#">Influence of different extraction methods on the quantification of selected flavonoids and phenolic acids from <i>Tilia cordata</i> inflorescence</a>	Oniszczyk, A.; Podgórski, R.	Ind. Crops Prod. <b>2015</b> , 76, 509-514.
Overview: Global market	<a href="#">Separation of Phenolic Acids and Flavonoids From <i>Trollius chinensis</i> Bunge by High Speed Counter-Current Chromatography</a>	Qin, Y.; Liang, Y.; Ren, D.; Qiu, X.; Li, X.	J. Chromatogr. B. <b>2015</b> , 1001, 82-89.
Measurement and analysis	<a href="#">Optimal extraction and fingerprinting of carotenoids by accelerated solvent extraction and liquid chromatography with tandem mass spectrometry</a>	Saha, S.; Walia, S.; Kundu, A.; Sharma, K.; Paul, R. K.	Food Chem. <b>2015</b> , 177, 369-375.
Instrumentation	<a href="#">Simultaneous Determination of Chlorogenic Acid, Caffeic Acid, Alantolactone and Isoalantolactone in <i>Inula helenium</i> by HPLC</a>	Wang, J.; Zhao, Y-M.; Zhang, M-L.; Shi, Q-W.	J. Chromatogr. Sci. <b>2015</b> , 53, 526-530.
Sample preparation	<a href="#">HPLC-PDA-CAD Fingerprints of Salt <i>Anemarrhenae rhizoma</i></a>	Wu, Y.; Gao, H.; Song, Z-B.	Zhong Yao Cai <b>2015</b> , 38, 942-947.
Separation	<a href="#">Accelerated, Microwave-Assisted, and Conventional Solvent Extraction Methods Affect Anthocyanin Composition from Colored Grains</a>	Abdel-Aal, E-S. M.; Akhtar, H.; Rabalski, I.; Bryan, M.	Food Sci. <b>2014</b> , 79, C138-C146.
Detection	<a href="#">Characterization of Flavonoid Glycosides From Fenugreek (<i>Trigonella foenum-graecum</i>) Crude Seeds by HPLC-DAD-ESI/MS Analysis</a>	Benayad, Z.; Gómez-Cordovés, C.; Es-Safi, N. E.	Int. J. Mol. Sci. <b>2014</b> , 15, 20668-20685.
Authentication of supplements	<a href="#">Ginseng total saponins reverse corticosterone-induced changes in depression-like behavior and hippocampal plasticity-related proteins by interfering with GSK-3<math>\beta</math>-CREB signaling pathway</a>	Chen, L.; Dai, J.; Wang, Z.; Zhang, H.; Huang, Y.; Zhao Y.	Evidence-Based Compl. Alt. Med. [Online] <b>2014</b> , 506735.
Application examples	<a href="#">Comparison of ultraviolet detection and charged aerosol detection methods for liquid-chromatographic determination of protoescigenin</a>	Filip, K.; Gryniewicz, G.; Gruza, M.; Jatczak, K.; Zagrodzki, B.	Acta Pol. Pharm. <b>2014</b> , 71, 933-938.
Substances A-C	<a href="#">Determination of C-glucosidic Ellagitannins in <i>Lythri salicariaeherba</i> by Ultra-High Performance Liquid chromatography coupled with charged aerosol detector: Method development and validation</a>	Granica, S.; Piwowarski, J. P.; Kiss, A. K.	Phytochem. Anal. <b>2014</b> , 25, 201-206.
Substances D-G	<a href="#">Phytochemical investigations of <i>Polygonum aviculare</i>, <i>Agrimonia eupatoria</i> and <i>Lythrum salicaria</i> by HPLC-DAD-MS3-CAD method - application of corona charged aerosol detection (CAD) for analysis of plant phenolics</a>	Granica, S.; Piwowarski, J.; Klbowska, A.; Krupa, K.; Kiss, A.	Planta Med. <b>2014</b> , 80, P2061.
Substances H-M	<a href="#">Development and Validation of a HPLC-UV-ESI-MS Method for the Simultaneous Quantitation of Ten Bioactive Compounds in Dahuang Fuzi Tang</a>	Guo, H.; Li, H.; Liu, X.; Cai, H.; Wu, L.; Cai, B-C.	Chin. J. Nat. Med. <b>2014</b> , 12, 952-960.
Substances N-T	<a href="#">Dereplication of microbial extracts and related analytical technologies</a>	Ito, T.; Masubuchi, M.	J. Antibiot. <b>2014</b> , 67, 353-360.
Substances U-Z	<a href="#">Simultaneous determination of bufadienolides and phenolic compounds in sea squill (<i>Drimia maritima</i> (L.) Stearn) by HPLC-DAD-MSn as a means to differentiate individual plant parts and developmental stages</a>	Knittel, D. N.; Stintzing, F. C.; Kammerer, D. R.	Anal. Bioanal. Chem. <b>2014</b> , 406, 6035-6050.
Literature	<a href="#">Study of quality control and uncertainty in estimation of capsaicinoids content and pungency in real chili samples using RP-HPLC</a>	Nagarnaik, M.; Dhakulkar, A.; Pandya, G. H.	Int. Food Res. J. <b>2014</b> , 21, 1101-1106.
Glossary	<a href="#">LC-MS Metabolic profiling for identification of active compounds in <i>Lonicera</i> species</a>	Ortmann, S.; Monschein, M.; Hartler, J.; Zhao, Y. M.; Miao, J. H.; Thallinger, G. G.; Bauer, R.	Planta Med. <b>2014</b> , 80, P1M12.

# Peer reviewed journal references

	Title	Authors	Publication
Table of contents			
Summary	<a href="#">Determination of total ginsenosides in ginseng extracts using charged aerosol detection with post-column compensation of the gradient</a>	Ouyang, L. F.; Wang, Z. L.; Dai, J. G.; Chen, L.; Zhao, Y. N.	Chinese J. Nat. Med. <b>2014</b> , 12, 857-868.
Overview: Dietary supplements and botanical natural products	<a href="#">Quantification of individual phenolic compounds' contribution to antioxidant capacity in apple: A novel analytical tool based on liquid chromatography with diode array, electrochemical, and charged aerosol detection</a>	Plaza, M.; Kariuki, J.; Turner, C.	J. Agric. Food Chem. <b>2014</b> , 62, 409-418.
Overview: Global market	<a href="#">Determination of flibanserin and tadalafil in supplements for women sexual desire enhancement using high-performance liquid chromatography with tandem mass spectrometer, diode array detector and charged aerosol detector</a>	Poplawska, M.; Blazewicz, A.; Zolek, P.; Fijalek, Z.	J. Pharm. Biomed. Anal. <b>2014</b> , 94, 45-53.
Measurement and analysis	<a href="#">Structural characterization of the degradation products of a minor natural sweet diterpene glycoside rebaudioside M under acidic conditions</a>	Prakash, I.; Chaturvedula, V. S. P.; Markosyan, A.	Int. J. Mol. Sci. <b>2014</b> , 15, 1014-1025.
Instrumentation	<a href="#">Profiling and quantitating the constituents of red clover extracts using UHPLC/UV/CAD/HRMS: A component of the safety assessment process</a>	Price, J. M.; Little, J. G.; Baker, T. R.	Planta Med. <b>2014</b> , 80, CL12.
Sample preparation	<a href="#">Characterization of fennel extracts and quantification of estragole: Optimization and comparison of accelerated solvent extraction and Soxhlet techniques</a>	Rodríguez-Solana, R.; Salgado, J. M.; Domínguez, J. M.; Cortés-Diéguez, S.	Indust. Crops Prod. <b>2014</b> , 52, 528-536.
Separation	<a href="#">Fumonisin measurement from maize samples by high-performance liquid chromatography coupled with corona charged aerosol detector</a>	Szekeres, A.; Budai, A.; Bencsik, O.; Németh, L.; Bartók, T.; Szécsi, Á.; Vágvölgyi, C.	J. Chromatogr. Sci. <b>2014</b> , 52, 1181-1185.
Detection	<a href="#">Simultaneous HPLC Quantitative Analysis of Active Compounds in Leaves of <i>Moringa oleifera</i> Lam</a>	Vongsak, B.; Sithisarn, P.; Gritsanapan, W.	J. Chromatogr. Sci. <b>2014</b> , 52, 641-645.
Authentication of supplements	<a href="#">Bioactive components on immuno-enhancement effects in the traditional Chinese medicine Shengqi Fuzheng Injection based on relevance analysis between chemical HPLC fingerprints and <i>in vivo</i> biological effects</a>	Wang, J.; Tong, X.; Li, P.; Liu, M.; Peng, W.; Cao, H.; Su, W.	J. Ethnopharmacol. <b>2014</b> , 155, 405-415.
Application examples	<a href="#">Simultaneous Determination of Seven Constituents in Si-Ni-San Decoction and a Compatibility Comparison Study Using HPLC-UV</a>	Wen, J.; Qiao, Y.; Xiong, Z.; Li, F.	Nat. Prod. Res. <b>2014</b> , 28, 1025-1029.
Substances A-C	<a href="#">Quantitative Determination of Triterpenoid Glycosides in <i>Fatsia japonica</i> Decne. &amp; Planch. Using High Performance Liquid Chromatography</a>	Ye, X.; Yu, S.; Lian, X-Y.; Zhang, Z.	J. Pharm. Biomed. Anal. <b>2014</b> , 88, 472-476.
Substances D-G	<a href="#">Evaluation and prediction of the antioxidant activity of <i>Epimedium</i> from multi-wavelength chromatographic fingerprints and chemometrics</a>	Zhang, L.; Zhang, Z.; Luo, Q.; Lu, H.; Lianga, Y.	Anal. Meth. <b>2014</b> , 6, 1036-1043.
Substances H-M	<a href="#">Dual-gradient liquid chromatography-tandem Charged Aerosol Detector for detection of macromolecules in traditional Chinese medicine injections</a>	Zhang, T. T.; Wang, Y.; Gu, D. H.; Zhang, L.; Jin, Y.	Chin. J. Anal. Chem. [Online] <b>2014</b> , 12.
Substances N-T	<a href="#">Preparation and quality assessment of high-purity ginseng total saponins by ion exchange resin combined with macroporous adsorption resin separation</a>	Zhao, Y. N.; Wang, Z. L.; Dai, J. G.; Chen, L.; Huang, Y. F.	Chin. J. Nat. Med. <b>2014</b> , 12, 382-392.
Substances U-Z	<a href="#">Profiling Hoodia extracts by HPLC with Charged Aerosol Detection and electrochemical array detection and pattern recognition</a>	Acworth, I.; Zhang, Q.; Thomas, D.	Planta Med. <b>2013</b> , 79, PP 7.
Literature	<a href="#">Novel GHB-derived natural products from European mistletoe (<i>Viscum album</i>)</a>	Amer, B.; Juvik, O. J.; Francis, G. W.; Fossen, T.	Pharm. Biol. <b>2013</b> , 51, 981-986.
Glossary	<a href="#">Methods for Extraction and Determination of Phenolic Acids in Medicinal Plants: A Review</a>	Arceusz, A.; Wesolowski, M.; Konieczynski, P.	Nat. Prod. Comm. <b>2013</b> , Epub.
	<a href="#">High-anthocyanin strawberries through cultivar selection</a>	Fredericks, C. H.; Fanning, K. J.; Gidley, M. J.; Netzel, G.; Zabar, D.; Herrington, M.; Netzel, M.	J. Sci. Food Agric. <b>2013</b> , 93, 846-852.

# Peer reviewed journal references

	Title	Authors	Publication
Table of contents			
Summary	<a href="#">Optimal extraction and fingerprint analysis of <i>Cnidii fructus</i> by accelerated solvent extraction and high performance liquid chromatographic analysis with photodiode array and mass spectrometry detections</a>	Gao, F.; Hu, Y.; Ye, X.; Li, J.; Chen, Z.; Fan, G.	Food Chem. <b>2013</b> , 141, 1962-1971.
Overview: Dietary supplements and botanical natural products	<a href="#">A new application of Charged Aerosol Detection in liquid chromatography for the simultaneous determination of polar and less polar ginsenosides in ginseng products</a>	Jia, S.; Li, J.; Yunusova, N.; Park, JH.; Kwon, SW.; Lee, J.	Phytochem. Anal. <b>2013</b> , 24, 374-380.
Overview: Global market	<a href="#">Lipophilic stinging nettle extracts possess potent anti-inflammatory activity, are not cytotoxic and may be superior to traditional tinctures for treating inflammatory disorders</a>	Johnson, T. A.; Sohn, J.; Inman, W. D.; Bjeldanes, L. F.; Rayburn, K.	Phytomed. <b>2013</b> , 20, 143-147.
Measurement and analysis	<a href="#">The effective parameters for subcritical water extraction of SDG lignan from flaxseed (<i>Linum usitatissimum</i> L.) using accelerated solvent extractor</a>	Kanmaz, E. O.; Ova, G.	Eur. Food Res. Tech. <b>2013</b> , 237, 159-166.
Instrumentation	<a href="#">Quantitative Analysis of Gleditsia Saponins in the Fruits of <i>Gleditsia sinensis</i> Lam. By High Performance Liquid Chromatography</a>	Lian, X-Y.; Zhang, Z.	J. Pharm. Biomed. Anal. <b>2013</b> , 75, 41-46.
Sample preparation	<a href="#">Researches on the Fingerprint of Dry Roots of <i>Angelica polymorpha</i> Maxim</a>	Lu, J.; Qin, W.; Jiaqi, X.; Xiaojin, Y.; Yuling, L.; Liang, H.	Med. Plant <b>2013</b> , 4, 89-91.
Separation	<a href="#">Application of high-performance liquid chromatography with charged aerosol detection for universal quantitation of undeclared phosphodiesterase-5 inhibitors in herbal dietary supplements</a>	Poplawska, M.; Blazewicz, A.; Bukowinska, K.; Fijalek, Z.	J. Pharm. Biomed. Anal. <b>2013</b> , 84, 232-243.
Detection	<a href="#">Rapid purification method for fumonisin B1 using centrifugal partition chromatography</a>	Szekeres, A.; Lorántfy, L.; Bencsik, O.; Kecskeméti, A.; Szécsi, Á.; Mesterházy, Á.; Vágvölgyi, C.	Food Addit. Contam. <b>2013</b> , 30, 147-155.
Authentication of supplements	<a href="#">Ethnopharmacological <i>in vitro</i> studies on Austria's folk medicine—An unexplored lore <i>in vitro</i> anti-inflammatory activities of 71 Austrian traditional herbal drugs</a>	Vogl, S.; Picker, P.; Mihaly-Bison, J.; Fakhrudin, N.; Atanasov, A. G.; Heiss, E. H.; Wawrosch, C.; Reznicek, G.; Dirsch, V. M.; Saukel, J.; Kopp, B.	J. Ethnopharmacol. <b>2013</b> , 149, 750-771.
Application examples	<a href="#">Determination of saikosaponins in <i>Bupleurum radix</i> from different locations by HPLC-CAD method and its immunomodulation effects on mouse splenocytes</a>	Wang, L. N.; Chen, B.; Wang, W.; Xu, N.; Jia, T. Z.	Lat. Am. J. Pharm. <b>2013</b> , 32, 1189-1195.
Substances A-C	<a href="#">Determination of Dehydroabietic Acid and Abietic Acid in Aqueous Alkali Extract of <i>Liquidambaris resina</i> by HPLC</a>	Wang, Y-F.; Wei, X-Y.	Zhongguo Zhong Yao Za Zhi <b>2013</b> , 38, 57-59.
Substances D-G	<a href="#">Influence of Extraction Methodologies on the Analysis of Five Major Volatile Aromatic Compounds of Citronella Grass (<i>Cymbopogon nardus</i>) and Lemongrass (<i>Cymbopogon citratus</i>) Grown in Thailand</a>	Chanthai, S.; Prachakoll, S.; Ruangviriyachai, C.; Luthria, D. L.	J. AOAC Int. <b>2012</b> , 95, 763-772.
Substances H-M	<a href="#">Quantification of <math>\alpha</math>-, <math>\beta</math>- and <math>\gamma</math>-mangostin in <i>Garcinia mangostana</i> fruit rind extracts by a reverse phase high performance liquid chromatography</a>	Aisha, A. F.; Abu-Salah, K. M.; Siddiqui, M. J.; Ismail, Z.; Majid, A. M.	J. Med. Plant Res. <b>2012</b> , 6, 4526-4534.
Substances N-T	<a href="#">Recent methodology in ginseng analysis</a>	Baek, S.; Bae, O.; Park, J.	J. Ginseng Res. <b>2012</b> , 36, 119-134.
Substances U-Z	<a href="#">Spectral analysis and chemical studies of the sweet constituent, rebaudioside A</a>	Chaturvedula, V.; Prakash, I.	Eur. J. Med. Plants <b>2012</b> , 2, 57-65.
Literature	<a href="#">Simultaneous Determination of Caffeine and Aspartame in Diet Supplements and Non-Alcoholic Beverages Using Liquid-Chromatography Coupled to Corona CAD and UV-DAD Detectors</a>	Grembecka, M.; Szefer, P.	Food Anal. Meth. <b>2012</b> , 5, 1010-1017.
Glossary	<a href="#">Accelerated Solvent Extraction of Alkylresorcinols in Food Products Containing Uncooked and Cooked Wheat</a>	Holt, M. D.; Moreau, R. A.; DerMarderosian, A.; McKeown, N.; Jacques, P. F.	J. Agric. Food Chem. <b>2012</b> , 60, 4799-4802.

# Peer reviewed journal references

Table of contents	Title	Authors	Publication
Summary	<a href="#">Simultaneous Analysis of Steviol and Steviol Glycosides by Liquid Chromatography with Ultraviolet Detection on a Mixed-Mode Column: Application to Stevia Plant Material and Stevia-Containing Dietary Supplements</a>	Jaworska, K.; Krynskiy, A.; Rader, J.	J. AOAC Int. <b>2012</b> , 95, 1588-1596.
Overview: Dietary supplements and botanical natural products	<a href="#">HPLC Fingerprint of Compound Xueshuantong Capsule</a>	Liang, J-P.; Liu, Z-Z.; Peng, W.; Su, W-W.	Zhong Yao Cai <b>2012</b> , 35, 1854-1858.
Overview: Global market	<a href="#">Characterization of Secondary Volatile Profiles in <i>Nigella sativa</i> Seeds From Two Different Origins Using Accelerated Solvent Extraction and Gas Chromatography-Mass Spectrometry</a>	Liu, X.; Abd El-Aty, A. M.; Cho, S. K.; Yang, A.; Park, J-H.; Shim, J. H.	Biomed. Chromatogr. <b>2012</b> , 26, 1157-1162.
Measurement and analysis	<a href="#">Phenolic Composition and Nutraceutical Properties of Organic and Conventional Cinnamon and Peppermint</a>	Lv, J.; Huang, H.; Yu, L.; Whent, M.; Niu, Y.; Shi, H.; Wang, T. T. Y.; Luthria, D.; Charles, D.; Yu, L. L.	Food Chem. <b>2012</b> , 132, 1442-1450.
Instrumentation	<a href="#">Accelerated Solvent Extraction for Natural Products Isolation</a>	Mottaleb, M. A.; Sarker, S. D.	In Natural Products Isolation. Methods in Molecular Biology (Methods and Protocols), vol 864; Sarker S., Nahar L. (eds); Humana Press; <b>2012</b> ; pp 75-87
Sample preparation			
Separation	<a href="#">Application of accelerated solvent extraction in the analysis of organic contaminants, bioactive and nutritional compounds in food and feed</a>	Sun, H.; Ge, X.; Lv, Y.; Wang, A.	J. Chromatogr. A <b>2012</b> , 1237, 1-23.
Detection	<a href="#">A consecutive preparation method based upon accelerated solvent extraction and high-speed counter-current chromatography for isolation of aesculin from <i>Cortex fraxinus</i></a>	Tong, X.; Zhou, T.; Xiao, X.; Li, G.	J. Sep. Sci. <b>2012</b> , 35, 3609-3614.
Authentication of supplements	<a href="#">Subcritical water extraction of alkaloids in <i>Sophora flavescens</i> Ait. and determination by capillary electrophoresis with field-amplified sample stacking</a>	Wang, H.; Lu, Y.; Chen, J.; Li, J.; Liu, S.	J. Pharm. Biomed. Anal. <b>2012</b> , 58, 146-151.
Application examples	<a href="#">Simultaneous Quantification Of Eleven Chemical Components In Traditional Herbal Medicinal Formula Socheongryongtang By HPLC-DAD And LC-MS</a>	Weon, J. B.; Yang, H. J.; Lee, B.; Ma, J. Y.; Ma, C. J.	J. Liq. Chromatogr. Rel. Technol. <b>2012</b> , 35, 2243-2254.
Substances A-C	<a href="#">Simultaneous determination of six active components in traditional herbal medicine 'Oyaksungisan' by HPLC-DAD</a>	Weon, J. B.; Yang, H. J.; Ma, J. Y.; Ma, C. J.	J. Nat. Med. <b>2012</b> , 66, 510-515.
Substances D-G			
Substances H-M	<a href="#">Rapid Separation and Identification of Anthocyanins from Flowers of <i>Viola yedoensis</i> and <i>V. prionantha</i> by High-performance Liquid Chromatography-Photodiode Array Detection-Electrospray Ionisation Mass Spectrometry</a>	Zhang, J.; Wang, L-S.; Gao, J-M.; Xu, Y-J.; Li, L-F.; Li, C-H.	Phytochem. Anal. <b>2012</b> , 23, 16-22.
Substances N-T			
Substances U-Z	<a href="#">Response Surface Modeling and Optimization of Accelerated Solvent Extraction of Four Lignans from <i>Fructus schisandrae</i></a>	Zhao, L-C.; He, Y.; Deng, X.; Yang, G-L.; Li, W.; Liang, J.; Tang, Q-L.	Molecules <b>2012</b> , 17, 3618-3629.
Literature	<a href="#">Screening of Medicinal Plants From Iranian Traditional Medicine for Acetylcholinesterase Inhibition</a>	Adhami, H-R.; Farsam, H.; Krenn, L.	Phytother. Res. <b>2011</b> , 25, 1148-1152.
	<a href="#">Dilute-and-shoot triple parallel mass spectrometry method for analysis of vitamin D and triacylglycerols in dietary supplements</a>	Byrdwell, W. C.	Anal. Bioanal. Chem. <b>2011</b> , 401, 3317-3334.
Glossary	<a href="#">Utilization of RP-HPLC fingerprinting analysis for the identification of diterpene glycosides from <i>Stevia rebaudiana</i></a>	Chaturvedula, V.; Prakash, I.	Int. J. Res. Phytochem. Pharmacol. <b>2011</b> , 1, 88-92.

# Peer reviewed journal references

	Title	Authors	Publication
Table of contents			
Summary	<a href="#">Acid and alkaline hydrolysis studies of stevioside and rebaudioside A</a>	Chaturvedula, V.; Prakash, I.	J. Appl. Pharm. Sci. <b>2011</b> , 1, 104-108.
Overview: Dietary supplements and botanical natural products	<a href="#">Evaluation of polyphenol contents in differently processed apricots using accelerated solvent extraction followed by high-performance liquid chromatography–diode array detector</a>	Erdoğan, S.; Erdemoğlu, S.	J. Food Sci. Nutr. <b>2011</b> , 62, 729-739.
Overview: Global market	<a href="#">Phenylphenalenones and related natural products from <i>Wachendorfia thyrsoiflora</i> L.</a>	Fang, J.; Paetz, C.; Hölscher, D.; Munde, T.; Schneider, B.	Phytochem. Letts. <b>2011</b> , 4, 203-208.
Measurement and analysis	<a href="#">Characterisation of phenolic acid derivatives and flavonoids from different morphological parts of <i>Helichrysum obconicum</i> by a RP-HPLC–DAD(–)–ESI–MSn method</a>	Gouveia, S.; Castilho, P. C.	Food Chem. <b>2011</b> , 129, 333-344.
Instrumentation			
Sample preparation	<a href="#">Anxiolytic activity of a supercritical carbon dioxide extract of <i>Souroubea sympetala</i> (Margaritaceae)</a>	Mullally, M.; Kramp, K.; Cayer, C.; Saleem, A.; Ahmed, F.; McRae, C.; Baker, J.; Goulah, A.; Otorola, M.; Sanchez, P.; Garcia, M.; Poveda, L.; Merali, Z.; Durst, T.; Trudeau, V. L.; Arnason, J. T.	Phytother. Res. <b>2011</b> , 25, 264-270.
Separation			
Detection	<a href="#">High-performance liquid chromatography analysis of plant saponins: An update 2005-2010</a>	Negi, J. S.; Singh, P.; Pant, G. J.; Rawat, M. S.	Pharmacogn. Rev. <b>2011</b> , 5, 155-158.
Authentication of supplements	<a href="#">Isolation and analysis of ginseng: advances and challenges</a>	Qi, L.; Wang, C.; Yuan, C.	Nat. Prod. Rep. <b>2011</b> , 28, 467-495.
Application examples			
Substances A-C	<a href="#">HPLC Analysis of Kaempferol and Quercetin Derivatives Isolated by Different Extraction Techniques From Plant Matrix</a>	Skalicka-Woźniak, K.; Szykowski, J.; Głowniak, K.	J. AOAC Int. <b>2011</b> , 94, 17-21.
Substances D-G	<a href="#">Identification and Quantification of Coumarins in <i>Peucedanum ostruthium</i> (L.) Koch by HPLC-DAD and HPLC-DAD-MS</a>	Vogl, S.; Zehl, M.; Picker, P.; Urban, E.; Wawrosch, C.; Reznicek, G.; Saukel, J.; Kopp, B.	J. Agric. Food Chem. <b>2011</b> , 59, 4371-4377.
Substances H-M	<a href="#">A HPLC-DAD Method for the Simultaneous Determination of Five Marker Components in the Traditional Herbal Medicine Bangpungtongsung-san</a>	Won, J. B.; Yang, H. J.; Ma, J. Y.; Ma, C. J.	Pharmacogn. Mag. <b>2011</b> , 7, 60-64.
Substances N-T	<a href="#">Application of accelerated solvent extraction coupled with high-performance counter-current chromatography to extraction and online isolation of chemical constituents from <i>Hypericum perforatum</i> L.</a>	Zhang, Y.; Liu, C.; Yu, M.; Zhang, Z.; Qi, Y.; Wang, J.; Wu, G.; Li, S.; Yu, J.; Hu, Y.	J. Chromatogr. A <b>2011</b> , 1218, 2827-2834.
Substances U-Z			
Literature	<a href="#">Bioactive molecules in <i>Kalanchoe pinnata</i> leaves: extraction, purification, and identification</a>	Abdellaoui, S. E.; Destandau, E.; Toribio, A.; Elfakir, C.; Lafosse, M.; Renimel, I.; André, P.; Cancellieri, P. Landemarre, L.	Anal. Bionanal. Chem. <b>2010</b> , 398, 1329-1338.
Glossary	<a href="#">Assessment of microcystin purity using charged aerosol detection</a>	Edwards, C.; Lawton, L. A.	J. Chromatogr. A <b>2010</b> , 1217, 5233-5238.
	<a href="#">Comparison between evaporative light scattering detection and charged aerosol detection for the analysis of saikosaponins</a>	Eom, H. Y.; Park, S. Y.; Kim, M. K.; Suh, J. H.; Yeom, H.; Min, J. W.; Kim, U.; Lee, J.; Youm, J. R.; Han, S. B.	J. Chromatogr. A <b>2010</b> , 1217, 4347-4354.
	<a href="#">Simultaneous analysis of seven alkaloids in Coptis-Evodia herb couple and Zuojin pill by UPLC with accelerated solvent extraction</a>	Gao, X.; Yang, X-W.; Marriott, P. J.	J. Sep. Sci. <b>2010</b> , 33, 2714-2722.
	<a href="#">Application of Response Surface Methodology to Optimize Pressurized Liquid Extraction of Antioxidant Compounds From Sage (<i>Salvia officinalis</i> L.), Basil (<i>Ocimum basilicum</i> L.) and Thyme (<i>Thymus vulgaris</i> L.)</a>	Hossain, M.; Brunton, N.; Martin-Diana, A.; Barry-Ryan, C.	Food Funct. <b>2010</b> , 1, 269-277.

# Peer reviewed journal references

	Title	Authors	Publication
Table of contents			
Summary	<a href="#">A Novel Method for Analyzing Solanesyl Esters in Tobacco Leaves Using Atmospheric Pressure Chemical Ionization/Mass Spectrometer</a>	Ishida, N.	J. Chromatogr. A <b>2010</b> , 1217, 5794-5801.
Overview: Dietary supplements and botanical natural products	<a href="#">Analysis of terpene lactones in a Ginkgo leaf extract by high- performance liquid chromatography using charged aerosol detection</a>	Kakigi, Y.; Mochizuki, N.; Icho, T.; Hakamatsuka, T.; Goda, Y.	Biosci. Biotechnol. Biochem. <b>2010</b> , 74, 590-594.
Overview: Global market	<a href="#">Optimization of pressurized liquid extraction for spicatoside A in <i>Liriope platyphylla</i></a>	Kim, S. H.; Kim, H. K.; Yang, E. S.; Lee, K. Y.; Du Kim, S.; Kim, Y. C.; Sung, S. H.	Sep. Pur. Technol. <b>2010</b> , 71, 168-172.
Measurement and analysis	<a href="#">Accelerated Solvent Extraction of Lignin from <i>Aleurites moluccana</i> (Candlenut) Nutshells</a>	Klein, A. P.; Beach, E. S.; Emerson, J. W.; Zimmerman, J. B.	J. Agric. Food Chem. <b>2010</b> , 58, 10045-10048.
Instrumentation	<a href="#">Production of surfactin and iturin by <i>Bacillus licheniformis</i> N1 responsible for plant disease control activity</a>	Kong, H. G.; Kim, J. C.; Choi, G. J.; Lee, K. Y.; Kim, H. J.; Hwang, E. C.; Lee, S. W.	J. Plant Pathol. <b>2010</b> , 26, 170-177.
Sample preparation	<a href="#">Application of Accelerated Solvent Extraction to the Investigation of Saikosaponins From the Roots of <i>Bupleurum falcatum</i></a>	Li, W.; Liu, Z.; Wang, Z.; Chen, L.; Sun, Y.; Hou, J.; Zheng, Y.	J. Sep. Sci. <b>2010</b> , 33, 1870-1876.
Separation	<a href="#">Subcritical Solvent Extraction of Anthocyanins from Dried Red Grape Pomace</a>	Monrad, J. K.; Howard, L. R.; King, J. W.; Srinivas, K.; Mauromoustakos, A.	J. Agric. Food Chem. <b>2010</b> , 58, 5, 2862-2868
Detection	<a href="#">Subcritical Solvent Extraction of Procyanidins from Dried Red Grape Pomace</a>	Monrad, J. K.; Howard, L. R.; King, J. W.; Srinivas, K.; Mauromoustakos, A.	J. Agric. Food Chem. <b>2010</b> , 58, 7, 4014-4021
Authentication of supplements	<a href="#">Quantification of the Total Amount of Artemisinin in Leaf Samples by Thin Layer Chromatography</a>	Quennoz, M.; Bastian, C.; Simonnet, X.; Grogg, A. F.	CHIMIA Int. J. Chem. <b>2010</b> , 64, 755-757.
Application examples	<a href="#">A biosynthetic pathway for BE-7585A, a 2-thiosugar-containing angucycline-type natural product</a>	Sasaki, E.; Ogasawara, Y.; Liu, H. W.	J. Am. Chem. Soc. <b>2010</b> , 132, 7405-7417.
Substances A-C	<a href="#">Characterization of Anthocyanins and Anthocyanidins in Purple-Fleshed Sweetpotatoes by HPLC-DAD/ESI-MS/MS</a>	Truong, V-D.; Deighton, N.; Thompson, R. T.; McFeeters, R. F.; Dean, L. O.; Pecota, K. V.; Yencho, G. C.	J. Agric. Food Chem. <b>2010</b> , 58, 404-410.
Substances D-G	<a href="#">Simultaneous Quantification of Marker Components in Ojeok-San by HPLC-DAD</a>	Weon, J. B.; Park, H.; Yang, H. J.; Ma, J. Y.; Ma, C. J.	J. Nat. Med. <b>2011</b> , 65, 375-380.
Substances H-M	<a href="#">Simultaneous Determination of Five Marker Constituents in Ssanghwa Tang by HPLC/DAD</a>	Won, J. B.; Ma, J. Y.; Um, Y. R.; Ma, C. J.	Pharmacogn. Mag. <b>2010</b> , 6, 111-115.
Substances N-T	<a href="#">Simultaneous determination of triterpenoid saponins from <i>Pulsatilla koreana</i> using high performance liquid chromatography coupled with a charged aerosol detector (HPLC-CAD)</a>	Yeom, H.; Suh, J. H.; Youm, J. R.; Han, S. B.	Bull. Kor. Chem. Soc. <b>2010</b> , 31, 1159-1164.
Substances U-Z	<a href="#">Comparison of extraction techniques and modeling of accelerated solvent extraction for the authentication of natural vanilla flavors</a>	Cicchetti, E.; Chaintreau, A.	J. Sep. Sci. <b>2009</b> , 32, 1957-1964.
Literature	<a href="#">Extraction of Bitter Acids from Hops and Hop Products Using Pressurized Solvent Extraction (PSE)</a>	Čulík, J.; Jurková, M.; Horák, T.; Čejka, P.; Kellner, V.; Dvořák, J.; Karásek, P.; Roth, R.	J. Inst. Brewing <b>2009</b> , 115, 220-225.
Glossary	<a href="#">Antioxidant Oligostilbenoids from the Stem Wood of <i>Hopea hainanensis</i></a>	Ge, H-M.; Yang, W-H.; Zhang, J.; Tan, R. X.	J. Agric. Food Chem. <b>2009</b> , 57, 5756-5761.

# Peer reviewed journal references

	Title	Authors	Publication
Table of contents			
Summary	<a href="#">Analysis of Phenolic Compounds From Different Morphological Parts of <i>Helichrysum devium</i> by Liquid Chromatography With On-Line UV and Electrospray Ionization Mass Spectrometric Detection</a>	Gouveia, S. C.; Castilho, P. C.	Rap. Commun. Mass Spectrom. <b>2009</b> , 23, 3939-3953.
Overview: Dietary supplements and botanical natural products	<a href="#">Pentacyclic Triterpene Distribution in Various Plants – Rich Sources for a New Group of Multi-Potent Plant Extracts</a>	Jäger, S.; Trojan, H.; Kopp, T.; Laszczyk, M. N.; Scheffler, A.	Molecules <b>2009</b> , 14, 2016-2031.
Overview: Global market	<a href="#">Certification of a pure reference material for the ginsenoside Rg1</a>	Kim, D.; Chang, J.; Sohn, H.; Cho, B.; Ko, S.; Nho, K.; Jang, D.; Lee, S.	Accredit. Qual. Assur. <b>2009</b> , 15, 81-87.
Measurement and analysis	<a href="#">Quality Evaluation of Golden Saxifrage (<i>Chrysosplenium alternifolium</i> L.) Through Simultaneous Determination of Four Bioactive Flavonoids by High-Performance Liquid Chromatography With PDA Detection</a>	Olszewska, M. A.; Gudej, J.	J. Pharm. Biomed. Anal. <b>2009</b> , 50, 771-777.
Instrumentation	<a href="#">Methods for Preparing Phenolic Extracts From Olive Cake for Potential Application as Food Antioxidants</a>	Suárez, M.; Romero, M-P.; Ramo, T.; Macià, A.; Motilva, M-J.	J. Agric. Food Chem. <b>2009</b> , 57, 1463-1472.
Sample preparation	<a href="#">Performance evaluation of Charged Aerosol and Evaporative Light Scattering detection for the determination of ginsenosides by LC</a>	Wang, L.; He, W. S.; Yan, H. X.; Jiang, Y.; Bi, K. S.; Tu, P. F.	Chromatograph. <b>2009</b> , 70, 603-608.
Separation	<a href="#">HPLC in natural product analysis: The detection issue</a>	Wolfender, J. L.	Planta Med. <b>2009</b> , 75, 719-734.
Detection	<a href="#">Aqueous extract of <i>Astragali radix</i> induces human natriuresis through enhancement of renal response to atrial natriuretic peptide</a>	Ai, P.; Yong, G.; Dingkun, G.; Qiuyu, Z.; Kaiyuan, Z.; Shanyan, L.	J. Ethnopharmacol. <b>2008</b> , 116, 413-421.
Authentication of supplements	<a href="#">Sensitive Determination of saponins in Radix et Rhizoma Notoginseng by Charged Aerosol Detector coupled with HPLC</a>	Bai, C. C.; Han, S. Y.; Chai, X. Y.; Jiang, Y.; Li, P.; Tu, P. F.	J. Liq. Chromatogr. Rel. Technol. <b>2008</b> , 32, 242-260.
Application examples	<a href="#">Chemotaxonomic differentiation between <i>Cortinarius infractus</i> and <i>Cortinarius subtortus</i> by supercritical fluid chromatography connected to a multi-detection system</a>	Brondz, I.; Høiland, K.	Trends Chromatogr. <b>2008</b> , 4, 79-87.
Substances A-C	<a href="#">Analysis of alkaloids in <i>Coptis chinensis</i> Franch by accelerated solvent extraction combined with ultra performance liquid chromatographic analysis with photodiode array and tandem mass spectrometry detections</a>	Chen, J.; Wang, F.; Liu, J.; Lee, F. S-C.; Wang, X.; Yang, H.	Analytica Chimica Acta <b>2008</b> , 613, 184-195.
Substances D-G	<a href="#">Photostability of rebaudioside A and stevioside in beverages</a>	Clos, J. F.; DuBois, G. E.; Prakash, I.	J. Agric. Food Chem. <b>2008</b> , 56, 8507-8513.
Substances H-M	<a href="#">Polyketide analysis using mass spectrometry, evaporative light scattering, and charged aerosol detector systems</a>	Pistorino, M.; Pfeifer, B. A.	Anal. Bioanal. Chem. <b>2008</b> , 390, 1189-1193.
Substances N-T	<a href="#">Influence of Altitudinal Variation on the Content of Phenolic Compounds in Wild Populations of <i>Calluna vulgaris</i>, <i>Sambucus nigra</i>, and <i>Vaccinium myrtillus</i></a>	Rieger, G.; Müller, M.; Guttenberger, H.; Bucar, F.	J. Agric. Food Chem. <b>2008</b> , 56, 9080-9086.
Substances U-Z	<a href="#">Analysis of Volatile Components in Qingshanlvshui Tea Using Solid-Phase Microextraction/ Accelerated Solvent Extraction-Gas Chromatography-Mass Spectrometry</a>	Zhan, J.; Lu, S.; Meng, Z.; Xiang, N.; Cao, Q.; Miao, M.	Se Pu. <b>2008</b> , 26, 301-305.
Literature	<a href="#">The real nature of the indole alkaloids in <i>Cortinarius infractus</i>: Evaluation of artifact formation through solvent extraction method development</a>	Brondz, I.; Ekeberg, D.; Høiland, K.; Bell, D.; Annino, A.	J. Chromatogr., A <b>2007</b> , 1148, 1-7.
Glossary	<a href="#">Phytochemical and analytical studies of extracts from <i>Rhodiola rosea</i> and <i>Rhodiola quadrifida</i></a>	Wiedenfeld, H.; Dumaa, M.; Malinowski, M.; Furmanowa, M.; Naranantuya, S.	Die Pharmazie <b>2007</b> , 62, 308-311.
	<a href="#">Linear aglycones are the substrates for glycosyltransferase DesVII in methymycin biosynthesis: analysis and implications</a>	Kao, C. L.; Borisova, S. A.; Kim, H. J.; Liu, H. W.	J. Am. Chem. Soc. <b>2006</b> , 128,5606-5607.



# Peer reviewed journal references

	Title	Authors	Publication
Table of contents			
Summary	<a href="#">Influence of Sample Preparation on Assay of Phenolic Acids from Eggplant</a>	Luthria, D. L.; Mukhopadhyay, S.	J. Agric. Food Chem. <b>2006</b> , 54, 1, 41-47.
Overview: Dietary supplements and botanical natural products	<a href="#">Optimization of extraction process for phenolic acids from Black cohosh (<i>Cimicifuga racemosa</i>) by pressurized liquid extraction</a>	Mukhopadhyay, S.; Luthria, D. L.; Robbins, R. J.	J. Sci. Food Agric. <b>2006</b> , 86, 156-162.
Overview: Global market	<a href="#">Pressurized hot water extraction of bioactive or marker compounds in botanicals and medicinal plant materials</a>	Ong, E. S.; Cheong, H. J. S.; Goh, D.	J. Chromatogr. A <b>2006</b> , 1112, 92-102.
Measurement and analysis	<a href="#">Comparison of soxhlet, ultrasound-assisted and pressurized liquid extraction of terpenes, fatty acids and Vitamin E from <i>Piper gaudichaudianum</i> Kunth</a>	Péres, V. F.; Saffi, J.; Inês, M.; Melechchi, S.; Abad, F. C.; de Assis Jacques, R.; Martinez, M. M.; Oliveira, E. C.; Caramão, E. B.	J. Chromatogr. A <b>2006</b> , 1105, 115-118.
Instrumentation	<a href="#">Comparison of the Chemical Composition of Extracts from <i>Scutellaria lateriflora</i> Using Accelerated Solvent Extraction and Supercritical Fluid Extraction versus Standard Hot Water or 70% Ethanol Extraction</a>	Bergeron, C.; Gafner, S.; Clausen, E.; Carrier, D. J.	J. Agric. Food Chem. <b>2005</b> , 53, 3076-3080.
Sample preparation	<a href="#">Determination of isoflavones in soy bits by fast column high-performance liquid chromatography coupled with UV-visible diode-array detection</a>	Klejduš, B.; Mikelová, R.; Petřlová, J.; Potěšil, D.; Adam, V.; Stiborová, M.; Hodek, P.; Vacek, J.; Kizek, R.; Kubáň, V.	J. Chromatogr. A <b>2005</b> , 1084, 71-79.
Separation	<a href="#">Identification and quantitation of eleven sesquiterpenes in three species of <i>Curcuma</i> rhizomes by pressurized liquid extraction and gas chromatography-mass spectrometry</a>	Yang, F. Q.; Li, S. P.; Chen, Y.; Lao, S. C.; Wang, Y. T.; Dong, T. T. X.; Tsim, K. W. K.	J. Pharm. Biomed. Anal. <b>2005</b> , 39, 552-558.
Detection	<a href="#">Free and bound phenolic compounds in barley (<i>Hordeum vulgare</i> L.) flours: Evaluation of the extraction capability of different solvent mixtures and pressurized liquid methods by micellar electrokinetic chromatography and spectrophotometry</a>	Bonoli, M.; Marconi, E.; Caboni, M. F.	J. Chromatogr. A <b>2004</b> , 1057, 1-12.
Authentication of supplements	<a href="#">Pressurized liquid extraction prior to liquid chromatography with electrochemical detection for the analysis of vitamin E isomers in seeds and nuts</a>	Delgado-Zamarreño, M. M.; Bustamante-Rangel, M.; Sánchez-Pérez, A.; Carabias-Martínez, R.	J. Chromatogr. A <b>2004</b> , 1056, 249-252.
Application examples	<a href="#">Pressurized fluid extraction of carotenoids from <i>Haematococcus pluvialis</i> and <i>Dunaliella salina</i> and kavalactones from <i>Piper methysticum</i></a>	Denery, J. R.; Dragull, K.; Tang, C. S.; Li, Q. X.	Anal. Chim. Acta <b>2004</b> , 501, 175-181.
Substances A-C	<a href="#">Simultaneous Determination of Ergosterol, Nucleosides and Their Bases From Natural and Cultured Cordyceps by Pressurized Solvent Extraction and High Performance Liquid Chromatography</a>	Li, P.; Li, S-P.; Gong, Y-X.; Wang, Y-T.	Yao Xue Xue Bao. <b>2004</b> , 39, 917-920.
Substances D-G	<a href="#">Determination of catechins by means of extraction with pressurized liquids</a>	Piñeiro, Z.; Palma, M.; Barroso, C. G.	J. Chromatogr. A <b>2004</b> , 1026, 19-23.
Substances H-M	<a href="#">Pressurized liquid extraction of isoflavones from soybeans</a>	Rostagno, M. A.; Palma, M.; Barroso, C. G.	Anal. Chim. Acta <b>2004</b> , 522, 169-177.
Substances N-T	<a href="#">Effects of Solvent and Temperature on Pressurized Liquid Extraction of Anthocyanins and Total Phenolics from Dried Red Grape Skin</a>	Ju, Z. Y.; Howard, L. R.	J. Agric. Food Chem. <b>2003</b> , 51, 5207-5213.
Substances U-Z	<a href="#">Determination of zearalenone from wheat and corn by pressurized liquid extraction and liquid chromatography-electrospray mass spectrometry</a>	Pallaroni, L.; von Holst, C.	J. Chromatogr. A <b>2003</b> , 993, 39-45.
Literature	<a href="#">Simultaneous determination of 13 quinolones from feeds using accelerated solvent extraction and liquid chromatography</a>	Pecorelli, I.; Galarini, R.; Bibi, R.; Floridi, A. I.; Casciarri, E.; Floridi, A.	Anal. Chim. Acta <b>2003</b> , 483, 81-89.
Glossary			

# Peer reviewed journal references

Table of contents
Summary
Overview: Dietary supplements and botanical natural products
Overview: Global market
Measurement and analysis
Instrumentation
Sample preparation
Separation
Detection
Authentication of supplements
Application examples
Substances A-C
Substances D-G
Substances H-M
Substances N-T
Substances U-Z
Literature
Glossary

Title	Authors	Publication
<a href="#">A review of modern sample-preparation techniques for the extraction and analysis of medicinal plants</a>	Huie, C. W.	Anal. Bioanal. Chem. <b>2002</b> , 373, 23-30.
<a href="#">Automated sample preparation by pressurized liquid extraction–solid-phase extraction for the liquid chromatographic–mass spectrometric investigation of polyphenols in the brewing process</a>	Papagiannopoulos, M.; Mellenthin, A.	J. Chromatogr. A <b>2002</b> , 976, 345-348.
<a href="#">Online coupling of pressurized liquid extraction, solid-phase extraction and high-performance liquid chromatography for automated analysis of proanthocyanidins in malt</a>	Papagiannopoulos, M.; Zimmermann, B.; Mellenthin, A.; Krappe, M.; Maio, G.; Galensa, R.	J. Chromatogr. A <b>2002</b> , 958, 9-16.
<a href="#">Pressurized liquid extraction for the determination of polyphenols in apple</a>	Alonso-Salces, R. M.; Korta, E.; Barranco, A.; Berrueta, L. A.; Gallo, B.; Vicente, F.	J. Chromatogr. A <b>2001</b> , 933, 37-43.
<a href="#">Determination of Polyphenolic Profiles of Basque Cider Apple Varieties Using Accelerated Solvent Extraction</a>	Alonso-Salces, R. M.; Korta, E.; Barranco, A.; Berrueta, L.A.; Gallo, B.; Vicente, F.	J. Agric. Food Chem. <b>2001</b> , 49, 3761-3767.
<a href="#">Stability of phenolic compounds during extraction with superheated solvents</a>	Palma, M.; Piñeiro, Z.; Barroso, C. G.	J. Chromatogr. A <b>2001</b> , 921, 169-174.
<a href="#">Indole alkaloids from the seeds of <i>Centaurea cyanus</i> (Asteraceae)</a>	Sarker, S. D.; Laird, A.; Nahar, L.; Kumarasamy, Y.; Jaspars, M.	Phytochem. <b>2001</b> , 57, 1273-1276.
<a href="#">Pressurized liquid extraction of medicinal plants</a>	Benthin, B.; Danz, H.; Hamburger, M.	J. Chromatogr. A <b>1999</b> , 837, 211-219.
<a href="#">Comparison of methods for extraction of flavanones and xanthenes from the root bark of the osage orange tree using liquid chromatography</a>	da Costa, C. T.; Margolis, S. A.; Benner Jr., B. A.; Horton, D.	J. Chromatogr. A <b>1999</b> , 831, 167-178.
<a href="#">Accelerated Solvent Extraction of Paclitaxel and Related Compounds from the Bark of <i>Taxus cuspidata</i></a>	Kawamura, F.; Kikuchi, Y.; Ohira, T.; Yatagai, M.	J. Nat. Prod. <b>1999</b> , 62, 244-247.
<a href="#">Accelerated Solvent Extraction: A Technique for Sample Preparation</a>	Richter, B. E.; Jones, B. A.; Ezzell, J. L.; Porter, N. L.; Avdalovic, N.; Pohl, C.	Anal. Chem. <b>1996</b> , 68, 1033-1039.

## Table of contents

Summary

Overview:  
Dietary supplements and  
botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of  
supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature


Glossary

## Aglycone


The compound that is left when the carbohydrate is removed from a glycoside. Examples include anthocyanidins (from anthocyanins), genistein (from genistin), and hesperetin (from hesperidin).

## Alkaloid


Alkaloids are a large group of structurally diverse, nitrogen-containing, secondary metabolites found in numerous flowering plants and to a lesser extent in animals (e.g., insects, amphibians, reptiles and mammals), fungi and bacteria. Their exact role in plants remains elusive but they may act as protectants preventing damage by herbivores and other pests. Many alkaloids are used in medicine; others are used as recreational drugs. Examples of alkaloids include atropine, caffeine, morphine, nicotine, quinine and strychnine.




Atropine  
(a tropane)




Caffeine  
(a methylxanthine)




Morphine  
(an isoquinoline)



Nicotine  
(a pyridine)



Quinine  
(a quinoline)



Strychnine  
(an indole)

## Table of contents

Summary

Overview:  
Dietary supplements and  
botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of  
supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z


Literature

Glossary

## Anthraquinone

Anthraquinone is a polycyclic aromatic compound with a 9,10-dioxanthracene core:

Numerous anthraquinone derivatives are found in plants including cascarin, catenarin, chrysophanol, emodin, physcion, and rhein. Naturally occurring anthraquinones and their derivatives are reported to possess anticancer, anti-inflammatory, antimicrobial and antimalarial properties.



Anthraquinone

## Ayurvedic Medicine

Ayurvedic medicine is a traditional approach used in India for thousands of years. It is a holistic in nature and uses herbs, diet, exercise, meditation, breathing exercises and other methods to treat illness by attempting to restore the balance between mind, body and spirit.

## Botanicals


Derived from plants. Includes fresh or dried parts of plants (e.g., roots, stems, leaves, flowers, fruits, and seeds) and their extracts.

## Chromophore

A chromophore is a chemical structure that can absorb light at a particular wavelength. The presence of a chromophore in a molecule is essential for it to be measured by UV/Vis absorbance detection.

## Coumestans

Coumestan (1-Benzoxolo[3,2-c]chromen-6-one) is the polycyclic aromatic backbone of a group of plant secondary metabolites called the coumestans. One such compound, coumestrol, is a phytoestrogen and can bind to mammalian estrogen receptors, although its activity is much weaker than the endogenous ligand, estradiol.



Coumestan

## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z


Literature

Glossary


## Flavonoids

Flavonoids, sometimes referred to as bioflavonoids are a class of plant secondary metabolites. Originally, the term flavonoid was used to describe only ketone containing compounds, but it now includes many non-ketone polyhydroxy-polyphenolic compounds (e.g., the flavanoids). It has been estimated that several thousand flavonoids occur naturally, and these can be subdivided based upon their chemical structure, as shown below.


Although flavonoids are reported to possess a wide range of biological and pharmacological activities (e.g., anti-allergic, anti-inflammatory, antioxidant, anti-microbial and anti-cancer) in *in vitro* studies, their potential health benefits *in vivo* are far less clear.




**Flavone**  
(Anthoxanthins  
- apigenin,  
luteolin,  
tangeritin)




**Flavonol**  
(Anthoxanthins  
- kaempferol,  
myricetin,  
quercetin)




**Flavan**  
(Flavan-3-ol: catechin, epicatechin  
gallate, epigallocatechin,  
epigallocatechin gallate,  
proanthocyanidins, theaflavins,  
thearubigins)




**Flavylium**  
(Anthocyanins and their  
aglycone anthocyanidins -  
cyanidin, delphinidin,  
malvidin, pelargonidin,  
peonidin, petunidin)



**Isoflavone**  
(Phytoestrogens such as  
daidzein and genistein and  
their glucosides  
daidzin and genistin)



**Flavanone**  
(Hesperetin,  
naringenin)



**Flavanonol**  
(Taxifolin)

## Flos

Flower blossoms.

## Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

## Glycoside


A glycoside is a secondary metabolite formed when a carbohydrate molecule is linked to another molecule (an aglycone) via a glycosidic bond. A wide variety of glycosides are found in plants – resulting from differences in the aglycone moiety, carbohydrate structure and type of glycosidic bond present.

## Herbal Medicine


As defined by The World Health Organization, herbal medicines include [“herbs, herbal materials, herbal preparations and finished herbal products, that contain as active ingredients parts of plants, or other plant materials, or combinations”](#).

## Mammalian Lignans

Mammalian lignans are phytoestrogens formed by the action of gut bacteria on inactive precursors contained within foods. For example, secoisolariciresinol and matairesinol, abundant in flax seed and whole grains, respectively, are converted in the intestine to the mammalian lignans, enterolactone and enterodiol. Flavonolignans are molecules composed of both a flavonoid and a lignan, for example, silibinin from milk thistle.



Secoisolariciresinol



Matairesinol

## Mass Spectrometry (MS)

LC-MS is a powerful analytical approach often used for the [identification and quantification](#) of analytes in botanical natural products. Additionally, LC-MS with [metabolomic workflows](#) can be used for authentication and adulteration studies.

ThermoFisher Scientific offers a complete range of MS detectors including [single quadrupole](#), [triple quadrupole](#), [ion trap](#) and [orbitrap](#) mass spectrometers.

To learn more about LC-MS see the ThermoFisher Scientific [Mass Spectrometry Learning Center](#).

Table of contents

Summary

Overview:  
Dietary supplements and botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

## Medhya effect


In Ayurvedic medicine the term Medhya Rasayanas refers to a group of several medicinal plants that can be used, singly or in combination, to improve memory and intellect – the medhya effect. *Bacopa monnieri*, called brahmi, is one of these plants.

## Metabolomics

Metabolomic approaches can be either untargeted or targeted. Untargeted metabolomics can be defined as the comprehensive measurement of all low-molecular-weight metabolites in a biological specimen. By contrast, targeted metabolomics focuses on the identification and quantitation of a defined subset of metabolites. Analytes are identified and quantified using a wide range of analytical technologies e.g., UHPLC-MS, UHPLC-absorbance detection, and UHPLC-CAD. Multivariate statistical methods can then be used for data interpretation and extraction of information (see PCA). Metabolomic approaches can be used, for example, to authenticate botanical supplements or to help identify possible adulteration.

## Oxypregnane Glycosides

These glycosated steroidal secondary plant metabolites contain a pregnane backbone and are particularly abundant in *Caralluma* and *Hoodia* species. Oxypregnane glycosides are proposed to be the compounds responsible for the purported biological activities of *Caralluma* and *Hoodia* supplements.



Pregnane

## Pharmacopeia

A pharmacopeia is a book containing detailed methods (monographs) for the determination and quality control of drugs, supplements, traditional medicines, food components etc. It is usually authorized by governments or an appropriate scientific society. Examples include the [United States Pharmacopeia](#), [The European Pharmacopeia](#) and [The Pharmacopeia of the People's Republic of China](#).

## Principal Component Analysis (PCA)

PCA is a mathematical approach used to look for patterns in large data sets. PCA software increases the ability to interpret complex data but without information loss. See "[Principal component analysis explained simply](#)" for greater details. PCA has a wide range of applications e.g., authentication of [botanicals](#), [supplements](#) and [foods](#).

## Table of contents

Summary

Overview:  
Dietary supplements and  
botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of  
supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary

## Phytoestrogen

Phytoestrogens, or dietary estrogens, are a group of plant-derived secondary metabolites that mimic the effects of endogenous estrogen. These xenoestrogens are structurally diverse and include coumestans, isoflavones, and mammalian lignans.

## Rattan

The dried stem of some climbing plants (e.g., honeysuckle and climbing palms).

## Saponins

Saponins are a group of terpene and steroid glycosides that are widely distributed in plants. They are characterized by their ability to form emulsions, produce a soapy lather, and acting as detergents.

## Secondary Metabolites


Primary metabolites consist of a wide range of biochemicals that are involved in normal growth and development of an organism. Secondary metabolites are compounds that are not essential for life but for example in plants, may play an adaptive role produced during periods of environmental stress or by acting as defense molecules formed in response to damage by pathogens or herbivores. Plant secondary metabolites include alkaloids, glycosides, phenolics and terpenes.

## Standards

The quantification of key analytes in a botanical natural product can be challenging. Method development may be hampered by the availability of authentic standards. Some standards are commercially available or can be obtained through collaboration. However, in many instances, the researcher may be required to [isolate](#), purify, [identify](#) and quantify analytes, themselves. The CAD is an extremely useful tool for analyte quantification as its response is [independent of chemical structure](#), and so can be used to measure the amount of an analyte even though [authentic standards are unavailable](#).

## Terpenes

Terpenes are a large group of structurally diverse plant secondary metabolites formed from units of isopentenyl pyrophosphate:



Isopentenyl Pyrophosphate

Terpenes include many compounds discussed in this document and include diterpenes e.g., taxane, triterpenes found in bacopa, Boswellia, black cohosh and ginseng, and the sesquiterpene, ursolic acid found in ginkgo. For a more in-depth discussion of the classification of the various terpenes found in plants, see: [Terpenes and Terpenoids by S. Perveen](#).



## Table of contents

Summary

Overview:  
Dietary supplements and  
botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of  
supplements

Application examples

Substances A-C

Substances D-G

Substances H-M

Substances N-T

Substances U-Z

Literature


Glossary

## Traditional Medicine

Traditional medicine (TM) (sometimes called indigenous medicine) can be defined as [“the sum total of knowledge, skills and practices based on the theories, beliefs and experiences indigenous to different cultures that are used to maintain health, as well as to prevent, diagnose, improve or treat physical and mental illnesses”](#). TM is used by many cultures to treat sickness, prevent illness and maintain good health. It typically integrates herbal, animal and mineral medicines, along with spiritual therapies and exercises. Examples of TM include Ayurveda, Siddha medicine, Traditional African medicine and Traditional Chinese medicine.

## Xanthones

Xanthone (9H-Xanthen-9-one) is the central core of many secondary metabolite phytochemicals. Examples include mangostin, a xanthone found in the pericarp of the mangosteen fruit, and mangiferin found in mangoes.



Xanthone

## Xenoestrogens

The xenoestrogens are a group of compounds, sometimes called endocrine disruptors, that can mimic the effects of endogenous estrogen. These “foreign” compounds include natural compounds (e.g., phytoestrogens), as well as synthetic chemicals (e.g., bisphenols).

## Table of contents

Summary

Overview:  
Dietary supplements and  
botanical natural products

Overview: Global market

Measurement and analysis

Instrumentation

Sample preparation

Separation

Detection

Authentication of  
supplements

Application examples

Substances A-C

Substances D-G


Substances H-M

Substances N-T

Substances U-Z

Literature

Glossary



[Download Charged Aerosol Detection Bibliography.](#)

Find out more at [www.thermofisher.com/CAD](http://www.thermofisher.com/CAD)

**For Research Use Only. Not for use in diagnostic procedures.** ©2020 Thermo Fisher Scientific Inc. All rights reserved. AddaVax and related products are trademarked and owned by Invivogen. ABISCO-100 and related products are trademarked and owned by ISCONOVA. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details. **EB73582-EN 0120S**

**ThermoFisher**  
SCIENTIFIC