

# LC Capillaries and Fittings

Golnar Javadi  
Applications Engineer  
Chemistries and Supplies Technical Support  
December 12, 2024



# Agenda

- Why is it important to use the right connecting capillaries?
- Dispersion contribution from capillaries and fittings
- Types of capillaries and fittings
- Connecting capillaries for each LC system



# Why Is It Important to Use the Right Capillaries and Fittings?

Problems caused by poor/wrong connections

- Leaks
- Dispersion due to the volume of the capillaries (extracolumn dispersion)
- Mixing/dispersion due to use of the incorrect fitting
- Secondary interactions between the analytes and the body of the capillaries and fittings

# Dispersion Reduces HPLC Performance

What is dispersion?


- It is the original sample concentration being diluted as it is carried through the system plumbing (extracolumn volume)

What increases dispersion?



- Connecting tubing that is too long
- Connecting tubing that is too large in diameter
- Connections that have gaps and form small mixing chambers

# Extracolumn Volume

- Flow cell
- Heat exchanger
- Switching valve
- Needle seat
- **Connecting capillaries**



# Extracolumn Volume



Column: StableBond SB-C18, 4.6 x 30 mm, 3.5  $\mu$ m, mobile phase: 85% H<sub>2</sub>O with 0.1% TFA: 15% ACN, flow rate: 1.0 mL/min  
Temperature: 35 °C, sample: 1. phenylalanine 2. 5-benzyl-3,6-dioxo-2-piperazine acetic acid 3. Asp-Phe 4. aspartame.

# Taylor-Aris Equation: Peak Dispersion in Cylindrical Tubing

$$\sigma_{v,\text{ext}}^2 = \frac{\pi d^4 L_{\text{cap}} u}{96 D_m}$$

$\sigma_{v,\text{ext}}^2$  is the volume variance  
 $d$  is the tubing diameter  
 $L$  is the tubing length  
 $u$  is the linear velocity of the liquid  
 $D_m$  is the molecular diffusion coefficient


# Tubing Volume

Tubing Length	10 mm	50 mm	100 mm	150 mm
Tubing id	Volume	Volume	Volume	Volume
0.17 mm (green)	0.227 $\mu\text{L}$	1.1 $\mu\text{L}$	2.27 $\mu\text{L}$	3.3 $\mu\text{L}$
0.12 mm (red)	0.113 $\mu\text{L}$	0.55 $\mu\text{L}$	1.13 $\mu\text{L}$	1.65 $\mu\text{L}$







# As Column Efficiency Increases, Peak Width Decreases, and System Dispersion Becomes More of a Factor



Columns: **Eclipse Plus C18**, as described above. Mobile phase: A: water, B: MeOH, (15:85), injection volume: 6  $\mu$ L  
Temperature: 25  $^{\circ}$ C, flow rate: 1 mL/min. Detection: 310, 4 nm, 0.5 s response time, semimicro flow cell, sample: sunscreens.



# Optimizing Connecting Tubing Volume For UHPLC columns



# Types of Capillaries and Fittings

## What to consider

- Material
- Inner diameter and length
- Rigid or flexible
- Socket
- Compatible fitting types
- LC instrument



# Types of Capillaries and Fittings

## Material of capillaries

Material	Internal Diameter	Application	Maximum Pressure	Instrument Model
Stainless steel (SST)	0.075 – 0.94 mm	Standard HPLC, UHPLC, GPC, SFC, preparative HPLC	1300 bar	1100 through 1290 Infinity II
PEEK	0.12 – 0.5 mm		250 bar	All
PEEK-coated fused silica (FS/PEEK)	25 – 100 µm	Microflow HPLC	400 bar	1100 through 1260 capillary/nano LC series
PEEK-lined stainless steel (PEEK/SST)	0.17 mm	Bio-inert	400 bar	1260/1260 Infinity II bio-inert system
Titanium (Ti)	0.17 mm	Bio-inert*	1300 bar	1260/1260 Infinity II bio-inert system
MP35N	0.075 – 0.17 mm	Bio-compatible	1300 bar	1290 Infinity II bio LC

\*Only from pump to sampler.

# Types of Capillaries and Fittings

## Color coding of capillary inner diameters

Internal diameter in mm		Color code
0.015		Orange
0.025		Yellow
0.05		Beige
0.075		Black
0.075	MP35N	Black with orange stripe
0.1		Purple
0.12		Red
0.12	MP35N	Red with orange stripe
0.17		Green
0.17	MP35N	Green with orange stripe
0.20/0.25		Blue
0.20/0.25	MP35N	Blue with orange stripe
0.3		Grey
0.50		Bone White

# Types of Capillaries and Fittings

## Inner diameter for standard HPLC/UHPLC

### 0.17 mm:

- Standard for 1100, 1200, and 1260 series
- On 1290 Infinity and 1260/1290 Infinity II only used for pump-to-sampler connection

### 0.12 mm:

- Low-dispersion alternative for 1100, 1200, and 1260 series
- Standard on 1200 Series rapid resolution LC system
- Used as standard on the 1290 Infinity and 1260/1290 Infinity II for connections beyond the point of injection

### 0.075 mm:

- Ultralow dispersion kit, optional to optimize system for low dispersion




# Types of Flexible Capillaries and Fittings

## Socket

Standard socket: 1.6 mm od, 17 mm long



Long socket: 1.6 mm od, 21 mm long



InfinityLab Quick Turn capillary



InfinityLab Quick Connect capillary

No socket: 0.8 mm od



# Types of Fittings

## Nonremovable (permanent) fittings

### Advantages

- Simple and familiar handling
- Robust
- Economical

### Disadvantages

- Nonremovable
- Can't be re-adjusted
- Require a wrench

S 

SI 

SL 

SX 

M4 



Stainless steel fittings (S),  
5062-2418



Stainless steel long fittings (SL),  
5065-4454



Stainless steel  
intermediate (SI)

M4



Typical preswaged (pilot) distance: 2.1 mm (0.090 in.)



Stainless steel extra long fitting (SX),  
5065-9967



# Types of Fittings

## Removable (re-adjustable) fittings



Finger-tight PEEK fitting (SPF),  
0100-1516



PEEK long fittings (SPFL), 5062-8541



Finger-tight PEEK fittings (SPF),  
5065-4426



Double winged fitting (SPF),  
5042-6500



Finger-tight polyketone fitting (SPF),  
5042-8957



1200 bar removable fitting (SV),  
5067-4733



1200 bar removable long fitting (SLV),  
5067-4738



1200 bar removable extra long fitting  
(SXV), 5067-4739



UHP-FF fitting, bio-inert (RLO), 5067-5695



InfinityLab Quick Turn fitting, 5067-5966


# Structure of a Standard Capillary

## Compatible fittings

### Suitable fittings:

- 1.6 mm (1/16 inch) id
- 10-32 coned thread (Swagelock)

### Examples:



Stainless steel fittings (S),  
5062-2418



Finger-tight PEEK fitting (SPF),  
0100-1516



Stainless steel long fittings (SL),  
5065-4454



PEEK long fittings (SPFL), 5062-8541



Stainless steel extra long fitting (SX),  
5065-9967



Finger-tight PEEK fittings (SPF),  
5065-4426




Finger-tight polyketone fitting (SPF),  
5042-8957



Double winged fitting (SPF),  
5042-6500

# Structure of a Quick Turn Capillary

## Compatible fittings





### Suitable fittings:

- InfinityLab Quick Turn fitting
- All "standard" capillary compatible fittings

# Structure of a Quick Connect Capillary

## Compatible fittings




One side is always intended for use with a Quick Connect fitting. The opposite side is preconfigured with various fittings, typically with a nonpreswaged steel fitting. SX for capillary lengths of up to 150 mm, above with SL fitting.

# Structure of a Bio-inert Capillary (PEEK/ST)

## Compatible fittings

### Suitable fitting types:

- UHP-FF fitting (typical)
- InfinityLab Quick Turn fitting (alternative)
- Quick Connect, only at InfinityLab bio-inert Quick Connect heat exchanger (G7116-60041)




The socket with PEEK cap is comparatively soft compared to the conventional steel variant and requires a more careful approach to connectivity. The use of the fitting tool (p/n 5043-0915) is recommended for both fittings.

# Structure of Titanium Capillaries

## Compatible fittings

A few Titanium (Ti) capillaries are for use only with the 1260 Infinity and 1260 Infinity II bio-inert series, as they are only for connections within the pump and between the pump and the autosampler.

These capillaries are used with SwageLok fittings with a gold ferrule, and with removable (SV) fittings.



5005-0067: Gold-plated front ferrule, 1/16 in, 10/pk



# Structure of MP35N Capillaries

## Compatible fittings


- System capillaries (predefined with fittings)
- Quick Turn capillaries (capillaries without fittings)
- Quick Connect capillaries

### Suitable fittings:

- UHP-FF fittings
- "Standard" stainless steel fittings with gold ferrule
- Quick Turn fitting
- Quick Connect fitting with use of MP35N Quick Connect capillaries



MP35N long socket Quick Turn capillary



MP35N Quick Connect capillary

Stainless steel  
fittings with gold  
ferrule

All MP35N free-end capillaries can easily be used with a Quick Turn fitting, unless they are defined as an "M4" fitting.

### Part numbers:

5005-0067: Gold-plated front ferrule, 1/16 in, 10/pk  
5005-0442: Stainless steel fitting set, includes stainless steel nut, stainless steel back ferrule, and gold-plated front ferrule, 1/16 in, bio, 10/pk



# PEEK-coated Fused Silica Capillaries

## Compatible fittings

- Especially for micro and nano LC systems
- Most capillaries have a dedicated connection and come with the necessary fittings.
- Outer diameter: 0.8 mm
- Sensitive handling – risk of breakage when bending



### Color code

15 µm id : black  
25 µm id: yellow  
50 µm id: green  
75 µm id: blue  
100 µm id: black  
125 µm id: red

### Suitable fittings:

- 0.8 mm id PEEK fittings
- Depending on the connection goal, a stainless steel 10-32 or M4 thread fitting is necessary



Stainless steel fittings, male (G),  
5063-6593



Ferrule and stainless steel lock ring (W),  
5065-4423



PEEK fittings, plugs (MP),  
5065-4410



Double winged PEEK nut and ferrule  
(WPF), 5065-4422



PEEK fitting, long (WPFL),  
5022-6536



# Capillaries with M4 Fitting Connections

- To connect to the corresponding multicolumn selection valves or microvalves
- Outer diameter at M4 end: 0.8 mm
- Stainless steel and MP35N available

## M4 fitting end

- No socket - **0.8 mm od**
- **No standard fitting possible**



## Standard (10-32) fitting end

- Socket: 1.6 mm od
- **No M4 fitting possible**

Capillaries for M4 fittings are marked accordingly with an M, for example, 5067-5110 capillary SST 0.17 x 90 SV/M

# Labeling of Capillaries and Fittings

Type The type gives some indication on the primary function, like a loop or a connection capillary.		Material The material indicates which raw material is used.		Fitting left/fitting right The fitting left/right indicate which fitting is used on both ends of the capillary.	
Key	Description	Key	Description	Key	Description
Capillary	Connection capillaries	ST	Stainless steel	W	Swagelok + 0.8 mm Port id
Loop	Loop capillaries	Ti	Titanium	S	Swagelok + 1.6 mm Port id
Seat	Autosampler needle seats	PK	PEEK	M	Metric M4 + 0.8 mm Port id
Tube	Tubing	FS/PK	PEEK-coated fused silica <sup>1</sup>	E	Metric M3 + 1.6 mm Port id
Heat exchanger	Heat exchanger	PK <sub>1</sub> /ST	Stainless steel-coated PEEK <sup>2</sup>	U	Swagelok union
			PTFE	L	Long
		FS	Fused silica	X	Extra long
		MP35N	Nickel-cobalt-chromium-molybdenum alloy	H	Long head
				G	Small head SW 4
				N	Small head SW 5
				F	Finger-tight
				V	1200 bar
				B	Bio
				P	PEEK
				I	Intermediate

<sup>1</sup> Fused silica in contact with solvent

<sup>2</sup> Stainless steel-coated PEEK

## Examples:



Capillary ST 0.12 x X mm SL/M ps/ps



Capillary ST 0.17 x X mm S/S ps/ns


ps: preswaged  
ns: nonswaged

# Styles of Fittings



## Swagelok

- Two-piece ferrule
- Used on Agilent LCs
- Short nut
- Also available with a long nut



## Parker

- One-piece ferrule
- Short nut
- Very similar to Swagelok




## Waters


- Longer nut
- Used on Alliance systems

# Styles of Fittings


## Swagelok




## Waters




## Valco




## Parker




## Rheodyne




## Upchurch




# Potential Issues with Fittings



- Leak




- Peak shape problem




- No dead volume


# Peak Shape



# Swaging the SwageLok Fittings




# Swaging Fittings






# Swaging Fittings




Finger-tighten the assembly and ensure the tubing is completely seated in the bottom end of fitting.

# Swaging Fittings






# Inspect the Position of the Ferrule



[Video](#): How to properly swage stainless steel fittings onto stainless steel capillary tubing

# Tools






# Tightening Fittings into a Column




# Overtightened Fittings



# Agilent InfinityLab Fittings: Unique Spring-loaded Design

## Quick Connect and Quick Turn



InfinityLab Quick Connect fitting, 5067-5965

- The unique spring-loaded design applies a constant force to eliminate dead volume
- Finger tight to 1300 bar



InfinityLab Quick Turn fitting, 5067-5966

400 bar (finger-tight, user-dependent)  
1300 bar (tightened with a mounting tool  
(5043-0915))



Flyer: [5991-5164EN](#)

# InfinityLab Quick Connect and Quick Turn Fittings




	Quick Connect Fitting	Quick Turn Fitting
Connects to	Columns (or inline filters)	Column, various receiving ports with 10-32 port geometry
Maximum pressure	1300 bar (finger-tight, by turning the lever)	To 400 bar (finger-tight, user-dependent) 1300 bar (with mounting tool, 5043-0915)
Features	<ul style="list-style-type: none"> <li>• Spring-loaded function for dead volume-free connections (special capillaries)</li> <li>• Replaceable ferrule and capillary</li> <li>• Capillaries in various lengths and diameters are available</li> </ul>	<ul style="list-style-type: none"> <li>• Spring-loaded function for dead volume-free connections</li> <li>• Replaceable ferrule and capillary</li> <li>• Capillaries in various lengths and diameters are available</li> </ul>
Wetted material	PEEK (ferrule)	PEEK (ferrule)






# InfinityLab Fittings Last Longer

Chromatogram overlap before and after 200 reconnections




No visible change of chromatogram after 200 reconnections.

# InfinityLab Fittings Last Longer




- Tailing factors and theoretical plates stayed constant within the experimental allowance through the reconnection procedures.


# Installing the InfinityLab Quick Connect Fitting



Finger tighten the fitting until you feel resistance, then close the lever.




Leak tight to 1300 bar




[Video](#): Making Great Connections – Less stress, more reliable fittings

# Installing InfinityLab Quick Turn Fitting



Finger tighten for 400 bar (user-dependent)




Tighten using a wrench/mounting tool for 1300 bar



Mounting tool, p/n 5043-0915

# Not Just for Agilent Columns



# Connecting Capillaries for Each LC System





# 1100/1200/1260 Series System

Connection	p/n	Description
Solvent bottle to vacuum degasser	G1311-60003	Bottle head assembly for screw bottle (GL45), with glass filter, 20 µm (5041-2168)
Degasser to pump	G1322-67300	Tubing kit degasser, 300 mm tubing, 4/pk
Pump to autosampler	G1312-87303	Capillary, 0.17 mm x 400 mm
Pump (purge valve) to waste	5062-2461	PTFE tube, 5000 mm
Autosampler to column compartment	G1313-87305	Capillary, 0.17 mm x 180 mm
	G1313-87304	Capillary, 0.12 mm x 180 mm
Thermostatted ALS to column compartment	01090-87306	Capillary, 0.17 mm x 380 mm
	01090-87610	Capillary, 0.12 mm x 280 mm
Column compartment to column	G1316-87300	Capillary, 0.17 mm x 90 mm
	01090-87611	Capillary, 0.12 mm x 105 mm
Column to VWD (std flow cell)	5062-8522	Inlet tubing assembly PEEK, 0.17 mm x 600 mm
Column to DAD/MWD	G1315-87311	Capillary, 0.17 mm x 380 mm (S/S, ps/ns)
	G1315-87312	Capillary, 0.12 mm x 150 mm
VWD to waste	5062-8535	Waste accessory kit
DAD to waste	5062-2462	PTFE tubing 0.7 mm id, 1.6 mm od, 5 m

0.17 mm id capillaries	Standard setup
0.12 mm id capillaries	Rapid Resolution LC setup

Solvent cabinet

Vacuum degasser


Pump (Iso/quat/binary)

Autosampler

Sampler thermostat

Column compartment

UV-detector (DAD/MWD/VWD)




# 1100/1200/1260 Infinity Series System

Dual – stack configuration with cooled ALS

Connection	p/n	Description
Solvent bottle to vacuum degasser	G1311-60003	Bottle head assembly for screw bottle (GL45), with glass filter, 20 µm (5041-2168)
Degasser to pump	G1322-67300	Tubing kit degasser, 300 mm tubing, 4/pk
Pump to autosampler	G1312-87304	Capillary, 0.17 mm x 700 mm
Pump (purge valve) to waste	5062-2461	PTFE tube, 5000 mm
Thermostatted ALS to column compartment	01090-87309	Capillary, 0.17 mm x 380 mm
	01090-87610	Capillary, 0.12 mm x 280 mm
Column compartment to column	G1316-87300	Capillary, 0.17 mm x 90 mm
	01090-87611	Capillary, 0.12 mm x 105 mm
Column to VWD (std flow cell)	5062-8522	Inlet tubing assembly PEEK, 0.17 mm x 600 mm
Column to DAD/MWD	G1315-87311	Capillary, 0.17 mm x 380 mm
	G1315-87312	Capillary, 0.12 mm x 150 mm
VWD to waste	5062-8535	Waste accessory kit
DAD to waste	5062-2462	PTFE tubing, 0.7 mm id, 1.6 mm od, 5 m


0.17 mm id capillaries	Standard setup
0.12 mm id capillaries	Rapid Resolution LC setup





# 1290 Infinity Series System


Connection	p/n	Description
Solvent bottle to pump	G7120-60007	Bottle head assembly for screw bottle (GL45), with glass filter 20 µm (5041-2168)
Pump to autosampler	5067-4657	Capillary, 0.17 mm x 300 mm
Autosampler to column compartment	5067-4659	Capillary, 0.12 mm x 340 mm
Column compartment to column	5500-1188	Capillary, 0.12 mm x 105 mm (capillary comes without fittings, use Quick Turn fittings or SST fittings)
Column to DAD	5067-4660	Capillary, 0.12 mm x 280 mm
DAD to waste	5062-2462	PTFE tubing, 0.7 mm id, 1.6 mm od, 5 m



# 1290 Infinity Series System






## Dual – stack configuration with cooled ALS

Connection	p/n	Description
Solvent bottle to vacuum degasser	G7120-60007	Bottle head assembly for screw bottle (GL45), with glass filter, 20 µm (5041-2168)
Pump to autosampler	5500-1217	Capillary, 0.17 mm x 900 mm
Autosampler to column compartment	5067-4659	Capillary, 0.12 mm x 340 mm
Column compartment to column (low dispersion heat exchanger double)	5500-1188	Capillary, 0.12 mm x 105 mm (capillary comes without fittings, use Quick Turn fittings or SST fittings)
Column to DAD	5067-4660	Capillary, 0.12 mm x 280 mm
DAD to waste	5062-2462	PTFE tubing, 0.7 mm id, 1.6 mm od, 5 m



# 1260/1290 Infinity II Series System

## Single stack with multisampler and MCT

		Connection	p/n	Description
S		Solvent bottle to pump	G7120-60007	Bottle head assembly for screw bottle (GL45), with glass filter, 20 µm (5041-2168)
SI		Pump to multisampler	5500-1246	Capillary, 0.17 mm x 500 mm (SI/SI, ps/ps)
SL		Multisampler to MCT	5500-1157	Capillary, 0.12 mm x 500 mm, (SI/SI, ps/ns)
SX		Heat exchanger to column	5067-5957	Quick Connect assembly, 0.12 mm x 105 mm
			5500-1173	Quick Connect capillary, 0.12 mm x 105 mm
			5067-5965	Quick Connect fitting
M4		Column to DAD	5500-1191	Quick Turn capillary, 0.12 mm x 280 mm (comes without fitting)
			5067-5966	Quick Turn fitting
		DAD to waste	5062-2462	PTFE tubing, 0.7 mm id, 1.6 mm od, 5 m


Solvent cabinet

DAD

Multicolumn thermostat (MCT)

Multisampler

Pump (flex/high speed)



# 1260 Infinity II Series System

Single stack with vialsampler and internal column compartment

	Connection	p/n	Description
	Solvent bottle to pump	G7120-60007	Bottle head assembly for screw bottle (GL45), with glass filter, 20 µm (5041-2168)
	Pump to vialsampler	5500-1217	Capillary, 0.17 mm x 900 mm (SI/SX, ps/ps)
S	Sampler to 3 µL heat exchanger	5500-1249	Capillary, 0.12 mm x 120 mm, (SL/SL, ps,ns)
SI		5500-1250	Capillary, 0.17 mm x 120 mm (SL/SL, ps,ns)
SL	Heat exchanger to column		
	3 µL heat exchanger	5500-1238	Capillary, 0.12 mm x 105 mm (SL/SL, ps,ps)
SX	6 µL heat exchanger	5500-1240	Capillary, 0.17 mm x 105 mm (SL/SL, ps,ps)
M4	Column to DAD	5500-1191	Quick Turn capillary, 0.12 mm x 280 mm (comes without fitting)
	DAD to waste	5062-2462	PTFE tubing, 0.7 mm id, 1.6 mm od, 5 m


Solvent cabinet

Vial sampler

Internal column compartment (ICC)






UV detector (DAD/MWD/VWD)

Pump (quat/binary)



# 1260/1290 Infinity II Series System

Single stack with vialsampler and multicolumn thermostat

		Connection	p/n	Description
		Solvent bottle to pump	G7120-60007	Bottle head assembly for screw bottle (GL45), with glass filter, 20 µm (5041-2168)
S		Pump to vialsampler	5500-1245	Capillary, 0.17 mm x 400 mm (SI/SX, ps/ps)
SI		Sampler to heat exchanger	5500-1157	Capillary, 0.12 mm x 500 mm, (SL/SL, ps,ns)
SL		Heat exchanger to column	5067-5957	Quick Connect assembly, 0.12 mm x 105 mm
SX			5500-1173	Quick Connect capillary, 0.12 mm x 105 mm
			5067-5965	Quick Connect fitting
M4		Column to DAD	5500-1191	Quick Turn capillary, 0.12 mm x 280 mm (comes without fitting)
			5067-5966	InfinityLab Quick Turn fitting
		DAD to waste	5062-2462	PTFE tubing, 0.7 mm id, 1.6 mm od, 5 m


Solvent cabinet

UV detector (DAD/MWD/VWD)



Multicolumn compartment

Vialsampler

Pump (flex/high speed)



# InfinityLab Flex Benches



Flyer: [5991-5163EN](#)

# Resources on the Web

## [InfinityLab fittings on Agilent.com](#)

### [InfinityLab LC Supplies catalog](#)

Product catalog

### [Less Stress. More Reliable Connections.](#)


Product flyer

### [Agilent Bio-inert Capillaries and Fittings flyer](#)

Product overview


### [Agilent InfinityLab UHPLC Fittings](#)

Technical overview





# Contact Agilent Chemistries and Supplies Technical Support



1-800-227-9770 option 3, option 3:

Option 1 for GC and GC/MS columns and supplies

Option 2 for LC and LC/MS columns and supplies

Option 3 for sample preparation

Option 4 for spectroscopy supplies

Option 5 for chemical standards

Option 6 for former Prozyme products

Available in the U.S. and Canada, 8–5 all time zones

[gc-column-support@agilent.com](mailto:gc-column-support@agilent.com)

[lc-column-support@agilent.com](mailto:lc-column-support@agilent.com)

[spp-support@agilent.com](mailto:spp-support@agilent.com)

[spectro-supplies-support@agilent.com](mailto:spectro-supplies-support@agilent.com)

[chem-standards-support@agilent.com](mailto:chem-standards-support@agilent.com)

[advancebio.glycan@agilent.com](mailto:advancebio.glycan@agilent.com)

Web chat: Product pages of [agilent.com](https://www.agilent.com)



Thank You!