

Application Note No. 078

The Thermal Desorption / Pyrolysis of Glass Fibre Air Filters

Diane Nicholas

- **Direct desorption of analytes from sample matrix to column**
- **Multi-step thermal desorption followed by pyrolysis of same sample in injector**
- **Can be automated using the CombiPAL and LINEX**

Keywords:

Pyrolysis, thermal desorption, in-injector pyrolysis, multi mode inlet pyrolysis

Instrumentation

- ATAS GL Optic 2-200 programmable injector
- Agilent 5890 gas chromatograph with 5971 mass selective detector
- SGE CO₂ cryotrap

Sample analysed

Borosilicate glass air filters impregnated with 5 % phenolic and epoxy resins and some fluorocarbons.

Principles

- Place 2-4 mg of filter paper in a fritted liner, place in injector and start run
- Liner is firstly swept and the cryotrap turned-on
- The analytes are thermally desorbed under static flow conditions
- The analytes are swept onto the column with a small split flow and trapped
- The cryotrap is turned off and the separation and analysis is performed
- At the end of the run the injector and oven are cooled and the sample is pyrolysed, following the same principles

Chromatograms




Figure 1: Filter A: Thermal desorption at 200 °C




Figure 2: Filter A: Pyrolysis at 600 °C

ATAS GL INTERNATIONAL B.V.

De Sleutel 9, 5652 AS, Eindhoven, The Netherlands
 Tel. +31 (0)40 254 95 31 Fax. +31 (0)40 254 97 79
 E-mail: info@atasgl.com Internet: www.atasgl.com



We would like to thank Dr Anthony Lawson from Hollingworth & Vose Co Ltd. for his kind permission to publish this data.

ATAS GL INTERNATIONAL B.V.

De Sleutel 9, 5652 AS, Eindhoven, The Netherlands
Tel. +31 (0)40 254 95 31 Fax. +31 (0)40 254 97 79
E-mail: info@atasgl.com Internet: www.atasgl.com




Figure 3: Filter B: Thermal desorption at 200 °C




Figure 4: Filter B: Pyrolysis at 600 °C




Figure 5: Filter D: Thermal desorption at 200 °C




Figure 6: Filter D: Pyrolysis at 600 °C




Figure 7: Filter H: Thermal desorption at 200 °C




Figure 8: Filter H: Pyrolysis at 600 °C




Figure 9: Filter L: Thermal desorption at 200 °C




Figure 10: Filter L: Pyrolysis at 600 °C

ATAS GL INTERNATIONAL B.V.

De Sleutel 9, 5652 AS, Eindhoven, The Netherlands
Tel. +31 (0)40 254 95 31 Fax. +31 (0)40 254 97 79
E-mail: info@atasgl.com Internet: www.atasgl.com

Optic Parameters:

Liner: ATASGL Fritted
Mode: Expert
Gas Flows: Vent: 75 ml/min
Split: 10 ml/min
Initial temperature: 35 °C
Isothermal time: 1 min
Ramp rate: 16 °C/s
Final temperature: 200 °C (thermal desorption)
Final temperature: 600 °C (pyrolysis)
End time: 43.5 mins
Sweep pressure: 8 psi
Sweep time: 0.5 mins
Split open time: 0.5 mins
Desorption pressure: 0 psi
Desorption time: 2.5 mins
Transfer pressure: 7.1 psi
Transfer time: 2 mins
Initial pressure: 7.1 psi
Final pressure: 24.8 psi

Cryotrap Parameters:

Cryo on: 0.25 mins
Cryo off: 4 mins

GC Parameters:

Column: HP5-MS 30m x 0.25mm i.d. x 0.25um film
Initial temperature: 45 °C
Initial time: 5 mins
Ramp rate: 10 °C/min
Final temperature: 330 °C
Final time: 10 mins

MS Parameters:

Acquisition mode: Scan
Low mass: 50 m/z
High mass: 300 m/z
Sampling number: 2
Threshold: 500
Transfer line: 330 °C
Solvent delay: 2 mins